

Pierce™ Zinc Reversible Stain Kit

24582

Pub. No. MAN0011374

Rev. B.0

Pub Part No. 2160797.4

Number	Description
24582	Pierce Zinc Reversible Stain Kit , contains sufficient reagents to stain up to 20 mini gels

Kit Contents:**Zinc Stain**, 500mL**Zinc Developer**, 500mL**Zinc Eraser**, 500mL**Storage:** Upon receipt store kit at room temperature. Product shipped at ambient temperature.**Introduction**

The Thermo Scientific Pierce Zinc Reversible Stain Kit provides rapid detection of proteins on SDS-PAGE with sensitivity equivalent to silver stain (0.25ng). The kit uses a simple, two-step procedure that requires less than 15 minutes to complete. The result is a negative stain with an opaque white background and clear protein bands. The bands are visualized by placing the gel against a dark background. The white background can be removed with the eraser solution and the gel can be restained with a permanent stain such as Invitrogen SimplyBlue SafeStain (Product No. LC6060) or Thermo Scientific Pierce Silver Stain Kit (Product No. 24612). Alternatively, the gel can be destained with a Tris-glycine buffer, allowing proteins to be eluted or transferred for Western blotting. The eluted protein is free of dye contamination, making it ideal for downstream applications such as antibody generation or immunological detection. All reagents are provided ready to use, making this an ideal stain for multiple applications.

Procedure for Gel Staining and Destaining

1. Place gel in a container with 25-50mL of Zinc Stain. The gel must be completely immersed in the solution.
Note: The amount of stain required depends on the tray size being used. Typically, 25mL of stain is sufficient for one 8 × 10cm gel.
2. Place container on a shaker and agitate for 10 minutes
3. Remove the stain solution and add 25-50mL of Zinc Developer. Allow 1-2 minutes for the stain to develop. Check stain development against a dark background.
4. Once protein bands are clearly discernible against a dark background, stop stain development by removing the developer solution and adding water. Rinse gel twice for 1 minute with water. The gel may be kept in water for up to two weeks.
5. To destain (see Note below), place gel in a container with 25-50mL of Zinc Eraser. Agitate gel until stain is removed. Rinse gel several times with water.

Note: For Western transfer or protein elution, destain gel with 25-50mL of Tris-glycine buffer at pH 8 (e.g., 25mM Tris, 192mM glycine; Product No. 28380).

Troubleshooting

Problem	Possible Cause	Solution
Protein bands are indiscernible after staining	The stain was overdeveloped	Partially destain gel until bands become discernible, then stop destaining with a water rinse
		To prevent overdevelopment, monitor staining against a dark background
	No protein was present in the sample	Check protein concentration in the original sample

Related Products

89871	In-Gel Tryptic Digestion Kit , sufficient reagents for approximately 150 in-gel digestions
NW04120BOX	Bolt™ Bis-Tris Plus protein gels (see thermofisher.com/proteingels for a complete listing)
24615	Imperial™ Protein Stain , 1L, sufficient reagent to stain up to 50 mini gels
LC6060	SimplyBlue™ SafeStain
24612	Pierce™ Silver Stain Kit , sufficient reagents to stain up to 20 mini gels
24597	Pierce Color Silver Stain Kit , sufficient reagents to stain 25 (18cm ²) or 40 (10 × 13cm) 2D gels

For research use only. Not for use in diagnostic procedures.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at thermofisher.com. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2017 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.