

User Guide

Nitrate
Ion Selective
Electrode



ROSS and the COIL trade dress are trademarks of Thermo Fisher Scientific Inc.

AQUAfast, Cahn, ionplus, KNiPHE, No Cal, ORION, perpHect, PerpHecT, PerpHecTion, pHISA, pHuture, Pure Water, Sage, Sensing the Future, SensorLink, ROSS, ROSS Ultra, Sure-Flow, Titrator PLUS and TURBO2 are registered trademarks of Thermo Fisher.

1-888-pHAX-ION, A+, All in One, Aplus, AQUAsnap, AssuredAccuracy, AUTO-BAR, AUTO-CAL, AUTO DISPENSER, Auto-ID, AUTO-LOG, AUTO-READ, AUTO-STIR, Auto-Test, BOD AutoEZ, Cable-Free, CERTI-CAL, CISA, DataCOLLECT, DataPLUS, digital LogR, DirectCal, DuraProbe, Environmental Product Authority, Extra Easy/Extra Value, FAST QC, GAP, GLPcal, GLPcheck, GLPdoc, ISEasy, KAP, LabConnect, LogR, Low Maintenance Triode, Minimum Stir Requirement, MSR, NISS, One-Touch, One-Touch Calibration, One-Touch Measurement, Optimum Results, Orion Star, Pentrode, pHuture MMS, pHuture Pentrode, pHuture Quatrode, pHuture Triode, Quatrode, QuiKcheK, rf link, ROSS Resolution, SAOB, SMART AVERAGING, Smart CheK, SMART STABILITY, Stacked, Star Navigator 21, Stat Face, The Enhanced Lab, ThermaSense, Triode, TRIUMpH, Unbreakable pH, Universal Access are trademarks of Thermo Fisher.

© 2008 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

The specifications, descriptions, drawings, ordering information and part numbers within this document are subject to change without notice.

This publication supersedes all previous publications on this subject.

Table of Contents

| | |
|--|-----------|
| Introduction | 1 |
| Required Equipment | 2 |
| Serial Dilutions | 4 |
| Electrode Setup | 5 |
| 9307BNWP Nitrate Half-Cell Electrode Preparation | 5 |
| 900200 Double Junction Reference Electrode Preparation | 5 |
| 9707BNWP Nitrate Combination Electrode Preparation | 6 |
| Checking Electrode Operation (Slope) | 8 |
| Measurement Units | 9 |
| Sample Requirements | 9 |
| Measuring Hints | 10 |
| Electrode Storage | 11 |
| Electrode Maintenance | 12 |
| Analytical Techniques | 14 |
| Direct Calibration Technique | 16 |
| Small Volume Direct Calibration Technique | 20 |
| Low Level Calibration Technique | 24 |
| Known Addition Technique | 26 |
| Electrode Characteristics | 34 |
| Electrode Response | 34 |
| Reproducibility | 34 |
| Limits of Detection | 35 |
| Electrode Life | 35 |
| Temperature Effects | 36 |
| Interferences | 37 |
| Theory of Operation | 40 |
| Troubleshooting | 42 |
| Assistance | 43 |
| Warranty | 43 |
| Troubleshooting Checklist | 44 |
| Ordering Information | 45 |
| Specifications | 46 |

Introduction

This user guide contains information on the preparation, operation and maintenance for the nitrate ion selective electrode (ISE). General analytical procedures, electrode characteristics and electrode theory are also included in this user guide. Nitrate electrodes measure free nitrate ions in aqueous solutions quickly, simply, accurately and economically.

Technical Support Chemists can be consulted for assistance and troubleshooting advice. Within the United States call 1.800.225.1480 and outside the United States call 978.232.6000 or fax 978.232.6031. In Europe, the Middle East and Africa, contact your local authorized dealer. For the most current contact information, visit www.thermo.com/contactwater.

For the latest application and technical resources for Thermo Scientific Orion products, visit www.thermo.com/waterapps.

Nitrate ionplus[®] Sure-Flow[®] Plastic Membrane Combination ISE, Cat. No. 9707BNWP

The nitrate combination electrode has the sensing and reference half-cells built into one electrode, which decreases the amount of required solutions and reduces waste. The built-in Sure-Flow reference junction prevents electrode clogging and provides fast and stable readings. The nitrate combination electrode is available with a waterproof BNC connector, Cat. No. 9707BNWP. Electrodes with a waterproof BNC connector can be used on any ISE or mV meter with a BNC connection.

Nitrate Plastic Membrane Half-Cell ISE, Cat. No. 9307BNWP

The nitrate half-cell electrode must be used with the double junction reference electrode, Cat. No. 900200. The nitrate half-cell electrode is available with a waterproof BNC connector, Cat. No. 9307BNWP. Electrodes with a waterproof BNC connector can be used on any ISE or mV meter with a BNC connection.

Required Equipment

1. Thermo Scientific Orion ISE meter, such as the 4-Star pH/ISE meter or 5-Star pH/ISE/DO/conductivity meter; equivalent ISE meter; or mV meter with a 0.1 mV resolution.

Nitrate electrodes can be used on any ISE or mV meter with a BNC connection. The electrodes can also be used on meters with a variety of inputs when an adapter cable is used. Visit www.thermo.com/water for details.

2. Thermo Scientific Orion nitrate electrode.

The 9307BNWP nitrate half-cell electrode requires a separate reference electrode, Cat. No. 900200.

3. Magnetic stirrer or Thermo Scientific Orion stirrer probe, Cat. No. 096019. The stirrer probe can be used with 3-Star, 4-Star and 5-Star benchtop meters.
4. Volumetric flasks, graduated cylinders and beakers. Plastic labware is required for low level nitrate analysis.
5. Distilled or deionized water.
6. Nitrate electrode filling solution.

Use Optimum Results™ F filling solution, Cat. No. 900046, for the 9707BNWP nitrate combination electrode.

Use inner chamber filling solution, Cat. No. 900002, and nitrate ISA, Cat. No. 930711, or Optimum Results F filling solution, Cat. No. 900046, for the double junction reference electrode that is used with the 9307BNWP nitrate half-cell electrode. Either the nitrate ISA or Optimum Results F filling solution can be used as the outer chamber filling solution.

7. Nitrate calibration standards.

| Cat. No. | Description |
|----------|--|
| 920706 | 0.1 M NaNO ₃ nitrate calibration standard |
| 920707 | 1000 ppm as N nitrate calibration standard |
| 930707 | 100 ppm as N nitrate calibration standard |

8. Nitrate ionic strength adjuster (ISA), Cat. No. 930711. ISA provides a constant background ionic strength for samples and standards.

Nitrate interference suppressor solution (NISS), Cat. No. 930710, can be used in place of the nitrate ISA to remove a variety of interfering anions, including chloride ions, present in samples such as drinking water, wastewater and soils. Refer to the **Interferences** section for details.

9. Preservative solution – add 1 mL of preservative solution to every 100 mL of standards and samples to prevent biological degradation of the solutions.

Prepare a 1 M boric acid preservative solution by dissolving 6.2 grams of reagent-grade boric acid in 100 mL of boiling water. Let the solution cool.

Serial Dilutions

Serial dilution is the best method for the preparation of standards. Serial dilution means that an initial standard is diluted, using volumetric glassware, to prepare a second standard solution. The second standard is similarly diluted to prepare a third standard, and so on, until the desired range of standards has been prepared.

1. **To prepare a 10^{-2} M standard (140 ppm as N) –**
Pipet 10 mL of the 0.1 M standard into a 100 mL volumetric flask. Dilute to the mark with deionized water and mix well.
2. **To prepare a 10^{-3} M standard (14.0 ppm as N) –**
Pipet 10 mL of the 10^{-2} M standard into a 100 mL volumetric flask. Dilute to the mark with deionized water and mix well.
3. **To prepare a 10^{-4} M standard (1.40 ppm as N) –**
Pipet 10 mL of the 10^{-3} M standard into a 100 mL volumetric flask. Dilute to the mark with deionized water and mix well.

To prepare standards with a different concentration use the following formula:

$$C_1 * V_1 = C_2 * V_2$$

C_1 = concentration of original standard

V_1 = volume of original standard

C_2 = concentration of standard after dilution

V_2 = volume of standard after dilution

For example, to prepare 100 mL of a 100 ppm nitrate standard from a 1400 ppm nitrate standard:

$$C_1 = 1400 \text{ ppm nitrate}$$

$$V_1 = \text{unknown}$$

$$C_2 = 100 \text{ ppm nitrate}$$

$$V_2 = 100 \text{ mL}$$

$$1400 \text{ ppm} * V_1 = 100 \text{ ppm} * 100 \text{ mL}$$

$$V_1 = (100 \text{ ppm} * 100 \text{ mL}) / 1400 \text{ ppm} = 7.14 \text{ mL}$$

Electrode Setup

9307BNWP Nitrate Half-Cell Electrode Preparation

Remove the sensing module from the vial and save the vial for electrode storage. Make sure that the rubber electrode washer on the sensing module is in place. See **Figure 1**. Screw the sensing module into the electrode body until the module is finger-tight. To ensure electrical continuity, shake the electrode down like a clinical thermometer. Rinse the nitrate electrode with distilled water and then soak it in a 100 ppm or 10^{-2} M nitrate standard for 1 to 2 hours prior to use.


Note: Do not immerse the electrode past the rubber electrode washer.

900200 Double Junction Reference Electrode Preparation

Prepare the reference electrode according to the reference electrode user guide. Fill the reference electrode with inner chamber filling solution, Cat. No. 900002, and either nitrate ISA, Cat. No. 930711, or Optimum Results F filling solution, Cat. No. 900046 as the outer chamber filling solution.

Note: Do not use the outer chamber filling solution that ships with the 900200 double junction reference electrode because it contains interferences for nitrate measurements.

Figure 1
9307BNWP Nitrate Half-Cell Electrode



9707BNWP Nitrate Combination Electrode Preparation


Note: Do not touch the sensing membrane or reference pellet during the electrode assembly.

1. Remove the sensing module from the vial and save the vial for storage. Make sure that both O-rings are in place on the module. Remove the electrode handle from the box.
2. Unscrew the electrode cap. Slide the cap and spring down the electrode cable.
3. Hold the outer body sleeve and gently push the inner stem through the outer body. Slide the outer body sleeve down the electrode cable until it is beyond the inner stem.
4. Grasp the middle of the inner stem without touching the reference pellet. If a red storage tip is connected to the inner stem, unscrew it and save it for storage.
5. Screw the sensing module into the stem until it stops and the module is flush against the stem. Tighten the module an additional one-quarter turn. The module should be firmly attached to the stem. Do not overtighten the module.
6. Hold the electrode cable and slide the outer body, spring and cap over the inner stem.
7. Grasp the outer body sleeve, without touching the sensing membrane, and gently screw the cap onto the inner stem while pulling on the cable. Stop when an opposite force is felt. Do not overtighten or continue to turn the cap. The cap will not completely stop. If the inner body turns at all, the cap is too tight. Remove the cap and reassemble.
8. Press on the top of the cap with your thumb to make sure that the electrode has a smooth flushing motion and the outer body sleeve returns to its original position.
9. Install the flip spout cap onto the Optimum Results F filling solution bottle and lift the flip spout to a vertical position. Insert the spout into the electrode fill hole and add a small amount of filling solution to the reference chamber.
10. Hold the electrode body and use your thumb to push down on the electrode cap to allow a few drops of filling solution to drain out of the electrode. Release the electrode cap.
11. If the sleeve does not return to its original position, add filling solution and repeat step 10 until the sleeve returns to its original position.

12. Add filling solution to the electrode up to the fill hole.
13. Rinse the electrode with distilled water and soak it in a 100 ppm or 10^{-2} M nitrate standard for 1 to 2 hours prior to use.

Note: Add filling solution each day before using the electrode. The filling solution level should be at least one inch above the level of sample in the beaker to ensure a proper flow rate. The fill hole should always be open when taking measurements.

Figure 2
9707BNWP Nitrate Combination Electrode



Checking Electrode Operation (Slope)

These are general instructions that can be used with most meters to check the electrode operation. Refer to the meter user guide for more specific information.

This procedure measures electrode slope. Slope is defined as the change in millivolts observed with every tenfold change in concentration. Obtaining the slope value provides the best means for checking electrode operation.

1. If the electrode has been stored dry, prepare the electrode as described in the **Electrode Preparation** section.
2. Connect the electrode to a meter with a mV mode. Set the meter to the mV mode.
3. Add 100 mL of distilled water and 2 mL of ISA, Cat. No. 930711, into a 150 mL beaker. Stir the solution thoroughly.
4. Rinse the electrode with distilled water and place the electrode into the solution prepared in step 3.
5. Select either a 0.1 M or 1000 ppm nitrate standard. Pipet 1 mL of the standard into the beaker and stir the solution thoroughly. When a stable reading is displayed, record the electrode potential in millivolts.
6. Pipet 10 mL of the same standard into the same beaker and stir the solution thoroughly. When a stable reading is displayed, record the electrode potential in millivolts.
7. There should be a -54 to -60 mV difference between the two millivolt readings when the solution temperature is between 20 to 25 °C. If the millivolt potential is not within this range, refer to the **Troubleshooting** section.

Measurement Units

Nitrate concentration can be measured in moles per liter (M), parts per million (ppm) or any convenient concentration unit.

Table 1
Concentration Unit Conversion Factors

| Moles/Liter (M) | ppm as NO ₃ ⁻ | ppm as N |
|------------------|-------------------------------------|----------|
| 1.0 | 62000 | 14000 |
| 10 ⁻¹ | 6200 | 1400 |
| 10 ⁻² | 620 | 140 |
| 10 ⁻³ | 62.0 | 14.0 |
| 10 ⁻⁴ | 6.20 | 1.40 |

Sample Requirements

All samples must be aqueous and must not contain organic solvents. Contact Technical Support for information on using the electrode for specific applications.

The solution temperature must be less than 40 °C.

Samples and standards should be at the same temperature. A 1 °C difference in temperature for a 10⁻³ M nitrate solution will give rise to about a 1.5% error.

Interferences should be absent from all samples. See the **Interferences** section for a list of possible interferences. If interferences are present in the sample and cannot be removed, use the nitrate interference suppressor solution (NISS), Cat. No. 930710, in a 1:1 ratio of solution to NISS. Do not use ISA when using the nitrate interference suppressor solution.

In all analytical procedures, ISA or NISS must be added to all samples and standards before measurements are taken.

Measuring Hints

- Stir all standards and samples at a uniform, moderate rate. Place a piece of insulating material, such as Styrofoam or cardboard, between the magnetic stir plate and beaker to prevent measurement errors from the transfer of heat to the sample.
- Always use freshly prepared standards for calibration.
- Always rinse the electrode with distilled water between measurements and shake the electrode to remove the water and prevent sample carryover. Do not wipe or rub the electrode sensing module.
- Store the nitrate electrode in a 10^{-2} M or 100 ppm nitrate standard between measurements.
- The 9307BNWP nitrate half-cell electrode should be immersed in standards and samples to approximately half the length of the nitrate module. Do not immerse the nitrate electrode past the electrode washer. Immerse the reference electrode to the same depth as the nitrate electrode.
- Allow all standards and samples to reach the same temperature for precise measurements.
- Verify the electrode calibration every two hours by placing the electrode in a fresh aliquot of the least concentrated standard used for calibration. If the value has changed by more than 2%, recalibrate the electrode.
- After immersing the electrode in a solution, check the electrode sensing surface for air bubbles and remove air bubbles by reimmersing the electrode in the solution and gently tapping it.
- For high ionic strength samples, prepare standards with a background composition similar to the sample.
- The fill hole cover must be open during measurements to ensure a uniform flow of filling solution.
- If the combination electrode is used and the electrode is used in dirty or viscous samples or the electrode response becomes sluggish, empty the electrode completely, hold the junction open and flush the junction with distilled water. Empty any water from the electrode and refill it with fresh filling solution. Press down on the electrode cap to let a few drops of the filling solution flow out of the electrode and then replenish any lost solution.

Electrode Storage

Nitrate Combination Electrode Storage, Cat. No. 9707BNWP

For storage between measurements and up to three days, store the electrode in a 10^{-2} M or 100 ppm nitrate standard. The filling solution inside the electrode should not be allowed to evaporate, as crystallization will result.

For storage longer than one week, drain the electrode, flush the reference chamber with distilled water, disassemble the electrode and store the sensing module in the glass vial.

1. Grasp the outer body sleeve and unscrew the electrode cap. Slide the cap and spring assembly down the electrode cable.
2. Push the inner stem of the electrode handle out through the outer electrode sleeve, exposing the sensing module.
3. Rinse the inner stem and module well with distilled water. Gently blot dry to prevent damaging the sensing module.
4. Carefully unscrew the sensing module from the inner stem, taking care not to touch the sensing membrane.
5. Place the nitrate sensing module in the glass vial until it is needed again. Gently blot dry the inside of the inner stem and O-ring area, reassemble the electrode handle without the module and store it dry.

Nitrate Half-Cell Electrode Storage, Cat. No. 9307BNWP

The nitrate half-cell electrode should be rinsed thoroughly with distilled water and stored in a 10^{-2} M or 100 ppm nitrate standard. When storing the electrode for more than three days, rinse the nitrate half-cell electrode thoroughly with distilled water, shake the electrode dry, disassemble the electrode and store the sensing module in the glass vial.

Double Junction Reference Electrode Storage, Cat. No. 900200

The double junction reference electrode may be stored in the nitrate ISA or Optimum Results F filling solution between sample measurements and up to one week. The filling solution inside the electrode should not be allowed to evaporate, as crystallization will result.

For storage longer than one week, drain the reference electrode, flush the inside with distilled water and store the electrode dry.

Electrode Maintenance

Cleaning the Nitrate Sensing Module

If the electrode is exposed to high levels of interfering ions, it may drift and become sluggish in response. When this happens, restore normal performance by soaking the electrode for an hour in distilled water, emptying the old filling solution, filling the electrode with fresh filling solution and then soaking the electrode for a few hours in a 10^{-2} M or 100 ppm nitrate standard. If soaking the electrode does not restore normal electrode performance, replace the nitrate sensing module.

Nitrate Combination Electrode and Double Junction Reference Electrode Flushing

If the area between the electrode outer body and inner cone becomes clogged with sample or precipitate, flush the area with filling solution or distilled water.

1. Hold the electrode body with one hand and use your thumb to push down on the electrode cap to drain all of the filling solution out of the electrode.
2. Fill the electrode with distilled water and then push down on the cap until all the water is drained from the chamber. Repeat this procedure until all of the sample or precipitate is removed from the electrode.
3. Fill the electrode with fresh filling solution up to the fill hole. Push down on the cap to allow a few drops of filling solution to drain out of the electrode and then replenish the lost filling solution.
4. Rinse the electrode with distilled water and soak it in a 10^{-2} M or 100 ppm nitrate standard for 1 to 2 hours.

Replacing the Nitrate Sensing Module

The sensing membrane of plastic membrane electrodes will wear over time, indicated by low slope values, drift, poor reproducibility and loss of response in low level samples. The electrode response can be restored by replacing the sensing module. Each sensing module will last about three months with normal laboratory use, but the actual lifespan of the sensing module will depend on the type of samples that are measured.

For the 9707BNWP nitrate combination electrode, use the 97 series nitrate module, Cat. No. 970701.

Drain the electrode and flush the reference chamber with distilled water. Hold the outer body sleeve and unscrew the electrode cap. Slide the cap and spring assembly down the electrode cable. Push the inner stem of the electrode handle out through the outer electrode sleeve, exposing the sensing module. Rinse the inner stem and module well with distilled water. Gently blot dry to prevent damaging the sensing module. Carefully unscrew the sensing module from the inner stem and dispose of the old sensing module. Obtain a new 97 series nitrate module, Cat. No. 970701, and refer to the **9707BNWP Nitrate Combination Electrode Preparation** section for detailed instructions on assembling the electrode.

For the 9307BNWP nitrate half-cell electrode, use the 93 series nitrate module, Cat. No. 930702.

Rinse the electrode with distilled water. Carefully unscrew the sensing module from the electrode and dispose of the old sensing module. Obtain a new 93 series nitrate module, Cat. No. 930702, and refer to the **9307BNWP Nitrate Half-Cell Electrode Preparation** section for detailed instructions on assembling the electrode.

Analytical Techniques

A variety of analytical techniques are available to the analyst. The following is a description of these techniques.

Direct Calibration is a simple procedure for measuring a large number of samples. Only one meter reading is required for each sample. Calibration is performed using a series of standards. The concentration of the samples is determined by comparison to the standards. ISA is added to all solutions to ensure that samples and standards have similar ionic strength.

Low Level Calibration is similar to the direct calibration technique. This method is recommended when the expected sample concentration is less than 10^{-4} M or 1.4 ppm nitrate as nitrogen (N). A minimum three point calibration is recommended to compensate for the electrode's non-linear response at these concentrations. A special calibration standard preparation procedure is the best means of preparing low level calibration standards.

Incremental Techniques provide a useful method for measuring samples, since a calibration is not required. The different incremental techniques are described below. They can be used to measure the total concentration of a specific ion in the presence of a large (50 to 100 times) excess of complexing agents. As in direct calibration, any convenient concentration unit can be used.

Known Addition is useful for measuring dilute samples, checking the results of direct calibration (when no complexing agents are present), or measuring the total concentration of an ion in the presence of an excess complexing agent. The electrode is immersed in the sample solution and an aliquot of a standard solution containing the measured species is added to the sample. From the change in potential before and after the addition, the original sample concentration is determined.

Table 2
Recommended Measuring Techniques


| | Direct | Small Volume Direct | Low Level | Known Addition |
|-----------------------------------|--------|---------------------|-----------|----------------|
| [N] < 1.4 ppm | | | ✓ | |
| [N] > 1.4 ppm | ✓ | ✓ | | ✓ |
| Occasional Sampling | | | | ✓ |
| Small sample volume | | ✓ | | ✓ |
| Large number of samples | ✓ | | ✓ | ✓ |
| Reduce chemical usage | | ✓ | | |
| Field measurement | | ✓ | | |
| Ionic strength greater than 0.1 M | ✓ | | | ✓ |

Direct Calibration Technique

Typical Direct Calibration Curve

In the direct calibration procedure, a calibration curve is constructed either in the meter memory or on semi-logarithmic paper. Electrode potentials of standard solutions are measured and plotted on the linear axis against their concentrations on the log axis. In the linear regions of the curves, only two standards are needed to determine a calibration curve. In non-linear regions, more points must be taken. These direct calibration procedures are given for concentrations in the region of linear electrode response. Low level measurement procedures are given in a following section for measurements in the non-linear electrode region.

Figure 3
Typical Direct Calibration Curve



Direct Calibration Overview

The following direct measurement procedures are recommended for moderate to high level measurements. Samples must be in the linear range of the electrode – greater than 10^{-4} M or 1.4 ppm nitrate as N. A two point calibration is sufficient, although more points can be used. When using an ISE meter, sample concentrations can be read directly from the meter. When using a mV meter, a calibration curve can be prepared on semi-logarithmic graph paper, or a linear regression (against logarithmic concentration values) can be performed using a spreadsheet or graphing program.

Calibration Hints

- Standard concentrations should bracket the expected sample concentrations.
- Always add 2 mL of ISA, Cat. No. 930711, per 100 mL of standard or sample. If interferences are present in the sample and cannot be removed, add 50 mL of NISS, Cat. No. 930710, per 50 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution.
- For high ionic strength samples that have an ionic strength of 0.1 M or greater, prepare standards with a background composition similar to that of the samples, or measure the samples using the known addition method.
- During calibration, measure the least concentrated standard first, and work up to the most concentrated standard.

Direct Calibration Setup

1. Prepare the electrode as described in the **Electrode Preparation** section. If using the 9707BNWP combination nitrate electrode, fill the electrode with Cat. No. 900046. If using the 9307BNWP half-cell nitrate electrode with the 900200 reference electrode, fill the reference electrode with inner chamber filling solution, Cat. No. 900002, and nitrate ISA, Cat. No. 930711, or Optimum Results F filling solution, Cat. No. 900046 as the outer chamber filling solution.
2. Connect the electrode to the meter.
3. Prepare at least two standards that bracket the expected sample range and differ in concentration by a factor of ten. Standards can be prepared in any concentration unit to suit the analysis requirement. See the **Serial Dilution** section for instructions preparing standards. All standards should be at the same temperature as the samples.

Direct Calibration Procedure Using a Meter with an ISE Mode

Note: See the meter user guide for more specific information.

1. Add 100 mL of the less concentrated standard and 2 mL of ISA to a 150 mL beaker and stir the solution thoroughly.
2. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the less concentrated standard. Wait for a stable reading and adjust the meter to display the value of the standard, as described in the meter user guide.
3. Add 100 mL of the more concentrated standard and 2 mL of ISA to a second 150 mL beaker and stir the solution thoroughly.
4. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the more concentrated standard. Wait for a stable reading and adjust the meter to display the value of the second standard, as described in the meter user guide.
5. Record the resulting slope value. The slope should be between -54 and -60 mV when the standards are between 20 and 25 °C.
6. Add 100 mL of sample and 2 mL of ISA to a clean 150 mL beaker and stir the solution thoroughly.
7. Rinse the electrode with distilled water, blot it dry and place it into the sample. The concentration of the sample will be displayed on the meter.

Note: Other solution volumes may be used, as long as the ratio of solution to ISA remains 50:1.

Note: If interferences are present in the sample and cannot be removed, add 50 mL of NISS, Cat. No. 930710, per 50 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution. Other solution volumes may be used, as long as the ratio of solution to NISS remains 1:1.

Direct Calibration Procedure Using a Meter with a mV Mode

Note: See the meter user guide for more specific information.

1. Set the meter to the mV mode.
2. Add 100 mL of the less concentrated standard and 2 mL of ISA to a 150 mL beaker and stir the solution thoroughly.
3. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the less concentrated standard. When a stable reading is displayed, record the mV value and corresponding standard concentration.
4. Add 100 mL of the more concentrated standard and 2 mL of ISA to a second 150 mL beaker and stir the solution thoroughly.
5. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the more concentrated standard. When a stable reading is displayed, record the mV value and corresponding standard concentration.
6. Using semi-logarithmic graph paper, prepare a calibration curve by plotting the millivolt values on the linear axis and the standard concentration values on the logarithmic axis.
7. Add 100 mL of sample and 2 mL of ISA to a clean 150 mL beaker and stir the solution thoroughly.
8. Rinse the electrode with distilled water, blot it dry and place it into the beaker. When a stable reading is displayed, record the mV value.
9. Using the calibration curve prepared in step 6, determine the unknown concentration of the sample.

Note: Other solution volumes may be used, as long as the ratio of solution to ISA remains 50:1.

Note: If interferences are present in the sample and cannot be removed, add 50 mL of NISS, Cat. No. 930710, per 50 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution. Other solution volumes may be used, as long as the ratio of solution to NISS remains 1:1.

Small Volume Direct Calibration Technique

Take advantage of special design features available with the 9707BNWP ionplus combination nitrate electrode to meet your measuring needs. Due to the Sure-Flow reference, this electrode is able to measure sample volumes as small as 5 mL using a modified direct measurement procedure. Because less solution volume is required, the chemical usage of nitrate standards and ISA is reduced. This method is also convenient when making field measurements, since the 9707BNWP combination nitrate electrode does not require a separate reference electrode. All samples should have a concentration greater than 10^{-4} M or 1.4 ppm nitrate as N. A two point calibration is sufficient, although more points can be used. The following procedure recommends using 25 mL of sample. Smaller sample volumes can be used, as long as the final volume of solution is sufficient to cover the bottom of the electrode.

Calibration Hints

- Use the 9707BNWP ionplus combination nitrate electrode.
- Standard concentrations should bracket the expected sample concentrations.
- Always add 0.5 mL of ISA, Cat. No. 930711, per 25 mL of standard or sample. Always keep the ratio of standard or sample to ISA at 50:1.
- If interferences are present in the sample and cannot be removed, add 25 mL of NISS, Cat. No. 930710, per 25 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution. Always keep the ratio of standard or sample to NISS at 1:1.
- For high ionic strength samples that have an ionic strength of 0.1 M or greater, prepare standards with a background composition similar to that of the samples, or measure the samples using the known addition method.
- During calibration, measure the least concentrated standard first, and work up to the most concentrated standard.
- Calibrate with the same volume of standard as the volume of sample that is available for analysis.

Small Volume Direct Calibration Setup

1. Prepare the 9707BNWP combination nitrate electrode as described in the **Electrode Preparation** section and fill the electrode with Optimum Results F filling solution, Cat. No. 900046.
2. Connect the electrode to the meter.
3. Prepare at least two standards that bracket the expected sample range and differ in concentration by a factor of ten. Standards can be prepared in any concentration unit to suit the particular analysis requirement. See the **Serial Dilution** section for instructions on how to prepare standards. All standards should be at the same temperature as the samples. For details on temperature effects on electrode performance, refer to the **Temperature Effects** section.

Small Volume Direct Calibration Procedure Using a Meter with an ISE Mode

Note: See the meter user guide for more specific information.

1. Add 25 mL of the less concentrated standard and 0.5 mL of ISA to a 50 mL beaker and swirl the solution to mix.
2. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the less concentrated standard. Wait for a stable reading and adjust the meter to display the value of the standard, as described in the meter user guide.
3. Add 25 mL of the more concentrated standard and 0.5 mL of ISA to a second 50 mL beaker and swirl the solution to mix.
4. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the more concentrated standard. Wait for a stable reading and adjust the meter to display the value of the second standard, as described in the meter user guide.
5. Record the resulting slope value. The slope should be between -54 and -60 mV when the standards are between 20 and 25 °C.
6. Add 25 mL of sample and 0.5 mL of ISA to a clean 50 mL beaker and swirl the solution to mix.
7. Rinse the electrode with distilled water, blot it dry and place it into the sample. The concentration of the sample will be displayed on the meter.

Note: Other solution volumes may be used, as long as the ratio of solution to ISA remains 50:1.

Note: If interferences are present in the sample and cannot be removed, add 25 mL of NISS, Cat. No. 930710, per 25 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution. Other solution volumes may be used, as long as the ratio of solution to NISS remains 1:1.

Small Volume Direct Calibration Procedure Using a Meter with a mV Mode

Note: See the meter user guide for more specific information.

1. Set the meter to the mV mode.
2. Add 25 mL of the less concentrated standard and 0.5 mL of ISA to a 50 mL beaker and swirl the solution to mix.
3. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the less concentrated standard. When a stable reading is displayed, record the mV value and corresponding standard concentration.
4. Add 25 mL of the more concentrated standard and 0.5 mL of ISA to a second 50 mL beaker and swirl the solution to mix.
5. Rinse the electrode with distilled water, blot it dry and place it into the beaker with the more concentrated standard. When a stable reading is displayed, record the mV value and corresponding standard concentration.
6. Using semi-logarithmic graph paper, prepare a calibration curve by plotting the millivolt values on the linear axis and the standard concentration values on the logarithmic axis.
7. Add 25 mL of sample and 0.5 mL of ISA to a clean 50 mL beaker and swirl the solution to mix.
8. Rinse the electrode with distilled water, blot it dry and place it into the beaker. When a stable reading is displayed, record the mV value.
9. Using the calibration curve prepared in step 6, determine the unknown concentration of the sample.

Note: Other solution volumes may be used, as long as the ratio of solution to ISA remains 50:1.

Note: If interferences are present in the sample and cannot be removed, add 25 mL of NISS, Cat. No. 930710, per 25 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution. Other solution volumes may be used, as long as the ratio of solution to NISS remains 1:1.

Low Level Calibration Technique

These procedures are for solutions that have a nitrate concentration of less than 10^{-4} M or 1.4 ppm nitrate as N. For solutions low in nitrate but high in total ionic strength (greater than 10^{-1} M), perform the same procedure by preparing a calibrating solution with a composition similar to the sample.

Accurate results require that the following conditions be met:

- Prepare at least three calibration standards that bracket the expected sample concentration.
- Always use low level ISA for standards and samples. If interferences are present in the sample and cannot be removed, use NISS instead of low level ISA.
- Plastic labware must be used for all low level nitrate measurements.
- Adequate time must be allowed for electrode stabilization. Longer response time will be needed at low level measurements.
- Stir all standards and samples at a uniform rate.

Low Level Setup

1. Prepare the electrode as described in the **Electrode Preparation** section.
2. Connect the electrode to the meter. Set the meter to the mV mode.
3. Prepare the low level ISA by pipetting 20 mL of the nitrate ISA, Cat. No. 930711, into a 100 mL volumetric flask and diluting to the mark with distilled water. Use low level ISA for low level measurements only.

If interferences are present in the sample and cannot be removed, use NISS instead of low level ISA. Add 10.1 mL of NISS, Cat. No. 930710, per 90.9 mL of distilled water or sample.

4. Select a standard solution. Use either a 100 ppm nitrate as N or 10^{-3} M nitrate standard.

Low Level Calibration and Measurement

1. Add 100 mL of distilled water and 1 mL of low level ISA to a 150 mL beaker.
2. Rinse the electrode with distilled water, blot it dry and place it into the beaker. Stir the solution thoroughly.
3. Add increments of the 100 ppm or 10^{-3} M nitrate standard to the beaker using the steps outlined in **Table 3**. Record the stable millivolt reading after each increment.
4. On semi-logarithmic paper, plot the concentration (log axis) against the millivolt potential (linear axis). Prepare a new calibration curve with fresh standards each day.
5. Measure 100 mL of sample and 1 mL of low level ISA and pour the solutions into a clean 150 mL beaker. Rinse the electrode with distilled water, blot it dry and place the electrode into the sample.
6. Stir the solution thoroughly. When a stable reading is displayed, record the mV value.
7. Determine the sample concentration corresponding to the measured potential from the low level calibration curve.

Table 3

Calibration Curve For Low Level Calibrations

Additions of standard to 100 mL distilled water and 1 mL low level ISA solution.

| Step | Pipet Size | Volume Added | Concentration ppm as N | M |
|------|------------|--------------|------------------------|----------------------|
| 1 | 1 mL | 0.1 mL | 0.1 | 1.0×10^{-6} |
| 2 | 1 mL | 0.1 mL | 0.2 | 2.0×10^{-6} |
| 3 | 1 mL | 0.2 mL | 0.4 | 3.9×10^{-6} |
| 4 | 1 mL | 0.2 mL | 0.6 | 5.9×10^{-6} |
| 5 | 1 mL | 0.4 mL | 1.0 | 9.8×10^{-6} |
| 6 | 2 mL | 2.0 mL | 2.9 | 2.9×10^{-5} |
| 7 | 2 mL | 2.0 mL | 4.7 | 4.7×10^{-5} |

Known Addition Technique

Known addition is a convenient technique for measuring samples in the linear range of the electrode (greater than 10^{-4} M or 1.4 ppm nitrate as N) because no calibration curve is required. It can be used to verify the results of a direct calibration or to measure the total concentration of an ion in the presence of a large excess of a complexing agent. The sample potential is measured before and after addition of a standard solution.

Accurate results require that the following conditions be met:

- Concentration should approximately double as a result of the addition.
- Sample concentration should be known to within a factor of three.
- Either no complexing agent or a large excess of the complexing agent may be present.
- The ratio of the uncomplexed ion to complexed ion must not be changed by addition of the standard.
- All samples and standards should be at the same temperature.
- With double or multiple known addition, the final addition should be 10 to 100 times the sample concentration.
- Add 2 mL of ISA to every 100 mL of sample before analysis. If interferences are present in the sample and cannot be removed, add 50 mL of NISS, Cat. No. 930710, per 50 mL of standard or sample. Do not use ISA when using the nitrate interference suppressor solution.

Known Addition Setup

1. Prepare the electrode as described in the **Electrode Preparation** section.
2. Connect the electrode to the meter.
3. Prepare a standard solution that will cause the nitrate concentration of the sample to double when added to the sample solution. Refer to **Table 4** for guidelines.
4. Determine the electrode slope by performing the procedure in the **Checking Electrode Operation (Slope)** section.
5. Rinse the electrode with distilled water.

Table 4
Guideline For Known Addition

| Volume of Addition | Concentration of Standard |
|---------------------------|----------------------------------|
| 1 mL | 100 times sample concentration |
| 5 mL | 20 times sample concentration |
| 10 mL* | 10 times sample concentration |

* Most convenient volume to use

Known Addition Using a Meter with a Known Addition Mode

Note: See the meter user guide for more specific information.

1. Set the meter to measure in the known addition mode.
2. Measure 100 mL of the sample and 2 mL of ISA and pour the solutions into a beaker. Rinse the electrode with distilled water and place it into the sample solution. Stir the solution thoroughly.
3. When a stable reading is displayed, set the meter as described in the meter user guide, if required.
4. Pipet the appropriate amount of the standard solution into the beaker. Stir the solution thoroughly.
5. When a stable reading is displayed, record the sample concentration.

Known Addition Using a Meter with a Millivolt Mode

1. Set the meter to the relative millivolt mode. If a relative millivolt mode is not available, use the millivolt mode.
2. Measure 100 mL of sample and 2 mL of ISA and pour the solutions into a 150 mL beaker. Stir the solution thoroughly.
3. Rinse the electrode with distilled water, blot it dry and place the electrode into the beaker. When a stable reading is displayed, set the meter to read 0.0 mV. If the reading cannot be adjusted to 0.0 mV, record the actual mV value.
4. Pipet the appropriate amount of standard solution into the beaker. Stir the solution thoroughly.
5. When a stable reading is displayed, record the mV value. If the meter could not be set to 0.0 mV in step 3, subtract the first reading from the second reading to calculate ΔE .
6. Use **Table 6** to find the Q value that corresponds to the change in potential, ΔE . To determine the original sample concentration, multiply Q by the concentration of the added standard:

$$C_{\text{sample}} = Q * C_{\text{standard}}$$

C_{standard} = standard concentration

C_{sample} = sample concentration

Q = value from **Table 6**

The table of Q values is calculated for a 10% volume change. The equation for the calculation of Q for different slopes and volume changes is given below.

$$Q = (p * r) / \{(1 + p) * 10^{\Delta E/S} - 1\}$$

Q = value from **Table 6**

$\Delta E = E_2 - E_1$

S = slope of the electrode

p = volume of standard / volume of sample and ISA

r = volume of sample and ISA / volume of sample

Calculating Known Addition for Samples using Lotus, Excel, or Quattro Spreadsheets

If it is more convenient, a simple spreadsheet can be set up to calculate the known addition results, using any ratio of sample to addition. A typical worksheet is shown in **Table 5**. The numbers shown are examples, but the formulas and their locations should be copied exactly.

Table 5
Known Addition Calculations using Lotus, Excel, or Quattro Spreadsheets

| A | B | C |
|----|---|--------------------------------|
| 1 | | Enter Value |
| 2 | Volume of sample and ISA (mL) | 102 |
| 3 | Volume of addition (mL) | 10 |
| 4 | Concentration of addition | 10 |
| 5 | Volume of sample | 100 |
| 6 | Initial mV reading | 45.3 |
| 7 | Final mV reading | 63.7 |
| 8 | Electrode slope | -59.2 |
| 9 | | |
| 10 | | Derived Values |
| 11 | Delta E | +C7 - C6 |
| 12 | Solution volume ratio | +C3/C2 |
| 13 | Antilog term | +10 [^] (C11/C8) |
| 14 | Sample volume ratio | +C2/C5 |
| 15 | Q term | +C12*C14/ (((1+C12)*C13)-1) |
| 16 | Calculated initial concentration in same units as addition | +C15*C4 |

Note: For Excel, use = instead of + at start of formulas.

Table 6
Q Values for a 10% volume change,
slopes (in column heading) are in units of mV/decade

| ΔE | Q Concentration Ratio | | | |
|-----------|------------------------------|--------------|--------------|--------------|
| | -57.2 | -58.2 | -59.2 | -60.1 |
| 5.0 | 0.2917 | 0.2957 | 0.2996 | 0.3031 |
| 5.2 | 0.2827 | 0.2867 | 0.2906 | 0.2940 |
| 5.4 | 0.2742 | 0.2781 | 0.2820 | 0.2854 |
| 5.6 | 0.2662 | 0.2700 | 0.2738 | 0.2772 |
| 5.8 | 0.2585 | 0.2623 | 0.2660 | 0.2693 |
| 6.0 | 0.2512 | 0.2550 | 0.2586 | 0.2619 |
| 6.2 | 0.2443 | 0.2480 | 0.2516 | 0.2548 |
| 6.4 | 0.2377 | 0.2413 | 0.2449 | 0.2480 |
| 6.6 | 0.2314 | 0.2349 | 0.2384 | 0.2416 |
| 6.8 | 0.2253 | 0.2288 | 0.2323 | 0.2354 |
| 7.0 | 0.2196 | 0.2230 | 0.2264 | 0.2295 |
| 7.2 | 0.2140 | 0.2174 | 0.2208 | 0.2238 |
| 7.4 | 0.2087 | 0.2121 | 0.2154 | 0.2184 |
| 7.6 | 0.2037 | 0.2070 | 0.2102 | 0.2131 |
| 7.8 | 0.1988 | 0.2020 | 0.2052 | 0.2081 |
| 8.0 | 0.1941 | 0.1973 | 0.2005 | 0.2033 |
| 8.2 | 0.1896 | 0.1927 | 0.1959 | 0.1987 |
| 8.4 | 0.1852 | 0.1884 | 0.1914 | 0.1942 |
| 8.6 | 0.1811 | 0.1841 | 0.1872 | 0.1899 |
| 8.8 | 0.1770 | 0.1801 | 0.1831 | 0.1858 |
| 9.0 | 0.1732 | 0.1762 | 0.1791 | 0.1818 |
| 9.2 | 0.1694 | 0.1724 | 0.1753 | 0.1779 |
| 9.4 | 0.1658 | 0.1687 | 0.1716 | 0.1742 |
| 9.6 | 0.1623 | 0.1652 | 0.1680 | 0.1706 |
| 9.8 | 0.1590 | 0.1618 | 0.1646 | 0.1671 |
| 10.0 | 0.1557 | 0.1585 | 0.1613 | 0.1638 |
| 10.2 | 0.1525 | 0.1553 | 0.1580 | 0.1605 |
| 10.4 | 0.1495 | 0.1522 | 0.1549 | 0.1573 |
| 10.6 | 0.1465 | 0.1492 | 0.1519 | 0.1543 |
| 10.8 | 0.1437 | 0.1463 | 0.1490 | 0.1513 |
| 11.0 | 0.1409 | 0.1435 | 0.1461 | 0.1485 |
| 11.2 | 0.1382 | 0.1408 | 0.1434 | 0.1457 |
| 11.4 | 0.1356 | 0.1382 | 0.1407 | 0.1430 |
| 11.6 | 0.1331 | 0.1356 | 0.1381 | 0.1404 |
| 11.8 | 0.1306 | 0.1331 | 0.1356 | 0.1378 |
| 12.0 | 0.1282 | 0.1307 | 0.1331 | 0.1353 |
| 12.2 | 0.1259 | 0.1283 | 0.1308 | 0.1329 |
| 12.4 | 0.1236 | 0.1260 | 0.1284 | 0.1306 |
| 12.6 | 0.1214 | 0.1238 | 0.1262 | 0.1283 |
| 12.8 | 0.1193 | 0.1217 | 0.1240 | 0.1261 |
| 13.0 | 0.1172 | 0.1195 | 0.1219 | 0.1239 |
| 13.2 | 0.1152 | 0.1175 | 0.1198 | 0.1218 |
| 13.4 | 0.1132 | 0.1155 | 0.1178 | 0.1198 |
| 13.6 | 0.1113 | 0.1136 | 0.1158 | 0.1178 |
| 13.8 | 0.1094 | 0.1117 | 0.1139 | 0.1159 |

| ΔE | Q Concentration Ratio | | | |
|------------------------------|------------------------------|--------------|--------------|--------------|
| | -57.2 | -58.2 | -59.2 | -60.1 |
| 14.0 | 0.1076 | 0.1098 | 0.1120 | 0.1140 |
| 14.2 | 0.1058 | 0.1080 | 0.1102 | 0.1121 |
| 14.4 | 0.1041 | 0.1063 | 0.1084 | 0.1103 |
| 14.6 | 0.1024 | 0.1045 | 0.1067 | 0.1086 |
| 14.8 | 0.1008 | 0.1029 | 0.1050 | 0.1069 |
| 15.0 | 0.0992 | 0.1012 | 0.1033 | 0.1052 |
| 15.5 | 0.0953 | 0.0973 | 0.0994 | 0.1012 |
| 16.0 | 0.0917 | 0.0936 | 0.0956 | 0.0974 |
| 16.5 | 0.0882 | 0.0902 | 0.0921 | 0.0938 |
| 17.0 | 0.0850 | 0.0869 | 0.0887 | 0.0904 |
| 17.5 | 0.0819 | 0.0837 | 0.0856 | 0.0872 |
| 18.0 | 0.0790 | 0.0808 | 0.0825 | 0.0841 |
| 18.5 | 0.0762 | 0.0779 | 0.0797 | 0.0813 |
| 19.0 | 0.0736 | 0.0753 | 0.0770 | 0.0785 |
| 19.5 | 0.0711 | 0.0727 | 0.0744 | 0.0759 |
| 20.0 | 0.0687 | 0.0703 | 0.0719 | 0.0734 |
| 20.5 | 0.0664 | 0.0680 | 0.0696 | 0.0710 |
| 21.0 | 0.0642 | 0.0658 | 0.0673 | 0.0687 |
| 21.5 | 0.0621 | 0.0637 | 0.0652 | 0.0666 |
| 22.0 | 0.0602 | 0.0617 | 0.0631 | 0.0645 |
| 22.5 | 0.0583 | 0.0597 | 0.0612 | 0.0625 |
| 23.0 | 0.0564 | 0.0579 | 0.0593 | 0.0606 |
| 23.5 | 0.0547 | 0.0561 | 0.0575 | 0.0588 |
| 24.0 | 0.0530 | 0.0544 | 0.0558 | 0.0570 |
| 24.5 | 0.0514 | 0.0528 | 0.0541 | 0.0553 |
| 25.0 | 0.0499 | 0.0512 | 0.0525 | 0.0537 |
| 25.5 | 0.0484 | 0.0497 | 0.0510 | 0.0522 |
| 26.0 | 0.0470 | 0.0483 | 0.0495 | 0.0507 |
| 26.5 | 0.0456 | 0.0469 | 0.0481 | 0.0492 |
| 27.0 | 0.0443 | 0.0455 | 0.0468 | 0.0479 |
| 27.5 | 0.0431 | 0.0443 | 0.0455 | 0.0465 |
| 28.0 | 0.0419 | 0.0430 | 0.0442 | 0.0452 |
| 28.5 | 0.0407 | 0.0418 | 0.0430 | 0.0440 |
| 29.0 | 0.0395 | 0.0407 | 0.0418 | 0.0428 |
| 29.5 | 0.0385 | 0.0396 | 0.0407 | 0.0417 |
| 30.0 | 0.0374 | 0.0385 | 0.0396 | 0.0406 |
| 30.5 | 0.0364 | 0.0375 | 0.0385 | 0.0395 |
| 31.0 | 0.0354 | 0.0365 | 0.0375 | 0.0384 |
| 31.5 | 0.0345 | 0.0355 | 0.0365 | 0.0374 |
| 32.0 | 0.0335 | 0.0345 | 0.0356 | 0.0365 |
| 32.5 | 0.0327 | 0.0336 | 0.0346 | 0.0355 |
| 33.0 | 0.0318 | 0.0328 | 0.0337 | 0.0346 |
| 33.5 | 0.0310 | 0.0319 | 0.0329 | 0.0337 |
| 34.0 | 0.0302 | 0.0311 | 0.0320 | 0.0329 |
| 34.5 | 0.0294 | 0.0303 | 0.0312 | 0.0321 |
| 35.0 | 0.0286 | 0.0295 | 0.0305 | 0.0313 |
| 35.5 | 0.0279 | 0.0288 | 0.0297 | 0.0305 |
| 36.0 | 0.0272 | 0.0281 | 0.0290 | 0.0298 |

| ΔE | Q Concentration Ratio | | | |
|-----------|------------------------------|--------------|--------------|--------------|
| | -57.2 | -58.2 | -59.2 | -60.1 |
| 36.5 | 0.0265 | 0.0274 | 0.0282 | 0.0290 |
| 37.0 | 0.0258 | 0.0267 | 0.0275 | 0.0283 |
| 37.5 | 0.0252 | 0.0260 | 0.0269 | 0.0276 |
| 38.0 | 0.0246 | 0.0254 | 0.0262 | 0.0270 |
| 38.5 | 0.0240 | 0.0248 | 0.0256 | 0.0263 |
| 39.0 | 0.0234 | 0.0242 | 0.0250 | 0.0257 |
| 39.5 | 0.0228 | 0.0236 | 0.0244 | 0.0251 |
| 40.0 | 0.0223 | 0.0230 | 0.0238 | 0.0245 |
| 40.5 | 0.0217 | 0.0225 | 0.0232 | 0.0239 |
| 41.0 | 0.0212 | 0.0219 | 0.0227 | 0.0234 |
| 41.5 | 0.0207 | 0.0214 | 0.0221 | 0.0228 |
| 42.0 | 0.0202 | 0.0209 | 0.0216 | 0.0223 |
| 42.5 | 0.0197 | 0.0204 | 0.0211 | 0.0218 |
| 43.0 | 0.0192 | 0.0199 | 0.0206 | 0.0213 |
| 43.5 | 0.0188 | 0.0195 | 0.0202 | 0.0208 |
| 44.0 | 0.0183 | 0.0190 | 0.0197 | 0.0203 |
| 44.5 | 0.0179 | 0.0186 | 0.0192 | 0.0198 |
| 45.0 | 0.0175 | 0.0181 | 0.0188 | 0.0194 |
| 45.5 | 0.0171 | 0.0177 | 0.0184 | 0.0190 |
| 46.0 | 0.0167 | 0.0173 | 0.0179 | 0.0185 |
| 46.5 | 0.0163 | 0.0169 | 0.0175 | 0.0181 |
| 47.0 | 0.0159 | 0.0165 | 0.0171 | 0.0177 |
| 47.5 | 0.0156 | 0.0162 | 0.0168 | 0.0173 |
| 48.0 | 0.0152 | 0.0158 | 0.0164 | 0.0169 |
| 48.5 | 0.0148 | 0.0154 | 0.0160 | 0.0166 |
| 49.0 | 0.0145 | 0.0151 | 0.0157 | 0.0162 |
| 49.5 | 0.0142 | 0.0147 | 0.0153 | 0.0158 |
| 50.0 | 0.0139 | 0.0144 | 0.0150 | 0.0155 |
| 50.5 | 0.0135 | 0.0141 | 0.0146 | 0.0151 |
| 51.0 | 0.0132 | 0.0138 | 0.0143 | 0.0148 |
| 51.5 | 0.0129 | 0.0135 | 0.0140 | 0.0145 |
| 52.0 | 0.0126 | 0.0132 | 0.0137 | 0.0142 |
| 52.5 | 0.0124 | 0.0129 | 0.0134 | 0.0139 |
| 53.0 | 0.0121 | 0.0126 | 0.0131 | 0.0136 |
| 53.5 | 0.0118 | 0.0123 | 0.0128 | 0.0133 |
| 54.0 | 0.0116 | 0.0120 | 0.0125 | 0.0130 |
| 54.5 | 0.0113 | 0.0118 | 0.0123 | 0.0127 |
| 55.0 | 0.0110 | 0.0115 | 0.0120 | 0.0125 |
| 55.5 | 0.0108 | 0.0113 | 0.0118 | 0.0122 |
| 56.0 | 0.0106 | 0.0110 | 0.0115 | 0.0119 |
| 56.5 | 0.0103 | 0.0108 | 0.0113 | 0.0117 |
| 57.0 | 0.0101 | 0.0106 | 0.0110 | 0.0114 |
| 57.5 | 0.0099 | 0.0103 | 0.0108 | 0.0112 |
| 58.0 | 0.0097 | 0.0101 | 0.0105 | 0.0110 |
| 58.5 | 0.0095 | 0.0099 | 0.0103 | 0.0107 |
| 59.0 | 0.0093 | 0.0097 | 0.0101 | 0.0105 |
| 59.5 | 0.0091 | 0.0095 | 0.0099 | 0.0103 |
| 60.0 | 0.0089 | 0.0093 | 0.0097 | 0.0101 |


Electrode Characteristics

Electrode Response

The electrode potential plotted against concentration on semi-logarithmic paper results in a straight line with a slope of about -54 to -60 mV per decade change in concentration.

The time response of the electrode (the time required to reach 99% of the stable potential reading) varies from several seconds in concentrated solutions to several minutes near the limit of detection.

Figure 4
Typical Electrode Response to Nitrate Concentration



Reproducibility

Reproducibility is limited by factors such as temperature fluctuations, drift and noise. Within the operating range of the electrode, reproducibility is independent of concentration. With hourly calibrations, direct electrode measurements reproducible to $\pm 2\%$ can be obtained.

Limits of Detection

In pure nitrate solutions, the upper limit of detection is 1 M. When possible, dilute the sample into the linear range of the electrode. If samples are not diluted, the possibility of a liquid reference junction potential and the salt extraction effect, need to be considered. At high salt concentrations, salts may be extracted into the electrode membrane, causing deviation from theoretical response. To measure samples between 10^{-1} and 1 M, calibrate the electrode at 4 or 5 intermediate points or dilute the sample.

The lower limit of detection is determined by the slight water solubility of the ion exchanger, which causes deviation from theoretical response. **Figure 3** shows the theoretical response at low levels of nitrate compared to the actual response. If nitrate measurements are made below 10^{-4} M or 1.4 ppm nitrate as N, a low level measurement procedure is recommended.

Electrode Life

Each sensing module will last approximately three months with normal laboratory use, but the actual lifespan of the sensing module will depend on the type of samples that the electrode is used in. Refer to the **Electrode Maintenance** section for instructions on changing the sensing module. In time, the electrode slope will decrease and readings will start to drift, indicating that the module should be changed. Before replacement, refer to the **Troubleshooting** section to make sure that the difficulties are caused by the sensing module.

Temperature Effects

Since electrode potentials are affected by changes in temperature, samples and standard solutions should be within ± 1 °C (± 2 °F) of each other. At the 10^{-3} M level, a 1 °C difference in temperature results in errors greater than 1.5 %. The absolute potential of the reference electrode changes slowly with temperature because of the solubility equilibria on which the electrode depends. The slope of the electrode also varies with temperature, as indicated by the factor S in the Nernst equation. Theoretical values of the slope at different temperatures are given in **Table 7**. If the temperature changes, the meter and electrode should be recalibrated.

The electrode can be used at temperatures from 0 to 40 °C, provided that temperature equilibrium has occurred. For use at temperatures substantially different from room temperature, calibration standards should be at the same temperature as samples.

Table 7
Theoretical Slope vs. Temperature Values

| Temperature (°C) | Slope (mV) |
|------------------|------------|
| 0 | -54.20 |
| 10 | -56.18 |
| 20 | -58.16 |
| 25 | -59.16 |
| 30 | -60.15 |
| 40 | -62.13 |

If sample temperatures vary, use of the 9707BNWP combination nitrate electrode is recommended. The Optimum Results F filling solution that is included with the electrode will minimize junction potentials and provide optimum temperature and time response. Optimum Results F filling solution produces an isopotential point of 3.2×10^{-3} M nitrate. The isopotential point is the concentration at which the potential of the electrode does not vary with temperature. Since the isopotential point of this electrode is known, the combination nitrate electrode may be used on meters that allow automatic temperature compensation for ISE measurements. By programming in the isopotential point and placing an ATC probe into the sample, any time the temperature changes the meter will automatically adjust the slope of the calibration curve, resulting in more accurate measurements.

Interferences

Some anions, if present at high enough levels, are electrode interferences and will cause measurement errors. **Table 8** indicates levels of common ions that will cause 10% errors at different concentrations of nitrate.

The nitrate interference suppressor solution (NISS), Cat. No. 930710, is recommended for the removal of a variety of interfering anions present in samples such as soils, drinking water, wastewater and plant tissues. The nitrate interference suppressor solution is mixed in an equal volume with samples and standards. For example, add 50 mL of NISS per 50 mL of standard or sample. This procedure ensures that samples and standards have a similar background and that no correction factor is needed for the dilution. Do not use ISA when using the nitrate interference suppressor solution.

If the electrode is exposed to high levels of interfering ions, it may drift and become sluggish in response. When this happens, restore normal performance by soaking the electrode for an hour in distilled water, emptying the old filling solution, filling the electrode with fresh filling solution and then soaking the electrode for a few hours in a 10^{-2} M or 100 ppm nitrate standard. If soaking the electrode does not restore normal electrode performance, refer to the **Electrode Maintenance** section for instructions on how to replace the sensing module.

When the level of interferences in samples is constant, it is sometimes possible to measure nitrate accurately when interference levels are higher than those in **Table 8**. For example, nitrate can be measured in sea water by using synthetic ocean water for calibration. Contact our Technical Support Chemists for more information.

Table 8
Nitrate Electrode Interferences

| Interferences Moles/Liter | 10⁻⁴ M Nitrate | 10⁻³ M Nitrate | 10⁻² M Nitrate |
|---|----------------------------------|----------------------------------|----------------------------------|
| (d) ClO ₄ ⁻ | 1 × 10 ⁻⁸ | 1 × 10 ⁻⁷ | 1 × 10 ⁻⁶ |
| (b) I ⁻ | 5 × 10 ⁻⁷ | 5 × 10 ⁻⁶ | 5 × 10 ⁻⁵ |
| (d) ClO ₃ ⁻ | 5 × 10 ⁻⁶ | 5 × 10 ⁻⁵ | 5 × 10 ⁻⁴ |
| (b) CN ⁻ | 1 × 10 ⁻⁵ | 1 × 10 ⁻⁴ | 1 × 10 ⁻³ |
| (b) Br | 7 × 10 ⁻⁵ | 7 × 10 ⁻⁴ | 7 × 10 ⁻³ |
| (c) NO ₂ ⁻ | 7 × 10 ⁻⁵ | 7 × 10 ⁻⁴ | 7 × 10 ⁻³ |
| (b) HS ⁻ | 1 × 10 ⁻⁴ | 1 × 10 ⁻³ | 1 × 10 ⁻² |
| (a) HCO ₃ ⁻ | 1 × 10 ⁻³ | 1 × 10 ⁻² | 0.1 |
| (a) CO ₃ ⁻² | 2 × 10 ⁻³ | 2 × 10 ⁻² | 0.2 |
| (b) Cl ⁻ | 3 × 10 ⁻³ | 3 × 10 ⁻² | 0.3 |
| (b) H ₂ PO ₄ ⁻ | 5 × 10 ⁻³ | 5 × 10 ⁻² | 0.5 |
| (b) HPO ₄ ⁻² | 5 × 10 ⁻³ | 5 × 10 ⁻² | 0.5 |
| (b) PO ₄ ⁻³ | 5 × 10 ⁻³ | 5 × 10 ⁻² | 0.5 |
| (e) OAc ⁻ | 2 × 10 ⁻² | 0.2 | 2 |
| F ⁻ | 6 × 10 ⁻² | 0.6 | 6 |
| SO ₄ ⁻² | 0.1 | 1.0 | 10 |

| Interferences ppm | 1 ppm Nitrate as N | 10 ppm Nitrate as N | 100 ppm Nitrate as N |
|---|---------------------------|----------------------------|-----------------------------|
| (d) ClO ₄ ⁻ | 7 × 10 ⁻⁴ | 7 × 10 ⁻³ | 7 × 10 ⁻² |
| (b) I ⁻ | 4 × 10 ⁻² | 0.4 | 4 |
| (d) ClO ₃ ⁻ | 0.3 | 3 | 30 |
| (b) CN ⁻ | 0.2 | 2 | 20 |
| (b) Br | 4 | 40 | 400 |
| (c) NO ₂ ⁻ | 2 | 23 | 230 |
| (b) HS ⁻ | 2 | 23 | 230 |
| (a) HCO ₃ ⁻ | 44 | 440 | 4400 |
| (a) CO ₃ ⁻² | 86 | 860 | 8600 |
| (b) Cl ⁻ | 76 | 760 | 7600 |
| (b) H ₂ PO ₄ ⁻ | 346 | 3464 | 34640 |
| (b) HPO ₄ ⁻² | 343 | 3430 | 34300 |
| (b) PO ₄ ⁻³ | 339 | 3390 | 33900 |
| (e) OAc ⁻ | 1042 | 10420 | 104200 |
| F ⁻ | 814 | 8140 | 81400 |
| SO ₄ ⁻² | 6857 | 68570 | 685700 |


- (a) Carbonate and bicarbonate can be removed by acidifying the sample to pH 4.5 with sulfuric acid, converting the ions to carbon dioxide.
- (b) These interferences can be minimized by precipitation with silver. Dissolve solid silver sulfate in samples to remove.
- (c) Nitrite can be removed by adding sufficient sulfamic acid to samples.
- (d) These interferences cannot be removed. Use the Thermo Scientific Orion nitrate test kit, Cat. No. 700005, to convert nitrate to ammonia. Measure samples using the ammonia electrode, Cat. No. 9512HPBNWP.
As an alternate method, convert nitrate to nitrite with a reduction column and measure nitrite levels with a nitrite electrode, Cat. No. 9746BNWP.
For more information contact our Technical Support Chemists.
- (e) Many organic (carboxylic) anions also interfere with the nitrate electrode. These anions can be removed by using a 1 M ISA solution.

Note: *The use of any of the above procedures require similar treatment of standards as well as samples.*

Theory of Operation

The nitrate electrode consists of a replaceable, pretested sensing module connected to an epoxy body. The sensing module contains a liquid internal filling solution in contact with a gelled organophilic membrane that contains a nitrate selective ion exchanger.

Figure 5
Example of an Ion Sensing Module



When the module is in contact with a solution containing nitrate ions, an electrode potential develops across the module. This potential, which depends on the level of free nitrate ion in solution, is measured against a constant reference potential with a digital pH/mV meter or ISE (concentration) meter. The measured potential corresponding to the level of nitrate ion in solution is described by the Nernst equation.

$$E = E_o + S * \log (A)$$

E = measured electrode potential

E_o = reference potential (a constant)

A = nitrate ion activity level in solution

S = electrode slope (about -57 mV per decade)

$S = (2.3 RT) / nF$

R and F are constants, T = temperature in degrees K and

n = ionic charge

The level of nitrate ions, A, is the activity or “effective concentration” of free nitrate ions in solution. The nitrate ion activity is related to free nitrate ion concentration, C_f , by the activity coefficient, γ .

$$A = \gamma * C_f$$

Ionic activity coefficients are variable and largely depend on total ionic strength. The ionic strength of a solution is determined by all of the ions present. It is calculated by multiplying the concentration of each individual ion by the square of its charge, adding all these values up and then dividing by two.

$$\text{Ionic strength} = 1/2 \sum (C_i Z_i^2)$$

C_i = concentration of ion i

Z_i = charge of ion i

\sum symbolizes the sum of all the types of ions in solutions

If background ionic strength is high and constant relative to the sensed ion concentration, the activity coefficient is constant and activity is directly proportional to concentration. Ionic strength adjustor (ISA) is added to all nitrate standards and samples so that the background ionic strength is high and constant relative to variable concentrations of nitrate. For nitrate, the recommended ISA is $(\text{NH}_4)_2\text{SO}_4$. Nitrate interference suppression solution (NISS), a specific solution for removal of nitrate-interfering ions, is recommended for samples with competing ions. Other solutions can be used as long as they do not contain ions that would interfere with the electrode response to nitrate.

If samples have a high ionic strength (above 0.1 M), standards should be prepared with a composition similar to the samples.

Reference electrode conditions must also be considered. Liquid junction potentials arise any time when two solutions of different composition are brought into contact. The potential results from the interdiffusion of ions in the two solutions. Since ions diffuse at different rates, the electrode charge will be carried unequally across the solution boundary resulting in a potential difference between the two solutions. In making electrode measurements, it is important that this potential is the same when the reference is in the standardizing solution as well as in the same solution; otherwise, the change in liquid junction potential will appear as an error in the measured specific ion electrode potential.

The most important variable that analysts have under their control is the composition of the liquid junction filling solution. The filling solution should be equitransferent. That is, the speed with which the positive and negative ions in the filling solution diffuse into the sample should be nearly as equal as possible. If the rate at which positive and negative charge is carried into the sample solution is equal, then no junction potential can result. Optimum Results filling solutions are specifically designed to meet all reference electrode conditions.

Troubleshooting

Follow a systematic procedure to isolate the problem. The measuring system can be divided into four components for ease in troubleshooting: meter, electrode, sample/application and technique.

Meter

The meter is the easiest component to eliminate as a possible cause of error. Thermo Scientific Orion meters include an instrument checkout procedure and shorting cap for convenience in troubleshooting. Consult the meter user guide for directions.

Electrode

1. Rinse the electrode thoroughly with distilled water.
2. Verify the electrode performance by performing the procedure in the **Checking Electrode Operation (Slope)** section.
3. If the electrode fails this procedure, review the **Measuring Hints** section. Clean the electrode thoroughly as directed in the **Electrode Maintenance** section. Drain and refill the electrode with fresh filling solution.
4. Repeat the procedure in the **Checking Electrode Operation (Slope)** section.
5. If the electrode fails this procedure again and the half-cell nitrate electrode is being used, determine whether the nitrate or reference electrode is at fault. To do this, substitute a known working electrode for the electrode in question and repeat the procedure in the **Checking Electrode Operation (Slope)** section.
6. If the electrode passes the procedure, but measurement problems persist, the sample may contain interferences or complexing agents, or the technique may be in error.
7. Before replacing a faulty electrode, review this user guide and be sure to thoroughly clean the electrode; correctly prepare the electrode; use the proper filling solution, ISA or NISS and standards; correctly measure the samples and review the **Troubleshooting Checklist** section.

Sample/Application

The quality of results depends greatly upon the quality of the standards. Always prepare fresh standards when problems arise, it could save hours of frustrating troubleshooting! Errors may result from contamination of prepared standards, accuracy of dilution, quality of distilled water, or a mathematical error in calculating the concentrations.

The best method for preparation of standards is serial dilution. Refer to the **Serial Dilution** section. The electrode and meter may operate with standards, but not with the sample. In this case, check the sample composition for interferences, incompatibilities or temperature effects. Refer to the **Sample Requirements, Temperature Effects** and **Interferences** sections.

Technique

If trouble persists, review operating procedures. Review calibration and measurement sections to be sure proper technique has been followed. Verify that the expected concentration of the ion of interest is within the limit of detection of the electrode.

Check the method of analysis for compatibility with your sample. Direct measurement may not always be the method of choice. If a large amount of complexing agents are present, known addition may be the best method. If working with low level samples, follow the procedure in the **Low Level Calibration** section.

Assistance

After troubleshooting all components of your measurement system, contact Technical Support. Within the United States call 1.800.225.1480 and outside the United States call 978.232.6000 or fax 978.232.6031. In Europe, the Middle East and Africa, contact your local authorized dealer. For the most current contact information, visit www.thermo.com/contactwater.

For the latest application and technical resources for Thermo Scientific Orion products, visit www.thermo.com/waterapps.

Warranty

For the most current warranty information, visit www.thermo.com/water.

Troubleshooting Checklist

- No electrode filling solution added –
Fill the electrode with filling solution up to the fill hole. Refer to the **Electrode Preparation** section for details.
- Incorrect electrode filling solution used –
Refer to the **Electrode Preparation** section to verify that the correct electrode filling solution was used.
- Electrode junction is dry –
Push down on the electrode cap to allow a few drops of filling solution to drain out of the electrode.
- No reference electrode present –
The 9307BNWP nitrate half-cell electrode require a separate reference electrode, Cat. No. 900200.
- Electrode is clogged or dirty –
Refer to the **Electrode Maintenance** section for electrode cleaning and flushing instructions.
- Sensing module is not installed properly, dirty or defective –
Refer to the **Electrode Preparation** section and verify that the electrode was assembled correctly. Refer to the **Electrode Maintenance** section for instructions on installing a new sensing module.
- Standards are contaminated or made incorrectly –
Prepare fresh standards. Refer to the **Serial Dilution**, **Measurement Hints** and **Analytical Techniques** sections.
- ISA not used or incorrect ISA used –
ISA must be added to all standards and samples. Refer to the **Required Equipment** section for information on the ISA.
- Interferences present –
Use the nitrate interference suppressor solution (NISS), Cat. No. 930710, instead of ISA.
- Samples and standards at different temperatures –
Allow solutions to reach the same temperature.
- Air bubble on sensing module –
Remove air bubble by reimmersing the electrode in solution.
- Electrode not properly connected to meter –
Unplug and reconnect the electrode to the meter.
- Meter or stir plate not properly grounded –
Check the meter and stir plate for proper grounding.
- Static electricity present –
Wipe plastic parts on the meter with a detergent solution.
- Defective meter –
Check the meter performance. See the meter user guide.

Ordering Information

| Cat. No. | Description |
|-----------------|--|
| 9707BNWP | Nitrate ionplus Sure-Flow combination electrode, waterproof BNC connector |
| 900046 | Optimum Results F electrode filling solution, 5 x 60 mL bottles |
| 9307BNWP | Nitrate half-cell electrode, waterproof BNC connector (requires separate reference electrode) |
| 900200 | Double junction reference electrode, pin tip connector |
| 900002 | Inner chamber filling solution for the double junction reference electrode, 5 x 60 mL bottles |
| 920706 | 0.1 M NaNO ₃ nitrate calibration standard |
| 920707 | 1000 ppm as N nitrate calibration standard |
| 930707 | 100 ppm as N nitrate calibration standard |
| 930711 | ISA for nitrate measurements, 475 mL bottle |
| 930710 | Nitrate interference suppressor solution (NISS), 475 mL bottle |
| 970701 | 97 series sensing module for 9707BNWP nitrate combination electrode |
| 9700BNWP | Replacement body for 9707BNWP nitrate combination electrode (requires separate 97 series sensing module) |
| 930701 | 93 series sensing module for 9307BNWP nitrate half-cell electrode, pack of 3 |
| 930702 | 93 series sensing module for 9307BNWP nitrate half-cell electrode |
| 9300BNWP | Replacement body for 9307BNWP nitrate half-cell electrode (requires separate 93 series sensing module) |

Specifications

Concentration Range

7×10^{-6} M to 1 M (0.1 ppm to 14,000 ppm nitrate as N)

pH Range

2.5 to 11

Low level measurements may be influenced by hydrogen or hydroxide ion interferences

Temperature Range

0 to 40 °C

Electrode Resistance

0.1 to 5 megohms

Reproducibility

± 2%

Minimum Sample Size (9707BNWP)

5 mL in a 50 mL beaker

Size— 9707BNWP

Body Diameter: 13 mm

Body Length: 110 mm

Cap Diameter: 16 mm

Cable Length: 1 meter

Size— 9307BNWP

Body Diameter: 12 mm

Body Length: 110 mm

Cap Diameter: 16 mm

Cable Length: 1 meter

** Specifications are subject to change without notice*

258510-001 Rev. A 12-08

Thermo Fisher Scientific

Environmental Instruments
Water Analysis Instruments



North America

166 Cummings Center
Beverly, MA 01915 USA
Toll Free: 1-800-225-1480
Tel: 1-978-232-6000
Dom. Fax: 1-978-232-6015
Int'l Fax: 978-232-6031

Europe

P.O. Box 254, 3860 AG Nijkerk
Wallerstraat 125K, 3862 BN
Nijkerk, Netherlands
Tel: (31) 033-2463887
Fax: (31) 033-2460832

Asia Pacific

Blk 55, Ayer Rajah Crescent
#04-16/24, Singapore 139949
Tel: 65-6778-6876
Fax: 65-6773-0836

www.thermo.com/water

Thermo
SCIENTIFIC