

Pro-Detect™ Rapid AVI Competitive Assay Kit

Catalog Numbers A38511

Pub. No. MAN0017619 Rev. A.0

WARNING! Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves. Safety Data Sheets (SDSs) are available from thermofisher.com/support.

Contents and storage

Product	Cat. No.	Contents	Storage
Pro-Detect™ Rapid AVI Competitive Assay Kit	A38511	Pro-Detect™ Rapid AVI Competitive Assay Strips, 10 strips Pro-Detect™ Rapid Assay Dilution Buffer, 15 mL	Store at 4°C. Do not freeze. After opening, store unused strips in the enclosed container containing desiccant.

Product description

The Pro-Detect™ Rapid AVI Competitive Assay Kit is a 10-minute dipstick lateral-flow assay to detect AVI-tagged proteins in cell culture and lysates during protein expression or in purified protein preparations. The assay is performed by simply applying the lateral flow strip into the sample of properly diluted tissue culture supernatant or lysate and visualizing via a loss of red bands in the test line section. A red control line will appear in both positive and negative results. In the competitive lateral flow assay, AVI-tagged protein antigens immobilized on the membrane as 3 distinct test lines. A control antibody is immobilized as the control line. The gold-conjugated capture antibodies, specific to AVI, are embedded in the sample pad (Fig. A).

In a negative result where no AVI-tagged protein is present or concentration of AVI-tagged protein is below detectable levels, the gold-conjugated capture antibodies will bind to the AVI-tagged protein antigen embedded on the test lines and form 3 visible red lines. In both positive and negative tests, the gold-conjugated anti-AVI capture antibodies will bind to the control antibodies at the control line (Fig. B).

In a positive test where the sample contains AVI-tagged proteins, the gold-conjugated antibodies embedded in the sample pad will bind to the available AVI-tagged proteins in the sample and therefore not bind to the AVI-tagged protein antigens immobilized on the test lines. As the concentration increases, the number of test lines will decrease until all test lines disappear. The concentration of the AVI-tagged proteins is inversely related to the number of test lines appearing on the strip (Fig. C).

Material preparation


Because the competitive lateral flow assays structure has an inverse relationship with signal loss, there is no true upper limit for the assay; however, dilution is recommended in efforts to preserve samples. The recommended working range of the AVI competitive lateral flow assay is 4 µg/mL to 20 µg/mL. The competitive lateral flow assays can detect down to 1 µg/mL, but loss of band intensity may be less distinguished. As there is no upper limit, the lateral flow strip will continue to show positive results at concentrations greater than 20 µg/mL, visualized by complete loss of test lines (test lines closest to the sample will disappear first). It is not recommended to use at concentrations less than 1 µg/mL.

Note: Proper sample dilution is essential for optimal results. Concentration ranges are based upon the concentration of the tagged protein of interest in the sample.

Sample dilution

If sample concentration is known, dilute the sample with Pro-Detect™ Rapid Assay Dilution Buffer to a concentration of 10 µg/mL. For unknown starting sample concentrations, dilute mammalian cell lysates at 1:4 dilution and bacterial lysates at 1:20.

For example, to make 150 µL at a 1:20 dilution, add 7.5 µL of sample (cell supernatant or lysate) to 142.5 µL of Pro-Detect™ Rapid Assay Dilution Buffer and vortex or pipette up and down to mix.



Perform test

Note: Perform all test at room temperature. To avoid condensation on the strips, allow the package to warm to room temperature for 15 minutes prior to removing the strips from the bag.

1. For each strip, add 150 µL of diluted sample to a microtiter plate or test tube.
2. Insert the lateral flow strip with arrow facing downward into the sample and wait 10-15 minutes for the color bands to appear.
3. Remove the strip(s) from the sample after 10-15 minutes (20 mins max.) of processing time.

A positive test will result in a red band at the control location and at the test location. Images may be acquired by photograph (camera, cell phone) or imaging equipment. Lateral flow assays may be further saved by placing in notebooks.

After results are obtained, if signal appears weak, perform a lower dilution. A weak signal can result from a high concentration, resulting in a weaker signal. The additional dilution will bring the concentration into the recommended working range. If upon detection with the second dilution no test line is observed, the original test performed was at a concentration below the recommended working range. When the concentration of actual tagged protein is not known in mammalian cell lysate, dilute 1:2; for a bacterial lysate, dilute 1:5.

Troubleshooting

Observation	Possible cause	Recommended action
All 3 test lines detected after applying sample.	Sample did not contain protein tag.	Verify correct assay strip is used.
		Verify presence of protein tag via alternative method (e.g., ELISA or Western blot).
Low intensity test lines.	Sample was below recommended concentration range.	Dilute samples as indicated in the Perform test section and re-test using a new lateral flow strip.
No control line detected.	Lateral flow strip was not sufficiently submerged in the sample.	Insert test strip fully into sample well and ensure enough volume is present in the sample well to fully cover the white application tip of the lateral flow strip.
	Sample contained an interfering substance.	Confirm additional lysis and extraction reagents are within the recommended ranges (see the Additional Information section).
High background detected.	Sample contained an interfering substance.	Confirm additional lysis and extraction reagents are within the recommended ranges (see the Additional Information section).
	Lateral flow strip was left in sample longer than 15 minutes.	Over-processing can result in higher background. Evaluate the results within 10-15 minutes and remove the strip from the sample after 15 minutes.

Additional information

Lateral flow assays are highly robust assays that can withstand many commonly used detergents, buffers, salts, and other lysis reagents. A list of commonly used reagents and the effective compatible concentrations are provided in the table below.

Table 1 Assay reagent compatibility.

Salts/Buffers		Detergents	
Substance	Compatible Concentration	Substance	Compatible Concentration
NaCl	0.25 M	SDS	0.2%
Urea	0.4 M	Triton™ X-100 Detergent	1%
RIPA Buffer	Undiluted	CHAPS Detergent	1%
B-PER™ Bacterial Protein Extraction Reagent	Undiluted	Misc. Reagents and Solvents	
M-PER™ Mammalian Protein Extraction Reagent	Undiluted	Glycerol	10%
KCl	0.25 M	NP40 Detergent	1%
		EDTA	5 mM



Manufacturer: Thermo Fisher Scientific | 3747 N. Meridian Road | Rockford, Illinois 61101 USA

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.