

HisPur™ Cobalt Purification Kit

90090 90091 90092

2060.3

Number	Description
90090	HisPur Cobalt Purification Kit, 0.2mL Kit Contents: HisPur Cobalt Spin Columns , 25 each, pre-filled with 0.2mL resin bed Equilibration/Wash Buffer , 100mL Elution Buffer , 50mL Collection Tubes , 80 each Tip Plugs , 25 each
90091	HisPur Cobalt Purification Kit, 1mL Kit Contents: HisPur Cobalt Spin Columns , 5 each, pre-filled with 1mL resin bed Equilibration/Wash Buffer , 100mL Elution Buffer , 50mL Tip Plugs , 5 each
90092	HisPur Cobalt Purification Kit, 3mL Kit Contents: HisPur Cobalt Spin Columns , 5 each, pre-filled with 3mL resin bed Equilibration/Wash Buffer , 225mL Elution Buffer , 75mL Tip Plugs , 5 each

Binding Capacity: ≥ 10 mg at $> 90\%$ purity of a 28kDa His-tagged protein from a bacterial source per milliliter of resin bed

Resin: Crosslinked 6% agarose in a 20% ethanol solution

Equilibration/Wash Buffer: 50mM sodium phosphate, 300mM sodium chloride, 10mM imidazole; pH 7.4.

Elution Buffer: 50mM sodium phosphate, 300mM sodium chloride, 150mM imidazole; pH 7.4

Note: Buffers contain 0.01% sodium azide as an antimicrobial.

Storage: Upon receipt store kit at 4°C. Product shipped at ambient temperature.

Introduction

The Thermo Scientific HisPur Cobalt Purification Kit contains HisPur Cobalt Spin Columns and pre-formulated Equilibration/Wash and Elution Buffers for the efficient purification of recombinant polyhistidine-tagged proteins from bacterial, mammalian and baculovirus-infected cells. His-tagged proteins are purified from total soluble protein extract using a cobalt-charged tetradentate chelator immobilized onto 6% crosslinked agarose. Many immobilized metal affinity chromatography (IMAC) resins contain nickel (Ni^{2+}) for purifying His-tagged proteins. Although Ni^{2+} chelate resins achieve high protein yields, purity is often suboptimal, requiring additional clean-up steps. Cobalt achieves both high yield and purity with minimal optimization. The HisPur Cobalt Resin binds fewer nonspecific proteins, displays less metal leaching and enables less stringent elution conditions compared to Ni^{2+} resins.

Important Product Information

- Typical binding capacity is ≥ 10 mg protein/mL resin. Depending upon expression of the protein, the maximum total protein (lysate) loading amount recommended is 8mg for the 0.2mL spin column, 40mg for the 1mL spin column and 120mg for the 3mL spin column. Typical yields are 10-25% of the total protein loaded onto the column. For optimal results, do not exceed the capacity of the resin.
- Avoid using chelator-containing protease inhibitors or other additives. EDTA and strong reducing agents, such as DTT and β -mercaptoethanol, disrupt the function of the cobalt resin.
- Optimization of the lysis procedure is critical for maximizing protein yield. Some methods for protein extraction include using commercially available detergent-based reagents, such as Thermo Scientific B-PER Bacterial Protein Extraction Reagent in Phosphate Buffer (Product No. 78266), and mechanical methods, such as freeze/thaw cycles, sonication or French press.

Additional Materials Required

- Regeneration MES Buffer: 20mM 2-(*N*-morpholine)-ethanesulfonic acid, 0.1M sodium chloride; pH 5.0
- EDTA-free protease inhibitors such as Thermo Scientific Halt Protease Inhibitor Single Use Cocktail, EDTA-Free (Product No 78425)

Procedure for Spin Purification of His-tagged Protein

Note: The total volume of the 0.2mL, 1mL and 3mL column devices are 1.0mL, 8mL and 22mL, respectively. If a sample volume is greater than the column device, perform multiple applications and centrifugations until the entire sample has been processed. Be careful not to exceed the working capacity of the resin. The HisPur Cobalt Spin Columns also can be used for gravity-flow purifications.

1. Equilibrate column(s) to working temperature. Perform purifications at room temperature or at 4°C.
2. Remove the bottom tab from the HisPur Cobalt Spin Column. Place column into a centrifuge tube and loosen top cap.
Note: Use 1.5mL, 15mL or 50mL centrifuge tubes for the 0.2mL, 1mL and 3mL spin columns, respectively.
3. Centrifuge column at $700 \times g$ for 2 minutes to remove storage buffer.
4. Equilibrate column by adding two resin-bed volumes of Equilibration/Wash Buffer. Invert column to mix.
5. Centrifuge column at $700 \times g$ for 2 minutes to remove buffer.
6. Prepare sample by mixing 1 part protein extract with 1 part Equilibration/Wash Buffer.
7. Cap bottom of column. Add the prepared protein extract to the tube and cap top of column. Invert column to mix.
Note: For maximal binding, incubate for 30 minutes at 4°C on an end-over-end rocking platform.
8. Centrifuge column at $700 \times g$ for 2 minutes and collect the flow-through in a new centrifuge tube.
9. Wash resin with two resin-bed volumes of Equilibration/Wash Buffer. Centrifuge at $700 \times g$ for 2 minutes and collect fraction in a centrifuge tube. Repeat this step two more times collecting each fraction in a separate centrifuge tube.
Note: If desired, perform additional washes. Monitor washes by measuring their absorbance at 280nm.
10. Elute His-tagged proteins from the resin by adding one resin-bed volume of Elution Buffer. Centrifuge at $700 \times g$ for 2 minutes. Repeat this step two more times, collecting each fraction in a separate tube.
Note: If performing gravity-flow add two resin bed volumes of Elution Buffer to achieve proper flow characteristics.
11. Monitor protein elution by measuring the absorbance of the fractions at 280nm or by Coomassie Plus (Bradford) Protein Assay (Product No. 23236). The eluted protein can be directly analyzed by SDS-PAGE. To remove excess imidazole for downstream applications, use gel filtration (e.g., Thermo Scientific Zeba Spin Desalting Columns) or dialysis (e.g., Thermo Scientific Slide-A-Lyzer Dialysis Cassettes).

Procedure for Column Regeneration

The cobalt spin column can be used multiple times without affecting protein yield or purity. After each use and before storing, perform the procedure as described below to remove residual imidazole and any nonspecifically adsorbed protein. To prevent cross contamination of samples, designate a given column to one specific fusion protein.

1. Wash column with 10 column volumes of Regeneration MES Buffer.
2. Wash column with 10 column volumes of ultrapure water. For storage, apply 20% ethanol to the column.
3. Before reuse, re-equilibrate the column with Equilibration/Binding Buffer until the pH returns to the buffer value.

Troubleshooting

Problem	Possible Cause	Solution
Low protein yield	Poor expression of soluble protein	Optimize bacterial expression conditions
	His-tagged protein forms inclusion bodies	Alter bacterial growth conditions to minimize inclusion body formation and maximize soluble protein yield
		Solubilize inclusion bodies and perform the purification with a compatible denaturant (e.g., Inclusion Body Solubilization Reagent, Product No. 78115)
	Insufficient cell lysis and extraction	Optimize cell lysis protocol
	His tag is absent	Verify the sequence or perform an ELISA or Western blot using an antibody against the His tag
	His tag is inaccessible using native conditions	See the Additional Information Section for denaturing conditions
	His-tagged protein has a low affinity to the cobalt column	Optimize the equilibration/wash buffer by decreasing the concentration of imidazole (i.e., 50mM sodium phosphate, 300mM sodium chloride, 0-10mM imidazole; pH 7.4)
Poor protein purity	Insufficient column washing	Wash column additional times or modify the imidazole concentration; alternatively, adjust the Equilibration/Wash Buffer to pH 7.0 to decrease nonspecific protein binding
Slow column flow	Column is overloaded	Apply less protein extract to the column and make sure the extract is not too viscous or contaminated with cell debris

Additional Information

A. Fusion Proteins Expressed in Inclusion Bodies

Over-expressed proteins are sometimes sequestered in inclusion bodies. Inclusion bodies can be solubilized in 8M urea, 6M guanidine or the Thermo Scientific Inclusion Body Solubilization Reagent (Product No. 78115); however, a denaturant must be added to the buffers to ensure the protein remains soluble throughout the procedure. (Follow the Procedure for Spin Purification of His-tagged Protein.)

If using 8M urea, you may proceed directly to SDS-PAGE without sample clean-up. If using 6M guanidine, perform one of the following steps before SDS-PAGE: 1) Dilute sample 1:6 in ultrapure water; 2) Dialyze against a compatible buffer; 3) Perform TCA precipitation; 4) Use the Thermo Scientific Pierce SDS-PAGE Sample Prep Kit (Product No. 89888).

For denaturing conditions prepare the following buffers:

- Equilibration/Wash Buffer: 50mM sodium phosphate, 300mM sodium chloride, 6M guanidine•HCl, 0-10mM imidazole; pH 7.4
- Elution Buffer: 50mM sodium phosphate, 300mM sodium chloride, 6M guanidine•HCl, 150mM imidazole; pH 7.4

B. Tech Tips Available from our Website

- Tech Tip #43: Protein stability and storage
- Tech Tip #6: Extinction coefficients guide
- Tech Tip #40: Convert between times gravity ($\times g$) and centrifuge rotor speed (RPM)

Related Thermo Scientific Products

88270	Pierce [®] High Capacity Endotoxin Removal Resin, 10mL
88282	Pierce LAL Chromogenic Endotoxin Quantitation Kit
16100	Pierce Glutathione Agarose, 10mL
78425	Halt [™] Protease Inhibitor Single Use Cocktail, EDTA Free, 24 \times 100 μ L microtubes, each 100 μ L microtube contains sufficient cocktail to treat 10mL of lysate
88221	HisPur Ni-NTA Resin, 10mL
89967	HisPur Cobalt Spin Columns, 0.2mL resin bed, 25/pkg
89968	HisPur Cobalt Spin Columns, 1.0mL resin bed, 5/pkg
89969	HisPur Cobalt Spin Columns, 3.0mL resin bed, 5/pkg
89964	HisPur Cobalt Resin, 10mL settled resin
78266	B-PER [®] Bacterial Protein Extraction Reagent (in Phosphate Buffer), 500mL
78248	B-PER Bacterial Protein Extraction Reagent, 500mL
78260	B-PER II Bacterial Protein Extraction Reagent, 250mL
89802	I-PER [®] Insect Cell Protein Extraction Reagent, 250mL
78115	Inclusion Body Solubilization Reagent, 100mL
89835	DNase I, 5000 units
23236	Coomassie Plus (Bradford) Assay Kit
89890	Zeba [™] Spin Desalting Columns, 7K MWCO, 2mL, 25/pkg
89892	Zeba Spin Desalting Columns, 7K MWCO, 5mL, 25/pkg
89894	Zeba Spin Desalting Columns, 7K MWCO, 10mL, 25/pkg
66385	Slide-A-Lyzer [®] Dialysis Cassettes Kit, 10K MWCO, 0.1-0.5mL
66382	Slide-A-Lyzer Dialysis Cassettes Kit, 10K MWCO, 0.5-3mL
66807	Slide-A-Lyzer Dialysis Cassettes Kit, 10K MWCO, 3-12mL

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.