Thermo Scientific™ Eutech™ 6+ Handheld Meters

Conductivity/TDS/Salinity Operation Instructions

Models:

COND 6+:
Conductivity/Temp


TDS 6+:
TDS/Temp

Salt 6+:
Salinity/Temp

Getting Started/Connections

Remove rubber armour & install (4) AAA batteries (included). Connect the desired sensor to the corresponding ports.

Utilize our **Grip-Clip™** to attach a 12 mm or 16 mm sensor to a beaker. Extend the stand as shown for bench use.

Keypad Functions

	Press to power ON in the mode that was previously used, or OFF. Meter will automatically power off in 20 minutes.
	Toggle between measurement and calibration modes. In SETUP mode, BACK serves to return to the previous menu option or setting.
	Confirm calibration values in CAL mode. Confirm selections in SETUP mode. Freeze ("HOLD" indicator) or release the measured reading.
	Toggle between available measurement types. Press simultaneously with ON/OFF to begin SETUP (See Setup Options).
	Increase value or scroll up in SETUP or manual calibration.
	Decrease value or scroll down in SETUP or manual calibration.

Setup Options

	With meter off, press MODE and ON/OFF keys together to enter "SETUP" mode. Use ▲/▼ keys to select the desired option. Press ENTER to accept the desired setting. Press CAL/BACK to exit SETUP. The meter will retain all settings unless they are changed – even upon battery removal.
	Cell Constant. Choose (0.1 / 1.0 / 10.0). Default is 1.0. (COND 6+, TDS 6+ only)
	View calibration data for each calibration point stored and the associated calibration range.
	View electrode efficiency of each calibration point stored and the associated calibration range.
	Reset to default settings. "n0": Exits option without action. "CAL": Erases all calibration data except temperature while retaining other settings. "FEL": Returns all settings to factory default after ENTER is pressed then restarts meter.
	Select °Celsius or °Fahrenheit. Default is °C.
	Select Single-Point Calibration (SPER "YES" to apply a single calibration value across all ranges or "n0" to use multiple calibration points to calibrate individual range(s). Default is "YES".
	Normalization Temperature. Select "200" or "250". Default is 25.0 °C. (COND 6+, TDS 6+ only). Salt 6+ is fixed at 15 °C.
	Adjust Temperature Compensation from 0.0 (none) to 10.0 % per °C. Default is 2.0 % per °C. (COND 6+, TDS 6+ only)
	Choose "YES" for automatic calibration standards or "n0" for manual calibration adjustment. Default is "YES". (COND 6+ only)
	Adjust TDS factor from 0.30 to 1.00. Default is 0.5. (TDS 6+ only)
	Select Measurement Unit, TDS (ppt) or percentage (%). (SALT 6+ only)

Conductivity/TDS/Salinity Calibration

For best results, periodic calibration with known accurate standards is recommended. Calibrate with standard(s) near your intended measuring range. Provide stirring for best results. If multi-point calibration is desired, select "n0" in setup option SPER.

Range #	Conductivity Range*	Automatic Calibration Values	
		25 °C	20 °C
r 1	0.00 - 19.99 µS/cm	None	None
r 2	20.00 - 199.9 µS/cm	84 µS/cm	76 µS/cm
r 3	200.0 - 1999 µS/cm	1413 µS/cm	1278 µS/cm
r 4	2.00 - 19.99 mS/cm	12.88 mS/cm	11.67 mS/cm
r 5	20.0 - 200.0 mS/cm	111.8 mS/cm	102.1 mS/cm

Table 1 * Based on 1.0 cell constant

Calibration Using Automatic Conductivity Recognition (COND 6+ only)

- From Conductivity measurement mode, dip the sensor into your lowest standard – either 84 µS/cm, 1413 µS/cm, 12.88 mS/cm, or 111.8 mS/cm. Note: if different standards are used, manual adjustment must be used by selecting "n0" in Setup option REAL.
- Wait for the µS or mS unit to stop blinking to indicate that reading is stable. Press **CAL**. "r R" will appear briefly. The meter will automatically search for the appropriate standard value (84, 1413, 1288, 1118).
- The standard value will begin blinking. Press **ENTER** to accept. A successful calibration will display "done" and return to live measurement.
- To calibrate another standard, rinse your electrode then dip into your next highest standard. Repeat steps 2 and 3 as needed.

Using Manual Adjustment

- From Conductivity/TDS/Salinity measurement mode, dip the sensor into your lowest standard.
- Wait for the µS or mS unit to stop blinking to indicate that reading is stable. Press **CAL**. "r R" will appear briefly.
- The standard value will begin blinking. Use **▲/▼** to adjust the value to match the value of your standard. Press **ENTER** to accept. A successful calibration will display "done" and return to live measurement.

- To calibrate another standard, rinse your electrode then dip into your next standard. Repeat steps 2 and 3. Note that only one calibration point can be retained in each range (see Table 1). Calibrations that already exist in a range will be replaced if calibration is performed in the same range.

Temperature Calibration

Temperature calibration is recommended prior to first use, after ATC sensor replacement, and periodically as needed.

- Dip the temperature sensor into a solution with a known accurate temperature.
- Press **MODE** to measure temperature, then press **CAL**. The value without offset will display briefly, before displaying the actual value with any offset applied. The temperature reading will begin to blink.
- Use **▲/▼** to adjust to the known temperature. When the temperature reading is stable, press **ENTER** to accept.

Error Messages

- "LO": Low battery condition / replace batteries
- Unstable readings → air bubbles present in sample / clean probe & recalibrate / make sure sample covers the probe sensors / move or switch off interfering motor / replace probe
- "not stable": Enter was pressed while reading was not stable / wait for reading to stabilize
- Inaccurate readings or can't calibrate → install probe guard / replace calibration solutions
- Er1: Calibration error, reading is out of range → ensure that probe is connected / check temperature / check calibration solution
- "Or": Over range

Precautions


- Do not use this equipment in potentially explosive atmospheres.
- Refer to the electrode instructions for use, storage and cleaning.
- Ensure that no liquid enters the instrument.
- There are no user serviceable parts inside. Attempts to service internal parts may void the warranty.
- **WARNING:** No modification of this equipment is allowed.

Instrument Operating Conditions

Operating Ambient Temp.	5 to 45 °C
Operating Relative Humidity	5 to 85 %, non-condensing
Storage Temp.	-20 to +60 °C
Storage Relative Humidity	5 to 85 %, non-condensing
Weight (Meter Only)	205 g
Size (L x W x H) (Meter Only)	155 x 80 x 33 mm
Regulatory & Safety	CE, FCC Class A
Battery Requirement	4 x AAA, 1.5 V batteries
Vibration	Shipping/handling per ISTA #1A
Shock	Drop test in packaging per ISTA #1A

Ordering Information

ECCON603PLUSK	COND 6+ with electrode & carrying kit set
ECCON603PLUS	COND 6+ with electrode
ECTDS603PLUSK	TDS 6+ with electrode & carrying kit set
ECTDS603PLUS	TDS 6+ with electrode
ECSALT603PLUSK	Salt 6+ with electrode & carrying kit set
ECSALT603PLUS	Salt 6+ with electrode
CONSEN91B	Replacement Ultem electrode with ATC, k=1.0
CONSEN72B	Epoxy electrode with ATC, k=0.1
GRIPCLIP	Grip-Clip™ electrode holder
ECCON10BS	10 µS/cm calibration solution, 20 x 20 mL
ECCON84BT	84 µS/cm calibration solution, 480 mL
ECCON447BS	447 µS/cm calibration solution, 20 x 20 mL
ECCON1413BT	1413 µS/cm calibration solution, 480 mL
ECCON1413BS	1413 µS/cm calibration solution, 20 x 20 mL
ECCON1288BT	12.88 mS/cm calibration solution, 480 mL
ECCON1118BT	11.18 mS/cm calibration solution, 480 mL
ECCON3000BS	3000 ppm salinity calibration solution, 20 x 20 mL

Stops blinking when measurement is stable

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Eutech Instruments Pte Ltd