Premium Series Heated Immersion Circulators

PC200

PC201

PC300

ARCTIC Series Refrigerated/Heated Bath Circulators

A25

A45HC

SAHARA Series Heated Bath Circulators

S7 S13

S45 S49

Thermo Scientific Laboratory Temperature Control Products

Manual Part Number U01047
Rev. 03/03/2021

Multilingual Quick Start Guides

Multilingual Essential Safety Instructions

Installation

Operation

Preventive Maintenance

Troubleshooting

Visit our Web site at:

<http://www.thermofisher.com/tc>
Product Service Information, Applications
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

Premium Series Heated Immersion Circulators

PC200

PC201

PC300

ARCTIC Series Refrigerated/Heated Bath Circulators

A25

A45HC

SAHARA Series Heated Bath Circulators

S7 S13

S45 S49

Thermo Scientific Laboratory Temperature Control Products

Manual Part Number U01047
Rev. 03/03/2021

Multilingual Quick Start Guides

Multilingual Essential Safety Instructions

Installation

Operation

Preventive Maintenance

Troubleshooting

Visit our Web site at:

<http://www.thermofisher.com/tc>
Product Service Information, Applications
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

Thermo Fisher Scientific

25 Nimble Hill Road
Newington, NH 03801
Tel : (800) 258-0830 or
(603) 436-9444
Fax : (603) 436-8411
www.thermofisher.com/tc

Sales, Service, and Customer Support

25 Nimble Hill Road
Newington, NH 03801
Tel: (800) 258-0830
Sales: 8:00 am to 5:00 pm
Service and Support: 8:00 am to 6:00 pm Monday
through Friday (Eastern Time)
Fax: (603) 436-8411
service.tc.us@thermofisher.com

Dieselstrasse 4
D-76227 Karlsruhe, Germany
Tel : +49 (0) 721 4094 444
Fax : +49 (0) 721 4094 300
info.tc.de@thermofisher.com

Building 6, No. 27
Xin Jinqiao Rd., Shanghai 201206
Tel : +86(21) 68654588
Fax : +86(21) 64457830
info.china@thermofisher.com

Statement of Copyright

Copyright © 2021 Thermo Fisher Scientific. All rights reserved.

This manual is copyrighted by Thermo Fisher Scientific.

Users are forbidden to reproduce, republish, redistribute, or resell any materials from this manual in either machine-readable form or any other form.

Table of Contents

Quick Starts

Preface	i
Compliance	i
After-Sale Support	i
Feedback.....	ii
Unpacking.....	ii
Section 1 Safety	1-1
Safety Warnings.....	1-1
Personal Protective Equipment	1-3
Training.....	1-3
Section 2 General Information	2-1
Description and Intended Use.....	2-1
Immersion Circulator Specifications.....	2-1
Bath/Circulator Specifications.....	2-3
Wetted Materials.....	2-7
Section 3 Installation	3-1
Ambient Conditions	3-1
Immersion Circulator Only	3-1
Bath Circulator	3-2
Grounding Strap and Nut	3-2
Ventilation	3-3
Electrical Requirements	3-3
External Circulation	3-7
Remote Temperature Sensor.....	3-8
Multifunction Port.....	3-8
USB Port	3-8
Tubing Requirements	3-8
Tubing.....	3-9
Approved Fluids.....	3-10
Additional Fluid Precautions	3-12
Filling Requirements	3-15
Draining.....	3-15
Section 4 Operation	4-1
PREMIUM Immersion Circulator.....	4-1
Setup	4-2
Initial Start Up.....	4-2
Daily Start Up.....	4-3
Status Display	4-4
Changing the Setpoint.....	4-5
Changing a Value Display	4-6
Menu Displays	4-7
Menu Tree	4-9

	Settings Display Options	4-10
	Settings Setpoint.....	4-11
	Settings Temperature Alarm	4-13
	Settings Audible Alarms	4-13
	Settings Features.....	4-14
	MultiFunction Port Standard I/O DB-15HD	4-16
	Settings Fluid Type	4-17
	Settings Pump Adjustment.....	4-18
	Settings On/Off Timer	4-18
	Advanced Charting.....	4-19
	Advanced Ramping	4-20
	Running a Ramp Program.....	4-22
	Advanced Interfaces.....	4-24
	System Configuration.....	4-25
	System Run Time.....	4-25
	System Messages	4-25
	System Password.....	4-25
	High Temperature Cutout	4-26
	Stopping the Unit.....	4-27
	Power Down.....	4-27
	Shut Down.....	4-27
	Restarting	4-27
	Decommissioning/Disposal	4-28
	Storage	4-28
Section 5	Preventive Maintenance.....	5-1
	Cleaning	5-1
	Condenser Fins	5-1
	Grounding Strap and Nut	5-2
	Testing the Safety Features	5-2
Section 6	Troubleshooting.....	6-1
	Error Displays	6-1
	Checklist	6-3
Section 7	Password	7-1
	Operator Mode	7-1
	Reset	7-2
	Calibration.....	7-2
	PID Tuning.....	7-5
	User Mode.....	7-6
Appendix	Communications.....	A-1
Warranty		



⚠️ This quick start guide is intended for initial start up only. For all other procedures you must refer to the manual. Also, if any of these steps are not clear refer to the manual before proceeding.

Safety:

- The bath is designed for indoor use only. Never place the bath in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present.
- Leave refrigerated baths in an upright position at ~25°C for 24 hours before starting.
- Connect the bath to a properly grounded outlet.
- For Sahara baths, and stand-alone circulators, secure the grounding strap to the top of the metal tank.
- Never operate the equipment with a damaged line cord.
- The refrigerants are heavier than air and will replace the oxygen causing loss of consciousness. Contact with leaking refrigerant will cause skin burns. Refer to the bath's nameplate and the manufacturer's most current MSDS for handling precautions and disposal.
- Move the bath with care, sudden jolts or drops can damage its components. Always turn the equipment off and disconnect it from its supply voltage before moving it.
- Never operate damaged or leaking equipment.

What you need to get started:

- An adjustable wrench
- Appropriate hose or plumbing
- Appropriate size hose clamps
- Flat head screwdriver

- Ensure none of the tubing comes in contact with the power cord.
- Extreme operating temperatures will lead to extreme temperatures on the tube surface, this is even more critical with metal nozzles.
- Ensure the tubing you select meets your maximum temperature and pressure requirements.
- Do not subject tubing to mechanical strain and ensure any specified bend radius is not exceeded.
- Always turn off the circulator and disconnect the power cord from the power source before installing the optional platform or bridge.
- Limit all acrylic bath's maximum high temperature setting to the temperature indicated on the label on the front of the bath, 80 °C.
- Use only the approved fluids shown in Table 1. Before using any fluid where contact with the fluid is likely, refer to the manufacturer's MSDS for handling precautions and disposal. Also, refer to the MSDS for venting requirements.
- Avoid overfilling, oil-based fluids expand when heated.

The circuit protection is designed to protect the circulator. The circulator's line cord is designed to act as a disconnecting device, position the circulator so it is not difficult to access the cord.

Refer to the bath nameplate on the rear of the bath for specific electrical requirements. Voltage deviations of $\pm 10\%$ are permissible. The outlet must be rated as suitable for total power consumption.

- Ensure the cords do not come in contact with any of the plumbing connections or tubing.

1A. For refrigerated baths make all supplied communication and electrical connections prior to starting.

- Never connect controller power inlet, B, to a power outlet. Never connect power outlet, C, to anything but the circulator.

Install the supplied communications cable, A, between the circulator and the bath RJ45 connectors.

Install the power cord from the connector, B, to the connector on the bath, C. Connect the bath's power cord, D, to a grounded power outlet, E.

1B. For non-refrigerated baths the power supply, B, goes directly to a grounded power outlet, E.

2. Plumbing connections for external circulation are on the rear of the circulator.

→ is the return flow from the external application. ● is the outlet flow to the external application. The connections are 16 mm O.D.

If desired, remove the union nuts and plates to install the supplied 8 or 12 mm hose barsbs and clamps.

- To prevent damage to the circulator's plumbing, use a 19 mm backing wrench when removing/installing the external connections.

Ensure the drain port on the front of the bath is closed and that all plumbing connections are secure.

To avoid spilling, place any application containers into the bath before filling. Fill the bath work area from 2.0 cm (3/4") to 4.0 cm (1 1/2") below the top.

Table 1. Approved Fluids:

Do not start the circulator until fluid is added to the bath reservoir. Have extra fluid on hand.

Filtered/ single distilled water (pH 7-8)		
Deionized water (1-3 MΩ-cm, compensated)		
Distilled water with Nalco biocide and inhibitor		
Distilled water with chlorine (5 ppm)		
0 to 75% Laboratory Grade Glycol/Water		
SIL 100	SIL 300	SYNTH 60
SIL 180		SYNTH 260

- When using water above 80 °C monitor the fluid level, frequent top-offs will be required. It also creates steam.
- Water/glycol mixtures require top-offs with pure water, otherwise the percentage of glycol will increase resulting in high viscosity and poor performance.

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

For refrigerated baths, place its circuit protector located on the rear of the bath to the | position.

For all circulators, place the circuit protector located on the rear to the | position. The blue LED on the front panel illuminates.

3. Press the Start Display appears.

Press to the left of the start symbol, . The circulator will start and the start symbol will turn into a stop symbol, . The pump starts immediately but the compressor takes up to 30 seconds.

4. If desired, press to bring up the Status Display. Press to toggle between the Start/Status Displays.

After starting check all plumbing connections for leaks. Adjust the High Temperature Cutout (HTC) safety device, refer to the manual.

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Preface

WEEE

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2002/96/EC. It is marked with this symbol.

Thermo Fisher Scientific has contracted with one or more recycling/disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at:

www.thermoscientific.com/WEEERoHS

After-sale Support

Thermo Fisher Scientific is committed to customer service both during and after the sale. If you have questions concerning operation, or questions concerning spare parts or Service Contracts, call our Sales, Service and Customer Support phone number, see this manual's inside cover for contact information.


Sample Nameplate

There can be up to three nameplates located on the rear of the bath. Before calling, please obtain the serial number printed on the complete system nameplate located on the upper rear of the bath.

Nameplate

Refer to nameplate when calling for after-sale support

Nameplates (Typical Locations)

Feedback

We appreciate any feedback you can give us on this manual. Please e-mail us at tcmanuals@thermofisher.com. Be sure to include the manual part number and the revision date listed on the front cover.

Warranty

Thermo Scientific Laboratory Temperature Control Products have a warranty against defective parts and workmanship for 36 months from date of shipment. See back page of this manual for more details.

Unpacking

Retain all cartons and packing material until the circulator is operated and found to be in good condition. If the circulator shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

If an immersion circulator was purchased without a bath, carefully remove the label and small piece of Styrofoam securing the float switch before using.

Refrigerated baths should be left in an upright position for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Section 1 Safety

Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your circulator. If you have any questions concerning the operation of your circulator or the information in this manual, please contact us. See inside cover for contact information.

DANGER indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

WARNING indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

CAUTION indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It also alerts against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the circulator's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates the presence of hot surfaces.

This label indicates read the manual.

This label indicates a flammable refrigerant is used.

Observe all warning labels. ▲

Never remove warning labels. ▲

Leave refrigerated baths in an upright position at room temperature (~25°C) for 24 hours before starting. This ensures the lubrication oil drains back into the compressor. ▲

The circuit protector located on the rear of the circulator is not intended to act as a disconnecting means. ▲

Operate the circulator using only the supplied line cord. The circulator's power cord is used as the disconnecting device, it must be easily accessible at all times. ▲

Never operate the bath with the immersion circulator removed. ▲

Do not mount the immersion circulator backwards; the line cord could contact the reservoir fluid. Ensure the electrical cords do not contact any of the plumbing connections or tubing. ▲

Never place the circulator in a location or atmosphere where excessive heat, moisture, or corrosive materials are present. ▲

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Never operate the circulator without fluid in the reservoir. ▲

Oil-based fluids expand when heated. Fill the reservoir only to the minimum level. ▲

Use only the approved fluids listed in this manual. Using other fluids voids the warranty. ▲

Other than water, before using any approved fluid, or when performing maintenance where fluid contact is likely, refer to the manufacturer's SDS and EC Safety Data Sheet for handling precautions. ▲

When using ethylene glycol and water, check the fluid concentration and pH on a regular basis. Changes in concentration and pH can impact system performance. ▲

Do not clean the circulator with solvents, use a soft cloth and water. ▲

Drain the reservoir before it is transported, and/or stored. ▲

Always turn the circulator off and disconnect the supply voltage from its power source before moving or before performing any service or maintenance procedures. ▲

Transport the circulator with care. Sudden jolts or drops can damage its components. ▲

Refer service and repairs to a qualified technician. ▲

The user is responsible for decontamination if hazardous materials are spilled. Consult the manufacturer regarding decontamination and or cleaning agents compatibility. ▲

Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and voids the manufacturer's warranty. ▲

The circulators construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the users responsibility to assure a proper ground connection is provided. ▲

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Many refrigerants which may be undetectable by human senses are heavier than air and will replace the oxygen in an enclosed area causing loss of consciousness. Contact with leaking refrigerant will cause skin burns. Refer to the chiller's nameplate for the type of refrigerant used and then the refrigerant's SDS for additional information. ▲

Ensure, that no toxic gases can be generated by the fluid. Flammable gases can build up over the fluid during usage. ▲

Ensure the fluid is at a safe temperature (~40°C) before handling or draining. ▲

Never operate damaged or leaking equipment, or with any damaged cords. ▲

Never operate the circulator or add fluid to the reservoir with panels removed. ▲

Personal Protective Equipment

There are no special personal protective equipment requirements needed to perform normal operation, unless required by the manufacturer of the fluid you are using. Consult their Safety Data Sheet and/or EC Safety Data Sheet. We do recommend wearing eye protection and gloves.

Training

The user must review and understand all the sections in this manual before operating the bath.

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
Section 2 General Information

Description and Intended Use

The Thermo Scientific PREMIUM Series of Heated Immersion Circulators are used with refrigerated and heated baths. All circulators can pump to an external system. All circulators have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms, and offer adjustable high temperature protection.

The circulator is designed for use in a clean laboratory environment and in accordance with the Letter of Compliance located at the end of this manual.

PREMIUM Heated Immersion Circulator


	PC200	PC201	PC300
Temperature Range °C* °F	Ambient +13 to +200 Ambient +23 to +392	Ambient +13 to +200 Ambient +23 to +392	Ambient +13 to +300 Ambient +23 to +572
Temperature Stability °C	±0.01	±0.01	±0.01
Heater Capacity 230V/115V watts	2000/1200	3000/NA	3000/NA
Circulator Dimensions (H x W x D) mm inches	421 x 189 x 233 16.6 x 7.4 x 9.2	421 x 189 x 233 16.6 x 7.4 x 9.2	421 x 189 x 233 16.6 x 7.4 x 9.2
Reservoir Depth Requirement mm inches	200 7.9	200 7.9	200 7.9
Net Weight kg lb	4.7 10.4	4.7 10.4	4.7 10.4
Pumping Capacity Max flow rate lpm/gpm Max pressure (mbar/psi) Max suction (mbar/psi) Pump speed	24/6.3 560/8.1 380/5.5 Variable (40% - 100%)	24/6.3 560/8.1 380/5.5 Variable (40% - 100%)	24/6.3 560/8.1 380/5.5 Variable (40% - 100%)
Electrical Requirements (Voltage ±10%)	100 V/50 Hz 100 V/60 Hz or 115 V/60 Hz or 200-230 V/50..60 Hz	200-230 V/50..60 Hz	200-230 V/50..60 Hz
Connectivity Remote sensor port USB port Multi function port RS232 RS485 Ethernet/LAN Analog I/O	Yes Yes Yes Optional Optional Optional Optional	Yes Yes Yes Optional Optional Optional Optional	Yes Yes Yes Optional Optional Optional Optional

- Performance specifications established in accordance with DIN 12 876 (using water at 70°C).
- The maximum wall thickness for circulators that have a factory installed clamp is 26 mm.
- Communication options require an accessory box.
- Thermo Fisher Scientific reserves the right to change specifications without notice.


*The lowest usable temperature may be 13°C (23°F) warmer than room temperature, lower temperature ranges require supplemental cooling.

ARCTIC Refrigerated/Heated Bath Circulator Specifications

Stainless Steel Refrigerated/Heated Bath Circulators		
	A25	A45HC
PC200 Temperature Range °C °F	-25 to 200 -13 to 392	-40 to 200 -40 to 392
Bath Volume liters gallons	7 - 12 1.8 - 3.2	7 - 12 1.8 - 3.2
Cooling Capacity watts	500	900
Refrigerant	R134a	R290
Overall Dimensions* (H x W x D) mm inches	749 x 273 x 483 29.5 x 10.7 x 19.0	787 x 385 x 580 31.0 x 15.1 x 22.8
Work Area Dimensions (D x W x L) mm inches	200 x 173 x 184 8.0 x 6.8 x 7.2	200 x 173 x 184 8.0 x 6.8 x 7.2
Net Weight kg lb	36.1 79.5	55.2 115
Electrical Requirements** (Voltage ±10%)	100 V/50 Hz 100 V/60 Hz or 115 V/60 Hz or 230 V/50 Hz	115 V / 60 230 V / 50

Overall Dimensions




Work Area Dimensions

*Add ~26 mm (1 inch) to D for drain fitting.

**See Section 3 for additional information.

- Thermo Fisher Scientific reserves the right to change specifications without notice.

Cooling Capacity





Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. The bath was operating with the lid closed and no heat load lid or any external circulation.

Minimum temperatures are only achieved with the above conditions. If your conditions vary from above, minimum temperature will most likely not be achieved.

Specifications are for reference only and are subject to change.

Time to Temperature

Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. Pump specifications are nominal values of ±10%. Specifications are for reference only and are subject to change.

SAHARA Heated Bath Circulator Specifications

Stainless Steel Heated Bath Circulators		
	S7	S13
PC200/201 Temperature Range* °C °F	Ambient +13 to 200 Ambient +23 to 392	Ambient +13 to 200 Ambient +23 to 392
PC300 Temperature Range* °C °F	Ambient +13 to 300 Ambient +23 to 572	Ambient +13 to 300 Ambient +23 to 572
Bath Volume liters gallons	4 - 7 1.1 - 1.8	7 - 11 1.8 - 2.9
Overall Dimensions** (H x W x D) mm inches	494 x 235 x 428 19.5 x 9.2 x 16.7	494 x 321 x 428 19.5 x 12.6 x 16.7
Work Area Dimensions (D x W x L) mm inches	200 x 154 x 112 7.9 x 6.1 x 4.4	200 x 293 x 112 7.9 x 9.4 x 4.4
Net Weight kg lb	10.6 23.4	12.3 27.0
PC200 Electrical Requirements*** (Voltage ±10%)	100-115V/50-60Hz or 200-230V/50-60Hz	
PC201/300 Electrical Requirements*** (Voltage ±10%)	200-230V/50-60Hz	

Stainless Steel Heated Bath Circulators		
	S45	S49
PC200/PC201 Temperature Range* °C °F	Ambient +13 to 150 Ambient +23 to 302	Ambient +13 to 200 Ambient +23 to 392
Bath Volume liters gallons	30 - 41 7.9 - 10.8	29 - 51 7.7 - 13.5
Overall Dimensions** (H x W x D) mm inches	594 x 381 x 628 23.4 x 15.0 x 24.7	494 x 579 x 746 19.5 x 22.8 x 29.4
Work Area Dimensions (D x W x L) mm inches	300 x 298 x 312 11.8 x 11.7 x 12.3	200 x 498 x 430 7.9 x 19.6 x 16.9
Net Weight kg lb	20.3 44.7	24.3 53.4
PC200 Electrical Requirements*** (Voltage ±10%)	100-115V/50-60Hz or 200-230V/50-60Hz	
PC201 Electrical Requirements*** (Voltage ±10%)	200-230V/50-60Hz	

Specifications achieved with reservoir cover removed.


*The lowest usable temperature may be 13°C (23°F) warmer than room temperature, lower temperature ranges require supplemental cooling. **See page 2-2. Add ~26 mm (1 inch) to D for drain fitting.

***See Section 3 for additional information.


- Thermo Fisher Scientific reserves the right to change specifications without notice.

Time to Temperature


S7, S13 (115 V - 1.2 kW)


S45, S49 (115 V - 1.2 kW)


S7, S13 (230 V - 2 kW)


S45, S49 (230 V - 2 kW)

S7, S13 (230 V - 3 kW)

S45, S49 (230 V - 3 kW)

Specifications obtained at sea level using water (above +5°C to +90°C) or a fluid with a specific heat of 2.3 kJ/kg-K or 0.55 Btu/lb-F (less than 5°C) as the recirculating fluid at a +20°C ambient condition, at nominal operating voltage. Other fluids, process temperatures, ambient temperatures, altitude or operating voltage will affect performance. Pump specifications are nominal values of ±10%. Specifications are for reference only and are subject to change. Heat-up rates for the 100V baths will take approximately 25% longer than the 115V.

Wetted Materials

PREMIUM Immersion Circulator

Viton

EPDM

Ryton (PC200/201 only)

Ceramic

Ultem

Vectra

Stainless Steel

Stainless Steel Baths/Circulators

Stainless Steel 316

Stainless Steel 304

EPDM (drain fitting)

Ryton

Zotek-N (cover seal)

Section 3 Installation

Ambient Conditions

Ambient Temperature Range	5°C to 40°C (41°F to 104°F)
Maximum Relative Humidity	80% for temperatures up to 31°C (88°F)*
Operating Altitude	Sea Level to 2000 meters (6560 feet)
Overvoltage Category	II
Pollution Degree	2
Degree of Protection	IP 20

*decreasing linearly to 50% relative humidity at 40°C (104°F)

Ensure the label and small piece of Styrofoam securing the float switch are removed before using.

Immersion Circulator Only

The circulator is designed for continuous operation and for indoor use.

For stand-alone immersion circulators use the supplied nut and washer to secure the grounding strap to the top of the metal tank. See next page. ▲

Never place the immersion circulator in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲

Carefully install the immersion circulator to ensure it does not fall into the container or that its line cord does not make contact with the container contents. ▲

For immersion circulators equipped with a bridge:

- The bridge is designed to fit containers that are 400 mm to 800 mm wide.
- Slide the bridge support rods to the desired length and secure them in place with the supplied eight 3 mm Phillips Head screws.
- If possible, secure the bridge to your container using the two supplied thumbscrews.
- The immersion depth is 75 to 145 mm (~ 3 to $5 \frac{3}{4}$ ").

Your container must be sturdy enough to support the weight of the assembly, approximately 3.8 kilograms (8.5 pounds). The circulator is designed for continuous operation and for indoor use only.

Bath Circulator


The circulators is designed for continuous operation and for indoor use.

The equipment normally ships with the immersion circulator mounted facing the reservoir. You may change the position $\pm 90^\circ$ by removing thumb screws, no tools are required.

Do not mount it backwards; the line cord could contact the reservoir fluid. ▲

Rotate immersion circulator only $\pm 90^\circ$ from shown position

Never place the circulator in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲


Leave refrigerated baths in an upright position at room temperature ($\sim 25^\circ\text{C}$) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Grounding Strap and Nut

Sahara baths, and stand-alone circulators, require immersion circulator grounding.


For Sahara baths, and stand-alone immersion circulators, use the supplied nut and washer to secure the grounding strap to the top of the metal tank. ▲

A45HC units utilize a flammable refrigerant, risk of fire or explosion, do not use sharp objects to clean evaporator coil in the tank. ▲

Ventilation

All circulators except A45HC can operate with 0 clearance on two exhaust sides as long as the third exhaust side has unrestricted air flow. A45HC must have at least 25mm [1"] of clearance on both sides to comply with UL leak mitigation safety standards. Blocked ventilation increases the circulator's temperature, reduce its cooling capacity and, on refrigerated baths, eventually lead to premature compressor failure..

Ventilation Options

Electrical Requirements

The circulator's construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

EN 61000-3-11 Compliance: This equipment complies with EN/ IEC 61000-3-11 provided that main power system impedance (Z_{max}) is less than or equal to $0.074\Omega + j0.046\Omega$ [$0.074\Omega + 146 \mu H$] at the interfac point between the user's supply connection and the public system. Consult with the local supply authority if necessary for system Z_{max} determination.

For Sahara baths, and stand-alone immersion circulators, use the supplied nut and washer to secure the grounding strap to the top of the metal tank. See previous page. ▲

The circulator is intended for use on a dedicated outlet. All circulators are equipped with automatic thermally-triggered 20 Amp circuit protector.

The circuit protector is designed to protect the circulator.

Note If the circuit protector activates allow the circulator to cool before resetting. Restart the circulator. Contact us if it activates again. ▲


The circulator's power cord is the disconnecting device, it must be easily accessible at all times. ▲

Refer to the nameplate on the rear, upper-left-hand corner of the bath for specific electrical requirements. Voltage deviations of $\pm 10\%$ are permissible. The outlet must be rated as suitable for the total power consumption, see next page.


The following power options are available:

Bath	Volts ¹ /Hertz/Phase	Amps ²	Total Wattage	Plug Type
A25	115/60/1	11.7	1185	N5-15
	100/50-60/1	11.5	1135	N5-15
	230/50/1	10.4	2395	Country Specific
A45HC	115/60/1	14.4	1660	N5-20
	230/50/1	11.3	2600	Country Specific
PC 200 All Heated Baths/Circulators	115/60/1	11.3	1300	N5-15
	100/50-60/1	10.0	1300	N5-15
	230/50/1	9.3	2135	Country Specific
PC 201/300 All Heated Baths/Circulators	230/50/1	13.6	3140	Country Specific

1. volts ± 10%
2. Maximum amp dr

20 Amp Outlet
(16 Amp)


15 Amp Outlet
(12 Amp)

For refrigerated baths (PC200 only):

Ensure all communication and electrical connections are made prior to starting the circulator. ▲

- Install the supplied RJ45 shielded cable between the immersion circulator and the bath RJ45 connectors (similar to Ethernet). **This is required for proper operation.**


- Install the power cord from the connector on the rear of the circulator, A, to the connector on the rear of the refrigerated bath, B.
- Connect the bath's power cord, C, to a grounded power outlet.


For refrigerated baths, never connect power inlet A to a power outlet. Never connect power outlet B to anything but the circulator. ▲

Ensure the electrical cords do not come in contact with any of the plumbing connections or tubing. ▲


A45HC Only


For non-refrigerated baths:

- For PC201 install the power cord from the connector on the rear of the circulator, A, to a grounded power outlet.


- For PC300 install the power cord from the extended cord on the rear of the circulator, A, to a grounded power outlet.

Ensure the electrical cords do not come in contact with any of the plumbing connections or tubing. ▲

External Circulation

The plumbing connections for external circulation are located on the rear of the immersion circulator.  is the return flow from the external application.  is the outlet flow to the external application (supply side). The connections are 16 mm O.D. Remove the union nuts and plates to install the 8 mm or 12 mm hose barbs and clamps supplied with the circulator.

To prevent damage to the circulator's plumbing, use a 19 mm backing wrench on the nut when removing/installing the external connections. ▲

Remote Temperature Sensor

The remote temperature sensor on the rear of the immersion circulator requires a 4-pin connector that must mate to a LEMO # ECP.1S.304.CLL. The immersion circulator uses a 3 wire sensor, but a 4 wire sensor can be used (pins 3 and 4 are interconnected in the control head). The pin-out is:

Pin 1 and 2 = Pt100 + Pin 3 and 4 = Pt100 -

See Section 4 for instructions to enable the remote sensor.

Multifunction Port

The Multifunction port on the rear of the immersion circulator is used in conjunction with an *optional* accessory box to select and activate additional features. See Section 4 for additional information.

USB Port

If your computer does not automatically recognize the USB driver, installation instructions are provided in the Appendix.

Tubing Requirements

Ensure none of the tubing comes in contact with the power cord. ▲

Tubing is normally used to connect the pump to an external application.

Note The maximum allowable length of tube depends largely on the size, form and material of the external reservoir. The length of tube and its diameter, combined with the circulating capacity, have a large affect on the temperature stability. Whenever possible, use a wider tube diameter and place the application as close as possible to the circulator. ▲

Extreme operating temperatures will lead to extreme temperatures on the tube surface, this is even more critical with metal nozzles. ▲

- the required tube material depends on the heat transfer liquid used
- tubes must not be folded or bent
- after prolonged use, tubes may become brittle or they may get very soft, check them on a regular basis and replace if necessary
- secure all tube connections using clamps

When using the internal reservoir only, the plumbing connections can be closed with the supplied plate and union nuts.

Tubing

Tubing for Thermo Scientific temperature control systems is optional.

Ensure the tubing you select meets your maximum temperature and pressure requirements. ▲

Plastic and rubber tubing

If other plastic and rubber tubes are used, ensure that the tubes selected are fully suitable for the particular application, i.e., that they will not split, crack or become disengaged from their connections.

Connect the tubing using the supplied tube fittings for 8 or 12 mm i.d. They are attached to the plumbing connections with a supplied coupling nut.

We highly recommend using foam rubber insulation on the tubing and the fittings.

Metal tubing

Thermo Scientific metal tubing (stainless steel insulated) offers a particularly high degree of safety and is suitable for both low and high temperatures/liquids.

The metal tubing is attached directly to the plumbing connections, gaskets are not required.

Do not subject tubing to mechanical strain and ensure any specified bend radius is not exceeded. ▲

Tubing is available in lengths of 0.5, 1.0 and 1.5 meters. Couplings for connecting tubes are also available.

The smallest opening inside the metal tubes is 10 mm. The metal tubing is provided with coupling nuts (M16 x 1, DIN 12 879, part 2) at either end.

Approved Fluids

Only use the approved fluids listed below. ▲

All immersion circulators:

Water

Ethylene Glycol-Water

Propylene Glycol-Water

SIL 100

SIL 180

SIL 300

SYNTH 60

SYNTH 260

Thermo Fisher Scientific takes no responsibility for damages caused by the selection of an unapproved fluid.

Other than water, handle and dispose of approved fluids in accordance with the fluid manufacturer's specification and/or the fluid's SDS. ▲

To set the fluid temperature alarms, use the circulator's menu selection to identify the fluid, see Section 4. ▲

When using water above 80°C closely monitor the fluid level, frequent top-offs will be required. It will also create steam. ▲

Water/glycol mixtures require top-offs with pure water, otherwise the percentage of glycol will increase resulting in high viscosity and poor performance. ▲

Fire point

Flammable thermal liquids can ignite when a specified temperature is surpassed. The liquid is limited to a temperature level 25°C below the fire point as defined by the EN 61010-2-010.

Viscosity

For optimum temperature accuracy, it is important that heat transfer liquids have a low viscosity.

Working temperature range

This is the recommended long-term operating range. The maximum viscosity is approximately 5 mPas.

Operating temperature range


Long-term operation is recommended only under certain conditions. The viscosity may rise to a maximum of 30 mPas. The pump capacity will not match specifications.

Heating-up range

Long-term operation is not recommended, the pump motor's excess temperature protection may switch off the pump.

Range of Application

	Sil 100	Sil 180	Sil 300	Synth 60	Synth 260
Fire Point °C	>100	>225	>325	70	275
Flash Point °C	57	170	300	59	260
Viscosity at 20°C (mPas)	3	11	200	2	140
Density at 20°C (kg/dm ³)	0.89	0.93	1.08	0.76	1.03
Specific heat capacity (kJ/Kg*K)	1.67	1.51	1.56	2.10	2.00

Order Number 10 liter container	999-0202	999-0204	999-0206	999-0210	999-0214
Order Number 5 liter container	999-0201	999-0203	999-0205	999-0209	999-0213
Color	trans-parent, colorless	trans-parent, colorless	trans-parent, colorless	trans-parent, colorless	trans-parent, yellow
Reacts with	Silicone	Silicone	Silicone	Rubber Silicone	Copper Light metals Bronze

EC-Safety Data Sheets will be delivered together with each container of liquid.

Chlorine

Short term usage of tap water may not cause any adverse effects on the circulator or your application, but in the long term problems may arise. To help alleviate these problems Thermo Fisher Scientific recommends the use of chlorine.

The duration of time that chlorine remains in solution depends on factors such as water temperature, pH and availability of direct sunlight. We recommend maintaining chlorine levels at proper levels using chlorine test strips, generally 1 to 5 ppm is adequate.

For best results, the pH of the fluid should be maintained between 6.5 and 7.5. Additional chlorine should not be added without first determining the concentration ratio that already exists in the fluid supply. Corrosion and degradation of the circulation components can result from concentration ratios that are too high. Contact our customer support for additional information.

5°C to 95°C — Distilled Water or Deionized Water (up to 3 MΩ-cm)

Normal tap water leads to calcareous deposits necessitating frequent circulator decalcification, see table page 3-14.

Calcium tends to deposit itself on the heating element. The heating capacity is reduced and service life shortened.

-30°C to 80°C — Water with Glycol

Below 5°C water has to be mixed with a glycol. The amount of glycol added should cover a temperature range 5°C lower than the operating temperature of the particular application. This will prevent the water/glycol from gelling (freezing) near the evaporating coil.

Excess glycol deteriorates the temperature accuracy due to its high viscosity.

-40°C to 200°C — SIL180

SIL180 is suitable for covering nearly the entire range with just one liquid, especially when used with the refrigerated baths.

Unfortunately SIL180 has a wetting tendency necessitating the occasional cleaning of the reservoir cover.

other temperatures:

Thermo Fisher Scientific offers a range of heat transfer fluids for these temperature control applications.

SYNTH 60 and SYNTH 260:

Synthetic thermal liquid with a medium life span (several months) and little smell annoyance.

SIL 100, SIL 180, SIL 300:

Silicone oil with a very long life span (over 1 year) and negligible smell.

Thermo Fisher heat transfer fluids are supplied with an EC Safety Data Sheet.

Ensure, when selecting the heat transfer fluid, that no toxic gases can be generated. Inflammable gases can build up over the fluid during usage. ▲

Ensure the over temperature cut-off point is set lower than the fire point for the heat transfer fluid selected. ▲

The highest working temperature, as defined by the EN 61010-2-010, must be limited to 25°C below the fire point of the fluid. ▲

Additional Fluid Precautions

When working with approved fluids other than water:

- Do not use any approved fluid until you have read and understood the label and the Material Safety Data Sheet (SDS).
- Ensure any fluid residue or any other material is thoroughly removed before filling the reservoir with a different fluid.
- Always wear protective clothing, especially a face shield and gloves.
- Avoid spattering on any of the circulator's components, always *slowly* add fluid. When adding, point the opening of a container away from yourself.
- Use fume hoods.
- Do not allow any ignition sources in the vicinity.

Water Quality and Standards

Process Fluid	Permissible (PPM)	Desirable (PPM)
Microbiologicals		
(algae, bacteria, fungi)	0	0
Inorganic Chemicals		
Calcium	<25	<0.6
Chloride	<25	<10
Copper	<1.3	<1.0
	0.020 ppm if fluid in contact with aluminum	
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates/Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<25	<1
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
Other Parameters		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*

* MΩ-cm (compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting. Eventually, the pitting becomes so extensive that refrigerant will leak into the reservoir.

As an example, raw water in the United States averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the reservoir with distilled or deionized water. Do not use untreated tap water as the total ionized solids level may be too high. This reduces the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

Filling Requirements

Ensure the reservoir drain port on the front of the bath is *closed* and that all plumbing connections are secure. Also ensure any residue is thoroughly removed before filling the reservoir.

Before using any fluid refer to the manufacturer’s SDS and EC Safety Data Sheets for handling precautions. ▲

To avoid spilling, place your containers into the reservoir before filling.

With a low level WARNING the circulator continues to run, with a FAULT the circulator shuts down the refrigeration, pump and heater, see Section 6. The low level warning is at approximately 43 mm (1 3/4”) below the top, the low level fault is at approximately 57 mm (2 1/4”).

Oil-based fluids expand when heated. Fill the reservoir only to the minimum level. Monitor the fluid level whenever heating the fluid. ▲

When pumping to an external system, keep extra fluid on hand to maintain the proper level in the circulating lines and the external system.

Draining

Before draining any fluid refer to the manufacturer’s SDS and EC Safety Data Sheet for handling precautions. ▲

Ensure the fluid is at a safe handling temperature, ~40°C. Wear protective clothing and gloves. ▲

- place a suitable container underneath the drain. If desired, attach an 8 mm i.d. tube on the drain.
- *slowly* turn the drain plug until flow is observed.

Turning the drain cap more than 1 1/2 turns will result in the drain cap and fitting coming off the bath. ▲

In this case screw the drain fitting back into the bath. Attaching the cap onto the fitting aids in installation. If required, contact us for additional information.


Installed Drain Fitting with Cap Removed

Section 4 Operation

PREMIUM Heated Immersion Circulator

The Thermo Scientific PREMIUM Heated Immersion Circulators have a digital display and easy-to-use touch pad, five programmable setpoint temperatures, acoustic and optical alarms and adjustable high temperature protection.

This label indicates read the instruction manual before starting the circulator.

Use this button to place the circulator in and out of stand by mode. The blue LED is on when stand by mode is enabled.

Use these arrows to move through the displays and adjust values.

Use these arrows to adjust values.


Pressing this button to save changes made on the circulator's displays.

Use this button to cancel any changes and to return the circulator to its previous display. Canceling a change can only be made before the change is saved. In some cases, it is also used to save changes.

Note Holding this button depressed for five seconds brings up the language selection display. ▲

Used for adjusting and resetting the High Temperature Cutout. Details are explained in this Section.

Setup


Leave refrigerated baths in an upright position at room temperature (~25°C) for 24 hours before starting. This will ensure the lubrication oil has drained back into the compressor. ▲

Before starting the circulator, double check all communication, electrical and plumbing connections. ▲

Initial Start Up

Ensure the HTC is in the full clockwise position.


Do not run the circulator until fluid is added to the reservoir. Have extra fluid on hand. If it does not start refer to Section 6 Troubleshooting.

- For refrigerated baths, place the circuit protector located on the rear of the bath to the **I** position.
- For all circulators, place the circuit protector located on the rear of the immersion circulator to the **I** position. The blue LED illuminates.
- Press . The circulator momentarily displays

Thermo Scientific

Start Display

- And then the Start Display appears. The blue LED goes out.
- Press to the left of the start symbol, . The circulator starts and the start symbol turns into a stop symbol, .

Note After start up, check all the plumbing connections for leaks. ▲

The **Menu** and **SP1** portions on the top of the display are used to view and/or change the circulator's settings. The reservoir fluid is also displayed. These features are explained in detail later in this Section.

indicates the circulator is using its internal sensor for temperature control. is displayed when the external sensor is selected for temperature control.

Daily Start Up

Before starting the circulator, double check all communication, electrical and plumbing connections. ▲

The blue LED on the immersion circulator illuminates to indicate it is in stand by mode.

If so, press . The blue LED goes out. If Auto Restart is enabled the circulator starts. If Auto Restart is not enabled, or if the blue LED is not illuminated, use the **Initial Start Up** procedure described on the previous page.

The Auto Restart feature is described later in this Section.

Note When using oil as a reservoir fluid, we recommend running the circulator at 95°C for 15 minutes to remove any moisture in the fluid. ▲

Status Display

If desired, press  to toggle between the Start and Status Displays.

Status Display

This icon is illuminated when the pump is running.

This icon is illuminated when the heater is on. The icon flashes when the reservoir fluid temperature is close to the desired setpoint.


This icon is illuminated when the refrigeration is on. The icon flashes when the circulator is in the energy saving mode, the mode is explained later in this section.

One of these icons is on to indicate which sensor is selected for temperature control, internal or external.

This icon is illuminated when the serial communications or analog option is enabled.


The icon is illuminated when both the circulator and timer are running. The icon flashes when the timer is enabled but not running.

Note If no operator inputs are being made, the circulator automatically switches to the Status Display after 60 seconds. If desired, change the time or disable this feature using the **Display Options** Menu. ▲


Changing the Setpoint

The Setpoint is the desired fluid temperature. The circulator can store up to five setpoints, **SP1** through **SP5**. The procedure for changing the stored setpoint values is discussed later in this Section.


- On the Start Display press  to the right of **SP1**.

The following display appears:

- Press  to the right of the desired setpoint (in this example **SP3**) and then press , or  to the left of the home symbol , to return to the Start Display.


Change a Value Display

This display automatically appears when you request a change in value.


Note The circulator does not need to be running to change a value. ▲

- A cursor appears on the left side of the entry window as shown below.

- Press  or  to enter the corresponding digit. Use the arrow keys on the circular keypad to move the cursor back to the left , enter a negative value  or insert a decimal point .

Change a Value Display

- Press  when the desired value is displayed.


Note The circulator beeps if you try to enter a value outside the minimum or maximum value or by trying to enter more than one negative dash or decimal point. ▲

Menu Displays


The circulator uses menus to view/change its settings.

Note Unless noted, the circulator does not need to be running to view/change settings. ▲


- Press  to the left of **Menu** to bring up the Main Menu Display.


- Press  to the left of the desired line to bring up submenus.

See page 4-10.


See page 4-18.


See page 4-23.

Menu Tree

Settings - Display Options

Display Options is used to view/adjust the circulator's temperature scale, select the temperature resolution and set the display delay.


Language is used to select which language is displayed, English, German, French, Spanish, Italian, Chinese or Japanese.

Temp. Units is used to select the displayed temperature scale, °C, °F or K.

Temp. Resolution is used to select the displayed temperature resolution, **0.01** or **0.1** degree.

Display Delay is used to enable/disable and to select the Start Display display delay. The delay range is **5** to **900** seconds. If no arrows are pressed, the Start Display changes to the Status Display after the delay expires, see pages 4-2 and 4-4.

Press  next to the **Delay** to bring up the Change a Value Display.

Settings - Setpoints

Setpoints is used to view/adjust the circulator's five Setpoints (**SP**) and Real Temperature Adjustments (**RTA**). The display shows the temperature measured by the internal or external sensor, if this temperature does not accurately reflect the actual temperature in the reservoir an RTA is required. The RTA can be set $\pm 10^{\circ}\text{C}$ ($\pm 18^{\circ}\text{F}$).

As an example, if the temperature is stabilized and displaying 20°C but a calibrated reference thermometer reads 20.5°C , set the RTA to -0.5°C . After you enter a RTA value allow display to stabilize before verifying the temperature in the bath. **Note** If display accuracy is required, we recommend repeating this procedure at various setpoint temperatures and on a regular basis. ▲

Note For refrigerated baths, the setpoint limits are determined by the combination of the circulator, bath and fluid. The highest setpoint is the lesser of the highest limit of the three while the lowest setpoint is the greater of the lowest limit of the three. For heated baths, only the circulator and fluid are used. ▲

Circulator Limits $^{\circ}\text{C}$:

	Min	Max
PC 200	-90	200
PC 201	-90	200
PC 300	-50	300


Bath Limits $^{\circ}\text{C}$:

	Min	Max
A25	-25	200
A45HC	-40	200


Fluid Limits $^{\circ}\text{C}$:

	Min	Max
SYNTH 60	-50	45
SYNTH 260	45	250
SIL 100	-75	75
SIL 180	-40	200
SIL 200	20	200
SIL 300	80	300
Water	5	95
EG/Water	-30	100
Other	-90	300

When selecting a different fluid, the limits shrink if they are outside of the new system limits. However, they do not expand if they are well within the system limits.


Press  next to the desired setpoint to bring up:

Press  next to the desired setting to bring up the Change a Value Display to make changes.

Settings - Temperature Alarms


Temperature Alarms is used to view/adjust the temperature warning and fault limits. If a limit is exceeded the alarm sounds, if enabled. In the case of a fault the circulator also shuts down. **Note** The alarm range limits depend on the fluid, circulator and bath used. ▲

Press  next to the desired fault/warning to make changes .

Settings - Audible Alarms


Audible Alarms is used to enable/disable the audible alarms.

Press  to enable/disable the desired alarm(s).

Settings - Features

Features is used to enable/disable immersion circulator features.

Press  next to the desired feature to enable/disable.

The **Remote Sensor** connector on the rear of the immersion circulator requires a 4-pin connector that must mate to a LEMO # ECP.1S.304. CLL. The immersion circulator uses a 3 wire sensor, but a 4 wire sensor can be used (pins 3 and 4 are interconnected in the control head). The pin-out is:

Pin 1 and 2 = Pt100 + Pin 3 and 4 = Pt100 -


When **Low Level Warning** is enabled, the circulator sounds an alarm if a low level condition occurs.

When **Auto Restart** is enabled, the circulator automatically restarts, if it was running, after a power interruption condition. If a ramp was running when power failed, the ramp program resumes where it left off. **Note** Consider any possible risks before enabling this mode of operation. ▲

The **Energy Saving** mode is primarily designed for applications running under a stable load. Enabling the mode saves energy by reducing the circulator's heater power and cooling requirements. This can result in substantial energy savings over the life of the circulator. The default setting is disable (enable for A45HC refrigerated circulators).

The refrigeration icon ❄️ blinks when the **Energy Saving** mode is enabled.

Multi Function Port


Booster Pump (pins 4 and 14) is on whenever the option is enabled.

Boost Heater (pins 5 and 10) is on when enabled and the setpoint is $>2^{\circ}\text{C}$ above the bath temperature. The boost heater is normally off when the bath temperature is within 2°C of setpoint.

Cooling Coil (pins 4 and 14) is on when the bath temperature is $>2^{\circ}\text{C}$ above setpoint or when the heater power drops to 0%, it's off when the heater power exceeds 80%. This feature is used in conjunction with an *optional* accessory box to select and activate additional features. The accessory box connects to the multifunction port.

Auto Refill (pins 5 and 10) is on when the fluid level drops below 47 mm ($\sim 1\ 7/8''$) from the top and is off when the fluid level is at 27 mm ($\sim 1''$) from the top. This feature is used in conjunction with an *optional* accessory box to select and activate additional features. The accessory box connects to the multifunction port.

Multifunction Port Pin Out


See next page for additional information.

MultiFunction Port Standard I/O DB-15HD


Alarm Output Pins 13(+) and 3 (-) Open Collector Output

This output turns on (conducts) when there is an alarm. Any alarm that turns off the circulator will also set this output. An alarm message is displayed. (Maximum current on the output cannot exceed 25mA and 24 VDC.) A resistor must be supplied to limit the current, see samples below.


External Alarm Input Pins 12(+) and 2 (-)

To enable this input, apply 10 - 20 mA to pins 12 and 2. This input turns off the circulator. The alarm continues to sound until the input is removed. See samples below.


External On/Off Input Pins 11(+) and 1 (-)

To enable this input, apply customer supplied 10 - 20 mA to bath pins 11 and 1. This input turns the circulator on and off. A message is displayed when the input turns off. The circulator logic is reversed for this input: conducting = circulator off, not conducting = circulator on. See samples below.


Samples:

ON STATE Current = 10 - 20 mA

Settings - Fluid Type

Fluid Type is used to identify the type of fluid used. The circulator uses the fluid type to automatically set certain operating parameters.

Press  or  next to the desired fluid or fluid type.


Note For high viscosity fluids with a wide "Heating-up range" (e.g., SIL300 - see Section 3), the circulator automatically slows the heating from ambient to the "Operating temperature range" to help prevent fluid degradation that is caused by rapid heating while still at high viscosities. ▲

Fluid Limits °C:

	Min	Max
SYNTH 60	-50	45
SYNTH 260	45	250
SIL 100	-75	75
SIL 180	-40	200
SIL 200	20	200
SIL 300	80	300
Water	5	95
EG/Water	-30	100
Other	-90	300

Settings - Pump Adjustment


Pump Adjustment is used to set the pump speed.

Press  next to **Pump Speed** to bring up the Change a Value Display to make changes. **Note** The minimum pump speed is 40%. ▲


Settings - On/Off Timer

On/Off Timer is used to set/verify the timer date and time as well as its on/off settings. It is also used to enable/disable the timer.

Advanced - Charting

Charting is used to view a graphical display of your program.

Enable the green box to display the external temperature sensor chart. The current external temperature sensor is also displayed in the upper left hand corner of the display.

Enable the yellow box to display the internal temperature sensor chart. The current internal temperature sensor is also displayed in the upper middle of the display.

Enable the blue box to display the setpoint chart. The current setpoint is also displayed in the upper right hand corner of the display.

The white horizontal axis on the bottom of the display represents time. The right end of the display always shows the current time. Use the white zoom tools to change the displayed time increments. There are eight incremental ranges available, from one minute to two hours. Units are minutes/division or hours/division.

The red vertical axis on the left hand side of the display represents temperature. The temperature scale depends on the units selected in **Menu - Display Options - Temp. Units**. Use the red zoom tools to change the displayed temperature increments. There are 12 incremental ranges available, from 0.1 of a degree to 500 degrees.

The red dash is used to horizontally center the setpoint chart on the display. If desired, use the triangles to move the desired setpoint chart up or down the display. Use the red dash again to center it.

Advanced - Ramping

Ramping is used to build/edit a ramp program.

▶	Ramp 1	Ramping	Step 1	◀
▶	Steps 2		Start Temp 40.00°C	◀
▶	Cycles 1		End Temp 20.00°C	◀
▶	Variance 0.05°C		Duration 1min	◀
▶	End State		<input type="checkbox"/> Assured Soak	◀

Define your program as a series of setpoints with a known period of time interval between each. Each interval is one step of the program. Pay careful attention to the first part of your program. What conditions must exist at the beginning of your process? For example, at the starting setpoint you may need to program an initial period of constant temperature to allow for thermal stabilization.

Note Consider the circulator's limitations when designing programs. Temperature or time parameters which exceed the performance capabilities of the circulator will result in unsatisfactory operation. If reaching the ramp setpoint temperatures is important, you will have to operate the bath between the desired setpoints and note the duration before programming the ramp. ▲

It is possible to create a program calling for very rapid changes in temperature. Although the circulator may not be capable of producing such changes, it may be practical to program such steps as a way to cause the fastest possible temperature change.

There can be up to 10 **Ramps**, each with up to 30 **Steps**. Each **Ramp** can be **Cycled** up to 100 times.

Cycles sets the number of times the entire ramp program is repeated after the last step is completed. For example, selecting 3 **Cycles** runs the entire ramp program a total of 3 times.

The **Variance** is used to set a temperature range for starting a program, the program will not start until the fluid temperature is within this range. For example, if the desired **Start Temp** is 25°C and the **Variance** is set to +5°C, the program automatically starts only when the fluid temperature is between 20°C to 30°C. **Variance** can only be viewed/changed when **Step 1** is displayed.

Note By default, the **End Temp** for any step becomes the **Start Temp** for the following step. ▲


The maximum **Duration** for a step is 1000 minutes.

End State configures the circulator to either **Shut Down** or **Maintain** (continue running) when the program is over.

The ramp program has an optional **Assured Soak** feature that can be enabled for each step independently. When enabled this feature pauses the ramp timer while the temperature reaches setpoint, $\pm 0.01^\circ\text{C}$. This assures the temperature reaches setpoint before the ramp program continues to the next step.


Running a Ramp Program

You can enable an alarm to sound when each step and/or the program is complete, see **Settings - Audible Alarms** in this section.


- Press  to the right of **SP1**.


The following display appears:

- Select the desired ramp and then press , or  to the left of the home symbol , to return to the Start Display.

- Press  to the left of  to start the program. The display changes to the Program Display

The program does not start until the process fluid temperature is at the **Start Temp** \pm the **Variance**.

- Press  to the left of  to stop the program or  to the right of  to pause the program at the current temperature.

When the program is complete the circulator maintains the last setpoint.

When  is pressed after the ramp program is complete the refrigeration, pump and circulator shut off.

Note Other than the current ramp step, you can edit any other step of the ramp during the run or in the pause mode. ▲

The numbers at the bottom indicates the cycle number, step number and total time remaining for the program (hours, minutes and seconds).

Note If assured soak is enabled the time stops counting down at the end of the step until the desired temperature $\pm 0.01^\circ$ is reached. ▲

When the program is complete the alarm, if enabled, sounds.

Note If you do not stop/end the ramp program before going back to the setpoint (**SP**) mode, changing the setpoint has no affect. For example on a program which ends with the bath at 25°C the bath continues to hold this setpoint. Changing from ramp mode to setpoint mode and then changing the setpoint to 20°C has no affect, the bath will continue to hold 25°C. ▲




Advanced - Interfaces

Interfaces is used to view/change communication and communication settings.

Note Keypad operation is still available with serial communications enabled. ▲


- Press  on the right to select Off, RS-232 or RS-485, or Analog IO. See Appendix for additional information.

System - Configuration

Configuration is used to view the circulator's firmware version(s).

System - Run Time

Run Time is used to view the operating hours for the circulator and its components.


System - Messages

Messages is used to view and, if desired, delete any stored messages. The circulator can store and display up to the last 10 Warning and last 10 Fault messages. Messages can be then deleted individually or as a group.

System - Password

Password is used if the circulator needs a calibration or if the circulator PID values require changing. Changing the Password also allows you to restore the circulator to factory preset values. See Section 7.

High Temperature Cutout

HTC (Temperature range varies with type of immersion circulator.)

To protect your application, the adjustable High Temperature Cutout (HTC) ensures the heater does not exceed temperatures that can cause serious damage. A temperature sensor is located in the reservoir. A HTC fault occurs when the temperature of the sensor exceeds the set temperature limit.

In the event of a fault the circulator shuts down and displays a fault message, see Section 6. The cause of the fault must be identified and corrected before the circulator can be restarted. A low reservoir fluid level is a primary reason for the HTC to trip.

The HTC is factory preset fully clockwise to the highest possible setting. To set the cutout start the circulator and adjust the setpoint a few degrees higher than the highest desired fluid temperature. Allow the circulator to stabilize at the temperature setpoint. Using a small diameter (less than 1/8") plastic/metal rod, turn the arrow counterclockwise until the circulator shuts down and the fault message appears. Press  to clear the message.

Note If using plastic rod ensure the fluid temperature is low enough to prevent possibly melting the plastic. ▲

Before you can restart the circulator it has to cool down a few degrees. To restart the circulator press the black reset ring surrounding the red dial - and then press  again. If Auto Restart is enabled the circulator restarts, if disabled use the Start Up procedure.

Note: We recommend periodically rechecking operation or if the circulator is moved. ▲

Stopping the Circulator

- Press  to the left of the stop symbol, . The circulator stops and the stop symbol turns into a start symbol, .

Power Down

- Press . The circulator's display goes blank and the blue LED illuminates indicating it is in the stand by mode.

Shut Down

- On refrigerated baths, place the circuit protector on the rear of the bath to the **O** position.
- Place the circuit protector on the rear of the immersion circulator to the **O** position. The blue LED extinguishes.

Using any other means to shut the circulator down can reduce the life of the compressor. ▲

Always turn the circulator off and disconnect it from its supply voltage before moving. ▲

The circuit protector(s) located on the rear of the component(s) is not intended to act as a disconnecting means. ▲

Restarting

Note When quickly restarting refrigerated baths, the compressor may take up to 10 minutes before it starts to operate. ▲

Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current SDS for handling precautions. ▲

Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲

A45HC units utilize flammable refrigerant (R290), risk of fire or explosion. Dispose of properly in accordance with federal or local regulations.

Consider decommissioning the circulator when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant (R134A) and oil (POE or Ester) used in all other refrigerated baths must be recovered from equipment before disposal.

Note Keep in mind any impact your application may have had on the circulator. ▲

Direct questions about chiller decommissioning or disposal to our Sales, Service and Customer Support.

Handling and disposal should be done in accordance with the manufacturer's specification and/or the SDS for the material used. ▲

Storage

If the circulator is to be transported and/or stored it needs to be drained and then flushed with a 50/50 laboratory grade glycol/water mixture. ▲

The circulator can be stored for up to 90 days inside the temperature range of -25°C to +60°C (-13°F to +140°F).

If necessary, when removing the circulator from storage allow it time to warm up and dry out in order to prevent any condensation issues.

Section 5 Preventive Maintenance

Disconnect the power cord prior to performing any maintenance. ▲

Handle the circulator with care. Sudden jolts or drops can damage its components. ▲

Cleaning

After time, the circulator's stainless steel surfaces may show spots and become tarnished. Clean the reservoir and built-in components at least every time the fluid is changed. Use only water and a soft cloth.

Do not use scouring powder. ▲

A45HC units utilize a flammable refrigerant (R290), risk of fire or explosion, do not use sharp objects to clean evaporator coil in the tank. In order to ensure a long service life keep the reservoir clean. Quickly remove substances containing acidic or alkaline substances and metal shavings that could harm the surfaces causing corrosion. If corrosion (e.g., small rust marks) occurs in spite of this, cleaning with stainless steel caustic agents is suitable. Apply these substances according to the manufacturer's recommendations.

Condenser Fins

When cleaning the reservoir do not use any substances which contain solvents. ▲

For refrigerated baths, in order to maintain the cooling capacity clean the fins two to four times per year, depending on the operating environment.

Switch off the circulator and unplug the power cord.

For ARCTIC A45HC

- 1 Remove the condenser panel.
- 2 Clean fins with brush or similar tool.
- 3 Replace the panel.

For all other refrigerated baths:
Clean the fins with compressed air.

For extreme cases have a qualified technician remove the cooling compressor casing.

Tubing


Inspect and tighten the tubing and clamps daily.

Grounding Strap and Nut

Sahara baths, and stand-alone circulators, require immersion circulator grounding.

If the circulator is removed for any reason you must ensure the nut and washer secures to the grounding strap located on the rear of the reservoir top when the circulator is reinstalled. ▲

Testing the Safety Features

Check the safety features for high temperature protection and low liquid level protection at regular intervals. The frequency depends on the circulator's designated application and the heat transfer fluid used.

High temperature protection

Using a small diameter (less than 1/8") plastic/metal rod, turn the arrow to the desired temperature.

Note If using plastic rod ensure the fluid temperature is low enough to prevent possibly melting the plastic. ▲

Set a cut-off temperature that is lower than the desired setpoint temperature.

Switch on the circulator and ensure the circulator shuts down at the set cut-off temperature.

If the circulator does not shut down have it checked by a qualified technician.

Reset the safety to the desired temperature.

Low liquid level protection

With the circulator on, use a screwdriver and slowly push down on each level sensor, one at a time, until an error message appears. See Section 6 for details on error messages.

If not, have the circulator checked by a qualified technician.

Error Displays

Section 6 Troubleshooting

The circulator can display error messages and, if enabled, sound an alarm. Error messages are cleared by pressing the enter key (↵). Restart the circulator once the cause of the error message is identified and corrected. If the cause was not corrected the error code will reappear, contact our Sales, Service and Customer Support. If **Auto start** is enabled the circulator will restart, if disabled use the Start Up procedure.

FAULT:
High Temperature
PRESS ENTER
to clear message

Error Display (Typical)

Fault Displays

The heating element, pump and, if applicable, refrigeration shut down with a fault. A fault also sounds the alarm, if enabled.

Message	Cause	Actions
Bath not found	<ul style="list-style-type: none"> a loss of communication between refrigerated bath and circulator 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check cable connections on the rear of the bath and circulator
High Fixed Temp.	<ul style="list-style-type: none"> circulator's nonadjustable high temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check fluid selection check environmental conditions
High Temperature	<ul style="list-style-type: none"> adjustable high temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check limit setting check fluid selection ensure circulator has adequate ventilation
High Temperature Refrigeration	<ul style="list-style-type: none"> high refrigeration temperature 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check voltage supply the refrigeration may need servicing, contact us
HPC High Press. Cutout	<ul style="list-style-type: none"> the high refrigeration pressure cutout activated 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check for obstructions to air flow the refrigeration may need servicing, contact us

Message	Cause	Actions
HTC-LLC	<ul style="list-style-type: none"> if the fluid level is greater than approximately 57 mm (2 1/4") below the reservoir top it is a low level fault. if not, high temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) if required, check for leaks and fill fluid to proper level on smaller baths, 2 - 4 liters, if the fluid is near the minimum level reducing the pump speed may eliminate nuisance messages. if fluid is not the issue, allow circulator to cool down to at least 10°C below the HTC setting turn the red knob on the HTC fully clockwise  °C press the HTC's black reset ring press the enter key (↵) again reset HTC to desired setting, see Section 4 if message reappears recycle power to circulator and repeat the procedure increasing the pump speed to maximum when using fluids with viscosity greater than 50Cst may eliminate nuisance messages
High Level	<ul style="list-style-type: none"> reservoir fluid level too high for safe operation 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check fluid level, drain excess fluid if required verify optional auto refill operation
Low Fixed Temp.	<ul style="list-style-type: none"> circulator's nonadjustable low temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check fluid selection
Low Temperature	<ul style="list-style-type: none"> adjustable high temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check limit setting check fluid selection
Motor Fault	<ul style="list-style-type: none"> pump is not turning due to possible motor overload 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) it can take over 10 minutes for the motor temperature to get low enough before the circulator can be restarted if the motor does not start, cycle the power if multiple attempts do no clear the issue the motor probably needs to be replaced, contact us

Message	Cause	Actions
MOL Motor Overload	<ul style="list-style-type: none"> high motor overload temperature 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) allow circulator to cool down
Open RTD1 Internal	<ul style="list-style-type: none"> open internal temperature sensor 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) contact us
Open RTD2 External	<ul style="list-style-type: none"> open external temperature sensor 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) contact us
Shorted RTD1 Internal	<ul style="list-style-type: none"> shorted internal temperature sensor 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) contact us
Shorted RTD2 External	<ul style="list-style-type: none"> shorted external temperature sensor 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) contact us

Warning Displays

The circulator will continue to run with a warning. A warning also sounds the alarm, if enabled.

Message	Cause	Actions
Bad Calibration	<ul style="list-style-type: none"> bad temperature probe calibration 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) redo calibration
High Temperature	<ul style="list-style-type: none"> adjustable high temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check limit setting check fluid selection
Low Level	<ul style="list-style-type: none"> reservoir fluid level too low for safe operation 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check fluid level on smaller baths, 2 - 4 liters, if the fluid is near the minimum level reducing the pump speed may eliminate nuisance messages.
Low Temperature	<ul style="list-style-type: none"> adjustable low temperature protection limit exceeded 	<ul style="list-style-type: none"> clear message by pressing the enter key (↵) check limit setting check fluid selection

Messages

Message	Actions
Bad Calibration	<ul style="list-style-type: none"> redo calibration
Bath not found	
Non-refrigerated?	
No	<ul style="list-style-type: none"> Note before connecting ensure the power is off. ▲
Yes	<ul style="list-style-type: none"> for non-refrigerated baths, select Yes
Calibration Saved	<ul style="list-style-type: none"> press enter (↵)
Calibration Restored	<ul style="list-style-type: none"> press enter (↵)
RS232 Enabled	<ul style="list-style-type: none"> press enter (↵)
RS485 Enabled	<ul style="list-style-type: none"> press enter (↵)
Remote Start Enabled	<ul style="list-style-type: none"> press enter(↵)

Checklist

Circulator will not start

Check the display for error codes, see Error Codes in this section.

Ensure the circuit protector(s) is in the on (I) position.

Make sure supply voltage is connected and matches the circulator's nameplate rating $\pm 10\%$.

No display

Recycle the circuit protector on the rear of the circulator.

Display reads 239°C

External temperature probe selected but no probe attached to circulator.

Attach external probe or select internal temperature sensor.

Circulator will not circulate process fluid

Check the reservoir level. Fill, if necessary.

Check the application for restrictions in the cooling lines.

The pump motor overloaded. The pump's internal overtemperature overcurrent device will shut off the pump causing the flow to stop. This can be caused by low fluid, debris in system, operating circulator in a high ambient temperature condition or excessively confined space. Allow time for the motor to cool down.

Ensure supply voltage matches the circulator's nameplate rating $\pm 10\%$.

Inadequate temperature control

Verify the setpoint.

For refrigerated baths, ensure the condenser is free of dust and debris.

Check the fluid concentration.

Ensure circulator installation complies with the site requirements in Section 3.

Ensure supply voltage matches the nameplate rating $\pm 10\%$.

If the temperature continues to rise, ensure your application's heat load does not exceed the rated specifications.

Enter the controller menu and ensure the ENERGY SAVER mode is on in order for the system to maintain a stable temperature.

Check for high thermal gradients (e.g., the application load is being turned on and off or rapidly changing).

Inadequate pump pressure

Ensure any user installed in-line valves are in the desired position.

Ensure the circulator's process fluid outlet is connected to the application's fluid inlet and not the outlet, see Section 3.

Ensure all connections are secure and that the proper sealant/lubricant for the fitting material is used.

Keep the distance between the circulator and the application as short as possible.

Ensure tubing is straight and without bends. If diameter reductions are required, make them at the inlet and outlet of your application, not at the circulator.

Circulator shuts down

Ensure  button wasn't accidentally pressed.

Ensure the circuit protector(s) is in the on (I) position.

Check the display for error codes.

Make sure supply voltage is connected and matches the nameplate rating $\pm 10\%$.

Restart the circulator.

Circulator fan not running

The fan comes on when internal head temp is greater than 30°C and shuts off when the temperature drops below 25°C.

USB Driver Not Recognized

If your operating system does not automatically recognize the optional driver log on to:


<http://www.ftdichip.com/FTDrivers.htm>

for instructions.

Please contact Thermo Fisher Scientific Sales Service and Customer Support if you need any additional information, see inside cover for contact instructions.


Section 7 Password

Password is used to perform a calibration or to change the circulator's PID values. It is also used to reset the circulator to factory preset values. To do this change the circulator's level from the **User** mode to the **Operator** mode.


Operator Mode

- Press  next to **Password** to display:


- Press  next to **Password** to display the Change a Value Display. Change the **Password** to **1** and press  to display:


- The circulator is now in the **Operator** mode.

- Press **esc** to display:

Reset


- If desired, press  next to **Reset** to display:

- Press  on the left of **User** adjustable settings, **PID** values or **Both** to reset factory values. **Note** The reset takes several seconds. ▲


Calibration

- Return to the Advanced menu screen to display:


Note Ensure RTAs are set to zero before performing a calibration. ▲

- To perform a calibration, press  on the right of **Calibration** to display:


- Press  on the left of the desired temperature sensor to display:

Note For ease of calibration, ensure the sensor selected is the sensor controlling the temperature. Do not calibrate the internal probe if the external probe is controlling temperature. ▲

- Press  on the right of **Calibrate** to display:

- Press  on the right of **SPx** to display the Change of Value Display. Enter the desired high of low calibration temperature. Press .


- Using the temperature display at the top right of the display, allow the temperature to stabilize for several minutes near the calibration temperature.

Internal RTD
xx.xx°C

- Once stabilized, press  on the right of either **High** or **Low**.
- Press  on the right of **SPx** to again display the Change of Value Display. Enter the other calibration temperature and repeat the procedure.
- Once both points are entered, press  on the right of **Calibrate** to complete the procedure. The calibration is stored in the circulator's non-volatile memory.


As a backup, once the calibration is completed you can store it into the circulator's memory by selecting **Save User Cal**. You can later restore that calibration by pressing **Restore User Cal**.

Another option is to restore the factory calibration values by pressing **Restore Factory Cal**.


PID Tuning


- To change the circulator's PID values press  to the right of **PID Tuning**.

- The circulator display:

- Press  to the left of either **Cool PID** or **Heat PID** to display:

- Press  on the right of the desired value to display the Change of Value Display.

Factory preset values are:

P = 0.60

I = 0.60

D = 0.00

User Mode

The circulator resets to the **User** mode by:

- changing the **Password** back to **0**.
- turning the circulator off

Note The circulator also reverts to the **User** mode whenever either the Start or Status Display is displayed continuously for 10 minutes. ▲

Appendix AC Serial Communications Protocol

Serial communication is accomplished either through the optional 9-pin Serial Communications Box or through the USB port on the immersion circulator. If your operating system does not automatically recognize the optional driver log on to: <http://www.ftdichip.com/FTDrivers.htm> for instructions.

Note This appendix assumes you have a basic understanding of communications protocols. Information on the NC, STANDARD and NAMUR protocols is available upon request.

Note Keypad operation is still available with serial communications enabled.

Note NC protocol is required to use RS485 device addressing.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Controller responses are either the requested data or an error message. The controller response must be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the controller's response. If the command is not understood, the controller responds with an error command. Otherwise, the controller responds with the requested data.

Commands are not case sensitive. Upper or lower case letters may be used. Commands are listed in the Commands Table, error responses are given in the Errors Table, and symbols are shown in the Key Table.

Key	
Symbol	Meaning
[B]	A binary value 0 or 1 (0 = Off, FALSE or Disable(d); 1 = On, TRUE or Enable(d)).
[CR]	Carriage return – used as the termination character.
[U]	Text representing the units associated with a value.
[V]	A value that can be requested in a read command or sent as part of a set command.
[V _{MAX}]	Maximum allowed value. Part of error message when set value is too high.
[V _{MIN}]	Minimum allowed value. Part of error message when set value is too low.

Value: Read commands return analog [V] or bit [B] values or settings, while set commands send analog or bit settings. Read commands return values in the same displayed precision Set command messages missing the space character between the command and the setting will be rejected, as the user's intent is unclear.

Units: A read command returning an analog [V] value or setting, will include the units [U] associated with that value or setting. A set command sending an analog value will not include the units. The units returned by the complementary read command are assumed.

Termination character: A carriage return [CR] is used to terminate command and response messages. (Typically the "Enter" key on the keyboard.)

Note The inter-character timeout (time between transmitted characters) is set to 30 seconds. Exceeding the timeout will clear the receiver buffer and require the message to be retransmitted.

Note Special characters (backspace, delete, insert, etc.) are not recognized and generate error responses.

Commands Table			
Commands All messages from master and slave are terminated with a carriage return [CR]			
Command Description	Master Sends	Sample Slave Response (echo off)	Alternate units
Read Temperature Internal	RT	[V]C	F K
Read Temperature 2 External (optional)	RT2	[V]C	F K
Read Displayed Setpoint	RS	[V]C	F K
Read Internal RTA1 – Internal RTA5	RIRTA1 – 5	[V]C	F K
Read External RTA1 – External RTA5	RERTA1 - 5	[V]C	F K
Read Setpoint X (X = 1 to 5)	RSX	[V]C	F K
Read High Temperature Fault	RHTF	[V]C	F K
Read High Temperature Warn	RHTW	[V]C	F K
Read Low Temperature Fault	RLTF	[V]C	F K
Read Low Temperature Warn	RLTW	[V]C	F K
Read Proportional Heat Band Setting	RPH	[V]%	
Read Proportional Cool Band Setting	RPC	[V]%	
Read Integral Heat Band Setting	RIH	[V]Repeats per minute	
Read Integral Cool Band Setting	RIC	[V]Repeats per minute	
Read Derivative Heat Band Setting	RDH	[V]Minutes	
Read Derivative Cool Band Setting	RDC	[V]Minutes	
Read Temperature Precision	RTP	[V]	
Read Temperature Units	RTU	[V] C,F,K	
Read Unit On	RO	[B]	
Read External Probe Enabled (optional)	RE	[B]	
Read Auto Restart Enabled	RAR	[B]	
Read Energy Saving Mode	REN	[B]	
Read Time	RCK	hh:mm:ss	
Read Date	RDT	mm/dd/yyyy or dd/mm/yyyy	
Read Date Format	RDF	mm/dd/yyyy or dd/mm/yyyy	
Read Ramp Status (optional)	RRS	Stopped, Running, Paused	
Read Firmware Version	RVER	[V]	
Read Firmware Checksum	RSUM	[V]	
Read Unit Fault Status	RUFS	[V1, V2 , V3, V4, V5] See page A-5	
Read Ramp Program (optional)	RRP	[V1, V2 , V3, V4, V5] See SRP on page A-4	

Commands All messages from master and slave are terminated with a carriage return [CR]		
Command Description	Master Sends	Sample Slave Response
Set Displayed Setpoint	SS [V]	OK
Set Internal RTA1 – Internal RTA5	SIRTA1 – SIRTA5 [V]	OK
Set External RTA1 – External RTA5	SERTA1 – SERTA5 [V]	OK
Set Setpoint X (X = 1 to 5)	SSX [V]	OK
Set High Temperature Fault	SHTF [V]	OK
Set High Temperature Warning	SHTW [V]	OK
Set Low Temperature Fault	SLTF [V]	OK
Set Low Temperature Warning	SLTW [V]	OK
Set Bath Faults Status (Clears all errors/faults)	SUFS	OK
Set Proportional Heat Band Setting	SPH [V]	OK
Set Proportional Cool Band Setting	SPC [V]	OK
Set Integral Heat Band Setting	SIH [V]	OK
Set Integral Cool Band Setting	SIC [V]	OK
Set Derivative Heat Band Setting	SDH [V]	OK
Set Derivative Cool Band Setting	SDC [V]	OK
Set Temperature Resolution	STR [V]	OK
Set Temperature Units	STU [V] C,F,K	OK
Set Unit On Status	SO [B]	OK
Set External Probe On Status (optional)	SE [B]	OK
Set Auto Restart Enabled	SAR [B]	OK
Set Energy Saving Mode	SEN [V]	OK
Set Pump Speed	SPS [V] L,M,H	OK
Set Ramp Number (optional)	SRN [V]	OK

Commands All messages from master and slave are terminated with a carriage return [CR]		
Command Description	Master Sends	Sample Slave Response
Set Ramp Program this will load defaults for all steps declared by V1	SRP [v1,v2,v3,v4,v5] V1 #of steps V2 variance V3 cycles V4 start temp V5 end state S=stop M=maintain	OK Note: If unit allows more than one ramp the ramp number must first be set.
Set Ramp Step	SRS[v1,v2,v3, v4] V1 step # V2 end temp V3 duration V4 assured soak	OK Note: If unit allows more than one ramp the ramp number must first be set and SRP defined. (0,1)
Set Ramp On Status	SRO [V] S,E,P Start, End, Pause/resume	OK Note: (P) pause toggles Pause/resume

Error Table	
Errors	
Error Description	Slave Responds
Not defined, not implemented or incorrectly formatted	? Unsupported command
Extra characters...	? Format error
Set value too high	? Maximum allowed is [V _{MAX}]
Set value too low	? Minimum allowed is [V _{MIN}]
Argument to binary set command not 0 or 1	? Value must be 0 or 1
Set command attempted while in read only mode	? Mode is read only
Set command failed (e.g. SO 1 with low level)	? Failed

RUFS Read Unit Fault Status

This command returns 5 values. These are decimal representations of hexadecimal values. Each individual bit of the value represents a different warning, fault or status.

decimal	hex	B7	B6	B5	B4	B3	B2	B1	B0
1	1	0	0	0	0	0	0	0	1
2	2	0	0	0	0	0	0	1	0
4	4	0	0	0	0	0	1	0	0
8	8	0	0	0	0	1	0	0	0
16	10	0	0	0	1	0	0	0	0
32	20	0	0	1	0	0	0	0	0
64	40	0	1	0	0	0	0	0	0
128	80	1	0	0	0	0	0	0	0

Value	Description of bits	Value	Description of bits
V1	B0 – B5 unused B6 rtd1 shorted B7 rtd1 open	V2	B0 HTC fault B1 high RA temperature fault B2 – B7 unused
V3	B0 low level warn B1 low temperature warn B2 high temperature warn B3 low level fault B4 low temperature fault B5 high temperature fault B6 low temperature fixed fault B7 high temperature fixed fault	V4	B0 PWM heat duty cycle > 0 B1 compressor on/off B2 pump on status B3 circulator on status B4 circulator stopping B5 circulator fault status B6 unused B7 beeper on status
V5	B0 pump speed fault B1 MOL fault B2 HPC fault B3 cool icon on steady (circulator is cooling at max capacity) B4 cool icon flashing (circulator is cooling) B5 heat icon on steady B6 heat icon flashing B7 external sensor controlling		

Refer to Key table on page 1 for explanation of symbols and their meanings.

Examples:

Read Temperature:

Host

R	T		CR
Command			[CR]

Controller:

2	0	.	0	C	CR
[V]			[U]	[CR]	

Set Setpoint:

Host

S	S		2	0	CR
Command			[V]	[CR]	

Controller:

O	K	CR
Command Accepted		[CR]

Read Temperature 2:

Host:

R	T	2	CR		
2	0	.	0	C	[CR]

Controller:

Set Setpoint to -22°C when minimum allowed is -20°C: Minimum allowed is $[V_{MIN}]$

Host:

S	S		-	2	2	CR															
?		M	i	n	i	m	u	m	a	l	l	o	w	e	d	i	s	-	2	0	CR

Controller:

Warranty

Thermo Fisher Scientific warrants for 36 months from date of shipment the Thermo Scientific PREMIUM series of Immersion Circulators, ARCTIC and GLACIER refrigerated bath circulators, and SAHARA heated bath circulators according to the following terms.

Any part of the unit manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The unit, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the unit to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any unit. With respect to units that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any unit that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the unit that is the result of improper installation or maintenance, or to any unit that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any unit that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the unit or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any unit without assuming any obligation to modify any unit previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE UNIT, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to units sold in the United States. Any units sold elsewhere are warranted by the affiliated marketing company of Thermo Fisher Scientific. This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

Thermo Fisher Scientific
168 Third Avenue
Waltham, Massachusetts 02451
United States

www.thermofisher.com