

# CO<sub>2</sub> Independent Medium

## Description

CO<sub>2</sub> Independent Medium is a non-HEPES proprietary medium used for supporting cell growth for a variety of suspension and adherent mammalian cells such as epithelial, fibroblast, and lymphoid cell lines without a CO<sub>2</sub> incubator. CO<sub>2</sub> Independent Medium is ideally suited for transporting cells or tissue, handling of mouse embryos under atmospheric conditions, and for use in toxicological and/or virological procedures where there is a risk of aerosol contamination or infection. CO<sub>2</sub> Independent Medium contains a unique buffering system composed of mono and dibasic sodium phosphate and β-glycerophosphate capable of maintaining long term pH stability under atmospheric CO<sub>2</sub> (0.04%). A small amount of sodium bicarbonate has been included in the formulation to meet essential bicarbonate dependent functions. No synthetic buffers are utilized, thus eliminating any cytotoxic effects associated with such buffering systems. Additionally, CO<sub>2</sub> Independent Medium has been formulated with components that enhance cellular production and utilization of CO<sub>2</sub> such that an exogenous source of CO<sub>2</sub> is not required for the maintenance of CO<sub>2</sub> dependent cellular functions.

Product	Catalog no.	Amount	Storage	Shelf life*
CO <sub>2</sub> Independent Medium	18045-088	500 mL	2°C to 8°C; Protect from light	12 months
	18045-054**	500 mL		
	18045-070**	10 × 500 mL		

\* Shelf Life duration is determined from Date of Manufacture.

\*\* Note: For European customers only.

## Product use

For Research Use Only. Not for use in diagnostic procedures.

## Important information

- CO<sub>2</sub> Independent Medium requires supplementation with GlutaMAX™-I or L-glutamine and when necessary, fetal bovine serum (FBS).
- Cultures can be placed in a 5% CO<sub>2</sub> incubator utilizing either an open or closed culture system without any deleterious effects to cellular growth.

## Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Caution: Human origin materials are non-reactive (donor level) for anti-HIV 1 & 2, anti-HCV, and HB<sub>s</sub>Ag. Handle in accordance with established bio-safety practices.

## Prepare medium

1. Aseptically add 20 mL GlutaMAX™-I (200 mM) or L-glutamine (200 mM) to 880 mL CO<sub>2</sub> Independent Medium before use.
2. Aseptically add 100 mL FBS to the medium before use.
3. Add antibiotics, if required.

## Culture conditions

**Media:** complete CO<sub>2</sub> Independent Medium

**Cells:** mammalian

**Culture type:** Adherent or Suspension

**Culture vessels:** T-Flasks

**Temperature range:** 36°C to 38°C

**Incubator atmosphere:** Humidified atmosphere of 0% CO<sub>2</sub> in air. Open or closed culture systems minimizing exposure of cultures to light.

## Adaptation of cells to CO<sub>2</sub> Independent Medium

For maximum growth performance, some cell lines may require either direct or sequential adaptation to CO<sub>2</sub> Independent Medium. In either case, the pre-adapted cell line should be in mid-logarithmic growth phase with high (>90%) viability. Success of the adaptation procedure will depend on the cell line being used and the culture conditions employed. It is recommended that the user first evaluate this product with unadapted cells since not all cell lines will require adaptation. If the growth assays employed are conducted in a closed culture system, stock cultures can be directly adapted and maintained in CO<sub>2</sub> Independent Medium using a closed, 0% CO<sub>2</sub> system.

### Direct adaptation

1. Inoculate cultures at normal seeding densities and incubate using a closed cap in a humidified (37°C) incubator with 0% CO<sub>2</sub>.
2. Monitor cell growth daily and subculture cells when they reach 80–90% confluency. Subsequent passages should utilize a humidified (37°C) 0% CO<sub>2</sub> atmosphere with an open cap vessel.

If the cell cultures fail to maintain acceptable growth and viability over 3–5 passages during direct adaptation, use the sequential adaptation method.

### Sequential adaptation

**Note:** Multiple passages at each step may be required.

1. Inoculate cells at normal seeding densities into a 50:50 ratio (v/v) of CO<sub>2</sub> Independent Medium and the currently utilized medium.
2. Maintain cultures under an open cap in a humidified (37°C) incubator with 0% CO<sub>2</sub>.
3. Monitor cell growth daily and subculture cells when they reach 80–90% confluency into a 75:25 ratio (v/v) of CO<sub>2</sub> Independent Medium and the currently utilized medium.
4. Monitor cell growth daily and subculture cells when they reach 80–90% confluency into 100% CO<sub>2</sub> Independent Medium.

Subsequent subculturing should use 100% CO<sub>2</sub> Independent Medium and maintained as described in the preceding adaptation procedure.

## Related products

Product	Catalog no.
L-Glutamine, 200 mM (100X), liquid	25030
GlutaMAX™-I, 200 mM (100X), liquid	35050
Certified FBS, Heat Inactivated, US	10082
Antibiotic-Antimycotic (100X), liquid	15240
Fungizone® Antimycotic, liquid	15290
Gentamicin	15750
Penicillin-Streptomycin, liquid	15140
Dulbecco's Phosphate Buffered Saline, without calcium and magnesium	14190
TrypLE™ Express (1X), liquid, without Phenol Red	12563
Trypsin-EDTA, 1X	25300
Trypan Blue Stain	15250
Countess® Automated Cell Counter	C10227

## Explanation of symbols and warnings

The symbols present on the product label are explained below:

Temperature Limitation	Manufacturer	Batch code	Use By:	Catalog number
Caution, consult accompanying documents	Consult instructions for use	Keep away from light	Sterilized using aseptic processing techniques	

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

## Important licensing information

This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

## References

Battista, P.J., DiSorbo, D.M. and Weiss, S.A. Development of a carbon dioxide independent medium. *In Vitro Cell. Dev. Biol.* 27:120A (1991).

For additional technical information such as Safety Data Sheets (SDS), Certificates of Analysis, visit [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support)  
For further assistance, email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

©2014 Thermo Fisher Scientific Inc. All rights reserved.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

DISCLAIMER - LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

[www.lifetechnologies.com](http://www.lifetechnologies.com)

*life*  
technologies™