

HIF-1A ELISA Kit

EHIF1A EHIF1A2 EHIF1A5

Number	Description
EHIF1A	HIF-1A ELISA Kit, sufficient reagents for 96 determinations
EHIF1A2	HIF-1A ELISA Kit, sufficient reagents for 2 × 96 determinations
EHIF1A5	HIF-1A ELISA Kit, sufficient reagents for 5 × 96 determinations

Kit Contents	EHIF1A	EHIF1A2	EHIF1A5
Antibody Coated Plate, 96-well plate	1 each	2 each	5 each
Lyophilized Recombinant HIF-1A Standard	2 vials	4 vials	10 vials
Standard Diluent	14mL	2 × 14mL	5 × 14mL
30X Wash Buffer	50mL	2 × 50mL	5 × 50mL
Biotinylated Antibody Reagent	12mL	2 × 12mL	5 × 12mL
Streptavidin-HRP Concentrate	75µL	2 × 75µL	5 × 75µL
HRP Dilution Buffer	14mL	2 × 14mL	5 × 14mL
TMB Substrate	13mL	2 × 13mL	5 × 13mL
Stop Solution, contains 0.16M sulfuric acid	13mL	2 × 13mL	5 × 13mL
Adhesive Plate Covers	6 each	12 each	30 each

For research use only. Not for use in diagnostic procedures.

Storage: Upon receipt store at 2-8°C.

Table of Contents












Introduction	2
Procedure Summary.....	2
Additional Materials Required.....	2
Precautions.....	3
Sample Preparation.....	3
Reagent Preparation.....	4
Assay Procedure	5
Performance Characteristics	6

Introduction

The Invitrogen™ HIF-1 Alpha ELISA Kit measures human hypoxia-inducible factor 1-alpha (HIF1A) protein in serum, plasma and cell culture lysates using the quantitative sandwich ELISA method. HIF1A standards and samples are captured by a polyclonal HIF1A antibody on the pre-coated plate and detected using a biotinylated monoclonal HIF1A antibody reactive to epitopes other than the capture antibody. The biotinylated detection antibody is then bound to streptavidin-HRP, which catalyzes the conversion of TMB to a colored derivative. Color development is linear for the assay's dynamic range and directly proportional to the amount of HIF1A present in the sample.

The HIF1A protein is the 92kDa alpha subunit of transcription factor hypoxia-inducible factor-1 (HIF-1), which is a heterodimer composed of an alpha and a beta subunit. HIF-1 functions as a master regulator of cellular and systemic homeostatic response to hypoxia by activating transcription of many genes, including those involved in energy metabolism, angiogenesis, apoptosis, and other genes whose protein products increase oxygen delivery or facilitate metabolic adaptation to hypoxia. HIF-1 thus plays an essential role in embryonic vascularization, tumor angiogenesis and pathophysiology of ischemic disease. Activation requires recruitment of transcriptional co-activators such as CREBBP and EP300. Activity is enhanced by interaction with either NCOA1 or NCOA2. Interaction with redox regulatory protein APEX seems to activate CTAD and potentiates activation by NCOA1 and CREBBP. The HIF-1 heterodimer binds to core DNA sequence 5'-[AG]CGTG-3' within the hypoxia response element (HRE) of target gene promoters.

Procedure Summary

-  1. Add 50µL of Standards or samples to each well.
-  2. Cover plate. Incubate at room temperature (20-25°C) for 2 hours on a plate shaker.
-  3. Wash plate THREE times. Add 50µL of Biotinylated detection antibody.
-  4. Cover plate. Incubate at room temperature (20-25°C) for 1 hour on a plate shaker.
-  5. Wash plate THREE times. Add 50µL of prepared Streptavidin-HRP to each well.
-  6. Cover plate. Incubate at room temperature (20-25°C) for 30 minutes on a plate shaker.
-  7. Wash plate THREE times. Add 100µL of TMB substrate.
-  8. Develop plate in the dark at room temperature for 30 minutes.
-  9. Stop the reaction by adding 100µL of Stop Solution to each well.
-  10. Measure absorbance on a plate reader at 450nm minus 550nm.
-  11. Calculate results using graph paper or curve-fitting statistical software.

Additional Materials Required

- Precision pipettors with disposable plastic tips to deliver 5-1000µL and plastic pipettes to deliver 5-15mL
- Glass or plastic 2L container to prepare Wash Buffer
- Squirt wash bottle or an automated 96-well plate washer
- 1.5mL polypropylene or polyethylene tubes to prepare standards – do not use polystyrene, polycarbonate or glass tubes
- Disposable reagent reservoirs (Product No. 15075)
- Microcentrifuge
- 15mL plastic tube to prepare Streptavidin-HRP Solution

-
- A standard ELISA reader for measuring absorbance at 450nm and 550nm. If a 550nm filter is not available, the absorbance can be measured at 450nm only. Refer to the instruction manual supplied with the instrument being used.
 - Graph paper or a computerized curve-fitting statistical software package

Precautions

- All samples and reagents must be at room temperature (20-25°C) before use in the ELISA.
- Review all instructions carefully and verify all components against the kit contents list (page 1) before beginning.
- Thaw samples at room temperature. Do not use a water bath to thaw samples.
- When preparing standard curve and sample dilutions in culture medium, use the same medium used to culture the cells. For example, if RPMI with 10% fetal calf serum (FCS) was used to culture cells, then use RPMI with 10% FCS to dilute the standard and samples. Do NOT use RPMI without serum supplement.
- To avoid cross-contamination, always use a new disposable reagent reservoir and new disposable pipette tips for each transfer. Also, use a new adhesive plate cover for each incubation step.
- Once reagents have been added to the plate, take care NOT to let the plate dry at any time during the assay.
- Avoid microbial contamination of reagents.
- Avoid exposing reagents to excessive heat or light during storage and incubation.
- Discard unused ELISA Kit components after assay completion. Do not mix reagents from different kit lots.
- Do not use glass pipettes to measure the TMB Substrate Solution. Do not contaminate the TMB Substrate Solution; if it is blue before use, DO NOT USE THE SOLUTION.
- Individual components can contain antibiotics or preservatives. Wear gloves while performing the assay to avoid contact with samples and reagents. Follow proper disposal procedures.

Sample Preparation

A. Sample Handling

- Serum; EDTA, heparin plasma; and culture lysates may be tested in this ELISA.
- 50µL per well of serum, plasma or culture lysate is required.
- Store samples at 2-8°C and assay within 24 hours. For long-term storage, aliquot and freeze samples at -70°C.
- Avoid repeated freeze-thaw cycles when storing samples.
- Test samples and standards must be assayed in duplicate or triplicate each time the ELISA is performed.
- Mix samples by gently inverting tubes.
- If samples are clotted, grossly hemolyzed, lipemic or microbially contaminated, or if there is any question about the integrity of a sample, make a note on the template and interpret the results with caution.

B. Sample Dilution

- If the HIF-1A concentration from the sample will exceed the highest point of the standard curve (i.e., 20,000pg/mL), prepare one or more five-fold dilutions of the test sample. When testing **culture lysates**, prepare the serial dilutions using the culture medium. When testing **serum or plasma**, prepare the serial dilutions using the provided Standard Diluent. For example, a five-fold dilution is prepared by adding 0.1mL (100µL) of test sample to 0.4mL (400µL) of appropriate diluent. Mix samples thoroughly between dilutions before assaying.

Reagent Preparation

For procedural differences when using partial plates, look for **(PP)** throughout these instructions.

Wash Buffer

1. Label a new glass or plastic 2L container "Wash Buffer." The 30X Wash Buffer may have a cloudy appearance.
2. Add entire contents of the 30X Wash Buffer (50mL) bottle to a 2L container and dilute to a final volume of 1.5L with ultrapure water. Mix thoroughly.

(PP) When using partial plates, store the reconstituted Wash Buffer at 2-8°C.

Note: Wash Buffer must be at room temperature before use in the assay. Do not use Wash Buffer if it becomes visibly contaminated during storage.

Standards

(PP) Reconstitute with Standard Diluent and use one vial of the lyophilized Standard per partial plate.


1. Prepare Standards just before use and use within one hour of reconstitution. Do not store reconstituted standards.
2. When testing **culture lysate samples**, reconstitute Standard with Standard Diluent. Reconstitution volume is stated on the Standard vial label. The Standard will dissolve in ~1 minute. Mix by gently inverting vial. Use the sample culture medium to prepare standard curve dilutions.

When testing **serum or plasma samples**, reconstitute Standard with Standard Diluent. Reconstitution volume is stated on the Standard vial label. The Standard will dissolve in ~1 minute. Mix by gently inverting vial. Use the Standard Diluent provided to prepare the standard curve serial dilutions.

When testing **serum, plasma and cell culture lysate samples on the same plate**, validate the media to establish if the same standard curve can be used for the different sample types. Prepare a standard curve (including a zero/blank) using culture medium to reconstitute and dilute the Standard. Use medium containing serum or other protein to maximize stability of the HIF-1A. Perform this curve in parallel with a standard curve reconstituted in Standard Diluent and diluted in the Standard Diluent provided. If the OD values of the two curves are within 10% of the mean for both curves, then the assay can be performed with Standard Diluent, whether you are testing culture lysate, plasma or serum samples.

3. Label eight tubes, one for each standard curve point: 20,000, 8,000, 3,200, 1,280, 512, 204.8, 81.92, and 0pg/mL. Prepare a 1:2 dilution of the reconstituted standard and then 1:2.5 serial dilutions for the remaining standards as follows:
4. Pipette 250µL of appropriate diluent into the first tube (i.e., 20,000pg/ml).
5. Pipette 300µL of appropriate diluent into all of the remaining tubes.
6. Pipette 250µL of the reconstituted standard into the first tube (i.e., 20,000pg/mL) and mix.
7. Pipette 200µL of this dilution into the second tube (i.e., 8,000pg/mL) and mix.
8. Repeat the serial dilutions (using 200µL) five more times to complete the standard curve points. These concentrations, 20,000, 8,000, 3,200, 1,280, 512, 204.8, 81.92, and 0pg/mL are the standard curve.

Serial Dilutions



Assay Procedure

A. Calibrator and Sample Incubation

- **(PP)** Determine the number of strips required. Leave these strips in the plate frame. Tightly seal remaining unused strips in the provided foil pouch with desiccant and store at 2-8°C. After completing the assay, retain the plate frame for second partial plate. When using the second partial plate, place strips securely in the plate frame.
- Use the Plate Template provided to record the locations of the zero standard (blank or negative control), HIF-1A standards and samples. Perform seven standard points and one blank in duplicate with each series of unknown samples.
- If using a multichannel pipettor, use a new reagent reservoir to add the Biotinylated Antibody Reagent. Remove from the vial only the amount required for the number of strips being used.

1. Add 50µL of diluted standards and test samples to each well. Mix well by gently tapping the plate several times.
2. Carefully cover the plate with an adhesive plate cover. Ensure all edges and strips are tightly sealed by running your thumb over the edges and down each strip. Incubate for 2 hours at room temperature (20-25°C) on a plate shaker.
3. Carefully remove the adhesive plate cover. Wash plate THREE times with Wash Buffer as described in the Plate Washing Section (Section B).

B. Plate Washing

1. Gently squeeze the long side of the plate frame before washing to ensure all strips securely remain in the frame.
2. Empty plate contents. Use a squirt bottle to vigorously fill each well completely with Wash Buffer, then empty plate contents. Repeat procedure two additional times for a total of THREE washes. Blot plate onto paper towels or other absorbent material.

Note: For automated washing, aspirate all wells and wash THREE times with Wash Buffer, overfilling wells with Wash Buffer. Blot plate onto paper towels or other absorbent material.

C. Biotinylated Antibody Reagent Incubation

1. Add 50µl of Biotinylated Antibody Reagent to each well.
2. Carefully attach a new adhesive plate cover, ensuring all edges and strips are tightly sealed. Incubate the plate for 1 hour at room temperature (20-25°C) on a plate shaker.
3. Carefully remove the adhesive plate cover, discard the plate contents and wash THREE times as described in the
4. Plate Washing Section.

D. Streptavidin-HRP Solution Preparation and Incubation

- Prepare Streptavidin-HRP Solution **just before use**. Do not prepare more solution than required.
 - Do not store the prepared Streptavidin-HRP Solution.
 - Use a 15mL plastic tube to prepare the Streptavidin-HRP Solution.
 - Use a new reagent reservoir and pipette tips when adding the prepared Streptavidin-HRP Solution.
1. Briefly centrifuge the Streptavidin-HRP Concentrate to force entire vial contents to the bottom.
 2. **(PP)** Use only the Streptavidin-HRP Solution amount required for the number of strips being used. For each strip, mix 2.5µL of Streptavidin-HRP Concentrate with 1mL of HRP Dilution Buffer. Store Streptavidin-HRP Concentrate reserved for additional strips at 2-8°C.
 3. For one complete 96-well plate, add 30µL of Streptavidin-HRP Concentrate to 12mL of HRP Dilution Buffer and mix gently.
 4. Add 50µL of prepared Streptavidin-HRP Solution to each well.
 5. Carefully attach a new adhesive plate cover, ensuring all edges and strips are tightly sealed. Incubate the plate for 30 minutes at room temperature (20-25°C) on a plate shaker.
 6. Carefully remove the adhesive plate cover, discard the plate contents and wash THREE times as described in the Plate Washing Section.

E. Substrate Incubation and Stop Solution Addition

- Use new disposable reagent reservoirs when adding the TMB Substrate Solution and Stop Solution.
 - From the bottle, dispense ONLY the required amount of 100µL per well for the number of wells being used. Do not use a glass pipette to measure the TMB Substrate Solution.
 - **(PP)** Do not combine leftover substrate with that reserved for the second partial plate. Take care not to contaminate the remaining TMB Substrate.
1. Pipette 100µL of TMB Substrate Solution into each well.
 2. Allow color reaction to develop at room temperature in the dark for 30 minutes. Do not cover plate with aluminum foil or a plate sealer. The substrate reaction yields a blue solution that turns yellow when Stop Solution is added.
 3. After 30 minutes, stop the reaction by adding 100µL of Stop Solution to each well.


F. Absorbance Measurement

Evaluate the plate within 30 minutes of stopping the reaction. Measure absorbance on an ELISA plate reader set at 450nm and 550nm. Subtract 550nm values from 450nm values to correct for optical imperfections in the microplate. If an absorbance at 550nm is not available, measure the absorbance at 450nm only.

Note: When the 550nm measurement is omitted, absorbance values will be higher.

G. Calculation of Results

- The standard curve is used to determine HIF-1A amount in an unknown sample. Generate the standard curve by plotting the average absorbance obtained for each Standard concentration on the vertical (Y) axis vs. the corresponding HIF-1A concentration (pg/mL) on the horizontal (X) axis.
- Calculate results using graph paper or curve-fitting statistical software. Determine the HIF-1A amount in each sample by interpolating from the absorbance value (Y-axis) to HIF-1A concentration (X- axis) using the standard curve.
- If the test sample was diluted, multiply the interpolated value obtained from the standard curve by the dilution factor to calculate pg/mL of HIF-1A in the sample.
- Absorbance values obtained for duplicates should be within 10% of the mean value. Carefully consider duplicate values that differ from the mean by greater than 10%.



Performance Characteristics

Sensitivity: ≤30pg/mL

The sensitivity, or Lower Limit of Detection (LLD), was determined by assaying replicates of zero and the standard curve. The mean signal of zero + two standard deviations read in dose from the standard curve is the LLD. This value is the smallest dose that is not zero with 95% confidence.

Assay Range: 81.92 – 20,000pg/mL

Suggested standard curve points are 20,000, 8,000, 3,200, 1,280, 512, 204.8, 81.92, and 0pg/mL

Reproducibility:

Intra-Assay CV: <10%

Inter-Assay CV: <10%

Specificity: This ELISA is specific for the measurement of natural and recombinant human HIF-1A. The kit shows cross reactivity with mouse HIF-1A and is known to cross react with rat HIF-1A.


Calibration: The standards in this assay were calibrated to an in-house standard.

Expected Values:

Cell Culture Lysates: NIH 3T3 cells were cultured in DMEM cell culture media with 10% fetal calf serum. Cultured NIH 3T3 cells were hypoxia induced via deferoxamine at a concentration of 100µM. Samples were collected after 24hr proliferation, cells were lysed by trypsinization, and 70pg/mL of HIF-1A was measured.

Precision: The standard curve's average coefficients of variation (CV) are plotted against HIF-1A concentration (pg/mL). The points represent samples evaluated in replicates of two in 10 representative pre-coated plates.

Recovery: Five samples of various matrices type: Human serum, EDTA plasma, and heparin plasma were initially diluted 1:10 with Sample Diluent, spiked with known concentrations of HIF-1A, and assayed to calculate mean % recovery.



Spiked concentration (pg/mL)	Mean % recovery		
	Serum (n=5)	Plasma EDTA (n=5)	Plasma heparin (n=5)
4000	101	93	89
2000	103	88	89
1000	88	87	92

Dilution Linearity: Five individual donor samples for each sample type: Human serum, EDTA plasma and heparin plasma were initially diluted 1:10 with Sample Diluent and spiked with a high concentration of HIF-1A. Samples were further serially diluted with appropriate Sample Diluent and assayed.

Dilution factor	Mean % recovery		
	Serum (n=5)	Plasma EDTA (n=5)	Plasma heparin (n=5)
1:2	104	103	98
1:4	105	93	95
1:8	110	88	90

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

Product label explanation of symbols and warnings

	Catalog Number		Batch code		Temperature limitation		Use by		Manufacturer		Consult instructions for use		Caution, consult accompanying documents
--	----------------	--	------------	--	------------------------	--	--------	--	--------------	--	------------------------------	--	---

Manufacturer's address: Bender MedSystems GmbH | Campus Vienna Biocenter 2 | 1030 Vienna, Austria
 The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Plate Templates

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												