

Click-iT[®] EdU Imaging Kit with Alexa Fluor[®] 488, 594, and 647 Azides

Cat. no. C10086

Click-iT[®] EdU cell proliferation assays eliminate DNA denaturation steps required by BrdU assays, steps that can damage sample morphology and integrity. This trial-size kit provides sufficient EdU to label up to 6 coverslips and enough Alexa Fluor[®] 488, 594, and 647 azide to detect new DNA synthesis on 2 coverslips with each fluorophore (Figure 1).

Kit contents:

- EdU (5-ethynyl-2'-deoxyuridine, purple cap)
- Alexa Fluor[®] 488 azide (green cap)
- Alexa Fluor[®] 594 azide (orange cap)
- Alexa Fluor[®] 647 azide (red cap)

Required but not included:

- PBS
- 0.5% Triton[®] X-100 in PBS
- 3% BSA in PBS
- 3.7% paraformaldehyde in PBS

- CuSO₄ (large bttle, clear cap)
- Click-iT[®] EdU reaction buffer concentrate (blue cap)
- Click-iT[®] EdU buffer additive (white cap)
- Deionized water or 18 megohm purified water
- Complete medium
- 6-well plate
- 18 × 18 mm sterile coverslips

Click-iT[®] EdU compatibility

Fluorescent molecule	Notes	
Qdot [®] nanocrystals	Use Qdot® nanocrystals after the Click-iT® detection reaction	
Fluorescent proteins (GFP)	Use anti-GFP antibodies* before the Click-iT [®] detection reaction or Use organic dye-based reagents for protein expression detection	
Organic dyes (i.e., Alexa Fluor® dyes, fluorescein (FITC))	These fluorescent molecules are completely compatible with $Click\text{-}iT^{\textcircled{B}}EdU$	
TC-FlAsH [™] /TC-ReAsH [™] reagents	Detect the tetracysteine (TC) tag before the Click-iT [®] detection reaction	

* Not all anti-GFP antibodies recognize the same antigen site. Rabbit and chicken anti-GFP antibodies perform well, whereas mouse monoclonal antibodies tested do not generate an acceptable amount of fluorescence, and are not recommended for this application.

Experimental protocol

Preparing cells

- Plate cells on coverslips at desired density and allow to recover overnight before additional treatment.
- 1.2 Treat cells as desired.

Labeling cells with 10 μM EdU

- 2.1 Transfer coverslips into the 6-well plate already containing medium so that each well contains a single coverslip.
- 2.2 Prepare a 20 μ M solution of EdU in prewarmed complete medium. The EdU is provided as a 10 mM solution in DMSO.
- 2.3 Add an equal volume of the 20 μM EdU solution to the volume of medium containing the cells to obtain a 10 μM solution.
- 2.4 Incubate for 15-60 minutes.


Cell fixation and permeabilization

- 3.1 Remove medium and add 1 mL 3.7% paraformaldehyde in PBS (fixative) to each well.
- 3.2 Incubate for 15 minutes at room temperature.
- Remove fixative and wash cells 2X with 1 mL 3% BSA in PBS (wash solution).
- 3.4 Remove wash solution. Add 1 mL 0.5% Triton[®] X-100 in PBS (permeabilization buffer) to each well.
- 3.5 Incubate for 20 minutes.

EdU detection for 2 coverslips with Alexa $\mathsf{Fluor}^{\circledast}$ dye of choice

- 4.1 Prepare 1X Click-iT $^{\circ}$ EdU buffer additive by adding 1.0 mL of dH $_2$ O to white-capped vial.
- 4.2 Prepare Click-iT[®] reaction cocktail by adding the following amounts to the Alexa Fluor[®] amber vial:
 - 100 µL Click-iT[®] reaction buffer (bluecapped vial)
 - 800 µL CuSO₄ (large clear bottle)
 - 100 µL 1X Click-iT[®] reaction buffer additive (white-capped vial; prepared in step 4.1)

Use this within 15 minutes after preparation.



- 4.3 Remove permeabilization buffer and wash cells 2X with 1 mL 3% BSA in PBS, then remove the wash solution.
- 4.4 Add 0.5 ml Click-iT[®] reaction cocktail to each well. Rock plate to evenly distribute solution.
- 4.5 Incubate for 30 minutes, protected from light.
- 4.6 Remove Click-iT[®] reaction cocktail and wash 2X with 1 mL 3% BSA in PBS.

Antibody and DNA staining (optional)

- Stain with primary and secondary antibodies.
- 5.2 Stain with DNA counterstain.

Imaging

6.1 Use filters appropriate for the Alexa Fluor[®] dye, live-cell stains, secondary antibodies, and DNA counterstain (Figures 2–5).




Figure 1-Detection of the incorporated EdU with the Alexa Fluor[®] 488 azide versus incorporated BrdU with an anti-BrdU antibody. The small size of the Alexa Fluor[®] 488 azide eliminates the need to denature DNA in order for the detection reagent to gain access to the modified base.




Figure 2-Click-iT[®] EdU: Amazing imagery following EdU administration *in vivo*. Mice were injected with EdU at 0.3 mg /10 g mouse weight. An intestinal section was stained with green fluorescent Click-iT[®] EdU Alexa Fluor[®] 488 Imaging Kit (C10337) and cells that incorporated EdU were visualized with green fluorescence. Nuclei are counterstained with blue fluorescent DAPI (D1306) and red-fluorescence is from autofluorescence.




Figure 3-Eliminate long detection procedures with Click-iT[®] EdU. Rat Ileum tissue section detected with red-fluorescent Click-iT[®] EdU Alexa Fluor[®] 594 Imaging Kit (C10339). EdU staining complete in 80 minutes, while BrdU protocols require harsh permeabilization and overnight anti-BrdU detection. Nuclei are counterstained with blue fluorescent Hoechst 33342 (H1399).




Figure 4—Click-iT[®] EdU compatibility with fluorescent proteins and other live-cell stains. HeLa cells transduced with Organelle Lights[™] NE GFP prior to being incubated with EdU and MitoTracker[®] Deep Red FM (M22426). GFP expressed in the nuclear envelope was detected with anti-green fluorescent protein, rabbit serum (A6455) and visualized with an Alexa Fluor[®] 488 goat anti-rabbit IgG antibody (A11034). Green fluorescence is also seen in the endoplasmic reticulum as it is formed from the nuclear envelope. Proliferating cells were detected with a Click-iT[®] EdU Alexa Fluor[®] 594 Imaging Kit (red-fluorescence, C10339). Mitochondria are pseudocolored pink.




Figure 5—Click-iT[®] EdU requires mild fixation; enabling truly spectacular multicolor imagery. Muntjac cells were treated with 10 μ M EdU. EdU incorporated into newly synthesized DNA was detected by the far redfluorescent Click-iT[®] EdU Alexa Fluor[®] 647 Imaging Assay (C10340). Tubulin was labeled with anti-tubulin antibody (A11126) and visualized with Alexa Fluor[®] 350 goat anti-mouse IgG antibody (A21049). The Golgi complex was stained with the green-fluorescent Alexa Fluor[®] 488 conjugate of lectin HPA from *Helix pomatia* (edible snail; L11271), and peroxisomes were labeled with an orange-fluorescent Alexa Fluor[®] 555 donkey anti-rabbit IgG antibody (A31572).

Product	Quantity 1 kit	Cat. no. C10086
<code>Click-iT$^{\circ}$ EdU Imaging Kit with Alexa Fluor$^{\circ}$ 488, 594, and 647 Azides</code>		
Click-iT® EdU Alexa Fluor® 488 Imaging Kit *for 50 coverslips*	1 kit	C10337
Click-iT® EdU Alexa Fluor® 555 Imaging Kit *for 50 coverslips*	1 kit	C10338
Click-iT® EdU Alexa Fluor® 594 Imaging Kit *for 50 coverslips*	1 kit	C10339
Click-iT® EdU Alexa Fluor® 647 Imaging Kit *for 50 coverslips*	1 kit	C10340
EdU (5-ethynyl-2'-deoxyuridine)	50 mg	A10044

Ordering information

For more information, please visit www.lifetechnologies.com/edu.

Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at **www.lifetechnologies.com/termsandconditions**.

If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Research Use Only. Not for use in diagnostic procedures.

DISCLAIMER - LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLIDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITUR, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

©2012 Life Technologies Corporation. All rights reserved.

The trademarks mentioned herein are the property of Life Technologies Corporation and/or their affiliate(s) or their respective owners.

Triton is a registered trademark of Union Carbide Corporation.

www.lifetechnologies.com

3 December 2012

