

IL-12 + p40 ELISA Kit

Catalog Number KAC1561 (96 tests)

Pub. No. MAN0019458 Rev. 1.0

CAUTION! This kit contains materials with small quantities of sodium azide. Sodium azide reacts with lead and copper plumbing to form explosive metal azides. Upon disposal, flush drains with a large volume of water to prevent azide accumulation. Avoid ingestion and contact with eyes, skin and mucous membranes. In case of contact, rinse affected area with plenty of water. Observe all federal, state, and local regulations for disposal.

Note: For safety and biohazard guidelines, see the “Safety” appendix in the ELISA Technical Guide (Pub. no. MAN0006706). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Product description

The Invitrogen™ IL-12 + p40 ELISA Kit is a solid-phase sandwich Enzyme-Linked Immunosorbent Assay (ELISA). This assay is designed to detect and quantify the level of IL-12 + p40 in serum, plasma, and cell culture supernatant. The assay recognizes both natural and recombinant IL-12 + p40.

Contents and storage

Upon receipt, store the kit at 2°C to 8°C. Store the Wash Buffer Concentrate at room temperature. When stored as indicated, all reagents are stable until the expiration date.

Contents	Cat. No. KAC1561 (96 tests)
Anti-IL-12 + p40 Antibody-Coated Wells, 96-well strip-well plate	1 plate
Diluent Spec. (bovine plasma); for cell culture or urine	2 vials
Incubation Buffer; for serum/plasma	11 mL
Standard; 0 ng/mL in buffer with bovine serum	1 vial
Standards 1 to 5; in bovine serum; lyophilized. Refer to vial label for reconstitution volume and range	5 vials
Anti-IL-12-HRP Conjugate; in buffer with proteins	6 mL
Controls 1 and 2; in buffer with human plasma	2 vials
Wash Buffer Concentrate (200X)	10 mL
Chromogenic TMB (tetramethylbenzidine) in DMF	1 mL
Substrate Buffer; H ₂ O ₂ in acetate/citrate buffer	3 vials
Stop Solution (1 N HCl)	6 mL

Materials required but not supplied

- Distilled or deionized water
- Microtiter plate reader with software capable of measurement at 450 nm, 490 nm, and 650 nm (polychromatic reading)
- Plate washer—automated or manual (squirt bottle, manifold dispenser, or equivalent)
- Calibrated adjustable precision pipettes and glass or plastic tubes for diluting solutions; beakers, flask and cylinders for preparation of reagents
- Horizontal microplate shaker capable of 700 rpm ± 100 rpm
- Magnetic stirrer

Prepare 1X Wash Buffer

1. Dilute 2 mL of Wash Solution Concentrate (200X) with 398 mL of deionized or distilled water. Label as 1X Wash Buffer.
2. Use a magnetic stirrer to mix the solution.

Note: Use 1X Wash Buffer on the same day it is prepared. Discard unused 1X Wash Buffer at the end of the day.

Prepare Chromogenic TMB

Pipette 0.2 mL of the Chromogenic TMB into one of the vials of Substrate Buffer.

Note: Solution is stable for only 15 minutes. Avoid exposure to light.

Before you begin

IMPORTANT! Reagents are lot-specific. Do not mix or interchange different reagent lots from various kit lots.

- Review the **Procedural guidelines** and **Plate washing directions** in the *ELISA Technical Guide* available at thermofisher.com.
- Allow reagents to reach room temperature before use. Mix to redissolve any precipitated salts.

Sample preparation guidelines

- Refer to the *ELISA Technical Guide* at thermofisher.com for detailed sample preparation procedures.
- Collect samples in pyrogen/endotoxin-free tubes.
- Freeze samples after collection if samples will not be tested immediately. Avoid multiple freeze-thaw cycles of frozen samples. Thaw completely and mix well (do not vortex) prior to analysis.
- Avoid the use of hemolyzed or lipemic sera. If large amounts of particulate matter are present in the sample, centrifuge or filter sample prior to analysis.
- Avoid impurities contained in sampling materials that can stimulate IL-12 + p40 production by blood cells and produce a false increase in plasma values for IL-12 + p40.

Reconstitute controls

Note: Controls are stable for 4 days at 2–8°C. For longer term storage, make aliquots and store at –20°C for up to 2 months. Avoid successive freeze thaw cycles.

Reconstitute Controls 1 and 2 by adding 1 mL of distilled water to each vial.

If the results obtained for Control 1 and/or Control 2 are not within the range specified on the vial label, the results cannot be used unless a satisfactory explanation for the discrepancy can be determined.

Dilute samples


Sample concentrations should be within the range of the standard curve. Because conditions may vary, each investigator should determine the optimal dilution for each application.





- Dilute serum, plasma, and pleural fluid with Diluent Spec.
- Dilute cell culture supernatant and urine with Incubation Buffer or the cell culture medium used.
- If samples generate values higher than the highest standard, dilute samples further and repeat the assay.

Perform Assay (Total assay time: 4.5 hours)

IMPORTANT! Perform a standard curve with each assay.

- Allow all components to reach room temperature before use. Mix all liquid reagents prior to use.
- Determine the number of 8-well strips required for the assay. Insert the strips in the frames for use. Re-bag any unused strips and frames, and store at 2°C to 8°C for future use.



1	Bind antigen 	<ol style="list-style-type: none">Pipet 100 µL of Incubation Buffer into the appropriate wells for the Standards and Controls.Pipet 100 µL of Incubation Buffer into the appropriate wells for serum/plasma samples or pipet 100 µL of Diluent Spec. into the appropriate wells for cell culture supernatant/urine samples.Add 100 µL of standards, controls, or samples (see “Pre-dilute samples” on page 2) to the appropriate wells. Leave the wells for chromogen blanks empty.Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.Thoroughly aspirate the solution and wash wells 3 times with 1X Wash Buffer.
2	Add HRP Conjugate solution 	<ol style="list-style-type: none">Add 100 µL of Diluent Spec., then 50 µL of anti-IL-12 + p40 conjugate into all the wells.Incubate for 2 hours at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm.Thoroughly aspirate the solution from the wells and wash wells 3 times with 1X Wash Buffer.
3	Add Chromogenic TMB 	<ol style="list-style-type: none">Add 200 µL of Chromogenic TMB to each well. The substrate solution begins to turn blue.Incubate for 30 minutes at room temperature on a horizontal shaker set at 700 rpm ± 100 rpm in the dark. Note: TMB should not touch aluminum foil or other metals.
4	Add stop solution 	Add 50 µL Stop Solution to each well. Tap the side of the plate to mix. The solution in the wells changes from blue to yellow.

Read the plate and generate the standard curve

1. Read the absorbance at 450 nm. Read the plate within 3 hours after adding the Stop Solution.
2. Use curve-fitting software to generate the standard curve. A 4 parameter algorithm provides the best standard curve fit. Optimally, the background absorbance may be subtracted from all data points, including standards, unknowns and controls, prior to plotting.
3. Read the concentrations for unknown samples and controls from the standard curve. Multiply value(s) obtained for sample(s) by the appropriate factor to correct for the sample dilution.

Note: Dilute samples producing signals greater than the upper limit of the standard curve in Standard Diluent Buffer and reanalyze. Multiply the concentration by the appropriate dilution factor.

Performance characteristics

Standard curve example

The following data are for illustration only and should never be used in place of a real-time standard curve.

Concentration (pg/mL)	Optical Density (450 nm)
1,300	3.26
400	1.63
100	0.56
33	0.21
15	0.12
0	0.03

Inter-assay precision

Samples were assayed 11 times in multiple assays to determine precision between assays.

Parameters	Sample 1	Sample 2
Mean (pg/mL)	91	313
Standard Deviation	8	33
% Coefficient of Variation	9	11

Intra-assay precision

Samples of IL-12 + p40 were assayed in replicates of 16 to determine precision within an assay.

Parameters	Sample 1	Sample 2
Mean (ng/mL)	78	288
Standard Deviation	5	16
% Coefficient of Variation	6.5	5.5

Expected values

These values are given only for guidance and it is recommended that each laboratory establishes its own normal values.

Sample	Range (pg/mL)	Average (pg/mL)	Standard deviation
Plasma (n=20)	6–72	27	17
Serum (n=20)	10–89	38	21

Recovery

Sample	Added IL-12 + p40 (pg/mL)	Recovered IL-12 + p40 (pg/mL)	Recovery %
Serum	228	231	82
	99	139	96
	40	87	108
	0	44	—
Plasma	228	236	88
	99	115	80
	40	74	95
	0	36	—
Cell culture medium	265	286	108
	104	115	111
	58	59	102
	0	0	—
Urine	265	259	98
	104	97	93
	0	0	0

Sensitivity

The minimum detectable dose of IL-12 + p40 is 1.5 pg/mL. This was determined by adding two standard deviations to the mean O.D. obtained when the zero standard was assayed 20 times, and calculating the corresponding concentration.

Specificity

No significant cross-reaction was observed in the presence of 150 ng of IL-1a, IL-1β, IL-1ra, IL-2, IL-3, IL-4, IL-5, IL-6, IL-7, IL-8, IL-10, IFN-γ, MCP-3, IL-15, TNF-α, TNF-β, IP-10, SCF, GRO, OSM, LIF, MCP-1, MIP-1a, MIP-1β, RANTES, NAP-2, PDGF, G-CSF, GM-CSF, gp-130 and sIL-6R. The IL-12 + p40 ELISA cross-reacts with Rhesus, Chimpanzee and cynomolgus IL-12. The IL-12 + p40 EASIA has been designed for the measurement of the bioactive p35-p40 heterodimer as well as the antagonist p40-p40 homodimer, allowing the detection of total IL-12 levels. The free subunits p35 and p40 are not detected by the IL-12 + p40 EASIA.

Ratio* = [(sTNF-RI measured in the presence of analyte) x 100] / [sTNF-RII added in the absence of analyte]

No significant cross-reaction was observed in the presence of sTNF-RII or TNF-α. This kit is specific for human natural and recombinant IL-12 + p40.

Linearity of dilution

Dilution	Serum			Plasma		
	Measured conc. (ng/mL)	Theor. conc. (ng/mL)	Recovery (%)	Measured conc. (ng/mL)	Theor. conc. (IU/mL)	Recovery (%)
1/1	206	206	—	210	210	—
1/2	99	103	96	103	105	98
1/4	48	52	92	48	53	91
1/8	24	26	92	26	26	100
1/16	13	13	100	13	13	100

Dilution	Cell Culture Medium			Pleural Fluid		
	Measured conc. (ng/mL)	Theor. conc. (ng/mL)	Recovery (%)	Measured conc. (ng/mL)	Theor. conc. (IU/mL)	Recovery (%)
1/1	954	954	—	600	600	—
1/2	480	477	101	303	300	101
1/4	224	239	94	140	150	93
1/8	105	119	88	73	75	97
1/16	54	59	92	35	38	92

Important Licensing Information


This product may be covered by one or more Limited Use Label Licenses. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

Product label explanation of symbols and warnings

 REF	Catalog Number	 LOT	Batch code		Temperature limitation		Use by		Manufacturer		Consult instructions for use		Caution, consult accompanying documents
---	----------------	--	------------	--	------------------------	--	--------	--	--------------	--	------------------------------	--	---

 Life Technologies Corporation | 7335 Executive Way | Frederick, MD 21704 | USA
For descriptions of symbols on product labels or product documents, go to thermofisher.com/symbols-definition.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

©2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.