USER GUIDE





BioPrime® Total Genomic Labeling System

Catalog Numbers 18097-010, 18097-011, and 18097-012

Document Part Number 25-1006 Publication Number MAN0000596 Revision 3.0



For Research Use Only. Not for use in diagnostic procedures.

Contents

Kit Contents and Storage	
Overview	6
Methods	9
Before Starting	9
Labeling	
Purification	
Assessing the Efficiency of the Labeling Procedure	
Troubleshooting	
Appendix	18
Technical Support	
Additional Products	
References	

Kit Contents and Storage

Kit Sizes and Modules	This manual supports the following BioPrime [®] Total kit configurations:			
	Catalog No.	# of Reactions	Modules	
	18097-010	10	Labeling and Purification	
	18097-011	30	Labeling and Purification	
	18097-012	30	Labeling only	
Shipping and Storage	Ig and The BioPrime [®] Total Labeling Module is shipped on dr and should be stored at -80°C. The 2X Reaction Mixes is be stored at 4°C for up to 4 weeks (store at -80°C long-		dry ice es may ng-term).	
	The BioPrime [®] be stored at roo	Purification Modu m temperature.	le is shipped and s	hould
Labeling Module	The following c 2X Reaction Min (store at –80°C l	omponents should kes may be stored ong-term).	l be stored at -80°C at 4°C for up to 4 v	C. The veeks

Component	10-reaction kit	30-reaction kits
Alexa Fluor [®] 3 2X Reaction Mix	125 µL	$125 \mu L \times 3$
Alexa Fluor [®] 5 2X Reaction Mix	125 µL	$125 \mu L \times 3$
Exo- Klenow Fragment, 40 U/µL	30 µL	30 µL × 3
5 mM EDTA (pH 8.0)	1.2 mL	1.2 mL
TE Buffer (10/1, pH 8.0)	660 μL	660 μL
Control DNA (Salmon Sperm), (10 mg/mL)	10 µL	10 µL

Important

The fluorescently labeled nucleotides in the Alexa Fluor[®] 2X Reaction Mixes are sensitive to photobleaching. Store the mixes protected from light.

Kit Contents and Storage, Continued

Purification Module

The following components should be stored at room temperature.

Component	10-reaction kit	30-reaction kit
PureLink [®] Spin Columns with Collection Tubes	10 columns/ tubes	30 columns/ tubes
Binding Buffer (B2) (combine with 100% isopropanol; see Preparing Binding Buffer B2 with Isopropanol)	9 mL	9 mL
Wash Buffer (W1) (combine with 100% ethanol; see Preparing Wash Buffer W1 with Ethanol)	10 mL	10 mL
Elution Buffer (E1) (10 mM Tris-HCl, pH 8.5)	3 mL	3 mL
Amber Recovery Tubes	10 tubes	30 tubes

Preparing Binding	Binding Buffer B2 supplied with the must be mixed with 100% isopropa	e Purification Module nol prior to use.	
Buffer B2 with Isopropanol Preparing Wash Buffer W1 with Ethanol	Add the amount of isopropanol below directly to the bottle of Binding Buffer B2, and mark the checkbox on the bottle to indicate that you have added the isopropanol.		
	Binding Buffer B2 100% Isopropanol Final Volume Store the buffer with isopropanol a	<u>Amount</u> 9 mL (entire bottle) <u>6 mL</u> 15 mL t room temperature.	
	Wash Buffer W1 supplied with the Purification Module must be mixed with 100% ethanol prior to use.		
	Add the amount of ethanol below directly to the bottle of Wash Buffer W1, and mark the checkbox on the bottle to indicate that you have added the ethanol.		
	Wash Buffer W1 100% Ethanol Final Volume Store the buffer with ethanol at roo	<u>Amount</u> 10 mL (entire bottle) <u>40 mL</u> 50 mL m temperature.	

Overview

Introduction Array comparative genomic hybridization (aCGH) is a microarray-based method for analyzing genomic DNA to detect variations in gene copy number between samples (Pollack et al., 1999; Pollack et al., 2002). In aCGH, two genomic DNA samples are labeled with different fluorophores. The samples are hybridized to a microarray and the ratio of the fluorescent intensities of the fluorophores is measured for each feature on the array (Beheshti et al., 2003; Cai et al., 2002; Snijders et al., 2001). This ratio provides a relative measure of the difference in gene copy number between the samples. The BioPrime® Total Genomic Labeling System uses a mutant form of the Klenow fragment of DNA Polymerase I (Exo- Klenow) and nucleotides labeled with two novel, application-specific dyes (Alexa Fluor® 3 and 5) to differentially label genomic DNA samples for analysis by aCGH. The kit provides the dye-labeled nucleotides in an optimized master-mix formulation that includes unlabeled nucleotides and random primers, for ease of reaction setup. Labeled DNA generated using this system can detect differences in gene copy number from as little as 50 ng of genomic DNA, depending on the sample and array type. Advantages of Amplified products labeled with novel Alexa Fluor®3 the System and 5 dyes have greater yields and higher signal intensities on the array. Exo- Klenow polymerase incorporates fluorescently modified nucleotides more effectively and provides higher yields than standard Klenow, for greater reproducibility of results. PureLink[®] purification columns, included with Cat. no. 18097-010 and 18097-011, are designed to effectively remove all unincorporated nucleotides for the most accurate quantitation of labeled product and reduced background on the array. Dye-specific 2X Reaction Mixes include random primers and both labeled and unlabeled nucleotides for simplified reaction setup and workflow.

• Provides a complete solution for fluorescent labeling of genomic DNA.

Overview, Continued



Overview, Continued

Alexa Fluor [®] 3 and Alexa Fluor [®] 5 Dyes	The novel, applicat Fluor® 5 dyes used used microarray sc correlation (R ²) val Cy®5 dye pair, imp The table lists the e each dye:	ion-specific Alexa Fluor [®] 3 in the kit are compatible v anners, and provide great ues than the spectrally sim roving the resolution of tv xcitation/emission maxim	3 and Alexa vith commonly er signal uilar Cy®3 and vo-color aCGH. na and color of	
	Dye	Excitation/Emission	Color	
	Alexa Fluor [®] 3	555/565 nm	Pink	
	Alexa Fluor® 5	650/670 nm	Light blue	
Control DNA Materials Supplied by the User	Control DNA (Salm help you determine Equations for calcu procedure using th In addition to the k following items on Genomic Labeling	non Sperm DNA) is includ e the efficiency of the label lating the efficiency of the e Control DNA are provid it components, you should hand before using the Bio System.	led in the kit to ing procedure. labeling led on page 15. d have the Prime® Total	
	• 50 ng-3 μg ger	nomic DNA (amount is arr	ay-dependent)	
	• Vortex mixer			
	Microcentrifug	ge		
	 Heat block, air lid 	incubator, or thermocycle	er with a heated	
	• Ice			
	 1.7-mL DNase tubes 	-free capped tubes or thin-	walled PCR	
	Aerosol-resistant pipette tips			
	• 100% isopropa purification bu	nol and 100% ethanol (for Iffers; see page 5)	preparing the	
	• Optional (if qu	enching the labeling react	ion): 0.5 M EDTA	

Methods

Before Starting

Amount of Starting Material	This kit has been designed to generate $\approx 8 \ \mu g$ of labeled DNA from as little as 50 ng to 3 μg of input genomic DNA. In general, larger amounts of starting material result in better fold changes on the array. Use the amount of starting material recommended by your array manufacturer.
Isolating Genomic DNA	Isolate genomic DNA using your method of choice. The PureLink [®] Genomic DNA Purification Kit (Cat. no. K1810-01) is a complete kit for the isolation of genomic DNA. See page 19 for ordering information. A wide range of ChargeSwitch [®] Genomic DNA purification kits is also available from Life Technologies.
DNA Treatment	The genomic DNA may be either intact or treated by enzymatic digestion or sonication, depending on the requirements of your array manufacturer. Note that DNA that has been fragmented by enzymatic digestion or sonication generally yields better results in aCGH.
General Handling of DNA	When handling DNA, use sterile conditions to ensure that no DNases are introduced. All equipment that comes into contact with DNA should be sterile, including pipette tips, microcentrifuge tubes, snap-cap polypropylene tubes, and pipettes. Be sure pipettor barrels are clean and treated with ethanol.
Checking DNA Quantity and Quality	Genomic DNA may be run on an agarose gel to check for quantity and quality. Bufferless E-Gel® Pre-cast Agarose Gels are available from Life Technologies for fast and easy electrophoresis. See page 19 for ordering information.
Storing DNA	After isolating the DNA, we recommend that you proceed directly to Labeling , page 10. Otherwise, store the isolated genomic DNA at 4°C. Note that storage in TE Buffer is recommended for greater stability.

Labeling

Before Starting	In addition to the components of the Labeli following materials are supplied by the use	ng Module, the r:			
	• 50 ng–3 µg genomic DNA (amount is array-dependent)				
	Vortex mixer				
	Microcentrifuge				
	• Heat block, incubator, or thermocycler	with a heated lid			
	• Ice				
	• 1.7-mL capped tubes or thin-walled PCR tubes				
	Optional (if quenching the labeling rea	ction): 0.5 M EDTA			
Preparing DNA Samples in TE	If the genomic DNA is in TE Buffer, add the provided in the kit to a final volume of 22 µ recommend using the TE in the kit (10 mM 1 mM EDTA); not all laboratory TE is provi If the sample is in water, see Preparing DN Water .	e TE Buffer IL. We Tris-HCl pH 8.0, ded at this ratio. A Samples in			
Preparing DNA Samples in Water	If the genomic DNA sample is in water, you of 5 mM EDTA to the sample to ensure opti (The EDTA will be at a concentration of 0.4 50-µL labeling reaction.) Add 5 mM EDTA bring the final volume to 22 µL as follows:	u must add 4.4 μL imal labeling. 4 mM in the final to the sample and			
	DNA sample in water 5 mM EDTA (provided in the kit) Sterile, distilled water	X μL 4.4 μL to 22 μL			
Important	 The addition of 5 mM EDTA to DNA s as described in the previous section, is ensure optimal labeling conditions. Fluorescently labeled nucleotides are s photobleaching. During all steps of the careful to minimize exposure of the 2X and labeled DNA to light. 	amples in water, necessary to ensitive to procedure, be Reaction Mixes			

Labeling, Continued

Preparing the Control DNA	The 10 r To j 1.	Control DNA is provided at a conc ng/mL, and should be diluted in TF prepare 1 μ g of Control DNA for lab Dilute the Control DNA in TE Buff concentration of 1 μ g/ μ L:	entration of E Buffer pric peling: Fer to a final	for to use.
		Control DNA (10 mg/mL) TE Buffer (provided in the kit) Final volume		1 μL <u>9 μL</u> 10 μL
	2.	Add 1 μ L of the diluted Control DI Buffer, for a final volume of 22 μ L.	NA to 21 μL	of TE
Incubation Methods	The incu read	incubation steps may be performed ubator, or thermocycler with a heate ction protected from light.	l in a heat b d lid. Incub	lock, air ate the
Labeling Procedure	1.	If necessary, thaw the Alexa Fluor [®] 2X Reaction Mixes at room temperature, protected from light (mixes may be stored at 4°C for up to 4 weeks).		
	2.	Briefly vortex each 2X Reaction Mit collect the contents of the tubes. Pla	x and centri ace the tube	ifuge to s on ice.
	3.	Genomic DNA samples should be in water with EDTA in a volume of described on page 10. Add the foll DNase-free 1.7-mL capped tubes of tubes:	e prepared i of 22 μL, as owing to se r thin-walle	n TE or parate d PCR
		Component	Tube 1	Tube 2
		Alexa Fluor [®] 3 2X Reaction Mix	25 µL	_
		Alexa Fluor [®] 5 2X Reaction Mix	_	25 µL
		Genomic DNA Sample 1	22 µL	_
		Genomic DNA Sample 2	_	22 µL
		Total volume	47 µL	47 µL
	4.	Gently pipet up and down to mix a protected from light, for 5 minutes on ice for 5 minutes.	and incubat . Immediate	e at 95°C, ely cool
	5.	On ice, add 3 μ L of Exo- Klenow Fr for a final reaction volume of 50 μ L	agment to e	each tube,

Labeling, Continued

Labeling Procedure, continued	6.	Vortex tubes briefly and centrifuge to collect the contents.
	7.	Incubate at 37°C for 2 hours in a heat block, air incubator, or thermocycler with a heated lid, protected from light.
	8.	After incubation, if you are storing the reaction for any length of time prior to purification, add 5 μ L of 0.5 M EDTA to each tube to quench the reaction. If you are proceeding directly to purification, you can skip this step.

Proceed to **Purification**, page 13. The reaction can be stored at –20°C overnight if necessary (following the addition of EDTA).

Purification

BioPrime [®] Purification Module	Cat Moe Foll DN	no. 18097-010 and 18097-011 include a Purification dule developed for use with the BioPrime® Total system. ow the procedure in this section to purify your labeled A using this module.	
Other Methods of Purification	Cat. no. 18097-012 does not include a Purification Module. Use your preferred method of purification, and then proceed to page 15. When assessing the labeling efficiency using a spectrophotometer, be sure to blank the spectrophotometer using the elution buffer from your purification system.		
PureLink [®] PCR Purification	The PureLink [®] PCR Purification System (K3100-01 and K3100-02) has been tested with the BioPrime [™] Total system, and is recommended if you are using cat. no. 18097-012. Ordering information is provided on page 19.		
Purification Procedure	Foll Cat pro 1.	ow the steps below using the Purification Module from no. 18097-010 and 18097-011 to purify the labeled DNA bes. Add 200 µL of Binding Buffer B2 (prepared with isopropanol as described on page 5) to each tube from	
	2.	Step 8, page 12, and vortex to mix. Load each sample onto a PureLink [®] Spin Column,	
	3.	Centrifuge at $10,000 \times g$ for 1 minute. Discard the flow- through and place the column back in the collection tube.	
	4.	Add 650 µL of Wash Buffer W1 (prepared with ethanol as described on page 5) to the column.	
	5.	Centrifuge at $10,000 \times g$ for 1 minute. Discard the flow- through and place the column back in the collection tube.	
	6.	Spin at maximum speed for an additional 2–3 minutes to remove any residual wash buffer. Discard the flow through.	

Purification, Continued

Purification Procedure, continued	7.	Place the Spin Column in a new, sterile Amber Recovery Tube (supplied in the kit).
	8.	Add 55 μ L of Elution Buffer E1 to the center of the column and incubate at room temperature for 1 minute.
	9.	Centrifuge at maximum speed (~20,000 × g) for 2 minutes. The flow-through contains the purified labeled DNA probes. (Discard the column after use.)
	To to A pag	determine the efficiency of the labeling reaction, proceed Assessing the Efficiency of the Labeling Procedure , ge 15.

For a list of array hybridization reagents available from Life Technologies, see page 19.

Assessing the Efficiency of the Labeling Procedure

the Results	To calculate the amount of labeled DNA using a UV/visible spectrophotometer:
	1. Transfer an appropriate volume of purified, labeled DNA from step 9, page 14, to a clean cuvette. Use an appropriate volume for your spectrophotometer. Blank the spectrophotometer using 10 mM Tris-HCl, pH 8.5.
	Important: The labeled DNA must be purified as described on page 14 before scanning, as any unincorporated labeled nucleotides will interfere with the detection of labeled DNA.
	2. Measure the absorbance of the sample at A ₂₆₀ , A ₃₂₀ , A ₅₅₅ , A ₆₅₀ , and A ₇₅₀ . Wash each cuvette thoroughly between samples.
	Yield:1
	DNA (μ g) = (A ₂₆₀ -A ₃₂₀) × 50 μ g/mL × volume in mL
	Dye Incorporation: ²
	Alexa Fluor [®] 3 (pmole) = $(A_{555}-A_{650})/0.15 \times \text{volume in } \mu\text{L}$
	Alexa Fluor [®] 5 (pmole) = $(A_{650}-A_{750})/0.24 \times \text{volume in } \mu\text{L}$
	Degree of Labeling: ³
	Alexa Fluor [®] 3 base/dye ratio = $((A_{260} - A_{320}) - ((A_{555} - A_{650}) \times 0.04)) \times 150,000/(A_{555} - A_{650}) \times 6600$
	Alexa Fluor [®] 5 base/dye ratio =
	$((A_{260} - A_{320}) - ((A_{650} - A_{750}) \times 0)) \times 239,000/(A_{650} - A_{750}) \times 6600$
	Notes:
	¹ Subtracting A_{320} from A_{260} corrects for any silica particles that may leak from the purification columns and artificially increase the yield calculations.
	² Subtracting A ₆₅₀ from A ₅₅₅ and A ₇₅₀ from A ₆₅₀ corrects for any fluorescent background that might artificially increase the measure of dye incorporation.
	³ Absorbance at A_{555} has a very slight effect on the A_{260} reading, and the formula (($A_{555} - A_{650}$) × 0.04)) corrects for this. Conversely, there is no effect of A_{650} on the A_{260} reading; the multiplication by zero was added to the second formula to keep the formulas consistent

Assessing the Efficiency of the Labeling Procedure, Continued

Control DNA	Typically, if starting with 1 µg of Control DNA as specified on page 11, you should expect the following: Yield: ≥7 µg of amplified DNA		
	Dye incorporation:	≥175 pmol Alexa Fluor® 3 ≥300 pmol Alexa Fluor® 5	
	Degree of labeling:	≥0.7 for Alexa Fluor® 3 ≥1.2 for Alexa Fluor® 5	
Note	The Alexa Fluor [®] dye conditions of the Bio have been optimized intensity and signal/ correlate directly wit labeling when compa Labeling with the Bio System yields microa signal/background r labeled with other dy dye incorporation an	e-labeled nucleotides and reaction Prime [®] Total Genomic Labeling System for use on microarrays. Signal 'background on microarrays does not th dye incorporation or degree of aring different fluorescent dyes. oPrime [®] Total Genomic Labeling array signal intensities and ratios greater than or equal to DNA ye-labeled nucleotides, even with lower ud/or degree of labeling.	

Troubleshooting

Problem	Cause	Solution
Yield of labeled DNA from the control reaction is low	DNA has been lost in the purification step after labeling	Make sure that isopropanol has been added to the Binding Buffer and ethanol has been added to the Wash Buffer, as specified on page 5. Measure the amount of labeled DNA in the control reaction before and after purification. Repeat the labeling and purification procedures, following all steps without modifications
	EDTA concentration in the sample is too low	The labeling reaction has been optimized for a final EDTA concentration of 0.44 mM. If your genomic DNA sample is in water, you should add EDTA to the sample before labeling as specified on page 10.
	Starting amount of DNA is too low	Increase the amount of starting DNA.
Cannot detect labeled probes	DNA has been lost in the purification step after labeling	Make sure that isopropanol has been added to the Binding Buffer and ethanol has been added to the Wash Buffer, as specified on page 5. Measure the amount of labeled DNA in the control reaction before and after purification. Repeat the labeling and purification procedures following all steps without
		modifications.
Amount of incorporated labeled nucleotides is low or fluorescence is low	Starting amount of DNA is too low	Increase the amount of starting DNA.
	EDTA concentration in the sample is too low	The labeling reaction has been optimized for a final EDTA concentration of 0.44 mM. If your genomic DNA sample is in water, you should add EDTA to the sample before labeling as specified on page 10.
	Reaction tubes have been exposed to light	Avoid direct exposure of the reaction tubes to light. Repeat the labeling procedure.
	Fluorescent nucleotides have been exposed to light	Repeat the labeling reaction, being careful to avoid direct exposure to light.
	Inefficient labeling due to improper purification	Follow all the purification steps as described in the procedures.

Appendix

Technical Support

Obtaining Support	For the latest services and support information for all locations, go to www.lifetechnologies.com .		
	At the website, you can:		
	• Access worldwide telephone and fax numbers to contact Technical Support and Sales facilities		
	• Search through frequently asked questions (FAQs)		
	 Submit a question directly to Technical Support (techsupport@lifetech.com) 		
	• Search for user documents, SDSs, vector maps and sequences, application notes, formulations, handbooks, certificates of analysis, citations, and other product support documents		
	Obtain information about customer training		
	Download software updates and patches		
Safety Data Sheets (SDS)	Safety Data Sheets (SDSs) are available at www.lifetechnologies.com/support .		
Certificate of Analysis	The Certificate of Analysis provides detailed quality control and product qualification information for each product. Certificates of Analysis are available on our website. Go to www.lifetechnologies.com/support and search for the Certificate of Analysis by product lot number, which is printed on the box.		
Limited Product Warranty	Life Technologies and/or its affiliate(s) warrant their products as set forth in the Life Technologies General Terms and Conditions of Sale found on the Life Technologies web site at www.lifetechnologies.com/termsandconditions . If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support .		

Additional Products

AdditionalLife Technologies has additional reagents that may be used
to prepare labeled probes for hybridization. Ordering
information is provided below.

Product	Quantity	Catalog No.
PureLink [®] PCR Purification System	50 reactions	K3100-01
	250 reactions	K3100-02
E-Gel [®] 1.2% Starter Pak	6 gels and base	G6000-01
Human Cot-1 DNA®-Fluorometric QC	1 mg	15279-101
Human Cot-1 DNA®	500 µg	15279-011
Mouse Cot-1 DNA®	500 µg	18440-016
Yeast tRNA	25 mg	15401-011
	50 mg	15401-029

References

- Beheshti, B., Braude, I., Marrano, P., Thorner, P., Zielenska, M., and Squire, J. A. (2003) Chromosomal localization of DNA amplifications in neuroblastoma tumors using cDNA microarray comparative genomic hybridization. *Neoplasia*, 5, 53-62
- Cai, W. W., Mao, J. H., Chow, C. W., Damani, S., Balmain, A., and Bradley, A.
 (2002) Genome-wide detection of chromosomal imbalances in tumors using BAC microarrays. *Nat Biotechnol*, 20, 393-396
- Pollack, J. R., Perou, C. M., Alizadeh, A. A., Eisen, M. B., Pergamenschikov, A., Williams, C. F., Jeffrey, S. S., Botstein, D., and Brown, P. O. (1999) Genome-wide analysis of DNA copy-number changes using cDNA microarrays. *Nat Genet*, 23, 41-46
- Pollack, J. R., Sorlie, T., Perou, C. M., Rees, C. A., Jeffrey, S. S., Lonning, P. E., Tibshirani, R., Botstein, D., Borresen-Dale, A. L., and Brown, P. O. (2002) Microarray analysis reveals a major direct role of DNA copy number alteration in the transcriptional program of human breast tumors. *Proc. Natl. Acad. Sci. USA*, 99, 12963–12968
- Snijders, A. M., Nowak, N., Segraves, R., Blackwood, S., Brown, N., Conroy, J., Hamilton, G., Hindle, A. K., Huey, B., Kimura, K., Law, S., Myambo, K., Palmer, J., Ylstra, B., Yue, J. P., Gray, J. W., Jain, A. N., Pinkel, D., and Albertson, D. G. (2001) Assembly of microarrays for genome-wide measurement of DNA copy number. *Nat Genet*, 29, 263-264

©2012 Life Technologies Corporation. All rights reserved. The trademarks mentioned herein are the property of Life Technologies Corporation or their respective owners.

Cy is a registered trademark of GE Healthcare UK Limited.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Notes

Notes

Headquarters 5791 Van Allen Way | Carlsbad, CA 92008 USA Phone +1 760 603 7200 | Toll Free in USA 800 955 6288 For support visit

lifetechnologies.com/support or email techsupport@lifetech.com

lifetechnologies.com

1 October 2012

