

# Pre-miR™ hsa-miR-1 miRNA Precursor

Store at or below  $-20^{\circ}\text{C}$ .  
Do not store in a frost-free freezer.

|   |  |
|---|--|
| <b>Catalog # (P/N):</b>                 | AM17150  |
| <b>Product Description:</b>             | A double-stranded RNA oligonucleotide designed for use as a positive control in experiments using Pre-miR™ miRNA Precursors.   |
| <b>Amount:</b>                          | 5 nmol   |
| <b>Appearance:</b>                      | Powder   |
| <b>Additional Material(s) Included:</b> | 1.75 mL Nuclease-free Water  |
| <b>Product Mimics Precursor for:</b>    | <p><u>Sanger miRNA Registry ID#:</u> MI0000651 (human mir-1-1), MI0000437 (human mir-1-2), MI0000139 (mouse mir-1-1), MI0000652 (mouse mir-1-2)</p> <p><u>Description:</u> hsa-miR-1, found in <i>Homo sapiens</i> and <i>Mus musculus</i></p> <p><u>Mature miRNA Sequence:</u> UGGA AUGUAAAGAAGUAUGUA</p>   |
| <b>miRNA Target Information:</b>        | <p><u>Gene Symbol:</u> PTK9</p> <p><u>Full Gene Name:</u> Protein tyrosine kinase 9</p> <p><u>Organism(s):</u> Human and Mouse</p> <p><u>RefSeq Number(s):</u> NM_002822 (human) and NM_008971.3 (mouse)</p> <p><u>Entrez Gene ID(s):</u> 5756 (human) and 19230 (mouse)</p> <p><u>TaqMan® Assay(s):</u> PTK9 TaqMan® Gene Expression Assays Hs00702289_s1(Human) and Mm01598981_g1 (Mouse) recommended (Applied Biosystems; not included)</p> |
| <b>Storage Conditions:</b>              | Store at or below $-20^{\circ}\text{C}$ . <b>Do not store in a frost-free freezer.</b> (Dried oligonucleotides are shipped at ambient temperature.)  |

---

## USER INFORMATION

---

|                               |  |
|-------------------------------|--|
| <b>General Information:</b>   | <p>Ambion® Pre-miR™ miRNA Precursors (P/N AM17100, AM17101, AM17103; patent pending) are designed to mimic endogenous mature microRNAs (miRNAs). Each Pre-miR miRNA Precursor includes one strand that is identical to a known mature miRNA. The design of Pre-miR miRNA Precursors ensures that the portion of the miRNA precursor identical to the mature miRNA is selected for uptake and activation by the miRNA processing pathway.</p> <p>Pre-miR™ hsa-miR-1 miRNA Precursor is an miRNA precursor mimic designed to be processed to mature miR-1 miRNA after delivery to mammalian cells. miR-1 is associated with downregulation of many genes [1]. In particular, the gene for Protein Tyrosine Kinase 9 (PTK9) has two predicted binding sites for miR-1 in its 3' UTR, and miR-1-mediated downregulation of PTK9 gene expression has been shown to occur at the mRNA level [1].</p> <p>miRNAs have been shown to regulate gene expression through three mechanisms: inhibition of translation without concomitant decrease in mRNA levels [2, 3], mRNA degradation [1, 4–6], and transcriptional silencing mediated through changes in chromatin state [7]. miRNA-mediated gene silencing is dependent on complementarity of the target messenger RNAs with critical sequences in the miRNA.</p> <p>Scientists at Applied Biosystems, Austin, have confirmed the finding by Lim et al [1], showing that expression of PTK9 mRNA decreased by 80% following Pre-miR hsa-miR-1 miRNA Precursor transfection, compared to transfection with Pre-miR™ miRNA Precursors—Negative Control #1 (P/N AM17110), using real-time PCR and TaqMan® Gene Expression Assays indicated above. Thus Pre-miR hsa miR 1 miRNA Precursor provides a convenient positive control for miRNA-mediated gene silencing of PTK9, when monitored using the appropriate TaqMan® Gene Expression Assays and real-time PCR.</p> <p>Pre-miR™ miRNA Precursors can be delivered to mammalian cells by chemical transfection or electroporation.</p> |
| <b>Handling Instructions:</b> | <p>RNA oligonucleotides are susceptible to degradation by exogenous ribonucleases introduced during handling. Wear gloves when handling this product. Use RNase-free reagents, tubes, and barrier pipette tips. Upon receipt, store in a non-frost-free freezer at or below <math>-20^{\circ}\text{C}</math> (dried oligonucleotides are shipped at ambient temperature).</p> <p><b>Resuspension Instructions</b><br/>Briefly centrifuge the tube to ensure that the dried oligonucleotide is at the bottom of the tube. Resuspend the oligonucleotide at a convenient concentration. To minimize freeze-thaw cycles, we recommend preparing a concentrated stock, such as 50 <math>\mu\text{M}</math>, and then further diluting to a practical working stock concentration. (Resuspend 5 nmol of oligonucleotide in 100 <math>\mu\text{L}</math> of Nuclease-free Water to obtain a 50 <math>\mu\text{M}</math> solution.)</p> <p>Ambion provides an online calculator for suspension of dry oligonucleotides on its web site at <a href="http://www.ambion.com/techlib/append/oligo_dilution.html">www.ambion.com/techlib/append/oligo_dilution.html</a></p>  |

Once reconstituted in Nuclease-free Water, the oligonucleotide is ready to transfect and can be used at your choice of final concentration.

Store the resuspended oligonucleotide at or below  $-20^{\circ}\text{C}$ . **Do not store in a frost-free freezer.**

## Applications:

Pre-miR™ hsa-miR-1 miRNA Precursor is designed for use as a positive control in miRNA experiments to confirm that the transfection procedure and cell cultures support miRNA-mediated gene silencing. It can also be used in development and optimization of miRNA transfection conditions. Expression of miR-1 is typically restricted to muscle cells and, thus, is not expressed in most other cell lines. Analysis of miR-1 overexpression is straightforward, since expression of miR-1 causes a profound decrease in PTK9 mRNA.

Optimization of transfection conditions is critical to ensure that oligonucleotides are delivered inside the cells and are capable of having biological activity once inside. Transfection optimization involves identifying conditions for maximal activity of the transfected oligonucleotides with minimal amounts of cellular toxicity. The delivery procedure itself can induce changes in gene expression profiles, making a nontransfected control meaningless for evaluation of silencing in delivery optimization. We recommend transfection with nontargeting Pre-miR™ miRNA Precursors–Negative Controls (P/N AM17110, AM17111) as a baseline reference. Nontargeting Pre-miR miRNA Precursors–Negative Controls should be used at the same concentration as positive control and experimental Pre-miR miRNA Precursors, because nucleic acid concentrations within cells can affect the activity and specificity of miRNAs.

To calculate the amount of reduction in target gene expression during transfection optimization or to analyze knockdown of a target gene by real-time PCR, Applied Biosystems Austin scientists typically use the  $\Delta\Delta\text{C}$ , method [8], normalizing to 18S rRNA, and comparing to Pre-miR miRNA Precursor–Negative Control-transfected samples.

### Transfection Starting Points for Mammalian Cells

As with other small nucleic acids, such as siRNAs and antisense oligonucleotides, the efficiency with which mammalian cells are transfected with Pre-miR™ miRNA Precursors will vary according to cell type and the transfection agent used. The optimal concentration used for transfections should be determined empirically. We have found that Pre-miR miRNA Precursors typically work best when transfected at a final concentration of 3–30 nM. However, a more extensive concentration range from 1–100 nM can be analyzed in optimization experiments. The following chart provides general starting points for transfection of Pre-miR miRNA Precursors into cultured mammalian cells.

### General Transfection Starting Points for Pre-miR miRNA Precursors in Cultured Mammalian Cells

| Plate Format                         | 96 wells              | 24 wells          | 12 wells          | 6 wells            |
|--------------------------------------|-----------------------|-------------------|-------------------|--------------------|
| Transfection Agent <sup>a</sup>      | 0.3–1.0 $\mu\text{L}$ | 1–3 $\mu\text{L}$ | 2–4 $\mu\text{L}$ | 3–6 $\mu\text{L}$  |
| Pre-miR miRNA Precursor <sup>b</sup> | 3 pmol                | 15 pmol           | 30 pmol           | 75 pmol            |
| Cell Density <sup>c</sup>            | 6,000 cells/well      | 40,000 cells/well | 80,000 cells/well | 200,000 cells/well |
| Final Volume per Well                | 0.1 mL                | 0.5 mL            | 1.0 mL            | 2.5 mL             |

<sup>a</sup> Refer to the instructions provided with your transfection agent for the recommended volume.

<sup>b</sup> The amount shown results in a final Pre-miR miRNA Precursor concentration of 30 nM. The amount of Pre-miR miRNA Precursor required for maximal Pre-miR miRNA Precursor activity will vary among cell types. For a 96-well plate and 100  $\mu\text{L}$  final transfection volume, 3 pmol of a 5  $\mu\text{M}$  oligonucleotide solution is 0.6  $\mu\text{L}$ . Robotic pipettors may require volumes of 2–5  $\mu\text{L}$  for accurate pipetting. To increase pipetting volumes and accuracy when preparing transfection complexes, we recommend first making a plate with a dilution of your stock oligonucleotide.

<sup>c</sup> Optimal cell density will vary among cell types, depending on cell size and growth characteristics. In general, we recommend 30–70% confluency.

### Transfection Optimization

Optimizing transfection efficiency is crucial for maximizing Pre-miR miRNA Precursor activity while minimizing cytotoxicity. Optimal transfection efficiencies are achieved by identifying an effective transfection agent for each cell type and by adjusting (in order of importance):

- Amount of transfection agent
- Amount and type of RNA oligonucleotide
- Cell density at the time of transfection
- Order of transfection (pre-plating cells or plating cells/transfecting in tandem)
- Length of exposure of cells to transfection agent/Pre-miR miRNA Precursor complexes

Most protocols recommend maintaining mammalian cells in the medium used for transfection; this avoids dilution or removal of RNA oligonucleotides from the cells by adding medium or washing the cells with new medium too soon after transfection. We have found that cells typically exhibit greater viability when existing medium is replaced with fresh medium 24 hours after transfection. Replacing medium after 24 hours generally does not change the activity of the transfected Pre-miR miRNA Precursor.

Once the conditions for optimal transfection efficiency are determined, they should be kept constant from experiment to experiment for a given cell type.

For additional information about small RNA transfection, including transfection conditions for many cell types and optimization protocols, see the Ambion siRNA Delivery Resource at:  
[www.ambion.com/techlib/resources/delivery](http://www.ambion.com/techlib/resources/delivery)

**Additional Information:**

For protocols, background information, a reference list, and miRNA research tools, see the Ambion miRNA information resource and product guide:  
[www.ambion.com/techlib/resources/miRNA/index.html](http://www.ambion.com/techlib/resources/miRNA/index.html).

**References:**

1. Lim LP, Lau NC, Garrett-Engle P, Grimson A, Schelter JM, Castle J, Bartel DP, Linsley PS, Johnson JM. (2005) Microarray analysis shows that some microRNAs downregulate large numbers of target mRNAs. *Nature* **433(7027)**:769–773. Epub 2005 Jan 30.
2. Lee RC, Feinbaum RL, Ambros V (1993) The *C. elegans* heterochronic gene *lin-4* encodes small RNAs with antisense complementarity to *lin-14*. *Cell* **75(5)**:843–854.
3. Wightman B, Ha I, Ruvkun G (1993) Posttranscriptional regulation of the heterochronic gene *lin-14* by *lin-4* mediates temporal pattern formation in *C. elegans*. *Cell* **75(5)**:855–862.
4. Yekta S, Shih IH, Bartel DP. 2004. MicroRNA-directed cleavage of HOXB8 mRNA. *Science* **304(5670)**:594–596.
5. Bagga S, Bracht J, Hunter S, Massirer K, Holtz J, Eachus R, Pasquinelli AE (2005) Regulation by *let-7* and *lin-4* miRNAs results in target mRNA degradation. *Cell* **122(4)**:553–563.
6. Jing Q, Huang S, Guth S, Zarubin T, Motoyama A, Chen J, Di Padova F, Lin SC, Gram H, Han J (2005) Involvement of microRNA in AU-rich element-mediated mRNA instability. *Cell* **120(5)**:623–634.
7. Schramke V, Sheedy DM, Denli AM, Bonila C, Ekwall K, Hannon GJ, Allshire RC. 2005. RNA-interference-directed chromatin modification coupled to RNA polymerase II transcription. *Nature* **435(7046)**:1275–1279.
8. Bustin SA (2002) Quantification of mRNA using real-time reverse transcription PCR (RT-PCR): trends and problems. *J Mol Endocrinol.* **29(1)**:23–39.

---

**RELATED PRODUCTS**

---

**Pre-miR™ miRNA Precursors (patent pending)**

P/N AM17100, AM17101, AM17103

Chemically modified and optimized nucleic acids designed to mimic the microRNA (miRNA) molecules in cells.

**Pre-miR™ miRNA Precursors—Negative Controls**

P/N AM17110, AM17111

Designed to serve as a negative control for experiments involving Pre-miR™ miRNA Precursors.

**Pre-miR™ miRNA Precursor Molecule Library**

P/N 4385830

An extensive collection of Pre-miR miRNA Precursors mimicking precursors to human miRNAs in miRBase Sequence Database Version 9.2.

**Anti-miR™ miRNA Inhibitors and Controls**

P/N AM17000, AM17001, AM17003, AM17010, AM17120, AM17121

Chemically modified and optimized nucleic acids designed to specifically inhibit microRNA (miRNA) molecules in cells.

**siPORT™ NeoFX™ Transfection Agent**

P/N AM4510 and AM4511

A versatile lipid-based agent for efficient and reproducible transfection of adherent cells while subculturing, without increased cytotoxicity.

**siPORT™ Amine Transfection Agent**

P/N AM4502 and AM4503

An easy-to-use blend of polyamines that delivers siRNA into mammalian cells with minimal cytotoxicity.

**TaqMan® Gene Expression Assays**

[www.allgenes.com](http://www.allgenes.com)

A comprehensive collection of over 700,000 probe and primer sets for quantitative gene expression analysis using real-time PCR.

---

**QUALITY CONTROL**

---

**Identity:**

The mass of a sample of each single-stranded RNA oligonucleotide is analyzed using MALDI-TOF mass spectrometry and compared to the calculated mass.

**Purity:**

Analytical HPLC of a sample of each purified single-stranded RNA oligonucleotide is used to confirm ≥90% purity.

**Annealing:**

The annealed Pre-miR miRNA Precursor is analyzed by nondenaturing gel electrophoresis.

---

**OTHER INFORMATION**

---

**Material Safety Data Sheets:**

Material Safety Data Sheets (MSDSs) can be printed or downloaded from product-specific links on our website at the following address: [www.ambion.com/techlib/msds](http://www.ambion.com/techlib/msds). Alternatively, e-mail your request to [MSDS\\_Inquiry\\_CCRM@appliedbiosystems.com](mailto:MSDS_Inquiry_CCRM@appliedbiosystems.com). Specify the catalog or part number(s) of the product(s), and we will e-mail the associated MSDSs unless you specify a preference for fax delivery. For customers without access to the internet or fax, our technical service department can fulfill MSDS requests placed by telephone or postal mail. (Requests for postal delivery require 1–2 weeks for processing.)

**Warranty and Liability:**

*For research use only. Not for use in diagnostic procedures.*

Applied Biosystems is committed to delivering superior product quality and performance, supported by industry-leading global service and technical support teams. Warranty information for the accompanying consumable product is available at [www.ambion.com/info/warranty](http://www.ambion.com/info/warranty) in "Limited Warranty for Consumables," which is subject to the exclusions, conditions, exceptions, and limitations set forth under the caption "EXCLUSIONS, CONDITIONS, EXCEPTIONS, AND LIMITATIONS" in the full warranty statement. Please contact Applied Biosystems if you have any questions about our warranties or would like information about post-warranty support.

Information in this document is subject to change without notice. Applied Biosystems assumes no responsibility for any errors that may appear in this document.

Applied Biosystems disclaims all warranties with respect to this document, expressed or implied, including but not limited to those of merchantability or fitness for a particular purpose. In no event shall Applied Biosystems be liable, whether in contract, tort, warranty, or under any statute or on any other basis for special, incidental, indirect, punitive, multiple or consequential damages in connection with or arising from this document, including but not limited to the use thereof.

**Trademarks, Patents, and Licensing:**

Applied Biosystems, AB (Design), and Ambion are registered trademarks, and Anti-miR, *NeoFX*, Pre-miR and siPORT are trademarks of Applied Biosystems or its subsidiaries in the US and/or certain other countries. TaqMan is a registered trademark of Roche Molecular Systems, Inc. All other trademarks are the sole property of their respective owners.

A license to perform the patented 5' Nuclease Process for research is obtained by the purchase of (i) both Licensed Probe and Authorized 5' Nuclease Core Kit, (ii) a Licensed 5' Nuclease Kit, or (iii) license rights from Applied Biosystems.

The TaqMan<sup>®</sup> Gene Expression Assay contains Licensed Probe. Use of the TaqMan Gene Expression Assays is covered by US patent claims and patent claims outside the US. The purchase of TaqMan Gene Expression Assays includes a limited, non-transferable immunity from suit under patent claims for using only this amount of product for the purchaser's own internal research. The right to use this product in the 5' Nuclease Process under the applicable claims of US patents, and patent claims outside the United States, can be obtained through purchase of an Authorized 5' Nuclease Core Kit. Except under separate license rights available from Applied Biosystems, no right under any other patent claim, or to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, or to sublicense, repackaging with other products, or resell in any form, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained from the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

© 2008, 2010 Ambion, Inc. All rights reserved. 4382576D