

HisPur™ Ni-NTA Purification Kit

88227 88228 88229

2206.1

Number	Description
88227	HisPur Ni-NTA Purification Kit, 0.2mL Kit Contents: HisPur Ni-NTA Spin Columns , 0.2mL resin bed, 25 each Phosphate-Buffered Saline (10X) , 30mL, 200mM sodium phosphate, 3M sodium chloride, pH 7.4 Imidazole (2M) , 8mL, pH 7.4 Collection Tubes , 80 each Column Plugs , 25 each
88228	HisPur Ni-NTA Purification Kit, 1mL Kit Contents: HisPur Ni-NTA Spin Columns , 1mL resin bed, 5 each Phosphate-Buffered Saline (10X) , 30mL, 200mM sodium phosphate, 3M sodium chloride, pH 7.4 Imidazole (2M) , 8mL, pH 7.4 Column Plugs , 5 each
88229	HisPur Ni-NTA Purification Kit, 3mL Kit Contents: HisPur Ni-NTA Spin Columns , 3mL resin bed, 5 each Phosphate-Buffered Saline (10X) , 30mL, 200mM sodium phosphate, 3M sodium chloride, pH 7.4 Imidazole (2M) , 25mL, pH 7.4 Column Plugs , 5 each

Binding Capacity: ≤ 60mg of a 28kDa 6xHis-tagged protein from a bacterial source per milliliter of settled resin

Resin: Crosslinked 6% agarose

Supplied: 50% slurry in 20% ethanol

Storage: Upon receipt store at 4°C. Product shipped at ambient temperature.

Introduction

The Thermo Scientific HisPur Ni-NTA Purification Kit enables effective immobilized metal affinity chromatography (IMAC) purification of polyhistidine-tagged proteins from a soluble protein extract. The Ni-NTA resin is composed of nickel-charged nitrilotriacetic acid (NTA) chelate immobilized onto 6% crosslinked agarose. The Ni-NTA resin is compatible with native or denaturing conditions and can be used in multiple formats, including conventional gravity-flow chromatography, spin column and FPLC. Ni-NTA resins are commonly chosen for His-tagged-protein purification because of the four metal-binding sites on the chelate, which allow for high-binding capacity and low metal ion leaching.

Important Product Information

- Protein yield and purity are dependent upon the expression level, conformation and solubility characteristics of the recombinant fusion protein. Therefore, it is important to optimize these parameters before attempting a large-scale purification. For best results, perform a small-scale test to estimate the expression level and determine the solubility of each His-tagged protein.
- Optimization of the lysis procedure is critical for maximizing protein yield. Some methods for protein extraction include using commercially available detergent-based reagents, such as Thermo Scientific B-PER Bacterial Protein Extraction Reagent with Enzymes (Product No. 90078), and mechanical methods, such as freeze/thaw cycles, sonication or French press. Add EDTA-free protease inhibitors, such as Thermo Scientific Halt Protease and Phosphatase Inhibitor Cocktail, EDTA-free (Product No. 78441), to protect proteins from degradation.
- Avoid using protease inhibitors or other additives that contain chelators, such as EDTA, or strong reducing agents, such as DTT or β -mercaptoethanol, which will disrupt the function of the nickel resin.
- When using the Thermo Scientific Coomassie Plus (Bradford) Assay (Product No. 23238) or Thermo Scientific Pierce 660 nm Protein Assay (Product No. 22660) to monitor protein concentration in the elution fractions, dilute the samples at least 1:2 before performing the protein assay.

Additional Materials Required

- Regeneration MES Buffer: 20mM 2-(*N*-morpholine)-ethanesulfonic acid, 0.1M sodium chloride; pH 5.0

Material Preparation

Depending on the specific protein, buffers might require some optimization. Use the table below to make buffers with different imidazole concentrations. Adjust the total volume depending on the resin volume used. For most proteins, the following imidazole concentrations are recommended:

- Equilibration Buffer: 10mM imidazole
- Wash Buffer: 25mM imidazole
- Elution Buffer: 250mM imidazole

<u>Imidazole Final</u> <u>Conc. (mM)</u>	<u>10X PBS*</u> <u>(mL)</u>	<u>2M Imidazole</u> <u>(μL)</u>	<u>Water</u> <u>(mL)</u>
10	1	50	8.95
25	1	125	8.875
40	1	200	8.8
60	1	300	8.7
75	1	375	8.625
150	1	750	8.25
200	1	1,000	8
250	1	1,250	7.75
500	1	2,500	6.5

*Phosphate-buffered saline.

Procedure for Spin Purification of His-Tagged Proteins

Note: The total volume of the 0.2, 1 and 3mL columns are 1, 8 and 22mL, respectively. If a sample volume is greater than the column, perform multiple applications and centrifugations until the entire sample has been processed. Be careful not to exceed the resin's binding capacity. The HisPur Ni-NTA Spin Columns also may be used for gravity-flow purifications.

1. Equilibrate column(s) to working temperature. Perform purifications at room temperature or at 4°C.
2. Prepare sample by mixing protein extract with Equilibration Buffer so the total volume equals two resin-bed volumes.
3. Remove the bottom tab from the HisPur Ni-NTA Spin Column by gently twisting. Place column into a centrifuge tube.

Note: Use 1.5, 15 or 50mL centrifuge tubes for the 0.2, 1 and 3mL spin columns, respectively.

4. Centrifuge column at $700 \times g$ for 2 minutes to remove storage buffer.

5. Equilibrate column with two resin-bed volumes of Equilibration Buffer. Allow buffer to enter the resin bed.
6. Centrifuge column at $700 \times g$ for 2 minutes to remove buffer.
7. Add the prepared protein extract to the column and allow it to enter the resin bed.
Note: For maximal binding, the sample can be incubated for 30 minutes at room temperature or 4°C on an end-over-end rocking platform.
8. Centrifuge column at $700 \times g$ for 2 minutes and collect the flow-through in a centrifuge tube.
9. Wash resin with two resin-bed volumes of Wash Buffer. Centrifuge at $700 \times g$ for 2 minutes and collect fraction in a centrifuge tube. Repeat this step two more times collecting each fraction in a separate centrifuge tube.
Note: If desired, perform additional washes. Monitor washes by measuring their absorbance at 280nm.
10. Elute His-tagged proteins from the resin by adding one resin-bed volume of Elution Buffer. Centrifuge at $700 \times g$ for 2 minutes. Repeat this step two more times, collecting each fraction in a separate tube.
Note: If performing gravity-flow add two resin-bed volumes of Elution Buffer to achieve proper flow characteristics.
11. Monitor protein elution by measuring the absorbance of the fractions at 280nm or by Coomassie Plus (Bradford) Assay Reagent (Product No. 23238). The eluted protein can be directly analyzed by SDS-PAGE.
Note: To remove imidazole for downstream applications, use gel filtration (e.g., Thermo Scientific Zeba Spin Desalting Columns) or dialysis (e.g., Thermo Scientific Slide-A-Lyzer Dialysis Cassettes).

Procedure for Ni-NTA Resin Regeneration

The Ni-NTA resin may be used at least five times without affecting protein yield or purity. Between each use, perform the procedure as described below to remove residual imidazole and any nonspecifically adsorbed protein. To prevent cross-contamination of samples, designate a given column to one specific fusion protein.

1. Wash resin with 10 resin-bed volumes of MES Buffer.
2. Wash resin with 10 resin-bed volumes of ultrapure water.
3. Store resin as a 50% slurry in 20% ethanol.

Troubleshooting

Problem	Possible Cause	Solution
Low protein yield	Poor expression of soluble protein	Optimize expression conditions
	His-tagged protein forms inclusion bodies	Alter growth conditions to minimize inclusion body formation and maximize soluble protein yield
		Solubilize inclusion bodies and perform the purification with a compatible denaturant (e.g., Thermo Scientific Inclusion Body Solubilization Reagent, Product No. 78115)
	Insufficient cell lysis and extraction	Optimize cell lysis protocol
	His tag is absent	Verify the sequence or perform an ELISA or Western blot using an antibody against the His tag
	His tag is inaccessible using native conditions	See the Additional Information Section for denaturing conditions
	His-tagged protein has a low affinity to the column	Optimize the Equilibration or Wash buffer by decreasing the concentration of imidazole (see Table in the Material Preparation Section)
Poor protein purity	Insufficient column washing	Wash column additional times or modify the imidazole concentration (see Table in Material Preparation Section)
Slow column flow	Column is overloaded	Apply less protein extract to the column and make sure the extract is not too viscous or contaminated with cell debris

Additional Information

A. Fusion Proteins Expressed in Inclusion Bodies

Over-expressed proteins are sometimes sequestered in inclusion bodies. Inclusion bodies can be solubilized in 8M urea, 6M guanidine or the Inclusion Body Solubilization Reagent (Product No. 78115); however, a denaturant must be added to the buffers to ensure the protein remains soluble throughout the procedure. (Follow the spin purification procedure.)

If using 8 M urea, you may proceed directly to SDS-PAGE without sample clean-up. If using 6 M guanidine, perform one of the following steps before SDS-PAGE: 1) Dilute sample 1:6 in ultrapure water; 2) Perform buffer exchange with dialysis or desalting; 3) Perform TCA precipitation; 4) Use the Thermo Scientific Pierce SDS-PAGE Sample Prep Kit (Product No. 89888).

For denaturing conditions prepare the following buffers:

- Equilibration Buffer: PBS with 6M guanidine•HCl, 10mM imidazole; pH 7.4
- Wash Buffer: PBS with 6M guanidine•HCl, 25mM imidazole; pH 7.4
- Elution Buffer: PBS with 6M guanidine•HCl, 250mM imidazole; pH 7.4

B. Visit the website for additional information relating to this product including the following:

- Tech Tip # 43: Protein Stability and Storage
- Tech Tip # 40: Convert Between Times Gravity ($\times g$) and Centrifuge Rotor Speed (RPM)
- Tech Tip # 6: Extinction Coefficients Guide

Related Thermo Scientific Products

88270	Pierce [®] High Capacity Endotoxin Removal Resin, 10mL
88282	Pierce LAL Chromogenic Endotoxin Quantitation Kit
88221	HisPur Ni-NTA Resin, 10mL
89967	HisPur Cobalt Spin Columns, 0.2mL, 25 each
89968	HisPur Cobalt Spin Columns, 1mL, 5 each
89969	HisPur Cobalt Spin Columns, 3mL, 5 each
24110	Guanidine•HCl, 500g
90078	B-PER [®] Bacterial Protein Extraction Reagent with Enzymes, 250mL
87785	Halt Protease Inhibitor Cocktail (100X), EDTA-free, 1mL
78441	Halt Protease and Phosphatase Inhibitor Cocktail, EDTA-free (100X), 1mL
23238	Coomassie Plus (Bradford) Assay Reagent, 300mL
22660	Pierce 660 nm Protein Assay Reagent, 750mL
78115	Inclusion Body Solubilization Reagent, 100mL
89890	Zeba Spin Desalting Columns, 7K MWCO, 2mL, 25 columns, for 200-700 μ L samples
89892	Zeba Spin Desalting Columns, 7K MWCO, 5mL, 25 columns, for 500-2,000 μ L samples
89894	Zeba Spin Desalting Columns, 7K MWCO, 10mL, 25 columns, for 1,500-4,000 μ L samples
87769	Zeba Spin Desalting Columns, 40K MWCO, 2mL, 25 columns, for 200-900 μ L samples
87771	Zeba Spin Desalting Columns, 40K MWCO, 5mL, 25 columns, for 300-2,000 μ L samples
87773	Zeba Spin Desalting Columns, 40K MWCO, 10mL, 25 columns, for 1,000-4,000 μ L samples
87730	Slide-A-Lyzer G2 Dialysis Cassettes, 10K MWCO, 3mL, 10 cassettes
87731	Slide-A-Lyzer G2 Dialysis Cassettes, 10K MWCO, 15mL, 8 cassettes

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.