




Bismaleimide Crosslinkers (BMOE, BMB and BMH)

22323 22331 22330

0649.3

Number	Description
22323	BMOE , bis(maleimido)ethane, 50mg Molecular Weight: 220.18 Spacer Arm: 8.0Å Net Mass Added: 220.05
22331	BMB , 1,4-bis(maleimido)butane, 50mg Molecular Weight: 248.23 Spacer Arm: 10.9Å Net Mass Added: 248.08
22330	BMH , bis(maleimido)hexane, 50mg Molecular Weight: 276.29 Spacer Arm: 13.0Å Net Mass Added: 276.11 CAS #: 4856-78-5



Storage: Upon receipt store desiccated at 4°C. Product is shipped at ambient temperature.

Introduction

Thermo Scientific BMOE, BMB and BMH are homobifunctional, maleimide crosslinkers for conjugation between sulfhydryl groups (-SH). Such bismaleimide crosslinkers are commonly used to explore and characterize protein structure (i.e., oligomerization) or protein interactions. Because BMOE, BMB and BMH have the same reactivity but differ in length, the relative success of these three reagents in forming crosslinks between sites in a protein oligomer or interaction can assist in determining intra- and intermolecular distances.

Reaction of a sulfhydryl to the maleimide group results in formation of a stable thioether linkage (Figure 1), which cannot be cleaved by reducing agents or physiological buffer conditions. Reaction of maleimides is very specific to sulfhydryls at pH 6.5-7.5.¹ Although maleimides will react to primary amines at pH > 8, the rate is 1000 times slower than the reaction to sulfhydryls at pH 7. Unlike iodoacetamides, maleimides do not react with tyrosines, histidines or methionines.

The maleimide moiety is temporarily stable in aqueous solutions devoid of reactive sulfhydryl targets, but hydrolysis to a nonreactive maleamic acid can occur during storage, especially at pH > 8 (Figure 1). For this reason, dissolved reagents are best used promptly and the remainder discarded. Hydrolysis of the ring structure also can occur following conjugation, resulting in an open-ring linkage² (Figure 1).




Figure 1. Reaction of maleimide-activated compounds to sulfhydryls.

Important Product Information

- Molecules to be reacted with maleimide compounds must have free (reduced) sulfhydryls. Reduce peptide disulfide bonds with Thermo Scientific Immobilized TCEP Disulfide Reducing Gel (Product No. 77712). Reduce disulfide bonds in high molecular weight proteins using 5mM TCEP (1:100 dilution of Thermo Scientific Bond-Breaker TCEP Solution, Product No. 77720) for 30 minutes at room temperature, followed by TCEP removal using a desalting column (e.g., Thermo Scientific Zeba Spin Desalting Columns). Proteins (e.g., antibodies) can be inactivated by complete reduction of their disulfide bonds. Selective reduction of hinge-region disulfide bonds in IgG can be accomplished with 2-Mercaptoethylamine•HCl (2-MEA, Product No. 20408). Sulfhydryls can be added to molecules using *N*-succinimidyl *S*-acetylthioacetate (SATA, Product No. 26102) or 2-iminothiolane•HCl (Traut's Reagent, Product No. 26101), which modify primary amines.
- Avoid extraneous sulfhydryl-containing components in the reaction buffers during conjugation (e.g., DTT), as they react with the maleimide portion of the reagent, inhibiting and reducing conjugation efficiency of the intended target.
- The maleimide group reacts predominantly with free sulfhydryls at pH 6.5-7.5, forming stable thioether bonds. At pH values > 7.5, reactivity toward primary amines and hydrolysis of the maleimide groups can occur. At pH 7, the maleimide group is ~1000 times more reactive toward a free sulfhydryl than to an amine.

Procedure for Crosslinking Proteins in Solution

Generally, a two- or three-fold molar excess of crosslinker over the amount of sulfhydryl-containing protein(s) results in sufficient conjugation between proximal sulfhydryl groups. Empirical testing of reagent and protein concentrations is necessary to determine optimal conditions for the experiment.

A. Material Preparation

- Conjugation Buffer: Phosphate buffered saline (PBS, pH 7.2; e.g., Product No. 28372) or other sulfhydryl-free buffer at pH 6.5-7.5. Include 5-10mM EDTA to help prevent the reoxidation of disulfides by trace divalent metals.
- Crosslinker Stock Solution: Immediately before use, weigh a small quantity of crosslinker and dissolve it in dimethylsulfoxide (DMSO) or dimethylformamide (DMF) at a 5-20mM concentration. For example, make a 20mM solution of each reagent as follows:
 - For BMOE, dissolve 2.2mg reagent in 0.5mL DMSO or DMF
 - For BMB, dissolve 2.5mg reagent in 0.5mL DMSO or DMF
 - For BMH, dissolve 2.75mg reagent in 0.5mL DMSO or DMF
- Sulfhydryl-containing protein, prepared as described the Important Product Information section.
- (Optional): Quenching Solution: concentrated (0.5-1M) cysteine, DTT, or other thiol-containing reducing agent.
- (Optional): Desalting column (e.g., Zeba™ Spin Desalting Columns) or dialysis unit (e.g., Thermo Scientific Slide-A-Lyzer Dialysis Cassettes) to separate crosslinked proteins from excess nonreacted crosslinker.

B. Procedure for Protein Crosslinking

1. Dissolve protein(s) in Conjugation Buffer at 0.1mM (e.g., 5mg in 1mL for a 50kDa protein).
2. Add crosslinker to the dissolved protein(s) at 0.2mM final concentration (= two-fold molar excess for 0.1mM protein solution) by adding 10µL of Crosslinker Stock Solution per milliliter of protein solution.

Note: The reaction solution may appear cloudy as a result of the low aqueous solubility of the crosslinker; usually, such solutions become clearer as the reaction proceeds. However, initial solubility can be increased by gentle heating and sonication. Other concentrations of Crosslinker Stock Solution can be used, as well as other final molar fold excesses of crosslinker. Most proteins remain soluble when the DMSO concentration does not exceed 10-15% of the final reaction volume; if protein solubility is not an issue, there is no limit to the DMSO concentration that may be used.

3. Incubate reaction mixture for 1 hour at room temperature or for 2 hours at 4°C.
4. Quench reaction by adding Quenching Solution at 10-50mM final and incubating for 15 minutes at room temperature. Alternatively (or in addition) remove the excess nonreacted reagent by desalting or dialysis.

Related Thermo Scientific Products

Table 1. Bismaleimide Crosslinkers.

Crosslinker Name	Spacer Arm Length (Å)	Spacer Arm Composition (between maleimide groups)	Product No.
BMOE	8.0	Alkane	22323
BMDB	10.2	Cis-diol (periodate cleavable)	22332
BMB	10.9	Alkane	22331
BMH	13.0	Alkane	22330
DTME	13.3	Disulfide (reducing agent cleavable)	22335
BM(PEO) ₂	14.7	Polyethylene glycol (PEG)	22336
BM(PEO) ₃	17.8	Polyethylene glycol (PEG)	22337

Cited References

1. Smyth, D.G., Blumenfeld, O.O. and Konigsberg, W. (1964). Reaction of *N*-ethylmaleimide with peptides and amino acids. *Biochem J* **91**:589.
2. Partis, M.D., *et al.* (1983). Crosslinking of protein by ω -maleimido alkanoyl *N*-hydroxysuccinimido esters. *J Protein Chem* **2**(3):263-77.

Product References

- Alisio, A. & Mueckler, M. (2004). Relative Proximity and Orientation of Helices 4 and 8 of the GLUT1 Glucose Transporter. *J Biol Chem* **279**:26540-45.
- Chen, L.L., *et al.* (1991). Production of Multimeric Forms of CD4 through a Sugar-based Crosslinking Strategy. *J Biol Chem* **226**(7):18237-43.
- Giron-Monzon, L., *et al.* (2004). Mapping Protein-Protein Interactions between MutL and MutH by Crosslinking. *J Biol Chem* **279**:49338-45.
- Green, N.S., Reisler, E. and Houk, K.N. (2001). Quantitative evaluation of the lengths of homobifunctional protein crosslinking reagents used as molecular rulers. *Protein Sci* **10**:1293-304.
- Ishmael, F.T., Trakselis, M.A. and Benkovic, S.J. (2003). Protein-Protein Interactions in the Bacteriophage T4 Replisome. The leading strand holoenzyme is physically linked to the lagging strand holoenzyme and the primosome. *J Biol Chem* **278**:3145-52.
- Kovalenko, O.V., *et al.* (2004). Evidence for specific tetraspanin homodimers: inhibition of palmitoylation makes cysteine residues available for crosslinking. *Biochem J* **344**:407-17.
- Kwon, Y-K., *et al.* (2003). Role of the pleckstrin homology domain of PLC γ 1 in its interaction with the insulin receptor. *J Cell Biol* **163**:375-84.
- Stalteri, M.A. and Mather, S.J. (1995). A crosslinked monoclonal antibody fragment for improved tumor targeting. *Bioconjugate Chem* **6**:179-86.
- Takebe, K., *et al.* (2003). Epimorphin acts to induce hair follicle anagen in C57BL/6 mice. *FASEB Journal* **17**:2037-47.
- Xie, Z., Turk, E. and Wright, E.M. (2000). Characterization of the *Vibrio parahaemolyticus* Na⁺/Glucose Cotransporter. A bacterial member of the sodium/glucose transporter (SGLT) family. *J Biol Chem* **275**:25959-64.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.