

resDNASEQ™ Quantitative DNA Kits

USER GUIDE

for use with:

(Genomic DNA quantitation)

- resDNASEQ™ Quantitative CHO DNA Kit
- resDNASEQ™ Quantitative *E. coli* DNA Kit
- resDNASEQ™ Quantitative HEK293 DNA Kit
- resDNASEQ™ Quantitative Human DNA Kit
- resDNASEQ™ Quantitative MDCK DNA Kit
- resDNASEQ™ Quantitative NS0 DNA Kit
- resDNASEQ™ Quantitative *Pichia* DNA Kit
- resDNASEQ™ Quantitative Sf9 and Baculovirus DNA Kit
- resDNASEQ™ Quantitative Vero DNA Kit
- resDNASEQ™ Quantitative Synthetic Vero DNA Kit

(Plasmid DNA quantitation)

- resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit

Catalog Numbers 4402085, 4458435, A46014, A26366, 4464335, 4458441, 4464336, A46066, A41797, A53242, and A50337

Publication Number 4469836

Revision K



Life Technologies Ltd | 7 Kingsland Grange | Woolston, Warrington WA1 4SR | United Kingdom

For descriptions of symbols on product labels or product documents, go to [thermofisher.com/symbols-definition](https://www.thermofisher.com/symbols-definition).

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4469836

Revision	Date	Description
K	8 November 2021	Update to include the resDNASEQ™ Quantitative Synthetic Vero DNA Kit (Cat. No. A53242).
J	28 April 2021	Update to the control serial dilutions required for the resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit (Cat. No. A50337).

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

TRADEMARKS: All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

©2021 Thermo Fisher Scientific Inc. All rights reserved.

Contents

■ CHAPTER 1	Product information	5
	Product description	5
	Contents	6
	Required materials not supplied	12
	Workflow	13
■ CHAPTER 2	Genomic DNA quantitation	14
	Prepare the control DNA serial dilutions for the standard curve	14
	Guidelines for standard dilutions	14
	Prepare the control serial dilutions (Genomic DNA)	14
	Prepare the PCR reaction mix	16
	Prepare the PCR plate	17
	Setup, run, and analyze samples with AccuSEQ™ Software on the QuantStudio™	
	5 Real-Time PCR Instrument	19
	Create a resDNASEQ™ experiment	19
	Start the run	22
	Analyze the results	23
■ CHAPTER 3	Plasmid DNA quantitation	25
	Prepare the control DNA serial dilutions for the standard curve	25
	Guidelines for standard dilutions	25
	Prepare the control serial dilutions (Plasmid DNA - Kanamycin Resistance)	25
	Prepare the samples (Kanamycin resistance and synthetic Vero assays only)	27
	Prepare the PCR reaction mix	28
	Prepare the PCR plate	29
	Setup, run, and analyze samples with AccuSEQ™ Software on the QuantStudio™	
	5 Real-Time PCR Instrument	30
	Create a resDNASEQ™ template (Plasmid DNA)	30
	Create a resDNASEQ™ experiment (Plasmid DNA)	31
	Start the run	35
	Analyze the results	36

- **APPENDIX A** Troubleshooting 38
- **APPENDIX B** Use the kit with the 7500 Fast Real-Time PCR Instrument and AccuSEQ™ software v2.x 40
 - Required materials not supplied 40
 - Create a plate document in the AccuSEQ™ software 41
 - Run the plate 42
 - Analyze the results 43
- **APPENDIX C** Use the kit with 7500 System SDS Software v1.5.1 44
 - Create the plate document, run the plate, and analyze the results with 7500 Fast SDS software 44
 - Create a plate document 44
 - Run the plate 48
 - Analyze the results 48
- **APPENDIX D** Good laboratory practices 49
 - Work area setup and lab design 49
 - Good laboratory practices for PCR and RT-PCR 49
 - Avoiding false positives due to cross-contamination 49
- **APPENDIX E** Safety 50
 - Chemical safety 51
 - Biological hazard safety 52
- Documentation and support** 53
 - Related documentation 53
 - Customer and technical support 53
 - Limited product warranty 54

References



Product information

IMPORTANT! Before using this product, read and understand the information in the “Safety” appendix in this document.

Product description

The resDNASEQ™ Quantitative DNA Kits are used to quantitate residual DNA from CHO, *E. coli*, HEK293, Human, MDCK, NS0, *Pichia*, Sf9 and Baculovirus, and Vero, or plasmid DNA for Kanamycin resistance, in cell lines which are used for production of biopharmaceutical products. Use the kit after you extract host-cell DNA from test samples. For extraction information, see the *PrepSEQ™ Sample Preparation Kits User Guide* (Pub. No. 4469838).

The resDNASEQ™ Quantitative DNA Kits use TaqMan™ quantitative PCR to perform rapid, specific quantitation of sub-picogram levels of residual host-cell or plasmid DNA. The assay is accurate and reliable across a broad range of sample types, from in-process samples to final product.

To generate the standard curve used to quantitate the DNA in test samples, the CHO, MDCK, NS0, Vero, synthetic Vero, and Plasmid DNA - Kanamycin resistance assays require six dilutions and the *E. coli*, HEK293, Human, *Pichia*, and Sf9 and Baculovirus assays require five dilutions. Control DNA for standard curve generation is included in the kits. In addition, the kits use an internal positive control (IPC) to evaluate the performance of each PCR reaction.

Contents

Table 1 resDNASEQ™ Quantitative CHO DNA Kit (Cat. No. [4402085](#))

Contents	Amount	Storage
resDNASEQ™ CHO Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	–25°C to –15°C before first use, protect from light 2–8°C after first use, protect from light
10X CHO DNA Real-Time PCR Assay Mix	300 µL	–25°C to –15°C, protect from light
Negative Control (Water)	1.0 mL	–25°C to –15°C before first use 2–8°C after first use
resDNASEQ™ CHO DNA Control		
CHO DNA Control, 30 ng/µL	40 µL	–25°C to –15°C
DNA Dilution Buffer	7 mL	2–8°C

Table 2 resDNASEQ™ Quantitative *E. coli* DNA Kit (Cat. No. [4458435](#))

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	–25°C to –15°C before first use, protect from light 2–8°C after first use, protect from light
Negative Control (Water)	1.0 mL	–25°C to –15°C before first use 2–8°C after first use
10X <i>E. coli</i> DNA Assay Mix	300 µL	–25°C to –15°C, protect from light
resDNASEQ™ <i>E. coli</i> DNA Control		
<i>E. coli</i> DNA Control, 30 ng/µL	40 µL	–25°C to –15°C
DNA Dilution Buffer	7 mL	2–8°C

Table 3 resDNASEQ™ Quantitative HEK293 DNA Kit (Cat. No. [A46014](#))

Contents	Amount	Storage
resDNASEQ™ HEK293 DNA Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X HEK293 Assay Mix	300 µL	-25°C to -15°C, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
resDNASEQ™ HEK293 DNA Control		
HEK293 DNA Control, 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 4 resDNASEQ™ Quantitative Human DNA Kit (Cat. No. [A26366](#))

Contents	Amount	Storage
resDNASEQ™ Human Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X Human DNA Assay Mix	300 µL	-25°C to -15°C, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
resDNASEQ™ Human DNA Control		
Human DNA Control, 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 5 resDNASEQ™ Quantitative MDCK DNA Kit (Cat. No. 4464335)

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X MDCK DNA Assay Mix	300 µL	-25°C to -15°C, protect from light
resDNASEQ™ MDCK DNA Control		
MDCK DNA Control, 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 6 resDNASEQ™ Quantitative NS0 DNA Kit (Cat. No. 4458441)

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X NS0 DNA Assay Mix	300 µL	-25°C to -15°C, protect from light
resDNASEQ™ NS0 DNA Control		
NS0 DNA Control, 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 7 resDNASEQ™ Quantitative *Pichia* DNA Kit (Cat. No. 4464336)

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	–25°C to –15°C before first use, protect from light 2–8°C after first use, protect from light
Negative Control (Water)	1.0 mL	–25°C to –15°C before first use 2–8°C after first use
10X <i>Pichia</i> DNA Assay Mix	300 µL	–25°C to –15°C, protect from light
resDNASEQ™ <i>Pichia</i> DNA Control		
<i>Pichia</i> DNA Control, 30 ng/µL	40 µL	–25°C to –15°C
DNA Dilution Buffer	7 mL	2–8°C

Table 8 resDNASEQ™ Quantitative Sf9 and Baculovirus DNA Kit (Cat. No. A46066)

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	–25°C to –15°C before first use, protect from light 2–8°C after first use, protect from light
Negative Control (Water)	1.0 mL	–25°C to –15°C before first use 2–8°C after first use
10X Sf9 + Baculovirus DNA Assay Mix	300 µL	–25°C to –15°C, protect from light
resDNASEQ™ Sf9 and Baculovirus DNA Control		
Multiplex DNA Control with: <ul style="list-style-type: none"> • Sf9 DNA Control, 30 ng/µL • Baculovirus DNA Control, 30 ng/µL 	40 µL	–25°C to –15°C
DNA Dilution Buffer	7 mL	2–8°C

Table 9 resDNASEQ™ Quantitative Vero DNA Kit (Cat. No. [A41797](#))

Contents	Amount	Storage
resDNASEQ™ Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X Vero DNA Assay Mix	300 µL	-25°C to -15°C, protect from light
resDNASEQ™ Vero DNA Control		
Vero DNA Control, 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 10 resDNASEQ™ Quantitative Synthetic Vero DNA Kit (Cat. No. [A53242](#))

Contents	Amount	Storage
resDNASEQ™ Vero PCR Reagents		
TaqMan™ Environmental Master Mix 2.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
10X Vero DNA Assay Mix	300 µL	-25°C to -15°C, protect from light
resDNASEQ™ Synthetic Vero DNA Control		
Vero PLASMID DNA Control, equivalent to 30 ng/µL	40 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Table 11 resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit (Cat. No. A50337)

Contents	Amount	Storage
resDNASEQ™ Kanamycin DNA Real-Time PCR Reagents		
TaqMan™ Environmental Master Mix 3.0	2 × 0.75 mL	-25°C to -15°C before first use, protect from light 2-8°C after first use, protect from light
10X KanR Assay Mix	300 µL	-25°C to -15°C, protect from light
Yeast tRNA (10 mg/mL)	500 µL	-25°C to -15°C
Negative Control (Water)	1.0 mL	-25°C to -15°C before first use 2-8°C after first use
resDNASEQ™ Kanamycin DNA Control		
KanR DNA Control, 3 × 10 ⁷ copies/µL	44 µL	-25°C to -15°C
DNA Dilution Buffer	7 mL	2-8°C

Required materials not supplied

Unless otherwise indicated, all materials are available through thermofisher.com. "MLS" indicates that the material is available from fisherscientific.com or another major laboratory supplier.

Item	Source
Instrument	
QuantStudio™ 5 Real-Time PCR System with AccuSEQ™ Real-Time PCR Software v3.1 or later	Contact your local sales representative.
Consumables	
MicroAmp™ Optical 96-Well Reaction Plate with Barcode, 20 plates, 0.1-mL wells; for use with Applied Biosystems™ 7500 Fast Real-Time PCR System	4346906
MicroAmp™ 96-Well Base	N8010531 , 10 bases
MicroAmp™ Optical Adhesive Film	4311971 , 100 covers 4360954 , 25 covers
MicroAmp™ Adhesive Film Applicator	4333183 , 5 applicators
Miscellaneous items	
Disposable gloves	Major lab supplier (MLS)
Pipettes	MLS
Aerosol-resistant micropipette tips	MLS
For the PCR plate: Fisher Scientific™ Mini Plate Spinner Centrifuge, 120- or 230-volt	14-100-143 (120-volt), 14-100-141 (230-volt)
Nonstick, RNase-free Microfuge Tubes, 1.5 mL (1 box; 250 tubes/box)	AM12450

Workflow

Genomic DNA

Prepare the control serial dilutions (Genomic DNA)
(page 14)



Prepare the PCR reaction mix
(page 16)



Prepare the PCR plate
(page 17)



AccuSEQ™ Software v3.1 for the QuantStudio™ 5 Instrument

Create a resDNASEQ™ experiment
(page 19)



Start the run
(page 22)



Analyze the results
(page 36)

Plasmid DNA

Prepare the control serial dilutions (Plasmid DNA - Kanamycin Resistance)
(page 25)



Prepare the PCR reaction mix
(page 28)



Prepare the PCR plate
(page 29)



AccuSEQ™ Software v3.1 for the QuantStudio™ 5 Instrument

Create a resDNASEQ™ experiment
(page 19)



Start the run
(page 22)



Analyze the results
(page 36)



Genomic DNA quantitation

Prepare the control DNA serial dilutions for the standard curve

Guidelines for standard dilutions

- Prepare the standard curve and the test samples in different areas of the lab.
- Use different sets of pipettors for test sample preparation and for standard curve preparation and aliquoting to avoid cross-contamination of test samples.
- Vortex each tube to mix the contents thoroughly before each dilution step.
- Briefly centrifuge to collect all the liquid at the bottom before making the next dilution.

Prepare the control serial dilutions (Genomic DNA)

1. Label nonstick 1.5-mL microfuge tubes: **NTC, SD1, SD2, SD3, SD4, SD5**.
For CHO, Vero, synthetic Vero, MDCK, and NS0 kits, label an additional tube with **SD6**.
2. Add 50 μL of DNA Dilution Buffer to tube NTC. Put aside.
3. Add 990 μL of DNA Dilution Buffer to tube SD1.
4. Add 450 μL of DNA Dilution Buffer to tubes SD2, SD3, SD4, SD5, and (for CHO, Vero, synthetic Vero, MDCK, and NS0 only) SD6.
5. Remove the tube of DNA control (30 $\text{ng}/\mu\text{L}$) from the freezer.
6. After the DNA thaws, vortex it gently for 2 seconds, then briefly centrifuge.
7. Perform the serial dilutions:

Note: The Vero PLASMID DNA Control, supplied in the resDNASEQ™ Quantitative Synthetic Vero DNA Kit, should be vortexed for 15–30 seconds after each dilution.

- a. Add 10 μL of the DNA control to the tube that is labeled SD1, then vortex thoroughly and briefly centrifuge.
- b. Transfer 50 μL of the DNA from tube SD1 to tube SD2, then vortex thoroughly and briefly centrifuge.

- c. Continue to transfer 50 μ L of DNA from the previous dilution tube to the next dilution tube until you add DNA to tube SD5 (*E. coli*, HEK293, Human, *Pichia*, and Sf9 and Baculovirus) or SD6 (CHO, Vero, synthetic Vero, MDCK, and NS0). Final dilutions are shown in the table. After each transfer, vortex thoroughly, then centrifuge briefly.

Serial dilution (SD) tube	Dilution	pg DNA/reaction (10 μ L of the diluted DNA used in final 30 μ L of PCR reaction)
Control	DNA control tube	300,000 pg
SD 1	10 μ L DNA control + 990 μ L DDB	3,000 pg
SD 2	50 μ L SD1 + 450 μ L DDB	300 pg
SD 3	50 μ L SD2 + 450 μ L DDB	30 pg
SD 4	50 μ L SD3 + 450 μ L DDB	3 pg
SD 5	50 μ L SD4 + 450 μ L DDB	0.3 pg
SD 6 (for CHO, Vero, synthetic Vero, MDCK, and NS0 only)	50 μ L SD5 + 450 μ L DDB	0.03 pg

8. Store the DNA dilution tubes:

Temperature	For use
4°C	Same day
-20°C	\leq 1 week
-20°C	SD1 in single-use aliquots \leq 6 months

Prepare the PCR reaction mix

1. Determine the number of reactions needed for the controls and test samples that you will quantify.
2. Thaw all kit reagents completely at room temperature, thoroughly mix reagents, and briefly centrifuge.
3. Prepare a PCR reaction mix using the reagents and volumes shown in the table below.
 - Multiply the PCR reaction volume for one reaction (30 μ L) by the number of reactions that you need to run.
 - Use 10% excess volume to compensate for pipetting losses.

Note: Use reagents from the same lot for all reactions.

Kit reagents	Volume for 1 30- μ L reaction	Volume for 36 30- μ L reactions (includes 10% overage)
Negative Control (Water)	2 μ L	79.2 μ L
10X DNA assay mix appropriate for the cell line being tested	3 μ L	118.8 μ L
TaqMan™ Environmental Master Mix 2.0	15 μ L	594 μ L
DNA template	10 μ L	Add DNA template to each well separately, not as part of Master Mix
Total	30 μL	792 μL of Master Mix

Prepare the PCR plate

Plate setup differs slightly for each AccuSEQ™ System. See your software user guide for specific instructions. Place samples, NTCs, and standards in different quadrants of the plate.




Figure 1 Default plate setup in the AccuSEQ™ Real-Time PCR Software v3.1

- ① No template controls
- ② Samples
- ③ Standard curve

1. Add 20 μ L PCR reaction mix to each well.
2. Add 10 μ L of PCR NTC to the appropriate wells.
3. Add 10 μ L each of extracted sample DNA to the appropriate wells.

Note: If you prepared samples with the automated protocol, use a multichannel pipette to transfer the extracted sample.

4. Add 10 μ L of standard dilutions to the appropriate wells.


Note: Use different sets of pipettors to dispense test sample and standard curve dilutions to avoid cross-contamination of test samples.

5. Use a film applicator to seal the plate with optical film, then briefly centrifuge with a miniplate centrifuge that is compatible with 96-well plates.

Setup, run, and analyze samples with AccuSEQ™ Software on the QuantStudio™ 5 Real-Time PCR Instrument

Create a resDNASEQ™ experiment

1. In the **Home** screen, click the **Factory default/Admin Defined Template** tab, then select a **resDNASEQ** template.



- ① **Factory default/Admin Defined Template** tab
- ② **resDNASEQ** template (ResDNA_5Std or ResDNA_6Std)

Serial Dilutions (Standards)	Template	Assays
5	_5Std	<i>E. coli</i> , HEK293, Human, <i>Pichia</i> , Sf9 and Baculovirus
6	_6Std	CHO, Vero, synthetic Vero, MDCK, NS0

2. In the **Experiment Properties** pane of the **Setup** tab:
 - a. (Optional) Change the system-generated name of the experiment.
 - b. (Optional) Enter the plate **Barcode**, then add **Comments**.
Default resDNASEQ™ settings (cannot be changed)
 - **Experiment Type** is **Quantitation-Standard Curve**
 - **Chemistry** is **TaqMan™ Reagents**
 - **Ramp Speed** is **Standard - 2hrs**
 - c. Click **Next**.

Note: Experiment names cannot be changed after this step.

- In the **qPCR Method** pane of the **Setup** tab, view the default volume and cycling conditions (cannot be changed).




Figure 2 resDNASEQ™ template default cycling conditions

- Click **Next**.
- In the **Samples** pane of the **Setup** tab, enter the sample **Name** and **Dilution**. **Sample Volume** is not applicable. Add additional **Samples** if needed.

IMPORTANT! Do not change the **Targets**.

① — Samples (24) Add

Color	Name	Dilution	Sample Volume
Orange	1A	1.00	0
Green	1B	1.00	0.00

- Samples** pane
- Add**—adds additional samples

6. In the **Samples** pane of the **Setup** tab, scroll to the right, then enter the spike information.

For more information on plate setup, see the *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

- **Sample Volume**—not applicable; leave as default (0).
- **Spike Volume**—volume of DNA added to the PCR (set to 10).
- **Spike Standard Concentration**—expected spike amount per reaction (for example, 10 pg/reaction or 30 copies/reaction).
- **Reference**—the non-spiked sample; the mean quantity of reference is subtracted during % recovery calculation.
- **Spike Input**—automatically calculated (double check if correct).

Note: If incorrect, be sure **Spike Volume** is set to 10 and **Spike Standard Concentration** is the expected pg spike per PCR reaction.

- (Optional) **Comments**
- **Protein Concentration**—Sample protein concentration (if Total DNA in pg DNA/mg Protein is required).

Samples (24) Add

Spike Volume	Spike Standard Concentration
0	0.00
0.00	0.00
0.00	0.00
0.00	0.00

- ① Textbox—type in the value, or use the up and down arrows
- ② Scroll bar—scroll to find the spike parameter

7. Click **Next**.

The **Run** tab is displayed.

8. Experiments are auto-saved in the software. To save, exit the experiment. The software prompts you to save changes. Click **Yes**.

Note: Clicking **Save As** will create a copy of the experiment.

Start the run

Start the run in the AccuSEQ™ Software.

Option	Description
If the experiment is open	Click Start Run .
If the experiment is closed	<ol style="list-style-type: none"> 1. Open the experiment. 2. Click the Run tab. 3. Click Start Run.

The screenshot shows the AccuSEQ software interface. At the top, there is a navigation bar with tabs: Setup, Run (highlighted with a circled '1'), Result, Report, and Audit & E-Sign. Below the navigation bar, there are several utility buttons: Analysis Settings, Save As, Save As Template, and Export Experiment. The main area displays the 'Run Control' section, which includes a 'Plots' tab and a 'Pool Run Summary' tab. Under 'Plots', there is a 'Select Plot' dropdown menu with 'Amplification' selected. To the left of the 'Amplification Plot' is a graph showing fluorescence intensity versus cycle number, with a legend for 'IPC' (orange) and 'E. Coli' (red). To the right of the graph is a 'Well Plate' view showing a 96-well plate layout with columns 1-12 and rows A-G. The plate contains various samples, including NTC (No Template Control) and TS (Target Sequence) samples. At the bottom right of the interface, there is a 'Start Run' button, which is circled with a '2'.


- 1 Run tab
- 2 Start Run button

A message stating **Run has been started successfully** is displayed when the run has started.

Analyze the results

After the qPCR run is finished, use the following general procedure to analyze the results. For more detailed instructions see the *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

1. In the AccuSEQ™ Real-Time PCR Software, open your experiment, then navigate to the **Result** tab.



- ① **Result** tab
- ② **Analysis Settings**
- ③ Plot horizontal scrollbar
- ④ **Analyse** button

2. In the **Result Analysis** tab, review the Amplification Curve plots for amplification profiles in the controls, samples, and the standard curve.
3. In the **Result Analysis** tab, review the **QC Summary** for any flags in wells.

4. In the **Result Analysis** tab, review the **Standard Curve** plot. Verify the values for the Slope, Y-intercept, R^2 , and Efficiency.

Note: The **Standard Curve** efficiency should be between 90-110% and the $R^2 > 0.99$. If these criteria are not met, up to two points, not in the same triplicate, can be removed from the standard curve data, and the analysis repeated.

5. (Optional) Navigate to the **Report** tab to generate a report of the experiment, or to export results.

3

Plasmid DNA quantitation

Sample residual plasmid DNA is measured in copy number, not concentration. To manually convert the copy number output to a mass measurement, multiply the copy number given by the AccuSEQ™ Real-Time PCR Software to the average molecular weight of your plasmid.

Prepare the control DNA serial dilutions for the standard curve

Guidelines for standard dilutions

- Prepare the standard curve and the test samples in different areas of the lab.
- Use different sets of pipettors for test sample preparation and for standard curve preparation and aliquoting to avoid cross-contamination of test samples.
- Vortex each tube for 15-30 seconds to mix the contents thoroughly before each dilution step.
- Vortex for 15-30 seconds, then tap down standards, before adding the standards to the PCR plate (plasmid standards are more fragile than genomic DNA standards).

Prepare the control serial dilutions (Plasmid DNA - Kanamycin Resistance)

1. Label nonstick 1.5-mL microfuge tubes: **NTC**, **SD1**, **SD2**, **SD3**, **SD4**, **SD5**, and **SD6**.
The dilution **SD6** will not be used for the standard curve. It will be used to confirm the limit of detection (LOD).
2. Add 70 µL of DNA Dilution Buffer to tube NTC. Put aside.
3. Add 990 µL of DNA Dilution Buffer to tube Dilution 1.
4. Add 180 µL of DNA Dilution Buffer to tubes SD1, SD2, SD3, SD4, and SD5.
5. Add 100 µL of DNA Dilution Buffer to tube SD6.
6. Remove the tube of KanR DNA control (3.0×10^7 copies/µL) from the freezer.
7. After the DNA thaws, vortex it thoroughly for 15-30 seconds, then tap gently to bring contents to the bottom of the tube. Do not centrifuge.
8. Perform the serial dilutions:
 - a. Add 10 µL of the KanR DNA control (3.0×10^7 copies/µL) to the tube that is labeled Dilution 1, vortex thoroughly for 15-30 seconds, then tap gently to bring contents to the bottom of the tube.

- b. Transfer 20 μL of the DNA from tube Dilution 1 to tube SD1, vortex thoroughly for 15-30 seconds, then tap gently to bring contents to the bottom of the tube..
- c. Continue to transfer 20 μL of DNA from the previous dilution tube to the next dilution tube until you add DNA to tube SD5.
- d. Transfer 100 μL of DNA from SD5 to SD6, then vortex thoroughly. Final dilutions are shown in the table. After each transfer, vortex thoroughly for 15-30 seconds, then tap gently to bring contents to the bottom of the tube. Do not centrifuge.

Serial dilution (SD) tube	Dilution	Concentration (copy number/ μL)	Copy number/ PCR reaction
Control	DNA control tube	3.0×10^7	N/A
Dilution 1	10 μL DNA control + 990 μL DDB	300,000	N/A
SD 1	20 μL Dilution 1 + 180 μL DDB	30,000	300,000
SD 2	20 μL SD1 + 180 μL DDB	3,000	30,000
SD 3	20 μL SD2 + 180 μL DDB	300	3,000
SD 4	20 μL SD3 + 180 μL DDB	30	300
SD 5	20 μL SD4 + 180 μL DDB	3	30
SD 6 (LOD)	100 μL SD5 + 100 μL DDB	1.5	15

9. Store the DNA dilution tubes:

Temperature	For use
4°C	≤ 2 days
-20°C	Dilution 1 ≤ 1 week
-20°C	Dilution 1 in single-use aliquots ≤ 6 months

Prepare the samples (Kanamycin resistance and synthetic Vero assays only)

Extract samples with the PrepSEQ™ Residual DNA Sample Preparation Kit (Cat. No. 4413686).

(Recommended) Use **Yeast tRNA**, supplied in the PrepSEQ™ Residual DNA Sample Preparation Kit or the resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit for the kanamycin resistance and synthetic Vero assays.

1. Dilute the **Yeast tRNA**.

Table 12 Diluted Yeast tRNA

Component	Volume
Yeast tRNA (10mg/mL)	5 µL
PBS (1X), pH 7.2	245 µL
Total	250 µL

2. Add 5 µL **Diluted Yeast tRNA** to 370 µL of each test sample. This is sufficient for triplicate 100 µL extractions.

Prepare the PCR reaction mix

1. Determine the number of reactions needed for the controls and test samples that you will quantify.
2. Thaw all kit reagents completely at room temperature, thoroughly mix reagents, and briefly centrifuge.
3. Prepare a PCR reaction mix using the reagents and volumes shown in the table below.
 - Multiply the PCR reaction volume for one reaction (30 μ L) by the number of reactions that you need to run.
 - Use 10% excess volume to compensate for pipetting losses.

Note: Use reagents from the same lot for all reactions.

Kit reagents	Volume for 1 30- μ L reaction	Volume for 48 30- μ L reactions (includes 10% overage)
Negative Control (Water)	2 μ L	106 μ L
10X DNA assay mix appropriate for the cell line being tested	3 μ L	159 μ L
TaqMan™ Environmental Master Mix 2.0	15 μ L	795 μ L
DNA template	10 μ L	Add DNA template to each well separately, not as part of Master Mix
Total	30 μL	1060 μL of Master Mix

Prepare the PCR plate

Plate setup differs slightly for each AccuSEQ™ System. See your software user guide for specific instructions. Place samples, NTCs, and standards in different quadrants of the plate.




Figure 3 Default plate setup in the AccuSEQ™ Real-Time PCR Software v3.1

- ① No template controls
- ② Samples
- ③ Standard curve

1. Add 20 μ L PCR reaction mix to each well.
2. Add 10 μ L of PCR NTC to the appropriate wells.
3. Add 10 μ L each of sample DNA (with tRNA added) to the appropriate wells.

IMPORTANT! Vortex samples for 15-30 seconds before adding the samples to the **PCR Plate**, then tap to bring contents to the bottom of the wells.

If you prepared samples with the automated protocol, use a multichannel pipette to transfer the extracted sample.

4. Add 10 μ L of standard dilutions to the appropriate wells.

IMPORTANT! Vortex samples for 15-30 seconds before adding the samples to the **PCR Plate**, then tap to bring contents to the bottom of the wells.

If you prepared samples with the automated protocol, use a multichannel pipette to transfer the extracted sample.


5. Use a film applicator to seal the plate with optical film, then briefly centrifuge with a miniplate centrifuge that is compatible with 96-well plates.

Setup, run, and analyze samples with AccuSEQ™ Software on the QuantStudio™ 5 Real-Time PCR Instrument

Create a resDNASEQ™ template (Plasmid DNA)

Plasmid DNA resDNASEQ™ assays do not use a **Factory default/Admin Defined Template**. Create a new template before the first use of these assays.


1. In the **Home** screen, click **Templates** in the left navigation pane.



Template Name	Template Location	Publish	Created Date	Created By User	Modified Date	Modified By User	Actions
2018-11-14_ProA 2020-06-03-...	My Template		03 Jun 2020 20:01:23 GMT-0500	Administrator	03 Jun 2020 20:01:46 GMT-0500	Administrator	+ ✎ 📄 🔒
Mycoseq	Factory Default	☑	01 Jun 2020 06:44:57 GMT-0500	system	01 Jun 2020 06:44:57 GMT-0500	system	+ ✎ 📄 🔒
ResDNA_5Std	Factory Default	☑	01 Jun 2020 06:44:57 GMT-0500	system	01 Jun 2020 06:44:57 GMT-0500	system	+ ✎ 📄 🔒
ProteinSeq HCP	Factory Default	☑	01 Jun 2020 06:44:57 GMT-0500	system	01 Jun 2020 06:44:57 GMT-0500	system	+ ✎ 📄 🔒
ProteinSeq Pro A	Factory Default	☑	01 Jun 2020 06:44:57 GMT-0500	system	01 Jun 2020 06:44:57 GMT-0500	system	+ ✎ 📄 🔒
ResDNA_6Std	Factory Default	☑	01 Jun 2020 06:44:57 GMT-0500	system	01 Jun 2020 06:44:57 GMT-0500	system	+ ✎ 📄 🔒

1. **Templates** icon

2. Click **+ Create New** next to the **ResDNA_5Std** factory default template.
3. Click **Next** to move to the **qPCR Method** screen.
4. Click **Next** to move to the **Plate Setup** screen.
5. In the **Plate Setup** screen, add the **Targets** and **Reporters**.



<input type="checkbox"/>	Color	Name	Reporter	Quencher	Task	Qua
<input type="checkbox"/>	■	Kan	FAM	NFQ-MGB		
<input type="checkbox"/>	■	IPC	NED	NFQ-MGB		

For the resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit, this is the FAM™ dye for the kanamycin (Kan) target and the NED™ dye for the IPC.

6. Click **Save as Template**.
7. Enter a **Template Name** and description, then select **Admin Defined** and **Locked**. Click **Save**.
The template is saved, and can be accessed from **Templates** in the **Home** screen.
8. Click **Templates** in the **Home** screen, then open the new template.

Note: The template must be saved prior to editing the **Analysis Settings**




9. Click  **Analysis Settings**, then deselect **Default Settings**.
10. Enter new thresholds for the targets.


Table 13 resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit

Target	Threshold
Kanamycin	0.04
IPC	0.02

11. Click **Apply**, then close the template.
12. In the  **Home** screen, click  **Templates** in the left navigation pane, then **Publish** the template. See *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287). The template is listed in the **Factory default/Admin Defined Templates**.

Create a resDNASEQ™ experiment (Plasmid DNA)

Plasmid DNA resDNASEQ™ assays do not use a **Factory default/Admin Defined Template**. Create a new template before the first use of these assays.

1. In the  **Home** screen, click **Factory default/Admin Defined Templates**, then select the custom plasmid DNA **resDNASEQ** template created in “Create a resDNASEQ™ template (Plasmid DNA)” on page 30.

Note: To create an experiment from an existing resDNASEQ™ experiment, see *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

2. In the **Experiment Properties** pane of the **Setup** tab:
 - a. (Optional) Change the system-generated name of the experiment.

Note: Names must be unique. Deleted experiment names can not be reused.

- b. (Optional) Enter the plate **Barcode**, then add **Comments**.

Note: Comments are not editable post analysis.

Default resDNASEQ™ settings (cannot be changed)

- **Experiment Type** is **Quantitation-Standard Curve**
- **Chemistry** is **TaqMan™ Reagents**
- **Ramp Speed** is **Standard - 2hrs**

- c. Click **Next**.

Note: Experiment names cannot be changed after this step.

- In the **qPCR Method** pane of the **Setup** tab, view the default volume and cycling conditions (cannot be changed).




Figure 4 resDNASEQ™ template default cycling conditions

- Click **Next**.
- In the **Samples** pane of the **Setup** tab, enter the sample **Name** and **Dilution**. **Sample Volume** is not applicable. Add additional **Samples** if needed.

IMPORTANT! Do not change the **Targets**.

① — Samples (24) Add

Color	Name ↕	Dilution	Sample Volume
Orange	1A	1.00	0
Green	1B	1.00	0.00

- Samples** pane
- Add**—adds additional samples

6. In the **Samples** pane of the **Setup** tab, scroll to the right, then enter the spike information.

For more information on plate setup, see *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

- **Sample Volume**—not applicable; leave as default (0).
- **Spike Volume**—volume of DNA added to the PCR (set to 10).
- **Spike Standard Concentration**—expected spike amount per reaction.
- **Reference**—the non-spiked sample; the mean quantity of reference is subtracted during % recovery calculation.
- **Spike Input**—automatically calculated (double check if correct).

Note: If incorrect, be sure **Spike Volume** is set to 10 and **Spike Standard Concentration** is the expected pg spike per PCR reaction.

- (Optional) **Comments**
- **Protein Concentration**—Drug substance protein concentration (if Total DNA in pg DNA/mg Protein is required).

Samples (24) Add

	Spike Volume	Spike Standard Concentration
①	0	0.00
	0.00	0.00
	0.00	0.00
	0.00	0.00
②		

① Textbox—type in the value, or use the up and down arrows

② Scroll bar—scroll to find the spike parameter

For more information on plate setup, see *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

Figure 5 shows the default sample plate layout for the resDNASEQ™_5Std template. The layout is a 96-well plate with 8 rows (A-H) and 12 columns. The wells are populated as follows:

Row	Col 1	Col 2	Col 3	Col 4	Col 5	Col 6	Col 7	Col 8	Col 9	Col 10	Col 11	Col 12
A	NTC	NTC	NTC	Sample	Sample	Sample	Sample	Sample	Sample	STD 3000	STD 3000	STD 3000
B	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	STD 300	STD 300	STD 300
C	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	STD 30	STD 30	STD 30
D	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	STD 3	STD 3	STD 3
E	TS 1	TS 1	TS 1	Sample	Sample	Sample	Sample	Sample	Sample	STD 0.3	STD 0.3	STD 0.3
F	TS 2	TS 2	TS 2	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample
G	TS 3	TS 3	TS 3	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample
H	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample	Sample

Well counts: U 27, N 3, S 15, 69 (Empty). Progress: 0%.

Figure 5 resDNASEQ™_5Std template default sample plate layout

- ① Toolbar (in order: Undo, Redo, Copy, Paste, Delete, View)
- ② 3 No Template Control (NTC) samples
- ③ 3 default **Samples**
- ④ Standard curve dilutions (S) in triplicate

7. Click **Next**.

The **Run** tab is displayed.

8. Experiments are auto-saved in the software. To save, exit the experiment. The software prompts you to save changes. Click **Yes**.

Note: Clicking **Save As** will create a copy of the experiment.

9. Assemble the PCR reactions following the manufacturer's instructions for the reagents and following the plate layout set up in the template.

Start the run

Start the run in the AccuSEQ™ Software.

Option	Description
If the experiment is open	Click Start Run .
If the experiment is closed	<ol style="list-style-type: none"> 1. Open the experiment. 2. Click the Run tab. 3. Click Start Run.

The screenshot shows the AccuSEQ software interface. At the top, there is a navigation bar with tabs: Setup, Run, Result, Report, and Audit & E-Sign. The 'Run' tab is selected and highlighted with a circled '1'. Below the navigation bar, there are buttons for Analysis Settings, Save As, Save As Template, and Export Experiment. The main area displays 'Run Control' with sub-tabs for Plots and Pool Run Summary. Under 'Plots', 'Amplification' is selected. An 'Amplification Plot' is shown on the left, and a 'Sample' grid is on the right. At the bottom right, a blue 'Start Run' button is highlighted with a circled '2'.

- ① Run tab
- ② Start Run button

A message stating **Run has been started successfully** is displayed when the run has started.

Analyze the results

After the qPCR run is finished, use the following general procedure to analyze the results. For more detailed instructions see the *AccuSEQ™ Real-Time PCR Software v3.1 User Guide* (Pub. No. 100094287).

1. In the AccuSEQ™ Real-Time PCR Software, open your experiment, then navigate to the **Result** tab.

The screenshot displays the AccuSEQ™ Real-Time PCR Software interface. The top navigation bar shows the 'Result' tab selected. The left sidebar contains navigation icons. The main content area is divided into 'Result Summary' and 'Result Analysis' sections. The 'Amplification Plot' shows a graph of fluorescence (y-axis, 1e-3 to 2e+0) versus Cycle (x-axis, 5 to 40). The plot displays multiple curves for CHO (red) and IPC (orange) targets. The 'Results' section shows a grid of wells with various sample types and their corresponding analysis results. Callouts 1-4 indicate the 'Result' tab, 'Analysis Settings' button, plot horizontal scrollbar, and 'Analyse' button.

- ① **Result** tab
- ② **Analysis Settings**
- ③ Plot horizontal scrollbar
- ④ **Analyse** button

2. In the **Result Analysis** tab, review the Amplification Curve plots for amplification profiles in the controls, samples, and the standard curve. Ensure that auto threshold is selected.

Note: (Plasmid- KanR assay only) Ensure that the threshold was set to 0.04 for Kanamycin target and 0.02 for IPC.

3. In the **Result Analysis** tab, review the **QC Summary** for any flags in wells.

4. In the **Result Analysis** tab, review the **Standard Curve** plot. Verify the values for the Slope, Y-intercept, R², and Efficiency.

Note: The **Standard Curve** efficiency should be between 90-110% and the R²>0.99. If these criteria are not met, up to two points, not in the same triplicate, can be removed from the standard curve data, and the analysis repeated.

5. (Optional) Navigate to the **Report** tab to generate a report of the experiment, or to export results.
6. (Optional) To manually convert the copy number output to a mass measurement, multiply the copy number given by the AccuSEQ™ Real-Time PCR Software to the average molecular weight of the plasmid, then divide by the Avogadro constant.

$$\text{Mass (g)} = \frac{\text{Copy number} \times \text{Molecular weight (gmol}^{-1}\text{)}}{6.0221 \times 10^{23} \text{ mol}^{-1}}$$



Troubleshooting

Observation	Possible cause	Action
<p>Slope for the standard curve is outside the typical range, or R_2 value is significantly less than 0.99.</p>	<p>When applying detectors for standards, the Task and Quantity were applied to the wrong detector.</p> <p><i>or</i></p> <p>The incorrect Quantity was entered.</p> <p><i>or</i></p> <p>Adjust baseline settings.</p> <p><i>or</i></p> <p>Poor standard curve preparation technique (forgot to mix, inaccurate pipetting).</p>	<ol style="list-style-type: none"> 1. In the SDS software, from the plate document, double-click a well containing a DNA standard to view the Well Inspector. 2. Ensure that the correct Task and Quantity are applied to the correct detector, then reanalyze. 3. Compare std curve statistics using autobaseline or manual baseline. The upper limit of the manual baseline setting must be 2 cycles before uptick in amplification. Verify in R_n vs C_t linear view.
<p>ΔR_n and C_t values are inconsistent with replicates</p>	<p>Evaporation of reaction mixture from some wells occurred because the optical adhesive cover was not correctly sealed to the reaction plate or due to over-drying the eluates in PrepSEQ™.</p>	<ol style="list-style-type: none"> 1. Select the Component tab. Confirm that affected wells generated significantly less fluorescence than unaffected replicates. 2. Check the amount of solution in each well of the reaction plate. Confirm that the wells affected by evaporation contained less solution than unaffected wells, and corresponded with the inconsistent results. 3. For subsequent runs, ensure that the optical adhesive cover is correctly sealed to the reaction plate.

(continued)

Observation	Possible cause	Action
ΔR_n and C_t values are inconsistent with replicates	Incorrect volume of PCR reaction mix was added to some reactions.	<ol style="list-style-type: none"> 1. Select the Component tab. Confirm that affected wells generated significantly less fluorescence than unaffected replicates. 2. Select the Spectra tab. Confirm that the wells with the incorrect volume of PCR reaction mix generated significantly different amounts of fluorescence than the unaffected wells. 3. For subsequent runs, ensure the correct volume of PCR reaction mix.
Jagged amplification plots	Weak lamp or incorrect replacement.	Replace the lamp or make sure that the existing replacement is correct.
No defined amplification plots	<p>An incorrect detector was selected on the amplification plot.</p> <p>or</p> <p>An incorrect detector was applied to the reactions when setting up the plate document.</p>	<ol style="list-style-type: none"> 1. Confirm that the correct detector was selected on the amplification plot. 2. If the correct detector was not selected, then in the plate document, double-click a well to view the Well Inspector, verify that the detector settings are correct, and reanalyze.
Abnormal ΔR_n values or negative ΔR_n values.	Incorrect passive reference was selected when setting up the plate document.	<ol style="list-style-type: none"> 1. From the plate document, double-click a well to view the Well Inspector. 2. Ensure that ROX™ dye was selected as the Passive Reference.
Standard curve for plasmid DNA assays is outside of the 90–110% efficiency range.	Incomplete vortexing of low level standards.	Repeat reactions, ensuring that samples and standards are vortexed for 15-30 seconds.
Wide variance of C_t values of plasmid DNA samples.	Incomplete vortexing of samples.	



Use the kit with the 7500 Fast Real-Time PCR Instrument and AccuSEQ™ software v2.x

Required materials not supplied

Unless otherwise indicated, all materials are available through thermofisher.com. "MLS" indicates that the material is available from fisherscientific.com or another major laboratory supplier.

Item	Source
Instrument	
7500 Fast Real-Time PCR System with AccuSEQ™ software v2.x	Contact your local sales representative
Generic consumables	
Disposable gloves	MLS
Aerosol-resistant pipette tips	MLS
Pipettors: <ul style="list-style-type: none">• Positive-displacement• Air-displacement• Multichannel	MLS
Consumables for the 7500 Fast Real-Time PCR System	
MicroAmp™ Fast Optical 96-Well Reaction Plate with Barcode, 0.1 mL, 20 plates; for use with 7500 Fast Real-Time PCR System	4346906
MicroAmp™ Optical 96-Well Reaction Plate with Barcode & Optical Adhesive Films, 100 plates with covers; for use with 7300 and 7500 Fast Real-Time PCR Systems	4314320
MicroAmp™ Optical 8-Cap Strips, 300 strips	4323032
MicroAmp™ Optical Adhesive Film Kit, 20 covers, 1 applicator, 1 optical cover compression pad	4313663
MicroAmp™ Optical Adhesive Film	4360954

Create a plate document in the AccuSEQ™ software

1. In the home screen, select **Create Custom Experiment**.
2. In the **Experiment name** field, enter a unique name for the experiment.
3. Specify experiment information.
 - a. Select experiment type **Quantitation – Standard Curve**.
 - b. Select reagents **TaqMan™ Reagents**.
 - c. Select ramp speed **Standard**.
4. In the **Plate Setup** screen, select the **Define Targets and Samples** tab.
5. Specify target information.
 - a. Click **Add New Target**.
 - b. Enter a host cell name in the target name field.
 - c. Select reporter **FAM™** dye and quencher **NFQ_MGB** dye.

Note: For the Sf9 and Baculovirus assay, select reporter **VIC™** dye for the Sf9 target and quencher **NFQ_MGB** dye. Add an additional target, and select reporter **FAM™** dye for the baculovirus target and quencher **NFQ_MGB** dye.

 - d. Select a color for this target.
6. Specify IPC target information.
 - a. Click **Add New Target**.
 - b. Enter **IPC** in the target name field.
 - c. Select reporter **VIC™** dye and quencher **NFQ_MGB** dye.

Note: For the Sf9 and Baculovirus assay and the Plasmid DNA - Kanamycin Resistance assay, select the reporter **NED™** dye and quencher **NFQ_MGB** dye.

 - d. Select a color for this target.
7. Define new samples.
 - a. Click **Add New Sample**.
 - b. In **Sample Name**, add the names of the samples you want to define.
 - c. Click **Next**, or select the **Assign Targets and Samples** tab.
8. In the **Assign Targets and Samples** tab, define new targets.
 - a. Follow the instructions in the top of the tab to set up the standards, unknowns, and negative controls.


3. Make any necessary changes to the test sample labels.
 - **Sample Volume**—not applicable; leave as default (0).
 - **Spike Volume**—volume of DNA added to the PCR (set to 10).
 - **Spike Standard Concentration**—expected spike amount per reaction (for example, 10pg).
 - **Reference**—the non-spiked sample; the mean quantity of reference is subtracted during % recovery calculation.
 - **Spike Input**—automatically calculated (double check if correct).

Note: If incorrect, be sure **Spike Volume** is set to 10 and **Spike Standard Concentration** is the expected pg spike per PCR reaction.

4. Select **Save As** to save the new experiment as an EDS experiment file with the same name as entered in the **Experiment Name** field.
5. Load the plate into the instrument.
6. Click **Start Run**.
7. Select a run screen (**Amplification** plot, **Temperature** plot, or **Run method**) to monitor the progress of the run.

Analyze the results

After the qPCR run is finished, use the following general procedure to analyze the results:

1. In the toolbar, select **Analysis ▶ Analysis Settings**.
2. Click  (**Analyze**).
3. Select **Analysis ▶ QC Summary** in the left panel of the screen. Review the flag summary.
4. In the left panel, select **Analysis ▶ Standard Curve**. Verify the values for the Slope, Y-Intercept, R2, and Efficiency.
5. Select **File ▶ Print Report** to generate a hardcopy of the experiment, or click **Print Preview** to view and save the report as a *.pdf or *.html file.
6. Optional: Select **File ▶ Export**. In the **Export Data** menu, select file type *.xls. Click **Start Export**.



Use the kit with 7500 System SDS Software v1.5.1

Create the plate document, run the plate, and analyze the results with 7500 Fast SDS software

The following instructions apply only to the Applied Biosystems™ 7500 Fast instrument with SDS v1.x software. If you use a different instrument or software, refer to the applicable instrument or software documentation. Genomic residual DNA is measured in concentration, while plasmid residual DNA is measured in copy number.

Create a plate document

Residual DNA assays are duplex assays, containing sample DNA and Internal Positive Control (IPC).

Plate document: settings

If you have run the assay frequently, you can use the table below to enter settings in Plate Document fields. If you are a new user, follow the detailed instructions in the following sections.

Summary of settings for the Plate Document		
In this field...		Use these settings
Detector	resDNASEQ™ kit target cell lines	<ul style="list-style-type: none">• Single target assays: FAM™ dye• Sf9 and Baculovirus assay: VIC™ and FAM™ dyes (Select NFQ_MGB dye for quencher)
	IPC	<ul style="list-style-type: none">• Single target assays: VIC™ dye• Sf9 and Baculovirus assay and Plasmid DNA - Kanamycin Resistance assay: NED™ dye (Select NFQ_MGB for quencher)
PCR	Hold	Temp: 95°C Time: 10:00



(continued)

Summary of settings for the Plate Document		
In this field...		Use these settings
PCR	Cycling (Standard Mode)	Cycles: 40 Temp: 95°C Time: 0:15 Temp: 60°C Time: 1:00
Analysis	CHO, <i>E. coli</i> , HEK293, Human, MDCK, NS0, <i>Pichia</i> , Plasmid DNA - Kanamycin resistance, Sf9 and Baculovirus, and Vero	Automatic Baseline or Manual Baseline ^[1] Threshold: 0.2 Note: For CHO, the upper limit for manual baseline analysis is 12. For the Plasmid DNA - Kanamycin Resistance assay, change the threshold to 0.04 for KanR and 0.02 for IPC.

^[1] You can analyze the assay using Automatic or Manual Baseline, use the setting that yields the best standard curve.

Plate document: procedure

In the SDS software:

1. In the template **Assay** drop-down list, select **Absolute Quantification**.
2. In the **Run Mode** drop-down list, select **Standard 7500**.
3. Enter **resDNA_Template** in the **Plate** name field, then click **Next**.
4. Click **New Detector**:
 - a. Enter the name of the target cell line in the **Name** field.
 - b. Select reporter **FAM™** dye and quencher **NFQ_MGB** dye.

Note: For the Sf9 and Baculovirus assay, select reporter **VIC™** dye for the Sf9 target and quencher **NFQ_MGB** dye. Add an additional target, and select reporter **FAM™** dye for the baculovirus target and quencher **NFQ_MGB** dye.

- c. Select a color for the detector, then click **Create Another**.
5. Click **New Detector**:
 - a. Enter **IPC** in the **Name** field.
 - b. Select reporter **VIC™** dye and quencher **NFQ_MGB** dye.

Note: For the Sf9 and Baculovirus assay and the Plasmid DNA - Kanamycin Resistance assay, select the reporter **NED™** dye and the quencher **NFQ_MGB** dye.

 - c. Select a color for the detector, then click **OK**.
 - d. Select the detectors, then click **Add>>** to add the detectors to the document (plate).



6. Select **ROX™** dye as the passive reference dye, then click **Next**.
7. Select the applicable set of wells for the samples, then select the target cell line and **IPC** detectors for each well. The following figure shows an example plate layout:

										Standard Curve (pg)		
	1	2	3	4	5	6	7	8	9	10	11	12
A	TS1	TS1	TS1		TS1 ERC	TS1 ERC	TS1 ERC			NTC	NTC	NTC
B	TS2	TS2	TS2		TS2 ERC	TS2 ERC	TS2 ERC					
C	TS3	TS3	TS3		TS3 ERC	TS3 ERC	TS3 ERC			0.03 pg	0.03 pg	0.03 pg
D										0.3 pg	0.3 pg	0.3 pg
E										3 pg	3 pg	3 pg
F	NEG	NEG	NEG							30 pg	30 pg	30 pg
G										300 pg	300 pg	300 pg
H										3,000 pg	3,000 pg	3,000 pg

8. Set tasks for each sample type by clicking on the **Task Column** drop-down list:
 - a. NTC: target cell line detector task = **NTC**
 - b. NEG, test samples, and ERC wells: target DNA detector task = **Unknown**
 - c. IPC = **Unknown** for all wells
9. Set up the standard curve:
 - a. Select the wells.



- b. Assign the tasks (target DNA = **Standard**) and enter the appropriate Quantity for each set of triplicates.

Tube label	Row-wells	Task	Quantity	Label (pg)
SD 1	H-10, 11, 12	Standard	3,000	3,000 pg
SD 2	G-10, 11, 12	Standard	300	300 pg
SD 3	F-10, 11, 12	Standard	30	30 pg
SD 4	E-10, 11, 12	Standard	3	3 pg
SD 5	D-10, 11, 12	Standard	0.3	0.3 pg
SD 6 (for CHO, Vero, MDCK, and NS0 only)	C-10, 11, 12	Standard	0.03	0.03 pg

Note: The Plasmid DNA - Kanamycin Resistance assay creates a standard curve based on copy number (300,000 to 30 copies per reaction).

10. Select the **Instrument** tab, then set thermal-cycling conditions:

- Set the **thermal cycling reaction volume** to **30 µL**.
- Set the reaction to **Standard Mode**.
- Set the temperature and the time as shown in the following table:

Step	AmpliTaq Gold™ enzyme activation	PCR	
		Denature	Anneal/extend
	Hold	Cycle (40 Cycles)	
Temp (°C)	95	95	60
Time (mm:sec)	10:00	0:15	1:00

Refer to the applicable 7500 Fast Real-Time PCR Systems instrument manual for additional information.

11. In the **Analysis Settings** window, enter the following settings, then click **OK**:

- Select **Manual Ct**.
- In **Threshold**, enter **0.2**.

Note: For the Plasmid DNA - Kanamycin Resistance assay, change the threshold to 0.04 for KanR and 0.02 for IPC.

- Select **Automatic Baseline** or **Manual Baseline**.

Note: You can analyze the assays using Automatic or Manual Baseline, use the setting that yields the best standard curve. For CHO, the upper limit threshold for manual baseline analysis is 12.



12. Select **File** ▶ **Save as**, confirm that the file is named “resDNA_Template”, then select **SDS Templates (*.sdt)** in the **Save as type** drop-down list and close the template plate document.

Note: You can reuse the plate template document whenever you run the assay.

13. Close the saved template file.

Run the plate

1. In the SDS software, select **File** ▶ **New**, navigate to the **resDNA_Template** file (created in “Plate document: procedure” on page 45), then click **Open**.
2. In **Plate Name**, enter an appropriate experiment name, then click **Finish**.
3. Make any necessary changes to the test sample labels.
4. Select **Save As** to save the new experiment as an SDS experiment file.
5. Load the plate on the instrument.
6. Select the **Instrument** tab, save the document, then click **Start** to start the real-time qPCR run.

Analyze the results

After the qPCR run is finished, use the following general procedure to analyze the results.

1. Select the **Results** tab.
2. Select the **Amplification Plot** tab.
3. Verify the analysis settings, change as appropriate, then click **Analyze**.
4. Select the **Results tab** ▶ **Standard Curve** tab, then verify the Slope, Intercept, and R2 values.
5. Right-click the **Standard Curve**, select **Export as JPEG**, then click **OK**. Alternatively, press **PrintScreen**, then paste the image in a WordPad file.
6. Select the **Report tab** ▶ **Report**, then review the mean quantity and standard deviation for each sample.
7. Optional: Select **File** ▶ **Export** ▶ **Results**. In the **Save as type** drop-down list, select **Results Export Files (*.csv)**, then click **Save**.



Good laboratory practices

Work area setup and lab design

The sensitivity of this kit (and other PCR-based tests) enables amplification of minute quantities of DNA, necessitating precautions to avoid contamination of samples yet to be amplified (Kwok and Higuchi, 1989).

Process samples carefully to prevent contamination by human DNA. Wear gloves at all times and change them frequently. Close sample tubes when not in use. Limit aerosol dispersal by handling sample tubes and reagents carefully.

Good laboratory practices for PCR and RT-PCR

- Wear clean gloves and a clean lab coat.
 - Do not wear the same gloves and lab coat that you have previously used when handling amplified products or preparing samples.
- Change gloves if you suspect that they are contaminated.
- Maintain separate areas and dedicated equipment and supplies for:
 - Sample preparation and reaction setup.
 - Amplification and analysis of products.
- Do not bring amplified products into the reaction setup area.
- Open and close all sample tubes carefully. Avoid splashing or spraying samples.
- Keep reactions and components capped as much as possible.
- Use a positive-displacement pipettor or aerosol-resistant barrier pipette tips.
- Clean lab benches and equipment periodically with 10% bleach solution or DNA decontamination solution.

Avoiding false positives due to cross-contamination

To avoid false positives due to cross-contamination:

- Do not open tubes after amplification.
- Use different sets of pipettors when pipetting negative control, unknown, and positive control samples.

Note: Refer to “Prepare the PCR plate” on page 17 for best practice.



Safety



WARNING! GENERAL SAFETY. Using this product in a manner not specified in the user documentation may result in personal injury or damage to the instrument or device. Ensure that anyone using this product has received instructions in general safety practices for laboratories and the safety information provided in this document.

- Before using an instrument or device, read and understand the safety information provided in the user documentation provided by the manufacturer of the instrument or device.
- Before handling chemicals, read and understand all applicable Safety Data Sheets (SDSs) and use appropriate personal protective equipment (gloves, gowns, eye protection, and so on). To obtain SDSs, see the “Documentation and Support” section in this document.

Chemical safety



WARNING! GENERAL CHEMICAL HANDLING. To minimize hazards, ensure laboratory personnel read and practice the general safety guidelines for chemical usage, storage, and waste provided below. Consult the relevant SDS for specific precautions and instructions:

- Read and understand the Safety Data Sheets (SDSs) provided by the chemical manufacturer before you store, handle, or work with any chemicals or hazardous materials. To obtain SDSs, see the "Documentation and Support" section in this document.
- Minimize contact with chemicals. Wear appropriate personal protective equipment when handling chemicals (for example, safety glasses, gloves, or protective clothing).
- Minimize the inhalation of chemicals. Do not leave chemical containers open. Use only with sufficient ventilation (for example, fume hood).
- Check regularly for chemical leaks or spills. If a leak or spill occurs, follow the manufacturer cleanup procedures as recommended in the SDS.
- Handle chemical wastes in a fume hood.
- Ensure use of primary and secondary waste containers. (A primary waste container holds the immediate waste. A secondary container contains spills or leaks from the primary container. Both containers must be compatible with the waste material and meet federal, state, and local requirements for container storage.)
- After emptying a waste container, seal it with the cap provided.
- Characterize (by analysis if needed) the waste generated by the particular applications, reagents, and substrates used in your laboratory.
- Ensure that the waste is stored, transferred, transported, and disposed of according to all local, state/provincial, and/or national regulations.
- **IMPORTANT!** Radioactive or biohazardous materials may require special handling, and disposal limitations may apply.



WARNING! HAZARDOUS WASTE (from instruments). Waste produced by the instrument is potentially hazardous. Follow the guidelines noted in the preceding General Chemical Handling warning.



WARNING! 4L Reagent and Waste Bottle Safety. Four-liter reagent and waste bottles can crack and leak. Each 4-liter bottle should be secured in a low-density polyethylene safety container with the cover fastened and the handles locked in the upright position.

Biological hazard safety



WARNING! Potential Biohazard. Depending on the samples used on this instrument, the surface may be considered a biohazard. Use appropriate decontamination methods when working with biohazards.



WARNING! BIOHAZARD. Biological samples such as tissues, body fluids, infectious agents, and blood of humans and other animals have the potential to transmit infectious diseases. Conduct all work in properly equipped facilities with the appropriate safety equipment (for example, physical containment devices). Safety equipment can also include items for personal protection, such as gloves, coats, gowns, shoe covers, boots, respirators, face shields, safety glasses, or goggles. Individuals should be trained according to applicable regulatory and company/ institution requirements before working with potentially biohazardous materials. Follow all applicable local, state/provincial, and/or national regulations. The following references provide general guidelines when handling biological samples in laboratory environment.

- U.S. Department of Health and Human Services, *Biosafety in Microbiological and Biomedical Laboratories (BMBL)*, 6th Edition, HHS Publication No. (CDC) 300859, Revised June 2020
<https://www.cdc.gov/labs/pdf/CDC-BiosafetymicrobiologicalBiomedicalLaboratories-2020-P.pdf>
- Laboratory biosafety manual, fourth edition. Geneva: World Health Organization; 2020 (Laboratory biosafety manual, fourth edition and associated monographs)
www.who.int/publications/i/item/9789240011311

Documentation and support

Related documentation

Document	Publication number	Description
<i>resDNASEQ™ Quantitative DNA Kits Quick Reference</i>	4469837	For brief instructions on using the resDNASEQ™ Quantitative DNA Kits.
<i>resDNASEQ™ Quantitative E1A DNA Fragment Length Kit User Guide</i>	MAN0025643	For instructions on running the resDNASEQ™ Quantitative E1A DNA Fragment Length Kit (Cat. No. A51969).
<i>PrepSEQ™ Sample Preparation Kits User Guide</i>	4469838	For information on preparing samples for extraction.
<i>PrepSEQ™ Residual DNA Sample Preparation Kit Quick Reference</i>	4469839	For brief instructions on preparing samples for extraction.
<i>AccuSEQ™ Real-Time PCR Software v3.1 User Guide</i>	100094287	For information on AccuSEQ™ Real-Time PCR Software v3.1 with the QuantStudio™ 5 Real-Time PCR System
<i>AccuSEQ™ Real-Time PCR Software v3.1 Quick Reference</i>	100094288	For basic information on AccuSEQ™ Real-Time PCR Software v3.1 with the QuantStudio™ 5 Real-Time PCR System
<i>Applied Biosystems™ 7300/7500/7500 Fast Real-Time PCR System Getting Started Guide: Absolute Quantitation using Standard Curve</i>	4347825	For information on the 7500 Fast instrument.
<i>AccuSEQ™ software: Custom Quick Reference Card</i>	4425585	For information on AccuSEQ™ software with the 7500 Fast instrument.

Customer and technical support

Visit thermofisher.com/support for the latest service and support information.

- Worldwide contact telephone numbers
- Product support information
 - Product FAQs
 - Software, patches, and updates
 - Training for many applications and instruments

- Order and web support
- Product documentation
 - User guides, manuals, and protocols
 - Certificates of Analysis
 - Safety Data Sheets (SDSs; also known as MSDSs)

Note: For SDSs for reagents and chemicals from other manufacturers, contact the manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

References

Kwok, S., and Higuchi, R. 1989. Avoiding false positives with PCR. *Nature* 339:237–238.

