**Thermo Scientific**  
**ThermoFlex™**  
**Recirculating Chillers**  
**(Deluxe Controller)**

Thermo Scientific Manual P/N U00939  
Rev. 06/21/2019

**Multilingual Quick  
Start Guides**

**Multilingual Essential  
Safety Instructions**

**Installation**

**Operation**

**Preventive  
Maintenance**

**Troubleshooting**

Visit our Web site at:

<http://www.thermofisher.com/tc>  
Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

**Thermo Scientific**  
**ThermoFlex™**  
**Recirculating Chillers**  
**(Deluxe Controller)**

Thermo Scientific Manual P/N U00939  
Rev. 06/21/2019

Label 1

---

**Multilingual Quick Start  
Guides**

**Multilingual Essential  
Safety Instructions**

**Installation**

**Operation**

**Preventive  
Maintenance**

**Troubleshooting**

Label 2

---

Visit our Web site at:

**<http://www.thermofisher.com/tc>**

Product Service Information, Applications  
Notes, SDS Forms, e-mail.

Voice Info: (800) 258-0830

## **Thermo Fisher Scientific**

25 Nimble Hill Road  
Newington, NH 03801  
Tel : (800) 258-0830 or  
(603) 436-9444  
Fax : (603) 436-8411  
[www.thermofisher.com/tc](http://www.thermofisher.com/tc)

## **Sales, Service, and Customer Support**

25 Nimble Hill Road  
Newington, NH 03801  
Tel: (800) 258-0830  
Sales: 8:00 am to 5:00 pm  
Service and Support: 8:00 am to 6:00 pm Monday  
through Friday (Eastern Time)  
Fax: (603) 436-8411  
[service.tc.us@thermofisher.com](mailto:service.tc.us@thermofisher.com)

Dieselstrasse 4  
D-76227 Karlsruhe, Germany  
Tel : +49 (0) 721 4094 444  
Fax : +49 (0) 721 4094 300  
[info.tc.de@thermofisher.com](mailto:info.tc.de@thermofisher.com)

Building 6, No. 27  
Xin Jinqiao Rd., Shanghai 201206  
Tel : +86(21) 68654588  
Fax : +86(21) 64457830  
[info.china@thermofisher.com](mailto:info.china@thermofisher.com)

## **Statement of Copyright**

Copyright © 2019 Thermo Fisher Scientific. All rights reserved.

This manual is copyrighted by Thermo Fisher Scientific.

Users are forbidden to reproduce, republish, redistribute, or resell any materials from this manual in either machine-readable form or any other form.

# Contents

## Quick Start

<b>Preface</b>	<b>i</b>
Compliance .....	i
WEEE .....	i
After-Sale Support .....	i
Feedback.....	ii
Warranty .....	ii
Unpacking.....	ii
<b>Section 1 Safety</b>	<b>1-1</b>
Warnings.....	1-1
<b>Section 2 General Information</b>	<b>2-1</b>
Description .....	2-1
Specifications.....	2-1
<b>Section 3 Installation</b>	<b>3-1</b>
Site Requirements .....	3-1
Electrical Requirements .....	3-2
Hard Wire Installation.....	3-9
Plumbing Requirements.....	3-10
Process Fluid Requirements.....	3-12
Compatibility with Approved Fluids .....	3-13
Additional Fluid Information .....	3-15
Process Water Quality - Standards and Recommendations .....	3-15
Facility Water Requirements - Standards and Recommendations.....	3-16
Facility Water Requirements (water-cooled chillers).....	3-17
Fluid Filter Bag.....	3-19
Priming .....	3-19
Initial Filling Requirements .....	3-20
Fluid Top Off.....	3-21
<b>Section 4 Operation</b>	<b>4-1</b>
Deluxe Controller .....	4-1
Initial Start Up.....	4-2
Daily Start Up.....	4-3
Initial Fill Mode.....	4-4
Process Fluid Temperature Display .....	4-5
SETPOINT .....	4-6
Status Display .....	4-7
MENU Display .....	4-8
Main Menu Tree.....	4-9
Shut Down.....	4-22

<b>Section 5</b>	<b>Additional Options/Accessories</b> .....	<b>5-1</b>
	Auto Refill .....	5-1
	Internal DI Cartridge .....	5-2
	P1 P2 T1 Pump Pressure Relief Valve (Internal Configuration) .....	5-3
	P1 P2 T1 Pump Pressure Relief Valve (External Configuration) .....	5-4
	Flow Control with Flow Readout.....	5-5
	P1 P2 T1 Pump Pressure Relief with Flow Readout.....	5-5
	T5 Pump Flow Control.....	5-6
	Anti Drainback.....	5-6
	SEMI.....	5-6
	Other Accessories.....	5-10
<b>Section 6</b>	<b>Preventive Maintenance</b> .....	<b>6-1</b>
	Preventive Maintenance Timer .....	6-1
	DIAGNOSTIC.....	6-2
	Fluid Bag Filter.....	6-3
	Fluid Diffuser .....	6-3
	Fluid Maintenance .....	6-4
	Reservoir Cleaning.....	6-4
	Condenser Filter.....	6-5
	Chiller Surface .....	6-6
	Hoses .....	6-6
	Testing the Alarm Features .....	6-7
	DI Filter (Optional) .....	6-7
	Preventive Maintenance Messages .....	6-8
<b>Section 7</b>	<b>Troubleshooting</b> .....	<b>7-1</b>
	Messages.....	7-1
	Faults and Warnings .....	7-1
	Safety Alarm .....	7-9
	Error .....	7-10
	MESSAGES.....	7-11
	Checklist.....	7-12
	TUNING .....	7-14
<b>Section 8</b>	<b>Additional Information</b> .....	<b>8-1</b>
	Draining.....	8-1
	CALIBRATION .....	8-3
	Wetted Materials.....	8-4
	Decommissioning/Disposal .....	8-5
	Shipment/Storage.....	8-5
<b>Appendix A</b>	Country Specific 230 VAC, 50 Hz, 1 Ø Requirements	
<b>Appendix B</b>	Voltage Configuration Instructions	
<b>Appendix C</b>	Serial Communications	
<b>Appendix D</b>	Analog I/O and Remote Sensor	

**Declaration of Conformity**

**WARRANTY**

**This quick start guide is intended for initial start up only. For all other procedures you must refer to the manual. Also, if any of these steps are not clear download the manual before proceeding.**

- The chiller is designed for indoor use only. Never place the chiller in a location with excessive heat, moisture, inadequate ventilation, or corrosive materials.
- Connect the chiller to a properly grounded outlet.
- Refrigerants used are heavier than air and will replace the oxygen causing loss of consciousness. Contact with leaking refrigerant will cause skin burns. Refer to the chiller's nameplate and the manufacturer's most current MSDS.
- Move the chiller with care. Sudden jolts or drops can damage its components. Always turn the equipment off and disconnect it from its supply voltage before moving it.
- Never operate damaged or leaking equipment.
- If your chiller is equipped with a positive displacement pump (P1 or P2), ensure your application plumbing lines and fittings are rated to withstand a minimum of 185 psi.
- Do not use a Deionization (DI) filter cartridge with inhibited EG or inhibited PG. A DI filter will remove inhibitors from the solution rendering the fluid ineffective against corrosion protection.
- Use only the approved fluids shown in Table 1. Before using any fluid or performing maintenance where contact with the fluid is likely, refer to the manufacturer's MSDS for handling precautions.
- To prevent freezing/glaizing of the plate exchanger, ThermoFlex7500-10000 chillers require the use of 50/50 EG/water or 50/50 PG/water below 10°C process temperature.


Figure A


Figure B

- What you need to get started:**
- An adjustable wrench
  - Facility water supply and return (water-cooled chillers)
  - Appropriate hose or plumbing
  - Appropriate size clamps or connection type
  - Teflon® Tape or appropriate sealant

Facility Water Connections (ENPT)  
 ThermoFlex1400 - 5000 Inlet/Outlet 1/2" cast bronze  
 ThermoFlex7500 - 10000 Inlet/Outlet 3/4" cast bronze

**Table 1 - Approved Fluids:**

Use of any other fluid will void the manufacturer's warranty.

**Standard Temp Chillers**

- Filtered/single distilled water (pH 7-8)
- Deionized water (1-3 MΩ-cm, compensated)
- Distilled water with Nalco biocide and inhibitor
- Distilled water with chlorine (5 ppm)
- 0 - 75% Laboratory Grade Ethylene Glycol/Water
- 0 - 75% Laboratory Grade Propylene Glycol/Water

**High Temp Chillers**

- Filtered/single distilled water (pH 7-8)\*
- 0 - 50% Laboratory Grade Ethylene Glycol/Water
- 0 - 50% Laboratory Grade Propylene Glycol/Water
- \* to 88°C for chillers with P1 and P2 pumps
- \* to 90°C for all other chillers

**Process Fluid Connections (FNPT)**

- ThermoFlex900 - 10000 P1 P2 T0 T1 1/2" cast bronze
- ThermoFlex3500 - 5000 P3 P4 3/4" cast bronze
- ThermoFlex7500 - 10000 P3 P5 T5 1" wrought copper
- Inlet - Same size as outlet
- all chillers stainless steel
- Supplied Adapters
- P1 P2 T0 T1 1/2" x 3/8" Polyethylene and 1/2" x 1/2" Nylon
- P3 P4 3/4 MPT x 1/2 barb PVC
- P3 P5 T5 1" MPT x 1" barb PVC and 1" MPT x 3/4" barb PVC

- 1** Remove all the plastic shipping plugs (2 or 4).
- 
- Locations vary with chiller size and selected options. The labels identify the exact location.

See Figure B.

- 2** Connect the ThermoFlex PROCESS OUTLET (A) to the fluid inlet on your application. Connect the ThermoFlex PROCESS INLET (B) to the fluid outlet on your application. Ensure the connections are sealed and secure. For air-cooled chillers skip to Step 4.
- 
- Never connect process fluid lines to your facility water supply or to any pressurized liquid source.**

See Figure B.

- 3** Connect the ThermoFlex FACILITY OUTLET (A) to a facility water return or drain. Connect the ThermoFlex FACILITY INLET (B) to a facility water supply. Ensure the connections are sealed and secure.
- 
- For water-cooled only.**
- Facility Water Maximum Inlet Pressure must not exceed 150 PSIG. Facility Water Maximum Pressure Differential must not exceed 50 PSID.**

See Figure B.

- 4** Refer to the name plate on the rear of the chiller and verify the appropriate voltage. For chillers supplied with a line cord, insert female end of line cord into chiller and then insert male end of line cord into power outlet. (The line cord is located under the shipping crate's lid. Do not discard the lid until the cord is located.)
- 
- Never operate the chiller with a damaged line cord.**
- Note:** ThermoFlex900-5000 chillers equipped with the Variable Voltage or Global Voltage option have a voltage configuration panel. Refer to the Voltage Instruction Sheet shipped with the chiller, or see manual Appendix B.
- Note:** For chillers requiring hard wiring see Section 3 in the manual.

See Figure B.

- 5** If equipped, place the optional GFCI breaker located on the rear to the up position. For ThermoFlex900 through 10000 chillers, place the circuit protector to the on ( I ) position.
- 
- The circuit protector is not intended to act as a disconnecting means.**

See Figure B.

- 6** Never operate the chiller without process fluid in the reservoir or without the fluid filter bag installed.
- 
- Gently pull up on the plastic funnel housing to remove it and install the supplied filter bag. Reinstall the housing.
- Remove the reservoir cap from the housing by unscrewing it counterclockwise.

See Figure A.

- 7** Slowly fill reservoir with clean process fluid (see Table 1), utilizing the reservoir cap, hand tight. Since the reservoir capacity may be small compared to your application and air may need to be purged from the lines, have extra fluid on hand to keep the system topped off when external circulation is started if the fluid level drops too low the chiller will shut down to prevent the pump from running dry. Acknowledge the alarm and re-start the priming process.
- 
- Note:** Be careful not to fill the reservoir above MAX LEVEL fill line. This will result in a unit over flow error (Over Flow) which will cause the chiller to shut down.

See Figure A.

- 8** Press
- 
- The controller will display **QUICK START**.
- Note:** Before proceeding, if the chiller is equipped with a deionization filter cartridge, refer to Section 5 in the manual for installation procedure.
- Please see reverse side for additional information.**
- 

See Figure A.

QUICK START	
<b>UNITS</b>	
Setpoint	20
TEMPERATURE	
PRESSURE	
FLUID LEVEL	
AUTO REFILL*	
FLOW*	
Line Frequency*	60HZ
<input type="checkbox"/> Auto Restart	
<input type="checkbox"/> Audible Alarms	
RA FAN SPEED MODE*	
Care Level	off
SERIAL COMM - DCOM*	
ANALOG COMM - ACOM*	
RESISTIVITY*	
mm/dd/yy	hh:mm:ss
◀ Quick Start Done	
*Displays only on chillers equipped with the option.	

**NOTE** Any setting can be changed after the chiller is started.

The controller can display only six lines of the QUICK START menu at a time.

**Use the down arrow to scroll through and highlight each line. If a change to any of the default settings is needed, with the desired line highlighted press the enter button.**

If the text on any setting is all capital letters, the setting has a sub-menu. Pressing enter will bring up the sub-menu. The sub-menus allow you to view/change the applicable values.

After pressing enter the line will start to flash indicating the value can be changed using the up/down arrows. Once the desired value is displayed, press enter again to stop the flashing and accept the new value. When all the desired changes are made press the left arrow or esc key to return to the QUICK START menu.

If the chiller exceeds the fault value it will shut down, the controller will display a fault message and, if enabled, sound the alarm. If the chiller exceeds the warning value the chiller will continue to run, the controller will display a warning message and, if enabled, sound the alarm.

Lines that are not all capital letters indicate the changes can be made directly on the QUICK START menu (e.g., Setpoint and Line Frequency). Use the same procedure to change these values.

If the line has a blank box (e.g., Auto Restart) pressing enter with that line highlighted will turn the box black. A black box indicates the feature is enabled. Press enter again to turn the box blank. A blank box indicates the feature is disabled.

Display	Indication	Range	Default
UNITS	Controller temperature, pressure and flow display scales. (Flow display is optional.)	C or F psi, bar or kPa gpm or lpm	C psi gpm
Setpoint	Setpoint value.	+5°C to +40°C (+5°C to +90° for high temp chillers)	+20°C
TEMPERATURE	Temperature fault and warning values.	+2°C to +43°C (+93°C for high temp chillers)	High +42°C (+92°C for high temp chillers) Low +3°C
PRESSURE	Pressure fault, warning and time delay values.  (The delay sets the length of time needed after a pressure fault before the chiller shuts down.)	Pump dependent - refer to Table 1  Time Delay: 0 to 30 seconds	Pump dependent - refer to Table 1  High Time Delay: 0 seconds (60 seconds for P3 - P5 pumps) Low Time Delay: 10 seconds
FLUID LEVEL	Reservoir low level fault and warning values.	Heater dependent - refer to Table 3	Heater dependent - refer to Table 3
AUTO REFILL	Optional auto refill values. (On is the % of fluid level in the reservoir needed to turn refill on. Time Out is the maximum time the option will operate. Setting the time to 0 disables the option.)	On: 70%- 100%  Time Out: 0 - 900 seconds	Heater dependent, see manual
FLOW	Process fluid flow fault and warning values.	Pump dependent - refer to Table 2	Pump dependent - refer to Table 2
Line Frequency	The incoming frequency (Chillers with a P3 - P5 pump and the capability to run on either 50 Hz or 60 Hz only. The selected frequency automatically adjusts the firmware's fixed high pressure default setting.)	50 Hz or 60 Hz	60 Hz
Auto Restart	Enables the auto restart.	<input type="checkbox"/> or <input checked="" type="checkbox"/>	<input type="checkbox"/>
Audible Alarms	Enables the audible alarm.	<input type="checkbox"/> or <input checked="" type="checkbox"/>	<input type="checkbox"/>
RA FAN SPEED MODE	Fan speed. ThermoFlex2500 air-cooled chillers only (Auto allows the fan to run under the conditions listed in the manual, see Section 3. On allows the fan to run at high speed all the time.)	On or Auto	Auto
Care Level	The preventive care cleaning frequency reminder for the chiller's air and fluid filters.	off, 1 (1000 hours), 2 (2000 hours) 3 (3000 hours)	off off
SERIAL COMM - DCOM	Optional serial communications. (See manual for additional information.)	off, RS232, RS485	off
ANALOG COMM - ACOM	Optional analog communications. (See manual for additional information.)	Refer to manual	
RESISTIVITY	Enables/configures the resistivity option. (See manual for additional information.)	<input type="checkbox"/> or <input checked="" type="checkbox"/> Setpoint: 0.2 to 3.0 MΩ-cm Interval: 0.1 to 0.5 MΩ-cm Warning High: 0.0 to 3.5 MΩ-cm Warning Low: 0.0 to 3.5 MΩ-cm	<input type="checkbox"/> Setpoint: 1.0 MΩ-cm Interval: 0.1 MΩ-cm Warning High: 3.0 MΩ-cm Warning Low: 0.5 MΩ-cm
mm/dd/yy	Sets the date. Some error messages display the date of occurrence.		
hh:mm:ss	Sets the time. Some error messages display the time of occurrence.		
Quick Start Done	To end the quick start procedure and save changes press the enter button. To leave the quick start and not save press the left arrow or esc button. In either case the screen will go blank.		

Table 1	Fault Range	High Default	Low Default
T0 T1 P1 P2 Pumps:	3 to 105 PSI	105 PSI	3 PSI
T5 Pumps:	2 to 105 PSI	105 PSI	4 PSI
P3 P4 P5 Pumps:	See Manual	See Manual	4 PSI

Table 2	Range	High Default	Low Default
T0 T1 P1 P2 Pumps:	0.0 to 10.5 GPM	0.0 GPM	0.0 GPM
T5 Pumps:	0.0 to 15.0 GPM	0.0 GPM	0.0 GPM
P3 P4 P5 Pumps:	0.0 to 30.0 GPM	0.0 GPM	0.0 GPM

Table 3	Warning Range	Default
Heater		
None	6 - 100%	6%
1 kW:	58 - 100%	58%
2.3 kW:	58 or 69 - 100%	58 or 69%
5.0 kW:	87 - 100%	87%
4.6 kW:	87 - 100%	87%
Heater		
None	0 - 100%	0%
1 kW:	52 - 100%	52%
2.3 kW:	52 or 63 - 100%	52 or 63%
5.0 kW:	81 - 100%	81%
4.6 kW:	81 - 100%	81%

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# Preface

## Compliance Third Party:

Designed to comply and listed/certified to:

UL 61010-1:2012, Third Edition

CAN/CSA-C22.2 No.61010-1-12, Third Edition

## European Union ( EU )

The Declaration of Conformity is located in the back of this manual.

## WEEE

This product is required to comply with the European Union's Waste Electrical & Electronic Equipment (WEEE) Directive 2012/19/EU. It is marked with 'wheelie bin' symbol.

Thermo Fisher Scientific has contracted with one or more recycling/ disposal companies in each EU Member State, dispose of or recycle this product through them. Further information on Thermo Fisher Scientific's compliance with these Directives is available at [www.thermofisher.com/WEEERoHS](http://www.thermofisher.com/WEEERoHS)

## After-sale Support

Thermo Fisher Scientific is committed to customer service both during and after the sale. If you have questions concerning the chiller operation, or questions concerning spare parts or Service Contracts, call our Sales, Service and Customer Support phone number, see this manual's inside cover for contact information.

When calling, please refer to the labels on the inside cover. These labels list all the necessary information needed to properly identify your chiller.

## Feedback

We appreciate any feedback you can give us on this manual. Please e-mail us at [tcmanuals@thermofisher.com](mailto:tcmanuals@thermofisher.com). Be sure to include the manual part number and the revision date listed on the front cover.

## Warranty

Thermo Scientific ThermoFlex chillers have a warranty against defective parts and workmanship for 24 months (**excluding MD1/MD2 Magnetic Drive and P1/P2 Positive Displacement pumps which are warranted for 12 months**) from date of shipment. See back page for more details.

## Unpacking

If the chiller has a line cord it is located under the shipping crate's lid. Do not discard the lid until the cord is located.

Locate the reservoir fluid filter bag and ensure it installed before the chiller is operated. See Section 3.

Retain all cartons and packing material until the chiller is operated and found to be in good condition. If it shows external or internal damage contact the transportation company and file a damage claim. Under ICC regulations, this is your responsibility.

### Out of Box Failure

An Out of Box Failure is defined as any product that fails to operate in conformance with sellers published specifications at initial power up. Install the chiller in accordance with manufacturer's recommended operating conditions within 30 days of shipment from the seller.

Any Temperature Control product meeting the definition of an Out of Box Failure must be packed and shipped back in the original packaging to Thermo Fisher Scientific for replacement with a new chiller; seller to pay the cost of shipping. Customer must receive a Return Material Authorization (RMA) from Thermo Fisher prior to shipping.

# Section 1 Safety

## Safety Warnings

Make sure you read and understand all instructions and safety precautions listed in this manual before installing or operating your chiller. If you have any questions concerning the operation or the information in this manual, please contact us. See inside cover for contact information.

**DANGER** indicates an imminently hazardous situation which, if not avoided, *will* result in death or serious injury.

**WARNING** indicates a potentially hazardous situation which, if not avoided, *could* result in death or serious injury.

**CAUTION** indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It is also be used to alert against unsafe practices.

The lightning flash with arrow symbol, within an equilateral triangle, is intended to alert the user to the presence of non-insulated "dangerous voltage" within the chiller's enclosure. The voltage magnitude is significant enough to constitute a risk of electrical shock.

This label indicates read the manual.

**Performance of installation, operation, or maintenance procedures other than those described in this manual may result in a hazardous situation and may void the manufacturer's warranty and safety compliance. ▲**

**Observe and never remove warning labels. ▲**

**Never place the chiller in a location where excessive heat, moisture, or corrosive materials are present. ▲**

**The chiller's construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲**

**Never operate equipment with damaged power cords. ▲**

**The circuit protector located on the rear is not intended to act as a disconnecting means. ▲**

**Never operate the chiller with panels removed. ▲**

**Never operate the chiller without process fluid in the reservoir. ▲**

**Never connect the process fluid inlet or outlet fittings to your building water supply or any water pressure source. ▲**

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**To prevent freezing/glazing of the plate exchanger, ThermoFlex7500 through ThermoFlex10000 chillers require the use of 50/50 EG/water or 50/50 PG/water below 10°C process temperature. ▲**

**When using a process fluid mixture of ethylene glycol and water or propylene glycol and water, check the fluid concentration and pH on a regular basis. Changes in concentration and pH can impact system performance. See Section 3. ▲**

**Do not use automotive antifreeze. Commercial antifreeze contains silicates that can damage the pump seals. Use of automotive antifreeze will void the manufacturer's warranty. ▲**

**Many refrigerants which may be undetectable by human senses are heavier than air and will replace the oxygen in an enclosed area causing loss of consciousness. Contact with leaking refrigerant will cause skin burns. Refer to the chiller's nameplate for the type of refrigerant used and then the refrigerant's SDS for additional information. ▲**

**Drain the chiller before it is transported and/or stored, see Draining in Section 8. Store the chiller in the temperature range -25°C to 60°C (with packaging), and <80% relative humidity. ▲**

**Always turn off the chiller and disconnect the power cord from the power source before performing any service or maintenance procedures, or before moving. ▲**

**Transport the chiller with care. Sudden jolts or drops can damage its components. ▲**

**Never operate damaged or leaking equipment. ▲**

**Refer service and repairs to a qualified technician. ▲**

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## **Fluorinated Greenhouse Gases**

The following information is included to comply with REGULATION (EU) No 517/2014 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 April 2014 on fluorinated greenhouse gases:

**This product contains fluorinated greenhouse gases in a hermetically sealed system.**

**If a leak in the sealed system is detected, the operator shall repair without undue delay.**

Refer to the F-Gas Declaration of Conformity for additional information.

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**Declaration of conformity with Article 14 of Regulation (EU) No 517/2014 of the European Parliament and of the Council**

We, Thermo Fisher Scientific (Asheville) LLC, F-Gas Portal Registration Number 23643, declare under our sole responsibility that when placing on the market pre-charged equipment, which we import to or manufacture in the Union, the hydrofluorocarbons contained in that equipment are accounted for within the quota system referred to in Chapter IV of Regulation (EU) No 517/2014 as:

We hold authorisation(s) issued in accordance with Article 18(2) of Regulation (EU) No 517/2014 and registered in the registry referred to in Article 17 of that Regulation, at the time of release for free circulation to use the quota of a producer or importer of hydrofluorocarbons subject to Article 15 of Regulation (EU) No 517/2014 that cover(s) the quantity of hydrofluorocarbons contained in the equipment.

The hydrofluorocarbons contained in the equipment have been placed on the market in the Union, subsequently exported and charged into the equipment outside the Union, and the undertaking that placed the hydrofluorocarbons on the market made a declaration stating that the quantity of hydrofluorocarbons has been or will be reported as placed on the market in the Union and that it has not been and will not be reported as direct supply for export in the meaning of Article 15(2)(c) of Regulation (EU) No 517/2014 pursuant to Article 19 of Regulation (EU) No 517/2014 and Section 5C of the Annex to Commission Implementing Regulation (EU) No 1191/2014.

January 1, 2017

\_\_\_\_\_  
Date

Mark Pearson, Director Global Regulatory Affairs

\_\_\_\_\_  
[name and position of legal representative]

\_\_\_\_\_  
[signature of legal representative]

Manufacturer:

Thermo Fisher Scientific (Asheville) LLC  
275 Aiken Road  
Asheville, NC 28804  
U.S.A.

EU Only Representative:

Thermo Electron LED GmbH  
Robert-Bosch-Strasse 1  
D-63505 Langenselbold  
Germany  
**VAT ID Number: DE 812 403 137**

## Section 2 General Information

### Description

The Thermo Scientific ThermoFlex™ recirculating chillers are designed to provide a continuous supply of fluid at a constant temperature and flow rate. The chiller consists of an air-cooled or water-cooled refrigeration system, heat exchanger, recirculating pump, reservoir and a microprocessor controller.

### Specifications

	ThermoFlex900	ThermoFlex1400	ThermoFlex2500
<b>Standard Chiller Process Fluid Temperature/Setpoint Range</b>	+5°C to +40°C +41°F to +104°F	+5°C to +40°C +41°F to +104°F	+5°C to +40°C +41°F to +104°F
<b>High-Temperature Chiller Process Fluid Temperature/Setpoint Range</b>	Not Available Not Available	+5°C to +90°C +41°F to +194°F	+5°C to +90°C +41°F to +194°F
<b>Ambient Temperature Range All Chillers</b>	+10°C to +40°C +50°F to +104°F	+10°C to +40°C +50°F to +104°F	+10°C to +40°C +50°F to +104°F
<b>Temperature Stability</b>	±0.1°C	±0.1°C	±0.1°C
<b>Cooling Capacity at 20°C</b> 60Hz 50Hz	900W (3074BTU) 750W (2561 BTU)	1400W (4781 BTU) 1170W (3996BTU)	2500W (8538BTU)* 2200W (7513BTU)
*To meet this specification, the ThermoFlex2500 air-cooled chillers require the fan to be operating in the high-speed mode, see Section 3			
<b>Heater Size</b> 208V/230V	Not Available	1.0kW/1.2kW or 2.3kW/2.8kW	2.3kW/2.8kW
<b>Refrigerant</b>	R134A	R134A	R134A
<b>Reservoir Volume</b> Gallons Liters	1.9 7.2	1.9 7.2	1.9 7.2
<b>Footprint or Dimensions (H x W x D)</b> Inches Centimeters	27.3x 14.2x 24.6 69.2x 36.0x 62.4	27.3x 14.2x 24.6 69.2x 36.0x 62.4	29.0x 17.2x 26.7 73.6x 43.6x 68.0
<b>Weight P2 Pump (empty)</b> lb kg	130.5 59.2	130.5 59.2	175.5 79.6
<b>Pumping Capacity</b>			
P1/MD1 - Positive Displacement 60Hz 50Hz		2.1 gpm @ 60 psig (7.9 lpm @ 4.1 bar) 1.7 gpm @ 60 psig (6.4 lpm @ 4.1 bar)	
P2/MD2 - Positive Displacement 60Hz 50Hz		4.0 gpm @ 60 psig (15.1 lpm @ 4.1 bar) 3.3 gpm @ 60 psig (12.5 lpm @ 4.1 bar)	
T0 - Turbine 60Hz* 50Hz*		2.0 gpm @ 60 psid (7.6 lpm @ 4.1 bar) 1.3 gpm @ 60 psid (4.9 lpm @ 4.1 bar)	
T1 - Turbine 60Hz* 50Hz*		3.5 gpm @ 60 psid (13.3 lpm @ 4.1 bar) 2.5 gpm @ 60 psid (9.5 lpm @ 4.1 bar)	

\* Pumping capacity pressure values for turbine pumps are differential pressures between the inlet and the outlet of the chiller.

- Cooling capacity based on P2 pumps with no backpressure. Heat input from the pump will result in a reduction in cooling capacity. The cooling capacity reduction will vary based on the pump chosen as well as pump backpressure and flow.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Additional dimensions are at the end of this section, add 1/8" (3mm) to height for SEMI.
- Add 5 pounds (2 kilograms) for global voltage chillers.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

# Specifications

	ThermoFlex3500	ThermoFlex5000
<b>Standard Chiller Process Fluid Temperature/Setpoint Range</b>	+5°C to +40°C +41°F to +104°F	+5°C to +40°C +41°F to +104°F
<b>High-Temperature Chiller Process Fluid Temperature/Setpoint Range</b>	+5°C to +90°C +41°F to +194°F	+5°C to +90°C +41°F to +194°F
<b>Ambient Temperature Range All Chillers</b>	+10°C to +40°C +50°F to +104°F	+10°C to +40°C +50°F to +104°F
<b>Temperature Stability</b>	±0.1°C	±0.1°C
<b>Cooling Capacity at 20°C</b> 60Hz 50Hz	3500W (11953BTU) 3050W (10416BTU)	5000W (17076BTU) 4400W (15027BTU)
<b>Heater Size</b> 208V/230V	2.3kW/2.8kW	2.3kW/2.8kW
<b>Refrigerant</b>	R407C	R407C
<b>Reservoir Volume</b> Gallons Liters	1.9 7.2	1.9 7.2
<b>Footprint or Dimensions (H x W x D)</b> Inches Centimeters	38.9 x 19.3 x 30.9 98.7 x 48.8 x 78.4	38.9 x 19.3 x 30.9 98.7 x 48.8 x 78.4
<b>Weight P1/ P2/P3/P4 (empty)</b> lb kg	264/264/270/303 120/120/123/138	NA/264/270/303 NA/120/123/138
<b>Pumping Capacity</b>		
P1/MD1 - Positive Displacement 60Hz 50Hz	2.1 gpm @ 60 psig (7.9 lpm @ 4.1 bar) 1.7 gpm @ 60 psig (6.4 lpm @ 4.1 bar)	Not Available Not Available
P2/MD2 - Positive Displacement 60Hz 50Hz	4.0 gpm @ 60 psig (15.1 lpm @ 4.1 bar) 3.3 gpm @ 60 psig (12.5 lpm @ 4.1 bar)	4.0 gpm @ 60 psig (15.1 lpm @ 4.1 bar) 3.3 gpm @ 60 psig (12.5 lpm @ 4.1 bar)
T1 - Turbine 60Hz* 50Hz*	3.5 gpm @ 60 psid (13.3 lpm @ 4.1 bar) 2.5 gpm @ 60 psid (9.5 lpm @ 4.1 bar)	3.5 gpm @ 60 psid (13.3 lpm @ 4.1 bar) 2.5 gpm @ 60 psid (9.5 lpm @ 4.1 bar)
P3 - Centrifugal Pump 60Hz* 50Hz*	10 gpm @ 32 psid (37.9 lpm @ 2.2 bar) 10 gpm @ 20 psid (37.9 lpm @ 1.4 bar)	10 gpm @ 32 psid (37.9 lpm @ 2.2 bar) 10 gpm @ 20 psid (37.9 lpm @ 1.4 bar)
P4 - Centrifugal Pump 60Hz* 50Hz*	15 gpm @ 57 psid (56.8 lpm @ 3.9 bar) 15 gpm @ 34 psid (56.8 lpm @ 2.3 bar)	15 gpm @ 57 psid (56.8 lpm @ 3.9 bar) 15 gpm @ 34 psid (56.8 lpm @ 2.3 bar)

\* Pumping capacity pressure values for turbine and centrifugal pumps are differential pressures between the inlet and the outlet of the chiller.

- Cooling capacity based on P2 pumps with no backpressure. Heat input from the pump will result in a reduction in cooling capacity. The cooling capacity reduction will vary based on the pump chosen as well as pump backpressure and flow.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Additional dimensions are at the end of this section, add 1/8" (3 cm) to height for SEMI.
- Add 30 pounds (14 kilograms) for global voltage chillers.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

# Specifications



	ThermoFlex7500	ThermoFlex10000
<b>Standard Chiller Process Fluid Temperature/Setpoint Range</b>	+5°C to +40°C +41°F to +104°F	+5°C to +40°C +41°F to +104°F
<b>High-Temperature Chiller Process Fluid Temperature/Setpoint Range</b>	+5°C to +90°C +41°F to +194°F	+5°C to +90°C +41°F to +194°F
<b>Ambient Temperature Range All Chillers</b>	+10°C to +40°C +50°F to +104°F	+10°C to +40°C +50°F to +104°F
<b>Temperature Stability</b>	± 0.1°C	± 0.1°C
<b>Cooling Capacity at 20°C</b>		
60 Hz	7500 W (25575 BTU)	10000 W (34100 BTU)
50 Hz	6425 W (21910 BTU)	8500 W (28985 BTU)
<b>Heater Size</b>	5.0kW at 208V/6.1kW at 230V 4.6kW at 400V/6.1kW at 460V	5.0kW at 208V/6.1kW at 230V 4.6kW at 400V/6.1kW at 460V
<b>Refrigerant</b>	R407C	R407C
<b>Reservoir Volume</b>		
Gallons	4.75	4.75
Liters	17.9	17.9
<b>Footprint or Dimensions (H x W x D)</b>		
Air-Cooled Inches	52.3 x 25.2 x 33.8	52.7 x 25.2 x 33.8
Centimeters	132.7 x 63.9 x 85.6	134.0 x 63.9 x 85.5
Water-Cooled Inches	45.9 x 25.2 x 33.8	46.4 x 25.2 x 33.8
Centimeters	116.6 x 63.9 x 85.6	117.8 x 63.9 x 85.6
<b>Weight P2/P3/P5 (empty)</b>		
Air-Cooled lb	356/372.5/405.5	356/372.5/405.5
kg	161.5/169/184	161.5/169/184
Water-Cooled lb	315/331.5/364.5	315/331.5/364.5
kg	143/150/165	143/150/165
<b>Pumping Capacity</b>		
P2/MD2 - Positive Displacement 60 Hz	4.0 gpm @ 60 psig (15.1 lpm @ 4.1 bar)	4.0 gpm @ 60 psig (15.1 lpm @ 4.1 bar)
50 Hz	3.3 gpm @ 60 psig (12.5 lpm @ 4.1 bar)	3.3 gpm @ 60 psig (12.5 lpm @ 4.1 bar)
P3 - Centrifugal Pump 60 Hz*	10 gpm @ 32 psid (37.9 lpm @ 2.2 bar)	10 gpm @ 32 psid (37.9 lpm @ 2.2 bar)
50 Hz*	10 gpm @ 20 psid (37.9 lpm @ 1.4 bar)	10 gpm @ 20 psid (37.9 lpm @ 1.4 bar)
P5 - Centrifugal Pump 60 Hz*	20 gpm @ 60 psid (75.7 lpm @ 4.1 bar)	20 gpm @ 60 psid (75.7 lpm @ 4.1 bar)
50 Hz*	20 gpm @ 35 psid (75.7 lpm @ 2.4 bar)	20 gpm @ 35 psid (75.7 lpm @ 2.4 bar)
T5 - Turbine Pump 60 Hz*	7.3 gpm @ 60 psid (27.6 lpm @ 4.1 bar)	7.3 gpm @ 60 psid (27.6 lpm @ 4.1 bar)
50 Hz*	6.2 gpm @ 60 psid (23.5 lpm @ 4.1 bar)	6.2 gpm @ 60 psid (23.5 lpm @ 4.1 bar)

\* Pumping capacity pressure values for centrifugal and turbine pumps are differential pressures between the inlet and the outlet of the chiller. P 5 pumping capacity in high temperature chillers is slightly lower, see P 5 Pumping Capacity curves in this Section.


- Cooling capacity based on P 2 pumps with no backpressure. Heat input from the pump will result in a reduction in cooling capacity. The cooling capacity reduction will vary based on the pump chosen as well as pump backpressure and flow.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Additional dimensions are at the end of this section.
- For global voltage chillers with a P2 pump add 30 pounds (14 kilograms). Add 10 pounds (4.5 kilograms) for chillers with a P3 or P5 pump.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

**Cooling Capacity**


- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage, on chillers with P2 pumps with no back pressure. Other fluids, fluid temperatures, ambient temperatures, altitude, operating voltages or pumps will affect performance. See Section 3
- Thermo Fisher Scientific reserves the right to change specifications without notice.


**Cooling Capacity**

- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage, on chillers with P2 pumps with no back pressure. Other fluids, fluid temperatures, ambient temperatures, altitude, operating voltages or pumps will affect performance. See Section 3.
- Chillers require the use of 50/50 EG/water or 50/50 PG/water below 10°C process temperature to prevent freezing/glazing of the plate exchanger.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

**Pumping Capacity  
Positive Displacement Pump P1/P2  
Magnetic Drive Pumps MD1/MD2**


**Pumping Capacity  
Turbine Pump T1/T5**


- Pump curves are nominal values. Pressure values for turbine pumps are differential pressures between the inlet and the outlet of the chiller.
- Pump performance results were obtained with no restrictions on the return to the system or with any options installed. For example, utilizing the DI option will result in a 0.5 gpm flow reduction.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

**Pumping Capacity  
Centrifugal Pump P3/P4**


**Pumping Capacity  
Centrifugal Pump P5**

- Pump curves are nominal values. Pressure values for centrifugal pumps are differential pressures between the inlet and the outlet of the chiller.
- Pump performance results were obtained with no restrictions on the return to the system or with any options installed. For example, utilizing the DI option will result in a 0.5 gpm flow reduction.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20° C process setpoint, 25° C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Thermo Fisher Scientific reserves the right to change specifications without notice.


**Pumping Capacity  
Turbine Pump T0**


- Pump curves are nominal values. Pressure values for turbine pumps are differential pressures between the inlet and the outlet of the chiller.
- Pump performance results were obtained with no restrictions on the return to the system or with any options installed. For example, utilizing the DI option will result in a 0.5 gpm flow reduction.
- Specifications obtained at sea level using water as the recirculating fluid, at a 20°C process setpoint, 25°C ambient condition, at nominal operating voltage. Other fluids, fluid temperatures, ambient temperatures, altitude or operating voltages will affect performance. See Section 3.
- Thermo Fisher Scientific reserves the right to change specifications without notice.

**ThermoFlex900/1400**  
**Dimensions**  
 (inches/centimeters)

Front View

Side View

**ThermoFlex900/1400**

Process discharge for chillers with optional flow transducer or Internal pressure regulator adjustment (Optional)  
1/2" FNPT Stainless Steel

Rear View


Process discharge fluid connection  
1/2" FNPT Cast Bronze


Process fluid return connection  
1/2" FNPT Stainless Steel

See Section 3 for additional plumbing information.

Top View

Shipping crate dimensions (approximate):

21" (53 cm) wide


35" (89 cm) tall

40" (102 cm) deep

• Thermo Fisher Scientific reserves the right to change specifications without notice.


**ThermoFlex2500**  
**Dimensions**  
 (inches/centimeters)

Front View

\* Add 1/8" (3 mm) for SEMI chillers , see Section 5.

Side View

**ThermoFlex2500**

Process discharge for chillers with optional flow transducer or Internal pressure regulator adjustment (Optional)  
1/2" FNPT Stainless Steel

Process discharge fluid connection  
1/2" FNPT Cast Bronze

Process fluid return connection  
1/2" FNPT Stainless Steel


Optional 1/2" Auxiliary and 1/4" Auto-Refill Ports

Rear View

See Section 3 for additional plumbing information.


Process fluid drain (1/4" FNPT) Riton

\* Add 1/8" (3 mm) for SEMI chillers, see Section 5.

Top View


Shipping crate dimensions (approximate):  
23" (58 cm) wide  
36" (91 cm) tall  
40" (102 cm) deep

• Thermo Fisher Scientific reserves the right to change specifications without notice.


**ThermoFlex3500/5000**  
**Dimensions**  
 (inches/centimeters)

Front View

\* Add 1/8" (3 mm) for SEMI chillers, see Section 5.

Side View

**ThermoFlex3500/5000**

Process discharge for chillers with optional flow transducer and P1, P2 & T1 pumps  
or  
Internal pressure regulator adjustment (Optional P1/MD1, P2/MD2 & T1 only)

1/2' FNPT Stainless Steel


Process discharge connection  
Cast Bronze

A P3, P4 pumps 3/4" FNPT  
B P1/MD1, P2/MD2, T1 pumps  
1/2" FNPT


Process return connection  
Stainless Steel  
P3, P4 pumps 3/4" FNPT  
P1/MD1, P2/MD2, T1 pumps  
1/2" FNPT

See Section 3 for additional plumbing information.

Rear View


\* Add 1/8' (3 mm) for SEMI chillers, see Section 5

1/4" FPT Process Drain  
Stainless Steel with  
1/4" Brass plug  
(P3, P4 pumps only)


Top View


• Thermo Fisher Scientific reserves the right to change specifications without notice.

**ThermoFlex7500/10000**  
**Dimensions**  
 (inches/centimeters)

Side View

Top View

Air-cooled shipping crate dimensions (approximate):  
 35 3/4" (91 cm) wide  
 61 1/2" (156 cm) tall  
 46 3/8" (118 cm) deep

The applicable options fit within this envelope, see Section 5

- Thermo Fisher Scientific reserves the right to change specifications without notice.

**ThermoFlex7500/10000**

Rear View  
(Air-Cooled)

Process Discharge   
 P2/MD2 = 1/2" FNPT  
 Cast Bronze  
 P3 - P5, T5= 1" FNPT  
 Wrought Copper


Process Return   
 Stainless Steel  
 P2/MD2 = 1/2" FNPT  
 P3 - P5, T5 = 1" FNPT

See Section 3 for  
 additional plumbing  
 information.



Process fluid drain (1/4" FNPT)  
 Stainless Steel with Brass plug or  
 a Riton connector

Rear View  
(Water-Cooled)

T5 flow control valve is used to adjust the flow rate. The handle is designed to identify the valve's position, from full flow to full bypass. When handle in the horizontal position (in line with the discharge line) the application is receiving full flow. With the handle is vertical the valve is in full bypass.

T 5 Flow Control - see above.

Facility water connections  
 Cast Bronze  
 3/4" FNPT

## Section 3 Installation

### Site Requirements

Ambient Temperature Range*	10°C to 40°C (50°F to 104°F)
Relative Humidity Range	10% to 80% (non-condensing)
Operating Altitude*	Sea Level to 6562 feet (2000 meters)
Overvoltage Category	II
Pollution Degree	2
Degree of Protection	IP 20

---

\*Because of the decrease in air density, maximum temperature for the air entering an air-cooled ThermoFlex is reduced by 1°C per 1,000 feet above sea level. In addition, cooling capacity is reduced 1.2% per 1,000 feet above sea level.

**Never place the chiller in a location where excessive heat, moisture, inadequate ventilation, or corrosive materials are present. ▲**

**Note** Refer to the nameplate information on the rear of the chiller. ▲

Air-cooled chillers retain their full rated capacity at 20°C setpoint in ambient temperatures up to 25°C (77°F). For ambient temperatures above 25°C please de-rate the cooling capacity 3% for every 1°C above 25°C (77°F), up to a maximum ambient temperature of 40°C (104°F). Note that when operating at a process temperature lower than 20°C the de-rate percentage may increase due to additional gains from losses to ambient.

**Note** Depending on the setpoint and ambient temperatures, there may be a heat gain or loss through the plumbing resulting in a variation from setpoint temperature at the application inlet. Applications with large temperature variations between ambient and setpoint temperatures, and/or long plumbing lengths, may require additional insulation. ▲

ThermoFlex2500 air-cooled chillers are equipped with a two-speed fan. Should the chiller's internal ambient temperature reach 50°C for 30 seconds, or reach 53°C, the fan speed switches from slow speed to high speed to maintain internal temperatures within acceptable limits. When the temperature reaches 44°C or below for at least 15 minutes the speed returns to low. When in high speed the chiller's decibel level increases significantly.

**Note** High speed is required for the chiller to achieve its 2500 watt cooling capacity. At high-end operating conditions the fan can be set to run at high speed all the time using the controller's **SETTINGS** display, see Section 4. ▲

Chillers installed below the end-user application may enable system fluid to drain back into the chiller and cause spillage. Thermo Fisher offers an anti-drainback kit to prevent any spillage, see Section 5.


Air-cooled chillers can be installed with both sides blocked, or one side and the rear. See Figure 3-1. The front needs a minimum clearance of 24". Air enters the front of the system and exits through the sides and rear.

Having two sides blocked can impact the chiller's performance due to changes in air flow. If your installation requires two blocked sides please ensure that the following requirements are met:

Process Setpoint Temperature: Below 30°C (86°F)

Ambient: Below 40°C (104°F)

Before operating the chiller in conditions outside any of those listed in this manual please contact Thermo Fisher Scientific's Sales, Service and Customer Support to review your installation.

**Figure 3-1** Minimum Clearance

## Electrical Requirements

The chiller's construction provides protection against the risk of electrical shock by grounding appropriate metal parts. The protection will not function unless the power cord is connected to a properly grounded outlet. It is the user's responsibility to assure a proper ground connection is provided. ▲

The chiller must be installed in accordance with the National Electrical Code and the with reference to the information on the chiller's nameplate located on the rear.

Locate the chiller so it is near, and has easy access to, its disconnecting device.

The user is responsible to ensure that the line cord provided meets local electrical codes. If not, contact qualified installation personnel.

The chiller is intended for use on a dedicated outlet. The ThermoFlex has an internal circuit protection that is equivalent (approximately) to the branch circuit rating. This is to protect the ThermoFlex, and is not intended as a substitute for branch circuit protection.

Electrical Service Requirements (Standard temperature chillers):

<b>ThermoFlex900</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	100 VAC	50Hz	1Ø	15A	5-15P
	115 VAC	60Hz	1Ø	15A	5-15P
	200 VAC	50Hz	1Ø	15A	6-15P
	208-230 VAC	60Hz	1Ø	15A	6-15P
	230 VAC	50Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-
<b>ThermoFlex1400</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	100 VAC	50Hz	1Ø	20A	5-20P
	115 VAC	60Hz	1Ø	20A	5-20P
	200 VAC	50Hz	1Ø	15A	6-15P
	208-230 VAC	60Hz	1Ø	15A	6-15P
	230 VAC	50Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-
<b>ThermoFlex2500</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC P 1, P 2 Pump	50 Hz	1Ø	15A	6-15P
	208-230 VAC P 1, P 2 Pump	60Hz	1Ø	15A	6-15P
	200 VAC T 1 Pump	50Hz	1Ø	20A	6-20P
	208-230 VAC T 1 Pump	60Hz	1Ø	20A	6-20P
	230 VAC	50 Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-

\* Refer to Appendix A for country specific ratings.

Continued on next page.

## Electrical Service Requirements (Standard temperature chillers):

<b>ThermoFlex3500/5000</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC P 1, P 2 Pump	50 Hz	1Ø	15A	6-15P
	200 VAC T 1, P 3 Pump	50 Hz	1Ø	20A	6-20P
	200 VAC P 4 Pump	50 Hz	1Ø	30A	L6-30P
	208-230 VAC P 1, P 2 Pump	60 Hz	1Ø	15A	6-15P
	208-230 VAC T 1, P 3 Pump	60 Hz	1Ø	20A	6-20P
	208-230 VAC P 4 Pump	60 Hz	1Ø	30A	L6-30P
	230 VAC P 1 - P 4 Pump	50 Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-

<b>ThermoFlex7500/10000 (Air-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	200 VAC P 2 Pump	50 Hz	3Ø	16.5	30	L15-20P
	200 VAC P 3 Pump	50 Hz	3Ø	18.7	30	L15-30P
	200 VAC P 5 Pump	50 Hz	3Ø	22.3	35	L15-30P
	200 VAC T 5 Pump	50 Hz	3Ø	17.3	30	L15-20P
	208-230 VAC P 2 Pump	60 Hz	3Ø	16.5	30	L15-20P
	208-230 VAC P 3 Pump	60 Hz	3Ø	18.7	30	L15-30P
	208-230 VAC P 5 Pump	60 Hz	3Ø	22.3	35	L15-30P
	208-230 VAC T 5 Pump	60 Hz	3Ø	17.3	30	L15-20P
	400 VAC P 2 Pump	50 Hz	3Ø	10.9	20	IEC309
	400 VAC P 3 Pump	50 Hz	3Ø	9.6	15	IEC309
	400 VAC P 5 Pump	50 Hz	3Ø	11.8	15	IEC309
	400 VAC T 5 Pump	50 Hz	3Ø	8.7	15	IEC309

<b>ThermoFlex7500/10000 (Water-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	200 VAC P 2 Pump	50 Hz	3Ø	16.2	30	L15-20P
	200 VAC P 3 Pump	50 Hz	3Ø	18.4	30	L15-30P
	200 VAC P 5 Pump	50 Hz	3Ø	22.0	35	L15-30P
	200 VAC T 5 Pump	50 Hz	3Ø	17.0	30	L15-20P
	208-230 VAC P 2 Pump	60 Hz	3Ø	16.2	30	L15-20P
	208-230 VAC P 3 Pump	60 Hz	3Ø	18.4	30	L15-30P
	208-230 VAC P 5 Pump	60 Hz	3Ø	22.0	35	L15-30P
	208-230 VAC T 5 Pump	60 Hz	3Ø	17.0	30	L15-20P
	400 VAC P 2 Pump	50 Hz	3Ø	10.6	20	IEC309
	400 VAC P 3 Pump	50 Hz	3Ø	9.3	15	IEC309
	400 VAC P 5 Pump	50 Hz	3Ø	11.5	20	IEC309
	400 VAC T 5 Pump	50 Hz	3Ø	8.4	15	IEC309

**MCA** = Minimum Current Ampacity

**MOPD** = Maximum Overcurrent Protective Device

Values reflect those on the nameplate located on the rear of the chiller.

Continued on next page.

## Electrical Service Requirements (Variable voltage chillers):

<b>ThermoFlex900</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	115 VAC	60Hz	1Ø	15A	5-15P*
	100 VAC	50/60Hz	1Ø	15A	5-15P*
<b>ThermoFlex1400</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	115 VAC	60Hz	1Ø	20A	-
	100 VAC	50/60Hz	1Ø	20A	-

\* United States and Japan only. All other plugs are country specific.

## Electrical Service Requirements (Global voltage chillers):

<b>ThermoFlex900</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200/208/230 VAC	60Hz	1Ø	15A	-
	200/230 VAC	50Hz	1Ø	**16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-
<b>ThermoFlex1400</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200/208/230 VAC	60Hz	1Ø	15A	-
	200/230 VAC	50Hz	1Ø	**16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-
<b>ThermoFlex2500</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC T 1 Pump	60Hz	1Ø	15A	-
	208-230 VAC T 1 Pump	60Hz	1Ø	20A	-
	230 VAC	50 Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-
<b>ThermoFlex3500/5000</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200/208-230 VAC P 1 P 3 Pump	50/60Hz	1Ø	15A	-
	200/208-230 VAC T 1 P 3 Pump	50/60Hz	1Ø	20A	-
	200/208-230 VAC P 4 Pump	50/60Hz	1Ø	30A	Hard wired

\*\* Chillers selected for 230 VAC operation have a range of -10% to +7%. Refer to Appendix A for country specific ratings.

For installation information on Variable Voltage and Global Voltage chillers refer to Appendix B. Refer to the nameplate label located on the rear of the chiller for specific electrical requirements.

Continued on next page.

## Electrical Service Requirements (Global voltage chillers):

<b>ThermoFlex7500/10000 (Air-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	400 VAC P 2 Pump	50 Hz	3Ø	8.8	15	Hard wire
	400 VAC P 3 Pump	50 Hz	3Ø	10.1	20	Hard wire
	400 VAC P 5 Pump	50 Hz	3Ø	12.3	20	Hard wire
	400 VAC T 5 Pump	50 Hz	3Ø	9.1	15	Hard wire
	460 VAC P 2 Pump	60 Hz	3Ø	8.8	15	Hard wire
	460 VAC P 3 Pump	60 Hz	3Ø	10.1	20	Hard wire
	460 VAC P 5 Pump	60 Hz	3Ø	12.3	20	Hard wire
	460 VAC T 5 Pump	60 Hz	3Ø	9.1	15	Hard wire

<b>ThermoFlex7500/10000 (Water-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	400 VAC P 2 Pump	50 Hz	3Ø	8.4	15	Hard wire
	400 VAC P 3 Pump	50 Hz	3Ø	9.7	20	Hard wire
	400 VAC P 5 Pump	50 Hz	3Ø	11.9	20	Hard wire
	400 VAC T 5 Pump	50 Hz	3Ø	8.8	15	Hard wire
	460 VAC P 2 Pump	60 Hz	3Ø	8.4	15	Hard wire
	460 VAC P 3 Pump	60 Hz	3Ø	9.7	20	Hard wire
	460 VAC P 5 Pump	60 Hz	3Ø	11.6	20	Hard wire
	460 VAC T 5 Pump	60 Hz	3Ø	8.8	15	Hard wire

## Electrical Service Requirements (High temperature chillers):

<b>ThermoFlex1400</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC 1KW	50 Hz	1Ø	20A	6-20P
	200-230 VAC 1KW	60 Hz	1Ø	20A	6-20P
	200 VAC 2KW	50 Hz	1Ø	30A	6-30P
	208-230 VAC 2KW	60 Hz	1Ø	30A	6-30P
	230 VAC 2KW	50 Hz	1Ø	32A	-
	230 VAC 1KW	50 Hz	1Ø	*16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	-

<b>ThermoFlex2500</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC P 1 P 2 Pump	50 Hz	1Ø	30A	6-30P
	208-230 VAC P 1 P 2 Pump	60 Hz	1Ø	30A	6-30P
	208-230 VAC T 1 Pump	60 Hz	1Ø	30A	6-30P
	200 VAC T 1 Pump	50 Hz	1Ø	40A	Hard wire
	230 VAC	50 Hz	1Ø	32A	-

\* Refer to Appendix A for country specific ratings.

Continued on next page.

## Electrical Service Requirements (High temperature chillers):

<b>ThermoFlex3500/5000</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>
	200 VAC P1 P2 Pump	50 Hz	1Ø	30A	6-30P
	200 VAC T1 P3 P4 Pump	50 Hz	1Ø	40A	Hard wire
	208-230 VAC P1 P2 Pump	60 Hz	1Ø	30A	6-30P
	208-230 VAC T1 P3 P4 Pump	60 Hz	1Ø	40A	Hard wire
	230 VAC P1- P4, T1 Pump	50 Hz	1Ø	32A	IEC309

<b>ThermoFlex7500/10000 (Air-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	200 VAC P2 Pump	50 Hz	3Ø	35.7	45	Hard wire
	200 VAC P3 Pump	50 Hz	3Ø	37.9	45	Hard wire
	200 VAC P5 Pump	50 Hz	3Ø	41.5	50	Hard wire
	200 VAC T5 Pump	50 Hz	3Ø	36.4	45	Hard wire
	208-230 VAC P2 Pump	60 Hz	3Ø	35.7	45	Hard wire
	208-230 VAC P3 Pump	60 Hz	3Ø	37.9	45	Hard wire
	208-230 VAC P5 Pump	60 Hz	3Ø	41.5	50	Hard wire
	208-230 VAC T5 Pump	60 Hz	3Ø	36.4	45	Hard wire
	400 VAC P2 Pump	50 Hz	3Ø	19.2	25	Hard wire
	400 VAC P3 Pump	50 Hz	3Ø	17.9	25	Hard wire
	400 VAC P5 Pump	50 Hz	3Ø	20.1	25	Hard wire
	400 VAC T5 Pump	50 Hz	3Ø	17.0	25	Hard wire

<b>ThermoFlex7500/10000 (Water-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	200 VAC P2 Pump	50 Hz	3Ø	35.3	45	L15-20P
	200 VAC P3 Pump	50 Hz	3Ø	37.5	45	L15-20P
	200 VAC P5 Pump	50 Hz	3Ø	41.1	50	L15-30P
	200 VAC T5 Pump	50 Hz	3Ø	36.1	45	L15-30P
	208-230 VAC P2 Pump	60 Hz	3Ø	35.3	45	L15-20P
	208-230 VAC P3 Pump	60 Hz	3Ø	37.5	45	L15-20P
	208-230 VAC P5 Pump	60 Hz	3Ø	41.1	50	L15-30P
	208-230 VAC T5 Pump	60 Hz	3Ø	36.1	45	L15-30P
	400 VAC P2 Pump	50 Hz	3Ø	18.8	25	IEC309
	400 VAC P3 Pump	50 Hz	3Ø	17.5	25	IEC309
	400 VAC P5 Pump	50 Hz	3Ø	19.7	25	IEC309
	400 VAC T5 Pump	50 Hz	3Ø	16.7	25	IEC309

**MCA** = Minimum Current Ampacity

**MOPD** = Maximum Overcurrent Protective Device

Values reflect those on the nameplate located on the rear of the chiller.

Continued on next page.

## Electrical Service Requirements (High temperature global voltage chillers):

<b>ThermoFlex1400</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>	
	200/208/230 VAC 1KW	60 Hz	1Ø	20A	6-20P	
	200/208/230 VAC 2KW	60 Hz	1Ø	30A	Hard wire	
	200/230 VAC 2KW	50 Hz	1Ø	32A	Hard wire	
	200/230 VAC 1KW	50 Hz	1Ø	**16A <sup>1</sup> , 15A <sup>2</sup> , 13A <sup>3</sup>	Hard wire	
<b>ThermoFlex2500</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>	
	200/208/230 VAC P 1 P 2 Pump	60 Hz	1Ø	30A	Hard wire	
	200/208/230 VAC T 1 Pump	60 Hz	1Ø	40A	Hard wire	
	200/230 VAC	50 Hz	1Ø	32A	Hard wire	
<b>ThermoFlex3500/5000</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>Branch Circuit Requirements</b>	<b>Line Cord Plug</b>	
	200/208-230 VAC P 1 P 2 Pump	50/60 Hz	1Ø	30, 32A	Hard wire	
	200/208-230 VAC T 1 P 3 P 4 Pump	50/60 Hz	1Ø	32, 40A	Hard wire	
<b>ThermoFlex7500/10000 (Air-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	400 VAC P 2 Pump	50 Hz	3Ø	18.3	25	Hard wire
	400 VAC P 3 Pump	50 Hz	3Ø	19.6	25	Hard wire
	400 VAC P 5 Pump	50 Hz	3Ø	21.8	30	Hard wire
	400 VAC T 5 Pump	50 Hz	3Ø	18.7	25	Hard wire
	460 VAC P 2 Pump	60 Hz	3Ø	18.3	25	Hard wire
	460 VAC P 3 Pump	60 Hz	3Ø	19.6	25	Hard wire
	460 VAC P 5 Pump	60 Hz	3Ø	21.8	30	Hard wire
	460 VAC T 5 Pump	60 Hz	3Ø	18.7	25	Hard wire
<b>ThermoFlex7500/10000 (Water-cooled)</b>	<b>Voltage ±10%</b>	<b>Frequency</b>	<b>Phase</b>	<b>MCA</b>	<b>MOPD</b>	<b>Line Cord Plug</b>
	400 VAC P 2 Pump	50 Hz	3Ø	18.0	25	Hard wire
	400 VAC P 3 Pump	50 Hz	3Ø	19.3	25	Hard wire
	400 VAC P 5 Pump	50 Hz	3Ø	21.5	30	Hard wire
	400 VAC T 5 Pump	50 Hz	3Ø	18.4	25	Hard wire
	460 VAC P 2 Pump	60 Hz	3Ø	18.0	25	Hard wire
	460 VAC P 3 Pump	60 Hz	3Ø	19.3	25	Hard wire
	460 VAC P 5 Pump	60 Hz	3Ø	21.5	30	Hard wire
	460 VAC T 5 Pump	60 Hz	3Ø	18.4	25	Hard wire

\*\* Chillers selected for 230 VAC operation have a range of -10% to +7%. Refer to Appendix A for country specific ratings.

For installation information Global Voltage chillers refer to Appendix B.


Refer to the nameplate on the chiller's rear for specific electrical requirements.

## Hard Wire Installation


For personal safety and equipment reliability, only a qualified technician should perform the following procedure. ▲

**Note** The technician is responsible for installing circuit protection for incoming power. Before wiring consult the nameplate on the rear of the chiller. Ensure installation is in accordance with the National Electrical Code and any other applicable country and local codes. ▲

**Figure 3-2** Electrical Box

- Remove the six screws securing the electrical box cover to the rear of the chiller.
- Remove the double knock out ( $\frac{7}{8}$ " and  $1\frac{3}{32}$ ").
- Insert the cable through the hole.
- Refer to the label in the electrical box to configure your chiller, see Figure 3-3.
- Secure the cable's ground wire to the ground stud.
- Reinstall the cover.

**Figure 3-3** Sample Label

## Plumbing Requirements

Ensure that all shipping plugs are removed before installation.

Never connect the process fluid lines to your facility water supply or any pressurized liquid source. ▲

To prevent damage to the chiller's plate exchanger, centrifugal pumps require a 4.0 gpm (15.1 lpm) minimum flow rate. ▲

P1 and P2 pumps are capable of producing 185 psig. Ensure your plumbing is rated to withstand this pressure at your operating temperature. An external pressure relief valve is available, see Section 5. ▲

**Note** Ensure your plumbing installation will develop a back pressure to the ThermoFlex greater than 3 PSIG. Lower pressure will shut down the chiller. ▲

The process fluid connections are located on the rear of the chiller and are labeled  (PROCESS OUTLET) and  (PROCESS INLET).

### Process Fluid Connections (FNPT)

#### Outlet


ThermoFlex900 - 10000	P 1 P 2 T 0 T 1	1/2" cast bronze
ThermoFlex3500 - 5000	P 3 P 4	3/4" cast bronze
ThermoFlex7500 - 10000	P 3 P 5 T 5	1" wrought copper
Inlet - Same size as outlet		all chillers stainless steel

### Supplied Adapters Standard Temperature Chillers


P 1 P 2 T 0 T 1	1/2" x 3/8" barb Polyethylene and 1/2" x 1/2" barb Nylon
P 3 P 4	3/4 MPT x 1/2 barb PVC
P 3 P 5 T 5	1" MPT x 1" barb PVC and 1" MPT x 3/4" barb PVC

The supplied adapters for high temperature chillers are brass.

See Section 2 for the specific locations on your chiller.

**Figure 3-4** Typical Plumbing Connections (1 of 2)


**Figure 3-4** Typical Plumbing Connections (2 of 2)

Connect the PROCESS OUTLET  to the fluid inlet on your application. Connect the PROCESS INLET  to the fluid outlet on your application. Ensure all connections are secure and that the proper sealant/lubricant for the fitting material is used. (If Teflon<sup>®</sup> tape is used, ensure the tape does not overhang the first thread as it could shred and get into the fluid.) Keep the distance between the chiller and the instrument being cooled as short as possible. Ensure tubing is straight and without bends. If diameter reductions are required, make them at the inlet and outlet of your application, not at the ThermoFlex.

### Water-cooled Chillers

For water-cooled chillers the facility water plumbing connections are also located on the rear and are labeled  FACILITY INLET and  FACILITY OUTLET. The connections are 1/2" Female NPT for ThermoFlex900 - 5000, 3/4" Female NPT for ThermoFlex7500 - 10000. Both connections for ThermoFlex900 to 10000 are cast bronze.

Connect the  FACILITY INLET to your facility water supply. Connect the  FACILITY OUTLET to your facility water return or drain. Ensure all connections are secure and that the proper sealant/lubricant for the fitting material is used. (If Teflon<sup>®</sup> tape is used, ensure the tape does not overhang the first thread as it could shred and get into the fluid.)

See Section 2 for the specific locations on your chiller.

**Figure 3-5** Typical Plumbing Connections, Water-cooled Chillers

## Process Fluid Requirements

Use of any fluid not listed below will void the manufacturer’s warranty. ▲

Standard Temp Chillers	High Temp Chillers
Filtered/Single Distilled water  0 - 75% Laboratory Grade Ethylene Glycol/Water  0 - 75% Laboratory Grade Propylene Glycol/Water  Deionized water (1 - 3 MΩ- cm, compensated)* *For applications requiring resistivity greater than 1 MΩ- cm please call and speak to an applications engineer for additional information.	Filtered water*  0 - 50% Laboratory Grade Ethylene Glycol/Water*  0 - 50% Laboratory Grade Propylene Glycol/Water*  *to 88°C for chillers with P1/P2 pumps, to 90°C for chillers with other pumps (when fluid properly maintained)

Ethylene glycol (EG) is poisonous and flammable. Before using any fluid or performing maintenance where contact with the fluid is likely refer to the manufacturer’s SDS for handling precautions. ▲

EG is also hygroscopic, it absorbs water from its environment. This can affect the freezing point and boiling point of the fluid over time and may result in system failure. ▲

To prevent freezing/glazing of the plate exchanger, all ThermoFlex7500 and 10000s require the use of 50/50 EG/water or 50/50 PG/water below 10°C process temperature. ▲

When using a process fluid mixture of ethylene glycol and water or propylene glycol and water, check the fluid concentration and pH on a regular basis. Changes in concentration and pH can impact system performance. ▲

When using EG/water or PG/water, top-off with plain water. After top-off check the fluid concentration. ▲

Do not use a Deionization (DI) filter cartridge with Inhibited EG or Inhibited PG. A DI filter will remove inhibitors from the solution rendering the fluid ineffective against corrosion protection. Also, inhibitors increase fluid conductivity. ▲

## Compatibility with Approved Fluids

### Filtered Tap Water/Single Distilled Water

Filtered drinking water and single distilled water are good choices for use in a recirculating chiller because the filtering/distilling process used removes microorganisms that could create biological fouling as well as harmful particulates and excessive minerals that could cause harmful deposits and scaling.

### Chlorine

Short term usage of tap water may not cause any adverse affects on the chiller or your application, but in the long term problems may arise. To help alleviate these problems Thermo Fisher Scientific recommends the use of chlorine.

The duration of time that chlorine remains in solution depends on factors such as water temperature, pH and availability of direct sunlight. We recommend maintaining chlorine levels at proper levels using chlorine test strips, generally 1 to 5 ppm is adequate.

For best results, maintain the pH of the fluid between 6.5 and 7.5. Do not add additional chlorine without first determining the concentration ratio that already exists in the fluid supply. Corrosion and degradation of the circulation components can result from concentration ratios that are too high. Contact our customer support for additional information.

### Deionized Water

Deionized water is water that has had its mineral ions removed using ion exchange resins. The purpose of this process is to remove the ions that allow electrical current to flow more easily through water. This helps to prevent electrical leaks to ground through the recirculating fluid. Deionized water is classified by the electrical resistance of the water, usually measured in MΩ-cm, with pure water having a resistance of 18 MΩ-cm.

Deionized water is in an unbalanced state and will leach the missing ions from the materials it comes in contact with. The aggressive nature of this leaching can cause pitting on metal surfaces. Note that the deionizing process does not remove microorganisms. Because of this, we recommend deionized water only with applications that have it as a specified requirement.

In any case, only deionized water with a 1 - 3 MΩ-cm resistivity is approved for use in Thermo Fisher Scientific recirculating chillers. For applications requiring resistivity greater than 1 MΩ- cm please call and speak to an applications engineer for additional information.

### Recommended Biocides and Inhibitors

Thermo Fisher Scientific offers a biocide and inhibitor package Thermo 200 (NALCO) premixed with 5 gallons of water or as a kit to be added to water. No other biocide or inhibitor is recommended for use in our chillers.

**Biocides are corrosive and can cause irreversible eye damage and skin burns. They are harmful if inhaled, swallowed or absorbed through the skin. Refer to the manufacturer's most current SDS. ▲**

To prevent freezing/glazing of the plate exchanger, ThermoFlex7500 and 10000 chillers require the use of 50/50 EG/water or 50/50 PG/water below 10°C process temperature. ▲

#### **Uninhibited Ethylene Glycol/Water**

Ethylene glycol is used to depress the freezing point of water and should only be used at temperatures where freeze point suppression is required. Ethylene glycol does not improve heat transfer and is not recommended for use as a biocide. Because glycols lower the surface tension of water and do not evaporate as readily as water, they may cause visible weepage past the pump seals. If weepage cannot be tolerated use seal-less, magnetically driven pumps where available.

Uninhibited simply means that the glycol does not contain any additives to prevent corrosion. While uninhibited ethylene glycol is acceptable for use, the pH level must be closely monitored and the fluid may need to be replaced more often. Since all glycols produce acids in the presence of air and the fluid change the glycol if the pH falls below 8. Note that litmus paper will not work to test the pH of ethylene glycol/water.

#### **Inhibited Ethylene Glycol/Water**

Inhibited glycol can help protect the wetted metals within the cooling circuit from corrosion caused by poor water quality, ethylene glycol oxidation (low pH) and mixed metals (electrolysis). The inhibitor works by either leaving a barrier coating on metal surfaces to buffer them from the corrosive fluid or by creating an oxidized layer that protects the underlying metal (passivating).

Inhibited automotive glycols are never acceptable. They use either silicates or Organic Acid Technology (OAT) as the inhibitor and these components are not compatible with the polymers used in recirculating chillers including the pump seals and internal hoses.

Inhibitors may also accelerate pump seal wear and seal-less, use magnetically driven pumps where available.

#### **Uninhibited Propylene Glycol/Water**

Propylene glycol does not transfer heat as well as ethylene glycol, but can be used when freeze point suppression is required as well as lower toxicity.

Propylene glycol does not function as a biocide and the pH needs to be maintained the same as with ethylene glycol as it also produces acid when oxidized.

#### **Inhibited Propylene Glycol/Water**

Inhibited propylene glycol has the same properties as uninhibited propylene glycol and the same concerns as inhibited ethylene glycol.

## Additional Fluid Information

When using the ThermoFlex chiller to circulate through aluminum, use a compatible corrosion inhibitor to prevent galvanic corrosion.

Ensure fluid viscosity is 50 cSt or less at the lowest temperature used.

Visible pump weepage may occur when compatible glycols, oils or other additives are used. Pump weepage is considered as a normal operating condition of mechanical seal pumps.

## Process Water Quality and Standards

Process Fluid	Permissible (PPM)	Desirable (PPM)
<b>Microbiologicals</b> (algae, bacteria, fungi)	0	0
<b>Inorganic Chemicals</b>		
Calcium	<25	<0.6
Chloride	<25	<10
Copper	<1.3	<1.0
	0.020 ppm if fluid in contact with aluminum	
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates/Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<25	<1
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
<b>Other Parameters</b>		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*

\* M<sup>-1</sup>cm (compensated to 25°C)

Unfavorably high total ionized solids (TIS) can accelerate the rate of galvanic corrosion. These contaminants can function as electrolytes which increase the potential for galvanic cell corrosion and lead to localized corrosion such as pitting. Pitting may become so extensive that it causes a breach in the cooling system resulting in a refrigerant leak into the fluid system.

As an example, raw water in the United States averages 171 ppm (of NaCl). The recommended level for use in a water system is between 0.5 to 5.0 ppm (of NaCl).

Recommendation: Initially fill the reservoir with distilled or 1-3 M<sup>-1</sup>cm deionized water. (It is acceptable to have the fluid drop to the other levels over-time.) Do not use untreated tap water as the total ionized solids level may be too high. This will reduce the electrolytic potential of the water and prevent or reduce the galvanic corrosion observed.

## Facility Water Quality and Standards (water-cooled chillers)

Facility Water	Permissible (PPM)	Desirable (PPM)
<b>Microbiologicals</b> (algae, bacteria, fungi)		
	0	0
<b>Inorganic Chemicals</b>		
Calcium	<40	<0.6
Chloride	<250	<25
Copper	<1.3 0.020 ppm if fluid in contact with aluminum	<1.0
Iron	<0.3	<0.1
Lead	<0.015	0
Magnesium	<12	<0.1
Manganese	<0.05	<0.03
Nitrates\Nitrites	<10 as N	0
Potassium	<20	<0.3
Silicate	<25	<1.0
Sodium	<20	<0.3
Sulfate	<250	<50
Hardness	<17	<0.05
Total Dissolved Solids	<50	<10
<b>Other Parameters</b>		
pH	6.5-8.5	7-8
Resistivity	0.01*	0.05-0.1*

\* MΩ-cm (compensated to 25°C)

**Note** A corrosion inhibitor is recommended if mixed metals are in the facility water loop. ▲

# Facility Water Requirements (water-cooled chillers)


Facility Water Maximum Inlet Pressure must not exceed 150 PSIG.

Facility Water Maximum Pressure Differential must not exceed 50 PSID.  
(Pressure Differential = Inlet Pressure - Outlet Pressure)

Contact us before using facility water that is above 35°C.


The facility water must meet the following conditions for the chillers to maintain their full rated capacity.

**ThermoFlex1400**

Example: See next page.


**ThermoFlex2500**

Example: See next page.


**ThermoFlex3500/5000**

The facility water must meet the following conditions for the chillers to maintain their full rated capacity.

Example: See below.

**ThermoFlex7500/10000**

Example:

Follow the ..... lines.

Start with a known, e.g., facility water temperature.

A - go across to temperature curve

B - drop down to determine the minimum required facility flow.

C - Where B crosses the PSID curve, go across to determine the minimum required PSID.

## Fluid Filter Bag

The reservoir has a fluid bag filter designed to prevent the introduction of particulates into the system.

**Install the filter bag before starting the chiller. ▲**

Place your fingers below the front of the housing and push up on the housing to remove it.

Install the bag.

Replace the housing. Slide its back edge under the lip of the chiller's top panel and then press down until the housing snaps into place.

**Figure 3-6** Fluid Filter Bag

## Priming

If able, pre-fill the process fluid lines. The chiller is designed to shut down if not properly primed.

Ensure that there is enough fluid prepared to fill both the chiller and your application. If able, pre-fill the process fluid lines to reduce the setup time.

Fill the reservoir to the max fill line on the reservoir sight tube. To prevent the introduction of particulates into the system, fill the chiller with the reservoir bag filter in place.

Start the chiller using the power button  on the control panel.

As the pump fills your application add fluid to the reservoir to maintain the fluid level.

Repeat this process until the fluid level no longer drops in the reservoir.

If you need to pause priming use the power button to turn the chiller off.

**Not completely filling the chiller and process fluid lines could damage the chiller's pump. ▲**

## Initial Filling Requirements

Ensure the reservoir drain plug on the back of the chiller is in place, or the Riton fitting is closed, and that all plumbing connections are secure.

**Before using any fluid refer to the manufacturer's SDS for handling precautions. ▲**

Locate and remove the reservoir cap by unscrewing it counterclockwise.


To prevent the introduction of particulates into the system, fill with the reservoir bag filter in place.

**Figure 3-7** Reservoir Cap

The reservoir has a sight tube and ball for easy fluid level monitoring. *Slowly* fill the reservoir with clean process fluid through the funnel only, failure to comply may result in internal spillage.

**Note** Filling the reservoir above MAX LEVEL fill line will result in a over flow error (**Over Flow**) causing the chiller to shut down. Also, fluids expand when heated. ▲

Since the reservoir capacity may be small compared to your application and air may need to be purged from the lines, have extra cooling fluid on hand to keep the system topped off when external circulation is started.

**Figure 3-8** Reservoir Sight Tube & Ball

**Before replacing the reservoir cap ensure the reservoir sight tube ball stopper is securely in place. ▲**

Replace the reservoir cap by screwing it clockwise. Cap should be hand tight.

## Fluid Top Off

Ensure the reservoir cap is at a safe handling temperature before removing. ▲

Remove the reservoir cap by unscrewing it counterclockwise.

To prevent the introduction of particulates into the system, fill the chiller with the reservoir bag filter in place.


The reservoir has a sight tube and ball for easy fluid level monitoring. *Slowly* fill the reservoir with clean process fluid through the funnel only, failure to comply may result in internal spillage.

**Note** Filling the reservoir above MAX LEVEL fill line will result in a over flow error (**Over Flow**) causing the chiller to shut down. Also, fluids expand when heated. ▲

**Note** Adding fluid that has a temperature differential with the fluid already in the reservoir will temporarily affect the chiller's stability performance. ▲

Before replacing the reservoir cap ensure the reservoir sight tube ball stopper is securely in place. ▲


**Figure 3-9** Reservoir Ball Stopper

## Section 4 Operation

### Deluxe Controller

The controller changes temperature using a Proportional-Integral-Derivative (PID) algorithm, and is designed with easy to use operator interface.

**Figure 4-1** Deluxe Controller

Press this button to start and stop the chiller.

Press this button to navigate through the controller displays and to increase adjustable values.

Press this button to navigate through the displays and to decrease adjustable values.

Press these buttons to navigate through the displays.

This button has four functions. Pressing it once allows changes to be made to the display, pressing it again saves the change and allows you to continue to other displays. It is also pressed to clear messages.

Depressing and holding the enter button for two seconds before starting the chiller allows you to view and to make changes to the controller settings.

esc

Press this button to abort any changes and at the same time return the controller to its previous display. Aborting a change can only be made *before* the change is saved.

## Initial Start Up Only

- Place the circuit protector located on the rear of the chiller to the up position.
- Press the  button on the controller.

The controller will display the quick start screen that will allow you to configure the chiller. Refer to the Quick Start Guide supplied with the chiller or the copy located after the Table of Contents.

Press the arrows to scroll through and highlight each line of the display. If a change is needed from the factory default value press the enter button.

**If the auto restart is enabled and the chiller shuts down as a result of a power failure, when power is restored the chiller will automatically restart and operate at the saved values. Consider any possible risks before enabling this mode of operation. ▲**

**Note** The chiller will not automatically restart if the **Show Quick Start** was enabled, see page 4-20. In this case only the controller is powered. ▲

If satisfied with all the entries, press enter when the last line - **Quick Start Done** - is highlighted. If not satisfied, press the left arrow, or esc button, to leave the quick start display.

In either case the screen will go blank.

If enter was pressed, press  to bring up the **Daily Start Up** displays and start the chiller. See next page.

## Daily Start Up

Before starting the chiller, double check all electrical and plumbing connections. Have extra recirculating fluid on hand. If the chiller will not start refer to Section 7 Troubleshooting.


Do not run the chiller until fluid is added.

Ensure the chiller's castors are locked. ▲

- Ensure the circuit protector located on the rear of the chiller is in the up position.
- Press  on the controller, the pump and the refrigeration system will automatically start.

The controller automatically sequences through the following displays:

**Figure 4-2** Daily Start Up Displays

If **Level %** is not displayed on the bottom of the screen refer to Status Displays in this section.

If the fluid level is above the required minimum value, see next page, the pump and compressor will automatically start but only after a 30 second delay.

If the fluid level is at or below the value, see next page, fluid will need to be added before the compressor will start, see **Initial Fill Mode** on next page. After start up, check the plumbing connections for leaks.

## Initial Fill Mode

- Press  on the controller.

If the **Initial Fill Mode** message appears the chiller needs fluid before it will start.

Required minimum level values depend on the chiller's heater size, see below.

If the **Level** is below the required amount **Add Fluid** will be displayed and flash until the required level is achieved, then **Add Fluid** will extinguish.


For all chillers, the pump will not start until the level is above **15%**. Once the level is above the required amount for five minutes, as the pump purges the system, the chiller will start.

If desired, press  to exit the initial fill mode and start the compressor.

**Note** The chiller can be turned ON/OFF using the ACOM and DCOM option in the **Initial Fill Mode**.

Heater	Required Minimum Level
None	15%
1 kW:	58%
2.3 kW:	64 or 58% (depending on chiller configuration)
5.0 kW:	87%
4.6 kW:	87%

## Process Fluid Temperature Display

**Figure 4-3** Process Fluid Temperature Display

The **SETPOINT** and **MENU** portions on the top of the display are used to view and/or change the controller's settings. They are explained in detail later in this Section.

The bars on the left of the display indicate the chiller's operating status.

-  Indicates the chiller is either not running or the ThermoFlex is at the desired setpoint.
-  Indicates the chiller is heating up to the setpoint.
-  Indicates the chiller is cooling down to the setpoint.

The small circle indicates which sensor is controlling the chiller.

-  When displayed, indicates the ThermoFlex is responding to the chiller's internal temperature sensor.
-  When displayed, indicates the ThermoFlex is responding to an external temperature sensor. See **Analog Comm - ACOM** in this Section to enable the external sensor.

The chiller is always in internal or external control. The default setting is internal.

The temperature displayed is the current process fluid temperature.


The message on the bottom is used to display the chiller's operational status. The status displays are explained later in this Section.


# SETPOINT

**Note** Depressing and holding the enter button for two seconds before starting the chiller allows you to view, as well as make changes to, the setpoint. ▲

Press the arrows to highlight **SETPOINT**.


Press enter and the setpoint value appears. If needed, press enter again and the value starts flashing indicating the setpoint can be changed.

The setpoint range for standard temperature chillers is +5°C to +40°C, the range for high temperature chillers is +5°C to +90°C. Press the up and down arrows to change the setpoint to the desired value.

Once the desired value is displayed press the enter button again to accept the change and to stop the flashing.


Return to the Process Fluid Temperature Display by pressing the left arrow or esc button.

## Status Displays

The controller can show up to four different messages. Press the down arrow to highlight the bottom of the display.

Press enter to get the displayed message to flash.

Press the arrows to scroll through the available displays. Once the desired display is showing press enter again.

### Available Displays:

**Pressure:** Indicates pump discharge pressure

**Level %:** Indicates reservoir fluid level

**Flow:** Indicates the process fluid flow rate (Optional)

**Resistivity:** Indicates the process fluid resistivity level (Optional)


The arrow at the end of the display indicates which direction you can scroll to get to another display. A ▼ means you can only scroll down to the next display. A ▲ means you can only scroll up to the next display. ▲▼ means you can scroll up or down.

# MENU Displays


The ThermoFlex controller uses menus to view/change the controller's settings.

**Note** Depressing and holding the enter button for two seconds before starting the chiller allows you to view and to make changes to the settings. ▲

1. Press the arrows to highlight MENU.


2. The controller brings up the Main Menu Display.


Main Menu Display

3. Press the up and down arrows to highlight the desired option.



4. With the desired option highlighted, press enter to display sub-menus.



5. There are various ways to navigate through all the menu options.

6. In the example in step 5, pressing the left arrow, or the esc button, on the controller will return you to the previous display.


7. Or, highlight the home icon at the bottom of the display.

8. Then press enter to return to the Main Menu Display.

## Main Menu Tree

**Note** The controller can only display six lines at a time. Press the up and down arrows to scroll through and highlight each line. ▲

**SETTINGS** allows you to view/change settings, see page.

**DRAINING** opens/closes the optional anti drainback valve, see page 4-20.


**TUNING** allows you to view/change the controller's **PID** values, see Section 7.

**CALIBRATE** calibrates the chiller's temperature, pressure and optional fluid flow sensors, see Section 8.

**DIAGNOSTICS** is used to display alarm messages, see Section 7, and also to display the chiller, filter and optional DI cartridge operating hours, see Section 6. **MONITOR VALUES** and **PASSWORD** are used only by qualified technicians.

**FIRMWARE VERSION** is used to display the controller's firmware information. Press the esc button to return to the Main Menu Display.

# SETTINGS Menu Tree

## SETTINGS

**SETTINGS** allows you to view/change settings.

**Note** The controller can only display six lines at a time. Press the up and down arrows to scroll through and highlight each line. ▲

If a change to a setting is needed, highlight the desired line and then press the enter button. The highlight will start to flash.

If the text on a line is all capital letters, e.g., **UNITS** and **TEMPERATURE**, the setting has a sub-menu. Pressing enter will bring up the sub-menu. The sub-menus, shown on the following pages, allow you to view/change the applicable settings.

Lines that are not all capital, e.g., **Setpoint** and **Line Frequency**, indicate the changes can be made directly on the **SETTINGS** display.

If the line has a box, **Auto Restart**, pressing enter with that line highlighted will toggle between on or off. See **Auto Restart** in this section for precautions and additional details.

SETTINGS	
◀	<b>UNITS</b>
<b>Setpoint</b>	<b>xx</b>
<b>TEMPERATURE</b>	
<b>PRESSURE</b>	
	
-----	
<b>FLUID LEVEL</b>	
<b>AUTO REFILL*</b>	
<b>FLOW*</b>	
<b>Line Frequency*</b>	<b>60HZ</b>
<input type="checkbox"/> <b>Auto Restart</b>	
<input type="checkbox"/> <b>Audible Alarms</b>	
<b>RA FAN SPEED MODE*</b>	
<b>Care Level</b>	<b>off</b>
<b>SERIAL COMM - DCOM*</b>	
<b>ANALOG COMM - ACOM*</b>	
<b>RESISTIVITY*</b>	
<b>RESTORE DEFAULTS</b>	
<b>mm/dd/yy</b>	<b>hh:mm:ss</b>
<input type="checkbox"/> <b>Show Quick Start</b>	


\*Only displays on chillers equipped with the option.

# UNITS

**UNITS** is used to view/change the controller's temperature, pressure and fluid flow (only chillers with an optional flow transducer) scales.

If a change to the setting is needed, highlight the desired scale and then press the enter button. The highlight will flash. Press the up or down arrow button to bring up the desired scale. Once the desired scale is displayed press enter to accept the change and stop the flashing.

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.


Scales: °C or °F  
 PSI, Bar or kPa  
 GPM or LPM

Defaults: °C  
 PSI  
 GPM

# Setpoint

**Setpoint** can be used to view/change the controller's setpoint temperature.

If a change to the setpoint is needed, highlight **Setpoint** and then press the enter button. The highlight will flash. Press the up or down arrow to bring up the desired setpoint value. Once the desired value is displayed press enter to accept the change and stop the flashing.

**Setpoint Range:**  
 +5°C to +40°C for standard temperature chillers  
 +5°C to +90°C for high temperature chillers,  
 +5°C to +88°C for P1 and P2 high temp chillers

**Setpoint Default:**  
 +20°C for all chillers

## TEMPERATURE

**TEMPERATURE** is used to view/change the chiller's high and low temperature fault and warning settings. If the chiller exceeds the fault setting it will shut down, the controller will display a fault message and, if enabled, sound the alarm. If the chiller exceeds the warning setting the chiller will continue to run, the controller will display a warning message and, if enabled, sound the alarm.

If a change is needed, press the arrows to highlight the desired setting and then press the enter button. The highlight will flash. Press the up or down arrow to change the value. Once the desired setting is displayed press enter to accept the change and stop the flashing.

**Note** You cannot set the **High Warning** value higher than the **High Fault** value. You cannot set the **Low Warning** value lower than the **Low Fault** value. ▲

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

TEMPERATURE	
◀High Fault	xxC
High Warning	xxC
Low Warning	xxC
Low Fault	xxC

### High Fault/Warning Range:

**+2°C to +43°C** for standard temperature chillers

**+2°C to +93°C** for high temperature chillers

### High Fault/Warning Default:

**+42°C** for standard temperature chillers

**+92°C** for high temperature chillers

### Low Fault/Warning Range:

**+2°C to +43°C** for all chillers

### Low Fault/Warning Default:

**+3°C** for all chillers

# PRESSURE

**PRESSURE** is used to view/change the chiller's high and low pressure fault and warning settings and set a delay time. If the chiller exceeds the fault setting it will shut down, the controller will display a fault message and, if enabled, sound the alarm. If the chiller exceeds the warning setting the chiller will continue to run, the controller will display a warning message and, if enabled, sound the alarm. The delay sets the length of time needed after a pressure fault before the chiller shuts down.

**Note** Since the controller can only display six lines at a time, keep pressing the down arrow until the **High P Delay** and **Low P Delay** lines are visible. ▲

If a change to the setting is needed, press the arrows to highlight the desired pressure and then press the enter button. The highlight will flash. Press the up or down arrow to change the value. Once the desired setting is displayed press enter to accept the change and stop the flashing.

**Note** You cannot set the **High Warning** value higher than the **High Fault** value. You cannot set the **Low Warning** value lower than the **Low Fault** value. ▲

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

PRESSURE	
◀ High Fault	xxpsi
High Warning	xxpsi
Low Warning	xxpsi
Low Fault	xxpsi
High P Delay	xsec
Low P Delay	xsec

### High/Low Fault/Warning Range:

Pump dependent - see below

### High/Low Pressure Time Delay Range:

0 to 30 seconds (0 to 60 for P3 - P5 pumps)

### High Time Pressure Delay Default:

0 seconds (60 seconds for P3 - P5 pumps)

### Low Time Pressure Delay Default:

10 seconds

Pump	Fault Range	High Default	Low Default
P1 P2 T0 T1:	3 to 105 PSI	105 PSI	3 PSI
P3 60 Hz:	3 to 46 PSI	46 PSI	3 PSI
P3 50 Hz:	3 to 32 PSI	32 PSI	3 PSI
P4 60 Hz:	3 to 85 PSI	85 PSI	3 PSI
P4 50 Hz:	3 to 60 PSI	60 PSI	3 PSI
P5 60 Hz:	3 to 87 PSI	87 PSI	3 PSI
P5 50 Hz:	3 to 56 PSI	56 PSI	3 PSI
T5 60 Hz:	2 to 105 PSI	105 PSI	2 PSI
T5 50 Hz:	2 to 105 PSI	105 PSI	2 PSI

Pump	Warning Range	High Default	Low Default
P1 P2 T0 T1:	4 to 100 PSI	100 PSI	4 PSI
P3 60 Hz:	4 to 46 PSI	46 PSI	4 PSI
P3 50 Hz:	4 to 32 PSI	32 PSI	4 PSI
P4 60 Hz:	4 to 85 PSI	85 PSI	4 PSI
P4 50 Hz:	4 to 60 PSI	60 PSI	4 PSI
P5 60 Hz:	4 to 87 PSI	87 PSI	4 PSI
P5 50 Hz:	4 to 56 PSI	56 PSI	4 PSI
T5 60 Hz:	4 to 105 PSI	105 PSI	4 PSI
T5 50 Hz:	4 to 105 PSI	105 PSI	4 PSI

## FLUID LEVEL

**FLUID LEVEL** is used to view/change the chiller's reservoir level fault and warning settings. If the chiller exceeds the fault setting it will shut down, the controller will display a fault message and, if enabled, sound the alarm. If the chiller exceeds the warning setting the chiller will continue to run, the controller will display a warning message and, if enabled, sound the alarm.

If a change to the setting is needed, press the arrows to highlight the desired line and then press the enter button. The highlight will flash. Press the up or down arrow to change the value. Once the desired setting is displayed press enter to accept the change and stop the flashing.

**Note** You cannot set the **Low Warning** value lower than the **Low Fault** value. ▲

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

FLUID LEVEL	
◀ Low Warning	x%
Low Fault	x%

**Low Warning/Fault Range/Default:**  
Heater dependent - see below

Heater	Warning Range	Default
None	6 - 100%	6%
1 kW:	58 - 100%	58%
2.3 kW:	58 or 69 - 100%	58 or 69%
5.0 kW:	87 - 100%	87%
4.6 kW:	87 - 100%	87%
Heater	Fault Range	Default
None	0 - 100%	0%
1 kW:	52 - 100%	52%
2.3 kW:	63 or 52 - 100%	63 or 52%
5.0 kW:	81 - 100%	81%
4.6 kW:	81 - 100%	81%

## AUTO REFILL (Optional)

**AUTO REFILL** is used to view/change the chiller's optional auto refill settings. **On** is the % of fluid level in the reservoir needed to turn the option on. **Time Out** is the maximum time the option will operate. Setting the time to **0** disables the option.

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

AUTO REFILL	
◀ On	x%
Time Out	xxsec

**On % Default:** Heater dependent - see below

Heater	Range	Default
None	70 - 100%	70%
1 kW:	70 - 100%	70%
2.3 kW:	69 or 70 - 100%	69 or 70%
5.0 kW:	87 - 100%	87%
4.6 kW:	87 - 100%	87%

**Time Out Range: 0 - 900** seconds

**Time Out Default: 45** seconds

**80** seconds for ThermoFlex7500/10000

Setting the time to **0** disables the option.


## FLOW (Optional)

**FLOW** is used to view/change the chiller's process fluid fault and warning settings. If the chiller exceeds the fault setting it will shut down, the controller will display a fault message and, if enabled, sound the alarm. If the chiller exceeds the warning setting the chiller will continue to run, the controller will display a warning message and, if enabled, sound the alarm.

If a change to the setting is needed, press the arrows to highlight the desired line and then press the enter button. The highlight will flash. Press the up or down arrow button to change the value. Once the desired setting is displayed press enter to accept the change and stop the flashing.

**Note** You cannot set the **High Warning** value higher than the **High Fault** value. You cannot set the **Low Warning** value lower than the **Low Fault** value. ▲

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

### High/Low Flow Fault/Warning Range/Defaults:

Pump	Range	High/Low Defaults
P1 P2 T1 T0:	0.0 to 10.5 GPM	0.0 GPM
P3 P4 P5:	0.0 to 30.0 GPM	0.0 GPM
T5	0.0 to 15.0 GPM	0.0 GPM

This feature is not enabled until the value is changed to something other than 0.0. If the feature is not enabled and the low flow rate drops below the flow rate listed below the chiller will continue to run and the controller, if displaying flow, will display **FLow: Low FLOW**.


P1, T0 ,T1 and T5 Pump 0.3 GPM

P2 Pump 1.0 GPM

P3, P4 and P5 Pump 4.0 GPM

## Line Frequency (Optional)

**Line Frequency** is used to identify the incoming frequency for chillers with a P3 - P5 pump and the capability to run on either **50 Hz** or **60 Hz**. The selected frequency automatically adjusts the firmware's *fixed* high pressure default setting.

Frequency Range: **50 Hz** or **60 Hz**

Frequency Default: **60 Hz**

## Auto Restart

**Auto Restart** is used to turn the auto restart feature on/off. The factory default setting is disabled but it may have been enabled using the quick start procedure.

If the auto restart is enabled and the chiller shuts down as a result of a power failure, when power is restored the chiller will automatically restart and operate at the saved values. Consider any possible risks before enabling this mode of operation. ▲

**Note** The chiller will not automatically restart if the **Show Quick Start** was enabled, see page 4-20. In this case only the controller is powered. ▲

Line Frequency	60Hz
<input type="checkbox"/> Auto Restart	
<input type="checkbox"/> Audible Alarms	
RA FAN SPEED MODE	
Care Level	off
SERIAL COMM - DCOM	

The color of the box next to **Auto Restart** indicates the status,  is enabled,  is disabled. In this screen the auto restart is disabled.

<input checked="" type="checkbox"/> Auto Restart
<input type="checkbox"/> Audible Alarms

To change the status highlight the line. Because the box is also now highlighted  indicates disabled and  indicates enabled. In this screen the feature is still disabled.

<input type="checkbox"/> Auto Restart
<input type="checkbox"/> Audible Alarms

Press the enter key to change the status. In this screen **Auto Restart** is enabled.

<input checked="" type="checkbox"/> Auto Restart
<input checked="" type="checkbox"/> Audible Alarms

Once the highlight is removed confirm the status of the Auto Restart feature. In this screen the **Auto Restart** feature is enabled.

## Audible Alarms

**Audible Alarms** is used to turn the audible alarm on/off.

Line Frequency	60Hz
<input type="checkbox"/> Auto Restart	
<input type="checkbox"/> Audible Alarms	
RA FAN SPEED MODE	
Care Level	off
SERIAL COMM - DCOM	

Use the same procedure described above to enable/disable the alarms.


In this screen the alarms are disabled.

## RA FAN SPEED MODE (ThermoFlex 2500 air-cooled chillers only)

**RA FAN SPEED MODE** controls the fan speed. **Auto** allows the fan to run under the conditions listed in Section 3. Selecting **On** allows the fan to run at high speed all the time.

**Note** **On** is required for the chiller to achieve its full cooling capacity specification. ▲

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.


Fan Speed Range: **On** or **Auto**

Fan Speed Default: **Auto**

## Care Level

**Care Level** sets the preventive care cleaning frequency reminder for the chiller's air and fluid filters.

If a change to the level is needed, highlight **Care Level** and then press the enter button. The highlight will flash. Press the up or down arrow to bring up the desired level. Once the desired level is displayed press enter to accept the change and stop the flashing.

Care Level Range: **off**  
**1** (1000 hours)  
**2** (2000 hours)  
**3** (3000 hours)

Care Level Default: **1** (1000 hours)

See Section 6 for additional details.

## SERIAL COMM - DCOM (Optional)

**SERIAL COMM - DCOM** is used to configure/enable the chiller's optional serial communications feature. If a change to the setting is needed, press the arrows to highlight the desired line and then press the enter button. The highlight will flash. Press the up or down arrow to change the setting. Once the desired setting is displayed press enter to accept the change and stop the flashing.

**Note** None of the controller's other menu displays are available when serial comm is enabled. ▲

SERIAL COMM - DCOM	
Serial Type	off
Baud	xxx
Parity	xxx
Data Bits	8
Stop Bits	X
Address	xx

Serial Type: **RS232, RS485, off**

Baud Range: **9600, 4800, 2400, 1200, 600 or 300**

Parity: **even, odd or none**

Data Bits: Fixed at **8**

Stop Bits: **2 or 1**

Unit ID: **1 to 99** (RS 485 only)

If serial communications is enabled a general message will appear.

**Note** Keypad operation is still available with serial communications enabled. ▲

**MESSAGE!**  
**RS-232 Enabled**  
**Press ENTER**

Press enter to extinguish the message.

Refer to Appendix C for additional information.

## ANALOG COMM - ACOM (Optional)

**ANALOG COMM - ACOM** is used to configure/enable the chiller's optional analog communications feature. Press the arrows to highlight the desired option and then press the enter button. The  box will turn black  indicating that option is enabled. Press enter again to turn it off.

When operating a ThermoFlex7500-10000 with the remote sensor enabled ensure the chiller's response to lowering the setpoint does not result in operation below 10°C process temperature. Operation below 10°C requires the use of 50/50 EG/water or 50/50 PG/water. ▲

ANALOG COMM - ACOM	
<input checked="" type="checkbox"/> Remote Sensor	
<input type="checkbox"/> Remote Start	
<input type="checkbox"/> Remote Setpt	volts
<input type="checkbox"/> Temp Out	volts
RELAY 1	
RELAY 2	
USER RANGE SETTINGS	

When an option is enabled a general message will appear, for example:

**MESSAGE!**  
**Remote Sensor Enabled**  
**Press ENTER**



Press enter to extinguish the message.

Refer to Appendix D for additional information.

## RESISTIVITY (Optional)



**Note** This option is only available on standard-temp chillers. ▲

**RESISTIVITY** enables/configures the resistivity option. With **DE-I CONTROL** highlighted press enter to enable the feature. Next set the limits to the desired values.

Resistivity	Range	Default
Setpoint:	0.2 to 3.0	1.0 MΩ-cm
Interval	0.1 to 0.5	0.1 MΩ-cm
Warning High:	0.0 to 3.5	3.0 MΩ-cm
Warning Low:	0.0 to 3.5	0.5 MΩ-cm


With **TEMP. COMPENSATION** highlighted press enter to turn compensation off or on.

When the desired changes are made press the left arrow, or esc button, to return to the **SETTINGS** display.

## RESTORE DEFAULTS


**RESTORE DEFAULTS** restores the controller back to factory default values. A reset is automatically performed whenever new firmware is installed.

Press the left arrow, or esc button, to return to the **SETTINGS** display.

## mm/dd/yy hh:mm:ss


**mm/dd/yy hh:mm:ss** sets the date and time. Some error messages can display the date and time of occurrence, see Section 7.

Press the arrows to highlight the desired setting. Press enter to start the highlight flashing. Press the arrows to change the setting. Press enter to accept the change and stop the flashing.

## Show Quick Start

**Show Quick Start** reruns the initial start up menu. The Quick Start menu will appear the next time the chiller is started.



With **Show Quick Start** highlighted press enter. The Quick Start menu will appear the next time the chiller is started.

**Note** If **Show Quick Start** is enabled and power is lost, if auto restart was enabled the chiller will not automatically restart once power is restored. In this case only the controller is powered. ▲

## DRAINING (Optional)

**DRAINING** opens/closes the optional anti drainback valve. The chiller must be off to drain the valve. Press enter to start the draining. The display will change to **Close Valve**. When the draining is complete press enter again to close the valve.


## Shut Down

Press the  button.

The controller displays:

With **YES** highlighted, press enter.

To protect the compressor the chiller will enter a 5 to 20 second shut down cycle (colder process fluids take longer) before the refrigeration system and pump shut down.

The display then goes blank.

Using any other means to shut the chiller down can reduce the life of the compressor. ▲

Always turn the chiller off and disconnect it from its supply voltage before moving. ▲

The circuit protector located on the rear of the chiller is not intended to act as a disconnecting means. ▲

## Section 5 Options/Accessories

### Auto Refill

The auto refill provides makeup fluid to replace any fluid lost to evaporation, etc. It requires a pressurized fluid source connection to the ¼" Female Pipe Thread fitting on the rear of the chiller. (If Teflon® tape is used, ensure the tape does not cover the connection's starting-end thread.)

**Note** ThermoFlex7500s and 10000s with a P3, P5 or T5 pump have a ¼" Male brass plug installed in the connection, remove the plug before connecting the makeup fluid. ▲

**Figure 5-1** Auto Refill Fitting

The auto refill fluid must also meet water quality standards or the valve may fail to operate as designed, see Section 3.

The auto refill valve input pressure must be < 80 PSI to ensure the valve functions properly.

**Note** Adding fluid that has a temperature differential with the fluid already in the reservoir will temporarily affect temperature stability performance. ▲

The auto refill operates when all of the following conditions are met:

- Fluid is available
- The chiller is turned on
- The fluid reaches a low level condition.

The auto refill shuts off when:

- The fluid reaches the correct operating level.
- The delay timer exceeds user fill time entered in the Quick Start or **SETTINGS** menu.
- The chiller shuts down for any reason.

**Note** Setting the fill time to 0 disables auto refill. ▲

## Internal DI Cartridge

**Note** This option is not available on high-temp chillers. ▲

A partial flow DI filter cartridge is designed to maintain water resistivity between 1 and 3 MΩ-cm.

**Note** The DI option results in a 0.5 gpm reduction of available flow. ▲


**Do not use a Deionization (DI) filter cartridge with Inhibited EG or Inhibited PG. A DI filter will remove inhibitors from the solution rendering the fluid ineffective against corrosion protection. Also, inhibitors increase fluid conductivity. ▲**

The controller will display a **High Resistivity** or a **Low Resistivity** warning message when the process fluid resistivity exceeds the limit set in the controller's **RESISTIVITY SETTINGS** display, see Section 4.

Remove the two thumbscrews securing the DI access panel to the top of the chiller.


If there is a cartridge in place, first undo the hose fitting by pressing on the quick disconnect located on the top white connection.

**Figure 5-2** Internal DI Cartridge

**The DI Cartridge will overpressurize if it is removed from the chiller before removing the hose fitting. ▲**

**Figure 5-3** DI Fittings

Next rotate the cartridge ¼ turn counter-clockwise and then pull the cartridge straight up to remove it.

Remove the new cartridge from the shipping bag. The cartridge has a blue and a white connector. Lower the cartridge into the chiller with the blue connector facing downward. Press down on the cartridge lightly to engage and then rotate it ¼ turn clockwise (do not over rotate) or until you feel the filter click into place.


Push the hose fitting into the quick disconnect located on the white end of the cartridge. Replace the access panel and thumbscrews.

**Note** The cartridge can be changed with the chiller running, however, since the cartridge runs in a parallel arrangement, disconnecting the cartridge adds 0.5 gpm to the main flow. The additional flow will cause an increase in system pressure which may cause a high fluid pressure fault. ▲


## P1 P2 T0 T1 Pump Pressure Relief Valve (Internal Configuration)

The pressure relief valve, located on the top left rear of the chiller, is used to set the desired system back pressure to your application. The valve is factory preset to  $80 \pm 5$  psi ( $5.5 \pm 0.4$  bar).

*If the chiller is not plumbed to an application,* set the pressure by installing a loop of hose equipped with a shut-off valve between the supply and return fittings. Start the chiller and allow it to prime, then close the valve.

**Figure 5-4** Nut and Screw

**Figure 5-5** Pressure Status Display

Use the Status Display to display P 1, it should display  $80 \pm 5$  psi.

Use a screwdriver to turn the adjusting screw (counterclockwise to reduce pressure) until the controller displays the desired setting.

**Note** Due to internal back pressure, the minimum pressure setting for a deadheaded pump is 32 psi (2.2 bar) for a P2 pump, and 8 psi (0.6 bar) for a P1 (these settings prohibit external flow from the chiller). ▲

*If the chiller is plumbed to an application,* ensure the chiller is off. Then back out the adjusting screw counterclockwise to reduce pressure. Turn the chiller on. Ensure that there is back pressure in the system. Turn the adjusting screw until the controller displays the desired setting.

**Do not exceed 100 psi (6.9 bar).** ▲

When complete, inspect the area around the  $\frac{5}{8}$ " packing nut for fluid. If fluid is present, slightly tighten the nut and reinspect.

**Note** Should the chiller start to vibrate the valve setting may be the cause. Changing the pressure setting  $\pm 5$  psi (0.3 bar) will eliminate the vibration. ▲

## P1 P2 T0 T1 Pump Pressure Relief Valve (External Configuration)


The pressure relief valve is used to set the desired system back pressure (P1) to your application. The valve is factory preset to  $80 \pm 5$  psi ( $5.5 \pm 0.4$  bar).

The valve's inlet/outlet connections are  $\frac{1}{2}$ " FNPT.

*If the chiller is not plumbed to an application, set the pressure by installing a loop of hose equipped with a shut-off valve between the supply and return fittings. Start the chiller and allow it to prime, then close the valve.*


Adjusting Screw

Packing Nut

**Figure 5-6** Nut and Screw

Use the Status Display to display P1, it should display  $80 \pm 5$  psi.

**Figure 5-7** Pressure Status Display

Use a screwdriver to turn the adjusting screw (counterclockwise to reduce pressure) until the controller displays the desired setting.

**Note** Due to internal back pressure, the minimum pressure setting for a deadheaded pump is 40 psi (2.8 bar) for a P2 pump, and 22 psi (1.5 bar) for a P1 (these settings prohibit external flow from the chiller). ▲

*If the chiller is plumbed to an application, ensure the chiller is off. Then back out the adjusting screw counterclockwise to reduce pressure. Turn the chiller on. Ensure that there is back pressure in the system. Turn the adjusting screw until the controller displays the desired setting.*


**Do not exceed 100 psi (6.9 bar). ▲**

When complete, inspect the area around the  $\frac{5}{8}$ " packing nut for fluid. If fluid is present, slightly tighten the nut and reinspect.

## Flow Control with Flow Readout

Flow control for P1, P2, T0 and T1 pumps on ThermoFlex900s - 5000s is achieved using a 3-way valve plumbed between the standard process outlet and the process inlet on the rear of the chiller. Use the auxiliary process outlet at the top left of the rear of the chiller as a connection point. The connections are 1/2" FNPT. See Figure 5-8.


**Figure 5-8** Flow Control

ThermoFlex3500s and 5000s with P3 and P4 pumps use a 2-way valve located on the rear of the chiller. The connections are 3/4" FNPT. See Figure 5-9.

ThermoFlex7500s and 10000s with P2 - P5 and T5 (see next page) pumps use a valve located on the rear of the chiller. The connections are 1/2" for P2, 1" FNPT for P3 and P5. See Figure 5-9.

Use the controller's Status Display to view the flow rate. Turn the valve handle until the desired rate is displayed.

**Figure 5-9** Flow Control Handle (Typical)


**Note** The valve is sensitive to slight adjustments. ▲

## P1 P2 T0 T1 Pump Pressure Relief with Flow Readout

The pressure relief with flow readout works just like the pressure relief Valve discussed on the previous page. It allows you to control the pressure going to your application.

This valve is plumbed into the chiller's auxiliary port, allowing you to also monitor the flow rate to your application using the controller's Status Display.


The valve's outlet connection is 1/2" FNPT. See Figure 5-10.

**Figure 5-10** Pressure Relief

## T5 Pump Flow Control

The flow control valve is used to adjust the flow rate. The valve's handle is designed to identify the valve's position, from full flow to full bypass. When the handle is in the horizontal position (in line with the discharge line) the application is receiving full flow. With the handle is vertical the valve is in full bypass.

## Anti Drainback

Chillers installed below the end-user application may allow system fluid to drain back into the chiller and cause spillage. The anti-drainback valve is designed to prevent any such spillage.

The valve opens just before the pump is turned on and it closes just after the pump shuts off.

This option is required if your chiller is more than 24 feet below your application, or if there is a possibility of drain back due to the occasional opening of the process lines for either application swaps or chiller servicing.

## Semiconductor Equipment and Materials International (SEMI) Chillers

### Compliance

SEMI chillers are compliant with:

**SEMI S2-0703 Product Safety Assessment**

**SEMI S8-0705 Ergonomic Assessment**

**SEMI S14-0704 Fire Risk Assessment**

**SEMI F47-0706**

### Emergency Off (EMO)

A guarded red mushroom shaped push-button switch with twist-to-reset is provided on the chiller's front to turn it off in case of an emergency. The button head is engraved with "EMO" in large white filled letters.

**Note** The EMO is controlled by a safety circuit and is not influenced by the chiller's firmware/software. ▲

Activation of the EMO button will remove power from the main contactor coil stopping operation of the chiller. The controller will display **External EMO**.

Resetting the EMO button will not restart the chiller. After all hazards have been removed reset the chiller by pressing **enter** on the controller. In the local mode, the chiller will restart by pressing the START STOP button again. In the serial communications mode, send the appropriate start command. In the analog I/O mode, the chiller starts when the error is cleared.

### Chiller Circuit Breaker Interrupt Rating

The main power circuit breaker located on the rear of the chiller has an Interrupting Capacity (AIC) of 10,000 amps.

### Lockout/Tagout (LOTO)

Before performing Chiller maintenance, the energy sources associated with the chiller system must be lockedout and tagged out (LOTO). Hazard control features added to the system (e.g., safety interlocks, EMO) are not a substitute for turning off and locking out electrical or fluid energy.


For chillers rated 20 Amps or less, electrical LOTO is accomplished by removing the power cord on the rear of the chiller then closing and locking the power receptacle locking device. For other chillers, electrical LOTO is the responsibility of the user and can be provided by:

- Using the main disconnect (knife switch at system control cabinet).
- Disconnecting main power at the facility power source prior to the system controller cabinet.
- In addition, follow all OSHA and local facility LOTO directives.

### Drip Pan and Drain

The chiller is equipped with a secondary containment (drip pan) in case there is a leak. The drip pan drain is located on the rear of the chiller. Install the supplied nylon 1/4 turn quick disconnect (QD) fitting into the drain fitting. The QD is barbed for a 1/2" ID hose.

Since the drip pan will not hold more than 110% of the reservoir volume, connect the drain to guide the fluid to an appropriate spillage location.

1/4 Turn Quick Disconnect Drip Pan Drain Fitting


**Figure 5-11** Drip Pan Drain

	900/1400		2500		3500/5000		7500/10000	
A	3 1/2"	8.8 cm	4"	10.1 cm	3 3/8"	11.3 cm	4 1/4"	10.8 cm
B	2 3/4"	7.0 cm	2 11/16"	6.8 cm	2 3/4"	7.1 cm	2 5/8"	6.6 cm
C	6 15/16"	17.7 cm	6 9/16"	16.7 cm	9 9/16"	24.3 cm	7 11/16"	19.5 cm

### Seismic Tie-Downs

Install the seismic tie-downs to the chiller as shown below. Then secure the chiller to the floor with user-supplied hardware.


**Figure 5-12** Seismic Tie-Downs

	900/1400		2500		3500/5000		7500/10000	
A	2 11/16"	6.8 cm	2 11/16"	6.8 cm	2 11/16"	6.8 cm	2"	5.1 cm
B*	18 1/2"	47.0 cm	20 1/16"	51.0 cm	24 1/2"	62.2 cm	17"	43.1 cm
C*	19 11/16"	50.0 cm	22 3/4"	57.8 cm	24 3/4"	62.9 cm	27 7/16"	69.6
D	21 3/16"	53.8 cm	24 1/4"	61.5 cm	26 1/4"	66.7 cm	28 15/16"	73.4

\* Distance between Ø.53 Seismic mounting holes

## Center of Gravity

**Figure 5-13** Center of Gravity

**Center of Gravity**  $\pm \frac{1}{2}$ ", air-cooled chiller, no fluid in tank

	900/1400 P2 Pump		2500 P2 Pump		3500/5000 P2 Pump		7500/10000 P3 Pump	
A	10 <sup>3</sup> / <sub>4</sub> "	27.3 cm	12"	30.5 cm	13 <sup>3</sup> / <sub>8</sub> "	34.0 cm	14 <sup>7</sup> / <sub>8</sub> "	37.8 cm
B	6 <sup>3</sup> / <sub>4</sub> "	17.2 cm	8 <sup>3</sup> / <sub>8</sub> "	21.3 cm	9"	22.9 cm	13 <sup>1</sup> / <sub>8</sub> "	33.3 cm
C	13 <sup>1</sup> / <sub>2</sub> "	34.3 cm	13 <sup>1</sup> / <sub>2</sub> "	34.3 cm	17"	43.2 cm	26"	66.0 cm

### 3500/5000 P4 Pump Global Voltage

A	12 <sup>3</sup> / <sub>8</sub> "	31.4 cm
B	9 <sup>3</sup> / <sub>4</sub> "	24.8 cm
C	19 <sup>1</sup> / <sub>2</sub> "	49.5 cm

**Center of Gravity**  $\pm \frac{1}{2}$ ", water-cooled chiller, no fluid in tank

### 3500/5000 P2 Pump

A	13"	33.0 cm
B	9 <sup>1</sup> / <sub>2</sub> "	24.1 cm
C	16"	40.6 cm

**Weight Distribution**  $\pm 2$  lbs, air-cooled chillers

	900/1400 P2		2500 P2		3500/5000 P2		7500/10000 P3	
Left Front	27.1 lbs	12.3 kg	40.7 lbs	18.5 kg	62.0 lbs	28.1 kg	97.8 lbs	44.4 kg
Left Rear	29.8 lbs	13.5 kg	42.0 lbs	19.1 kg	63.7 lbs	28.9 kg	99.9 lbs	45.3 kg
Right Front	32.9 lbs	14.9 kg	45.7 lbs	20.7 kg	68.2 lbs	30.9 kg	89.2 lbs	40.5 kg
Right Rear	36.2 lbs	16.4 kg	47.1 lbs	21.4 kg	70.0 lbs	31.8 kg	91.1 lbs	41.3 kg

## Other Accessories

Installation kit - includes replacement air and fluid filters

Maintenance kit - includes a set of hoses, adaptor fittings and Teflon<sup>®</sup> tape

### Fluids

Fluid treatment kit designed to minimize the effects of corrosion, scale, fouling, and microbial contamination. It allows the system to continue providing reliable service with optimal efficiency for the life of the chiller.

The kit includes a biocide and corrosion inhibitor capable of treating up to ten gallons of application water and is designed to provide protection for a period of six months. This kit is compatible with the following fluids:

- Filtered/Single Distilled Water
- Uninhibited Ethylene Glycol/Water
- Uninhibited Propylene/Water
- Deionized (DI) Water\*
- Reverse Osmosis (RO) Water

\*Do not use the Thermo Fisher Water Treatment Kit with a DI filtered system; the filter will remove a portion of the reagent's active ingredients limiting its effectiveness.

Please contact Thermo Fisher Scientific's Sales, Service and Customer Support to assist you with questions that you may have regarding accessories for your ThermoFlex, see inside front cover for contact information.

## Section 6 Preventive Maintenance

### Preventive Maintenance Timer (Care Level)

Only Thermo Fisher should provide any required replacement parts.

The ThermoFlex chiller has an integrated preventive maintenance timer that can alert you when it is time to perform preventive maintenance. This unique feature can be set to remind you to change your air and fluid filters.

Based on the environment in which your chiller is located, you can choose from four levels of preventive maintenance off, L1, L2, and L3:

- off – Disables the alert
- L1 – 1,000 hours - default setting
  - Heavy manufacturing environment
  - Airborne particulate created during manufacturing process
- L2 – 2,000 hours
  - Typical production environment
- L3 – 3,000 hours
  - Clean environment – filtered air
  - Typically laboratory or research environment

Each time the chiller exceeds the chosen time, the controller will flash **Filter PM** and, if enabled, an audible alarm will sound.

Change/set the level using the controller's **DIAGNOSTICS** menu, see next page.


To clear this message press . This will temporarily clear the alarm until the chiller is restarted. To permanently clear the alarm the filter run time hours must be reset in the **DIAGNOSTICS** menu. Also, if you change your filters before the preventive timer trips, you can clear the timer by using the **DIAGNOSTICS** menu.

**Note** For air-cooled chillers, both the air and fluid filters in the ThermoFlex can be changed while the chiller is running. For water-cooled chillers, only the fluid filter can be changed while it is running. ▲


## DIAGNOSTIC - OPERATING HOURS

**OPERATING HOURS** displays the chiller (**Unit**), filter and optional DI cartridge operating hours. The display is also used to select the preventive maintenance schedule (**Care Level** and optional **De-I Period**).


1. Use the arrow buttons to highlight **OPERATING HOURS**.

2. Press enter to display:

3. Highlight the desired item and press enter to bring up one of the displays shown on the right:


4. With **Reset** highlighted, press enter to reset the filter and cartridge operating hours to zero. You cannot reset the chiller's operating hours.

5. Use the arrow buttons to highlight **Care Level**. Press enter and the display will flash. Use the arrow buttons to display the desired care level and then press enter again to stop the flashing and accept the new value.

6. The optional **De-I Period** is the operating time needed to display the **DI** message. Use the arrow buttons to highlight **De-I Period**. Press enter and the display will flash. Use the arrow buttons to display the desired time, from 0 to 9999 hours. Press enter again to stop the flashing and accept the new value.

(Optional Display)

## Fluid Filter Bag

The reservoir has a fluid bag filter designed to prevent the introduction of particulates into the system.

**Note** The fluid bag filter can be removed with the chiller operating. ▲

**Before using any fluid or performing maintenance where contact with the fluid is likely refer to the manufacturer's SDS for handling precautions.** ▲

**For high-temperature chillers, ensure the fluid is at a safe temperature (below 40°C) before handling.** ▲

Fluid Bag Filter

When it is time to replace the bag, gently pull up on the plastic funnel housing to remove it and simply pull the bag out of the chiller. Replacement bags are available from Thermo Fisher Scientific.

**Figure 6-1** Fluid Filter Bag

**Before replacing the reservoir housing ensure the reservoir sight tube ball stopper is securely in place, see next page.** ▲

## Fluid Diffuser

On ThermoFlex900-5000s, when you remove the bag you will notice a wire mesh fluid diffuser inside the reservoir supply line, see Figure 6-2. The diffuser is used to help streamline the flow into the reservoir. After several bag replacements turn the chiller off and remove the diffuser to inspect it for debris/damage.

**The fluid velocity into the reservoir will rapidly increase with the diffuser removed and cause splashing. Turn the chiller off before removing the diffuser. This is especially critical when using ethylene or propylene glycol.** ▲

**Note** To prevent particulates from entering the reservoir, ensure the fluid bag filter is in place before removing the diffuser. ▲

**Do not operate the chiller unless the diffuser is installed.** ▲

## Fluid Maintenance

An effective recommended maintenance plan would include changing the fluid every six months to optimize chiller reliability, see Section 3 for additional information.

## Reservoir Cleaning

The user is responsible for maintaining reservoir fluid quality. Check the fluid on a regular interval. Start with frequent checks until a regular interval (based on your application) is established.

If cleaning is necessary, flush the reservoir with a fluid compatible with the process fluid and the chiller's wetted parts, see Section 8.

**Before using any fluid or performing maintenance where contact with the fluid is likely refer to the manufacturer's SDS for handling precautions. ▲**

**For high-temperature chillers, ensure the fluid is at a safe temperature (below 40°C) before handling. ▲**


### Reservoir Sight Tube

Clean the sight tube by gently pulling up on the plastic funnel housing to remove it (see illustration on previous page) and then gently pulling out the black sight ball stopper from the tube. Use a long soft-bristle 1/4" brush. Use caution not to scratch the glass.

**Before replacing the reservoir housing ensure the reservoir sight tube ball stopper is securely in place. ▲**

For easier replacement, wet the stopper first and then use a twisting motion to install it in the sight tube.

**Figure 6-2** Reservoir Cleaning and Diffuser

## Condenser Filter

Failure to clean/replace the condenser filter will cause a loss of cooling capacity and lead to premature failure of the cooling system. ▲

### ThermoFlex900s - 5000s

Clean the filter through the grill using a vacuum with a soft-bristle brush.

When it is time for a more thorough cleaning, remove the one-piece grill assembly by first pulling the bottom of the assembly away from the chiller and then pulling it away from the top.

The condenser framing and fins located behind the grill assembly are very sharp. Use caution when removing the assembly. ▲

**Note** ThermoFlex900s - 5000s water-cooled chillers have an embedded screw(s) located at the top (and bottom) of the grill securing it to the chiller. Loosen the screw(s) to remove the grill. ▲


Water-cooled chillers also have a fan with sharp blades, ensure the chiller is off before removing the assembly. ▲

Shake off as much of the excess water as possible before reinstalling. Press the grill back into place.

For water-cooled chillers, tighten the screw(s) at the top (and bottom) of the grill.

Replacement grill assemblies are available from Thermo Fisher.

**Figure 6-3** ThermoFlex900s - 5000s Condenser Grill Removal

### **ThermoFlex7500s - 10000s**

For air-cooled chillers, remove the one-piece grill assembly by pulling the assembly away from the chiller.


Water-cooled chillers do not have a filter.

The filter goes over four studs and plastic "fast nuts" that hold it in place.

Replace it or vacuum the old filter with a soft-bristle brush, or wash it. Shake off as much of the excess water as possible before reinstalling.

Tuck the filter around the perimeter of the grill and over the four studs, use the plastic "fast nuts" to hold it in place.

Replacement grills are available from Thermo Fisher.

**Figure 6-4** Filter Removal/Replacement ThermoFlex7500s - 10000s Air-Cooled

### **Chiller Surface**

Clean the chiller's surface with a soft cloth and warm water only.

### **Hoses**

Inspect the chiller's external hoses and clamps on a daily basis.

## Testing the Alarm Features

Using the Setup Loop, adjust each temperature alarm limit towards the setpoint and ensure the chiller reacts accordingly. Reset each alarm limit to the desired value. See Section 4.

We recommend slowly draining the chiller to ensure the low level warning and fault messages activated. See Section 8.

For chillers equipped with auto refill switch we recommend slowly draining and ensure the auto refill activates.

## DI Filter (Optional)

Establish a preventive maintenance schedule for the DI filter cartridge based on your specific application.

The controller will display a **High Resistivity** or a **Low Resistivity** warning message when the process fluid resistivity exceeds the limit set in the controllers **RESISTIVITY SETTINGS** display, see Section 4.


The chiller also has an integrated alarm. The alarm is based on chiller run hours that will alert you when it is time to change your filter. The alarm is enabled using the controller's **DIAGNOSTICS** menu.

If you already know how often your DI filter needs changing, you can input the number of hours into the menu. When the time is reached, the controller will flash **DEI PM** and the audible alarm, if enabled, will sound.

To clear this message and stop the audible alarm press  .

## Preventive Maintenance Messages

The controller also displays Preventive Maintenance Messages. These messages are based on the component run time and are established to avoid unplanned failures. The chiller will continue to run.

**Figure 6-5** Sample PM Message

<b>Filter PM</b>	<b>Fluid Filter Bag and Condenser Filter</b>
<b>DEI PM</b>	<b>Optional DI Filter</b>
<b>XXXXX Run Hours</b>	<b>Chiller operating hours (appears every 20,000 hours)</b>

You can press the enter button to clear any message. Note the message because once the Filter or DEI PM message is cleared it will not reappear.

## Section 7 Troubleshooting

### Messages

WARNING and FAULT Messages are a result of exceeding one of the controller's SETTINGS, see Section 4, exceeding a sensor factory preset safety value, or a safety switch is activated. An ERROR Message indicates an unusual condition. In the case of an ERROR or WARNING Message the chiller, if running, will continue to run. Press **enter** to see if the message clears, a limit may have been only temporarily exceeded. In the case of a FAULT Message the chiller will shut down and the controller will continue to flash the message. Press **enter** to clear the display and, if enabled, silence the alarm. Once the cause of the shut down is identified and corrected, restart the chiller. If the cause was not corrected the message will reappear.

<p><b>WARNING!</b> Low Level</p> <p>Press ENTER</p>	<p><b>FAULT!</b> High Temp</p> <p>Press ENTER</p>	<p><b>ERROR!</b> CAN Bus</p> <p>Press ENTER</p>
---	---	---

Figure 7-1 Sample Messages

### Faults and Warnings

Message	Reaction	Cause	Actions
<b>Auto Refill</b>	The chiller will continue to run. Auto refill, if installed, will shut down. (Optional display - only chillers equipped with auto refill.)	The auto refill did not reach the minimum operating level within the time chosen for the customer adjustable <i>fill</i> setting.  The auto refill successfully filled within the time frame chosen for the customer adjustable <i>fill</i> setting, but the chiller tries to refill 5 times in 40 hours.	<ul style="list-style-type: none"> <li>• Check for leaks.</li> <li>• Check the supply pressure on the auto refill supply line. With low pressure the auto refill time span setting may be set too low and the reservoir does not have time to fill. Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Add fluid to the tank.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Bad Calibration Data</b>	Chiller continues to run.	Bad sensor calibration detected several seconds after performing a calibration.	<ul style="list-style-type: none"> <li>• Redo calibration, see Section 8.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>Drip Pan</b>	Chiller will shut down. (SEMI chillers only)	Fluid in the drip pan.	<ul style="list-style-type: none"> <li>• Check for leaks.</li> <li>• Remove the fluid from the drip pan and reset the fault.</li> <li>• Contact our Sales, Service and Customer Support.</li> <li>• When able, reset external EMO..</li> </ul>
<b>External EMO</b>	Chiller will shut down. (SEMI chillers only)	External EMO depressed.	<ul style="list-style-type: none"> <li>• Allow chiller to cool down.</li> <li>• For air-cooled chillers, clean air filter, see Section 6.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Fan Motor Overload</b>	Chiller will shut down. (3- $\Phi$ chillers only)	Fan motor overload activated.	<ul style="list-style-type: none"> <li>• Check application valves and ensure that they have not changed or been closed. <b>Note</b> If routine shut-off of the process flow is required then add an external pressure regulator accessory - contact Thermo Fisher. ▲</li> <li>• May occur as a result of changing the internal DI cartridge. Disconnecting the cartridge adds an additional 0.5 GPM to the main flow, see Section 5.</li> <li>• Check for debris in the application or clogged external filters.</li> <li>• Double check fluid lines. Excessive bends, long tubing and diameter reductions can affect the pump's discharge pressure. <b>Note</b> If diameter reductions must be made, make them at the inlet and outlet of your application, not at the chiller. ▲</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>High Fixed Pressure</b>	Chiller will shut down. This error code has priority over <b>High Pressure</b> error code.	Process pressure (P1) exceeded <i>factory preset</i> value for greater than 30 seconds. Preset Values: P1, P2 and T1- 105 psi P3 60 Hz - 48 psi P3 50 Hz - 32 psi P4 60 Hz - 85 psi P4 50 Hz - 60 psi P5 60 Hz - 87 psi P5 50 Hz - 56 psi	

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>High Fixed Temp</b>	<p>Chiller will shut down.</p> <p>This error code has priority over <b>High Temp</b> error code.</p> <p><b>Note</b> Chiller will not restart until process fluid temperature is below +43°C (+93°C for high-temp chillers). ▲</p>	<p>Reservoir fluid temperature exceeded the <i>factory preset</i> value of +43°C (+93°C for high-temp chillers).</p>	<ul style="list-style-type: none"> <li>• Ensure all environmental requirements are met, see Section 3.</li> <li>• Ensure chiller has adequate ventilation, see Section 3.</li> <li>• Clean air filter. Dirt and debris on filter can prevent the chiller from functioning at full capacity, see Section 6.</li> <li>• Ensure that the heat load being applied to the chiller is not too high. Contact Thermo Fisher for assistance on calculating heat loads.</li> <li>• Bring cooler air in from another area or exhaust the hot air into another location using an auxiliary fan.</li> <li>• Verify/adjust controller PID values, see <b>TUNING</b> in this Section.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>High Flow</b>	<p>With a warning message the chiller will continue to run.</p> <p>With a fault message the chiller will shut down.</p> <p>(Optional display - only chillers equipped with a flow transducer.)</p>	<p>The process fluid flow rate has exceeded the adjustable setting's high value.</p>	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the code clears, the limit may have been only temporarily exceeded.</li> <li>• Verify <b>SETTINGS</b>, see Section 4, and adjust setting if necessary.</li> <li>• Check all application and plumbing shut off valves for correct position.</li> <li>• Adjust flow, if chiller is equipped with an optional flow control valve, see Section 5.</li> <li>• If flow transducer was recently calibrated double check calibration, see Section 8.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>High Pressure</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down.	The pump's discharge pressure exceeded the adjustable high value.	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Check application valves and ensure that they have not changed or been closed. <b>Note</b> If routine shut-off of the process flow is required then add an external pressure relief valve, see Section 5. ▲</li> <li>• Check for debris in the application or external filters.</li> <li>• May occur as a result of changing the internal DI cartridge. Disconnecting the cartridge adds an additional 0.5 gpm to the main flow, see Section 5.</li> <li>• Double check fluid lines. Excessive bends, long tubing and diameter reductions can affect the pump's discharge pressure. <b>Note</b> If diameter reductions must be made, make them at the inlet and outlet of your application, not at the chiller. ▲</li> <li>• Contact our Sales, Service and Customer Support.</li> <li>• Make sure supply voltage matches the chiller's nameplate rating <math>\pm 10\%</math>.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>High RA Temperature</b>	Chiller will shut down.	Refrigeration suction gas temperature exceeded 50°C.	<ul style="list-style-type: none"> <li>• Press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Replace DI filter and/or process fluid.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>High Resistivity</b>	Chiller will continue to run. (Optional display)	The process fluid resistivity exceeded the upper adjustable value.	<ul style="list-style-type: none"> <li>• Press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Replace DI filter and/or process fluid.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>High Temp</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down. <b>Note</b> If the chiller does shut down it can be restarted provided the temperature is still within the factory-set high fixed temperature limit. However, the error will reoccur if the temperature goes below the adjustable setting and then again exceeds it. ▲	The process fluid temperature exceeded the adjustable high value.	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Ensure all environmental requirements are met, see Section 3.</li> <li>• Ensure chiller has adequate ventilation, see Section 3.</li> <li>• Clean air filter. Dirt and debris on filter can prevent the chiller from functioning at full capacity, see Section 6.</li> <li>• Ensure that the heat load being applied to the chiller is not too high. Contact Thermo Fisher for assistance on calculating heat loads.</li> <li>• Bring cooler air in from another area or exhaust the hot air into another location using an auxiliary fan.</li> <li>• Verify/adjust controller PID values, see <b>TUNING</b> in this Section.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>HPC</b>	Chiller will shut down.	High refrigeration pressure cutout activated.	<p><b>Air-cooled chillers</b></p> <ul style="list-style-type: none"> <li>• Ensure that the ambient temperature is not exceeding the recommended range, see Section 3.</li> <li>• Ensure chiller has adequate ventilation, see Section 3.</li> <li>• Clean air filter, see Section 6.</li> <li>• Bring cooler air in from another area or exhaust the hot air into another location using an auxiliary fan.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul> <p><b>Water-cooled chillers</b></p> <ul style="list-style-type: none"> <li>• Ensure facility water is on and connected.</li> <li>• Check facility water flow rate and pressure.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>HTC</b>	Chiller will shut down.	Low fluid level or High heater temperature.	<ul style="list-style-type: none"> <li>• Check reservoir fluid level.</li> <li>• Allow heater to cool and then remove the small black plastic cap(s) on the rear of the chiller, press the reset node(s) until you hear a “click.” See <b>Checklist</b> in this Section.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>Invalid Level</b>	Chiller will shut down.	Invalid level fault. Chiller sensed both a high level and low level reservoir fluid level.	<ul style="list-style-type: none"> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Invalid rem setpt</b>	Chiller will continue to run using the last valid setpoint received. (Optional display - only chillers equipped with Analog I/O).	Analog remote setpoint is enabled and the chiller receives a voltage or current level that is outside the chiller's set point range.	<ul style="list-style-type: none"> <li>The error can be cleared only after a valid set point is received, or the remote analog setpoint is turned off.</li> </ul>
<b>LLC</b>	Chiller will shut down.	Low Level Cutout activated.	<ul style="list-style-type: none"> <li>Check for leaks.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Local EMO</b>	Chiller will shut down. (Optional display)	Chiller's EMO button depressed.	<ul style="list-style-type: none"> <li>When able, reset chiller's EMO.</li> </ul>
<b>Low Fixed Flow</b>	Chiller will shut down. (Optional display - only chillers equipped with a flow transducer) (This error code has priority over <b>Low Flow</b> .)	Low flow fault. For chillers with a P1, P2 or T1 pump the flow dropped below 0.8 GPM for more than 15 seconds.  For chillers with a P3, P4 or P5 pump the flow dropped below 3.8 GPM for more than 15 seconds.	<ul style="list-style-type: none"> <li>Adjust flow, if chiller is equipped with an optional flow control valve, see Section 5.</li> <li>Check all valves in your application and plumbing lines to ensure that they have not changed or closed.</li> <li>If flow transducer has recently been calibrated, double check calibration to ensure it was done properly, see Section 4.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Low Fixed Pressure</b>	Chiller will shut down. (This error code has priority over <b>Low Pressure</b> .)	Process pressure (P1) below <i>factory preset</i> low limit of 2 psi for greater than 15 seconds.	<ul style="list-style-type: none"> <li>Ensure that the chiller reservoir is not empty.</li> <li>Chiller requires &gt;2 PSIG application pressure drop. If a bypass valve has been installed, some restriction may need to be added to the bypass line.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>Low Fixed Temp</b>	Chiller will shut down.  <b>Note</b> Chiller will not restart until process fluid temperature exceeds +2°C. ▲ (This error code has priority over <b>Low Temp</b> .)	Reservoir fluid temperature below the <i>factory preset</i> low value of +2°C.	<ul style="list-style-type: none"> <li>• Check ambient temperature. Chiller may not be able to reach setpoint at low ambient temperatures. If your load is constant, then turn your chiller on. Chiller will control setpoint when sufficient heat is added.</li> <li>• Verify/adjust controller PID values, see <b>TUNING</b> in this Section.</li> <li>• Add insulation to external plumbing lines to reduce the heat-loss to the environment.</li> <li>• Ensure that the ambient temperature is not exceeding the recommended range, see Section 3.</li> <li>• For water-cooled chillers check facility water temperature.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Low Flow</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down.  (Optional display - only chillers equipped with a flow transducer.)	The process fluid flow rate has gone below the adjustable setting's low value.	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the code clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Adjust flow, if chiller is equipped with an optional flow control valve, see Section 5.</li> <li>• Check all valves in your application and plumbing lines to ensure that they have not changed or closed.</li> <li>• If flow transducer has recently been calibrated, double check calibration to ensure it was done properly, see Section 4.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Low Level</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down.	Reservoir fluid level too low for normal operation.	<ul style="list-style-type: none"> <li>• Excessive evaporation. Ensure the chiller is operating with the funnel and cap in place.</li> <li>• Check for leaks.</li> <li>• Check auto refill operation, see Section 5.</li> <li>• Check the supply pressure on the auto refill supply line. With low pressure the auto refill time span setting may be set too low and the reservoir does not have time to fill. Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings

Message	Reaction	Cause	Actions
<b>Low Pressure</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down.	Pump's discharge pressure is below adjustable low setting.	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Ensure that chiller reservoir is not empty.</li> <li>• Chiller requires &gt;3 PSIG application pressure drop. If a bypass valve has been installed, some restriction may need to be added to the bypass line.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Low Resistivity</b>	Chiller will continue to run. (Optional display)	The process fluid resistivity exceeded the lower adjustable value.	<ul style="list-style-type: none"> <li>• Press enter to see if the message clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Replace DI filter and/or process fluid.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>Low Temp</b>	With a warning message the chiller will continue to run. With a fault message the chiller will shut down. <b>Note</b> If the chiller does shut down it can be restarted provided the temperature is still above the factory-set low fixed temperature limit. However, the error will reoccur if the temperature goes above the adjustable setting and then again drops below it. ▲	Process fluid temperature is below adjustable low value.	<ul style="list-style-type: none"> <li>• If the chiller is still running press enter to see if the code clears, the limit may have been only temporarily exceeded.</li> <li>• Verify controller's <b>SETTINGS</b>, see Section 4.</li> <li>• Ensure that the ambient temperature is not exceeding the recommended low-range, see Section 3. If your application load is constant and/or the lower temperature can be temporarily tolerated, then continue operation. (The ThermoFlex will control setpoint when sufficient heat is added.)</li> <li>• Add insulation to external plumbing lines to reduce the heat-loss to the environment.</li> <li>• For water-cooled chillers check facility water temperature.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>LPC</b>	Chiller will shut down.	Low refrigeration pressure cutout switch activated.	<ul style="list-style-type: none"> <li>• Check for refrigerant leak.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

## Faults and Warnings


Message	Reaction	Cause	Actions
<b>Open RTD</b>	Chiller will shut down.	Internal sensor open.	<ul style="list-style-type: none"> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Open Remote RTD</b>	Chiller will shut down.	Remote temperature sensor not connected or open.	<ul style="list-style-type: none"> <li>Check the connections on the rear of the chiller.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Over Flow</b>	Chiller will shut down.	There is an overflow condition in the reservoir.	<ul style="list-style-type: none"> <li>Ensure the reservoir was not filled above the MAX LEVEL line.</li> <li>Check for clogged reservoir filter.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Shorted Remote RTD</b>	Chiller will shut down.	Remote temperature sensor shorted.	<ul style="list-style-type: none"> <li>Check the connections on the rear of the chiller.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Shorted RTD</b>	Chiller will shut down.	Internal sensor shorted.	<ul style="list-style-type: none"> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Motor Overload</b>	Chiller will shut down. (Chiller equipped with 3- $\Phi$ pump motor overload)	Pump motor exposed to excessive current due to high pressure, flow or ambient temperature.	<ul style="list-style-type: none"> <li>Allow pump to cool down.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>Phase Monitor</b>	Chiller will shut down. (3- $\Phi$ chillers only)	Phase rotation is wrong.	<ul style="list-style-type: none"> <li>Disconnect chiller from power source and reverse any two line conductors on the line side of the main circuit breaker.</li> <li>Contact our Sales, Service and Customer Support.</li> </ul>
<b>UNIT FAULT</b>	Chiller will shut down.	Critical error.	<ul style="list-style-type: none"> <li>Contact our Sales, Service and Customer Support.</li> </ul>

**Error**


Message	Reaction	Cause	Actions
<b>Async Rec Err Int</b>	Chiller will continue to run.	Communication error between display and main control board.	<ul style="list-style-type: none"> <li>• Check the serial communication connection.</li> <li>• Cycle circuit protector on the rear of the chiller off and on.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>CAN Bus</b>	Chiller will continue to run.	Internal communications error.	<ul style="list-style-type: none"> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>I2C Bus</b>	Chiller will continue to run.	Internal communications error.	<ul style="list-style-type: none"> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>NO 5V CAL</b>	Chiller will continue to run.	No 5V calibration stored in controller memory.	<ul style="list-style-type: none"> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>NO PRES CAL</b>	Chiller will continue to run.	No pressure calibration stored in controller memory.	<ul style="list-style-type: none"> <li>• Perform a pressure calibration, see Section 8.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>NO RTD CAL</b>	Chiller will continue to run.	No temperature calibration stored in controller memory.	<ul style="list-style-type: none"> <li>• Perform a temperature calibration, see Section 8.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>NVS 3 4 RESET</b>	Chiller will not start.	Bad checksum when chiller is started. (Normal when new software is installed.)	<ul style="list-style-type: none"> <li>• Press enter to try and clear the message.</li> <li>• Try a restart.</li> <li>• Contact our Sales, Service and Customer Support.</li> </ul>
<b>SENSE 5V</b>	Chiller will continue to run.	Voltage out of tolerance.	<ul style="list-style-type: none"> <li>• Contact our Sales, Service and Customer Support.</li> </ul>

**MESSAGES** displays the ten most recent warnings and faults. The date and time each error occurred is shown. The alarm value is also displayed. Use the arrow buttons to scroll through the list of messages. Press enter to show the deletion options, **Single** or **All**. Choose **Single** to delete the selected error or **All** to delete all errors in the category. Press enter to show the delete confirmation screen. Finally, press enter again to make the deletion.


1. From the Main Menu display use the arrow buttons to highlight **DIAGNOSTICS**.


2. Press enter and use the arrow buttons to highlight **MESSAGES**.


3. Press enter to display:


4. Highlight the desired type message and press enter again to display:


5. If desired, press enter again to display:

6. If **Yes** is selected press enter again to display:

For either message type, if there aren't any messages the display will indicate:

## Checklist

### Chiller will not start

Check the electrical connections.

For first time use, please refer to the quick start instructions included with your chiller or the copy in this manual. The manual's copy follows the Table of Contents.

Check the controller for messages, see Messages in this Section.

Ensure the circuit protector is in the on ( **I** ) position.

Make sure supply voltage is connected and matches the chiller's nameplate rating  $\pm 10\%$

### Chiller shuts down

Check the electrical connections.

Ensure  button wasn't accidentally pressed.

Ensure the circuit protector is in the on ( **I** ) position.

Check the controller for messages, see Messages in this Section.

The chiller is designed to shut down if not properly primed, refer to Section 3 for priming instructions.

Make sure supply voltage is connected and matches the chiller's nameplate rating  $\pm 10\%$ .

Restart the chiller.

### Clearing Messages

Note the code in case it clears before you are done troubleshooting.

If desired, silence the audible alarm by pressing enter.

*If the chiller shut down* the controller will continue to display the message. Press **enter** to clear the display and silence any alarm. Refer to Messages in this section. Once the cause of the shut down is identified and corrected, start the chiller. If the cause was not corrected the message will reappear.

*If the chiller is still running* press **enter** to see if the message clears, a limit may have been only temporarily exceeded. If the message does not clear refer to Messages in this Section.

### Inadequate pump pressure

Ensure any user installed in-line valves are in the desired position.

Ensure the chiller's process fluid outlet is connected to the application's fluid inlet and not the application's fluid outlet, see Section 3.

Ensure all connections are secure and that the proper sealant/lubricant for the fitting material is used.

Keep the distance between the chiller and the instrument being cooled as short as possible.

Ensure tubing is straight and without bends. If diameter reductions are required, make them at the inlet and outlet of your application, not at the chiller.

### **Chiller will not circulate process fluid**

Check the reservoir level. Fill, if necessary.

Ensure the reservoir bag filter is not clogged.

Check the application for restrictions in the cooling lines.

Chiller requires >3 PSIG application pressure drop. If a bypass valve has been installed, a restriction may need to be added to the bypass line.

The pump motor overloaded. The pump's internal overtemperature overcurrent device will shut off the pump causing the flow to stop. This can be caused by low fluid, debris in system, operating chiller in a high ambient temperature condition or excessively confined space. Allow time for the motor to cool down.

Make sure supply voltage matches the chiller's nameplate rating  $\pm 10\%$ .

### **Inadequate temperature control**

Adding fluid that has a temperature differential with the fluid already in the reservoir will temporarily affect temperature stability performance.

Verify the setpoint.

For high-temperature chillers check the heater's high temperature safety. Remove the two small black plastic covers on the rear of the chiller and press on each of the black stems. If the safety activated you should hear a "click." If it activates again contact us.

If the chiller is over-cooling, recycle the power.

Make sure the condenser/air filter is free of dust and debris.

Check the fluid concentration, see Section 3.

Ensure your chiller's installation complies with the site requirements listed in Section 3.

Make sure supply voltage matches the chiller's nameplate rating  $\pm 10\%$ . For ThermoFlex900 - 5000 Global Voltage chillers ensure it is properly configured, see Appendix B.

If the temperature continues to rise, make sure your application's heat load does not exceed the rated specifications.

Check for high thermal gradients (e.g., the application load is being turned on and off or rapidly changing).

If operating at high altitude note that heat removal capacity decreases 1.2% per 1,000 feet above sea level. Also, reduce the maximum temperature for the air entering the ThermoFlex by 1°C per 1,000 feet above sea level.

Verify/adjust controller PID values, see next page. Ensure the chiller was shut down properly, see Section 4. If not the compressor may be damaged.

### **Chiller vibration**

The optional pressure relief valve setting may be the cause. If it is, change the pressure setting  $\pm 5$  psi to eliminate the vibration.

Please contact Thermo Fisher Scientific Sales Service and Customer Support if you need any additional information, see inside cover for contact instructions.


# TUNING

## Verifying/ Adjusting the Controller PID Values


The controller controls temperature using a Proportional-Integral-Derivative (PID) algorithm. Should your chiller experience temperature control issues, verifying/adjusting the controller's PID values may correct the condition.

**Note** Thermo Fisher recommends that only a qualified technician adjust the PID values. Incorrect values will hamper chiller performance. ▲

1. From the Main Menu display use the arrow buttons to highlight **TUNING**.


2. Press enter and highlight **COOL PID** or **HEAT PID**.

**Note** **HEAT PID** only appears on chillers equipped with a heater. ▲

3. Highlight and change the desired values as required.

### COOL

**P** proportional value, % of span (100°C)

**Range:** 0.0 to 99.9

**Factory Preset:** ThermoFlex900-5000 10.0  
ThermoFlex7500-10000 20.0

**I** integral value, repeats/minute


**Range:** 0.00 to 9.99

**Factory Preset:** ThermoFlex900-10000 0.50

**D** derivative value, minutes

**Range:** 0.0 to 9.9

**Factory Preset:** ThermoFlex900-10000 0.0

### HEAT

**P** proportional value, % of span (100°C)

**Range:** 0.0 to 99.9

**Factory Preset:** ThermoFlex1400-10000 5.0

**I** integral value, repeats/minute

**Range:** 0.00 to 9.99

**Factory Preset:** ThermoFlex1400-10000 0.10

**D** derivative value, minutes

**Range:** 0.0 to 9.9

**Factory Preset:** ThermoFlex1400-10000 0.0

## Section 8 Additional Information

### Draining

Before using any fluid or performing maintenance where contact with the fluid is likely refer to the manufacturer's SDS for handling precautions. ▲

Ensure the fluid is below the safe-handling temperature (below 40°C) before draining the chiller. ▲

The drain port is located on the rear. There are two different types, depending on the pump, a 1/4" stainless steel FPT with 1/4" brass pipe plug or a 1/4" MPT Riton fitting.

Position a suitable pan beneath the drain port. The drain pan must be shallow (under 3 1/2" in height) and have a volume of approximately 3 gallons (6 gallons for ThermoFlex7500 - 10000). Remove the 1/4" Male NPT pipe plug or, using a 9/16" wrench, open the Riton fitting by turning either counter clock wise. For ThermoFlex7500-10000, open the drain valve. This will drain the return line, reservoir, plate exchanger, and the suction side of the pump.

To drain the discharge side of the pump disconnect the Female NPT outlet connection on the rear of the chiller.

**Note** Internally the chiller does not contain a large quantity of fluid on the discharge side however take care to contain what fluid does drain, a wet-vac can be employed to minimize the potential for spillage. ▲

If the chiller is equipped with the flow control or pressure relief with flow control option, open the valve or remove the drain plug in order to drain the discharge line, see Section 5.


If the chiller is equipped with the anti drainback option, use the controller's **SETTINGS** menu to open the valve, see Section 4. Opening the valve allows the fluid to drain out of the chiller.

Reinstall 1/4" Male NPT pipe plug using a sealant suitable for the wetted materials or close the Riton filling prior to refilling the chiller.

**Do not overtighten the fitting.** ▲

For ThermoFlex7500-10000, close the drain valve.

**Figure 8-1** Drains

### Water-Cooled

Draining ThermoFlex1400 - 2500 water-cooled chillers is accomplished by removing the right side panel. Use a Phillips head screwdriver to remove the five screws indicated in the illustration below. Slide the panel back approximately one inch, then lift slightly from the rear to disengage the panel's two tabs from their slots.

**Figure 8-2** Water-Cooled

Install a  $\frac{7}{16}$ " ID tube on the drain petcock valve located on the lower end of the exchanger. Open the valve to allow fluid to drain into an external device. When draining is complete close the valve and replace the panel.


A wet-vac is needed on the facility water inlet connection to thoroughly drain any remaining fluid from the lines.

# CALIBRATION


The ThermoFlex has been designed to minimize the need for calibration. However, if calibration is desired or recommended by our Sales, Service and Customer Support, please use the following procedure.

**CALIBRATION** calibrates the chiller's temperature (**t1**), pressure (**p2**) and optional fluid flow (**flow1**) sensors. Each calibration requires a running chiller and a calibrated reference device. Typically, a 2-point calibration is used. Select which point to calibrate, **Low** or **High**. Press enter and enter the value as read by the reference device. Press enter. Select **Cal.** to the right of the value just entered by pressing the down button. Press enter and note that the **Current** value now equals the value just entered.


1. Highlight **CALIBRATION**.


2. Press enter and highlight **CALIBRATE**.

3. Press enter to display:

4. Highlight the desired sensor then press enter to display:

For **t1**, run the chiller to a suitable high-end calibration point. Place a calibrated reference thermometer in the reservoir. Ensure the fluid temperature is stabilized before performing the calibration. If it is more convenient, perform the low-end calibration before doing the high-end. Do not pick calibration points that are outside the safe operating limits of the fluid in your application. For example, with water, 40°C and 5°C are typical high and low calibration points (90°C and 5°C for high-temperature chillers).

For **p2**, connect a calibrated reference pressure gauge to the outlet line. Use an external flow control valve to adjust the pressure to suitable calibration points. Ensure the pressure is stabilized before calibrating.

For **flow1**, connect a calibrated reference flow meter to the outlet line. Use an external flow control valve to adjust the flow to suitable calibration points. Ensure the flow is stabilized before calibrating.

To restore a selected sensor to the factory values you have two options, **BACKUP** and **RESET**. **BACKUP** restores both the controller board and sensor calibration, **RESET** only restores the board calibration.

## **Wetted Materials**

### **P1, P2, MD1 and MD2 Pumps**

300 Series Stainless Steel

Bronze

Carbon Graphite

Ceramic

Fluorocarbon (Viton®)

Polysulfone

### **T0 and T1 Pumps**

Stainless Steel AISI 304

Bronze ASTM B62

Bronze ASTM B16

Buna N

Buna/Ceramic

Buna/Carbon

### **T5 Pumps**

Stainless Steel AISI 304

Bronze w/monel

Carbon

Buna N

Ceramic

### **Tank**

Polyethylene (standard temp chiller)

Polyvinylidene Difluoride (high temp)

Brass

EPDM

Pyrex®

Riton® (optional drain fitting)

Viton® (optional drain fitting o-ring)

### **P3, P4 and P5 Pumps**

316 Series Stainless Steel

Carbon

Silicon Carbide

Fluorocarbon (Viton®)

### **Plumbing**

300 Series Stainless Steel

Bronze

Fluorocarbon (Viton®)

Nickel

Polypropylene

EPDM

Brass

Copper

Teflon®

PPS (flow transducer)

Nitrile (Buna-n®)

### **Filter bag**

Polypropylene

Mono-filament nylon

### **Cap and Funnel**

Acetal Copolymer

## Decommissioning/ Disposal

Decommissioning prepares equipment for safe and secure transportation.

Laboratory Grade Ethylene glycol (EG) is poisonous and flammable. Before disposing refer to the manufacturer's most current SDS for handling precautions. ▲

Decommissioning must be performed only by qualified dealer using certified equipment. All prevailing regulations must be followed. ▲

Consider decommissioning the chiller when:

- It fails to maintain desired specifications
- It no longer meets safety standards
- It is beyond repair for its age and worth

Refrigerant and compressor oil must be recovered from equipment before disposal.

**Note** Keep in mind any impact your application may have had on the chiller. ▲

Direct questions about chiller decommissioning or disposal to our Sales, Service and Customer Support.

Handle and dispose in accordance with the manufacturers specification and/or the SDS for the material used. ▲

## Shipment/Storage

Follow the manufacturer's SDS instructions if decontamination is required. ▲

Transporting and/or storing the chiller requires draining, see Draining in this Section. Store the chiller in the temperature range of -25°C to 60°C (with packaging), and <80% relative humidity. ▲

If the chiller is stored for more than 90 days it must be flushed with clean fluid before operating. ▲

# Appendix A Country Specific

## 230 VAC, 50 Hz, 1Ø Requirements

Refer to the nameplate label located on the rear of the chiller for specific electrical requirements.

### 1. Chillers shipped to the following locations require a **16 Amp service**:

Afghanistan, Albania, Algeria, Andorra, Angola, Argentina, Armenia, Austria, Azerbaijan, Belarus, Belgium, Benin, Bolivia, Bosnia and Herzegovina, Brazil, Bulgaria, Burkina Faso, Burundi, Cambodia, Cameroon, Cape Verde, Central African Republic, Chad, Chile, Comoros, Congo, Croatia, Czech Republic, Denmark, Djibouti, DR Congo, Ecuador, Egypt, Eritrea, Estonia, Ethiopia, Finland, France, French Guiana, Gabon, Georgia, Germany, Greece, Guinea, Hungary, Iceland, Indonesia, Iran, Iraq, Israel, Italy, Ivory Coast, Jordan, Kazakhstan, Kyrgyzstan, Latvia, Lebanon, Liberia, Libya, Liechtenstein, Lithuania, Luxembourg, Madagascar, Mali, Mauritania, Moldova, Monaco, Mongolia, Morocco, Mozambique, Namibia, Nepal, Netherlands, Niger, North Korea, Norway, Paraguay, Peru, Poland, Portugal, Romania, Russia, Rwanda, Saint Vincent and the Grenadines, San Marino, Sao Tome and Principe, Saudi Arabia, Senegal, Serbia, Slovakia, Slovenia, Somalia, South Africa, South Korea, Spain, Sweden, Switzerland, Syria, Tajikistan, Thailand, Togo, Tunisia, Turkey, Turkmenistan, Ukraine, Uruguay, Uzbekistan, Vanuatu, Vatican City, Vietnam.

---

### 2. Chillers shipped to the following locations require a **15 Amp service**:

Australia, China, Fiji Islands, Nauru, New Zealand, Papua New Guinea, Solomon Island, Tonga, Tuvalu.

---

### 3. Chillers shipped to the following locations require a **13 Amp service**:

Abu Dhabi, Bahrain, Bangladesh, Botswana, Brunei, Cyprus, Dominica, Gambia, Ghana, Gibraltar, Grenada, Hong Kong, India, Ireland, Kenya, Kiribati, Kuwait, Lesotho, Malawi, Malaysia, Maldives, Malta, Mauritius, Myanmar, Nigeria, Oman, Pakistan, Qatar, Saint Lucia, Seychelles, Sierra Leone, Singapore, Sri Lanka, Sudan, Swaziland, Tanzania, Uganda, United Arab Emirates, United Kingdom, Yemen, Zambia, Zimbabwe.

---


# Appendix B Voltage Configuration Instructions

ThermoFlex 900 and 1400 chillers equipped with the 115V 60Hz, 100v 50/60Hz Variable Voltage option and ThermoFlex 900 to 5000 chillers equipped with 200-230V 50/60Hz Global Voltage option have a voltage configuration panel located on the rear of the chiller behind an access panel, see Figure B-1.

- Use a 1/4" socket to remove the four screws securing the access panel to the chiller.
- The configuration panel has two 3-position toggle switches, one for voltage and one for frequency. All chillers are shipped with the toggle switch in the center **SHIP** position. Place each switch to the settings that match the voltage/frequency supplied to the chiller.


**Note** For ThermoFlex900-2500 global voltage chillers, the compressor and fan will not operate when the switch is in the **SHIP** position. ▲

- Reinstall the access panel.

**Figure B-1** Variable/Global Voltage Chillers

# Appendix C NC Serial Communications Protocol

**Note** Appendix C assumes you have a basic understanding of communications protocols. ▲

Connect your PC to the applicable connector on the rear of the chiller. Use the controller, see Section 4, to enable serial communications.

**Note** Keypad operation is still available with serial communications enabled. ▲

**WARNING** Never apply line voltage to any of the connections. ▲

**Figure C-1** Connectors

All data is sent and received in binary form, do not use ASCII. In the following pages the binary data is represented in hexadecimal (hex) format.

The NC Serial Communications Protocol is based on a master-slave model. The master is a host computer, while the slave is the chiller's controller. Only the master can initiate a communications transaction (half-duplex). The slave ends the transaction by responding to the master's query. The protocol uses RS-232/RS-485 serial interface with the default parameters: 9600 baud, 8 data bits, 1 stop bit, and no parity. RS-485 offers a slave address selection, default parameter: 1.

The chiller can be controlled through your computer's serial port by using the chiller's standard female 9-pin connection.

### RS-232 COMM


Pin #	Function
1	No connection
2	TX
3	RX
4	No connection
5	GND = Signal ground
6 - 9	No connection

TX = Transmitted data from controller  
RX = Received data to controller.

### RS-485 COMM

Pin #	Function
1-7	No connection
8	T+
9	T-

**Hardware Mating Connector**  
**AMP Part# 745492-2 or equivalent**


Communication cables are available from Thermo Fisher. Contact us for additional information.

All commands must be entered in the exact format shown in the tables on the following pages. The tables show all commands available, their format and responses. Controller responses are either the requested data or an error message. The controller response *must* be received before the host sends the next command.

The host sends a command embedded in a single communications packet, then waits for the controller's response. If the command is not understood or the checksums do not agree, the controller responds with an error command. Otherwise, the controller responds with the requested data. If the controller fails to respond within 1 second, the host should resend the command.

**Note** All byte values are shown in hex, hex represents the binary values that must be sent to the chiller. **Do not use ASCII.** ▲

The framing of the communications packet in both directions is:

<i>Lead char</i>	<b>0xCA (RS-232) 0xCC (RS-485)</b>
	Device address is 1 (RS-232)
<i>Addr-msb</i>	Most significant byte of slave address (RS-232: 0)
<i>Addr-lsb</i>	Least significant byte of slave address (RS-232: 1)
<i>Command</i>	Command byte (see Table of Commands)
<i>n d-bytes</i>	Number of data bytes to follow
<i>d-byte 1</i>	1 <sup>st</sup> data byte (the qualifier byte is considered a data byte)
...	...
<i>d-byte n</i>	n <sup>th</sup> data byte.
<i>Checksum</i>	Bitwise inversion of the 1 byte sum of bytes beginning with the most significant address byte and ending with the byte preceding the checksum. (To perform a bitwise inversion, "exclusive OR" the one byte sum with FF hex.)

When a command has no value associated with it, (e.g., REQ ACK), "n d-bytes" will be set to 0. Values such as temperature and flow are sent as either 2 or 4 byte signed integers, depending on how they are stored in the controller RAM.

When the controller sends a value, a qualifier byte is sent first, followed by a 2 or 4 byte integer (the least significant byte is sent last). The qualifier indicates the precision and units of the value. The host does not send the qualifier byte; it must send the value using the correct precision, units and number of bytes. The host first inquires about a value it wants to change, then uses the number of data bytes and the qualifier byte it receives to generate the proper integer to send.

**Analog Values**

*Qualifier Byte	
b.7	Precision of measurement
b.6	
b.5	
b.4	
b.3	Unit of measure index
b.2	
b.1	
b.0	

Unit of Measure	
Index	Unit
0	NONE
1	Temperature in °C
2	Temperature in °F
3	Flow liters per minute
4	Flow in gallons per minute
5	Time in seconds
6	Pressure in PSI
7	Pressure in bars
8	Resistivity in MΩ-cm
9	%
10	Volts
11	Pressure in K Pascals
12	Conductivity in μs/cm

E.g., the integer 986 preceded by a qualifier byte of 0x12 is 98.6°F.

Example to set setpoint to 25°C:

If the temperature units are unknown, before changing the setpoint send a command to request setpoint. The response will include both the precision and units. Precision is fixed at 0.1 and units can be either °C or °F. If the units are already known skip to step 3.

1. Master sends: CA 00 01 70 00 8E (REQ SETPOINT1)
2. Slave responds: CA 00 01 70 03 11 00 C8 B2 (0.1°C x 200)  
Response indicates:  
uses a 2 byte integer (nn=03)  
precision and units are 0.1°C (d1=11)
3. Master sends: CA 00 01 F0 02 00 FA 12 (Set Setpoint 1 to 25.0°C)
4. Slave responds: CA 00 01 F0 03 11 00 FA 00 (0.1°C x 250)

See Additional Command Examples in this Appendix.

**Set Commands** – When a Set Command is received and the value is within the allowable limits, the new value will take affect immediately and the new value will be returned as part of the response to the command. If the value is outside of the allowable limits, the value will be rejected and the old value will be returned as part of the response to the command.

**Error Response F0** - The "Bad Data" and "Bad Checksum" error responses will not be used. Reject out of range values and return old settings. Do not reply to message frames with bad checksums. This behavior is compatible with existing implementations of the protocol.

## Table of Commands

Command	M: Master Sends S: Slave Responds	Notes
REQUEST STATUS		
REQ ACK	M: lc a1 a2 <b>00</b> 00 cs S: lc a1 a2 <b>00</b> 02 v1 v2 cs	protocol version v1=0; v2=1
REQ CONTROLLER SW VER	M: lc a1 a2 <b>02</b> 01 d1 cs S: lc a1 a2 <b>02</b> nn d1 ... dn cs	d1 = 0 Controller SW version in ASCII d1 = 1 Controller SW checksum
REQ STATUS	M: lc a1 a2 <b>09</b> 00 cs S: lc a1 a2 <b>09</b> nn d1 ... dn cs	see Request Status Table in this Appendix
ERROR	M: S: lc a1 a2 <b>0F</b> 02 en ed cs	Response Only! ed = Error Data en = Error Number 1: Bad Command See Error in this Appendix
REQUEST MEASUREMENTS		
REQ SETPT1	M: lc a1 a2 <b>70</b> 00 cs S: lc a1 a2 <b>70</b> 03 d1 d2 d3 cs	Process Fluid Setpoint
REQ FLOW1	M: lc a1 a2 <b>10</b> 00 cs S: lc a1 a2 <b>10</b> 03 d1 d2 d3 cs	Process Fluid Flow
REQ TEMP1	M: lc a1 a2 <b>20</b> 00 cs S: lc a1 a2 <b>20</b> 03 d1 d2 d3 cs	Process Fluid Supply Temperature (RTD1)
REQ TEMP2	M: lc a1 a2 <b>21</b> 00 cs S: lc a1 a2 <b>21</b> 03 d1 d2 d3 cs	Process Fluid Return Temperature (RTD2)
REQ TEMP4	M: lc a1 a2 <b>23</b> 00 cs S: lc a1 a2 <b>23</b> 03 d1 d2 d3 cs	Entering Air/Facility Water (RTD4)
REQ TEMP7	M: lc a1 a2 <b>26</b> 00 cs S: lc a1 a2 <b>26</b> 03 d1 d2 d3 cs	ThermoFlex 2500 Air-cooled Fan Speed
REQ ANALOG1	M: lc a1 a2 <b>28</b> 00 cs S: lc a1 a2 <b>28</b> 03 d1 d2 d3 cs	Process Fluid Supply Pressure (P2)
REQ ANALOG2	M: lc a1 a2 <b>29</b> 00 cs S: lc a1 a2 <b>29</b> 03 d1 d2 d3 cs	Refrigeration Suction Pressure (P5)
REQ ANALOG3	M: lc a1 a2 <b>2A</b> 00 cs S: lc a1 a2 <b>2A</b> 03 d1 d2 d3 cs	Process Fluid Return Pressure (P1)
REQ ANALOG4	M: lc a1 a2 <b>2B</b> 00 cs S: lc a1 a2 <b>2B</b> 03 d1 d2 d3 cs	Condensing Pressure (P6)

REQ ANALOG5	M: lc a1 a2 <b>1C</b> 00 cs S: lc a1 a2 <b>1C</b> 03 d1 d2 d3 cs	Facility Inlet Pressure (P7)
REQ ANALOG6	M: lc a1 a2 <b>1D</b> 00 cs S: lc a1 a2 <b>1D</b> 03 d1 d2 d3 cs	Facility Outlet Pressure (P8)
REQ ANALOG7	M: lc a1 a2 <b>1E</b> 00 cs S: lc a1 a2 <b>1E</b> 03 d1 d2 d3 cs	Analog Level (LEV4)
REQ ANALOG9	M: lc a1 a2 <b>2F</b> 00 cs S: lc a1 a2 <b>2F</b> 03 d1 d2 d3 cs	+5V Sense
REQ REMOTE RTD	M: lc a1 a2 <b>1B</b> 00 cs S: lc a1 a2 <b>1B</b> 03 d1 d2 d3 cs	Remote Temperature from Analog Board
REQ RES1	M: lc a1 a2 <b>2C</b> 00 cs S: lc a1 a2 <b>2C</b> 03 d1 d2 d3 cs	Process Fluid Resistivity
REQUEST LOW ALARM VALUES		
REQ LO FLOW1	M: lc a1 a2 <b>30</b> 00 cs S: lc a1 a2 <b>30</b> 03 d1 d2 d3 cs	Process Warning
REQ LO FLOW3	M: lc a1 a2 <b>32</b> 00 cs S: lc a1 a2 <b>32</b> 03 d1 d2 d3 cs	Process Fault
REQ LO ANALOG1	M: lc a1 a2 <b>48</b> 00 cs S: lc a1 a2 <b>48</b> 03 d1 d2 d3 cs	Pressure Process Supply Warning
REQ LO ANALOG2	M: lc a1 a2 <b>49</b> 00 cs S: lc a1 a2 <b>49</b> 03 d1 d2 d3 cs	Pressure Process Supply Fault
REQ LO ANALOG7	M: lc a1 a2 <b>3E</b> 00 cs S: lc a1 a2 <b>3E</b> 03 d1 d2 d3 cs	Level Warning
REQ LO ANALOG8	M: lc a1 a2 <b>3F</b> 00 cs S: lc a1 a2 <b>3F</b> 03 d1 d2 d3 cs	Level Fault
REQ LO TEMP1	M: lc a1 a2 <b>40</b> 00 cs S: lc a1 a2 <b>40</b> 03 d1 d2 d3 cs	Process Warning
REQ LO TEMP2	M: lc a1 a2 <b>41</b> 00 cs S: lc a1 a2 <b>41</b> 03 d1 d2 d3 cs	Process Fault
REQ AUTO REFILL ON	M: lc a1 a2 <b>45</b> 00 cs S: lc a1 a2 <b>45</b> 03 d1 d2 d3 cs	Auto refill On Setting
REQ LO RES1	M: lc a1 a2 <b>4C</b> 00 cs S: lc a1 a2 <b>4C</b> 03 d1 d2 d3 cs	Process Warning

## REQUEST HIGH ALARM VALUES

REQ HI FLOW1	M: lc a1 a2 <b>50</b> 00 cs S: lc a1 a2 <b>50</b> 03 d1 d2 d3 cs	Process Warning
REQ HI FLOW3	M: lc a1 a2 <b>52</b> 00 cs S: lc a1 a2 <b>52</b> 03 d1 d2 d3 cs	Process Fault
REQ HI TEMP1	M: lc a1 a2 <b>60</b> 00 cs S: lc a1 a2 <b>60</b> 03 d1 d2 d3 cs	Process Warning
REQ HI TEMP2	M: lc a1 a2 <b>61</b> 00 cs S: lc a1 a2 <b>61</b> 03 d1 d2 d3 cs	Process Fault
REQ HI ANALOG1	M: lc a1 a2 <b>68</b> 00 cs S: lc a1 a2 <b>68</b> 03 d1 d2 d3 cs	Pressure Process Supply Warning
REQ HI ANALOG2	M: lc a1 a2 <b>69</b> 00 cs S: lc a1 a2 <b>69</b> 03 d1 d2 d3 cs	Pressure Process Supply Fault
REQ HI RES1	M: lc a1 a2 <b>6C</b> 00 cs S: lc a1 a2 <b>6C</b> 03 d1 d2 d3 cs	Process Warning

## REQUEST PID SETTINGS

REQ COOL P	M: lc a1 a2 <b>74</b> 00 cs S: lc a1 a2 <b>74</b> 03 d1 d2 d3 cs
REQ COOL I	M: lc a1 a2 <b>75</b> 00 cs S: lc a1 a2 <b>75</b> 03 d1 d2 d3 cs
REQ COOL D	M: lc a1 a2 <b>76</b> 00 cs S: lc a1 a2 <b>76</b> 03 d1 d2 d3 cs
REQ HEAT P	M: lc a1 a2 <b>71</b> 00 cs S: lc a1 a2 <b>71</b> 03 d1 d2 d3 cs
REQ HEAT I	M: lc a1 a2 <b>72</b> 00 cs S: lc a1 a2 <b>72</b> 03 d1 d2 d3 cs
REQ HEAT D	M: lc a1 a2 <b>73</b> 00 cs S: lc a1 a2 <b>73</b> 03 d1 d2 d3 cs

## SET STATUS SETTINGS

SET KEYSTROKE	M: lc a1 a2 <b>80</b> 01 d1 cs S: lc a1 a2 <b>80</b> 01 d1 cs	See Set Keystroke in this Appendix
SET ON/OFF ARRAY	M: lc a1 a2 <b>81</b> nn d1 ... dn cs S: lc a1 a2 <b>81</b> nn d1 ... dn cs	See Set On/Off Array in this Appendix di: 0 = OFF, 1 = ON, 2 = no change

## SET MEASUREMENT

SET SETPT1	M: lc a1 a2 <b>F0</b> 02 d1 d2 cs S: lc a1 a2 <b>F0</b> 03 d1 d2 d3 cs	Process Fluid Setpoint
------------	---	------------------------

## SET LOW ALARM VALUES

SET LO FLOW1	M: lc a1 a2 <b>B0</b> 02 d1 d2 cs S: lc a1 a2 <b>B0</b> 03 d1 d2 d3 cs	Process Warning
SET LO FLOW3	M: lc a1 a2 <b>B2</b> 02 d1 d2 cs S: lc a1 a2 <b>B2</b> 03 d1 d2 d3 cs	Process Fault
SET LO TEMP1	M: lc a1 a2 <b>C0</b> 02 d1 d2 cs S: lc a1 a2 <b>C0</b> 03 d1 d2 d3 cs	Process Warning
SET LO TEMP2	M: lc a1 a2 <b>C1</b> 02 d1 d2 cs S: lc a1 a2 <b>C1</b> 03 d1 d2 d3 cs	Process Fault
SET LO ANALOG1	M: lc a1 a2 <b>C8</b> 02 d1 d2 cs S: lc a1 a2 <b>C8</b> 03 d1 d2 d3 cs	Pressure Process Supply Warning
SET LO ANALOG2	M: lc a1 a2 <b>C9</b> 02 d1 d2 cs S: lc a1 a2 <b>C9</b> 03 d1 d2 d3 cs	Pressure Process Supply Fault
SET LO ANALOG7	M: lc a1 a2 <b>BE</b> 02 d1 d2 cs S: lc a1 a2 <b>BE</b> 03 d1 d2 d3 cs	Level Warning
SET LO ANALOG8	M: lc a1 a2 <b>BF</b> 02 d1 d2 cs S: lc a1 a2 <b>BF</b> 03 d1 d2 d3 cs	Level Fault
SET AUTO REFILL ON	M: lc a1 a2 <b>C5</b> 02 d1 d2 cs S: lc a1 a2 <b>C5</b> 03 d1 d2 d3 cs	When level % drops below this, turn on auto refill
SET LO RES1	M: lc a1 a2 <b>CC</b> 02 d1 d2 cs S: lc a1 a2 <b>CC</b> 03 d1 d2 d3 cs	Process Warning
SET LO RES1	M: lc a1 a2 <b>CD</b> 02 d1 d2 cs S: lc a1 a2 <b>CD</b> 03 d1 d2 d3 cs	Process Fault

## SET HIGH ALARM VALUES

SET HI FLOW1	M: lc a1 a2 <b>D0</b> 02 d1 d2 cs S: lc a1 a2 <b>D0</b> 03 d1 d2 d3 cs	Process Warning
SET HI FLOW3	M: lc a1 a2 <b>D2</b> 02 d1 d2 cs S: lc a1 a2 <b>D2</b> 03 d1 d2 d3 cs	Process Fault
SET HI TEMP1	M: lc a1 a2 <b>E0</b> 02 d1 d2 cs S: lc a1 a2 <b>E0</b> 03 d1 d2 d3 cs	Process Warning
SET HI TEMP2	M: lc a1 a2 <b>E1</b> 02 d1 d2 cs S: lc a1 a2 <b>E1</b> 03 d1 d2 d3 cs	Process Fault
SET HI ANALOG1	M: lc a1 a2 <b>E8</b> 02 d1 d2 cs S: lc a1 a2 <b>E8</b> 03 d1 d2 d3 cs	Pressure Process Supply Warning
SET HI ANALOG2	M: lc a1 a2 <b>E9</b> 02 d1 d2 cs S: lc a1 a2 <b>E9</b> 03 d1 d2 d3 cs	Pressure Process Supply Fault
SET HI ANALOG7	M: lc a1 a2 <b>DE</b> 02 d1 d2 cs S: lc a1 a2 <b>DE</b> 03 d1 d2 d3 cs	Level Warning
SET HI ANALOG8	M: lc a1 a2 <b>DF</b> 02 d1 d2 cs S: lc a1 a2 <b>DF</b> 03 d1 d2 d3 cs	Level Fault
SET AUTO REFILL OFF	M: lc a1 a2 <b>E5</b> 02 d1 d2 cs S: lc a1 a2 <b>E5</b> 03 d1 d2 d3 cs	When level % drops below this, turn off auto refill
SET HI RES1	M: lc a1 a2 <b>EC</b> 02 d1 d2 cs S: lc a1 a2 <b>EC</b> 03 d1 d2 d3 cs	Process Warning
SET HI RES2	M: lc a1 a2 <b>ED</b> 02 d1 d2 cs S: lc a1 a2 <b>ED</b> 03 d1 d2 d3 cs	Process Fault
SET PID SETTINGS		
SET COOL P	M: lc a1 a2 <b>F4</b> 02 d1 d2 cs S: lc a1 a2 <b>F4</b> 03 d1 d2 d3 cs	Cool P Term
SET COOL I	M: lc a1 a2 <b>F5</b> 02 d1 d2 cs S: lc a1 a2 <b>F5</b> 03 d1 d2 d3 cs	Cool I Term
SET COOL D	M: lc a1 a2 <b>F6</b> 02 d1 d2 cs S: lc a1 a2 <b>F6</b> 03 d1 d2 d3 cs	Cool D Term
SET HEAT P	M: lc a1 a2 <b>F1</b> 02 d1 d2 cs S: lc a1 a2 <b>F1</b> 03 d1 d2 d3 cs	Heat P Term
SET HEAT I	M: lc a1 a2 <b>F2</b> 02 d1 d2 cs S: lc a1 a2 <b>F2</b> 03 d1 d2 d3 cs	Heat I Term
SET HEAT D	M: lc a1 a2 <b>F3</b> 02 d1 d2 cs S: lc a1 a2 <b>F3</b> 03 d1 d2 d3 cs	Heat D Term

## Request Status Table

nn	4		
	b0	Chiller Running	b0 HTC (High Temperature Cutout)
	b1	Chiller Faulted	b1 LLC (Low Level Cutout)
	b2	Process Supply RTD open or shorted	b2 MOL (Motor Overload)
d1	b3	Process Return RTD open or shorted	d3 b3 Phase Monitor
	b4	Suction RTD open or shorted	b4 HPC (High Pressure Cutout)
	b5	Entering Air or Facility Water RTD open or shorted	b5 LPC (Low Pressure Cutout)
	b6	High Temp Error	b6 EMO
	b7	Low Temp Error	b7 External EMO
	b0	High Pressure Error	b0 RA T_MAX (High Temperature)
	b1	Low Pressure Error	b1 Not used
	b2	High Flow Error (user set able)	b2 Auto Refill Valve Open
d2	b3	Low Flow Error (user set able)	d4 b3 Anti Drainback Valve Open
	b4	High Level Error	b4 Clogged Fluid Filter Fault
	b5	Low Level Error	b5 Temp Fault Startup Bypass
	b6	Drip Pan fault	b6 System Low Flow
	b7	Auto Refill fault	b7 Not used

**Error**

The slave detected an error in the message it received from the master, so it returns this command instead of echoing the command sent by the master. The slave returns the command it received from the master in the ed byte, and an error code in the en byte.

en	Error
1	Bad command – not recognized by slave
2	Reject value and return old setting
3	Do not respond at all

Some errors may not result in any response. The slave ignores incoming bytes until it sees the valid lead character and its slave address. Then it must receive the correct number of bytes (determined by the length byte) before it can respond. If an incomplete frame is received, the slave will time out and clear its input buffer without responding.

**Set On/Off Array**

This command is used to set the state of the chiller on or off. Sending a 0 in the array turns off the chiller while sending a 1 turns it on. Sending a 2 does not change the state of the feature. The array is returned showing the state of each feature after the command has been carried out. Sending a 2 effectively turns this command into a request status command.

nn	1
d1	Chiller On/Off

**Set Keystroke**

This command is used to affect a keystroke remotely as if someone pressed the button on the controller.

Value	
0	Null
1	Enter
2	Up/Yes
3	Down/No
4	Esc
5	On/Off
6	Left
7	Right

## Set Special Command

Used for product specific commands, in the ThermoFlex this command configures the analog options (DAC out, Analog in and turns on and off DAC output)

CA 00 01 8D nn d1...dn cs

### Set Special Command

Byte		Notes
d1	Command Byte	Indicates what command
d2	Description Byte	See description below
d3	Entered Value if necessary	

### Command Byte

Command	Description
00	Set Analog Option
01	Set Resistivity Control (for basic derivatives)
02	Unused

### Description Byte

Byte	Notes
b.7 - b.6	Unused
b.5 and b.4	Remote DAC output enable      00 = Off, 01 = On, 1x = No Change
b.3 - b.2	Configure DAC Output          00h = Volt, 01h = millivolt, 02h = milliamp, 03h = No Change
b.1 - b.2	Configure Analog Input        00h = Volt, 01h = millivolt, 02h = milliamp, 03h = No Change

The following command string example enables the DAC output, configures the DAC for millivolts and configures the analog input for voltage:

CA 00 01 8D 02 00 14 cs

# Appendix D Analog I/O and Remote Sensor

## Analog I/O Connector Pinout ♀

Install your analog input/output device to the 15-pin female connector on the rear of the chiller. Analog I/O is activated using the controller, see Section 4.

PIN	NAME	NOTES	DEFINITION
1	DIGITAL GROUND		Common round connection for pins 12, 13 and 14
2	RESISTIVITY OUT	Optional	Analog voltage output 0v = 0 MΩ - cm 10v = 20 MΩ - cm
3	LOW LEVEL (Only if option chosen)	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Closes if either level switch is in the "low" position for more than 1 second.
4	CONFIGURABLE RELAY 2	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Closes when any configured fault or warning occurs.
5	PUMP ON	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Closes when pump is turned on. Opens when pump is turned off.
6	ANALOG GROUND		Common for analog signals (pins 2, 7 and 15)
7	RESERVOIR TEMP OUT <b>OR</b> EXTERNAL SENSOR TEMPERATURE IF EXTERNAL SENSOR ENABLED	Note 2	Analog Voltage Output 0-10VDC, 10mV/°C, or 4-20mA: Reference to pin 6. This voltage output is proportional to the reservoir fluid temperature: Default scale= 0–10V (where: 0V = Low Temp Span, 10V = Hi Temp Span) Optional Range = 10mV/ °C. (Ex: 200mV = 20°C) (Max Load @ 10V = 5mA) or 4-20mA, 4mA = low temp span, 20 mA = high temp span (maximum output current = 5mA @10VDC.
8	LOW FLOW (Only if option chosen)	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Closes when a low flow occurs while the pump is on. Note: To allow the pump to get up to speed at startup, the pump runs for 3 - 5 seconds before the low flow sensor is read.
9	CONFIGURABLE RELAY 1 (Normally Open)	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Closes when any of the configured faults occur.
10	CONFIGURABLE RELAY 1 (Normally Closed)	Note 1	<u>Dry Relay Contact</u> : Reference to pin 11. Complement of pin 9 (open when pin 9 is closed).
11	RELAY COMMON		Common for all relay contacts (pins 3, 4, 5, 8, 9, 10).
12	REMOTE START ENABLE	Note 3	Connect to pin 1 to allow chiller to be remotely turned on/off through pin 14 REMOTE START.

Note 1: All relay contacts (except for Pin 10) are normally OPEN when power is off. Pin 10 contacts are normally CLOSED when power is off. Relay contacts are rated: 24V AC/DC, 2A, <= 0.08 Ohm maximum each or 5A total for all relays combined, 1mA minimum, switching capacity: 48VA/48W (Resistive load only).

Note 2: Default = 0-10VDC. These ranges are set by the user.

Note 3: Connect to digital ground (pin 1) using a low resistance connection (gold contact relay).


PIN	NAME	NOTES	DEFINITION
13	REMOTE SETPOINT ENABLE	Note 3	Connect to pin 1 to allow the setpoint to be changed remotely through pin 15 REMOTE SETPOINT.
14	REMOTE START	Note 3	Connect to pin 1 to turn chiller on. Disconnect to turn chiller off. Note: Pins 1 and 12 must be connected to allow operation from this pin.
15	REMOTE SETPOINT	Note 2, 4	Analog Voltage Input 0-10VDC, 10mV/°C, or 4-20mA: Reference to pin 6. Apply a DC voltage to this pin to adjust the setpoint: Default Range = 0 – 10V (where: 0V = Low Temp Span, 10V = Hi Temp Span) (Input Impedance > 600K) Optional Range = 10mV/ °C. (Ex: 200mV = 20°C) (Max Input Voltage = 10VDC, or 4-20mA, 4mA = low temp span, 20 mA = high temp span.

Note 1: All relay contacts (except for Pin 10) are normally OPEN when power is off. Pin 10 contacts are normally CLOSED when power is off. Relay contacts are rated: 24V AC/DC, 2A, <= 0.08 Ohm maximum each or 5A total for all relays combined, 1mA minimum, switching capacity: 48VA/48W (Resistive load only).

Note 2: Default = 0-10VDC. These ranges are set by the user.

Note 3: Connect to digital ground (pin 1) using a low resistance connection (gold contact relay).


Note 4: Remote setpoint must be enabled, pin 13

**WARNING** Never apply line voltage to any of the connections. ▲

When making your connection to the ThermoFlex Analog I/O connector, in order to comply with the EMC directive:

- Use a shielded I/O cable
- Connect the remote end of the cable shield to earth ground.
- Connect cable shield to ThermoFlex end connector.

**Figure D-1** Analog I/O Connector

## Analog Comm and Configurable Relays

Use the up/down arrows to scroll through the entire menu (only six lines are displayed at a time). When any desired feature is highlighted press enter to select it. Enabled features are indicated by .

ANALOG COMM - ACOM	
<input checked="" type="checkbox"/> Remote Sensor	
<input type="checkbox"/> Remote Start	
<input type="checkbox"/> Remote Setpt	volts
<input type="checkbox"/> Temp Out	volts
⏠	
-----	
Remote Start	
Remote Setpt	volts

**Note** Enabling analog I/O remote start/stop disables the chiller's local controller start/stop capability. Enabling analog I/O remote also overrides serial communications start/stop commands. ▲

To change the **Remote Setpoint** or **Temp Out** input type use the up and down arrows until the input is highlighted. (The left and right arrows have no affect in moving the highlight.)

ANALOG COMM - ACOM	
<input type="checkbox"/> Remote Sensor	
<input type="checkbox"/> Remote Start	
<input type="checkbox"/> Remote Setpt	volts
<input type="checkbox"/> Temp Out	volts
⏠	

Press enter to get the input type to flash. Use the up and down arrows to bring up the desired input: **volts**, **mv** or **ma**. With the desired input displayed press enter again to accept the change.

With **RELAY1** or **RELAY2** highlighted press enter to bring up the corresponding **RELAY** display. Highlight the desired error(s) and press enter. Enabled errors are indicated by are indicated by .


**Note** Except for **Warning**, errors will automatically clear once the condition is corrected. **Warning** also requires pressing enter once the condition is corrected.

RELAY 1	
<input checked="" type="checkbox"/> Low Level	
<input type="checkbox"/> Tank Overflow	
<input type="checkbox"/> Drip Pan Full	
<input type="checkbox"/> Low Temperature	
⏠	
-----	
<input type="checkbox"/> High Temperature	
<input checked="" type="checkbox"/> Low Flow	
<input type="checkbox"/> High Flow	
<input type="checkbox"/> Low Resistivity	
<input type="checkbox"/> High Resistivity	
<input type="checkbox"/> High Pressure	
<input type="checkbox"/> Low Pressure	
<input checked="" type="checkbox"/> Unit Fault	
<input type="checkbox"/> Pump/Unit Shut Off	
<input type="checkbox"/> Refrig. Shut Off	
<input checked="" type="checkbox"/> Limit Fault	
<input type="checkbox"/> Sensor Fault	

RELAY 2	
<input type="checkbox"/> Low Level	
<input type="checkbox"/> Tank Overflow	
<input type="checkbox"/> Drip Pan Full	
<input type="checkbox"/> Auto Refill Error	
⏠	
-----	
<input checked="" type="checkbox"/> Low Temperature	
<input checked="" type="checkbox"/> High Temperature	
<input type="checkbox"/> Low Flow	
<input type="checkbox"/> High Flow	
<input type="checkbox"/> Low Resistivity	
<input checked="" type="checkbox"/> High Resistivity	
<input type="checkbox"/> High Pressure	
<input type="checkbox"/> Low Pressure	
<input type="checkbox"/> Warning	
<input type="checkbox"/> PM Timer	
<input type="checkbox"/> Comm Error	
<input checked="" type="checkbox"/> Sensor Fault	

## USER RANGE SETTINGS

This menu allows the you to customize the scaling of the analog input and DAC output to meet the application needs. With **USER RANGE SETTINGS** highlighted press enter to display:

The ThermoFlex supports three standard analog interface types and defaults to the following scaling:

<b>Voltage</b>	0v – 10v = the operating range of the chiller (5° – 40° in standard chillers, 5° – 90° in high temp chillers)
<b>MilliVolts</b>	10mv/°C
<b>Current</b>	4ma – 20ma = the operating range of the chiller

The following values reflect the standard scaling and are the factory defaults:

**Voltage (0v – 10v = the operating range of the chiller)**

<b>High volts</b>	<b>10.0v</b>
<b>Low volts</b>	<b>0.0v</b>
<b>High temp</b>	<b>40.0°</b> ( <i>standard temp chillers</i> ), <b>90.0°</b> ( <i>high temp chillers</i> )
<b>Low temp</b>	<b>5.0°</b>

With **Voltage** highlighted press enter to display:

<b>USER VOLTAGE</b>	
◀ <b>High Volts</b>	<b>10.0</b>
<b>Low Volts</b>	<b>0.0</b>
<b>High Temp</b>	<b>40.0</b>
<b>Low Temp</b>	<b>5.0</b>

Highlight the desired value and then press enter, the high-light will flash. Use the up and down arrows to change the value. Once the desired value is displayed press enter to accept the change and stop the flashing.

These values are used to calculate a gain and offset and are applied to the setpoint input voltage.

$$\text{VoltageInputGain} = (\text{high temp} - \text{low temp}) / (\text{high volt} - \text{low volt})$$

$$\text{VoltageInputOffset} = \text{high temp} - (\text{high volt} * \text{VoltageInputGain})$$

The default settings for a high temp chiller will yield the following gain and offset for setpoint input:

$$\text{VoltageInputGain} = (90 - 5) / (10 - 0) = 8.5$$

$$\text{VoltageInputOffset} = 90 - (10 * 8.5) = 5$$

Now 10v will yield a setpoint of 90°.  $(10 * \text{gain}) + \text{offset}$

$$(10v * 8.5) + 5 = 90^\circ$$

It could just as easily be configured for other ranges such as:

1v – 9v = 5° - 150°. However the chiller will only recognize a set point between 5° - 90°.

A separate gain and offset is also be calculated and are applied to the DAC temperature output value.

$$\text{VoltageOutputGain} = (\text{high volt} - \text{low volt}) / (\text{high temp} - \text{low temp})$$

$$\text{VoltageOutputOffset} = \text{high volt} - (\text{high temp} * \text{VoltageOutputGain})$$

The default settings for a high temp chiller will yield the following gain and offset for temperature output:

$$\text{VoltageOutputGain} = (10 - 0) / (90 - 5) = 0.117647$$

$$\text{VoltageOutputOffset} = 10 - (90 * 0.117647) = -0.58823$$

Now 90° will output 10v.  $(90 * \text{gain}) + \text{offset}$

$$(90^\circ * 0.117647) + (-0.58823) = 10.0v$$

The following values are the factory defaults for the millivolt setting that will effectively yield 10mv/°C:

**MilliVolts (10mv/°C)**

<b>High mv</b>	<b>1000.0mv</b>
<b>Low mv</b>	<b>0.0mv</b>
<b>High temp</b>	<b>100.0°</b>
<b>Low temp</b>	<b>0.0°</b>

With **Millivolts** highlighted press enter to display:

<b>USER MV</b>	
<b>High MV</b>	<b>1000.0</b>
<b>Low MV</b>	<b>0.0</b>
<b>High Temp</b>	<b>100.0</b>
<b>Low Temp</b>	<b>0.0</b>

Highlight the desired value and then press enter, the highlight will flash. Use the up and down arrows to change the value. Once the desired value is displayed press enter to accept the change and stop the flashing.

These values are used to calculate a gain and offset and are applied to the set point input voltage.

$$\begin{aligned}\text{MilliVoltInputGain} &= (\text{high temp} - \text{low temp}) / (\text{high mv} - \text{low mv}) \\ \text{MilliVoltInputOffset} &= \text{high temp} - (\text{high mv} * \text{MilliVoltInputGain})\end{aligned}$$

So that the default settings will yield the following gain and offset for set point input:

$$\begin{aligned}\text{MilliVoltInputGain} &= (100 - 0) / (1000 - 0) = 0.1 \\ \text{MilliVoltInputOffset} &= 100 - (1000 * 0.1) = 0\end{aligned}$$

$$\begin{aligned}\text{Now } 100\text{mv will yield a set point of } 10^\circ & ((100 * \text{gain}) + \text{offset}) \\ (100\text{mv} * 0.1) + 0 &= 10.0^\circ = 10\text{mv}/^\circ\text{C}\end{aligned}$$

Just like voltage the mv range can be set for other ranges.

A separate gain and offset will also be calculated and are applied to the DAC temperature output value.

$$\begin{aligned}\text{MilliVoltOutputGain} &= (\text{high mv} - \text{low mv}) / (\text{high temp} - \text{low temp}) \\ \text{MilliVoltOutputOffset} &= \text{high mv} - (\text{high temp} * \text{MilliVoltOutputGain})\end{aligned}$$

So that the default settings will yield the following gain and offset for temperature output.

$$\begin{aligned}\text{MilliVoltOutputGain} &= (1000 - 0\text{mv}) / (100 - 0) = 10 \\ \text{MilliVoltOutputOffset} &= 1000 - (100 * 10) = 0\end{aligned}$$

$$\begin{aligned}\text{Now } 10^\circ \text{ will output } 100\text{mv} & ((10 * \text{gain}) + \text{offset}) \\ (10^\circ * 10) + 0 &= 100\text{mv} = 10\text{mv}/^\circ\text{C}\end{aligned}$$

The following values reflect the standard scaling and are the factory defaults:

**Current (4ma – 20ma = the operating range of the chiller)**

<b>High ma</b>	<b>20.0ma</b>
<b>Low ma</b>	<b>4.0ma</b>
<b>High temp</b>	<b>40.0° (standard temp chillers) 90.0° (high temp chillers)</b>
<b>Low temp</b>	<b>5.0°</b>

With **Current** highlighted press enter to display:

USER MA	
◀ High MA	20.0
Low MA	4.0
High Temp	40.0
Low Temp	5.0

Highlight the desired value and then press enter, the highlight will flash. Use the up and down arrows to change the value. Once the desired value is displayed press enter to accept the change and stop the flashing.

Calculating the default gain and offset for set point input for a high temp chiller:

$$\begin{aligned} \text{MilliAmpInputGain} &= (\text{high temp} - \text{low temp}) / (\text{high ma} - \text{low ma}) \\ \text{MilliAmpInputOffset} &= \text{high temp} - (\text{high ma} * \text{MilliAmpInputGain}) \end{aligned}$$

The default settings for a high temp chiller will yield the following gain and offset for setpoint input:

$$\begin{aligned} \text{MilliAmpInputGain} &= (90.0 - 5.0) / (20.0 - 4.0) = 5.3125 \\ \text{MilliAmpInputOffset} &= 90 - (20 * 0.188235) = -16.25 \end{aligned}$$

Now 20ma will represent a setpoint of 90°

$$(20.0 * 5.3125) + (-16.25) = 90.0^\circ$$

A separate gain and offset will also be calculated and are be applied to the DAC temperature output value.

$$\begin{aligned} \text{MilliAmpOutputGain} &= (\text{high ma} - \text{low ma}) / (\text{high temp} - \text{low temp}) \\ \text{MilliAmpOutputOffset} &= \text{high ma} - (\text{high temp} * \text{MilliAmpOutputGain}) \end{aligned}$$

The default settings for a high temp chiller will yield the following gain and offset for temperature output:

$$\begin{aligned} \text{MilliAmpOutputGain} &= (20.0 - 4.0) / (90 - 5) = 0.188235294 \\ \text{MilliAmpOutputOffset} &= 20 - (90 * 0.188235) = 3.058823529 \end{aligned}$$

Now 90° = 20ma temperature output.



$$(90 * 0.188235294) + 3.058823529 = 20\text{ma}$$

**Note:** Ensure that there is sufficient range in the voltage or current to provide enough resolution of temperature. For example a voltage range of 0v – 1v to represent a temperature range of 0° - 150° would not have enough resolution to provide stable readings.

# Remote Sensor Connector Pinout

**Table 1**

Pin	Color
1	White
2	NA
3	NA
4	White
5	NA
6	NA
7	Red
8	NA
9	Red (4th wire not connected to the control board)

Never apply line voltage to any of the connections. ▲

When operating a ThermoFlex7500-10000 with the remote sensor enabled ensure the chiller's response to lowering the setpoint does not result in operation below 10°C process temperature. Operation below 10°C requires the use of 50/50 EG/water or 50/50 PG/water. ▲

**Figure D-2** Remote Sensor Connector

**[1] Manufacturer**

Thermo Fisher Scientific  
25 Nimble Hill Rd  
Newington, NH 03801  
USA

**[2] Authorized Representative**

Thermo Electron LED GmbH  
Robert-Bosch-Str. 1  
63505 Langenselbold  
Germany

**[3] Product Type**

ThermoFlex product line of chillers and heat exchangers for process flow applications

**[4] Model Name(s)**

Unit has a 15-digit or 16-digit part numbers defined on pages 2-3 of this document.

**[5] Rated Voltage, Rated Frequency**

Voltage ratings are reflected as part of the Model Name, defined on the following pages of this document.

[6] The object of the declaration described above is in conformity with the relevant Union harmonization legislation:

**[7] EC Directives**

2014/30/EU	Electromagnetic Compatibility Directive
2014/35/EU	Low Voltage Directive
2011/65/EU	RoHS Directive

[8] The object of the declaration described above is in conformity with Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment.

**[9] Standards**

EN 61326-1:2013  
EN 61010-1:2010, EN 61010-2-010:2014, EN 61010-2-011:2017  
EN 50581:2012

[10] This declaration is valid from the product manufactured after the date.

**28 September, 2018**

---

Mark Pearson  
Global Regulatory Affairs Director  
Laboratory Products Division  
Thermo Fisher Scientific

**Model and Voltage Number Definition :**

Unit has a 15-digit part number consisting of UU C VV PP c XXXXXXX as follows:

UU = unit type can be:

10= TF900	11= TF1400	12= TF2500	13= TF3500	14= TF5000
15= TF7500	16= TF10000	17= TF15000	18= TF20000	19= TF24000

C = Cooling type and Temperature Range where: A/C= Air-cooled, W/C-Water-cooled

1 - A/C Std Temp 5-40C,	2 - A/C Hi Temp 5-90C	3 - W/C Std Temp 5-40C	4 - W/C Hi Temp 5-90C
5 - A/C Temp -5 to 90C	6 - W/C Temp -5 to 90C	7 - A/C Temp -5 to 40C	8 - W/C Temp -5 to 40C
A - Hi Temp No Cooling	B - A/C Temp -5 to 55C	C - W/C Temp -5 to 55C	

VV = Voltage Ratings:

Single Phase:		Three Phase:	
10 = 115V/60Hz,	100V/50Hz	17 = (200/208/230V)/60Hz,	200V/50Hz
11 = (100/115V)/60Hz,	(100/115V)/50Hz	18 = 400V/50Hz	
12 = (208/230V)/60Hz,	200V/50Hz	20 = (200/208/230V)/60Hz,	(200-230V)/50Hz
16 = 230V/50Hz		21 = 460V/60Hz,	400V/50Hz
20 = (200/208/230V)/60Hz,	(200/230V)/50Hz		

PP = Pump type, can be 10 through 40 inclusive.

c = Controller type. Can be any digit from 1-8, inclusive.

X = Any digit from 0-9, used as sequential numbering only.

16-digit part number consisting of UU C V P c XXXXXXXX defined as follows:

UU = unit type can be:

10= TF900	11= TF1400	12= TF2500	13= TF3500	14= TF5000
15= TF7500	16= TF10000	17= TF15000	18= TF20000	19= TF24000

C = Cooling type and Temperature Range where: A/C= Air-cooled, W/C-Water-cooled

1 - A/C Std Temp 5-40C	2 - A/C Hi Temp 5-90C	3 - W/C Std Temp 5-40C	4 - W/C Hi Temp 5-90C
5 - A/C Temp -5 to 90C	6 - W/C Temp -5 to 90C	7 - A/C Temp -5 to 40C	8 - W/C Temp -5 to 40C
A - Hi Temp No Cooling	B - A/C Temp -5 to 55C	C - W/C Temp -5 to 55C	D - W/C Temp 5 to 80C

V = Voltage Ratings:

Single Phase:		Three Phase:	
1 = 115V/60Hz,	100V/50Hz	5 = (200/208/230V)/60Hz,	200V/50Hz
2 = (100/115V)/60Hz,	(100/115V)/50Hz	6 = 400V/50Hz	
3 = (208/230V)/60Hz,	200V/50Hz	8 = 460V/60Hz,	400V/50Hz
4 = 230V/50Hz			
7 = (200-230V)/60Hz,	(200-230V)/50Hz		

P = Pump type, can be 0-9 or A-Z.

c = Controller type. Can be any digit from 1-6, inclusive.

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## WARRANTY

Thermo Fisher Scientific warrants for 24 months (**excluding MD1/MD2 Magnetic Drive and P1/P2 Positive Displacement pumps which are warranted for 12 months**) from date of shipment the Thermo Scientific ThermoFlex chiller according to the following terms.

Any part of the chiller manufactured or supplied by Thermo Fisher Scientific and found in the reasonable judgment of Thermo Fisher to be defective in material or workmanship will be repaired at an authorized Thermo Fisher Repair Depot without charge for parts or labor. The chiller, including any defective part must be returned to an authorized Thermo Fisher Repair Depot within the warranty period. The expense of returning the chiller to the authorized Thermo Fisher Repair Depot for warranty service will be paid for by the buyer. Our responsibility in respect to warranty claims is limited to performing the required repairs or replacements, and no claim of breach of warranty shall be cause for cancellation or rescission of the contract of sales of any chiller. With respect to chillers that qualify for field service repairs, Thermo Fisher Scientific's responsibility is limited to the component parts necessary for the repair and the labor that is required on site to perform the repair. Any travel labor or mileage charges are the financial responsibility of the buyer.

The buyer shall be responsible for any evaluation or warranty service call (including labor charges) if no defects are found with the Thermo Scientific product.

This warranty does not cover any chiller that has been subject to misuse, neglect, or accident. This warranty does not apply to any damage to the chiller that is the result of improper installation or maintenance, or to any chiller that has been operated or maintained in any way contrary to the operating or maintenance instructions specified in this Instruction and Operation Manual. This warranty does not cover any chiller that has been altered or modified so as to change its intended use.

In addition, this warranty does not extend to repairs made by the use of parts, accessories, or fluids which are either incompatible with the chiller or adversely affect its operation, performance, or durability.

Thermo Fisher Scientific reserves the right to change or improve the design of any chiller without assuming any obligation to modify any chiller previously manufactured.

THE FOREGOING EXPRESS WARRANTY IS IN LIEU OF ALL OTHER WARRANTIES, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO WARRANTIES OR MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE.

OUR OBLIGATION UNDER THIS WARRANTY IS STRICTLY AND EXCLUSIVELY LIMITED TO THE REPAIR OR REPLACEMENT OF DEFECTIVE COMPONENT PARTS AND Thermo Fisher Scientific DOES NOT ASSUME OR AUTHORIZE ANYONE TO ASSUME FOR IT ANY OTHER OBLIGATION.

Thermo Fisher Scientific ASSUMES NO RESPONSIBILITY FOR INCIDENTAL, CONSEQUENTIAL, OR OTHER DAMAGES INCLUDING, BUT NOT LIMITED TO LOSS OR DAMAGE TO PROPERTY, LOSS OF PROFITS OR REVENUE, LOSS OF THE CHILLER, LOSS OF TIME, OR INCONVENIENCE.

This warranty applies to chillers sold by Thermo Fisher Scientific. (Refer to the warranty for chillers sold by the affiliated marketing company of Thermo Fisher Scientific for any additional terms.) This warranty and all matters arising pursuant to it shall be governed by the law of the State of New Hampshire, United States. All legal actions brought in relation hereto shall be filed in the appropriate state or federal courts in New Hampshire, unless waived by Thermo Fisher Scientific.

Thermo Fisher Scientific  
168 Third Avenue  
Waltham, Massachusetts 02451  
United States

[www.thermofisher.com](http://www.thermofisher.com)