

Instruction Manual

## SuperScript<sup>™</sup> III CellsDirect cDNA Synthesis System

Catalog Nos. 18080-200 and 18080-300

**Version B** 18 April 2005 25-0731

## **Table of Contents**

Kit Contents and Storage	v
Additional Products	vi
Introduction	1
Lysing Cells	3
First-Strand cDNA Synthesis	8
PCR	9
Troubleshooting	13
Purchaser Notification	14
Technical Service	16
References	

## Kit Contents and Storage

Shipping and Storage	Kit components are shipped on dry ice and	should be stor	ed at –20°C.
Kit Components	Catalog no. 18080-200 provides reagents for Catalog no. 18080-300 provides reagents for		
•	Component	<u>25 Rxns</u>	<u>100 Rxns</u>
	Resuspension Buffer	250 µl	1 ml
	RNaseOUT <sup>™</sup> Recombinant		
	Ribonuclease Inhibitor (40 units/µl)	50 µl	200 µl
	DNase I (1 U/µl)	125 µl	500 µl
	10X DNase I Buffer	40 µl	160 µl
	25 mM EDTA	30 µl	120 µl
	Oligo(dT)20 (50 µM)	50 µl	120 µl
	10 mM dNTP Mix	25 µl	100 µl
	SuperScript <sup>™</sup> III RT (200 units/µl)	25 µl	100 µl
	5X RT Buffer*	150 µl	600 µl
	0.1 M DTT	50 µl	100 µl
	E. Coli RNase H (2 U/ $\mu$ l)	30 µl	100 µl
	HeLa Total RNA (10 ng/µl)	10 µl	10 µl
	Forward Control Primer (10 µM)	10 µl	10 µl
	Reverse Control Primer (10 µM)	10 µl	10 µl
	* 5X RT Buffer composition: 250 mM Tris-H 375 mM KCl; 15 mM MgCl <sub>2</sub>	Cl (pH 8.3, roo	om temp.),
Note	DNA polymerase is not included in this kit. polymerases and optional amplification pro starting on page 9.		
Quality Control	This kit was verified in an end-point RT-PC dilution of HeLa cells and GAPDH primers, product. No template controls and negative controls were also prepared. Results were co electrophoresis.	yielding a 1.1 reverse transc	8 kb PCR criptase

## **Additional Products**

## Additional

The following related products are available from Invitrogen.

### Products

Product	Size	Cat. No.
SuperScript <sup>™</sup> III Platinum <sup>®</sup> CellsDirect Two-Step qRT-PCR Kit	25 RT/100 qPCR 100 RT/500 qPCR	11737-030 11737-038
SuperScript <sup>™</sup> III Platinum <sup>®</sup> CellsDirect Two-Step qRT-PCR Kit with SYBR <sup>®</sup> Green	25 RT/100 qPCR 100 RT/500 qPCR	11738-060 11738-068
CellsDirect Resuspension Buffer and Lysis Enhancer	10 ml/1 ml	11739-010
Platinum <sup>®</sup> Taq DNA Polymerase	100 reactions 250 reactions 500 reactions	10966-018 10966-026 10966-034
Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase High Fidelity	100 reactions 500 reactions	11304-011 11304-029
Platinum <sup>®</sup> PCR SuperMix	100 reactions	11306-016
Platinum <sup>®</sup> PCR SuperMix High Fidelity	100 reactions	12532-016
E-Gel® Pre-cast Agarose Gels 0.8% Starter Pak 1.2% Starter Pak 2% Starter Pak 4% Starter Pak	9 gels and base 9 gels and base 9 gels and base 9 gels and base	G5000-08 G5000-01 G5000-02 G5000-04
100-bp DNA Ladder	50 µg	15628-019
1-Kb Plus DNA Ladder	250 μg 1,000 μg	10787-018 10787-026

#### Additional Materials Required

The following materials are provided by the user:

- Lysis Enhancer (optional; see page 4) (Catalog no. 11739-010)
- Mammalian cell cultures in growth media ٠
- Coulter Counter or hemacytometer •
- Centrifuge (for pelleting cells) •
- Incubator, water bath, or thermal cycler preheated to 75°C
- Trypsin (for adherent cell cultures only)
- 1X cold phosphate-buffered saline (PBS), without Ca++ or Mg++
- 0.2-ml thin-walled PCR tubes or 96-well PCR plates
- Ice
- Pipettes
- DNA polymerase and associated PCR reagents (see pages 9-11)
- Thermal cycler for PCR

## Introduction


# System The SuperScript<sup>™</sup> III CellsDirect cDNA Synthesis System is an optimized kit for synthesizing first-strand cDNA directly from mammalian cell lysate without first isolating the RNA. Lysis and reverse transcription are performed in the same tube, and the resulting first-strand cDNA is ready to use in cloning and PCR. For real-time quantitative RT-PCR, see the note on the following page.

In traditional RT-PCR, RNA is first isolated from cells in a time-consuming procedure that can lead to a loss of material. Using the SuperScript<sup>™</sup> III CellsDirect cDNA Synthesis System, the cells are lysed and the cDNA is generated from the lysate in a single tube with minimal handling and no sample loss. DNase I is added to eliminate genomic DNA prior to first-strand synthesis.

This kit has been optimized for small cell samples, ranging from 10,000 cells down to a single cell (as measured by serial dilution). The use of SuperScript<sup>™</sup> III Reverse Transcriptase ensures high specificity and high yields of cDNA from small amounts of starting material—as little as 10 pg total RNA.

After synthesis, the first-strand cDNA can be amplified with specific primers by PCR without intermediate organic extractions or ethanol precipitations.

The diagram below outlines the procedure:



## Introduction, continued

Note	For real-time quantitative RT-PCR (qRT-PCR) from cell lysate, we recommend the SuperScript <sup>™</sup> III Platinum <sup>®</sup> CellsDirect Two-Step qRT-PCR Kit (Catalog nos. 11737-030 and 11737-038) or the SuperScript <sup>™</sup> III Platinum <sup>®</sup> CellsDirect Two-Step qRT-PCR Kit with SYBR <sup>®</sup> Green (Catalog nos. 11738-060 and 11738-068). These kits include reagents and protocols that have been specifically optimized for real-time qRT-PCR.
Advantages of the Kit	<ul> <li>This kit has the following advantages:</li> <li>Compatible with a wide range of mammalian cell types grown under different treatment conditions</li> <li>Single-tube format minimizes reagent loss, sample loss, and</li> </ul>
	<ul> <li>handling time</li> <li>Total lysate volume is used in first-strand cDNA synthesis reaction, providing greater yields with a limited number of cells and allowing for detection of rare transcripts</li> </ul>
	<ul> <li>SuperScript<sup>™</sup> III Reverse Transcriptase, with reduced RNase H activity and higher thermal stability, produces high yields of cDNA in the first-strand synthesis reaction, for greater sensitivity and enhanced detection of rare transcripts</li> </ul>
	<ul> <li>Generates high-quality cDNA for use in a variety of applications, including cloning and PCR</li> <li>Simple protocol takes less than 2 hours</li> </ul>
SuperScript <sup>™</sup> III RT	SuperScript <sup>™</sup> III Reverse Transcriptase is an engineered version of M- MLV RT with reduced RNase H activity and increased thermal stability. The enzyme can be used to synthesize first-strand cDNA at temperatures up to 55°C, providing increased specificity, higher yields of cDNA, and more full-length product than other reverse transcriptases. Because SuperScript <sup>™</sup> III RT is not inhibited significantly by ribosomal and transfer RNA, it can effectively synthesize first-strand cDNA directly from total RNA. The concentration of SuperScript <sup>™</sup> III RT in this system has been optimized to synthesize first-strand cDNA from total RNA in cell lysate.
Control RNA and Primers	The control RNA provided with this system consists of HeLa Total RNA (10 ng/ $\mu$ l). The Forward Control Primer and Reverse Control Primer provided with this kit are designed from the human GAPDH gene and produce a 1.18-kb PCR product.

## Methods

Lysing Cells	
Introduction	In this step, you lyse your cells in Resuspension Buffer or a Resuspension Buffer/Lysis Enhancer solution and perform a DNase I digestion to remove genomic DNA from the sample.
Cell Types and Density	This kit has been optimized for small cell samples, ranging from 1 to 10,000 cells. The performance of this kit was verified using several different mammalian cell lines, including HeLa, COS-7, 293, Jurkat, CV1, and K562. Cells may be grown under variety of conditions and treatments. Any type of culture vessel can be used.
<b>Q</b> Important	<ul> <li>We recommend using a maximum of 10,000 cells per reaction. Higher numbers of cells may inhibit reverse transcription and result in reduced yields and/or truncated cDNA product.</li> <li>Make sure that all solutions and equipment that come in contact with the cells are sterile. Always use proper sterile technique and work in a laminar flow hood when handling cells.</li> </ul>
Required Materials	<ul> <li>The following materials are provided by the user:</li> <li>Lysis Enhancer (optional) (Catalog no. 11739-010)</li> <li>Mammalian cell cultures in growth media</li> <li>Coulter Counter or hemacytometer</li> <li>Centrifuge (for pelleting cells)</li> <li>Incubator, water bath, or thermal cycler preheated to 75°C</li> <li>Trypsin (for adherent cell cultures only)</li> <li>1X cold phosphate-buffered saline (PBS), without Ca<sup>++</sup> or Mg<sup>++</sup></li> <li>0.2-ml thin-walled PCR tubes or 96-well PCR plates</li> <li>Ice</li> <li>Pipettes</li> <li>The following materials are provided in the kit:</li> <li>Resuspension Buffer</li> <li>RNaseOUT<sup>™</sup> (40 U/µl)</li> <li>DNase I, Amplification Grade (1 U/µl)</li> <li>10X DNase I Buffer</li> <li>EDTA, 25 mM</li> <li>Optional: Control HeLa Total RNA</li> </ul>

Lysis Enhancer	We recommend using Lysis Enhancer in the following procedures. Lysis Enhancer has been specially formulated for use with this kit to facilitate cell lysis. Catalog no. 11739-010 provides 1 ml of Lysis Enhancer and 10 ml of Resuspension Buffer (the additional Resuspension Buffer may be required for larger-volume tissue- culture wells). <b>For Cells in Tissue-Culture Wells:</b> Prepare a 10:1 solution of Resuspension Buffer/Lysis Enhancer immediately before use (e.g., 10 μl of Resuspension Buffer to 1 μl of Lysis Enhancer).
Note	All steps should be performed on ice, and reagents should be chilled and/or thawed immediately prior to use. The incubator should be <b>preheated</b> to 75°C.
Control Reaction	For the control reaction, use 1 $\mu l$ of the HeLa Total RNA provided in the kit instead of cell lysate.
Lysing Adherent Cells or Cells in Suspension	<ul> <li>Use the following procedure for lysing adherent cell cultures in vessels larger than 24-well plate wells. For cells in suspension, skip Steps 1–4 and proceed to Step 5 below.</li> <li>1. Add enough trypsin to cover the adherent cells in your tissue culture dish, plate, or flask (e.g., for a 10-cm dish, use ~1 ml; for a T75 flask, use ~3 ml).</li> </ul>
	<ol> <li>Incubate for 5 minutes at room temperate or in a 37°C incubator.</li> <li>Check for cell detachment under a microscope. If cells have not detached, gently tap the disk or flask to dislodge the cells, or let the cells incubate longer, checking them every minute under a microscope.</li> </ol>
	<ol> <li>When all the cells have detached, add serum-containing media to a final volume of 10 ml (for 6- and 12-well plates, add a 1X–2X volume of media). Note that the media must contain serum to inactivate the trypsin.</li> </ol>
	5. Pipet the cells gently up and down to mix, and then transfer the cell suspension to a centrifuge tube.
	6. Spin the cells at $200 \times g$ for 5 minutes to pellet.
	<ol> <li>Aspirate the media and wash the cell pellet with 5–10 ml of 1X cold PBS.</li> </ol>
	8. Spin the cells at $200 \times g$ for 5 minutes to pellet.
	<ol> <li>Aspirate the PBS and resuspend the pellet in 500 µl to 1 ml of 1X cold PBS. Mix the cell solution gently.</li> </ol>
	Procedure continued on next page

Lysing	Proc	cedure continued from p	revious page
Adherent Cells or Cells in Suspension,	10.	concentration. Deter	ot to verify that the cells are at the desired mine cell density electronically using a nanually using a hemacytometer chamber.
continued	11.		ty using cold PBS so that it falls within the ls/μl. Count the cells again to verify cell
	12.	of Resuspension Buf	ed PCR tube or plate well <b>on ice,</b> add 10 µl fer. Then add 1 µl of Lysis Enhancer previous page) <i>or</i> 1 µl of RNaseOUT™
	13.	Transfer 1–2 µl of cel	ls (<10,000 cells) to the PCR tube/well.
			trol reaction, add 1 μl of the Control HeLa R tube or plate well instead of cell lysate.
	14.	cycler preheated to 7	ate to an incubator, water bath, or thermal 5°C and incubate for 10 minutes. <b>Control:</b> on, incubate for 3 minutes.
	15.	After incubation, spi proceed to <b>DNase I</b>	n briefly to collect the condensation and <b>Digestion,</b> page 7.
Amount of Resuspension Buffer in	Res	0	volumes of Resuspension Buffer or is Enhancer solution are required for lysing ls:
Tissue-Culture	Ν	umber of Wells in	Volume of Resuspension Buffer (or
Wells	Т	issue-Culture Plate	Buffer/ Lysis Enhancer) per Well
		24	100 µl
		48	50 µl
	1	96	10 µl

Higher volumes may be required. Cells should be completely covered by solution.



Extra Resuspension Buffer may be required for lysing cells in largervolume tissue-culture wells. Extra Resuspension Buffer is provided in catalog no. 11739-010, which includes 10 ml of Resuspension Buffer and 1 ml of Lysis Enhancer.

Lysing Cells in Tissue-Culture Wells	cell For	<b>te:</b> Seed cells in tissue-culture wells so that 10 μl of resuspended s will yield the desired concentration.
	folle 1.	owing lysis procedure. Aspirate the media in each well and wash each well with 1X cold
		PBS. Aspirate the PBS.
	2.	Add Resuspension Buffer or Resuspension Buffer/Lysis Enhancer solution (recommended; see page 4) to each well. See the table on page 5 for amounts. The buffer should cover the cells in the well.
	3.	Incubate the plates <b>on ice</b> for up to 10 minutes. During that period, tap the plate periodically and check the cells under a microscope every 2–3 minutes to see whether they have detached or burst.
	4.	After 10 minutes, gently pipet the cells up and down to dislodge the remaining attached cells. If the cells are difficult to detach, incubate the plates at room temperature or 37°C for an additional 5 minutes, checking the cells under a microscope periodically.
	5.	Count the cells or estimate their density based on the seeding density (10 $\mu$ l should contain <10,000 cells).
	6.	Transfer 10 $\mu$ l of the cell suspension to a 0.2-ml thin-walled PCR tube or plate well.
		<b>Control:</b> For the control reaction, add 10 $\mu$ l of Resuspension Buffer to a PCR tube or plate well, and then add 1 $\mu$ l of Control HeLa Total RNA.
	7.	Add 1 $\mu$ l of RNaseOUT <sup>M</sup> (40 U/ $\mu$ l) to the PCR tube/well.
	8.	Transfer the tube/plate to an incubator or thermal cycler preheated to 75°C and incubate for 10 minutes. <b>Control:</b> For the control reaction, incubate for 3 minutes.
	9.	After incubation, spin briefly to collect the condensation, and proceed to <b>DNase I Digestion</b> , page 7.
		Continued on next page

DNase I Digestion	In this step, you treat the cell lysate with DNase I to degrade any contaminating DNA.		degrade any
	1. Place each tube/plate from Step 15, page 5, c ice, and add the following:		Step 9, page 6, on
		<u>Component</u>	Amount
		DNase I, Amplification Grade (1 U/µl) 10X DNase I Buffer	5 μl 1.6 μl
	2.	Mix by gently pipetting up and down, and spin the contents.	n briefly to collect
	3.	Incubate for 5 minutes at room temperature. <b>N</b> incubation time (up to 10 minutes) may be user samples (>1,000 cells). However, incubation tim minutes can greatly reduce cDNA yield.	d for larger
	4.	Spin briefly, and add 1.2 $\mu$ l of 25 mM EDTA to on ice. Mix by gently pipetting up and down, a collect the contents.	
	5.	Incubate at 70°C for 5 minutes.	
	6.	Spin briefly and proceed to <b>First-Strand cDNA</b> page 8.	Synthesis,

## **First-Strand cDNA Synthesis**

#### The following materials are provided by the user: Required Materials Thermal cycler preheated to 70°C • Ice Pipettes . The following materials are provided in the kit: $Oligo(dT)_{20}$ (50 $\mu$ M) 10 mM dNTP Mix • 5X RT Buffer • RNaseOUT<sup>™</sup> (40 U/µl) SuperScript<sup>™</sup> III RT (200 U/µl) • 0.1 M DTT • RNase H (2 U/ $\mu$ l) • Place each tube from DNase I Digestion, Step 6, page 7, on ice, First-Strand 1. and add the following: **cDNA** Synthesis Component Amount $Oligo(dT)_{20}$ (50 mM) 2 µl 10 mM dNTP Mix 1 ul 2. Mix by gently pipetting up and down, and spin the tube briefly to collect the contents. Incubate the tube at 70°C for 5 minutes. Spin the tube briefly to 3. collect the condensation. 4. Place the tube on ice for 2 minutes, and then add the following: Component Amount 5X RT Buffer 6 µl RNaseOUT<sup>TM</sup> (40 U/ $\mu$ l) 1 ul SuperScript<sup>™</sup> III RT (200 U/µl)\* 1 µl 0.1 M DTT 1 ul \*For negative RT controls, use 1 µl of sterile, distilled water instead of SuperScript<sup>™</sup> III RT 5. Mix by gently pipetting up and down, and spin the tube briefly to collect the contents. Transfer the tube to a thermal cycler preheated to 50°C. Incubate 6. for 50 minutes. 7. Inactivate the reaction at 85°C for 5 minutes. 8. Add 1 µl of RNase H (2 U/µl) to each tube and incubate at $37^{\circ}$ C for 20 minutes. Note: This step is optional if you are amplifying short targets (<1.0 kb) in end-point PCR. 9 Chill the reaction on ice.

Introduction	The first-strand cDNA generated using this kit PCR without additional purification. This section protocols for PCR and high-fidelity PCR.		
PCR Products	The following PCR products are available separ	ately:	
	<u>Product</u>	Size	Cat. No.
	Platinum <sup>®</sup> Taq DNA Polymerase	250 rxns 500 rxns	10966-018 10966-026 10966-034 10966-083
	Platinum <sup>®</sup> Taq DNA Polymerase High Fidelity	500 rxns	11304-011 11304-029 11304-102
	Platinum <sup>®</sup> PCR SuperMix	100 rxns	11306-016
	Platinum <sup>®</sup> PCR SuperMix High Fidelity	100 rxns	12532-016
PCR Enzymes	For amplifying the first-strand cDNA generated using this kit, we recommend Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase for targets < 1.0 kb and Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase High Fidelity for targets > 1.0 kb. <b>Platinum<sup>®</sup> <i>Taq</i> DNA Polymerase</b> is recombinant <i>Taq</i> DNA polymerase complexed with proprietary Platinum <sup>®</sup> antibodies that block polymerase activity at ambient temperatures. Activity is restored after the denaturation step in PCR cycling at 94°C, providing an automatic "hot start" for <i>Taq</i> DNA polymerase in PCR. Hot starts in PCR provide increased sensitivity, specificity, and yield, while allowing assembly of reactions at room temperature. The use of Platinum <sup>®</sup> antibodies helps reduce PCR optimization requirements, reaction set-up and handling time, and contamination risk, thereby improving PCR results for templates up to 5 kb.		
	Platinum® <i>Taq</i> DNA Polymerase High Fidelity recombinant <i>Taq</i> DNA polymerase, <i>Pyrococcus</i> s polymerase, and Platinum® <i>Taq</i> antibody. Platir complexes with <i>Taq</i> DNA polymerase and inhib temperatures, allowing room-temperature setup after the PCR denaturation step at 94°C, provid start" for the enzyme and increasing specificity, <i>Pyrococcus</i> species <i>GB-D</i> polymerase is a proofm possesses a 3' to 5' exonuclease activity. The enzy Platinum® <i>Taq</i> DNA Polymerase High Fidelity of	y is a mixtu species <i>GB</i> - num <sup>®</sup> antibu pits activity p. Activity ing an auto sensitivity eading enzy zyme mixtu	re of D ody at ambient is restored omatic "hot and yield. yme that ure in
	increase in fidelity over <i>Taq</i> DNA polymerase a amplification of simple and complex DNA temprange of target sizes, up to 12 kb with no optimi	plates over	

## PCR, continued



- Since PCR is a powerful technique capable of amplifying trace amounts of DNA, take all appropriate precautions to avoid sample contamination.
- Annealing and extension conditions are dependent on primer T<sub>m</sub>, and should be determined independently for each reaction.
- If PCR efficiency is not optimal, repeat the reaction with a primer titration from 100 to 500 nM (final conc.) in 100-nM increments.

#### PCR — Targets Up to 1 KB

The following protocol uses Platinum<sup>®</sup> *Taq* DNA Polymerase in a standard PCR reaction. Adjust the reaction size as needed. Optimal reaction conditions—including incubation times and temperatures, and concentrations of enzyme, primers, and MgCl<sub>2</sub>—may vary.

**Note:** A concentration of  $1.5 \text{ mM MgCl}_2$  is sufficient for most targets. For further optimization, prepare a titration from 1.5 mM to 3 mM in 0.25-mM increments.

1. Add the following components to a sterile 0.2- or 0.5-ml PCR tube or plate well at room temperature or on ice. For multiple reactions, prepare a master mix of common components.

<u>Components</u>	Volume	Final Conc.
10X PCR Buffer, Minus Mg	5 µl	1X
10 mM dNTP mixture	1 µl	0.2 mM each
50-mM MgCl <sub>2</sub>	1.5 µl	1.5 mM
Sense primer (10 µM)	1 µl	0.2 µM
Antisense primer (10 µM)	1 µl	0.2 µM
cDNA from Step 10, page 8	2 µl	—
Platinum <sup>®</sup> Taq DNA Polymerase	0.4 µl	2.0 units*
Autoclaved, distilled water	to 50 µl	n/a

\*2.0 units are recommended for amplifying cDNA from the CellsDirect kit. In some cases, more enzyme may be required (up to 2.5 units).

- 2. Mix contents of the tubes and overlay with 50  $\mu l$  of mineral or silicone oil, if necessary.
- 3. Cap the tubes and centrifuge briefly to collect the contents.
- 4. Incubate tubes in a thermal cycler at 94°C for 30 seconds to 2 minutes to denature the template and activate the enzyme.
- 5. Perform 30–40 cycles of PCR amplification as follows:

Denature	94°C for 15–30 seconds
Anneal	55–65°C for 30 seconds
Extend	72°C for 1 minute per kb

- 6. Maintain the reaction at 4°C after cycling. The samples can be stored at -20°C until use.
- Analyze the products by agarose gel electrophoresis and visualize by ethidium bromide staining. Use appropriate molecular weight standards.

PCR — Targets Above 1 KB	The following protocol uses Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase High Fidelity. Adjust the reaction size as needed. Optimal reaction conditions—including incubation times and temperatures, and the concentrations of Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase High Fidelity, primers, MgSO <sub>4</sub> , and template DNA—may vary.					
	<b>Note:</b> A concentration of 2 mM MgSO <sub>4</sub> is sufficient for most targets. For further optimization, prepare a titration from 2 mM to 4 mM in $0.25$ -mM increments.					
	1.	Add the following components to a sterile 0.2- or 0.5-ml PCR tube or plate well at room temperature or on ice. For multiple reactions, prepare a master mix of common components.				
		<u>Component</u> 10X High Fidelity PCR Buffer 10-mM dNTP mixture 50-mM MgSO <sub>4</sub> Sense primer (10 $\mu$ M) Antisense primer (10 $\mu$ M) cDNA from Step 10, page 8 Platinum <sup>®</sup> <i>Taq</i> High Fidelity Autoclaved, distilled water	<u>Volume</u> 5 μl 1 μl 2 μl 1 μl ≥1 μl 0.2 μl to 50 μl	Final Conc 1X 0.2 mM each 2 mM 0.2 µM 0.2 µM — 1.0 unit* Not applicable		
	2.	*1.0 unit is sufficient for amplifying more enzyme may be required (up Mix contents of the tubes and over	to 2.5 units	).		
	<ol> <li>3.</li> <li>4.</li> <li>5.</li> <li>6.</li> </ol>	silicone oil, if necessary. Cap the tubes and centrifuge brieff Incubate tubes in a thermal cycler 2 minutes to denature the template Perform 30–40 cycles of PCR ampl Denature 94°C for 15–30 seco Anneal 55–65°C for 30 seco Extend 68°C for 1 minute p Maintain the reaction at 4°C after of stored at -20°C until use.	at 94°C for 3 a and activa ification as a onds onds oer kb	30 seconds to te the enzyme. follows:		
	7.	Analyze the products by agarose g visualize by ethidium bromide sta molecular weight standards.				
- Real-Time qRT- PCR	reco qRT Sup SYE The	real-time quantitative RT-PCR (qRT ommend the SuperScript <sup>™</sup> III Platinu Г-PCR Kit (Catalog nos. 11737-030 ar eerScript <sup>™</sup> III Platinum <sup>®</sup> CellsDirect ' BR <sup>®</sup> Green (Catalog nos. 11738-060 ar se kits include reagents and protoco imized for real-time qRT-PCR.	um <sup>®</sup> CellsDi nd 11737-038 Two-Step ql nd 11738-06	rect Two-Step 8) or the RT-PCR Kit with 8).		

## PCR, continued

PCR — Control Reaction	The following protocol uses Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase and the Control Primers provided in the kit.			
	1.	Prepare a PCR mixture for each control reaction from Step 10, page 8. For each control reaction, add the following to a sterile 0.2- or 0.5-ml PCR tube or plate well at either room temperature or on ice:		
		Component	Volume	
		DEPC-treated water	38.1 µl	
		10X PCR buffer minus Mg <sup>++</sup>	5 μl	
		50 mM MgCl <sub>2</sub> 10 mM dNTP mix	1.5 μl	
		Forward Control Primer (10 µM)	1 μl 1 μl	
		Reverse Control Primer (10 $\mu$ M)	1 μl	
		cDNA from control RNA/	1 µ1	
		negative RT control, Step 10, page 8	2 µl	
		Platinum <sup>®</sup> <i>Taq</i> DNA polymerase (5 units/µl)	<u>0.4 µl</u>	
		final volume	50 µl	
	2.	Mix the contents of the tube. Centrifuge briefl reaction components.	y to collect the	
	3.	Place reaction mixture in preheated (94°C) the Perform an initial denaturation step: 94°C for		
	4.	Perform 40 cycles of PCR:		
		Denature 94°C for 15 sec		
		Anneal 60°C for 30 sec		
		Extend 72°C for 1 min		
		<b>Note:</b> For slow-ramping thermal cyclers, follo directions.	w manufacturer's	
	5.	Upon completion, maintain reactions at 4°C.		
	6.	Analyze $10 \ \mu$ l of each sample using agarose gel electrophoresis and ethidium bromide staining. A 1.18-kb band corresponding to at least 25 ng of product should be visible for the control sample. No band should be visible for the negative RT control sample.		

Problem	Possible Cause	Suggested Solution
No bands after electrophoretic analysis of amplified products	Procedural error	Confirm that all steps were followed. Use the Control RNA to verify the efficiency of the first- strand reaction (see page 12 for troubleshooting with the Control RNA).
	RNA is degraded	Add control total HeLa RNA to sample to determine if RNase is present in the first-strand reaction.
		Confirm that RNaseOUT <sup>™</sup> was added at the appropriate steps in the protocol.
		A longer DNase I digestion can hydrolyze the RNA in the sample. Use a digestion time of <10 minutes.
		Maintain aseptic conditions to prevent RNase contamination.
	Target mRNA contains strong transcriptional pauses	Use random hexamers (Cat. no. 48190-011) instead of $oligo(dT)_{20}$ in the first-strand reaction.
		Maintain an elevated temperature after the annealing step.
		Increase the temperature of first-strand reaction (up to $55^{\circ}$ C).
		Use PCR primers closer to the 3' terminus of the target cDNA.
	Too much first-strand product was used in PCR	Use no more than 5 $\mu$ l of the first-strand product in PCR.
Unexpected bands after electrophoretic analysis	Contamination by genomic DNA	Do not omit the DNase Digestion step on page 7. For larger samples (>1,000 cells), use a longer DNase I incubation time, i.e., up to 10 minutes.
		Design primers that anneal to sequence in exons on both sides of an intron or exon/exon boundary of the mRNA to allow differentiation between amplification of cDNA and products potential contaminating genomic DNA.
		To test if products were derived from DNA, prepare a negative RT control.
	Nonspecific annealing of primers	Vary the annealing conditions. Use Platinum <sup>®</sup> <i>Taq</i> DNA Polymerase for automatic hot-start PCR. Optimize magnesium concentration for each
		template and primer combination.

## Troubleshooting

## **Purchaser Notification**

Limited Use Label License No: 4: Products for PCR that include no rights to perform PCR This product is optimized for use in the Polymerase Chain Reaction (PCR) covered by patents owned by Roche Molecular Systems, Inc. and F. Hoffmann-La Roche, Ltd. ("Roche"). No license under these patents to use the PCR process is conveyed expressly or by implication to the purchaser by the purchase of this product. A license to use the PCR process for certain research and development activities accompanies the purchase of certain reagents from licensed suppliers such as Invitrogen, when used in conjunction with an Authorized Thermal Cycler, or is available from Applied Biosystems. Further information on purchasing licenses to practice the PCR process may be obtained by contacting the Director of Licensing at Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404 or at Roche Molecular Systems, Inc., 1145 Atlantic Avenue, Alameda, California 94501.

#### Limited Use Label License No. 138: SuperScript<sup>™</sup> III Reverse Transcriptase

The purchase of this product conveys to the buyer the nontransferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) not to transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of patents owned by Invitrogen Corporation and claiming this product based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

Limited Use Label License No: 18: RNaseOUT<sup>™</sup> Ribonuclease Inhibitor This product is the subject of U.S. Patent No. 5,965,399 owned by Invitrogen Corporation. The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and components of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Licensing Department, Invitrogen Corporation, 1600 Faraday Avenue, Carlsbad, California 92008. Phone (760) 603-7200. Fax (760) 602-6500.

#### Trademarks

TaqMan® is a registered trademark of Roche Molecular Systems, Inc.

## **Technical Service**

World Wide Web



Visit the Invitrogen Web Resource using your World Wide Web browser. At the site, you can:

- Get the scoop on our hot new products and special product offers
- View and download vector maps and sequences
- Download manuals in Adobe<sup>®</sup> Acrobat<sup>®</sup> (PDF) format
- Explore our catalog with full color graphics
- Obtain citations for Invitrogen products
- Request catalog and product literature

Once connected to the Internet, launch your Web browser (Internet Explorer 5.0 or newer or Netscape 4.0 or newer), then enter the following location (or URL):

#### http://www.invitrogen.com

...and the program will connect directly. Click on underlined text or outlined graphics to explore. Don't forget to put a bookmark at our site for easy reference!

### **Contact Us**

For more information or technical assistance, call, write, fax, or email. Additional international offices are listed on our Web page (www.invitrogen.com).

#### **Corporate Headquarters:**

Invitrogen Corporation 1600 Faraday Avenue Carlsbad, CA 92008 USA Tel: 1 760 603 7200 Tel (Toll Free): 1 800 955 6288 Fax: 1 760 602 6500 E-mail: tech\_service@invitrogen.com Japanese Headquarters: Invitrogen Japan K.K. Nihonbashi Hama-Cho Park Bldg. 4F 2-35-4, Hama-Cho, Nihonbashi Tel: 81 3 3663 7972 Fax: 81 3 3663 8242 E-mail: jpinfo@invitrogen.com

#### European Headquarters:

Invitrogen Ltd Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK Tel: +44 (0) 141 814 6100 Tech Fax: +44 (0) 141 814 6117 E-mail: eurotech@invitrogen.com

#### 

Limited Warranty	Invitrogen is committed to providing our customers with high- quality goods and services. Our goal is to ensure that every customer is 100% satisfied with our products and our service. If you should have any questions or concerns about an Invitrogen product or service, contact our Technical Service Representatives.		
	Invitrogen warrants that all of its products will perform according to specifications stated on the certificate of analysis. The company will replace, free of charge, any product that does not meet those specifications. This warranty limits Invitrogen Corporation's liability only to the cost of the product. No warranty is granted for products beyond their listed expiration date. No warranty is applicable unless all product components are stored in accordance with instructions. Invitrogen reserves the right to select the method(s) used to analyze a product unless Invitrogen agrees to a specified method in writing prior to acceptance of the order.		
	Invitrogen makes every effort to ensure the accuracy of its publications, but realizes that the occasional typographical or other error is inevitable. Therefore Invitrogen makes no warranty of any kind regarding the contents of any publications or documentation. If you discover an error in any of our publications, please report it to our Technical Service Representatives.		
	Invitrogen assumes no responsibility or liability for any special, incidental, indirect or consequential loss or damage whatsoever. The above limited warranty is sole and exclusive. No other warranty is made, whether expressed or implied, including any warranty of merchantability or fitness for a particular purpose.		

## References

- Berger, S.L. and Kimmel, A.R. (1987) Methods Enzymol 152, 316.
- Bracete, A.M., Mertz, L.M., Fox, D.K. (1999) Focus® 21, 38.
- Chomczynski, P. (1993) Biotechniques Vol. 15, 532.
- Chomczynski, P. and Sacchi, N. (1987) Anal. Biochem. 162, 156.
- Compton, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 39, Academic Press, Inc.
- D'Alessio, J. M., Gruber, C. E., Cain, C., and Noon, M. C. (1990) Focus® 12, 47.
- Frohman, M.A., Dush, M.K, and Martin, G.R. (1988) Proc. Nat. Acad. Sci USA 85, 8998.
- Gerard, G.F. (1994) Focus® 16, 102.
- Gerard, G.F., D'Alessio, J.M., and Kotewicz, M.L. (1989) Focus® 11, 66.
- Gerard, G.F., Schmidt, B.J., Kotewicz, M.L., and Campbell, J.H. (1992) Focus® 14, 91.
- Hu, A.W., D'Alessio, J.M., Gerard, G.F., and Kullman, J. (1991) Focus® 13, 26.
- Lee, C.C. and Caskey, T. (1990) in *PCR Protocols: A Guide to Methods and Applications* (Innis, M., Gelfand, D., Sninsky, J., and White, T., eds.), p. 46, Academic Press, Inc.
- Sambrook J., Fritsch, E.F., and Maniatis, T. (1989) Molecular Cloning: A Laboratory Manual, Cold Spring Harbor Laboratory Press, Cold Spring Harbor.
- Simms, D., Cizdziel, P.E., and Chomczynski, P. (1993) Focus® 15, 99.
- Westfall, B., Sitaraman, K., Solus, J., Hughes, J., and Rashtchian, A. (1997) Focus® 19, 46.
- Westfall, B., Sitaraman, K., Berninger, M., and Mertz, L.M. (1995) Focus® 17, 62.
- Westfall, B., Sitaraman, K., Lee, J., Borman, J. and Rashtchian, A. (1999) Focus® 21, 49.
- Takagi, M., Nishioka, M., Kakihara, H., Kitabayashi, M., Inoue, H., Kawakami, B., Oka, M., and Imanaka, T. (1997) *Appl. Environ. Microbiol.* 63, 4504.
- Sitaraman, K., Darfler, M., and Westfall, B. (1999) Focus<sup>®</sup> 21, 10.
- Nathan, M., Mertz, L., Fox, D. (1995) Focus® 17, 78.
- Schwabe, W., Lee, J.E., Nathan, M., Xu, R.H., Sitaraman, K., Smith, M., Potter, R.J., Rosenthal, K., Rashtchian, A., Gerard, G.F. (1998) *Focus*<sup>®</sup> 20, 30.

©2004–2005 Invitrogen Corporation. All rights reserved. For research use only. Not intended for any animal or human therapeutic or diagnostic use.



Corporate Headquarters: Invitrogen Corporation 1600 Faraday Avenue Carlsbad, California 92008 Tel: 1 760 603 7200 Tel (Toll Free): 1 800 955 6288 Fax: 1 760 602 6500 Email: tech\_service@invitrogen.com

#### **European Headquarters:**

Invitrogen Ltd Inchinnan Business Park Inchinnan Business Park 3 Fountain Drive Paisley PA4 9RF, UK Tel (Free Phone Orders): 0800 269 210 Tel (General Enquiries): +44 (0) 141 814 6100 Fax: +44 (0) 141 814 6260 Email: wateche @initiates and enm Email: eurotech@invitrogen.com

## International Offices: Argentina 5411 4556 0844

Australia 1 800 331 627 Brazil 55 11 5051 7422 Canada 800 263 6236 China 10 6849 2578 Hong Kong 2407 8450 Japan 03 3663 7974 New Zealand 0800 600 200 Singapore 65 686 186 38 Taiwan 2 2651 6156

For other countries see our website

www.invitrogen.com