INSTRUCTIONS



1897.13

DyLight Antibody Labeling Kits

Number Description 62275 **DyLight 350 Antibody Labeling Kit** 53020 **DyLight 405 Antibody Labeling Kit** 53024 DyLight 488 Antibody Labeling Kit 84530 DyLight 550 Antibody Labeling Kit 53044 DyLight 594 Antibody Labeling Kit 53046 DyLight 633 Antibody Labeling Kit 84535 DyLight 650 Antibody Labeling Kit 53056 DyLight 680 Antibody Labeling Kit 84538 **DyLight 755 Antibody Labeling Kit** 53062 **DyLight 800 Antibody Labeling Kit Note:** Each kit contains sufficient reagents to label and purify 3×1 mg (2mg/mL) of IgG or similar amounts of other proteins. **Kit Contents: DyLight NHS Ester**, $3 \times 50 \mu g$ vials (except 350 and 594 Kits, which contain $3 \times 65 \mu g$ vials) Borate Buffer (0.67M), 1mL Purification Resin, 5mL Spin Columns, 6 each Microcentrifuge Collection Tubes, 12 each

Storage: Upon receipt, store the DyLight NHS Ester at -20°C. Store all other kit components at 4°C.

Introduction

Each Thermo ScientificTM DyLightTM Antibody Labeling Kit contains all the necessary components for three labeling reactions of any protein with a molecular weight of 50-150K and subsequent excess dye removal. The DyLight Dyes included in these kits are activated with *N*-hydroxysuccinimide (NHS) esters, which is the most commonly used reactive group for labeling proteins. NHS esters react with primary amines, forming stable, covalent amide bonds and releasing the NHS groups. The DyLight Dyes have absorption spectra ranging from 350nm to 770nm (Table 1). These reagents fluoresce in a broad pH range, are more intense than Alexa FluorTM or CyTM Dyes in many applications, and match the output wavelengths of common fluorescence instrumentation. Additionally, the water solubility of the DyLight Reagents allows conjugations using a high dye-to-protein ratio without precipitation.

The purification resin and spin columns eliminate equilibration steps and the need to collect and monitor gravity-flow fractions. This system enables efficient removal of excess dye and, therefore, accurate determination of the dye-to-protein ratio and exceptional protein recovery.



Table 1. Properties of the Thermo Scientific DyLight NHS-Ester Dyes.						
DyLight Dye	Ex/Em*	ε†	MW (g/mol)	Spectrally Similar Dyes		
350	353 / 432	15,000	874	Alexa Fluor 350, AMCA		
405	400 / 420	30,000	793	Alexa Fluor 405		
488	493 / 518	70,000	1011	Alexa Fluor 488, Cy2		
550	562 / 576	150,000	1040	Alexa Fluor 555, Cy3		
594	593 / 618	80,000	1078	Alexa Fluor 594, Texas Red		
633	638 / 658	170,000	1066	Alexa Fluor 633		
650	652 / 672	250,000	1066	Alexa Fluor 647, Cy5		
680	682 / 715	140,000	950	Alexa Fluor 680		
755	754 / 776	220,000	1092	Alexa Fluor 750		
800	770 / 794	270,000	1050	IRDye 800		

* Excitation and emission maxima in nanometers

†Molar extinction coefficient (M⁻¹ cm⁻¹)

Important Product Information

- NHS ester-activated dyes are moisture-sensitive. Prepare the DyLight Labeling Reagent immediately before use and discard any unused reconstituted reagent. Do not store reconstituted labeling reagent.
- Low concentrations of sodium azide (≤ 3 mM or 0.02%) or thimerosal (≤ 0.02 mM or 0.01%) will not significantly interfere with protein labeling; however, 20-50% glycerol will reduce labeling efficiency.
- Use the following fluorescent imagers:
 - 350 dye: UV argon-ion laser at 351-363nm
 - 405 dye: Spectral line of the blue diode laser
 - 488 dye: Green (526) laser
 - 550 and 594 dyes: Green (532) laser
 - 633 and 650 dyes: Red (633) laser
 - 680, 755 and 800 dyes: laser- and filter-based instruments that emit in the 700nm and 800nm region of the spectrum, respectively; these dyes are well-suited for the 700 and 800 channels of the LI-COR OdysseyTM and the LI-COR AeriusTM Infrared Imaging Systems.

Additional Materials Required

- Variable-speed centrifuge
- Phosphate-buffered saline (PBS; for measuring the dye-to-protein ratio)

Procedure for Labeling Proteins with DyLight Dyes

A. Protein Preparation

Note: When labeling with the DyLight 594 Dye, prepare the protein in phosphate-buffered saline.

Note: If Borate Buffer precipitates during storage, solubilize it by warming at 37-50°C and vigorously vortexing the vial.

- 1. The optimal labeling buffer is 50mM sodium borate, pH 8.5 (please see note above). For best results use 1mg of protein at ~2mg/mL. Prepare the protein as follows:
 - **Proteins Lyophilized in PBS:** Just before use, prepare the labeling buffer by diluting the Borate Buffer (0.67M) to 0.05M in PBS or ultrapure water. Prepare only enough labeling buffer required for the reaction (for example, to prepare 1mL, add 75µL of Borate Buffer (0.67M) to 925µL of ultrapure water or PBS). Reconstitute 1mg of protein with 0.5mL of labeling buffer.
 - Proteins in PBS Solution: Add 40µL of the Borate Buffer (0.67M) to 0.5mL of 2mg/mL protein in PBS. If the protein is > 2mg/mL, adjust the concentration to 2mg/mL with labeling buffer (for example, 0.05M sodium borate see the above bullet point: Proteins Lyophilized in PBS).
 - **Proteins in Other Buffers:** Protein must not be in a buffer containing ammonium ions or primary amines (for example, Tris or glycine). If necessary, replace buffer with 50mM sodium borate (Product No. 28384), pH 8.5 by dialysis or buffer exchange.



B. Protein Labeling

- Tap the bottom of the DyLight Reagent vial against a hard surface to ensure the dye is in the bottom of the tube. Add 1. 0.5mL of the prepared protein to the vial of DyLight Reagent, vortex gently and pipette up and down to mix.
- Briefly centrifuge the vial to collect the sample in the bottom of the tube. Incubate the reaction mixture for 60 minutes at 2. room temperature protected from light.

C. Protein Purification

- Place two spin columns in two microcentrifuge collection tubes. 1.
- Mix the Purification Resin to ensure uniform suspension and add 250µL of the suspension into both spin columns. 2. Centrifuge for 1 minute at $\sim 1000 \times g$ to remove the storage solution. Discard the used collection tubes and place the columns in new collection tubes.
- 3. Add 250-270µL of the labeling reaction to each spin column and mix the sample with the resin by briefly vortexing.
- Centrifuge columns for 1 minute at $\sim 1000 \times g$ to collect the purified proteins. Combine the samples from both columns 4. (~0.5mL total). Discard the used columns.
- Store the labeled protein protected from light at 4°C for up to one month. Alternatively, store labeled protein in single-5. use aliquots at -20°C. Avoid repeated freeze/thaw cycles. If the final concentration of conjugate is < 1mg/mL, add a stabilizing agent, such as bovine serum albumin at 1-10mg/mL.

D. Dye-to-Protein Ratio Estimation

- Dilute a small amount of labeled purified protein in PBS. 1.
- 2. Use a 1cm path length cuvette to measure absorbance at 280nm and the A_{max} of the specific dye (Table 2).

Table 2. I Toper ties of the Thermo Scientific Dylight Dyes.					
DyLight Dye A _{max} * ε† CF‡					
350 353 15,000 0.144					
405 405 30,000 0.564					
488 493 70,000 0.147					
550 557 150,000 0.081					
594 595 80,000 0.585					
633 627 170,000 0.110					
650 655 250,000 0.037					
680 684 140,000 0.128					
755 755 220,000 0.030					
800 777 270,000 0.045					

Table 2. Properties of the Thermo Scientific DyLight Dyes

* Excitation wavelength in nanometers - note that upon protein conjugation the absorption maximum shifts to the right of the spectra

†Molar extinction coefficient (M⁻¹ cm⁻¹) at A_{max}

‡Correction factor (A₂₈₀/A_{max})

3 Calculate protein concentration as follows:

Protein concentration (M) =
$$\frac{[A_{280} - (A_{max} \times CF)]}{\epsilon_{notein}} \times dilution factor$$

- $\varepsilon_{\text{protein}}$ = protein molar extinction coefficient (e.g., the molar extinction coefficient of IgG is ~210,000 M⁻¹ cm⁻¹)
- $CF = Correction \ factor = \frac{A_{280} \ of \ the \ dye}{A_{max} \ of \ the \ dye} \ (see \ Table \ 2)$
- Calculate the degree of labeling as follows: 4.

Moles dye per mole protein = $\frac{A_{max}}{\epsilon_{dye} \times \text{protein concentration (M)}}$

 ε_{dve} = dye (fluorophore) molar extinction coefficient (see Table 2)



Example calculations for DyLight 550 Dye conjugated to antibodies:

- Dilution factor = 10
- $A_{280} = 0.287$
- A_{max} at 557nm = 0.878

Protein concentration (M) = $\frac{[0.287 - (0.878 \times 0.081)]}{210,000} \times 10 = 0.00001028 \text{ M}$

Moles dye per mole protein = $\frac{0.878 \times 10}{150,000 \times 0.00001028} = 5.7$

Troubleshooting

Problem	Possible Cause	Solution	
Protein was not labeled	Protein buffer contained amines that interfered with labeling	Perform buffer exchange via dialysis or other method into 50mM sodium borate	
	The NHS ester is hydrolyzed and non-reactive	Prepare labeling reagent immediately before use – do not store reagent in aqueous solution	
The downstream application was unsuccessful	Protein was not labeled	Determine if the protein was labeled by calculating the dye-to-protein ratio	
Sample or buffer does not flow through resin	Centrifugation problem	Ensure that centrifuge is in proper working condition	
Low yield	Improper centrifugation	Make sure to use the indicated centrifugation speed	
	Unstable protein	Equilibrate the column with PBS or other suitable buffer before adding the labeled protein	
Problems with image capture	Photostability	DyLight 350 Dye is a relatively dim dye and quickly photobleaches under microscope light	

Visit our Website for Additional Information

- Tech Tip #43: Protein stability and storage
- Tech Tip #31: Calculate dye:protein (F/P) molar ratios

Related Thermo Scientific Products

28384

BupH[™] Borate Buffer Packs, 40 packs, each pack yields 500mL

28372 BupH Phosphate Buffered Saline Packs, 40 packs, each pack yields 500mL

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at <u>www.thermofisher.com/techresources</u>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2015 Thermo Fisher Scientific Inc. All rights reserved. Alexa Fluor is a trademark of Life Technologies, Inc. Cy is a trademark of GE Healthcare Limited

© 2015 Thermo Fisher Scientific Inc. All rights reserved. Alexa Fluor is a trademark of Life Technologies, Inc. Cy is a trademark of GE Healthcare Limited UK. Odyssey and Aerius are trademarks of LI-COR Biosciences. Unless otherwise indicated, all other trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.

Pierce Biotechnology	PO Box 117	(815) 968-0747	www.thermofisher.com/techresources
3747 N. Meridian Road	Rockford, IL 61105 USA	(815) 968-7316 fax	