

## Revised: 24–June–2004

# Alexa Fluor <sup>®</sup> Hydrazides

# Quick Facts

Storage upon receipt:

A10440, A10441, A10442:

- ≤-20°C
- Protect from light

A10436, A10437, A10438, A10439, A20501MP, A20502, A30634

- Room temperature
- Desiccate
- Protect from light

| Hydrazide  | Solid    | 10 mM Solution<br>for<br>Microinjection | Abs* | Em* |  |  |
|--|----------|---|------|-----|--|--|
| Alexa Fluor 350  | A10439   |   | 345  | 445 |  |  |
| Alexa Fluor 488  | A10436   | A10440                                  | 493  | 517 |  |  |
| Alexa Fluor 555  | A20501MP |   | 554  | 567 |  |  |
| Alexa Fluor 568  | A10437   | A10441                                  | 576  | 599 |  |  |
| Alexa Fluor 594  | A10438   | A10442                                  | 588  | 613 |  |  |
| Alexa Fluor 633  | A30634   |   | 624  | 643 |  |  |
| Alexa Fluor 647  | A20502   |   | 649  | 666 |  |  |
| * Approximate absorption (Abs) and fluorescence emission (Em) maxima, in nm. |          |   |      |     |  |  |

#### Table 1. Molecular Probes' Alexa Fluor hydrazides.

## Introduction

The Alexa Fluor<sup>®</sup> hydrazides offer a significant improvement over many conventional fluorophores for intracellular labeling experiments. The Alexa Fluor dyes are much brighter and more resistant to photobleaching, especially compared to fluorescein or lucifer yellow. The excitation spectra of the Alexa Fluor dyes closely matches the output of commonly used light sources



Figure 1. The APR motor neuron of a larval moth, Manduca sexta, was labeled by the intracellular injection of Alexa Fluor 488 hydrazide (A10436). This pseudocolored image was created by combining 21 optical sections obtained with a scanning confocal microscope equipped with a bandpass filter appropriate for fluorescein. Photo contributed by Jack Gray, Institute of Neuroscience, University of Oregon, and Walter K. Metcalfe, Molecular Probes, Inc.

(mercury-arc lamps and laser), they are readily visualized using common filter sets, they can be used with a variety of electrolytes and they can be iontophoretically injected using hyperpolarizing current. These properties make them the dye of choice for most intracellular applications.

## Materials

## Contents

Alexa Fluor hydrazides (Table 1) are now offered in solution, ready to inject into cells. Alexa Fluor 488 hydrazide (A10440), Alexa Fluor 568 hydrazide (A10441) and Alexa Fluor 594 hydrazide (A10442) are prepared as 10 mM solutions in 200 mM KCl and filter sterilized. Electrodes using these prepared solutions have conductances that are very nearly like those in 200 mM KCl alone, and are suitable for intracellular physiology, as well as for cell labeling.

Alexa Fluor hydrazides are also offered in solid form (Table 1) in a 1 mg unit size. Solutions may be prepared as needed if an experiment requires the use of electrolytes other than KCl (e.g., lithium chloride or potassium acetate). However, Alexa Fluor dyes may not be completely soluble in solutions of high ionic strength. For example, Alexa Fluor hydrazide 488 precipitates in 1 M potassium acetate at room temperature, but dissolves well in 200 mM potassium acetate.

## Storage and Handling

The Alexa Fluor hydrazides in solutions should be stored at  $\leq -20^{\circ}$ C, protected from light. The Alexa Fluor hydrazides as solids can be stored at room temperature, desiccated and protected from light. When stored properly, Alexa Fluor hydrazides, in either form, should remain stable for at least 6 months.

## Applications

#### **Cell Labeling**

Intracellular labeling with 10 mM Alexa Fluor hydrazide solutions in 200 mM KCl is accomplished by passing hyperpolarizing current through the electrode. For example, to label the huge accessory planta retractor (APR) motoneurons in larvae of *Manduca sexta*, 2-3 nA DC hyperpolarizing current for 15 minutes worked well for an electrode with 50 megohm resistance (Figure 1). Smaller cells will require much less time to fill. The amount of time required for dye diffusion to label fine processes is similar to the time required for other commonly used dyes and depends on the distance and diameter of the processes.

#### Fixation

These dyes may be fixed using standard aldehyde fixation (e.g. 4% paraformaldehyde in 100 mM phosphate buffer, pH 7). However, we have observed that the fluorescence of Alexa Fluor 488 hydrazide is largely quenched after overnight fixation. Very bright labeling with excellent histological preservation was obtained with only 30 minutes of fixation at room temperature in the APR motoneurons mentioned above. Following fixation, the preparation may be dehydrated in ethanol and cleared in xylene. For the best long-term observation and storage, the fixed specimen may be mounted in an antifade medium such as ProLong<sup>®</sup> reagent (P7481).

#### **Product List** Current prices may be obtained from our Web site or from our Customer Service Department.

| Cat #    | Product Name   | Unit Size |
|----------|--|-----------|
| A10439   | Alexa Fluor® 350 hydrazide, sodium salt  | 5 mg      |
| A10436   | Alexa Fluor® 488 hydrazide, sodium salt  | 1 mg      |
| A10440   | Alexa Fluor® 488 hydrazide, sodium salt *for microinjection* *10 mM in 200 mM KCI* | 125 µL    |
| A20501MP | Alexa Fluor® 555 hydrazide, tris(triethylammonium salt)                            | 1 mg      |
| A10437   | Alexa Fluor® 568 hydrazide, sodium salt  | 1 mg      |
| A10441   | Alexa Fluor® 568 hydrazide, sodium salt *for microinjection* *10 mM in 200 mM KCI* | 125 µL    |
| A10438   | Alexa Fluor® 594 hydrazide, sodium salt  | 1 mg      |
| A10442   | Alexa Fluor® 594 hydrazide, sodium salt *for microinjection* *10 mM in 200 mM KCI* | 125 µL    |
| A30634   | Alexa Fluor® 633 hydrazide, bis(triethylammonium salt)                             | 1 mg      |
| A20502   | Alexa Fluor® 647 hydrazide, tris(triethylammonium salt)                            | 1 mg      |

## **Contact Information**

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly from Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Leiden, the Netherlands. All others should contact our Technical Assistance Department in Eugene, Oregon.

Please visit our Web site - www.probes.com - for the most up-to-date information

Molecular Probes, Inc. 29851 Willow Creek Road, Eugene, OR 97402 Phone: (541) 465-8300 ● Fax: (541) 335-0504

Customer Service: 6:00 am to 4:30 pm (Pacific Time) Phone: (541) 335-0338 ● Fax: (541) 335-0305 ● order@probes.com

**Toll-Free Ordering for USA and Canada:** Order Phone: (800) 438-2209 • Order Fax: (800) 438-0228

Technical Assistance: 8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 ● Toll-Free: (800) 438-2209 Fax: (541) 335-0238 ● tech@probes.com Molecular Probes Europe BV Poortgebouw, Rijnsburgerweg 10

2333 AA Leiden, The Netherlands Phone: +31-71-5233378 ● Fax: +31-71-5233419

**Customer Service**: 9:00 to 16:30 (Central European Time) Phone: +31-71-5236850 • Fax: +31-71-5233419 eurorder@probes.nl

Technical Assistance: 9:00 to 16:30 (Central European Time) Phone: +31-71-5233431 ● Fax: +31-71-5241883 eurotech@probes.nl

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. Our products are not available for resale or other commercial uses without a specific agreement from Molecular Probes, Inc. We welcome inquiries about licensing the use of our dyes, trademarks or technologies. Please submit inquiries by e-mail to busdev@probes.com. All names containing the designation <sup>®</sup> are registered with the U.S. Patent and Trademark Office.

Copyright 2004, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.