

Pierce[®] Phosphoprotein Enrichment Kit

90003

2021.1

Number	Description
90003	Pierce Phosphoprotein Enrichment Kit Kit Contents: Phosphoprotein Enrichment Column , 10 each, 1mL resin bed Lysis/Binding/Wash Buffer , 325mL Elution Buffer , 60mL CHAPS , 1g White Column Caps , 10 each Pierce Concentrator, 7mL/9K , 10 devices

Storage: Upon receipt store at 4°C. Product shipped at ambient temperature.

Introduction

The Thermo Scientific Pierce Phosphoprotein Enrichment Kit enables efficient enrichment of phosphorylated proteins derived from mammalian cells in culture and tissue. The column contains a proprietary metal/buffer composition that produces superior yields with negligible nonspecific binding. For convenience the Pierce Phosphoprotein Enrichment Kit incorporates a spin format that produces enriched phosphoproteins within 2 hours. Once enriched, the phosphoproteins can be used for downstream proteomic analyses by Western blotting, mass spectrometry, 2D-PAGE and protein arrays.

Important Product Information

- Adding CHAPS (0.25%) to the Lysis/Binding/Wash Buffer is required and helps maintain protein solubility. However, adding CHAPS to the Elution Buffer is optional and NOT recommended if performing mass spectrometry downstream.
- Typical yields are 10-25% of the total cellular protein loaded onto the column. Factors that influence yield include cell type and cell treatment (pharmacological agents, growth factors, etc.).
- Specific binding of phosphoproteins is dependent upon protein isoelectric point (pI) and buffer pH. If required, adjust the pH of the Lysis/Binding/Wash Buffer or Elution Buffer to enrich for different populations of phosphoproteins.
- The Elution Buffer consists of 75mM sodium phosphate, 500mM sodium chloride; pH 7.5. Therefore, samples must be desalted before performing 2D gel analysis.
- The Lysis/Binding/Wash and Elution Buffers contain 0.01% sodium azide as an antimicrobial.
- Tissues must be homogenized in appropriate buffer, centrifuged and diluted 1:1 in Lysis/Binding/Wash Buffer before loading the sample onto a Phosphoprotein Enrichment Column in Section B.

Additional Materials

- EDTA-free protease inhibitors such as Thermo Scientific Halt Protease Inhibitor Single-Use Cocktail EDTA-Free (Product No. 78425)
- Phosphatase inhibitors such as Halt™ Phosphatase Inhibitor Cocktail (100X), (Product No. 78420)
- Phosphate-free buffer such as 50mM HEPES, pH 7.0
- 50mL conical tubes

- Centrifuge with swinging-bucket or fixed-angle rotor for 15mL and 50mL conical tubes, rated for $\geq 2000 \times g$
- Thermo Scientific T-PER Tissue Protein Extraction Reagent (Product No. 78510) (For tissue lysis protocol)
- Optional: DNase I (5000-10,000 units/mL) (Product No. 89835)
- Optional: PMSF (Product No. 36978)

Material Preparation

10% CHAPS Solution	Add 10mL of ultrapure water to the bottle containing 1g of CHAPS. Mix by inversion and vortex. Store solution at 4°C.
Lysis/Binding/Wash Buffer with CHAPS (0.25%)	Prepare enough Lysis/Binding/Wash Buffer with a final concentration of 0.25% CHAPS for each step, typically 24mL for a 2mg load.
Elution Buffer with CHAPS (0.25%)	If performing downstream mass spectrometry applications DO NOT add CHAPS to the Elution Buffer. Otherwise, to maintain protein solubility, prepare 5mL of Elution Buffer by adding 125 μ L of 10% CHAPS for each column being used.

Procedure for Enrichment of Phosphoproteins

Note: For cell lysis, use Section A and continue to Section C. For tissue lysis, proceed to Section B and continue with Section C, protocol step 2.

A. Cell Lysis

1. Wash cells 2X with a non-phosphate-based buffer (e.g., 50mM HEPES, pH 7.0). Remove HEPES wash buffer.
2. Lyse cells by adding 1mL (15cm plate) of Lysis/Binding/Wash Buffer with CHAPS, 1X Halt Protease Inhibitor EDTA-free and 1X Halt Phosphatase Inhibitor Cocktail to each cell culture plate.
3. Scrape cells into appropriate-sized collection tube and place cells on ice for 45 minutes, vortexing periodically.
4. Centrifuge lysed cells at $10,000 \times g$ for 20 minutes at 4°C to pellet cellular debris.
5. Collect supernatant and proceed to the next section, or freeze sample at -80°C for later use.

B. Tissue Lysis

1. Wash tissue two times with a sufficient volume of non-phosphate-based buffer (e.g., 50mM HEPES, pH 7.0) to completely cover the tissue. Remove wash buffer.
2. Weigh out tissue; for every 1gm of tissue, add 20mL of T-PER[®] Reagent or 10-20mL of RIPA Buffer containing 1X Halt Protease Inhibitor EDTA-free and 1X Halt Phosphatase Inhibitor Cocktail.
3. Using a polytron or tissue homogenizer, process tissue until homogeneous throughout.
4. Centrifuge the sample at $10,000 \times g$ for 5 minutes to pellet tissue debris.
5. Collect supernatant and proceed to the next step, or freeze sample at -80°C for later use.
6. Prior to Phosphoprotein Enrichment, determine the protein concentration of the supernatant by performing a protein assay.
7. Adjust the concentration of sample to 0.5mg/mL with Lysis/Binding/Wash Buffer. Up to 4mg of total protein may be applied to each column. Add a volume of Lysis/Binding/Wash Buffer that is at least equal to the volume of the sample; larger volumes can also be added if greater dilution is required.

C. Phosphoprotein Enrichment

1. Determine the protein concentration of the supernatant by performing a protein assay. Adjust the concentration to 0.5mg/mL with Lysis/Binding/Wash Buffer. Up to 4mg of total protein may be applied to each column.
2. Invert column to suspend the resin slurry.
3. Snap off bottom of column and loosen cap. Place column into a 50 ml conical tube.

4. To remove the storage solution, centrifuge column at $1000 \times g$ for 1 minute at 4°C and discard the flow-through.
5. Equilibrate resin by adding 5mL of Lysis/Binding/Wash Buffer with CHAPS. Place column into a 50mL conical tube, centrifuge at $1000 \times g$ for 1 minute at 4°C and discard the flow-through.
6. Plug bottom of column and add diluted lysate. Screw cap on top of column and invert several times to mix. Place on a platform rocker for 30 minutes at 4°C .
Note: The presence of protease and phosphatase inhibitors might make the resin slightly yellow; however, this color change has not adversely affected results.
7. Remove bottom plug and place column in a 50mL conical tube. Centrifuge column at $1000 \times g$ for 1 minute at 4°C . If desired, save the flow-through for downstream analysis.
8. Wash resin by adding 5mL Lysis/Binding/Wash Buffer with CHAPS. Place column into a 50mL conical tube and centrifuge at $1000 \times g$ for 1 minute at 4°C . If desired, save the wash fraction for downstream analysis.
9. Repeat Step 7 two times for a total of three washes.
10. Plug bottom of column and transfer to a new 50mL conical tube. Add 1mL of Elution Buffer (with optional CHAPS) and incubate at room temperature with occasional agitation for 3 minutes.
11. Remove plug and centrifuge at $1000 \times g$ for 1 minute at 4°C .
12. Repeat elution steps four times for a total of five elution fractions (5mL total). Pool the fractions and proceed to the next section or freeze sample at -80°C for later use.

D. Phosphoprotein Concentration

Note: The Pierce Concentrators can be used effectively at a relative centrifugal force of $2000\text{-}9000 \times g$; the recommended speeds are $2500\text{-}4500 \times g$ for swinging-bucket rotors and $5000\text{-}7000 \times g$ for fixed-angle rotors. Set centrifuge to the appropriate speed at 4°C .

1. Place sample into upper sample chamber of concentrator. Maximum sample volume is 7mL using a swinging bucket rotor and approximately 4.5-5.0mL using a fixed-angle rotor. The pooled elution fractions will be approximately 5mL.
2. To prevent over-concentrating, control the dead-stop volume by removing cap and upper chamber and placing 2.5mL of ultrapure water in the lower chamber. Replace the upper chamber and add the eluate (from Step B.11).
3. Cap and place concentrator assembly into rotor with proper counterbalance.
4. Centrifuge for approximately 30 minutes or until the sample is 150-200 μL . Remove concentrated protein promptly.
5. Use a 200 μL or similar pipette tip to gently aspirate concentrated sample from upper chamber.

Troubleshooting

Problem	Possible Cause	Solution
Low yield	Cell harvest conditions were not optimal	Before harvest, wash cells 2X with a non-phosphate-based buffer such as 50mM HEPES, pH 7.0
	Column was not equilibrated with Binding Buffer	Equilibrate column with 5mL (5X bed volume) of Lysis/Binding/Wash Buffer with CHAPS
	Sample was not equilibrated in Binding Buffer	Dilute the sample to 0.25-0.5mg/mL total protein in Lysis/Binding/Wash Buffer with CHAPS
	Sample was too concentrated	Use a sample total protein concentration from 0.25 to 0.5mg/mL
	Sample was too viscous	Add DNase to the Lysis/Binding/Wash Buffer
	Insufficient sample/resin mixing	Ensure that sample mixes well at 4°C on a rocking platform for 30 minutes
	Chelating/reducing agent was present in sample	Buffer exchange the sample to eliminate chelating agents such as EDTA, or reducing agents such as β -mercaptoethanol or DTT before column chromatography

Continued on next page

Continued from previous page

Nonspecific protein contamination	Column was overloaded	For best results, use $\leq 4\text{mg}$ of total cellular protein at $\leq 0.5\text{mg/mL}$
	Inadequate washing	Make sure to wash the column 3X for a total of 15mL using the Lysis/Binding/Wash Buffer with CHAPS
2D smearing	Salt was present	Perform a buffer exchange before 2-D analysis

Related Thermo Scientific Products

89884	Pierce Concentrator, 7mL/9K, 25 devices
89885	Pierce Concentrator, 20mL/9K, 25 devices
78425	Halt Protease Inhibitor Single-Use Cocktail EDTA-Free, 100 μl \times 24 microtubes
78420	Halt Phosphatase Inhibitor Cocktail (100X), 1mL
23236	Pierce Coomassie Plus (Bradford) Assay Kit
23227	Pierce BCA Protein Assay Kit
23208	Pre-Diluted Protein Assay Standards: Bovine Serum Albumin (BSA) Set
34080	SuperSignal[®] West Pico Chemiluminescent Substrate
46428	Restore[™] PLUS Western Blot Stripping Buffer, 30mL
34090	CL-XPosure[™] Film, 5 \times 7in (13 \times 18cm), 100/pkg
89983	Pierce SILAC Protein Quantitation Kit – DMEM
89982	Pierce SILAC Protein Quantitation Kit – RPMI 1640
89893	Zeba[™] Spin Desalting Columns, 10mL, 5 columns
78510	T-PER Tissue Protein Extraction Reagent, 500mL
89900	RIPA Buffer, 100mL
89901	RIPA Buffer, 250mL

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.