

Pierce™ Albumin Depletion Kit

85160

2518.0

Number**Description**

85160

Pierce Albumin Depletion Kit, contains sufficient resin and buffer for 24 spin columns**Kit Contents:****Pierce Albumin Depletion Resin**, 10mL, supplied as 50% resin slurry

Binding Capacity: 2.0mg of human serum albumin (HSA) for 200µL of settled resin

Binding/Wash Buffer (25mM Tris, 75mM NaCl; pH 7.5), 11mL**Pierce Spin Columns**, 24 columns**Storage:** Upon receipt store at 4°C. Product shipped with an ice pack.**Introduction**

The Thermo Scientific Pierce Albumin Depletion Kit improves serum component analysis by rapidly removing abundant albumin protein from serum samples. The Pierce Albumin Depletion Resin provided in the kit is a high capacity, immobilized Cibacron Blue dye agarose resin, which binds a variety of species-specific albumins. Each aliquot of 200µL settled resin can process 5-50µL of serum sample in less than 10 minutes by low-speed centrifugation using the Pierce Spin Columns provided in the kit. The kit has been optimized for HSA, but also effectively binds swine, sheep and rabbit serum albumin. With modification of the binding buffer, Pierce Albumin Depletion Resin can bind bovine, calf, goat and rat albumin; however, this product is not for use with mouse albumin.

Important Product Information

- Use a clarified (centrifuged or filtered) protein solution to ensure proper resin flow.
- Albumin binding occurs at pH 6.0-9.0. Samples must be free of excess salt and, therefore, methods for preparing serum samples involving cell clotting may be incompatible with this kit. For best results, use serum separators or filtration when preparing serum samples. To ensure proper ionic strength, desalt samples using a desalting column or dialyze samples vs. a low-salt buffer (25mM Tris, 75mM NaCl; pH 7.5 for human and swine albumin) (25mM Tris, 25mM NaCl; pH 7.5 for bovine and calf albumin).
- The Binding/Wash Buffer of 25mM Tris, 75mM NaCl; pH 7.5 provided in the kit has been optimized for binding HSA and reduced binding of non-albumin proteins. For other albumin species such as bovine, calf, rat or goat, an alternative low-salt Binding/Wash Buffer of 25mM Tris, 25mM NaCl; pH 7.5 must be used in all steps. Some species of albumin may not bind to Pierce Albumin Depletion Resin.
- Differing sample types contain a wide range of albumin concentrations. Take care not to exceed the resin's binding capacity. Each 200µL aliquot of Pierce Albumin Depletion Resin binds 2.0mg of HSA.
- Because of the high albumin concentration present in serum, each Pierce Albumin Depletion Resin aliquot of 200µL settled resin can bind sufficient albumin to process only 5-50µL of serum sample. To process larger serum samples volumes, use more resin volume or multiple spin columns. For serum samples containing high albumin concentrations (e.g., > 50mg/mL), dilute sample with the appropriate Binding/Wash Buffer (as described above). For best results, optimize sample conditions for each specific application.

Additional Materials Required

- 1.5-2.0mL microcentrifuge tubes for sample collection

Procedure for Albumin Depletion Using Centrifugation

The following protocol is an example application for this product. Specific applications may require optimization.

Note: Each 200 μ L aliquot of Pierce Albumin Depletion Resin binds 2.0mg of HSA. Binding capacity can be increased by doubling volume measurements in the procedure below.

1. Shake the resin bottle to resuspend the resin. Using a wide-bore micropipette tip, transfer 400 μ L of the slurry (corresponding to 200 μ L settled resin volume) into a spin column and loosely cap the column.

Note: The amount of resin to use depends on the volume and albumin concentration of the sample being processed.

2. Twist off bottom closure of the spin column and place spin column into a 1.5 to 2.0mL collection tube. Centrifuge at $\sim 12,000 \times g$ for 1 minute to remove excess liquid. Discard flow-through and place spin column back into the same collection tube.

Note: Do not snap off bottom. To remove, twist slightly in one direction followed by the other direction.

3. Add 200 μ L of Binding/Wash Buffer to the spin column.
4. Centrifuge at $12,000 \times g$ for 1 minute. Discard flow-through and place spin column into a new collection tube.
5. Apply 5-50 μ L of albumin-containing sample to resin and incubate for 1-2 minutes at room temperature.

Note: Sample must have < 100 mM salt for proper albumin binding. Dialyze or dilute sample as needed (see Important Product Information Section concerning sample preparation).

6. Centrifuge at $12,000 \times g$ for 1 minute. Re-apply flow-through to spin column and incubate for 1-2 minutes at room temperature to ensure maximal albumin binding.
7. Centrifuge at $12,000 \times g$ for 1 minute. Retain flow-through. Place spin column in a new collection tube.
8. Wash resin to release unbound proteins by adding 50 μ L of Binding/Wash Buffer for each 200 μ L of resin used. Use the appropriate Binding/Wash Buffer for your albumin type (see Important Product Information Section).
9. Centrifuge at $12,000 \times g$ for 1 minute. Retain flow-through. Place spin column in a new collection tube.
10. Repeat Steps 8 and 9 three to four additional times.

Note: Once the procedure is optimized for a particular application, the wash steps can be increased in volume and reduced in number to simplify sample processing.

11. Analyze retained fractions by SDS-PAGE analysis or by protein concentration determination. Combine desired fractions. Concentrate albumin-depleted sample as needed for further processing.

Note: For 2D PAGE or mass spectrometry analysis, albumin-depleted samples must be precipitated, dialyzed or desalted to remove interfering salts (see Related Thermo Scientific Products Section).

Optional Elution of Albumin from Resin

1. To elute bound albumin, wash the resin with 200 μ L of 20mM sodium phosphate, 250mM sodium thiocyanate; pH 7.2.
2. Centrifuge at $\sim 12,000 \times g$ for 1 minute. Retain flow-through. Place spin column in a new collection tube.
3. Repeat Steps 1 and 2 three to four additional times.
4. Analyze retained fractions by SDS-PAGE or protein concentration determination. Discard the used spin column.

Note: Stepwise elution of albumin and other bound proteins using NaCl concentrations between 300mM and 1.5M can be used in place of the 20mM sodium phosphate, 250mM sodium thiocyanate; pH 7.2 elution buffer.

Troubleshooting

| Problem | Possible Cause | Solution |
|---------------------------------|--|---|
| High residual albumin in sample | Albumin concentration exceeded resin binding capacity | Reduce sample amount loaded |
| | | Increase resin amount per column |
| | Salt concentration in sample was too high | Dialyze sample before use or dilute with low ionic strength buffer |
| | Incompatible serum preparation method was used | Dialyze sample before use with low ionic strength buffer |
| | | Prepare serum by alternative method such as filtration |
| | Salt concentration in Binding Buffer is too high for specific albumin type | Use low salt (e.g., < 25mM NaCl) or no salt 25mM Tris, pH 7.5 for Binding/Wash Buffer Note: This may increase nonspecific protein binding |
| Non-compatible albumin type | Ensure albumin type is compatible with the depletion kit | |

Related Thermo Scientific Products

| | |
|--------------|--|
| 69705 | Pierce Spin Columns – Screw Cap, 25 units |
| 89877 | Zeba™ Micro Spin Desalting Columns, 7K MWCO, 75µL |
| 89882 | Zeba Spin Desalting Columns, 7K MWCO, 0.5mL |
| 88512 | Protein Concentrators PES, 3K MWCO, 0.5mL |
| 88400 | Slide-A-Lyzer™ MINI Dialysis Device, 3.5K MWCO, 0.5mL |
| 24615 | Imperial™ Protein Stain, 1L |
| 24590 | GelCode™ Blue Stain Reagent, 500mL |
| 23227 | Pierce BCA Protein Assay Kit |
| 23235 | Micro BCA Assay Reagent Kit |
| 22662 | Pierce 660nm Protein Assay Kit |

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts (“Documentation”). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER’S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER’S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.