



Package contents

Catalog No.	Size
K8050-10	10 reaction
K8050-20	20 reaction

[Kit contents](#)



Storage conditions

Store Box 1 at -30°C to -10°C
 Store Box 2 at room temperature
 Store Box 3 at -85°C to -68°C



Required materials

- Template: DNA
- PCR thermocycler
- Microcentrifuge
- Agarose gel electrophoresis equipment
- 42°C water bath
- 37°C incubator



Timing

- Perform PCR: 30 minutes to 3 hours
- Purify PCR product: 15–20 minutes
- TOPO™ Cloning reaction: 5 minutes
- Transformation reaction: 2 hours



Product description

- The TOPO™ XL-2 Complete PCR Cloning Kit provides all the components needed for long fragment PCR cloning in a single kit format.
- Platinum™ SuperFi™ Green PCR Master Mix is used to generate long PCR products with high specificity and fidelity.
- The PureLink™ Quick Gel Extraction and PCR Purification Kit is used to obtain high yields of purified PCR product despite large fragment sizes.
- One Shot™ OmniMAX™ 2 T1^R Chemically Competent *E. coli* have superior transformation efficiency to ensure sufficient transformants for evaluation.



Online resources

Visit our product page for additional information and protocols. For support, visit thermofisher.com/support.

Protocol outline

- Produce blunt-ended PCR product
- Purify PCR product
- Perform TOPO™ Cloning reaction
- Transform TOPO™ Cloning reaction
- Select and analyze colonies

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.thermofisher.com/us/en/home/global/terms-and-conditions.html.

If you have any questions, please contact Life Technologies at www.thermofisher.com/support.

Important licensing information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.








Disclaimer: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Corporate entity: Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.






PCR protocol for TOPO™ XL-2 Complete PCR Cloning Kit

The following procedure is used to produce blunt-ended PCR products for TOPO™ Cloning. See the [Platinum™ SuperFi™ Green PCR Master Mix User Guide](#) for additional information on thermocycling conditions and performing PCR. For details on purification of PCR products, see the [TOPO™ XL-2 Complete PCR Cloning Kit User Guide](#).

Step	Action	Procedure details																								
1 	Prepare PCR mix	<p>Combine the following components in a sterile PCR tube.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>50-µL reaction</th> </tr> </thead> <tbody> <tr> <td>Water, nuclease-free</td> <td>to 50 µL</td> </tr> <tr> <td>Platinum™ SuperFi™ Green PCR Master Mix (2X)</td> <td>25 µL</td> </tr> <tr> <td>10 µM forward primer (0.5 µM final concentration)</td> <td>2.5 µL</td> </tr> <tr> <td>10 µM reverse primer (0.5 µM final concentration)</td> <td>2.5 µL</td> </tr> <tr> <td>Template DNA (5–50 ng of gDNA or 1 pg to 10 ng of plasmid DNA)</td> <td>varies</td> </tr> <tr> <td>(Optional) SuperFi™ GC Enhancer (5X)</td> <td>10 µL</td> </tr> </tbody> </table>	Component	50-µL reaction	Water, nuclease-free	to 50 µL	Platinum™ SuperFi™ Green PCR Master Mix (2X)	25 µL	10 µM forward primer (0.5 µM final concentration)	2.5 µL	10 µM reverse primer (0.5 µM final concentration)	2.5 µL	Template DNA (5–50 ng of gDNA or 1 pg to 10 ng of plasmid DNA)	varies	(Optional) SuperFi™ GC Enhancer (5X)	10 µL										
Component	50-µL reaction																									
Water, nuclease-free	to 50 µL																									
Platinum™ SuperFi™ Green PCR Master Mix (2X)	25 µL																									
10 µM forward primer (0.5 µM final concentration)	2.5 µL																									
10 µM reverse primer (0.5 µM final concentration)	2.5 µL																									
Template DNA (5–50 ng of gDNA or 1 pg to 10 ng of plasmid DNA)	varies																									
(Optional) SuperFi™ GC Enhancer (5X)	10 µL																									
2 	Amplify PCR product	<p>Use the following cycling parameters for amplicons <10 kb.</p> <table border="1"> <thead> <tr> <th>Step</th> <th>Temperature</th> <th>Time</th> <th>Cycles</th> </tr> </thead> <tbody> <tr> <td>Initial Denaturation</td> <td>98°C</td> <td>30 seconds</td> <td>1X</td> </tr> <tr> <td>Denaturation</td> <td>98°C</td> <td>5–10 seconds</td> <td rowspan="3">25–35X</td> </tr> <tr> <td>Anneal</td> <td> varies</td> <td>10 seconds</td> </tr> <tr> <td>Extend</td> <td>72°C</td> <td>15–30 seconds/kb</td> </tr> <tr> <td rowspan="2">Final Extension</td> <td>72°C</td> <td>5 minutes</td> <td rowspan="2">1X</td> </tr> <tr> <td>4°C</td> <td>Hold</td> </tr> </tbody> </table>	Step	Temperature	Time	Cycles	Initial Denaturation	98°C	30 seconds	1X	Denaturation	98°C	5–10 seconds	25–35X	Anneal	 varies	10 seconds	Extend	72°C	15–30 seconds/kb	Final Extension	72°C	5 minutes	1X	4°C	Hold
Step	Temperature	Time	Cycles																							
Initial Denaturation	98°C	30 seconds	1X																							
Denaturation	98°C	5–10 seconds	25–35X																							
Anneal	 varies	10 seconds																								
Extend	72°C	15–30 seconds/kb																								
Final Extension	72°C	5 minutes	1X																							
	4°C	Hold																								
3 	Analyze PCR product	Analyze 5–10 µL by agarose gel electrophoresis to verify the size, quality, and quantity of the PCR product.																								
4 	Purify PCR product	<ul style="list-style-type: none"> If you have a single discrete band, prepare the PCR product by column purification or gel purification. Note: Gel purification results in optimal cloning efficiency. If you do not have a single, discrete band, isolate the desired PCR product by gel purification. 																								

Cloning and transformation protocol for TOPO™ XL-2 Complete PCR Cloning Kit

The following procedure is used to perform TOPO™ Cloning and transformation into One Shot™ OmniMAX™ 2 T1^R Chemically Competent *E. coli*. For details on analyzing transformants, see the [TOPO™ XL-2 Complete PCR Cloning Kit User Guide](#).

Step	Action	Procedure details								
1 	Prepare TOPO™ reaction mix	<p>Combine the following components in a sterile PCR tube. Use a 1:1 molar ratio of insert to vector.</p> <table border="1"> <thead> <tr> <th>Component</th> <th>6-μL reaction</th> </tr> </thead> <tbody> <tr> <td>Column or gel purified PCR product</td> <td>up to 4 μL</td> </tr> <tr> <td>pCR-XL-2-TOPO™ Vector (10 ng/μL)</td> <td>1 μL</td> </tr> <tr> <td>Salt Solution</td> <td>1 μL</td> </tr> </tbody> </table>	Component	6-μL reaction	Column or gel purified PCR product	up to 4 μL	pCR-XL-2-TOPO™ Vector (10 ng/μL)	1 μL	Salt Solution	1 μL
Component	6-μL reaction									
Column or gel purified PCR product	up to 4 μL									
pCR-XL-2-TOPO™ Vector (10 ng/μL)	1 μL									
Salt Solution	1 μL									
2  5 min	Incubate TOPO™ reaction mix	<ol style="list-style-type: none"> Gently mix, then briefly centrifuge. Incubate at room temperature for 5 minutes. Place the tube on ice. 								
3  42°C 30 s	Transform competent cells	<ol style="list-style-type: none"> Add 2 μL of the TOPO™ Cloning reaction into a vial of competent cells. Incubate for 30 minutes on ice. Heat-shock the cells for 30 seconds in a 42°C water bath. Immediately place the tubes on ice and incubate for 2 minutes. Add 250 μL of room temperature S.O.C. medium. Cap the tube tightly and shake the tube horizontally at 225 rpm for 1 hour at 37°C. 								
4 	Plate transformed cells	<ol style="list-style-type: none"> Spread 50–150 μL from each TOPO™ Cloning transformation reaction on a pre-warmed LB plate containing kanamycin or ampicillin, and 1 mM IPTG. Incubate plates overnight at 37°C. 								
5 	Analyze transformants	<ul style="list-style-type: none"> Analyze transformants by colony PCR. Analyze transformants by restriction enzyme digestion. 								