

# Pierce Antibody Biotinylation Kit for IP

 Pub. No. MAN0016152  
 Pub. Part No. 2162590.1  
 Rev. A.0

**90407**


Number	Description
90407	<p><b>Pierce Antibody Biotinylation Kit for IP</b>, sufficient reagents for 8 labeling reactions, each containing 50-200µg of antibody in 100µL reaction volumes</p> <p><b>Kit Contents:</b></p> <p><b>EZ-Link NHS-PEG<sub>4</sub>-Biotin, No-Weigh Format</b>, 8 × 0.5mg microtubes            Molecular Weight: 588.67            Spacer Arm: 29Å</p> <p><b>20X Phosphate Buffered Saline</b>, 15mL</p> <p><b>Zeba Spin Desalting Column (7K MWCO), 0.5mL</b>, 8 columns, for 30-130µL samples</p> <p><b>Storage:</b> Upon receipt store at 4°C. Product is shipped with an ice pack.</p>

## Introduction

The Thermo Scientific™ Pierce™ Antibody Biotinylation Kit for IP provides optimized reagents for labeling antibodies and desalting columns for purifying the labeled molecule. Each reaction is sufficient for labeling 50-200µg of antibody in 100µL reaction volumes. The hydrophilic polyethylene oxide (PEO), also called polyethylene glycol (PEG), spacer arm imparts water solubility that is transferred to the biotinylated molecule. Consequently, antibodies labeled with NHS-PEG<sub>4</sub>-Biotin exhibit less aggregation when stored in solution compared to antibodies labeled with reagents having only hydrocarbon spacers. This kit is specifically optimized to label antibodies for downstream immunoprecipitation (IP) applications involving streptavidin magnetic or agarose supports while having little to no effect on antigen binding. The NHS-PEG<sub>4</sub>-Biotin is provided in Thermo Scientific™ No-Weigh™ Format as convenient single-use microtubes, eliminating difficulties associated with weighing small quantities of reagent.

Biotin is a small, naturally occurring vitamin that binds with high affinity to avidin and avidin-like proteins. Biotinylated antibodies typically retain biological activity because the biotin group is relatively small. An antibody conjugated with several biotin molecules can interact rapidly and tightly with streptavidin, making this a good platform for IP applications.

N-Hydroxysuccinimide (NHS) esters are the most popular biotinylation reagents. In pH 7-9 buffers, NHS esters react efficiently with primary amino groups (-NH<sub>2</sub>) by nucleophilic attack, forming an amide bond and releasing the NHS (Figure 1). Proteins typically have many sites for labeling, including the primary amine in the side chain of lysine (K) residues and the N-terminus of each polypeptide.



**Figure 1. Reaction of NHS-PEG<sub>4</sub>-Biotin with primary amine on an antibody.** The leaving group and any non-reacted biotin reagent are removed during the desalting step.

## Important Product Information

- Use reconstituted NHS-PEG<sub>4</sub>-Biotin immediately. The NHS-ester moiety readily hydrolyzes and becomes non-reactive; therefore, do not prepare solutions for storage. Discard any unused reconstituted reagent.
- NHS-PEG<sub>4</sub>-Biotin is moisture-sensitive. Immediately before use, puncture the microtube foil with a pipette tip, add solvent and mix by pipetting up and down. After use, cut off the used microtube and discard. Store the microtube strip at 4°C in the foil pouch provided.
- Avoid buffers containing primary amines (e.g., Tris or glycine), as these will compete with the intended reaction (Figure 1). If necessary, dialyze or desalt to exchange the antibody into phosphate-buffered saline (PBS).
- Antibodies must be carrier-free (i.e., free of BSA and/or gelatin). Any carrier will compete with the intended reaction. To remove carrier, Protein A, G, or A/G resins may be used, as well as other commercially available antibody clean-up kits.
- This kit has been successfully used to label multiple mouse IgG<sub>1</sub>, mouse IgG<sub>2</sub>, rabbit monoclonal and polyclonal antibodies without impairing antigen binding. The number of biotin groups per antibody molecule ranges from 3-7, as determined by the Thermo Scientific™ Fluorescence Biotin Quantitation Kit (Product No. 46610).

## Procedure for Biotinylating Antibodies

**Note: The maximum total reaction volume must not exceed 110µL; this is the volume limit on the Thermo Scientific™ Zeba™ Spin Desalting Column.**

### A. Calculations

1. Calculate millimoles of NHS-PEG<sub>4</sub>-Biotin to add to the reaction for a 40-fold molar excess:

$$\mu\text{g IgG} \times \frac{1\text{mg}}{1000\ \mu\text{g}} \times \frac{1\ \text{mmol IgG}}{150,000\ \text{mg IgG}} \times \frac{40\ \text{mmol Biotin}}{1\ \text{mmol IgG}} = \text{mmol Biotin}$$

- 40 = Recommended molar fold excess of biotin per protein sample
- 150,000 = Molecular weight of IgG

2. Calculate microliters of 8.5mM NHS-PEG<sub>4</sub>-Biotin (prepared in Step B.3) to add to the reaction:

$$\text{mmol Biotin} \times \frac{589\ \text{mg}}{1\ \text{mmol Biotin}} \times \frac{100\ \mu\text{L}}{0.5\ \text{mg}} = \mu\text{L Biotin Solution}$$

- 589 = Molecular weight of NHS-PEG<sub>4</sub>-Biotin
- 100 = Microliters of solvent in which 0.5mg of NHS-PEG<sub>4</sub>-Biotin is dissolved to make 8.5mM

**Example:** For 85µg of a 1mg/mL IgG (150,000 MW) solution, 2.7µL of 8.5mM NHS-PEG<sub>4</sub>-Biotin will be added.

$$85\ \mu\text{g IgG} \times \frac{1\text{mg}}{1000\ \mu\text{g}} \times \frac{1\ \text{mmol IgG}}{150,000\ \text{mg IgG}} \times \frac{40\ \text{mmol Biotin}}{1\ \text{mmol IgG}} = 0.000023\text{mmol Biotin}$$

$$0.000023\text{mmol Biotin} \times \frac{589\ \text{mg}}{1\text{mmol Biotin}} \times \frac{100\ \mu\text{L}}{0.5\ \text{mg}} = 2.7\ \mu\text{L Biotin Solution}$$

3. Calculate the volume of 1X PBS to bring the total reaction to 100µL:  
100µL — (volume of antibody + volume of biotin) = volume of 1X PBS

## B. Biotin-labeling Reaction

**Note: The antibody must be in an amine-free buffer (preferably PBS) without carrier to perform this labeling reaction (see Important Product Information).**

1. Dilute the 20X PBS to 1X (10mM sodium phosphate, 0.15M NaCl, pH 7.5) with ultrapure water. The 1X PBS will be used to equilibrate the Zeba Spin Desalting Column.

**Note:** If precipitate formed in the 20X PBS, warm solution to 37°C in a water bath and vortex before use. The PBS can be stored at 1X if desired.

2. Cut off one microtube of NHS-PEG<sub>4</sub>-Biotin from the No-Weigh Microtube Strip. Return the unused strip of microtubes to its pouch and store desiccated at 4°C.
3. With a pipette tip, puncture the foil top on the biotin reagent microtube. To prepare 8.5mM solution, add 100µL of 1X PBS to the tube and mix by pipetting up and down several times.
4. Combine antibody, 1X PBS and NHS-PEG<sub>4</sub>-Biotin (as calculated in Section A) and mix by gently pipetting up and down.
5. Incubate the reaction at room temperature for 30 minutes.

## C. Buffer Exchange and Excess Biotin Removal

1. Remove bottom closure of Zeba Spin Desalting Column and loosen cap (do not remove cap).
2. Place column in a 1.5-2.0mL collection tube. Centrifuge at 1,500 × g for 1 minute to remove storage solution.
3. Place a mark on the side of the column where the compacted resin is slanted upward. Place column in the microcentrifuge with the mark facing outward in all subsequent centrifugation steps.
4. Add 300µL of 1X PBS on top of the resin bed. Centrifuge at 1,500 × g for 1 minute to remove buffer.
5. Repeat Step 4 two additional times, discarding buffer from the collection tube.
6. Place column in a new collection tube, remove cap and apply labeling reaction sample (from Section B) to the top of the compact resin bed.
7. Centrifuge at 1,500 × g for 2 minutes to collect the sample. Discard desalting column after use.
8. Store biotinylated antibody at 4°C for < 1 month. For longer periods, store at -20°C or -80°C.

## Troubleshooting

Problem	Possible Cause	Solution
Low or no biotinylation	Buffer contained primary amines	Buffer exchange the antibody into a non-amine-containing buffer such as the PBS provided by desalting columns or dialysis
	NHS-PEG <sub>4</sub> -Biotin was hydrolyzed	Use reagent immediately upon reconstitution
	Carrier protein was present in the antibody solution	Remove carrier protein before biotinylation by using Protein A, G or A/G resin or an antibody clean-up kit. This will reduce competition for labeling

## Additional Information

Please visit our website for additional information including the following:

- Tech Tip #43: Protein stability and storage

---

**Related Thermo Scientific Products**

<b>90408</b>	<b>Pierce™ MS-Compatible Magnetic IP Kit (Streptavidin)</b>
<b>88816</b>	<b>Pierce™ Streptavidin Magnetic Beads, 1mL</b>
<b>89882</b>	<b>Zeba Spin Desalting Columns, 7K MWCO, 0.5mL</b>
<b>69570</b>	<b>Slide-A-Lyzer™ MINI Dialysis Device, 10K MWCO, 0.1mL</b>
<b>46610</b>	<b>Fluorescence Biotin Quantitation Kit</b>
<b>20423</b>	<b>Pierce™ Protein A/G Plus Agarose Resin</b>

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at [thermofisher.com](http://thermofisher.com). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2016 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.