

Pierce™ Quantitative Peroxide Assay Kits

23280 23285

0526.2

Number	Description
23280	Pierce Quantitative Peroxide Assay Kit, aqueous-compatible formulation Kit contents:
23281	Reagent A , 1mL, Composition: 25mM ammonium ferrous (II) sulfate, 2.5M H ₂ SO ₄
23282	Reagent B , 2 × 50mL, Composition: 100mM sorbitol, 125μM xylenol orange in water
23285	Pierce Quantitative Peroxide Assay Kit, lipid-compatible formulation Kit contents:
23281	Reagent A , 1mL, Composition: 25mM ammonium ferrous (II) sulfate, 2.5M H ₂ SO ₄
23283	Reagent C , 2 × 50mL, Composition: 4mM BHT, 125μM xylenol orange in methanol

Xylenol orange is (*o*-cresolsulfonephthalein-3'-3'-bis-[methyliminodiacetic acid sodium salt]).

CAS numbers for chemical components of formulations:

methyl alcohol: #67-56-1	sorbitol: #50-70-4	xylenol orange: #1611-35-4
BHT: #28-37-0	sulfuric acid: #7664-93-9	ferrous sulfate: #7783-85-9

Storage: Upon receipt store kit at 4° C. Product shipped at ambient temperature.

Introduction

The Thermo Scientific Pierce Quantitative Peroxide Assay Kits detect peroxide based on oxidation of ferrous to ferric ion in the presence of xylenol orange. One assay kit has an aqueous-compatible formulation that includes sorbitol, which provides sensitivity enhancement. The other assay kit has a lipid-compatible formulation that may be used without extraction; however, because it does not contain sorbitol, the assay may be less sensitive. Most proteins do not interfere with these assays, although some metal chelators may require use of a blank (i.e., control). Pierce Quantitative Peroxide Assays are suitable for H₂O₂ measurement in the following applications: quantitating detergent peroxide, monitoring cellular activity, determining protein glycation and measuring peroxide accumulation in lipid.

In these assays, hydroperoxides convert the Fe²⁺ to Fe³⁺ at acidic pH. With the aqueous-compatible formulation, peroxide first reacts with sorbitol, converting it to a peroxy radical, which in turn initiates Fe²⁺ oxidation to Fe³⁺. In the lipid-compatible formulation, the peroxide converts the Fe²⁺ to Fe³⁺ directly. In a sulfuric acid solution, the Fe³⁺ complexes with the xylenol orange dye to yield a purple product with maximum absorbance at 560nm. The molar extinction coefficient of the xylenol orange-Fe³⁺ complex is 1.5 × 10⁴ M⁻¹cm⁻¹ in 25mM H₂SO₄ at room temperature. The complex in both assays absorbs strongest between 540 and 580nm but absorbs measurably up to 620nm. When using microplates, the best wavelength for measurement is 595nm (best signal-to-noise). The maximum absorbance of the assay components before addition of peroxide is at 430nm.

Important Product Information

- Excessive H₂O₂ (above 1mM) can result in low absorbance measurements caused by a bleaching effect on the dye. Ensure that results are accurate by preparing a 1:100 dilution of the sample and perform the assay in the same manner as the undiluted sample. If the absorbance of the diluted sample is higher than or similar to the original reading, then there is excessive peroxide in the sample.
- The protocols use a 1:10 (sample:WR) ratio. This ratio can be changed to accommodate samples with relatively high levels of peroxide [e.g., use 1:100 (sample:WR)] as long as the dilution is accounted for when comparing to the assay standards.
- The working range of the assay is from 1μM to 1mM peroxide in a sample assayed using the Standard Procedure [1:10 (sample:WR) ratio]. Above 1mM, a bleaching effect occurs, causing a decrease in absorbance and loss of linearity in the assay.
- The standard curve is not linear over the entire assay range (see Figure 1). Use a quadratic or best-fit curve to the standard points. Alternatively, use a linear fit for individual 100μM increments (e.g., from 1-100μM, 101-200μM, etc.); each incremental standard curve will be more horizontal in slope than the previous one until a plateau is reached at about 1mM.

Reagent Preparation for Aqueous-Compatible Procedure

Working Reagent (WR) Mix 1 volume of Reagent A with 100 volumes of Reagent B.

Notes:

- Prepare at least 1mL of WR for each sample and standard replicate to be assayed in cuvettes. Prepare at least 200μL of WR for each sample and standard replicate to be assayed in a microplate.
- The WR is stable for at least 12 hours.

Peroxide Standards Serially dilute a 30% (8.8M) hydrogen peroxide stock solution to achieve 8-10 standards in the concentration range of 1-1000μM (Figure 1).

Notes:

- A 30% hydrogen peroxide solution can be standardized using $43.6\text{M}^{-1}\text{cm}^{-1}$ as the molar extinction coefficient for H₂O₂ at 240nm. Then, prepare the high standard (1000μM) by adding 100μL of 30% hydrogen peroxide to 880mL of ultrapure water or buffer.

Aqueous-Compatible Procedure

- Add 10 volumes of WR to 1 volume of sample. For example, in a microplate, add 20μL of sample (or standard) and 200μL of WR to each well.
- Mix and incubate assay reactions for 15-20 minutes at room temperature. This incubation step is necessary for the reaction to reach an endpoint. Once formed, the complex is relatively stable, but absorbance measurements are optimal if taken on the same day the reaction is performed.
- Measure absorbance at 560-620nm (560nm is optimal in a spectrophotometer; 595nm is optimal if using a plate reader).
- Calculate the concentration of peroxide in the sample by reference to its assay absorbance compared to the standard curve (see the Important Product Information Section for correct curve-fitting procedures for the standards).

Reagent Preparation for Lipid-Compatible Procedure

Working Reagent (WR) Mix 1 volume of Reagent A with 100 volumes of Reagent C.

Notes:

- Prepare at least 1mL of WR for each sample and standard replicate to be assayed in cuvettes. Prepare at least 200 μ L of WR for each sample and standard replicate to be assayed in a microplate.
- The WR is stable for at least 12 hours.

Peroxide Standards Serially dilute a 30% (8.8M) hydrogen peroxide stock solution in methanol to make 8-10 standards in the concentration range of 1-1000 μ M (Figure 1). Refer to the Aqueous-Compatible Procedure for more information on diluting the peroxide standard.


For the “TCEP Reference,” dilute Thermo Scientific Bond-Breaker TCEP Solution, Neutral pH (Product No. 77720; 0.5M) (see Related Thermo Scientific Products) to make a 10mM stock solution by diluting 1 part of the TCEP Solution with 49 parts of methanol.

Lipid-Compatible Procedure

- If samples may contain metals, then produce “test” or methanol-only replicates, as well as “TCEP Reference” replicates. Subtract the absorbance values of the “TCEP Reference” from that of the “test” replicates to eliminate the metal effects. Create both types of replicates for each peroxide concentration standard. Use the absorbance difference at each standard concentration as the value for creating the standard curve. If the absorbance of the “test” and “TCEP Reference” are equal, then metals are not present and the “TCEP Reference” steps do not need to be performed.
 - Plasma hydroperoxide levels may be determined using the lipid-compatible formulation of the Pierce Quantitative Peroxide Assay Kit (Prod. No. 23285). Since plasma is quite low in peroxide content, the assay absorbance measurement often will be less than about 0.100. Because of endogenous plasma iron, a “TCEP Reference” must be assayed as described. This protocol is based on the work by Nourooz-Zadeh, J., *et.al.*¹
 - If using the assay with lipid samples that are not extracted, prepare a blank that omits Reagent A. Subtracting the absorbance for this blank from the sample tested in WR controls endogenous iron interferences. A blank may also be necessary if the sample contains other transition metals or a protein having chelating properties or strong absorbance characteristics at the wavelengths used for measurement.
1. For each sample to be tested, pipette 90 μ L of sample into microcentrifuge tubes.
Note: For each plasma sample or other metal-containing sample to be tested, pipette 90 μ L of sample into microcentrifuge tubes labeled “test” and “TCEP Reference.” Use a minimum of two replicates for each “test” and “TCEP Reference.” In this case, each peroxide standard must also be created as “test” and “TCEP Reference” replicates.
 2. Add 10 μ L of methanol to each “test” tube and 10 μ L of 10mM TCEP in methanol to each “TCEP Reference” tube.
 3. Vortex tubes and incubate for 20-30 minutes at room temperature.
 4. To each tube, add 900 μ L of WR.
 5. Vortex tubes and incubate for 20-30 minutes.
 6. Centrifuge at 12,000 \times g for 5-10 minutes.
 7. Transfer supernatant to a cuvette (900 μ L) or microwell plate (300 μ L).
 8. Measure absorbance at 560-620nm (560nm is optimal in a spectrophotometer or 595nm if using a plate reader) as described in Aqueous Compatible Procedure.
 9. For samples containing metals, subtract the “TCEP Reference” absorbance from the “test” absorbance and generate the hydrogen peroxide standard curve using the difference at each peroxide concentration. Calculate the concentration of the samples against this standard curve.

Additional Information

Typical standard curve over the entire assay range



Standard Curve at low end of assay range




Figure 1. Thermo Scientific Pierce Quantitative Peroxide Assay absorbance response to hydrogen peroxide concentration. A 1mM (1000μM) solution of hydrogen peroxide was initially made by diluting a 30% H₂O₂ stock 1:9000 (11μL of 30% H₂O₂ into 100mL of ultrapure water). This sample was then serially diluted with ultrapure water 1:2 (100μL of ultrapure water + 100μL of the previous dilution) for a total of seven samples. Ultrapure water was used as the blank. In microplate wells, 200μL of the Working Reagent was added to 20μL of the diluted H₂O₂ standards. Samples were mixed and incubated for 15 minutes at room temperature. Absorbances were measured at 595nm using a plate reader. The blank value was subtracted from all sample measurements.

Related Thermo Scientific Products

77720	Bond-Breaker™ TCEP Solution, Neutral pH
20490	TCEP•HCl, Tris[2-carboxyethyl]phosphine hydrochloride, 1g
15041	96-Well EIA Plates, corner-notched, clear, 100 per package
15031	8-Well EIA Strip Plates, corner-notched, clear, 100 per package
15082	MicroTube Racked System

General References

- Jiang, Z.-Y., Hunt, J.V. and Wolff, S.P. (1992). Ferrous ion oxidation in the presence of xylenol orange for detection of lipid hydroperoxide in low-density lipoprotein. *Anal Biochem* **202**:384-9.
- Jiang, Z.-Y., Woollard, A.C.S. and Wolff, S.P. (1990). Hydrogen peroxide production during experimental protein glycation. *FEBS* **268(1)**:69-71.
- Jiang, A.-Y., Woollard, A.C.S. and Wolff, S.P. (1991). Lipid hydroperoxide measurement by oxidation of Fe²⁺ in the presence of xylenol orange. Comparison with the TBA assay and an iodometric method. *Lipids* **26**:853-6.

Cited Reference

- Nourooz-Zadeh, J., Tajaddini-Sarmadi, J. and Wolff, S.P. (1994). Measurement of plasma hydroperoxide concentrations by the ferrous oxidation - Xylenol Orange assay. *Anal Biochem* **220**:403-9.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.