

Active Rac1 Pull-Down and Detection Kit

16118

2271.1

Number	Description
16118	Active Rac1 Pull-Down and Detection Kit , contains sufficient reagents for 30 pull-down reactions

Kit Contents:**Box 16118X (these items ship together on dry ice; upon receipt store at -20°C):**

GST-human Pak1-PBD, 600µg, contains 1-2mg/mL in 25mM Tris•HCl, pH 7.2, 150mM NaCl and 10% glycerol; ~35kDa; GST-human Pak1-PBD interacts with Rac1 from human and mouse, and possibly from all mammalian species, store at -20°C

Anti-Rac1 Antibody, 50µL (5 units), mouse monoclonal IgG_{2b}; Anti-Rac1 antibody reacts with Rac1 of human, dog, rat, mouse and chicken; store at -20°C. Note: One unit of Anti-Rac1 antibody is defined as the amount of antibody required to efficiently detect Rac1 in 20µg NIH3T3 whole cell lysate by Western blotting (8.5 × 7.5cm membrane).

100X GTPγS, 50µL, 10mM in sterile water, store at -20°C

100X GDP, 50µL, 100mM in sterile water, store at -20°C

Box 16118Y (these items ship together with an ice pack; upon receipt store at 4°C):

Glutathione Resin, 3.0mL, supplied as 50% slurry containing 0.05% sodium azide, store at 4°C

1X Lysis/Binding/Wash Buffer, 100mL, 25mM Tris•HCl, pH 7.2, 150mM NaCl, 5mM MgCl₂, 1% NP-40 and 5% glycerol, store at 4°C

2X SDS Sample Buffer, 1.5mL, 125mM Tris•HCl, pH 6.8, 2% glycerol, 4% SDS (w/v) and 0.05% bromophenol blue, store at 4°C

Spin Cups, 30 each, maximum volume 850µL, store at room temperature or 4°C

Collection Tubes, 90 each, store at room temperature or 4°C

Table of Contents

Introduction	1
Important Product Information	2
Additional Materials Required.....	2
Procedure for the Active Rac1 Pull-Down and Detection	3
A. Cell Lysis.....	3
B. <i>In vitro</i> GTPγS or GDP Treatment (Optional).....	3
C. Affinity Precipitation of Activated Rac1	3
D. Western Blot Analysis.....	4
Troubleshooting.....	5
Additional Information	5
Related Thermo Scientific Products	6
General References.....	6

Introduction

The Thermo Scientific Active Rac1 Pull-Down and Detection Kit is a simple and fast tool to monitor Rac1 small GTPase activation. The kit provides a GST-fusion protein containing the p21-binding domain (PBD) of human p21-activated protein kinase 1 (Pak1) along with glutathione agarose resin to specifically pull down active Rac1 and an anti-Rac1 antibody for Western blot detection. Also included are two control nucleotides, GTP γ S and GDP, which can be used to generate positive and negative control lysates, respectively. Each kit is functionally tested to ensure component performance.

GTPase Rac1 plays an important role in the organization of actin filament networks and in membrane ruffling. Rac1 (~22kDa) has been implicated in cell transformation, apoptosis and migration. Like other small GTPases, Rac1 is active when bound to GTP and is inactive when bound to GDP.

Important Product Information

- Rac1-GTP is quickly hydrolyzed to Rac1-GDP; use fresh lysate for each assay.
- Lysis/Binding/Wash buffer is compatible with Thermo Scientific Pierce BCA (Product No. 23227) and Pierce 660nm (Product No. 22660) Protein Assays but not the Bradford Protein Assay.
- For best results always use protease inhibitors for cell lysis, and keep lysates on ice between steps.
- For optimal pilot experiments, use from 500 μ g to 1mg of total lysate per assay.
- For best results when performing the Western blotting procedure, use Pierce[®] Goat Anti-Mouse IgG (H+L), Peroxidase Conjugated (Product No. 31430) and Thermo Scientific SuperSignal West Pico Chemiluminescent Substrate (Product No. 34080). (Refer to Additional Information section, Figure 1.) If similar products from other vendors are used, the Western blotting procedure must be optimized.

Additional Materials Required

- Protease inhibitors (e.g., Thermo Scientific Halt Protease Inhibitor Single-Use Cocktail – EDTA Free, Product No. 78425)
- Pierce BCA Protein Assay Reagent (Product No. 23227) or Pierce 660nm Protein Assay (Product No. 22660)
- β -mercaptoethanol (Product No. 35602) or dithiothreitol (DTT) (Product No. 20291)
- Polyacrylamide gel, 12% or 4-20% (Thermo Scientific Precise Protein Gels; see catalog or website)
- Nitrocellulose (Product No. 88014) or PVDF (Product No. 88585) membrane
- Tris-buffered saline (TBS; 25mM Tris•HCl, pH 7.5, 150mM NaCl; Product No. 28379 or 28358)
- Tween[®]-20 Detergent (Product No. 28320)
- BSA, Fraction V
- Nonfat Dry Milk
- Pierce Goat Anti-Mouse IgG-Horseradish Peroxidase Conjugate (Product No. 31430)
- SuperSignal[®] West Pico Chemiluminescent Substrate (Product No. 34080)
- Thermo Scientific CL-XPosure X-ray Film (Product No. 34090 or 34091) or a CCD camera
- 0.5M EDTA, pH 8.0
- 1M MgCl₂
- Sodium azide (NaN₃)
- Electrophoresis apparatus
- Variable-speed Bench-top Microcentrifuge

Procedure for the Active Rac1 Pull-Down and Detection

A. Cell Lysis

Note: Add protease inhibitors to Lysis/Binding/Wash Buffer before use.

- **For adherent cells:**

1. Carefully remove the culture medium and gently rinse the cells once with ice-cold TBS.
2. Add 0.5-1.0mL Lysis/Binding/Wash Buffer per 75cm² flask or 0.3-0.5mL Lysis/Binding/Wash Buffer per 100mm plate with cells at 80-90% confluency.
3. Scrape the cells and transfer to a microcentrifuge tube. Vortex the tube briefly and incubate on ice for 5 minutes.
4. Centrifuge at 16,000 × *g* at 4°C for 15 minutes.
5. Transfer the supernatant (total lysate) to a new tube.

- **For non-adherent cells:**

1. Pellet cells from one 75cm² flask (approx. 1-2 × 10⁷ cells) at 100 × *g* for 5 minutes and then resuspend cells in 10mL ice-cold TBS.
2. Pellet the cells at 100 × *g* for 5 minutes and carefully remove TBS.
3. Add 0.5-1.0mL Lysis/Binding/Wash Buffer to the cell pellet and resuspend the pellet.
4. Transfer the sample to a microcentrifuge tube and incubate on ice for 5 minutes.
5. Centrifuge at 16,000 × *g* at 4°C for 15 minutes.
6. Transfer the supernatant (total lysate) to a new tube.

B. *In vitro* GTPγS or GDP Treatment (Optional)

Perform the following treatments, GTPγS (positive control) and GDP (negative control) to ensure the pull-down procedures are working properly. Use 500μg of cell lysate for each treatment. For best results, aliquot GTPγS and GDP at first use to minimize freeze/thaw cycles.

1. For 500μL lysate, add 10μL 0.5M EDTA pH 8.0 (for a final concentration of 10mM) and vortex the sample.
2. Add 5μL of 10mM GTPγS (for a final concentration of 0.1mM) or 5μL 100mM GDP (for a final concentration of 1mM) and vortex the sample.
3. Incubate the mixture at 30°C for 15 minutes with constant agitation.
4. Terminate the reaction by placing the sample on ice and adding 32μL of 1M MgCl₂ (for a final concentration of 60mM) and vortex the sample.

C. Affinity Precipitation of Activated Rac1

1. Save a sample of the cell lysate for protein assay using the Pierce BCA or 660nm Protein Assay.
2. Place a spin cup into a collection tube for each sample.
3. Swirl the bottle of Glutathione Resin to thoroughly resuspend the agarose beads. Add 100μL of the 50% resin slurry to the spin cup with collection tube. Centrifuge the tubes at 6000 × *g* for 10-30 seconds.
4. Discard the flow-through. Add 400μL of Lysis/Binding/Wash Buffer to each tube with resin. Invert the tubes gently several times. Centrifuge the tubes at 6000 × *g* for 10-30 seconds. Discard the flow-through.
5. Thaw the GST-human Pak1-PBD on ice and immediately make 20μg aliquots. Store aliquots for later use at -70°C.
6. Add 20μg of GST-human Pak1-PBD to the spin cup containing the glutathione resin.
7. Immediately transfer up to 700μL of the cell lysate (containing at least 500μg of total proteins) to the spin cup, close the cap and vortex the sample.

8. Seal cap of the collection tube with laboratory film to prevent leakage, which may result from the presence of detergent in the lysate, and vortex the sample.
9. Incubate the reaction mixture at 4°C for 1 hour with gentle rocking.
10. Centrifuge the spin cup with collection tube at 6000 × g for 10-30 seconds.
11. Remove the laboratory film and transfer the spin cup to a new collection tube.
12. To wash resin, add 400µL of Lysis/Binding/Wash Buffer, invert the tube three times, and centrifuge at 6000 × g for 10-30 seconds. Decant the buffer. Repeat this wash step two additional times.
13. Transfer the spin cup to a new collection tube.
14. Prepare 50µL of reducing sample buffer for each pull-down reaction by mixing 1 part β-mercaptoethanol to 20 parts 2X SDS Sample Buffer (e.g., mix 2.5µL of β-mercaptoethanol to 50µL of 2X SDS Sample Buffer), or by adding dithiothreitol (DTT) to a final concentration of 200mM.
15. Add 50µL 2X reducing sample buffer to the resin. Vortex the sample and incubate at room temperature for 2 minutes.
16. Centrifuge the tube at 6000 × g for 2 minutes. Remove and discard the spin cup containing the resin.
17. Heat the eluted samples for 5 minutes at 95-100°C. Samples may be electrophoresed on a gel or stored at -20°C until use.
18. Apply at least 25µL per lane for a 10 × 10cm mini-gel (12% or 4-20% acrylamide gel provides the best separation).

D. Western Blot Analysis

Notes:

- This procedure has been optimized for use with SuperSignal West Pico Chemiluminescent Substrate (see Important Product Information section).
 - Include unfractionated cell lysate as a control to verify that the Western blot analysis is functioning properly.
 - Perform all blocking, probing and washing incubation steps using constant agitation.
1. Separate the proteins by SDS-PAGE and transfer to nitrocellulose or PVDF membrane.
 2. Block the membrane in TBS containing 3% BSA at room temperature for 1-2 hours.
 3. Rinse the membrane with TBS containing 0.05% Tween-20 Detergent (TBST) for 5 minutes.
 4. Prepare a solution containing the Anti-Rac1 mouse monoclonal antibody (1:1,000 dilution) in 3% BSA and 0.1% NaN₃ in TBST. An example of a 1:1,000 dilution is to add 10µL of the stock antibody solution to 10mL of buffer.
 5. Incubate the membrane in the anti-Rac1 antibody solution at 4°C overnight.
Note: If the number of pull-down reactions per blot is low, the diluted anti-Rac1 antibody solution can be re-used up to three times with no performance loss. Store the diluted anti-Rac1 antibody solution at 4°C for up to two months.
 6. Wash the membrane five times for 5 minutes each with TBST.
 7. Dilute the anti-mouse IgG-HRP-conjugate in TBST containing 5% nonfat dry milk [e.g., if using Pierce Goat Anti-Mouse IgG (H+L) Peroxidase Conjugated (Product No. 31430) dilute within a range from 1:20,000 to 1:100,000].
Note: Ensure the dry milk is completely dissolved in TBST (e.g., mix the milk in TBST on a stir-plate for 30 minutes at room temperature), otherwise the milk residuals can cause background on the Western blot.
 8. Incubate membrane in the anti-mouse IgG-HRP Conjugate solution at room temperature for 1 hour.
 9. Wash the membrane five times for 5 minutes each with TBST.
 10. Incubate the membrane with chemiluminescent substrate (e.g., SuperSignal West Pico Chemiluminescent Substrate).
 11. Immediately expose the membrane to X-ray film or a CCD camera.

Note: The Rac1 band is located at ~22 kDa.

Troubleshooting

Problem	Cause	Solution
No activated Rac1 detected	Primary antibody requires optimization	Optimize the primary antibody concentration
	Incorrect secondary antibody used for detection	Use goat anti-mouse IgG
	No activated Rac1 present in lysates	Include GTP γ S-treated lysate as positive control for pull-down
	Insufficient activated Rac1	Increase the amount of lysate used for detection
	GST-Pak1-PBD was not added	Add GST-Pak1-PBD to the reactions
	Degraded GST-Pak1-PBD	Avoid multiple freeze/thaw cycles of GST-Pak1-PBD
	Degraded proteins	Add protease inhibitors to the Lysis/Binding/Wash Buffer before lysing the cells
Detection system is not functioning properly or requires optimization	Consult the instructions for the detection system being used	
No signal with GTP γ S or strong signal with GDP	GTP γ S or GDP are no longer functional	Aliquot GTP γ S or GDP after the first thaw and store at -20°C; avoid repeated freeze/thaw cycles of the resuspended solution
	Incorrect concentration of EDTA or MgCl ₂	Prepare new solutions with correct concentration
Western blot resulted in high background	Inadequate blocking and/or washing	Consult the instructions for the detection system being used
	Secondary antibody concentration is too high	

Additional Information

Rac1 is active when bound to GTP and is inactive when bound to GDP. Active Rac1 binds specifically to the p21-binding domain (PBD) of p21-activated protein kinase 1 (Pak1), leading its activation. Therefore, the PBD of Pak1 can be used as a probe to specifically isolate the active form of Rac1 (Figure 1).




Figure 1. Western blot of control reactions.

NIH3T3 cell lysates (500 μ g) were treated *in vitro* with GTP γ S or GDP to activate or inactivate Rac1 (refer to optional step B). The lysates were then incubated with 20 μ g of GST-Pak1-PBD and a Glutathione Resin. GTP γ S-treated lysate was also incubated with GST alone in the presence of a Glutathione Resin (negative control). Half of the volume of the eluted samples (25 μ L) and 20 μ g of cell lysate were separated by 4-20% SDS-PAGE, transferred to a nitrocellulose membrane and probed with Anti-Rac1 Antibody. Pierce Goat Anti-Mouse IgG (H+L), Peroxidase Conjugated (Product No. 31430; 1:20,000 dilution) was used as the secondary antibody. The detection was performed with SuperSignal West Pico Chemiluminescent Substrate (Product No. 34080) and followed by exposure to X-ray film. The exposure time was one second.

Related Thermo Scientific Products

25200-44	Precise™ Protein Gels (see catalog or website for a complete listing)
21065	Pierce Background Eliminator Kit, for eliminating background from overexposed X-ray film
23236	Pierce Coomassie Plus (Bradford) Protein Assay Reagent
23227	BCA Protein Assay Reagent Kit
22660	Pierce 660nm Protein Assay Reagent, 750mL
28320	Surfact-Amps® 20 (Active Ingredient: Tween-20), 6 × 10mL
28379	BupH™ Tris Buffered Saline Packs, 10 packs, each makes 500mL
28358	Tris Buffered Saline, 20X, 500mL
78425	Halt™ Protease Inhibitor Single-Use Cocktail, EDTA-free (100X), 24 × 100µL microtubes
31430	Pierce Goat Anti-Mouse IgG (H+L), Peroxidase Conjugated, 2mL
34079	SuperSignal West Pico Chemiluminescent Substrate, 500mL
34090	CL-XPosure™ Film (5" × 7" sheets), 100 sheets/pkg
34091	CL-XPosure Film (8" × 10"), 100 sheets/pkg
20291	Dithiothreitol (DTT), No-Weigh™, 7.7mg DTT/Tube × 48 tubes
88014	Nitrocellulose Membrane, 0.45µm, 7.9cm × 10.5cm
88585	PVDF Membrane, 0.45µm, 7.9cm × 10.5cm
21059	Restore® Western Blot Stripping Buffer, 500mL

General References

- Benard, V. and Bokoch, G. (2002). Assay of Cdc42, Rac, and Rho GTPase activation by affinity methods. *Methods. Enzymol.* **345**:349-59.
- Burbelo, P.D., Drechsel, D. and Hall A. (1995). A conserved binding motif defines numerous candidate target proteins for both Cdc42 and Rac GTPases. *J. Biol. Chem.* **270**:29071.
- Benard, V., Bohl, B. and Bokoch, G. (1999). Characterization of Rac and Cdc42 activation in chemoattractant-stimulated human neutrophils using a novel assay for active GTPases. *J Biol Chem.* **274(19)**:13198-204.

Tween is a trademark of Croda International Plc.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.