## iBind<sup>™</sup> Fluorescent Detection (FD) Solution Kit

novex life technologies

Pub. Part No. 100022233 Pub. No. MAN0009825 Rev. D.0

(

Package Contents

Size Cat. No. SLF1019 10 rxns

- iBind<sup>™</sup> FD 5X Buffer
- iBind<sup>™</sup> 100X Additive • iBind™ FD 10% SDS

- - Storage Conditions
- Store buffer and additive at 4°C (do not freeze)
- Store SDS at room temperature



Required Materials

- iBind<sup>™</sup> Window Cover (or aluminum foil)
- iBind™ Device iBind<sup>™</sup> Card
- Blotting roller
- Alexa Fluor® 680 or IRDye® 680LT Goat Anti-Rabbit
- Alexa Fluor® 790 or IRDye® 800CW Goat Anti-Mouse



'imina

Preparation time: ~10 minutes Run time: Allow well 4 to empty; ~2.5 hours or longer



Product

Important

Guidelines

- This iBind<sup>™</sup> Fluorescent Detection (FD) Solution Kit facilitates primary and secondary antibody binding
- of membrane-bound proteins in western detection. ■ This kit is compatible with Alexa Fluor® and/or IRDye® conjugated secondary antibodies.
- Wear gloves when handling the iBind<sup>™</sup> card. Do not open the iBind<sup>™</sup> device after closing it over
- an iBind<sup>™</sup> card. Make sure iBind<sup>™</sup> FD 10% SDS is completely in solution before use. SDS is added to the secondary antibody solution to decrease background for both PVDF and nitrocellulose membranes.
- Protect the device from light with the iBind<sup>™</sup> window cover or aluminum foil.



Online

Visit our product page for additional information and protocols. For support, visit www.lifetechnologies.com/support.



For Research Use Only. Not for use in diagnostic procedures.

### iBind™ Fluorescent Detection (FD) Solution Kit

#### Prepare Antibodies

Dilute antibodies with 1X iBind™ FD Solution. Each blot requires 2 mL of each antibody solution.

Component	Primary Antibody (1° Ab) Solution
1X iBind™ FD Solution	2 mL
Primary Antibody	Final antibody concentration equal to the manufacturer's recommended dilution
Component	Secondary Antibody (2° Ab) Solution
1X iBind™ FD Solution	2 mL
iBind™ FD 10% SDS	10 μL
■ Alexa Fluor® 680 <b>OR</b>	• 1 μL (1:2000 dilution)
■ IRDye® 680LT	■ 0.5 μL (1:4000 dilution)
■ Alexa Fluor® 790 <b>OR</b>	• 1 μL (1:2000 dilution)
■ IRDye® 800CW	■ 0.67 µL (1:3000 dilution)

#### Limited Product Warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

#### Important Licensing Information

These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

#### Disclaimer

23 October 2014

LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT, TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. LI-COR, Odyssey, and IRDve are registered trademarks of LI-COR, Inc.



# iBind™ Fluorescent Detection (FD) Solution Kit

### iBind™ Fluorescent Detection Solution Kit Protocol

Follow the diagram below to prepare solutions and membrane and perform antibody binding with the iBind™ device.

Timeli	ine	Steps	Procedure Details			
			Standard 1X iBind™ FD Solution is recommended for most primary antibodies.  Use the Optional solution only if initial results give low sensitivity.			
1	Prepare 1X FD Solution	Component	1X iBind™ FD Solution (for 1 mini blot)			
			Standard	Optional		
			iBind™ FD 5X Buffer	6 mL	1.5 mL	
			iBind™ 100X Additive	75 μL	300 μL	
		Distilled water	23.9 mL	28.2 mL		
2	10	Prepare diluted antibodies	See "Prepare Antibodies" for details.			
3		Load iBind™ card and membrane in iBind™ device	<ol> <li>Immerse the membrane in 5 mL of 1X iBind™ FD Solution.</li> <li>Place the iBind™ card on the iBind™ device, making sure the card is aligned with the stoppers.</li> <li>Evenly apply 5 mL of 1X iBind™ FD Solution to the card so that the Flow Region is completely wet.</li> <li>Add 1 mL 1X iBind™ FD Solution to the center of the Membrane Region.</li> <li>Place the blot on the card, protein side down, with the low MW region toward the stack.</li> <li>Remove any air bubbles between the iBind™ card and the membrane with a blotting roller.</li> </ol>			
4	П		iBind™ Device Well Number	Solution	1	
	Add prepared antibody and wash solutions	1 2	2 mL diluted primary antibody 2 mL 1X iBind™ FD Solution			
		3	2 mL diluted secondary antibody			
		4	6 mL 1X iBind™ FD Solution			
5		Run	Close the iBind™ well cover to minimize evaporation. Cover the viewing area with the iBind™ window cover or aluminum foil.			
			Leave the device undisturbed until the last well (Well 4) is empty, ~2.5 hours or longer.			
6	1111	Detect blot	Rinse the membrane in distilled water.  Scan wet or dry membranes on the LI-COR® Odyssey® CLx scanner according to manufacturer recommendations.			