

FluoSpheres® Europium Luminescent Microspheres

Catalog no. F20880, F20881, F20883, F20884

Table 1. Contents and storage information.

Material	Amount	Concentration	Storage*	Stability
Standard carboxylate-modified FluoSpheres® europium luminescent microspheres	2 mL	0.5% solids in deionized water containing 2 mM sodium azide	• 2–6°C • Protect from light • DO NOT FREEZE	When stored as directed this product is stable for 1 year.
NeutrAvidin™-labeled FluoSpheres® europium Iuminescent microspheres	0.4 mL	0.5% solids in 50 mM sodium phosphate, 50 mM NaCl, pH 7.5, containing 5 mM sodium azide		

Approximate fluorescence excitation/emission maxima: FluoSpheres® microspheres: 340-370/610 nm.

Introduction

In biological specimens, autofluorescence is a common source of background fluorescence. One approach to increasing detectability is the use of time-resolved luminescence reagents. FluoSpheres[®] europium luminescent microspheres incorporate Eu³⁺ in an organic coordination complex. This unique dye confers luminescence with a decay time of >500 µs, far longer than that of conventional fluorescent probes or autofluorescent samples, typically having decay times of <50 ns. Thus, time-resolved fluorometry can virtually eliminate autofluorescence. 1-6 In addition, these europium luminescent microspheres feature longwavelength emission (~610 nm) that is well separated from the excitation peak (~365 nm) (Figure 1). This unusually large Stokes shift permits the use of filter combinations that effectively isolate the desired luminescence signal. These FluoSpheres® europium luminescent microspheres, with nominal diameters of 0.04 μm or 0.2 μm, are available uncoated or conjugated to NeutrAvidin™ biotin-binding protein (Table 2). The NeutrAvidin™ proteinlabeled beads can be used for the indirect detection of antigens and DNA targets in many biotin/avidin reaction-based assays.

Guidelines for Use

Before sampling, mix well by sonication, shaking, or vortexing. Aliquot the NeutrAvidin™labeled microspheres before sonication.

Note: Avoid excessive vortexing/sonication of the beads as it may result in protein damage.

For all applications involving protein-labeled microspheres, such as the NeutrAvidin[™]-labeled




Figure 1. Fluorescence excitation and emission spectra of the FluoSpheres® europium luminescent microspheres.

Table 2. Summary of FluoSpheres® europium luminescent microspheres.

Missosubovo tuno (Ev/Em)*	Microsphere diameter†			
Microsphere type (Ex/Em)*	0.04 μm	0.2 μm		
Standard (365/610)	F20880	F20881		
NeutrAvidin™-labeled (365/610)	F20883	F20884		

^{*}The microspheres are available as standard "carboxylate-modified" microspheres or as microspheres conjugated with NeutrAvidin™ biotin-binding protein. Approximate excitation (Ex) and emission (Em) maxima, in nm, are indicated in parentheses.

FluoSpheres® europium luminescent microspheres, we recommend the use of BlockAid™ blocking solution (Cat. no. B10710) to reduce nonspecific binding.

References

1. J Membrane Biol 19, 1 (1974); 2. 1. J Photochem Photobiol B 27, 3 (1995); 2. Clin Chem 43, 1937 (1997); 3. Biophys J 74, 2210 (1998); 4. J Histochem Cytochem 44, 1091 (1996); 5. Histochem J 31, 45 (1999); 6. Biochemistry 37, 2372 (1998).

 $^{^{\}dagger}$ The sizes listed are nominal diameters; actual sizes are printed on the product labels.

Product List Current prices may be obtained from our website or from our Customer Service Department.

Cat. no.	Product Name	Unit Size
F20880	FluoSpheres® carboxylate-modified microspheres, 0.04 µm, europium luminescent (365/610) *0.5% solids*	2 mL
F20881	FluoSpheres® carboxylate-modified microspheres, 0.2 µm, europium luminescent (365/610) *0.5% solids*	2 mL
F20883	FluoSpheres® NeutrAvidin™ labeled microspheres, 0.04 µm, europium luminescent (365/610) *0.5% solids*	0.4 mL
F20884	FluoSpheres® NeutrAvidin™ labeled microspheres, 0.2 μm, europium luminescent (365/610) *0.5% solids*	0.4 mL
Related prod	luct	
B10710	BlockAid™ blocking solution *for use with microspheres*	50 mL

Contact Information

Molecular Probes, Inc.

29851 Willow Creek Road Eugene, OR 97402 Phone: (541) 465-8300 Fax: (541) 335-0504

Customer Service:

6:00 am to 4:30 pm (Pacific Time) Phone: (541) 335-0338 Fax: (541) 335-0305 probesorder@invitrogen.com

Toll-Free Ordering for USA:

Order Phone: (800) 438-2209 Order Fax: (800) 438-0228

Technical Service:

8:00 am to 4:00 pm (Pacific Time) Phone: (541) 335-0353 Toll-Free (800) 438-2209 Fax: (541) 335-0238 probestech@invitrogen.com

Invitrogen European Headquarters

Invitrogen, Ltd. 3 Fountain Drive Inchinnan Business Park Paisley PA4 9RF, UK Phone: +44 (0) 141 814 6100 Fax: +44 (0) 141 814 6260 Email: euroinfo@invitrogen.com Technical Services: eurotech@invitrogen.com

For country-specific contact information. visit www.invitrogen.com.

Further information on Molecular Probes products, including product bibliographies, is available from your local distributor or directly an extension of the contract of the product offrom Molecular Probes. Customers in Europe, Africa and the Middle East should contact our office in Paisley, United Kingdom. All others should contact our Technical Service Department in Eugene, Oregon.

Molecular Probes products are high-quality reagents and materials intended for research purposes only. These products must be used by, or directly under the supervision of, a technically qualified individual experienced in handling potentially hazardous chemicals. Please read the Material Safety Data Sheet provided for each product; other regulatory considerations may apply.

Limited Use Label License No. 223: Labeling and Detection Technology

The purchase of this product conveys to the buyer the non-transferable right to use the purchased amount of the product and composition of the product and composition of the product and composition of the product conveys to the buyer the non-transferable right to use the purchased amount of the product and composition of the product conveys to the buyer the non-transferable right to use the purchased amount of the product and composition of the product conveys to the buyer the non-transferable right to use the purchased amount of the product and composition of the product conveys to the purchased amount of thnents of the product in research conducted by the buyer (whether the buyer is an academic or for-profit entity). The buyer cannot sell or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (a) this product (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (b) its components to a third party or otherwise transfer (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (b) its components or (c) materials made using this product or its components to a third party or otherwise transfer (b) its components or (c) materials made using this product or its components or (c) materials made using the component of the c erwise use this product or its components or materials made using this product or its components for Commercial Purposes. The buyer may transfer information or materials made through the use of this product to a scientific collaborator, provided that such transfer is not for any Commercial Purpose, and that such collaborator agrees in writing (a) to not transfer such materials to any third party, and (b) to use such transferred materials and/or information solely for research and not for Commercial Purposes. Commercial Purposes means any activity by a party for consideration and may include, but is not limited to: (1) use of the product or its components in manufacturing; (2) use of the product or its components to provide a service, information, or data; (3) use of the product or its components for therapeutic, diagnostic or prophylactic purposes; or (4) resale of the product or its components, whether or not such product or its components are resold for use in research. Invitrogen Corporation will not assert a claim against the buyer of infringement of the above patents based upon the manufacture, use or sale of a therapeutic, clinical diagnostic, vaccine or prophylactic product developed in research by the buyer in which this product or its components was employed, provided that neither this product nor any of its components was used in the manufacture of such product. If the purchaser is not willing to accept the limitations of this limited use statement, Invitrogen is willing to accept return of the product with a full refund. For information on purchasing a license to this product for purposes other than research, contact Molecular Probes, Inc., Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

Several Molecular Probes products and product applications are covered by U.S. and foreign patents and patents pending. All names containing the designation ® are registered with the U.S. Patent and Trademark Office.

Copyright 2009, Molecular Probes, Inc. All rights reserved. This information is subject to change without notice.