

Pierce™ High Capacity Streptavidin Chromatography Cartridge

87739 87740

2165.1

Number	Description
87739	Pierce High Capacity Streptavidin Chromatography Cartridge, 2 × 1mL
87740	Pierce High Capacity Streptavidin Chromatography Cartridge, 1 × 5mL Support: 6% highly crosslinked fast-flow agarose supplied in 0.02% sodium azide Binding Capacity: ≥ 10mg biotinylated BSA/ml resin; ≥ 100µg biotin para-nitrophenyl ester/ml resin Storage: Upon receipt store at 4-8°C. Product shipped at ambient temperature. Do not freeze.

Introduction

The Thermo Scientific™ Pierce™ High Capacity Streptavidin Chromatography Cartridges are convenient, ready-to-use devices for separating biotinylated molecules from non-biotinylated molecules and for purifying antigens using biotinylated antibodies. This resin has 2-3 times higher binding capacity than the standard streptavidin agarose resin and is ideal when a high amount of recovery is required. The cartridges are effective for assay development and immunoprecipitation¹⁻³ and for physically separating DNA strands produced in a polymerase chain reaction by incorporating biotin in one of the amplification polymers.⁴

These cartridges enable fast, easy and reproducible chromatographic separations and are compatible with the major automated liquid-chromatography systems or manual syringe processing (see Table 1 for cartridge general properties). The cartridges attach directly to ÄKTA™ or FPLC Systems without additional connectors. The included accessory pack readily adapts cartridges for use with Luer-Lok™ Syringe Fittings or 1/16" tubing. These cartridges enable fast, easy and reproducible chromatographic separations and can be regenerated multiple times for affinity purification.

Table 1. Properties of the Thermo Scientific Pierce High Capacity Streptavidin Chromatography Cartridges.

Particle Size	45-165µm
Cartridge Dimensions	0.7 × 2.7cm (1mL column) 1.3 × 3.8cm (5mL column)
Void Volume	0.32mL (1mL column) 1.5mL (5mL column)
Recommended Flow Rate	0.1-1mL/min (1mL column); 0.5-2mL/min (5mL column)
Maximum Recommended Flow Rate*	4mL/min (1mL column); 5mL/min (5mL column)
pH Limits	2-11 (short-term); 4-9 (long-term)
Maximum Operating Pressure	0.3 MPa, 43.5 psi or 3 bar
Cartridge Material	Polypropylene
Frit	Polyethylene, 10µm
Accessory Pack	Luer-Lok Adapter to 10-32 male Finger-tight 10-32 connector fitting for 1/16" OD tubing Column plug for 10-32 coned port Cap 1/16 male

*At room temperature using water

Important Product Information

- The example protocols included in these instructions describe applications using liquid chromatography (LC) systems using 1/16" tubing or syringes. For best results, process all buffers through a 0.45µm filter before use in LC systems. Additional connectors and fittings are required to attach these cartridges to the Bio-Rad BioLogic™ System. Specific applications and systems will require optimization.
- For optimum binding of the biotinylated molecule, apply sample at the indicated low flow rates.
- All purification procedures listed can be performed using a syringe fitted with the Luer-Lok Adapter included in the accessory pack.
- To elute biotinylated molecules from the streptavidin cartridges, use 8M guanidine•HCl, pH 1.5 (Product No. 24115). For non-denaturing elution conditions, biotinylate the protein using NHS-Iminobiotin (Product No. 21117), which binds to streptavidin at pH 9.5 and dissociates at pH 4. Alternatively, use a thiol-cleavable biotinylation reagent such as NHS-SS-biotin (Product No. 21331) or NHS-SS-PEG₄-biotin (Product No. 21442).
- Guanidine•HCl is a harsh denaturant and may irreversibly damage the protein of interest and result in leaching of streptavidin subunits. The cartridge cannot be reused after eluting with guanidine•HCl.
- The cartridge was not designed to be refilled.

Procedure for Purifying Biotinylated Molecules

A. Additional Materials Required

- Binding Buffer: Phosphate-buffered saline (e.g., 0.1M phosphate, 0.15M sodium chloride; pH 7.2; Product No. 28372)
- Elution Buffer: 8M guanidine•HCl, pH 1.5 (Product No. 24115)

B. Procedure

Note: For syringe application, 30 drops per minute is equivalent to a flow rate of 1mL per minute.

1. Use approximately 10mg of biotinylated protein/ml of streptavidin resin. Adjust the sample to the ionic strength and pH of the Binding Buffer by dilution or buffer exchange before applying to the cartridge. If the sample contains insoluble matter, centrifuge or filter (0.45µm filter) before use.
2. Equilibrate the cartridge and all buffers to working temperature. Perform purifications at room temperature or 4°C. Ensure that all solutions are degassed.
3. Prepare the LC system by filling tubing with buffer. Remove top plug from cartridge and carefully snap off the end-tab (do not twist). To avoid introducing air into the system, let a few drops of buffer flow from tubing into the cartridge then connect cartridge top to the tubing; allow a few drops to emerge from the cartridge before connecting to the LC inlet port.
4. Equilibrate the cartridge with 10 column volumes of Binding Buffer at a flow rate of 0.1-1mL/minute for the 1mL cartridge or 0.5-2mL/minute for the 5mL cartridge.
5. Apply biotinylated sample to the cartridge. Use a 0.1-0.2mL/minute flow rate for a 1mL cartridge or 0.5-1mL/minute for a 5mL cartridge. Apply any volume that does not exceed cartridge capacity. Collect fractions.
6. Wash with at least 10 column volumes of Binding Buffer or until the absorbance of the flow-through approaches baseline. Use 0.1-1mL/minute or 0.5-2mL/minute flow rate for washing 1mL or 5mL columns, respectively.
7. Elute the bound biotinylated sample with approximately 10-20 column volumes of Elution Buffer and collect fractions. Immediately dialyze or buffer exchange the fractions of interest.
8. Monitor protein elution by measuring the absorbance of the fractions at 280nm or by performing a protein assay (e.g., Thermo Scientific BCA Protein Assay, Product No. 23227).

Procedure for Affinity Purification of Antigens

A. Additional Materials Required

- Biotinylated antibody: Prepare biotinylated antibody at 0.2-10mg/mL in Binding Buffer and use approximately 10 mg of biotinylated antibody/ml of streptavidin resin
- Antigen-containing sample
- Binding Buffer: Phosphate-buffered saline (e.g., 0.1M phosphate, 0.15M NaCl; pH 7.2; Product No. 28372)
- Elution Buffer: IgG Elution Buffer (Product No. 21004), Gentle Ag/Ab Elution Buffer (Product No. 21027) or 0.1M glycine•HCl, pH 2.8

B. Procedure

Note: For syringe application, 30 drops per minute is equivalent to a flow rate of 1mL per minute.

1. Equilibrate the cartridge and all buffers to working temperature. Perform purifications at room temperature or 4°C. Ensure that all solutions are degassed.
2. Prepare the LC system by filling tubing with buffer. Remove top plug from cartridge and carefully snap off the end-tab (do not twist). To avoid introducing air into the system, let a few drops of buffer flow from the tubing into the cartridge and then connect cartridge top to the tubing; allow a few drops to emerge from the cartridge before connecting to the LC inlet port.
3. Equilibrate the cartridge with 10 column volumes of Binding Buffer at a flow rate of 0.1-1 mL/minute for the 1mL cartridge or 0.5-2mL/minute for the 5mL cartridge.
4. Apply biotinylated antibody solution to the cartridge. Use a 0.1mL/minute flow rate for a 1mL cartridge or 0.5mL/minute for a 5mL cartridge. Alternatively or for syringe application allow one column volume of solution to enter the cartridge. Pause the flow and incubate for 10-15 minutes. If the entire solution volume cannot be added at once, repeat with more antibody solution and incubate. Apply any volume that does not exceed cartridge capacity.
5. Wash with at least 10 column volumes of Binding Buffer or until the flow-through absorbance approaches baseline. Use a flow rate of 0.1-1mL/minute for the 1mL cartridge or 0.5-2mL/minute for the 5mL cartridge.
6. Apply antigen sample to the cartridge at 0.1mL/minute flow rate for a 1mL cartridge or 0.5mL/minute for a 5mL cartridge. Alternatively or for syringe application allow the sample to enter the cartridge, pause the flow and incubate at room temperature for 30 minutes or overnight at 4°C.
7. Wash with at least 10 column volumes of Binding Buffer or until the flow-through absorbance approaches baseline. Use a flow rate of 0.1-1mL/minute for the 1mL cartridge or 0.5-2mL/minute for the 5mL cartridge.

Note: If using Gentle Ag/Ab Elution Buffer, wash cartridge with three column volumes of Tris-buffered saline before antigen elution. The Gentle Elution Buffer is not compatible with phosphate-based buffers.

8. Elute the bound antigen with approximately 10-20 column volumes of Elution Buffer and collect fractions.
9. If using IgG Elution Buffer or 0.1M glycine•HCl, pH 2.8, immediately adjust the pH to neutral by adding 100µL of 1M Tris; pH 7.5-9.0 per 1mL of eluate collected or by dialysis or buffer exchange.
10. Monitor protein elution by measuring the absorbance of the fractions at 280nm or by using the BCA Protein Assay. The eluted protein can be directly analyzed by SDS-PAGE.
11. Re-equilibrate the immobilized antibody cartridge with at least 10 column volumes of Binding Buffer until the flow-through absorbance returns to baseline. The ability to reuse the cartridge is dependent upon each specific immobilized antibody and the sample. For best results reuse cartridge with identical samples to minimize contamination. Store cartridge at 4°C in 0.02% sodium azide.

Procedure for Purifying Iminobiotinylated Molecules

A. Additional Materials Required

- NHS-Iminobiotin-labeled sample adjusted to the ionic strength and pH of the Binding Buffer
- Binding Buffer: 50mM ammonium carbonate, 0.5M NaCl, pH 10.0
- Elution Buffer: 50mM ammonium acetate, 0.5M NaCl, pH 4.0

B. Procedure

Note: For syringe application, 30 drops per minute is equivalent to a flow rate of 1mL per minute.

1. Equilibrate the cartridge and all buffers to working temperature. Perform purifications at room temperature or 4°C. Ensure that all solutions are degassed.
2. Prepare the LC system by filling tubing with buffer. Remove top plug from cartridge and carefully snap off the end-tab (do not twist). To avoid introducing air into the system, let a few drops of buffer flow from tubing into cartridge top then connect cartridge top to the tubing; allow a few drops to emerge from the cartridge before connecting to the LC inlet port.
3. Equilibrate the cartridge with 10 column volumes of Binding Buffer at a flow rate of 0.1-1mL/minute for the 1mL cartridge or 0.5-2mL/minute for the 5mL cartridge.
4. Apply iminobiotinylated sample to the cartridge at 0.1-0.2mL/minute flow rate for a 1mL cartridge or 0.5-1mL/minute for a 5mL cartridge. Apply any volume that does not exceed cartridge capacity.
5. Wash with at least 10 column volumes of Binding Buffer or until the absorbance of the flow-through approaches baseline. Use 0.1-1mL/minute or 0.5-2mL/minute flow rate for washing 1mL or 5mL columns, respectively.
6. Elute with approximately 10-20 column volumes of Elution Buffer and collect fractions.
7. Immediately adjust the pH of the eluate to neutral by buffer exchange.
8. Monitor protein elution by measuring the absorbance of the fractions at 280nm or using the BCA Protein Assay. The eluted protein can be directly analyzed by SDS-PAGE.
9. Typically the cartridge can be reused multiple times. For best results, reuse cartridge with identical samples to minimize contamination. Store cartridge at 4°C in 0.02% sodium azide. For reuse, equilibrate the cartridge with at least 10 column volumes of Binding Buffer until the flow-through absorbance returns to baseline.

Related Thermo Scientific Products

87741	Pierce High Capacity NeutrAvidin Chromatography Cartridges, 2 × 1mL
87742	Pierce High Capacity NeutrAvidin Chromatography Cartridges, 2 × 5mL
21117	NHS-Iminobiotin Trifluoroacetamide
21435	EZ-Link™ Sulfo-NHS-LC-Biotinylation Kit
21440	EZ-Link NHS-PEG Solid Phase Biotinylation Kit: <i>pre-packed column</i>
20227	Pierce Monomeric Avidin Kit
66425	Slide-A-Lyzer™ Dialysis Cassette, 10K MWCO, 3mL, 10 pack
89934	Pierce Desalting Chromatography Cartridges, 1 ml × 5 columns
89935	Pierce Desalting Chromatography Cartridges, 5 ml × 5 columns
89891	Zeba™ Spin Desalting Columns, 7K MWCO, 5mL, 5/pkg
89893	Zeba Spin Desalting Columns, 7K MWCO, 10mL, 5/pkg
89894	Zeba Spin Desalting Columns, 7K MWCO, 10mL, 25/pkg

Visit Our Website for Additional Information

- Tech Tip #27: Optimize elution conditions for immunoaffinity purification
- Tech Tip #29: Degas solutions for use in affinity and gel filtration columns

Cited References

1. Bruch, R.C. and White, H.B., III. (1982). Compositional and structural heterogeneity of avidin glycopeptides. *Biochemistry* **21**:5334-41.
2. Swack, J.A., *et al.* (1978). Use of avidin-sepharose to isolate and identify biotin polypeptides from crude extracts. *Anal Biochem* **87**:14-26.
3. Wilchek, M. and Bayer, E.A. (1989). A universal affinity column using avidin-biotin technology. In Protein Recognition of Immobilized Ligands. (Hutchins, T.W., ed.). Alan R. Liss, Inc. 83-90.
4. Pellegrini, M., *et al.* (1977). Application of the avidin-biotin method of gene enrichment to the isolation of long double-stranded DNA containing specific gene sequences. *Nucl Acids Res* **4**:2961-73.

General References

- Chaiet, L. and Wolf, F.J. (1964). The properties of streptavidin, a biotin-binding protein produced by streptomycetes. *Arch Biochem Biophys* **106**:1-5.
- Gitlin, G., *et al.* (1987). Studies of the biotin-binding site of avidin. *Biochem J* **242**:923-6.
- Green, N.M. (1975). Avidin. In Adv. in Protein Chemistry, Academic Press, New York, **29**:85-133.
- Manning, J., *et al.* (1977). A method for gene enrichment based on the avidin-biotin interaction. Application to the *Drosophila* ribosomal RNA genes. *Biochemistry* **16**:1364-70.
- Updyke, T.V. and Nicolson, G.L. (1984). Immunoaffinity isolation of membrane antigens with biotinylated monoclonal antibodies and immobilized streptavidin matrices. *J Immunol Meth* **73**:83-95.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2013 Thermo Fisher Scientific Inc. All rights reserved. ÄKTA is a trademark of GE Healthcare. Bio-Logic is a trademark of Bio-Logic SAS. Luer-Lok is a trademark of Becton Dickinson, Inc. All (other) trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.