INSTRUCTIONS

HABA 4´-hydroxyazobenzene-2-carboxylic acid

28010

Number

28010

Description

HABA, 10g CAS # 1634-82-8

Molecular Weight: 242.23



Storage: Upon receipt store product at room temperature. Product shipped at ambient temperature.

Introduction

Antibodies or other proteins are often biotinylated with reagents such as Thermo Scientific EZ-Link Sulfo-NHS-LC-Biotin (Product No. 21335). Thermo Scientific HABA (4'-hydroxyazobenzene-2-carboxylic acid) is a reagent that enables a quick estimation of the mole-to-mole ratio of biotin to protein. To quantify biotin label incorporation, a solution containing the biotinylated protein is added to a mixture of HABA and avidin. Because of its higher affinity for avidin, biotin displaces the HABA from its interaction with avidin and the absorption at 500nm decreases proportionately. By this method, an unknown amount of biotin present in a solution can be evaluated in a single cuvette by measuring the absorbance of the HABA-avidin solution before and after addition of the biotin-containing sample. The change in absorbance relates to the amount of biotin in the sample.

Important Product Information

- Following any biotin-labeling reaction, the biotinylated protein sample to be assayed must be dialyzed or desalted to remove nonreacted and hydrolyzed biotinylation reagent.
- Samples must be in one of the recommended buffers (PBS or TBS, see Reagent Preparation Section) for the assay. Avoid buffers containing potassium (such as Modified Dulbecco's PBS), which will cause precipitation in the assay. Other buffers may interfere and should not be used unless first validated by comparison to results using PBS or TBS.
- Slight color variation between the HABA/Avidin Premix microtubes does not affect product performance.

Procedure for Measuring the Level of Biotin Incorporation

Note: The biotin-labeled protein sample must be desalted or dialyzed to remove all traces of nonreacted and hydrolyzed biotinylation reagent before performing the HABA assay. (This assay detects total biotin in the solution.)

A. Additional Materials Required

- Avidin (Product No. 21121)
- 1N NaOH
- Phosphate Buffered Saline (PBS): 100mM sodium phosphate, 150mM sodium chloride; pH 7.2 (Product No. 28372)
- Cuvettes and spectrophotometer, or 96-well microplate and plate reader capable of measuring at 500nm.

B. Reagent Preparation

HABA Solution	Add 24.2mg HABA to 9.9mL of ultrapure water and then add 0.1mL of 1N NaOH. If the HABA does not completely dissolve, add another 0.1mL of 1N NaOH and filter solution before use. Store HABA solution at 4°C.
HABA/Avidin Solution	Add 10mg of avidin and 600 μL of the HABA Solution to 19.4mL of PBS. Store solution at 4°C for up to 2 weeks.

Thermo

0212.6



C. Procedure (Cuvette Format)

- 1. Pipette 900 μ L of HABA/Avidin Solution into a 1mL cuvette and measure the absorbance at 500nm. The A₅₀₀ of this solution should be 0.9-1.3. Record the value as A₅₀₀ HABA/Avidin.
- 2. Add 100µL of biotinylated protein sample to the cuvette containing HABA/Avidin and mix well.
- 3. Measure the absorbance of the solution at 500nm. Once the absorbance value remains constant for at least 15 seconds, record the value as A_{500} HABA/Avidin/Biotin Sample. If the A_{500} HABA/Avidin/Biotin Sample is ≤ 0.3 , dilute sample and repeat assay. Proceed to Calculations (Section E).

D. Procedure (Microplate Format)

- 1. Add 180µL of the HABA/Avidin Solution to a microplate well. Thoroughly mix plate on orbital shaker or equivalent.
- 2. Measure the absorbance of the solution in the well at 500nm and record the value as A_{500} HABA/avidin.
- 3. Add 20µL of biotinylated sample to the well containing the HABA/Avidin Solution. Mix as described above.
- Measure the absorbance of the solution in the well at 500nm and record the value as A₅₀₀ HABA/avidin/biotin sample once the value remains constant for at least 15 seconds. If the A₅₀₀ HABA/Avidin/Biotin Sample is ≤ 0.15, dilute sample and repeat assay.
- 5. Proceed to the Calculations (Section E).

E. Calculations for Moles of Biotin Per Mole of Protein

Note: The HABA Calculator, which is available from the Technical Resources menu from the website, will calculate the moles of biotin per mole of protein upon entering the required values.

These calculations are based on the Beer Lambert Law (Beer's Law): $A_{\lambda} = \varepsilon_{\lambda} bC$

Where:

A is the absorbance of the sample at a particular wavelength (λ). The wavelength for the HABA assay is 500nm. There are no units for absorbance.

ε is the absorptivity or extinction coefficient at the wavelength (λ). For HABA/avidin samples at 500nm, pH 7.0 extinction coefficient is equal to 34,000M⁻¹cm⁻¹.

b is the cell path length expressed in centimeters (cm). A 10mm square cuvette has a path length of 1.0cm. In the microplate format with the recommended volumes, the path length is typically 0.5cm.

C is the concentration of the sample expressed in molarity (= mol/L = mmol/mL).

The values needed for calculating the number of moles of biotin per mole of protein or sample are as follows:

- Concentration of the protein or sample used, expressed as mg/ml
- Molecular weight (MW) of the protein, expressed as grams per mole (e.g., HRP = 40,000; IgG = 150,000)
- Absorbance at 500nm for HABA/avidin reaction mixture $(A_{500} H | A)$
- Absorbance at 500nm for HABA/avidin/biotin reaction mixture (A_{500} H\A\B)
- Dilution factor, if the sample is diluted before adding to the HABA/avidin reaction mixture
- 1. Calculation #1 is for the concentration of biotinylated protein in mmol/mL (before any dilution for the assay procedure):

mmol protein per mL =
$$\frac{\text{protein concentration (mg/mL)}}{\text{MW of protein (mg/mmol)}} = Calc#1$$

- 2. Calculation #2 is for the change in absorbance at 500nm:
 - Cuvette:

```
\Delta A_{500} = (0.9 \times A_{500} \text{ H} \land \text{A}) - (A_{500} \text{ H} \land \text{A}) = Calc#2
```

Note: The cuvette format requires the 0.9 correction factor to adjust for dilution of the H\A mixture by the biotinylated protein sample. The microplate format does not require this correction factor because the



•	Microplate:
---	-------------

 $\Delta A_{500} = (A_{500} \text{ H} \land A) - (A_{500} \text{ H} \land A \land B) = Calc #2$

dilution effect is exactly offset by the increased height and light path length of solution in the well.

www.thermoscientific.com/pierce

3. Calculation #3 is for the concentration of biotin in mmol per ml of reaction mixture:

	mmol biotin mL reaction mixture	$=\frac{\Delta A500}{(34,000\times b)}=$	$=\frac{Calc#2}{(34,000\times b)}$	= <i>Calc</i> #3	Note: <i>b</i> is the light path length (cm) of the sample. Use $b = 1$ with the cuvette format. Use $b = 0.5$ with the microplate format when using a standard 96-well plate and the volumes specified in the procedure. The exact path length is the height of the solution through which the plate reader measures the absorbance.
4. Calculation #4 is for the mmol of biotin per mmol of protein:					
	$=\frac{\text{mmol biotin in original sample}}{\text{mmol protein in original sample}}$		Note: The original biotinylated protein sample was diluted 10-fold in the reaction mixture. Therefore, a multiplier of 10 is used in this step to convert the		
	(mmol per mL biotin	in in reaction mixture)(10)(dilution factor)			biotin concentration in the reaction mixture to the biotin concentration in the original sample. If the
	=		original sample was diluted before performing the assay, then the dilution factor must be used as well. Calculation #4 yields the biotin:protein molar ratio (average # of biotin molecules per protein molecule).		
	$=\frac{(Calc#3)\times10\times010}{Calc#1}$				

EXAMPLE: In this example, the labeled protein is IgG (MW 150,000) at 0.69 mg/mL. The absorbance measurements were A_{500} H\A = 0.904 and A_{500} H\A = 0.771.

- 1. mmol biotinylated protein per mL = $\frac{0.69 \text{ mg/mL}}{150,000 \text{ mg/mmol}} = 4.6 \times 10^{-6}$
- 2. $\Delta A_{500} = (0.9 \times 0.904) 0.771 = 0.0426$

3.
$$\frac{\text{mmol biotin}}{\text{mL reaction mixture}} = \frac{0.0426}{34,000} = 1.25 \times 10^{-60}$$

4. $\frac{\text{mmol biotin}}{\text{mmol protein}} = \frac{(1.25 \times 10^{-6}) \times 10}{4.6 \times 10^{-6}} = \frac{12.5 \times 10^{-6}}{4.6 \times 10^{-6}} = 2.72 \text{ biotin molecules per IgG molecule}$

Troubleshooting

Problem	Possible Cause	Solution
$\Delta A_{500} \text{ in HABA}$ assay is ≤ 0	The protein sample had no or a low level of biotinylation because of limited accessible functional groups on the protein	Repeat biotinylation with alternative chemistry (e.g., amine reactive rather than sulfhydryl reactive) or use a higher molar ratio of biotinylation reagent
	Incomplete mixing of reagent	Completely solubilize and mix HABA/Avidin before diluting
	Particulate in sample contributed to absorbance	Filter protein sample to remove particulate
High levels of biotinylation	Nonreacted biotin was not removed	Dialyze or desalt sample before performing the assay



Related Thermo Scientific Products

46610	Fluorescence Biotin Quantitation Kit
21329	EZ-Link [®] NHS-PEG₄-Biotin, No-Weigh [™] Format, 8 × 2 mg microtubes
21335	EZ-Link Sulfo-NHS-LC-Biotin, 100mg
69576	Slide-A-Lyzer [®] MINI Dialysis Device Kit, 10K MWCO, 0.1mL
66382	Slide-A-Lyzer Dialysis Cassette Kit, 10K MWCO, 3mL
20347	Streptavidin Agarose Resin, 2mL
21126	Streptavidin, Horseradish Peroxidase Conjugated, 1mg
28005	Pierce [®] Biotin Quantitation Kit

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at <u>www.thermoscientific.com/pierce</u>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.