

## WT Expression Kits

Catalog Numbers 4411973, 4411974, 4440536, and 4440537

Pub. No. 4425210 Rev. C

**Note:** For safety and biohazard guidelines, see the “Safety” appendix in the *WT Expression Kit User Guide* (Pub. No. 4425209). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

### Before you begin

1. Prepare the Control RNA.
2. Prepare your total RNA according to your laboratory’s procedure.
3. Determine your input RNA quantity.
4. Prepare the Poly-A RNA Controls.
5. Evaluate RNA quality by determining its  $A_{260}/A_{280}$  ratio. RNA of acceptable quality is in the range 1.7–2.1.
6. Evaluate RNA integrity by microfluidic analysis or denaturing agarose gel electrophoresis.
7. Program your thermal cycler.

Method	Heated lid temp.	Alternate Protocol <sup>[1]</sup>	Step 1	Step 2	Step 3	Step 4
First-Strand cDNA Synthesis	50°C	105°C	25°C, 60 minutes	42°C, 60 minutes	4°C, 2 minutes	—
Second-Strand cDNA Synthesis	RT or disable	Lid open	16°C, 60 minutes	65°C, 10 minutes	4°C, 2 minutes	—
In Vitro Transcription cRNA Synthesis	50°C	40°C oven	40°C, 16 hours	4°C, Hold	—	—
2nd-Cycle cRNA Denaturation	75°C	105°C	70°C, 5 minutes	25°C, 5 minutes	4°C, 2 minutes	—
2nd-Cycle cDNA Synthesis	75°C	105°C	25°C, 10 minutes	42°C, 90 minutes	70°C, 10 minutes	4°C, 2 minutes
RNase H Hydrolysis	75°C	105°C	37°C, 45 minutes	95°C, 5 minutes	4°C, 2 minutes	—

<sup>[1]</sup> For thermal cyclers that lack a programmable heated lid.

## Day 1 workflow

### 1 Synthesize first-strand cDNA

- a. At room temperature, prepare the First-Strand Master Mix in a nuclease-free tube.

First-Strand Master Mix component	Volume for one reaction
First-Strand Buffer Mix	4 $\mu$ L
First-Strand Enzyme Mix	1 $\mu$ L
Total volume	5 $\mu$ L

- b. Mix thoroughly by gently vortexing. Centrifuge briefly. Proceed immediately to the next step.  
c. Transfer 5  $\mu$ L of the First-Strand Master Mix to the supplied PCR tubes.  
d. Add 5  $\mu$ L of total RNA, mix thoroughly by gently vortexing, centrifuge briefly, then proceed immediately to the next step.  
e. Incubate for 1 hour at 25°C, then for 1 hour at 42°C, then for at least 2 minutes at 4°C in a thermal cycler. Immediately after the incubation, centrifuge briefly (~5 seconds), then place the sample on ice for 2 minutes.

### 2 Synthesize second-strand cDNA

- a. On ice, prepare the Second-Strand Master Mix in a nuclease-free tube.

Second-Strand Master Mix component	Volume for one reaction
Nuclease-free Water	32.5 $\mu$ L
Second-Strand Buffer Mix	12.5 $\mu$ L
Second-Strand Enzyme Mix	5 $\mu$ L
Total volume	50 $\mu$ L

- b. Mix thoroughly by gently vortexing. Centrifuge briefly and proceed immediately to the next step.  
c. Transfer 50  $\mu$ L of the Second-Strand Master Mix to each (10  $\mu$ L) first-strand synthesis cDNA sample. Mix thoroughly by gently vortexing or flicking the tube 3 or 4 times. Centrifuge briefly and proceed immediately to the next step.  
d. Incubate for 1 hour at 16°C, for 10 minutes at 65°C, then for at least 2 minutes at 4°C in a thermal cycler.

**IMPORTANT!** Disable the heated lid of the thermal cycler or keep the lid off during the second-strand cDNA synthesis.

- e. Immediately after the incubation, centrifuge briefly (~5 seconds). Place the sample on ice.

### 3 Synthesize cRNA by in vitro transcription

- a. At room temperature, prepare an IVT Master Mix in a nuclease-free tube.

IVT Master Mix component	Volume for one reaction
IVT Buffer Mix	24 $\mu$ L
IVT Enzyme Mix	6 $\mu$ L
Total volume	30 $\mu$ L

- b. Mix thoroughly by gently vortexing. Centrifuge briefly, proceed immediately to the next step.  
c. Transfer 30  $\mu$ L of the IVT Master Mix to each 60- $\mu$ L Second-Strand cDNA sample. Mix thoroughly by gently vortexing, then centrifuge briefly.  
d. Incubate the IVT reaction for 16 hours at 40°C, then overnight at 4°C in a thermal cycler.  
e. After the incubation, centrifuge briefly (~5 seconds).  
f. Place the cRNA on ice briefly, or freeze immediately.

## Day 2 workflow

### 1 Purify cRNA

Before beginning the cRNA purification:

- Preheat the bottle of Elution Solution to 50–58°C for at least 10 minutes.
- Add 100% ethanol to the bottle of Nucleic Acid Wash Solution Concentrate before use.
- Make sure that the Nucleic Acid Binding Buffer Concentrate is completely dissolved. If not, warm the solution to < 50°C until the concentrate is solubilized.
- Vortex the Nucleic Acid Binding Beads vigorously before use to ensure that they are fully dispersed.

To purify the cRNA:

- a. Prepare the cRNA Binding Mix.

cRNA Binding Mix component	Volume for one reaction
Nucleic Acid Binding Beads	10 µL
Nucleic Acid Binding Buffer Concentrate	50 µL

- b. Add 60 µL of cRNA Binding Mix to each sample. Transfer each sample to a well of a U-Bottom Plate.
- c. Add 60 µL of isopropanol to each sample, then shake gently for 2 minutes.
- d. Capture the Nucleic Acid Binding Beads and discard the supernatant.
1. Move the plate to a magnetic stand to capture the magnetic beads.
  2. Carefully aspirate and discard the supernatant without disturbing the magnetic beads, then remove the plate from the magnetic stand.
- e. Wash twice with 100 µL of Nucleic Acid Wash Solution.
1. Add 100 µL of Nucleic Acid Wash Solution to each sample, then shake at moderate speed for 1 minute (setting 7 on the Lab-Line Titer Plate Shaker).
  2. Move the plate to a magnetic stand and capture the Nucleic Acid Binding Beads.
  3. Carefully aspirate and discard the supernatant without disturbing the Nucleic Acid Binding Beads, then remove the plate from the magnetic stand.
  4. Repeat substep e1 to substep e3 to wash each sample again with 100 µL of Nucleic Acid Wash Solution.
  5. Move the plate to a shaker and shake the plate vigorously for 1 minute (setting 10 on the Lab-Line Titer Plate Shaker).
- f. Elute cRNA with 40 µL of preheated Elution Solution.
1. Add to each sample 40 µL of preheated (55–58°C) Elution Solution. Incubate without shaking for 2 minutes.
  2. Vigorously shake the plate for 3 minutes (setting 10 on the Lab-Line Titer Plate Shaker), then check to make sure that the Nucleic Acid Binding Beads are fully dispersed. If they are not, continue shaking until the beads are dispersed and/or pipette up/down 3 times.
  3. Move the plate to a magnetic stand to capture the Nucleic Acid Binding Beads.
  4. Transfer the supernatant, which contains the eluted cRNA, to a nuclease-free multiwell plate.
- g. Place the cRNA on ice briefly, or freeze immediately.

- 2 Assess cRNA yield and size distribution Determine cRNA yield by UV absorbance or by using Quant-iT™ RiboGreen™ RNA Reagent. Optionally, use a bioanalyzer to determine cRNA size distribution.

- 3 Synthesize 2nd-cycle cDNA
- a. On ice, prepare 455 ng/µL cRNA. This is equal to 10 µg cRNA in a volume of 22 µL. If necessary, use nuclease-free water to bring the cRNA sample to 22 µL.
  - b. On ice, using supplied PCR tubes or plate, combine:
    - 22 µL of cRNA (10 µg)
    - 2 µL of Random Primers
  - c. Mix thoroughly by gently vortexing. Centrifuge briefly (~5 seconds). Place on ice.
  - d. Incubate for 5 minutes at 70°C, then 5 minutes at 25°C, then 2 minutes at 4°C. After the incubation, centrifuge briefly (~5 seconds).
  - e. At room temperature, prepare the 2nd-Cycle Master Mix in a nuclease-free tube.

2nd-Cycle Master Mix component	Volume for one reaction
2nd-Cycle Buffer Mix	8 µL
2nd-Cycle Enzyme Mix	8 µL
Total volume	16 µL

- f. Mix thoroughly by gently vortexing. Centrifuge briefly (~5 seconds). Proceed immediately to the next step.
- g. Transfer 16 µL of 2nd-Cycle Master Mix to each (24-µL) cRNA/Random Primer sample. Mix thoroughly by gently vortexing. Centrifuge briefly. Proceed immediately to the next step.
- h. Incubate for 10 minutes at 25°C, then 90 minutes at 42°C, then 10 minutes at 70°C, then for at least 2 minutes at 4°C. Immediately after the incubation, centrifuge briefly (~5 seconds), then place the sample on ice.

- 4 Hydrolyze using RNase H**
- Add 2  $\mu\text{L}$  of RNase H to the 2nd-Cycle cDNA. Mix by gently vortexing, then centrifuge briefly. Proceed immediately to the next step.
  - Incubate for 45 minutes at 37°C, then 5 minutes at 95°C, then for at least 2 minutes at 4°C. After the incubation, centrifuge briefly (~5 seconds), then place the sample on ice.

- 5 Purify 2nd-cycle cDNA**
- Before beginning the cDNA purification:
- Preheat the bottle of Elution Solution to 50–58°C for at least 10 minutes.
  - Make sure to add ethanol to the bottle of Nucleic Acid Wash Solution Concentrate before use.
  - Make sure that the Nucleic Acid Binding Buffer Concentrate is completely dissolved. If not, warm the solution to < 50°C until the concentrate is solubilized.
  - Vortex the Nucleic Acid Binding Beads vigorously before use to ensure they are fully dispersed.

To purify the 2nd-cycle cDNA:

- Prepare the cDNA Binding Mix for the experiment.

cDNA Binding Mix component	Volume for one reaction
Nucleic Acid Binding Beads	10 $\mu\text{L}$
Nucleic Acid Binding Buffer Concentrate	50 $\mu\text{L}$

- Add 18  $\mu\text{L}$  of nuclease-free water and 60  $\mu\text{L}$  of cDNA Binding Mix to each sample. Transfer each sample to a well of a U-Bottom Plate.
- Add 120  $\mu\text{L}$  of ethanol to each sample, then shake gently for 2 minutes.
- Capture the Nucleic Acid Binding Beads, then discard the supernatant.
  - Move the plate to a magnetic stand to capture the magnetic beads.
  - Carefully aspirate and discard the supernatant without disturbing the magnetic beads, then remove the plate from the magnetic stand.
- Wash twice with 100  $\mu\text{L}$  of Nucleic Acid Wash Solution.
  - Add 100  $\mu\text{L}$  of Nucleic Acid Wash Solution to each sample, then shake the samples at moderate speed for 1 minute (setting 7 on the Lab-Line Titer Plate Shaker).
  - Move the plate to a magnetic stand to capture the Nucleic Acid Binding Beads.
  - Carefully aspirate and discard the supernatant without disturbing the Nucleic Acid Binding Beads, then remove the plate from the magnetic stand.
  - Repeat substep e1 to substep e3 to wash a second time with 100  $\mu\text{L}$  of Nucleic Acid Wash Solution.
  - Move the plate to a shaker, then shake the plate vigorously for 1 minute to evaporate residual ethanol from the beads (setting 10 on the Lab-Line Titer Plate Shaker).
- Elute cDNA with 30  $\mu\text{L}$  of preheated Elution Solution.
  - Elute the purified cDNA from the Nucleic Acid Binding Beads by adding 30  $\mu\text{L}$  of preheated (55–58°C) Elution Solution to each sample. Incubate for 2 minutes at room temperature without shaking.
  - Vigorously shake the plate for 3 minutes (setting 10 on the Lab-Line Titer Plate Shaker).
  - Move the plate to a magnetic stand to capture the Nucleic Acid Binding Beads.
  - Transfer the supernatant, which contains the eluted cDNA, to a nuclease-free multiwell plate.
- Place the cDNA on ice briefly, or freeze immediately.

- 6 Assess cDNA yield and size distribution**
- Determine the concentration of a cDNA solution by measuring its absorbance at 260 nm. We recommend evaluating the absorbance of 1.5  $\mu\text{L}$  of cDNA sample using a NanoDrop™ Spectrophotometer.
  - (Optional) Use Quant-iT™ PicoGreen™ RNA Reagent to Assess cRNA yield.

- 7 Fragment and label the single-stranded cDNA**
- For instructions on how to fragment and label the single-stranded cDNA, see the *Affymetrix™ GeneChip™ WT Terminal Labeling and Hybridization User Manual* (Pub. No. 702880).

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.thermofisher.com/us/en/home/global/terms-and-conditions.html](http://www.thermofisher.com/us/en/home/global/terms-and-conditions.html). If you have any questions, please contact Life Technologies at [www.thermofisher.com/support](http://www.thermofisher.com/support).

The information in this guide is subject to change without notice.

### DISCLAIMER

TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

Corporate entity: Life Technologies Corporation | Carlsbad, CA 92008 USA | Toll Free in USA 1 800 955 6288

©2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. GeneChip is a trademark of Affymetrix, Inc.

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)

thermofisher.com

**ThermoFisher**  
SCIENTIFIC