

Denaturing Cell Extraction Buffer

Catalog Number FNN0091

Pub. No. MAN0016093 Rev. 2.0 (30)

Product Description

Cat. No.	FNN0091
Lot No.	See product label.
Quantity	100 mL
Description	Denaturing Cell Extraction Buffer is suitable for use in ELISA and western blotting. This can be used for Invitrogen™ phosphoELISAs (i.e., ERK, JNK, and MEK pELISAs) that require Sample Treatment (Invitrogen™ Cat. No. ST001) or boiling. After lysing cells with Denaturing Cell Extraction Buffer, lysates will not have to be sample-treated or boiled.
Buffer Formulation	<p>Proprietary</p> <p>Note: This Denaturing Cell Extraction buffer must be supplemented with 1 mM PMSF and protease inhibitor cocktail (i.e., Sigma-Aldrich™ Cat. No. P-2714) just prior to use.</p> <ul style="list-style-type: none"> For the PMSF addition, we recommend making a 0.3 M stock in DMSO, and adding sufficient volume for a final concentration of 1 mM (i.e., 17 µL per 5 mL extraction buffer). PMSF is very unstable and must be added just prior to use, even if added previously. For the protease inhibitor cocktail addition, we recommend Sigma-Aldrich™ Cat. No. P-2714, reconstituted according to the manufacturer's instructions, added at 100 µL per 1 mL Denaturing Cell Extraction Buffer. The stability of protease inhibitor-supplemented cell extraction buffer is 24 hours at 4°C.
Instructions	<p>This protocol has been successfully applied to several cell lines. Some optimization may be required for each specific application.</p> <ol style="list-style-type: none"> 1. Collect cells in PBS by centrifugation (non-adherent) or scraping from culture flasks (adherent). 2. Wash cells twice with cold PBS. 3. Remove and discard the supernatant and collect the cell pellet. 4. Lyse the cell pellet in Denaturing Cell Extraction Buffer for 30 minutes, on ice, with vortexing at 10-minute intervals. The volume of extraction buffer depends on the cell number and expression of target protein and level of phosphorylation. A suitable starting concentration is 10⁷ cells per mL of extraction buffer. 5. Transfer the extract to microcentrifuge tubes and centrifuge at 13,000 rpm for 15 minutes at 4°C. 6. Aliquot the clear lysate into clean microcentrifuge tubes. These samples are ready to assay. Lysates can be stored at -80°C. Avoid multiple freeze/thaw cycles.
Storage	Store at ≤-20°C. Thaw this buffer on ice. This buffer is stable for 2-3 weeks at 2-8°C or for up to 1 year when apportioned into working aliquots and stored at ≤-20°C.

thermofisher.com/support | thermofisher.com/askaquestion



Manufacturer: Bender MedSystems GmbH | Campus Vienna Biocenter 2 | 1030 Vienna, Austria

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2017 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.