

Pierce[®] Recombinant Protein A Agarose

20365 20366

0104.4

| Number | Description |
|--------|---|
| 20365 | Pierce Recombinant Protein A Agarose, 5mL settled resin |
| 20366 | Pierce Recombinant Protein A Agarose, 25mL settled resin Support: Crosslinked 6% beaded agarose supplied as a 50% slurry (e.g., 5mL of settled resin is equivalent to 10mL of 50% slurry) containing 0.02% sodium azide Binding Capacity: 6-8mg mouse or 15-17mg human polyclonal IgG per mL of settled resin using Protein A IgG Purification Buffers |

Storage: Upon receipt store product at 4-8°C. Product is shipped at ambient temperature.

Table of Contents

| | |
|--|---|
| Introduction | 1 |
| Important Product Information | 2 |
| Column Procedure for Antibody Purification Using Protein A Agarose | 2 |
| Example Immunoprecipitation (IP) Procedure Using Protein A Agarose | 3 |
| Troubleshooting | 4 |
| Additional Information Available on Our Website..... | 4 |
| Related Thermo Scientific Products | 4 |

Introduction

The Thermo Scientific Pierce Recombinant Protein A Agarose is a genetically engineered form of Protein A that is produced from a nonpathogenic strain of *Bacillus*. Although some non-essential regions have been removed, the recombinant is similar to native Protein A in size (MW ~44,600; apparent MW by SDS-PAGE ~45,000) and IgG-binding sites. The protein has four high affinity binding sites ($K_a = 10^8/M$) capable of specifically binding to the Fc region of immunoglobulin molecules from several species. Protein A is heat-stable and retains its native conformation even when exposed to denaturing reagents such as 4M urea, 4M thiourea and 6M guanidine hydrochloride.

Covalently immobilized Protein A matrices have been extensively used to purify IgG from several species of mammals. However, the interaction between Protein A and IgG is not equivalent for all species. Even within a species, Protein A interacts with some subgroups of IgG and not others. For instance, human IgG₁, IgG₂ and IgG₄ bind strongly, whereas IgG₃ does not bind and mouse IgG₁ binds poorly to Protein A. There are also many instances in which monoclonal antibodies do not bind to Protein A, such as the majority of rat immunoglobulins. Despite its variable binding characteristics, Protein A possesses properties that make it ideal for isolation of IgG. When using samples from species not well referenced in the literature, conduct preliminary experiments to determine if a particular immunoglobulin binds.

The Pierce Recombinant Protein A Agarose is prepared using a coupling method that results in excellent resin stability and binding characteristics. The special Protein A Buffers enhance IgG binding, resulting in higher yields compared with conventional methods.

Important Product Information

- The recommended buffers provide the highest efficiency of IgG binding and elution for most species. Using other buffer formulations may significantly alter the binding capacity and the wash volumes required for efficient purification. Therefore, optimization may be necessary when using other buffers.
- For optimal recovery, use a sample size such that the expected IgG load on the column is less than 80% of the maximum binding capacity. The total IgG content of serum is approximately 10-15mg/mL. The concentration of antibody in tissue culture supernatant varies considerably among hybridoma clones. Be aware that antibodies from fetal bovine serum (FBS) culture media supplement will be purified along with the antibody of interest.
- The crosslinked 6% beaded agarose support can tolerate commonly used water-miscible solvents when they are added in a stepwise gradual manner. The agarose support will compress under pressure causing column flows to slow. Never freeze agarose supports, as this will cause irreversible damage to the bead structure.
- Serum samples, ascites fluid, plasma or tissue culture supernatant may be used with this product.

Column Procedure for Antibody Purification Using Protein A Agarose

Note: The following protocol is for using a gravity-flow column packed with 1mL of resin (i.e., 2mL of the 50% slurry). When using columns containing other resin volumes, reagent amounts must be adjusted accordingly. See the Additional Information Available on Our Website section for batch and spin cup methods.

A. Additional Materials Required

- Column capable of containing at least 1mL resin bed volume such as the Disposable Polypropylene Columns (Product No. 29922) or the Column Trial Pack (Product No. 29925) that contains two each of three column sizes.
- Binding Buffer: Protein A IgG Binding Buffer (Product No. 21001 or 21007)
- Elution Buffer: IgG Elution Buffer (Product No. 21004 or 21009) or 0.1M glycine, pH 2-3
- Neutralization Buffer: 1mL of high-ionic strength alkaline buffer such as 1M phosphate or 1M Tris (pH 7.5-9)
- (Optional): Thermo Scientific Slide-A-Lyzer Dialysis Cassette or Zeba Spin Desalting Columns (Product No. 89893) for buffer exchange

B. Antibody Purification Procedure

1. Equilibrate Immobilized Protein A and all buffers to room temperature.
2. Carefully pack the column with 2mL of resin slurry, following the instructions provided with the columns.
3. Equilibrate the column by adding 5mL of the Binding Buffer and allowing the solution to drain through the column.
Note: To avoid air bubbles being drawn into the resin, remove the top cap before the bottom cap when opening column.
4. Dilute sample at least 1:1 with Binding Buffer before application to the Protein A Column to maintain the proper ionic strength and pH for optimal binding.
Note: Plasma may become hazy upon dilution with the Binding Buffer because of lipoprotein precipitation. Centrifuge the diluted sample at $10,000 \times g$ for 20 minutes and apply the supernatant to the equilibrated Immobilized Protein A.
5. Apply the diluted sample to the column and allow it to flow completely into the resin. Do not allow the resin bed to run dry. Any volume may be applied provided the total amount of antibody is less than 80% of column capacity.
Note: If the sample contains more IgG than can bind to the Protein A column (or is an antibody type that does not bind to Protein A), the flow-through will contain excess antibody. By saving the flow-through, non-bound antibody can be recovered and examined by antibody-specific assays.
6. Wash the Protein A column with 15mL of the Binding Buffer.
Note: If desired, verify that all non-bound proteins are removed from the column by collecting separate 2mL fractions as the solution drains and measuring their absorbance at 280nm. The last fractions should have absorbances similar to Binding Buffer alone.

7. Elute antibodies with 5mL of Elution Buffer and collect 0.5-1mL fractions. Immediately adjust eluted fractions to physiologic pH by adding 100 μ L of the Neutralization Buffer per 1mL of eluate. Monitor the elution by measuring the absorbance at 280nm or by protein assay such as Thermo Scientific BCA Protein Assay Kit (Product No. 23225).
8. Pool the eluted IgG fractions that contain the highest absorbance. The purified antibodies may be used directly for SDS-PAGE, or the buffer may be exchanged by dialysis or desalting column to one that is compatible with the specific downstream application (see Related Thermo Scientific Products).
9. Regenerate column by washing with 12mL of Elution Buffer. Columns may be regenerated at least 10 times without significant loss of binding capacity.
10. For storage, wash column with 5mL of water containing 0.02% sodium azide. When approximately 3mL of solution remains, replace the bottom cap followed by the top cap on the column. Store columns upright at 4°C.

Example Immunoprecipitation (IP) Procedure Using Protein A Agarose

A. Additional Materials Required

- 1.5-2mL microcentrifuge tube
- IP Buffer: 25mM Tris, 150mM NaCl; pH 7.2 (Thermo Scientific BupH Tris Buffered Saline Pack, Product No. 28379)
- Antigen Sample: Antigen-containing lysate or sample prepared in IP Buffer or other buffer that is compatible with both the desired antibody binding interaction and the binding of antibody to Protein A
- Elution Buffer: IgG Elution Buffer (Product No. 21004) or 0.1-0.2M glycine•HCl buffer, pH 2.5-3.0
- Electrophoresis Loading Buffer: Lane Marker Reducing Sample Buffer (5X), (Product No. 39000)
- Neutralization Buffer (optional): 1mL of strong alkaline buffer, such as 1M phosphate or 1M Tris, (pH 7.5-9)

B. Immunoprecipitation Procedure

Note: This procedure uses 50 μ L of settled resin (100 μ L resin slurry). This amount of resin is sufficient to bind 25-250 μ g of antibody. Depending on the amount of antibody needed to immunoprecipitate the desired amount of antigen, scale the amount of resin and suggested wash and elution volumes accordingly. To allow for proper mixing, make sure the total reaction volume does not completely fill the microcentrifuge tube.

1. In a microcentrifuge tube, combine 50-1000 μ L of the Antigen Sample and the optimized amount of antibody. Incubate the reaction overnight at 4°C.
2. Add 100 μ L of Immobilized Protein A resin slurry to the antigen-antibody complex.
3. Incubate reaction with gentle mixing for 2 hours at room temperature.
4. Add 0.5mL of IP Buffer, centrifuge for 2-3 minutes at 2500 \times g and discard supernatant. Repeat this step several times.
5. To elute the immune complex, add 50 μ L of Elution Buffer and incubate for 5 minutes. Centrifuge tube for 1-3 minutes at 2500 \times g and collect the supernatant. Repeat this step and combine the two supernatant fractions.

Alternatively, wash the complex-bound resin with 0.5mL water, centrifuge for 2-3 minutes at 2500 \times g, and discard supernatant. Add Electrophoresis Loading Buffer to the complex-bound resin and incubate for 5 minutes at 95°C. Centrifuge the resin mixture at 2500 \times g, collect the supernatant and evaluate by SDS-PAGE.

6. Adjust eluate to physiological pH by adding ~10 μ L of the Neutralization Buffer per 100 μ L of eluate. The IP products may be used directly for SDS-PAGE, or the buffer may be exchanged by dialysis or desalting column to one that is compatible with the specific downstream application (see Related Thermo Scientific Products).

Troubleshooting

| Problem | Possible Cause | Solution |
|--|--|--|
| Flow of the column is exceedingly slow (i.e., < 0.5mL/minute) | Outgassing of buffers or sample on the column is blocking the resin pores with microscopic air bubbles | Degas buffers and remove air bubbles from column (see Additional Information Available on Our Website section for suggested Tech Tip protocol) |
| Considerable antibody purified but no specific antibody of interest detected | Antibody of interest is at very low concentration | Use serum-free medium for cell supernatant samples |
| | | Affinity purify the antibody using the specific antigen coupled to an affinity support such as Thermo Scientific AminoLink Plus Immobilization Kit (Product No. 44894) |
| Antibody was purified but is degraded, as determined by lack of function | Antibody is sensitive to low-pH Elution Buffer | Try Gentle Ag/Ab Elution Buffer (see Related Thermo Scientific Products) |
| | Downstream application is sensitive to neutralized Elution Buffer | Desalt or dialyze eluted sample into suitable buffer |
| No protein detected in any elution fraction | Sample devoid of antibody species or isotype that binds to Protein A | Refer to the Binding Characteristics Table for Protein A (see Additional Information Available on Our Website section) |
| Antibody bands interfere with antigen detection after immunoprecipitation | Antibody is eluted along with the antigen | Crosslink antibody to the Protein A by using the Thermo Scientific Pierce Crosslink IP Kit (Product No. 26147) |

Additional Information Available on Our Website

- Tech Tip #34: Binding characteristics of Protein A, Protein G, Protein A/G and Protein L
- Tech Tip #4: Batch and spin cup methods for affinity purification of proteins
- Tech Tip #7: Remove air bubbles from columns to restore flow rate
- Tech Tip #13: Pack beaded affinity resin into columns
- Tech Tip #29: Degas buffers for use in affinity and gel filtration columns
- Tech Tip #43: Protein stability and storage

Related Thermo Scientific Products

| | |
|--------------|---|
| 53139 | Protein A UltraLink Resin, 5mL |
| 20333 | Pierce Protein A Agarose, 5mL |
| 22810 | Pierce Protein A Plus Agarose, 1mL |
| 20397 | Pierce Protein G Agarose, 25mL |
| 22851 | Pierce Protein G Plus Agarose, 2mL |

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.