

PrepSEQ™ Residual DNA Sample Preparation Kit

For (Genomic DNA) CHO, *E. coli*, HEK293, Human, MDCK, NSO, *Pichia*, Sf9 and Baculovirus, and Vero
 For (Plasmid DNA) Kanamycin resistance

Catalog Numbers 4402085, 4458435, A46014, A26366, 4464335, 4458441, 4464336, A46066, A41797, A50337

Pub. No. 4469839 Rev. G

Note: For safety and biohazard guidelines, see the “Safety” appendix in the *PrepSEQ™ Sample Preparation Kits User Guide* (Pub. No. 4469838). Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

- Prepare the reagents: before first use of the kit 1
- Prepare reagents: before each use of the kit 1
- Manual protocol for DNA/RNA extraction 2
- Automated protocol for DNA/RNA extraction 4
- Limited product warranty 5

Prepare the reagents: before first use of the kit

Magnetic beads

1. Set a block heater to 37°C.
2. Incubate the Magnetic Particle suspension at 37°C for a minimum of 10 minutes with intermittent vortexing at setting #7, or until the particles are completely suspended.

Binding Solution

1. Add 30 mL of 100% isopropanol to the Binding Solution bottle.
2. Label the bottle to indicate that it contains isopropanol, then store the bottle at ambient temperature.

Wash Buffer Concentrate

1. Add 74 mL of 95% ethanol to one bottle of PrepSEQ™ Wash Buffer Concentrate, then mix completely.
2. Label the bottle to indicate that it contains ethanol, then store the bottle at room temperature.

Prepare reagents: before each use of the kit

Proteinase K (PK) mix

- Use Proteinase K (PK) Buffer II for new manual protocols and automated protocols.
Note: Proteinase K (PK) Buffer is also provided in the kit for use with existing manual protocols that have been internally validated with this buffer.
- Prepare a fresh mix before each use of the kit.
- Include a 10% overage to account for pipetting losses.

Component	Number of extractions				
	1	7	10	13	25
Proteinase K, 20 mg/mL	10 µL	70 µL	100 µL	130 µL	250 µL
Proteinase K (PK) Buffer II or Proteinase K (PK) Buffer	60 µL	420 µL	600 µL	780 µL	1,500 µL

Lysis solution

- Prepare a fresh mixture immediately before use or during Proteinase K incubation.
- Prepare 360 µL (amount required) of lysis solution mix per sample.

Reagent	Volume for ~20 extractions
Glycogen, 5 mg/mL	180 µL
Yeast tRNA, 10 mg/mL	4 µL
Lysis Buffer	7,600 µL
Total	7,784 µL

Manual protocol for DNA/RNA extraction

(Plasmid samples only) Add Yeast tRNA

Plasmid samples require the use of **Yeast tRNA** as a carrier during the extraction process.

1. Dilute the **Yeast tRNA**.

Table 1 Diluted Yeast tRNA

Component	Volume
Yeast tRNA (10mg/mL)	5 µL
PBS (1X), pH 7.2	245 µL
Total	250 µL

2. Add 5 µL **Diluted Yeast tRNA** to 370 µL of each test sample before extraction. This is sufficient for triplicate 100 µL extractions.

Digest the test samples and controls

1. Set a block heater to 56°C. If available, set a second block heater to 70°C.

2. Label 2-mL Safe-Lock tubes:
 - 3 for each sample
 - 3 for each sample + ERC
 - 3 for NEG

- 3 for each sample
- 3 for each sample + ERC
- 3 for NEG

3. Add 100 µL of sample, sample + ERC, or 1X PBS to into each tube.

Note: Ensure that **Diluted Yeast tRNA** was added to each plasmid sample.

4. Add 10 µL of 5 M NaCl and 70 µL Proteinase K/Proteinase K Buffer II mix.

5. Briefly vortex and centrifuge.

6. Incubate at 56°C for 30 minutes.

If only one block heater is available, after this incubation step is complete, reset the block heater to 70°C for the elution step.

Note: For samples with high protein concentration, extending the incubation time to 60 minutes can increase recovery.

7. Cool samples to room temperature.
8. Add 360 μL freshly made Lysis solution mix to each tube.

Bind the DNA/RNA

1. Vortex the Magnetic Particles to resuspend the particles.
Note: The appearance of the mixture should be homogeneous.
2. Add 30 μL of the Magnetic Particles to each tube.
3. Add 400 μL Binding Solution to each tube.
4. Mix and vortex the tubes:
 - a. Close the caps, **immediately** invert each tube twice to mix.
 - b. Vortex the tubes in the vortex adaptor for 5 minutes at setting #7.
5. Briefly centrifuge the tubes for 15 seconds at top speed ($>15,000 \times g$) to collect the Magnetic Particles at the bottom of the tubes.
6. Place the tubes in the magnetic stand with the pellet against the magnet, then let the tubes stand for 5 minutes or until the solution is clear.
7. Without disturbing the magnetic beads, remove the supernatant using a PIPETMAN™ pipette or by aspiration.

Wash the DNA/RNA

For aspiration of liquid supernatants and wash buffers during sample preparation, we recommend attaching the waste-collection bottle to the vacuum using tubing that can accommodate 200- μL pipette tips.

1. Remove the tube rack (with tubes) from the magnetic stand, then add 300 μL of Wash Solution to the tubes. Vortex the tubes for 5 seconds at room temperature at setting #7.
2. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
3. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA/RNA are magnetically captured after approximately 1 minute.
4. Without disturbing the Magnetic Particles, remove the supernatant using a PIPETMAN™ pipette or by aspiration.
5. Remove the tube rack (with tubes) from the magnetic stand, then add 300 μL of Wash Solution to each tube for a second wash. Vortex the tubes for 5 seconds at room temperature at setting #7.
6. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
7. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA/RNA are magnetically captured after approximately 1 minute.
8. Open all tubes, then start the 5-minute timer.
9. Without disturbing the Magnetic Particles, remove the supernatant using a PIPETMAN™ pipette or by aspiration.
Use a P200 to remove the remaining solution from the bottom of the tube.
10. With the tube lid open, air-dry the Magnetic Particles pellet in the magnetic stand for no more than 5 minutes at room temperature.

IMPORTANT! Air-dry to remove ethanol from the Wash Solution. After dried, the DNA/RNA stays bound to the magnetic beads. Do not over-dry; over-dried beads are not easily resuspended.

Elute the DNA/RNA

1. Add 50 μL of Elution Buffer to each tube.
2. Vortex the tubes for 20 seconds at high speed, then incubate the tubes at 70°C for 7 minutes. Vortex the tubes two to three times during the incubation to help resuspension.
Note: (*Optional*) If vortexing does not result in full resuspension, then wash the beads off the tube by pipetting.
3. Centrifuge the tubes in a microcentrifuge at top speed ($>15,000 \times g$) for a maximum of 20 seconds. Do not centrifuge for >20 seconds.
4. Place the tubes in the magnetic stand, then let the tubes stand for 1 minute.
Note: The Magnetic Particles with the bound DNA/RNA are magnetically captured after approximately 1 minute.
5. Without disturbing the Magnetic Particles, transfer the liquid phase containing the eluted DNA/RNA to a new nonstick 1.5-mL microcentrifuge tube.
6. Centrifuge the tube at top speed ($>15,000 \times g$) for 3 minutes to collect the Magnetic Particles at the bottom, then place the tubes in the magnetic stand for 1 minute.
7. Without disturbing the Magnetic Particles, transfer the liquid phase containing the eluted DNA/RNA to the nonstick 1.5-mL microcentrifuge tube.

Note: Magnetic Particles can inhibit PCR.

When you finish the sample extraction procedure, see the *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836) or the *ViralSEQ™ Real-Time PCR Kits User Guide* (Pub. No. 4445235) to set up PCR for DNA/RNA quantification.

Use 10 μL of the eluate in the real-time PCR.

Automated protocol for DNA/RNA extraction

You can use the KingFisher™ Flex or MagMAX™ Express 96-deep well automation platforms to automate the extraction of DNA/RNA. For all chemicals, read the Safety Data Sheet (SDS) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

(Plasmid samples only) Add Yeast tRNA

Plasmid samples require the use of **Yeast tRNA** as a carrier during the extraction process.

1. Dilute the **Yeast tRNA**.

Table 2 Diluted Yeast tRNA

Component	Volume
Yeast tRNA (10mg/mL)	5 μL
PBS (1X), pH 7.2	245 μL
Total	250 μL

2. Add 5 μL **Diluted Yeast tRNA** to 370 μL of each test sample before extraction. This is sufficient for triplicate 100 μL extractions.

Prepare the plates

Prepare the Wash 1, Wash 2, and Elution plates:

Plate name	Plate type	Volume of buffer to add
Wash 1	96 deep-well plate	300 μL of Wash buffer
Wash 2	96 deep-well plate	300 μL of Wash buffer
Elution	96 deep-well plate	200 μL of Elution buffer

Prepare the lysis plate

In all steps that require pipetting, dispense liquid at bottom center of the wells.

1. Add 100 µL to the appropriate wells of the 96 deep-well Lysis plate:
 - 3 wells for each sample
 - 3 wells for each sample + ERC
 - 3 wells for NEG

Note: Ensure that **Diluted Yeast tRNA** was added to each plasmid sample.

2. Add 10 µL of 5 M NaCl to each sample well.
3. Add 70 µL Proteinase K/Proteinase K (PK) Buffer II mix to each sample well.

Process samples on the instrument

1. Select the script or program for the instrument you are using:

Instrument	Select
KingFisher™ Flex	PrepSEQ_resDNA_v1 script
MagMAX™ Express-96	PrepSEQ_resDNA_2011 PrepSEQ_1hr_resDNA (if installed)

2. Load the plates into the instrument in the order listed below. After loading each plate, press **START** to move the turntable.
 - a. Comb loading plate
 - b. Elution plate with 200 µL of Elution Buffer
 - c. Wash 2 plate with 300 µL of wash buffer
 - d. Wash 1 plate with 300 µL of wash buffer
 - e. Lysis plate

3. Press **START** to begin the PK digestion process.

The instrument mixes the samples for 10 seconds at fast speed, then incubates the samples at 56°C for 30 minutes, mixing at slow speed. When digestion is complete, the instrument pauses and returns the Lysis plate to the loading position.

4. After the digestion step is complete, add additional components to the Lysis plate:
 - a. Remove the Lysis plate from the instrument.
 - b. Add 360 µL of Lysis Solution to each sample well.
 - c. Add 30 µL of Magnetic Particle suspension to each sample well.
 - d. Add 400 µL of Binding Solution to each sample well, then immediately pipet up-and-down three times to mix.
 - e. Place the plate back into the instrument loading position, then press **START** to begin binding.

5. When DNA/RNA extraction is finished, the instrument returns the Elution plate to the loading position.

When you finish the sample extraction procedure, see the *resDNASEQ™ Quantitative DNA Kits User Guide* (Pub. No. 4469836) or the *ViralSEQ™ Real-Time PCR Kits User Guide* (Pub. No. 4445235) to set up PCR for DNA/RNA quantitation.

Note: Use 10 µL of the eluate in the real-time PCR.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale at www.thermofisher.com/us/en/home/global/terms-and-conditions.html. If you have any questions, please contact Life Technologies at www.thermofisher.com/support.



The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, THERMO FISHER SCIENTIFIC INC. AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Revision history: Pub. No. 4469839

Revision	Date	Description
G	2 April 2021	Added the resDNASEQ™ Quantitative Plasmid DNA - Kanamycin Resistance Gene Kit (Cat. No. A50337).
F	5 August 2020	Added the resDNASEQ™ Quantitative Sf9 and Baculovirus DNA Kit (Cat. No. A46066).

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products, you accept the terms and conditions of all applicable Limited Use Label Licenses.

©2021 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.