

CaptureSelect™ FcXL Ligand Leakage ELISA

INSTRUCTIONS

Publication Number MAN0010266 Rev. A.0

■ Product description	1
■ Principle of the assay	1
■ Kit contents	1
■ Required materials and equipment (not provided)	1
■ Procedure 1: Samples without target protein	2
■ Procedure 2: Samples with target protein	3
■ Validate the assay	4
■ Ordering information	4
■ Support	4
■ Safety information	4
■ Limited product warranty	4

Product description

The CaptureSelect™ FcXL Ligand Leakage ELISA (Enzyme Linked Immuno-Sorbent Assay) is designed for the detection of FcXL affinity ligand as low as 1 ng/mL that may be present in product purified with CaptureSelect™ FcXL affinity media, which contains the FcXL affinity ligand as capturing agent. The assay is designed to minimize interference and to provide accurate quantitation in the presence of human IgG antibodies and other proteins. The FcXL Ligand Leakage ELISA can be used as a tool to aid in optimal purification process development and in routine quality control of in-process streams as well as final product.

Principle of the assay

The CaptureSelect™ leakage assays enable detection of the affinity ligand in solutions with and without the presence of the target protein. These sandwich assays involve the following steps:

- A microtiter plate is coated with affinity-purified antiaffinity ligand.
- Samples containing the affinity ligand are incubated in the coated plate wells.
- Bound affinity ligand is detected by biotinylated affinity ligand.
- Streptavidin-horseradish peroxidase conjugate is added to bind to the biotinylated antibody in the sandwich complex.
- Substrate reactive with horseradish peroxidase (tetramethylbenzidine-hydrogen peroxide) is added.
- The amount of hydrolyzed substrate is determined and is directly proportional to the concentration of affinity ligand present.

Kit contents

Note: After thawing and before use, spin the tubes to ensure that all reagents are at the bottom of the tube.

Item	Description	Storage
Coating reagent (green label)	Goat IgG anti-FcXL affinity ligand, 100 µL	-20°C (-4°F)
Standard solution (blue label)	CaptureSelect™ FcXL affinity ligand, 100 µL	
Biotinylated reagent (yellow label)	Biotinylated Goat IgG anti-FcXL affinity ligand, 100 µL	

Required materials and equipment (not provided)

- PBS: Phosphate buffered saline pH 7.4
- PBST: Phosphate buffered saline (PBS) pH 7.4 + 0.05 (v/v)% Tween® 20 Solution
- Bovine Serum Albumin (BSA), Fraction V 99% pure (Sigma-Aldrich A3059)
Note: Other qualities of Bovine Serum Albumin or other blocking proteins may result in higher background levels.
- Dilution Buffer A for assays *without* target protein:
 - Dilution Buffer A: 0.05 (v/v)% Tween® 20 Solution in PBS pH 7.4
 - 2X Dilution Buffer A: 0.1 (v/v)% Tween® 20 Solution in PBS pH 7.4
- Standard Dilution Buffer B for assays *with* target protein: 0.05 (v/v)% Tween® 20 Solution in PBS pH 7.4 plus IgG at a concentration that is half of the concentration of target in samples
- Blocking solution: 4 (w/v)% BSA in PBS pH 7.4
- Human IgG (for protocol for samples containing IgG)
- Streptavidin-horseradish peroxidase diluted immediately before use according to manufacturer guidelines
- Tetramethylbenzidine (TMB) and hydrogen peroxide (H₂O₂) substrate (prepare 1:1 solution immediately before use)
- 1 M H₂SO₄
- Microtiter plate (Maxisorp, Nunc)
- Microtiter plate shaker
- Microtiter plate reader (450 nm)
- Milli-Q® Water

Procedure 1: Samples without target protein

Coat the plate

1. Make a 1:100 dilution of the Coating reagent with PBS pH 7.4.
2. Add 100 μL diluted Coating reagent to each well in a microtiter plate and incubate overnight at 4°C (39°F). This step captures fragments of the leached ligand.

Prepare standards

1. Prepare a 6.4 $\mu\text{g}/\text{mL}$ stock Standard dilution series: Add 10 μL Standard solution to 770 μL Dilution Buffer A.
2. Using the stock Standard solution from step 1, prepare a standard dilution series according to the table below.

Tube	Concentration (ng/mL)	Standard	Dilution Buffer A
1	64.0	10 μL diluted Standard solution	990 μL
2	16.0	250 μL 64.0 ng/mL	750 μL
3	8.0	500 μL 16.0 ng/mL	500 μL
4	4.0	500 μL 8.0 ng/mL	500 μL
5	2.0	500 μL 4.0 ng/mL	500 μL
6	1.0	500 μL 2.0 ng/mL	500 μL
7	0.5	500 μL 1.0 ng/mL	500 μL
8	0.25	500 μL 0.5 ng/mL	500 μL
9	0	0	500 μL

Prepare assay samples

Dilute 75 μL sample with 75 μL 2X Dilution Buffer A.

ELISA assay procedure


1. Block the plate:
 - a. Wash the coated plate 5 times with PBST.
 - b. Add 250 μL /well of Blocking solution to the coated plate. Leave at room temperature for 30 minutes on a microtiter plate shaker.
 - c. Wash the plate 1 time with PBST.
2. Add samples and standards:
 - a. Add 100 μL of each concentration of the standard dilution series (0 to 64.0 ng/mL) or sample to appropriate wells.
 - b. Incubate the plate 1 hour at room temperature on a microtiter plate shaker.
 - c. Wash the plate 5 times with PBST.
3. Add Biotinylated reagents (detects leached ligand):
 - a. Make a 1:100 dilution of the Biotinylated reagents with Dilution Buffer A.

- b. Add 100 μL diluted Biotinylated reagents to each well containing sample or standard and incubate the plate 1 hour at room temperature.
 - c. Wash the plate 5 times with PBST.
4. Add streptavidin-horseradish peroxidase (colorimetric reagent that binds to the biotinylated reagents):
 - a. Dilute in Dilution Buffer A according to the manufacturer guidelines.
 - b. Add 100 μL diluted streptavidin-horseradish peroxidase to each well containing sample or standard.
 - c. Incubate the plate 1 hour at room temperature on a microtiter plate shaker.
 - d. Wash the plate 5 times with PBST.
 - e. Wash the plate 2 times with Milli-Q® water.
 5. Develop and read the plate:
 - a. Prepare a 1:1 solution of TMB:H₂O₂ substrate.
 - b. Add 100 μL to each well containing sample or standard.
 - c. Incubate the plate for approximately 10 minutes on a microtiter plate shaker.
 - d. When the background signal starts to develop, add 50 μL 1 M H₂SO₄ to stop the coloring reaction and achieve a maximal signal-to-noise ratio.
 - e. Measure the OD of the microtiter plate at 450 nm with a microtiter plate reader.

Calculate results

Create a standard curve using the OD values from the standards reported in ng/mL. Use curve-fitting routines such as 4-parameter logistic fit. Do not use linear regression analysis to interpolate values for samples, which may lead to significant inaccuracies.

Figure 1 Example calibration curve CaptureSelect™ FcXL Ligand Leakage ELISA without target protein



Procedure 2: Samples with target protein

Coat the plate

1. Make a 1:100 dilution of the Coating reagent with PBS pH 7.4.
2. Add 100 μL diluted Coating reagent to each well in a microtiter plate and incubate overnight at 4°C (39°F).

Prepare standards

1. Prepare a 6.4 $\mu\text{g}/\text{mL}$ stock Standard solution: Add 10 μL Standard solution to 770 μL Standard Dilution Buffer B.
2. Using the stock Standard solution from step 1, prepare a standard dilution series according to the table below.

Tube	Concentration (ng/mL)	Standard	Standard Dilution Buffer B
1	64.0	10 μL diluted Standard solution	990 μL
2	16.0	250 μL 64.0 ng/mL	750 μL
3	8.0	500 μL 16.0 ng/mL	500 μL
4	4.0	500 μL 8.0 ng/mL	500 μL
5	2.0	500 μL 4.0 ng/mL	500 μL
6	1.0	500 μL 2.0 ng/mL	500 μL
7	0.5	500 μL 1.0 ng/mL	500 μL
8	0.25	500 μL 0.5 ng/mL	500 μL
9	0	0	500 μL

Prepare assay samples

1. Dilute 75 μL sample with 75 μL PBS Dilution Buffer A.
2. Incubate the samples and standard dilution series for 15 minutes at 95°C (203°F).
3. Centrifuge the heat-treated samples and standard dilution series for 5 minutes at 20,000 $\times g$.
4. Transfer the supernatants to a clean tube.

ELISA assay procedure


1. Block the plate:
 - a. Wash the coated plate 5 times with PBST.
 - b. Add 250 μL /well of Blocking solution to the coated plate. Leave at room temperature for 30 minutes on a microtiter plate shaker.
 - c. Wash the plate 1 time with PBST.
2. Add samples and standards:
 - a. Add 100 μL of each concentration of the standard dilution series (0 to 64.0 ng/mL) or sample to appropriate wells.
 - b. Incubate the plate 1 hour at room temperature on a microtiter plate shaker.
 - c. Wash the plate 5 times with PBST.

3. Add Biotinylated reagents:
 - a. Make a 1:100 dilution of the Biotinylated reagents with Dilution Buffer A.
 - b. Add 100 μL diluted Biotinylated reagents to each well and incubate the plate 1 hour at room temperature.
 - c. Wash the plate 5 times with PBST.
4. Add diluted streptavidin-horseradish peroxidase:
 - a. Dilute in Dilution Buffer A according to the manufacturer guidelines.
 - b. Add 100 μL diluted streptavidin-horseradish peroxidase to each well containing sample or standard.
 - c. Incubate the plate 1 hour at room temperature on a microtiter plate shaker.
 - d. Wash the plate 5 times with PBST.
 - e. Wash the plate 2 times with Milli-Q® water.
5. Develop and read the plate:
 - a. Add 100 μL 1:1 mixed TMB/ H_2O_2 substrate per well.
 - b. Incubate the plate for approximately 25 minutes on a microtiter plate shaker.
 - c. When the background signal starts to develop, add 50 μL 1 M H_2SO_4 to stop the coloring reaction and achieve a maximal signal-to-noise ratio.
 - d. Measure the OD of the microtiter plate at 450 nm with a microtiter plate reader.

Calculate results

Construct a standard curve with values reported in ng/mL. Use curve-fitting routines such as 4-parameter logistic fit. Do not use linear regression analysis to interpolate values for samples, which may lead to significant inaccuracies.

Figure 2 Example calibration curve CaptureSelect™ FcXL Ligand Leakage ELISA in the presence of 1 mg/mL human IgG



Validate the assay

Perform validation studies that include at least the following experiments to validate this kit for your application: 1) Intra- and inter-assay precision experiments to establish reproducibility, 2) Recovery experiments using test samples with known amounts of the 500 µg/mL Standard solution, which is included in the kit.

Ordering information

CaptureSelect™ FcXL Ligand Leakage ELISA	Cat. no.
1 assay	810328001
10 assays	810328010

Support

For more information on CaptureSelect™ products, go to www.lifetechnologies.com/captureselect.

For the latest services and support information for all locations, go to www.lifetechnologies.com, then click the link for **Support**, or contact your local representative.

Safety information

Obtaining SDSs

Safety Data Sheets (SDSs) are available from www.lifetechnologies.com/support.

Note: For the SDSs of chemicals not distributed by Thermo Fisher Scientific, contact the chemical manufacturer.

Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at www.lifetechnologies.com/termsandconditions. If you have any questions, please contact Life Technologies at www.lifetechnologies.com/support.

For Research Use Only. Not intended for use in diagnostic procedures.

The information in this guide is subject to change without notice.

DISCLAIMER: LIFE TECHNOLOGIES CORPORATION AND/OR ITS AFFILIATE(S) DISCLAIM ALL WARRANTIES WITH RESPECT TO THIS DOCUMENT, EXPRESSED OR IMPLIED, INCLUDING BUT NOT LIMITED TO THOSE OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, OR NON-INFRINGEMENT. TO THE EXTENT ALLOWED BY LAW, IN NO EVENT SHALL LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) BE LIABLE, WHETHER IN CONTRACT, TORT, WARRANTY, OR UNDER ANY STATUTE OR ON ANY OTHER BASIS FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING BUT NOT LIMITED TO THE USE THEREOF.

Important Licensing Information: This product may be covered by one or more Limited Use Label License s. By use of this product, you accept the terms and conditions of all applicable Limited Use Label Licences.

All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

Tween is a registered trademark of Uniqema Americas. LLC. Milli-Q is a registered trademark of Merck KGaA.

© 2014 Thermo Fisher Scientific Inc. All rights reserved.

For support visit lifetechnologies.com/support or email techsupport@lifetech.com
lifetechnologies.com

10 April 2014

