invitrogen™

Lot #: 091503 Expiration : 30.06.2010

Catalog #: CMC2223

1.	Coating Antibody:	Anti-Mouse IL-13 (0.25 mg / 0.125 mL)
	Part Number:	5M.222.09
	Lot Number:	9C15/3
	Form:	Liquid, 2 vials, contains 0.1% sodium azide
	Storage:	Store at 2-8°C until expiration date.
	Recommended Dilution:	Dilute to 4 µg/mL with Coating Buffer A (Cat. # CB07100, or see Recommended Buffers). For example, to make
		10 mL (enough to coat 1 plate), add 20 µL coating antibody to 9.980 mL Coating Buffer A.
2.	Detection Antibody:	Anti- Mouse IL-13 Biotin (3.125 ug/ 0.125 mL)

2.	Detection Antibody:	Anti- Mouse IL-13 Blotin (3.125 μg/ 0.125 mL)
	Part Number:	5M.222.03
	Lot Number:	9C15/3
	Form:	Liquid, 2 vials, contains 0.1% sodium azide
	Storage:	Store at 2-8°C until expiration date.
	Recommended Dilution:	Dilute to 0.0625 µg/ mL with Assay Buffer (Cat. # DS98200, or see Recommended Buffers). For example, to
		make enough for 1 plate, add 13.75 µL detection antibody to 5.48625 mL Assay Buffer.

next tube and vortexing each tube. Assay Buffer should be used as the zero standard.

3. Standard:

Part Number:
Lot Number:
Form:
Storage:
Reconstitution:

Standard Curve:


4. Streptavidin-HRP:

Part Number:		
Lot Number:		
Form:		
Storage:		
Recommended Dilution:		

make enough for 1 plate, add 13.75 μL detection antibody to 5.48625 mL Assay Buffer.
Recombinant Mouse IL-13
5M.222.10 (additional vials of standard may be purchased using this part number)
3E1/1
Lyophilized, 3 vials
Store at 2-8°C.
Reconstitute with Assay Buffer (Cat. # DS98200 or see Recommended Buffers) to yield a stock of 5,000 pg/mL.
After 10 minutes of rehydratation, use the standard stock immediately or aliquot in polypropylene tubes and freeze at -80°C. *Do not store at room temperature or at 4°C and do not subject to more than one freeze-thaw cycle*.
Dilute standard stock to 500 pg/mL (60 μL stock plus 540 μL Assay Buffer) with Assay Buffer (Cat. # DS98200 or see Recommended Buffers). Add 300 μL Assay Buffer to 6 tubes and label as 250, 125, 62.5, 31.2, 15.6 and 7.8 pg/mL. Make serial dilutions starting with 500 pg/mL by transferring 300 μL of each standard to

0.25 mL
41.000.19
9D7/1
Liquid, 1vial, contains 0.05% thymol
Store at 2-8°C until expiration date.
Dilute to 1/4000. For example, to make

Dilute to 1/4000. For example, to make enough for 1 plate, add 2.5 μ L of streptavidin-HRP to 9.9975 mL of Assay Buffer (Cat. # DS98200 or see Recommended Buffers).



Representative standard curve was generated by following the recommended assay procedure, which includes the use of the **BioSource** CytoSet[™] Buffer Set (Cat. # CNB0011)

> biosource.com BioSource International, Inc., USA • 542 Flynn Road, Camarillo, CA 93012 • (800) 242-0607 • FAX (805) 987-3385 BioSource Europe, S.A. • Rue de l'Industrie 8, B-1400 Nivelles, Belgium • +32 67 88 99 99 • FAX +32 67 88 99 96 This product is for research use only. Not for use in diagnostic procedures.

Intended Use and Materials Provided

The CytoSetTM for Mouse IL-13 contains components required to construct an enzyme-linked immunoassay for the specific and quantitative measurement of IL-13. Sufficient quantities of all reagents are provided to yield 10 plates of 96 wells if the recommended assay procedure and recommended storage and handling of materials are followed as specified on this insert. The materials provided are **FOR RESEARCH USE ONLY**.

Recommended Buffers and Solutions

The BioSource CytoSetTM Buffer Set (Cat. # CNB0011) containing Coating Buffers A and B, Assay Buffer, Substrate Solution (TMB), Stop Solution, and Wash Buffer is recommended.

1.	Coating Buffer A:	Coating Buffer A (Cat. # CB07100) from BioSource is recommended. Alternate buffer choice listed below.
		8.0 g NaCl, 1.13 g Na ₂ HPO ₄ , 0.2 g KH ₂ PO ₄ , 0.2 g KCl, 0.1% ProClin TM ; q.s. to 1.0 L with distilled H ₂ O, pH to 7.4.
2.	Coating Buffer B:	Coating Buffer B (Cat. # CB01100) from BioSource is recommended. Alternate buffer choice listed below.
		4.3 g NaHCO ₃ , 5.3 g Na ₂ CO ₃ , 0.1% ProClin TM , q.s. to 1.0 L with distilled H ₂ O, pH to 9.4.
3.	Assay Buffer:	Assay Buffer (Cat. # DS98200) from BioSource is recommended. Alternate buffer choice listed below.
	-	8.0 g NaCl, 1.13 g Na ₂ HPO ₄ , 0.2 g KH ₂ PO ₄ , 0.2 g KCl, 5.0 g bovine serum albumin (fraction V), 1 mL Tween 20 and 0.5% ProClin TM as a preservative; q.s. to 1.0 L with distilled H ₂ O, pH to 7.4.
		and 0.5% ProClin TM as a preservative; q.s. to 1.0 L with distilled H ₂ O, pH to 7.4.
4.	Wash Buffer:	Wash Buffer 25x (Cat. # WB01) from BioSource is recommended. Alternate buffer choice listed below.
		0.2 g KH ₂ PO ₄ , 1.9 g, K ₂ HPO ₄ , 3H ₂ O, 0.4 g EDTA, 0.5 mL Tween 20; q.s. to 1.0 L with distilled H ₂ O, pH to 7.4.
5.	Substrate Solution:	TMB (Cat. # SB01) from BioSource is recommended. Alternate solution choice listed below.
		Tetramethylbenzidine (TMB) and Hydrogen Peroxide.
6.	Stop Solution:	Stop Solution (Cat.# SS01100) from BioSource is recommended. Alternate solution choice listed below.
		$1.8 \text{ N H}_2 \text{SO}_4.$

Assay Optimization

CytoSets[™] from BioSource are designed to be very flexible for your experiments. Consequently, the assay procedure contains only recommendations. The assay procedure has been optimized for use with tissue culture samples. However, serum and plasma samples may be used but may require that certain assay parameters be modified. Investigators are advised to determine optimal buffer formulations, concentrations and incubation times for individual applications.

Recommended Assay Procedure

- 1. Prepare coating solution by diluting the coating antibody. See "coating antibody" section for the recommended coating antibody dilution.
- 2. Coat plates with 100 µL per well of the coating solution. Cover plates and incubate overnight (12-18 hr.) at 4°C.
- 3. Aspirate wells and wash 1 time with > 400 μL of Wash Buffer per well. Following wash, invert and tap on absorbent paper to remove excess liquid.
- 4. Block plate with $300 \,\mu\text{L}$ per well of Assay Buffer for 1 hour at room temperature.
- 5. Aspirate, invert, and tap on absorbent paper to remove excess liquid.
- 6. Prepare standards and sample dilutions in Assay Buffer (or in a diluent that most closely matches the matrix of your sample). For recommended dilutions and storage of the standard, see "standard" section.
- 7. Pipette $100 \,\mu\text{L}$ of standards (in duplicate) and samples into designated wells.
- 8. Immediately following step 7, add 50 μL of the working detection antibody into each well. For recommended dilutions, see "detection antibody" section. *Incubate for 2 hours at room temperature with continual shaking (700 rpm)*.
- 9. Aspirate and wash 5 times using the method in step 3.
- 10. Add 100 μL of the working streptavidin-HRP solution into each well. For recommended dilutions, see "streptavidin-HRP conjugate" section. *Incubate for 30 minutes at room temperature with continual shaking (700 rpm).*
- 11. Aspirate and wash 5 times using the method in step 3.
- 12. Add 100 µL of the TMB substrate to each well. Incubate plate for 30 minutes at room temperature with continual shaking (700 rpm).
- 13. Add 100 µL of Stop Solution to each well.
- 14. Measure absorbance at 450 nm (reference absorbance: 650 nm) within 30 minutes of adding Stop Solution. Calculate results using a loglog or 4-parameter curve fit.

Additional Materials Required

- 96 well NUNC MaxiSorp microplates; NUNC Cat. # 434797.
- Pipettes, shaker and timer.
- Microplate reader with a detector that can measure absorbance at 450 nm.
- 1 L graduated cylinder; plate washer or wash bottle.
- Polypropylene tubes for standards and sample dilutions, if needed.

biosource.com

BioSource International, Inc., USA • 542 Flynn Road, Camarillo, CA 93012 • (800) 242-0607 • FAX (805) 987-3385 BioSource Europe, S.A. • Rue de l'Industrie 8, B-1400 Nivelles, Belgium • +32 67 88 99 99 • FAX +32 67 88 99 96 This product is for research use only. Not for use in diagnostic procedures.