CERTEST *SARS-CoV-2*


ONE STEP SARS-CoV-2 Ag CARD
TEST

CERTEST BIOTEC S.L.

The test should be used by professional trained in IVD devices

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ENGLISH

INTENDED USE

CerTest SARS-CoV-2 one step card test is a coloured chromatographic immunoassay for the qualitative detection of nucleoprotein antigen of SARS-CoV-2 from nasopharyngeal swab samples from patients suspected of COVID-19 infection.

CerTest SARS-CoV-2 card test offers a simple and a highly sensitive screening assay as an aid in the diagnosis of SARS-CoV-2 infection.

INTRODUCTION

An unseen pneumonia spread swiftly in Wuhan, China, in December 2019. Scientists sequenced and recognized a novel β -coronavirus.

The disease, caused by the Severe Acute Respiratory Syndrome Corona Virus 2 (SARS-CoV-2), is called as Coronavirus disease 2019 (COVID-19).

Coronaviruses (CoVs) are RNA, enveloped viruses, with a large genome (29.9 kb), amongst which β -CoVs and α -CoVs can infect mammals. Viruses belonging to the Coronaviridae family have the most abundant membrane protein (M) among other proteins, i.e. spike glycoprotein (S), nucleocapsid protein (N) and an envelope protein (E). Spike glycoprotein (S) is one of the targets of T cell response in the immune system. The S protein also promotes the binding of the virus envelope to the ACE2 receptor and the entry of the virus into the cell target. ACE2 receptors are present in cells of the arteries, veins, smooth muscles, small intestine, alveoli of the lungs, hair follicles, cardiac myofibroblasts, skin, brain, and kidney, thus SARS-CoV-2 could potentially infect these tissues.

Clinical forms of this disease comprise from mild to very severe symptoms: pneumonia, fever and respiratory symptoms are the most frequent. Other symptoms of the viral infection include a sore throat, aches, pains, and shortness of breath. In some cases, runny nose, nausea, and diarrhea also occur.

The outbreak of newly discovered coronavirus begins mainly via discharge from the nose or droplets of saliva, once an infected individual sneezes or coughs, and the estimated time for incubation is within 2 weeks. Due to the high infectious rate of SARS-CoV-2, detection of asymptomatic positive patients is probably one of the key points to controlling the outbreak.

TEST PRINCIPLE

CerTest SARS-CoV-2 is based on the principle of a qualitative immunochromatographic assay for the qualitative detection of SARS-CoV-2 from nasopharyngeal swab samples from patients suspected of COVID-19 infection.

The strip consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against SARS-CoV-2 and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with test label solution (mouse monoclonal antibodies anti-SARS-CoV-2) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-SARS-CoV-2 monoclonal antibodies-red polystyrene microspheres), which was previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-SARS-CoV-2 antibodies present on the membrane (test line) capture the coloured conjugate and the red line will be visible. This band is used to interpret the result.

If the sample is negative, there is no SARS-CoV-2 antigens presence and yet, the antigens may be present in a concentration lower than the detection limit value, for which the reaction will not take place with the red-coloured conjugate complex. The anti-SARS-CoV-2 antibodies present on the membrane (test line) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red line will not appear.

Whether the sample is positive or not, the mixture continues to move across the membrane to the immobilized specific antibodies placed in the control line. The anti-specific protein antibodies present on the membrane will capture control green-conjugate complex and the control line will always appear. The presence of this green line serves as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.

- All the specimens should be considered potentially hazardous (especially samples from patients suspected of SARS-CoV-2 infection) and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- Use proper infection control practices when collecting and handling specimens. These practices should include, but are not limited to, personal protective equipment (PPE), such as a laboratory coat, surgical or appropriate mask, or face shield, disposable gloves, and eye protection. Follow local or national regulations regarding collecting and handling specimens.
- The tests (swabs, vials, positive control swab and devices) should be discarded in a proper biohazard container after testing. These containers should be discarded in accordance with local or national laws or regulations.
- Clean up spills thoroughly using an appropriate disinfectant.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest SARS-CoV-2 test. Do not use any other commercial kit component.
- Sterile swabs provided in the kits should be only used for taking the nasopharyngeal sample collection. They cannot be reuse.
- Do not touch the head of the sterile swabs provided when opening their primary packaging to avoid contamination.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- All positive results should be processed following local or national laws and regulations.

SPECIMEN COLLECTION AND PREPARATION

Samples should be process as soon as possible after collection. Follow proper infection control practices. If this is not possible, the samples can be stored in the refrigerator (2-8 °C) for a maximum of 8 hours prior to testing. If the samples are preserved in validated transport media (VTM, UTM o Saline Buffer) could be preserved on it until 6 hours at room temperature or in the refrigerator (2-8°C).

Specimen preparation (see illustration):

Nasopharyngeal swab method:

1. Remove the sterile swab from its packing.
2. Collect specimen with a sterile swab from one nostril.
3. Insert the swab into the nostril to the nasopharynx, rotating against the wall.
(to ensure swab contains cells as well as mucus). Remove the swab from the nostril carefully.
4. Consider repeating procedure using other nostril only if the protocol of the professional taking the sample required do it.
5. Process the swab as soon as possible after collecting the specimen.


Follow the test procedure (samples).

Nasopharyngeal samples previously extracted by transport media:

Even while using direct nasopharyngeal swab is the preferred protocol, the device can be also used with samples from nasopharyngeal swabs previously diluted in transport media such as: VTM, UTM or Saline Buffer.

1. Use the minimum volume of transport media in order to avoid loss in sensitivity. Example: 1.0 mL. Using higher volumes will affect sensitivity of the system.
2. Dilute the extracted sample 1:1 in Reagent *Respiratory* (sample diluent) provided (0.5mL/vial).

Follow immediately with the point 3 of test procedure (samples).

Specimen collection

MATERIALS

MATERIALS PROVIDED

- CerTest SARS-CoV-2 card tests
- Sample extraction vials with diluent (sample and control diluent)
- Sterile swabs
- CerTest SARS-CoV-2 Positive Control
- Laboratory rack
- Instructions for use

MATERIALS REQUIRED BUT NO PROVIDED



- Disposable gloves (other personal protective equipment that will be considered necessary)
- Timer
- Vortex (optional)

TEST PROCEDURE (SAMPLES)

Follow proper infection control practices. Allow tests, samples to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

- Put the sample extraction vials in the laboratory rack (identify properly samples from the different patients).
- Open the blue cap (prefilled with sample diluent Reagent *Respiratory*) and immediately put the swab into the sample extraction vial (1). Mix and homogenize the solution by rotating the swab forcefully against the side of the extraction vial to allow the solution to mix for 15 seconds. Best results are obtained when it is mixed to extract as much liquid as possible from the swab avoiding splashes and aerosols. Rotate the swab against the side of the sample extraction vial as the swab is withdrawn to extract as much liquid as possible. Discard the swab. Close the extraction vial with his blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. The sample can be shake using a vortex (this is optional).
- Remove the CerTest SARS-CoV-2 card test from its sealed bag just before using it.
- Be sure that use a separate sample extraction vial and test for each sample. Open the white-translucent cap of the sample extraction vial (2) and dispense 3 drops (approx. 100-120µL) from the sample extraction vial into the circular window marked with the letter S (3). After adding the drops into the sample window close the white-translucent cap (4).
- Read the results at 10 minutes.** Do not read the test result later than 10 minutes.


The Sample diluent buffer (Reagent *Respiratory*) has SARS-CoV-2 inactivation capability: 99.68% after 1 minute and 99.98% after 10 minutes.

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(4) Close the white-translucent cap.

Wait 10 minutes

INTERPRETATION OF RESULTS (please refer to the illustration below)

	SARS-CoV-2		Interpretation of the results
1.	NEGATIVE		There is no SARS-CoV-2 antigen presence in the sample.
	GREEN		
2.	POSITIVE		There is SARS-CoV-2 antigen presence in the sample. The presence of these lines (control green line and test red line), independent of their intensity should be considered as a positive result.
	GREEN-RED		
3.	INVALID		Invalid result, we recommend repeating the assay using the same sample with another test.
	Any other result		

NEGATIVE: Only one GREEN line appears across the results window in the site marked with the letter C (control line).

POSITIVE: In addition to the GREEN control line, a RED line also appears in the site marked with the letter T (test line) in the results window.

The result would be considered as positive independent on the intensity of the RED test line.

INVALID: Total absence of the control coloured line (GREEN) regardless the appearance or not of the test line (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persists, discontinue using the test kit and contact your local distributor.

NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured band in the test line (T) in the results window will vary depending on the concentration of antigen presents in the specimen. However, neither the quantitative value nor the rate of increase in antigen can be determined by this qualitative test.

QUALITY CONTROL

The CerTest SARS-CoV-2 test presents two controls:

Internal control: internal procedural control is included in the test. A green line appearing in the control line (C) in the results window is an internal control, which confirms sufficient specimen volume and correct procedural technique.

External control: Positive control swab.

EXTERNAL CONTROL


The CerTest SARS-CoV-2 Positive Control is an external quality control for CerTest SARS-CoV-2 card test. Use this control swab to check that the extraction reagents and the test are working properly.

Store as packaged in the sealed pouch at 2-30°C. The positive control swab is stable through the expiration date printed on the sealed pouch. The positive control swab must remain in the sealed pouch until use. Do not freeze.

POSITIVE CONTROL PROCEDURE


Allow tests and positive control swabs to reach room temperature (15-30°C) prior to testing. Do not open the package until the performance of the assay.

1. Open the blue cap (prefilled with sample diluent Reagent *Respiratory*) and immediately put the **positive control swab** into the sample extraction vial (1)
2. Mix the solution by rotating the positive control swab forcefully against the side of the sample extraction vial at least 15 seconds. Extract as much liquid as possible from the swab rotating the swab against the side of the sample extraction vial as the swab is withdrawn (2). Discard the swab. Close the extraction vial with his blue cap, please be sure that both caps (blue and white-translucent caps) are properly closed. The positive control extraction can be shake using a vortex (this is optional).
3. Remove the CerTest SARS-CoV-2 card test from its sealed bag just before using it.
4. Use a separate sample extraction vial and test for each control. Open the white-translucent cap of the sample extraction vial (2) and dispense 3 drops (approx. 100-120µL) from the sample extraction vial into the circular window marked with the letter S (3). After adding the drops in the sample window, close the white-translucent cap.
5. Read the result at 10 minutes. Do not read the test result later than 10 minutes.

Positive controls could be tested once for each new test kit opened and as otherwise required by your laboratory's standard quality control procedures.

INTERPRETATION POSITIVE CONTROL PROCEDURE

SARS-CoV-2 Positive Control	Interpretation of the results
2. POSITIVE GREEN-RED	 <p>Control + SARS-CoV-2 POSITIVE: A GREEN line (control line (C)) and a RED line (test line (T)) appear across the results window during the test performance due to presence of SARS-CoV-2 antigens in the positive control swab.</p>
3. INVALID Any other result	 <p>INVALID: Total absence of the control coloured line (GREEN), regardless the appearance or not of the test line (RED), or only the control (C) coloured line (GREEN). Incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control line failure. Review the procedure and repeat the control with a new test and a new positive control swab.</p>

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RECOMMENDATIONS

Recommendations World Health Organization: Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays- 11 September 2020.

1. To optimize performance, testing should be conducted by trained operators in strict accordance with the test procedure and within the first 5-7 days following the onset of symptoms.
2. Where possible, all positive samples giving positive results should be transported to laboratories with NAAT (Nucleic Acid Amplification Test) capability for confirmatory testing.
3. This test could be used to screen at-risk individuals and rapidly isolate positive cases in NAAT- confirmed COVID-19 outbreaks.
4. To monitor trends in disease incidence in communities.
5. For early detection and isolation of positive cases in health facilities, where there is widespread community transmission.
6. A negative result cannot completely exclude an active COVID-19 infection, repeat testing or preferably confirmatory testing should be performed (NAAT) whenever possible, particularly in symptomatic patients.
7. Even CerTest SARS-CoV-2 test was not validated using samples from asymptomatic contacts of cases, asymptomatic cases have been demonstrated to have viral loads similar to symptomatic cases, so CerTest SARS-CoV-2 could detect as positive.

LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. CerTest SARS-CoV-2 should be used only with nasopharyngeal swab samples. The use of swab specimens taken from other sites or the use of other samples such as saliva, sputum or urine has not been established. The quality of the test depends on the quality of the sample; proper nasopharyngeal specimens must be obtained.
3. The intensity of test line may vary from very strong at high antigens concentration to faint when the antigens concentration is close to the detection limit value of the test.
4. Positive results determine the presence of SARS-CoV-2 antigens. A diagnosis of SARS-CoV-2 infection should only be made by a physician after all clinical and laboratory findings have been evaluated and must be based in the correlation of the results with further clinical observations.
5. Positive results do not rule out co-infections with other pathogens.
6. Only VTM, UTM and Saline buffer transport media has been validated with CerTest SARS-CoV-2 test. When using transport media, the test sensitivity can be reduced due to excessive dilution of sample.
7. A negative result is not meaningful because of it is possible the antigens concentration in the nasopharyngeal samples is lower than the detection limit value. If the symptoms or situation still persist, it is recommended that all negative results undergo confirmatory testing using other method and/or virus identification by cell culture, PCR or other technique.

EXPECTED VALUES

On January 30, 2020, the World Health Organization declared COVID-19 a public health emergency of international concern. Since the onset of the pandemic, SARS-CoV-2 has affected millions of people, causing thousands of victims worldwide.

Recent studies have shown that the primary viral infection found in symptomatic patients was identical to that seen in asymptomatic patients, suggesting that COVID-19 patients with asymptomatic infections are capable of transmitting SARS-CoV-2. During the incubation period, asymptomatic patients are also likely to transmit the virus to healthy people.

Although SARS-CoV-2 infection can affect people of any age, children and the elderly with asthma are in a higher risk category for the disease. Those with underlying health problems such as cancer, chronic respiratory diseases, cardiovascular diseases, and diabetes are also reported to have a higher chance of developing serious illnesses.

PERFORMANCE CHARACTERISTICS

Analytical sensitivity (detection limit)

Detection limit value (typical) of CerTest SARS-CoV-2 is 1.0 ng/mL of SARS-CoV-2 recombinant protein or 1×10^3 TCID₅₀/mL of 2019nCoV/USA-WA1/2020.

Clinical sensitivity and specificity

An external evaluation (Hospital Universitario Príncipe de Asturias, Madrid (Spain)), with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2, Certest) vs qPCR technique (VIASURE). The results were as follows:

		qPCR technique		
		+	-	Total
IC test: CerTest SARS-CoV-2.	+	79	0	79
	-	5	150	155
	Total	84	150	234

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	94.0%	86.7–98.0%
Specificity	100.0%	97.6–100.0%
PPV	100.0%	95.4–100.0%
NPV	96.8%	92.6–98.9%

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)* from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

An external evaluation (Laboratoire de Virologie CHU Amiens-Picardie, Amiens (Francia)), with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2, Certest) vs qPCR technique. The results were as follows:

		qPCR technique		
		+	-	Total
IC test: CerTest SARS-CoV-2.	+	88	0	88
	-	8	62	70
	Total	96	62	158

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	91.7%	84.2–96.3%
Specificity	100.0%	94.2–100.0%
PPV	100.0%	95.9–100.0%
NPV	88.6%	78.7–94.9%

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)* from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

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An evaluation, with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2, Certest) vs qPCR technique (VIASURE). The results were as follows:

		qPCR technique		
		+	-	Total
IC test: CerTest SARS-CoV-2.	+	26	1	27
	-	2	233	235
	Total	28	234	262

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	92.9%	76.5–99.1%
Specificity	99.6%	97.6–100.0%
PPV	96.3%	81.0–99.9%
NPV	99.1%	97.0–99.9%

(*) Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)* from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection).

An evaluation multi-center, with nasopharyngeal samples from people suspected of infection by SARS-CoV-2 virus, was performed comparing the results obtained by an immunochromatographic test (CerTest SARS-CoV-2, Certest) vs qPCR technique. The results were as follows:

		qPCR technique		
		+	-	Total
IC test: CerTest SARS-CoV-2.	+	93	1	94
	-	7	455	462
	Total	100	456	556

CerTest SARS-CoV-2 vs qPCR technique		
	Mean Value	95% confidence interval
Sensitivity (*)	93.0%	86.1–97.1%
Specificity	99.8%	98.8–100.0%
PPV	98.9%	94.2–100.0%
NPV	98.5%	96.9–99.4%

Taking into account the recommendations for *Antigen-detection in the diagnosis of SARS-CoV-2 infection using rapid immunoassays (11 September 2020)* from the WHO, the sensitivity of the test was calculated with nasopharyngeal samples with high viral load (high viral loads is expected in early symptomatic phases of the illness (with the first 5-7 days of illness) in the range of Ag-RDT test detection). (*) For this evaluation the positive samples correspond to patients with COVID-19 symptoms within seven days after onset of symptoms).

The results showed a high sensitivity and specificity to detect SARS-CoV-2 using CerTest SARS-CoV-2.

Hook Effect

CerTest SARS-CoV-2 test does not show inhibitory hook effect at concentration higher than 100,000xLoD.

Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest SARS-CoV-2; no cross reactivity against organism, pathogens, substances that could cause infections:

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CerTest SARS-CoV-2 showed some cross reaction with SARS and very low with MERS.

Interferences

An evaluation was performed to determine the possible interferences of CerTest SARS-CoV-2; no interferences against the substances tested were detected:

	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metronidazole	3.0 mg/mL	Negative	Loperamide hydrochloride (<i>Fortasec</i>)	0.15 mg/mL	Negative
	Ampicillin	3.0 mg/mL	Negative	Heparin (<i>Hibor</i>)	350.0 IU/mL	Negative
	Oseltamivir	3.0·10 ⁻² mg/mL	Negative	Almagato (<i>Almax</i>)	3.0 mg/mL	Negative
	Amantadine	0.3 mg/mL	Negative	Fosfamycin (<i>Monurol</i>)	3.0·10 ⁻³ mg/mL	Negative
	Ribavirin	3.0 mg/mL	Negative	Acetylcysteine (<i>Fluimucil</i>)	3.0·10 ⁻² mg/mL	Negative
	Codeine (<i>Toseina</i>)	0.2 mg/mL	Negative	Dexketoprofen trometamol (<i>Enantyum</i>)	0.3 mg/mL	Negative
	Benzocaine (<i>Angileptol</i>)	3.0·10 ⁻² mg/mL	Negative	Levofloxacin	3.0·10 ⁻² mg/mL	Negative
	Cloperastine (<i>Flutox</i>)	0.3 mg/mL	Negative	Ciprofloxacin	0.3 mg/mL	Negative
	Carbocisteine (<i>Iniston mucolítico</i>)	3.0·10 ⁻² mg/mL	Negative	Rifampicin (<i>Rifaldin</i>)	0.3 mg/mL	Negative
	Loratadine	0.3 mg/mL	Negative	Phenoxyethylpenicillin potassium	3.0 mg/mL	Negative
	Dexchlorpheniramine (<i>Polaramine</i>)	0.3 mg/mL	Negative	Ambroxol hydrochloride (<i>Mucosan</i>)	0.3 mg/mL	Negative
	Ebastine (<i>Ebastel</i>)	3.0 mg/mL	Negative	Macrogol 3350 (<i>Movicol</i>)	3.0 mg/mL	Negative
	Acetyl Salicylic (<i>Adiro</i>)	0.3 mg/mL	Negative	Lysine Carbocysteinate (<i>Pectox</i>)	3.0·10 ⁻² mg/mL	Negative
	Ibuprofen (<i>Espidifen</i>)	0.3 mg/mL	Negative	Hydroxyzine dihydrochloride	0.3 mg/mL	Negative
	Paracetamol (<i>Dolocatil</i>)	5.0 mg/mL	Negative	Lorazepam	3.0·10 ⁻³ mg/mL	Negative

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	Potencial interfering substance	Concentration	Results	Potencial interfering substance	Concentration	Results
Exogenous	Metamizole (<i>Nolotil</i>)	5.0 mg/mL	Negative	Amoxicillin	3.0 mg/mL	Negative
	Prednisone	0.3 mg/mL	Negative	Mercaptopurine	0.3 mg/mL	Negative
	Omeprazole	2.0·10 ⁻³ mg/mL	Negative	Biotine	100.0 µg/mL	Negative
	Naso GEL	0.9 mg/mL	Negative	Sore Throat Phenol spray	0.51 mg/mL	Negative
	CVS Nasal Spray (<i>Cromolyn</i>)	4 mg/mL	Negative	Tobramycin	0.3 mg/mL	Negative
	Afrin (<i>Oxymetazoline</i>)	0.05 mg/mL	Negative	Mupirocin	0.025 mg/mL	Negative
	CVS Nasal Drops (<i>Phenylephrine</i>)	10 mg/mL	Negative	Fluticasone Propionate	0.05 mg/mL	Negative
	ZICAM	0.1 mg/mL	Negative	Chloraseptic (<i>Menthol/Benzocaine</i>)	1.5 mg/mL	Negative
	Homeopathic	DIL 1/10	Negative			
Endogenous	Human haemoglobin	5.0 µg/mL	Negative	Human lactoferrin	5.0 µg/mL	Negative
	Human transferrin	0.5 µg/mL	Negative	Mucine	5.0 mg/mL	Negative
	Human calprotectin	5.0 µg/mL	Negative	Human blood	50.0 mg/mL	Negative
Others	HAMAs (Human anti-mouse antibodies)	10 µg/mL	Negative			

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