

# GlycanTune™ A+ Total Feed

## Description

GlycanTune™ A+ Total Feed is an AGT™ dry-format, single-part product for use in fed batch culture of multiple cell lines. When used, this total feeding supplement can shift glycan profiles from heavily G0F to mostly G1F and G2F glycans. GlycanTune™ A+ Total Feed is an animal-origin free, chemically defined formulation that contains no proteins, hydrolysates, or components of unknown composition. The product provides an easy-to-reconstitute, pH-neutral solution that can be concentrated from the suggested 73.0 g/L (1X) up to 219.0 g/L (3X). As a result, key nutrients can be delivered in smaller volumes.

Product	Catalog No.	Amount	Storage	Shelf life
GlycanTune™ A+ Total Feed (all in AGT™ dry-format)	A29719-04	1 × 1 L	2°C to 8°C, Store dark and dry.	12 months
	A29719-05	1 × 10 L		
	A29719-01	1 × 100 L		

## Product use

Caution: For use as a raw material in further manufacturing applications.

## Safety information

Read the Safety Data Sheets (SDSs) and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

## Reconstitute GlycanTune™ A+ Total Feed

- Measure 80% of the final volume of deionized or distilled water at room temperature (15°C to 30°C).
- Add GlycanTune™ A+ Total Feed (73.0 g/L) to water. Mix until completely dissolved and for a minimum of 20 minutes.

**Note:** Solution will be cloudy and become clear over time.

- Using a calibrated vessel, dilute to final volume with deionized or distilled water. Mix for an additional 10 minutes.
- Measure pH and osmolality. Confirm final pH of 6.5–7.2. Final osmolality should be 490–520 mOsm/kg.
- Sterilize immediately by membrane filtration (positive pressure recommended).  
**Note:** Once product is filtered, use immediately or store at 2°C to 8°C for up to 30 days. Protect from light.

## Reconstitution options

Depending on culture needs, GlycanTune™ A+ Total Feed can be concentrated following these recommendations:

Concentration	Dissolution time
1X (73.0 g/L)	20 min
3X (219.0 g/L)	75 min

## Use

To maximize protein glycosylation, GlycanTune™ A+ Total Feed can be added as a complete stand-alone feed every 1–3 days or as a continuous feed. To target glycosylation profiles, initiate feeding with EfficientFeed™ A+ AGT™ Supplement, and then transition to feeding with GlycanTune™ A+ Total Feed (both feeds can be added every 1–3 days or as continuous feeds). The timing of the transition to GlycanTune™ A+ will determine the glycosylation profile. We recommend that you test at least three different transition time points (early, middle and late culture) to assist in choosing a transition point that will result in the desired glycan profile. Total feed volume can range from 10–45% (1X) to 3–15% (3X) of the starting culture volume depending on the cell line and process. For best results, we recommend that you initiate feeding only after reaching mid to late exponential phase.

## Feeding examples

Example feeding strategy to maximize glycosylation using GlycanTune™ A+ Total Feed for either 1X or 3X concentration (Feed volume in percent of starting culture volume):

GlycanTune™ A+ Total Feed (Concentration)	Culture days	
	0–2	3–14
1X	—	3.8%
3X	—	1.3%

Example feeding strategy to target a glycosylation profile using EfficientFeed™ A+ AGT™ Supplement, and then transitioning to GlycanTune™ A+ Total Feed both at 3X concentration (Feed volume in percent of starting culture volume):











Concentration	Culture days		
	0–2	3–7	8–14
3X EfficientFeed™ A+ AGT™ Supplement	—	1.3%	—
3X GlycanTune™ A+ Total Feed	—	—	1.3%

## Related products

Product	Catalog No.
CD OptiCHO™ Medium, (1X) liquid	12681
CD CHO Medium, (1X) liquid	10743
CD FortiCHO™ Medium, (1X) liquid	A11483
Dynamis™ Medium AGT™ dry format	A26175
EfficientFeed™ A+ AGT™ Supplement	A25023
EfficientFeed™ B+ AGT™ Supplement	A25030
EfficientFeed™ C+ AGT™ Supplement	A25031
Trypan Blue Solution	15250
Water, Distilled	15230

## Explanation of symbols and warnings

The symbols present on the product label are explained below:

				
Caution, consult accompanying documents	Temperature Limitation	Keep away from light	Use By:	Consult instructions for use
				
Batch Code	Catalog number	Manufacturer	Sterilized using aseptic processing techniques	Read Safety Data Sheet

## Limited product warranty

Life Technologies Corporation and/or its affiliate(s) warrant their products as set forth in the Life Technologies' General Terms and Conditions of Sale found on Life Technologies' website at [www.lifetechnologies.com/termsandconditions](http://www.lifetechnologies.com/termsandconditions). If you have any questions, please contact Life Technologies at [www.lifetechnologies.com/support](http://www.lifetechnologies.com/support).

## Limited Use Label License No. 517: Internal Research and Bioproduction Use

Notice to Purchaser: The purchase of this product conveys to the purchaser the limited, non-transferable right to use the purchased amount of the product to (a) perform internal research for the sole benefit of the purchaser; (b) manufacture protein (or other biological material) for resale; and (c) perform research or manufacturing services conducted by the purchaser on a fee for service or contract basis for or on behalf of third parties. However, the purchaser may transfer this product, its components, or materials made using this product to a third party (including contract research/manufacturing organizations), provided that each such third party agrees in writing to use such product, components, or materials solely on behalf of the purchaser, and such third party is restricted from further transferring any such product, components, or materials to any individual or entity other than the purchaser. No additional rights are granted. By purchasing this product, the purchaser agrees not to: (1) resell the product in any form; (2) use the product as a therapeutic agent or diagnostics test component; (3) reverse engineer the product or cause the product to be reverse engineered; or (4) use the product for purposes other than what is indicated in this Limited Use Label License. Life Technologies is not aware of Intellectual Property ("IP") that would be infringed by the manufacture or sale of its products. Customers are urged to perform an IP search and analysis specific to their manufacturing processes and target biologic products. Should that analysis require information about Life Technologies products, Life Technologies will provide to the customer's IP counsel the information necessary to analyze any relevant IP in a manner that protects Life Technologies proprietary information. The purchaser is responsible for obtaining all regulatory approvals necessary for any therapeutic or diagnostic use of the protein (or biological material) manufactured using this product. For information on obtaining additional rights, please contact [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com) or Out Licensing, Life Technologies, 5791 Van Allen Way, Carlsbad, California 92008.

The information in this guide is subject to change without notice.

**DISCLAIMER:** TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Corporate entity

Life Technologies | Carlsbad, CA 92008 USA | Toll Free in USA 1.800.955.6288

© 2016 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit [thermofisher.com/support](http://thermofisher.com/support) or email [techsupport@lifetech.com](mailto:techsupport@lifetech.com)  
[thermofisher.com](http://thermofisher.com)

**ThermoFisher**  
S C I E N T I F I C