

Product P/N 4331182, 4351372
 Insert P/N 4335626 REV D
 Printed in (United States)

For Research Use Only.
 Not for use in diagnostic procedures.

For safety guidelines, refer to the "Safety" section in the TaqMan® Gene Expression Assays Protocol, P/N 4333458_b. For all chemicals in bold type below, read the MSDS, and follow the handling instructions. Wear appropriate protective eyewear, clothing, and gloves.

Overview TaqMan® Gene Expression Assays consist of a 20X mix of unlabeled PCR primers and TaqMan® MGB probe (FAM™ dye-labeled). These assays are designed for the detection and quantitation of specific human genetic sequences in RNA samples converted to cDNA. Gene expression quantitation using TaqMan® Gene Expression Assays is performed as the second step in a two-step reverse transcription-polymerase chain reaction (RT-PCR) protocol on any ABI PRISM® Sequence Detection System Instrument. All TaqMan® Gene Expression Assays are optimized to work with either TaqMan® Universal PCR Master Mix, No AmpErase® UNG (P/N 4324018) or TaqMan® Universal PCR Master Mix (P/N 4304437) and with complementary DNA (cDNA). These products utilize the universal thermal cycling parameters described below in Table 2.

Procedure

To prepare the reaction components for a single 20µL reaction (384-well plate) or a single 50µL reaction (96-well plate) refer to Table 1 below for singleplex reactions.

Table 1. Singleplex PCR Reaction Mix using TaqMan® Universal PCR Master Mix, No AmpErase® UNG (P/N 4324018)

Reaction Component	Volume/Well (20µL volume reaction 1)	Volume/Well (50µL volume reaction 1)	Final Concentration
TaqMan® Universal PCR Master Mix, No AmpErase® UNG (2X) 2	10	25	1X
20X TaqMan® Gene Expression Assay Mix	1	2.5	1X
cDNA diluted in RNase-free water	9	22.5	-
Total	20	50	

1. If different reaction volumes are used, amounts should be adjusted accordingly.
2. Volumes should be the same if using TaqMan® Universal PCR Master Mix (2X) (P/N 4304437).

Table 2. Thermal Cycler Conditions

Thermal Cycler	Times and Temperatures			
	Initial Setup		Each of 40 Cycles	
			Denature	Anneal/Extend
Sequence Detection Systems (7900HT, 7700, 7000) and Real-Time PCR Systems (7300, 7500)	HOLD 3	HOLD	CYCLE	
	UNG activation 2 min 50°C	10 min 95°C	15 sec 95°C	1 min 60°C

3. The two-minute, 50°C step is required for optimal AmpErase® UNG activity when using TaqMan® Universal PCR Master Mix (P/N 4304437). This step is not needed when using the TaqMan® Universal PCR Master Mix, No AmpErase® UNG (P/N 4324018).

For further information on the plate set-up procedure and data analysis refer to the User's Manual for the appropriate Sequence Detection System Instruments (7900HT, 7700, 7000, 7300, 7500). Gene expression using TaqMan® Gene Expression Assays should be performed in separate wells (singleplex assay). We recommend that the endogenous control of choice be run in separate wells (singleplex) as this does not require any validation experiments. If performing multiplex* experiments, we recommend running multiplex and singleplex assays in parallel to confirm that the C_T values are not affected by multiplex PCR amplification. For additional information regarding relative quantitation of gene expression experiments refer to the ABI PRISM® 7700 Sequence Detection System User Bulletin #2 (P/N 4303859). *Multiplex PCR is the use of more than one primer pair in the same tube. Refer to the ABI PRISM® Sequence Detection System User Bulletin #5 (P/N 4306236) for information regarding multiplex reactions.

Storage

Store between -15°C and -20°C; minimize freeze thaw cycles. The TaqMan® Gene Expression Assay Mix may be diluted in TE (final concentration of TE should be 10mM Tris-HCl/1mM EDTA pH 8.0, use RNase-free water)

Table 3. Product Information

Product	Part Number	Contents	Volume / # of Rxns
TaqMan® Gene Expression Assays	4331182	20X mix of unlabeled PCR primers and TaqMan® MGB probe (FAM™ dye-labeled)	250ul / 250rxns
TaqMan® Gene Expression Assays	4351372	20X mix of unlabeled PCR primers and TaqMan® MGB probe (FAM™ dye-labeled)	360ul / 360rxns

For Research Use Only. Not for use in diagnostic procedures.

Notice to purchaser: Limited License

TaqMan® probes are covered by U.S. patents and foreign counterparts and patents pending owned by Applied Biosystems.

