

Water-Soluble Bolton-Hunter Reagent

27712

0476.3


Number

27712

Description

Water-Soluble Bolton-Hunter Reagent (Sulfo-SHPP), 100mg (sulfosuccinimidyl-3-[4-hydroxyphenyl]propionate)

Molecular Weight: 365.29



Storage: Upon receipt store at -20°C. Product is shipped at ambient temperature.

Introduction

The Thermo Scientific™ Water-Soluble Bolton-Hunter Reagent (Sulfo-SHPP) conjugates tyrosine-like residues to primary amines for increasing the yield of subsequent iodination. The reagent is ideal for iodinating proteins that do not contain tyrosine residues or contain tyrosine residues with limited accessibility to iodination. Labeling with Sulfo-SHPP also can help to preserve tyrosine residues that are essential for protein function or immunogenicity.

Water-Soluble Bolton-Hunter Reagent contains an *N*-hydroxysuccinimide (NHS) ester reactive group. NHS esters react efficiently with primary amino groups (-NH₂) in pH 7-9 buffers to form stable amide bonds. The side chain of lysine (K) residues and the N-terminus of each polypeptide contain a primary amine available for labeling. The sulfonated NHS ester of Sulfo-SHPP provides increased reagent solubility in comparison to the non-sulfonated form of Bolton-Hunter Reagent, which is no longer commercially available.

Procedure for Iodinating Purified Protein

The following protocol describes modification of a protein (avidin is the example) before iodination.²

Materials Required

- Modification Buffer: 200mM borate buffer, pH 9.0
- Avidin (Product No. 21121)
- Phosphate-buffered saline (PBS): 0.1M sodium phosphate, 150mM sodium chloride, pH 7.2 (Product No. 28372)

Method

1. Dissolve 190mg of avidin in 2mL of Modification Buffer.
2. Immediately before use dissolve 5mg of Water-Soluble Bolton-Hunter Reagent in 1mL of Modification Buffer. (Alternatively, add 0.5mg of the reagent directly to the avidin solution and then skip Step 3.)
3. Add 100μL of the Water-Soluble Bolton-Hunter Reagent solution to the avidin solution.
4. Incubate on ice for 3 hours with periodic mixing.
5. Remove the non-reacted Water-Soluble Bolton-Hunter Reagent by dialyzing against PBS or other suitable buffer.
6. Radiolabel with ¹²⁵I using Thermo Scientific™ Pierce™ Iodinating Reagents.

Protocol for Iodinating Cell Surface Proteins

The following general protocol uses Water-Soluble Bolton-Hunter Reagent (Thompson, *et al.*).

Materials

- DMSO
- Sodium phosphate buffer: 50mM sodium phosphate, pH 7.5
- Chloramine-T: just before use, prepare 5mg/mL chloramine-T in sodium phosphate buffer
- Hydroxyphenyl acetic acid: just before use, prepare 1mg/mL in water
- Sodium metabisulfite: just before use, prepare 12mg/mL in sodium phosphate buffer
- Borate buffer: 0.1M sodium borate, pH 8.3
- Cells: 10^7 - 10^9 cells/mL in Borate buffer.
- Na^{125}I
- Glycine buffer: 0.1M glycine in sodium borate buffer
- Phosphate-buffered saline (PBS): 100mM sodium phosphate, 150mM sodium chloride, pH 7.2 (Product No. 28372)

Method

1. Immediately before use, prepare 0.1mg/mL Water-soluble Bolton-Hunter reagent in DMSO. Place 1-20 μL of this solution into a microcentrifuge tube.
2. The following reagents must be added in the specified order: quickly add 1-20 μL of 100mCi/mL of Na^{125}I , 10 μL of chloramine-T solution, 100 μL of hydroxyphenyl acetic acid, and 10 μL of sodium metabisulfite.
3. Add 100 μCi of ^{125}I -labeled Bolton-Hunter Reagent to 1mL of cells and incubate for 30 minutes on ice.
4. Stop the reaction with 100 μL of 0.1M glycine buffer.
5. Incubate cells for 15 minutes on ice.
6. Wash cells three times with PBS. This step removes any non-conjugated ^{125}I -labeled Bolton-Hunter Reagent and labeled internal components released by autolysis.

Related Thermo Scientific Products

28372	BupH™ Phosphate Buffered Saline Packs, 40 each
69576	Slide-A-Lyzer™ MINI Dialysis Device Kit, 10K MWCO, 0.1mL
43233	Pierce Dextran Desalting Columns, 5 × 10mL columns
28665	Pierce Iodination Beads, 50 beads
28601	Pierce Iodination Tubes, 10 tubes
28600	Pierce Iodination Reagent, 1g

General References

- Bolton, A. E. and Hunter, W.M. (1973). The labelling of proteins to high specific radioactivities by conjugation to a ^{125}I -containing acylating agent. *Biochem J* **133**:529-39.
- Feitelson, M.A., *et al.* (1981). Tryptic peptide mapping of picomolar quantities of protein labelled with the Bolton-Hunter Reagent. *Anal Biochem* **116**:473-9.
- McCarthy, R.C., *et al.* (1988). Human prostatic acid phosphatase: purification, characterization, and optimization conditions for the radioimmunoassay. *Clin Chimica Acta* **132**:287-99.
- Mock, D.M. and Dubois, D.B. (1986). A sequential, solid-phase assay for biotin in physiological fluids that correlates with expected biotin status. *Anal Biochem* **153**:272-8.
- Raja, R.H., *et al.* (1984). Preparation of alkylamine and ^{125}I -radiolabelled derivatives of hyaluronic acid uniquely modified at the reducing end. *Anal Biochem* **139**:168-77.
- Roll, F.J., *et al.* (1979). A new iodinating method of iodinating collagens for use in radioimmunoassay. *Anal Biochem* **96**:489-99.
- Thirkell, D., *et al.* (1989). Serotype 8- and serocluster-specific surface-expressed antigens of *Ureaplasma urealyticum*. *Infect Immunity* **57**:1697-701.
- Thompson, J.A., *et al.* (1987). Selective radiolabeling of cell surface proteins to a high specific activity. *Biochemistry* **26**:743-50.

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2013 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.