

| | Package Contents | <p>Catalog Number: K9000-20</p> <table border="1"> <thead> <tr> <th></th> <th>Amount</th> </tr> </thead> <tbody> <tr> <td>FreeStyle™ CHO-S™ Cells</td> <td>1 mL</td> </tr> <tr> <td>FreeStyle™ MAX Reagent</td> <td>1 mL</td> </tr> <tr> <td>FreeStyle™ CHO Expression Medium</td> <td>1 Liter</td> </tr> <tr> <td>OptiPRO™ SFM</td> <td>100 mL</td> </tr> <tr> <td>pCMV SPORT-βgal</td> <td>25 µg</td> </tr> </tbody> </table> | | Amount | FreeStyle™ CHO-S™ Cells | 1 mL | FreeStyle™ MAX Reagent | 1 mL | FreeStyle™ CHO Expression Medium | 1 Liter | OptiPRO™ SFM | 100 mL | pCMV SPORT-βgal | 25 µg |
|----------------------------------|-----------------------------|--|--|--------|-------------------------|------|------------------------|------|----------------------------------|---------|--------------|--------|-----------------|-------|
| | Amount | | | | | | | | | | | | | |
| FreeStyle™ CHO-S™ Cells | 1 mL | | | | | | | | | | | | | |
| FreeStyle™ MAX Reagent | 1 mL | | | | | | | | | | | | | |
| FreeStyle™ CHO Expression Medium | 1 Liter | | | | | | | | | | | | | |
| OptiPRO™ SFM | 100 mL | | | | | | | | | | | | | |
| pCMV SPORT-βgal | 25 µg | | | | | | | | | | | | | |
| | Storage Conditions | <ul style="list-style-type: none"> Store cells in liquid nitrogen. Store reagent, and media at 4°C. Protect media from light. Store the control vector at -20°C. | | | | | | | | | | | | |
| | Required Materials | <ul style="list-style-type: none"> 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells Orbital shaker in temperature and CO₂ controlled incubator | | | | | | | | | | | | |
| | Timing | <p>Thawing and Recovery: 2–3 days Subculturing: Every 2–3 days Transfection: 1–7 days</p> | | | | | | | | | | | | |
| | Selection Guide | <p>Protein Expression Systems Go online to view related products.</p> | | | | | | | | | | | | |
| | Product Description | <ul style="list-style-type: none"> The FreeStyle™ MAX CHO Expression System facilitates large-scale transient transfection of Chinese Hamster Ovary (CHO) cells, in a defined, serum-free medium, for expression of proteins and virus. Transfection and expression experiments may be performed directly in FreeStyle™ CHO Expression Medium without the need for media change. The kit provides enough reagents to perform 25 transfections and one control transfection in a 30-mL volume, but larger volume transfections may be performed using simple scale-up of reagents. All reagents are completely animal origin-free, including the defined, serum-free medium, which may be imperative for regulatory requirements. | | | | | | | | | | | | |
| | Important Guidelines | <ul style="list-style-type: none"> General Cell Handling Preparing Media | | | | | | | | | | | | |
| | Online Resources | <p>Visit our product page for additional information and protocols. For support, visit www.thermofisher.com/support.</p> | | | | | | | | | | | | |

Protocol Outline

- Thaw cells.
- Subculture cells.
- Transfect cells and generate protein or virus.

FreeStyle™ MAX CHO Expression System Kit Characteristics

- Expression system based on CHO-S™ Cells for best compatibility with downstream bioproduction cell lines
- High protein yields in 2 to 7 days
- Scalable from multi-well plates to liter scale









FreeStyle™ MAX CHO Expression System Individual Components

The FreeStyle™ MAX CHO Expression System includes the following major components:

Click the next to each product to go to its specific protocol.

- FreeStyle™ CHO-S™ Cells:** This cell line is adapted to high density, serum-free suspension culture in FreeStyle™ CHO Expression Medium and is capable of producing high levels of recombinant protein.
- FreeStyle™ CHO Expression Medium:** This is a defined, serum-free medium formulated specifically to allow growth and large-scale transfection of suspension FreeStyle™ CHO-S™ Cells.
- FreeStyle™ MAX Reagent:** This transfection reagent provides high transfection efficiency in suspension FreeStyle™ CHO-S™ Cells.

Limited Product Warranty and Disclaimer Details

| | | |
|---|-----------------------------|---|
|  | Package Contents | Catalog Number R800-07 Size 1 vial <ul style="list-style-type: none"> One vial containing 1×10^7 cells |
|  | Storage Conditions | <ul style="list-style-type: none"> Store in liquid nitrogen. |
|  | Required Materials | <ul style="list-style-type: none"> FreeStyle™ CHO Expression Medium 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells Orbital shaker in temperature and CO₂ controlled incubator Reagents and equipment to determine cell viability (e.g., hemocytometer with trypan blue or cell counter) |
|  | Timing | Thawing and Recovery: 2–3 days Subculturing: Every 2–3 days |
|  | Selection Guide | Protein Expression Systems Go online to view related products. |
|  | Product Description | <ul style="list-style-type: none"> FreeStyle™ CHO-S™ Cells are derived from the CHO cell line, and are adapted to suspension culture in FreeStyle™ CHO Expression Medium. Chinese Hamster Ovary (CHO) cells are among the most commonly used cell lines for transfection, expression, and large-scale production of recombinant proteins. |
|  | Important Guidelines | <ul style="list-style-type: none"> Subculture the FreeStyle™ CHO-S™ Cells a minimum of three times to allow them to recover from thawing before using them in transfection experiments. Keep cell densities between 1×10^6–3×10^6 cells/mL of culture for best performance. We recommend maintaining cells in a 125-mL or a 250-mL polycarbonate, disposable, sterile Erlenmeyer flask containing 25–40 mL or 50–80 mL total working volume of cell suspension, respectively. Glass flasks without baffles may be used, but clean them thoroughly after each use to avoid potential toxicity. |
|  | Online Resources | Visit our product page for additional information and protocols. For support, visit www.thermofisher/support . |

Protocol Outline

- A. Thaw cells.
- B. Passage cells every 2–3 days.

FreeStyle™ CHO-S™ Cells Protocol

 See page 2 to view a typical procedure for thawing and subculturing the cells.

FreeStyle™ CHO-S™ Cells Characteristics

Growth properties: Suspension

Doubling time: 18 hours. Doubling times may vary based on cell health, handling, and passage number.

Viability during log phase culture: >95%

Subculture conditions: Grow to 1×10^6 – 3×10^6 cells/mL. Passage by splitting back to 0.2×10^6 – 0.5×10^6 cells/mL (every 2–3 days). Discard cells when they reach passage number 30.

Scaling Up FreeStyle™ CHO-S™ Cell Culture

You can scale up the FreeStyle™ CHO-S™ cultures in spinner flasks or bioreactors. Determine the optimal spinner or impeller speed and seeding density for your culture system. We recommend that the cells be seeded at 0.5×10^6 viable cells/mL. Optimum spinner speed is approximately 100–130 rpm, and optimum impeller speed in Celligen™ stirred tank bioreactors is 70–100 rpm.






If the split ratio of cells to fresh media is less than 1:2, centrifuge the cell suspension and re-suspend the cell pellet in fresh medium before inoculating the culture.









 **Cryopreserving FreeStyle™ CHO-S™ Cells**

 **Limited Product Warranty and Disclaimer Details**

Thawing and Passaging FreeStyle™ CHO-S™ Cells in FreeStyle™ CHO Expression Medium

Follow the procedure below to recover and subculture FreeStyle™ CHO-S™ Cells.


| | Timeline | Steps | Procedure details | | |
|----------|---|--|--|---|---|
| Day 1 | 1  | Thaw cells | Rapidly thaw the cells in a water bath, decontaminate the vial using 70% ethanol, and open the cryovial in a class II biological cabinet. | | |
| | 2  | Add cells to medium | Add cells to 29 mL of pre-warmed medium in a 125-mL shake flask. | | |
| | 3  | Count cells and determine viability | Within 1–2 hours post-thaw, count cells and determine viability. Use hemocytometer and trypan blue exclusion method or automated cell counter. Cell density should be approximately 0.3×10^6 cells/mL and cell viability >95%. | | |
| | 4  | Incubate | Temperature 37°C | Humidified Atmosphere 8% CO ₂ in air | Orbital Shaker Platform 125 rpm |
| Days 3–4 | 5  | Subculture cells | <p>First passage: When cell density reaches $>1 \times 10^6$ cells/mL at $\geq 90\%$ viability (typically 2–3 days post-thaw), split the culture to 0.2×10^6–0.5×10^6 cells/mL in FreeStyle™ CHO Expression Medium.</p> <p>Subsequent passages: Every 2–3 days, cells should reach 1×10^6–3×10^6 cells/mL. Split to 0.2×10^6–0.5×10^6 cells/mL. Do not grow above 3×10^6 cells/mL. We recommend using a 125-mL or a 250-mL flask containing 25–40 mL or 50–80 mL of medium, respectively.</p> | | |

|  | Package Contents | <table border="1"> <thead> <tr> <th>Catalog Number</th> <th>Size</th> </tr> </thead> <tbody> <tr> <td>12651-014</td> <td>1000 mL</td> </tr> <tr> <td>12651-022</td> <td>6 × 1000 mL</td> </tr> </tbody> </table> | Catalog Number | Size | 12651-014 | 1000 mL | 12651-022 | 6 × 1000 mL |
|--|-----------------------------|---|----------------|------|-----------|---------|-----------|-------------|
| Catalog Number | Size | | | | | | | |
| 12651-014 | 1000 mL | | | | | | | |
| 12651-022 | 6 × 1000 mL | | | | | | | |
|  | Storage Conditions | <ul style="list-style-type: none"> Store at 4°C for a 12-month shelf life. | | | | | | |
|  | Required Materials | <ul style="list-style-type: none"> FreeStyle™ CHO-S™ Cells or other CHO cells 125-mL polycarbonate, disposable, sterile, vent-cap Erlenmeyer shaker flask or other appropriate vessel for culturing suspension cells Orbital shaker in temperature and CO₂ controlled incubator L-Glutamine-200 mM | | | | | | |
|  | Timing | <p>Thawing and Recovery: 2–3 days Subculturing: Every 2–3 days</p> | | | | | | |
|  | Selection Guide | <p>Protein Expression Systems Go online to view related products.</p> | | | | | | |
|  | Product Description | <ul style="list-style-type: none"> FreeStyle™ CHO Expression medium is a serum-free, protein-free, chemically-defined medium for the growth of Chinese Hamster Ovary (CHO) cells and expression of recombinant proteins in suspension culture. | | | | | | |
|  | Important Guidelines | <ul style="list-style-type: none"> FreeStyle™ CHO Expression Medium requires supplementation with L-glutamine. Aseptically add 8 mM to the medium before use. Subculture FreeStyle™ CHO-S™ Cells when they reach a density of approximately 1 to 3 × 10⁶ viable cells/mL, typically every 2–3 days. Split the culture to between 0.2 and 0.5 × 10⁶ cells/mL. Keep cell densities between 1–3 × 10⁶ cells/mL of culture for best performance. Do not add anti-clumping agent to the culture prior to transfection. It can be added post-transfection. | | | | | | |
|  | Online Resources | <p>Visit our product page for additional information and protocols. For support, visit www.thermofisher.com/support.</p> | | | | | | |

Protocol Outline

- Thaw cells.
- Passage cells every 2–3 days.

FreeStyle™ CHO-S™ Cell Culturing

 See page 2 to view a typical procedure for thawing and culturing CHO Cells.

Scaling Up FreeStyle™ CHO-S™ Cell Culture

You can scale up FreeStyle™ CHO-S™ cultures in spinner flasks or bioreactors. Determine the optimal spinner or impeller speed and seeding density for your culture system. We recommend that the cells be seeded at 0.2 to 0.5 × 10⁶ viable cells/mL. Optimum spinner speed is approximately 100–130 rpm, and optimum impeller speed in Celligen™ stirred tank bioreactors is 70–100 rpm.

If the split ratio of cells to fresh media is less than 1:2, centrifuge the cell suspension and resuspend in fresh, prewarmed FreeStyle™ CHO Expression Medium before inoculating the culture.

At high stirring speeds (i.e. >130 rpm) and/or depending on the impeller design, you may need to supplement the FreeStyle™ CHO Expression Medium with additional Pluronic™ F-68 (2.5–5 mL/L of 10% Pluronic™ F-68) to avoid shear stress in the culture.






Adapting Other CHO Cells to FreeStyle™ CHO Expression Medium









Cryopreserving FreeStyle™ CHO-S™ Cells

Limited Product Warranty and Disclaimer Details

Thawing and Culturing FreeStyle™ CHO-S™ Cells in FreeStyle™ CHO Expression Medium

Follow the procedure below to recover and passage CHO Cells in FreeStyle™ CHO Expression Medium.

| | Timeline | Steps | Procedure Details | | |
|----------|---|--|--|---|---|
| Day 1 | 1  | Thaw cells | Rapidly thaw the cells in a water bath, decontaminate the vial using 70% ethanol, and open the cryovial in a class II biological cabinet. | | |
| | 2  | Add cells to medium | Add cells to 29 mL of pre-warmed medium in 125-mL shake flask. | | |
| | 3  | Count cells and determine viability | Within 1–2 hours post-thaw, count cells and determine viability. Use hemocytometer and trypan blue exclusion method or automated cell counter. Cell density should be approximately 0.3×10^6 cells/mL and cell viability >90%. | | |
| | 4  | Incubate | Temperature 37°C | Humidified Atmosphere 8% CO ₂ in air | Orbital Shaker Platform 125 rpm |
| Days 3–4 | 5  | Subculture cells | <p>First passage: When cell density reaches $>1 \times 10^6$ cells/mL at $\geq 90\%$ viability (typically 2–3 days post-thaw), split cells to $0.2\text{--}0.5 \times 10^6$ cells/mL in FreeStyle™ CHO Expression Medium.</p> <p>Subsequent passages: Every 2–3 days, cells should reach $1\text{--}3 \times 10^6$. Split to $0.2\text{--}0.5 \times 10^6$ cells/mL. Do not grow above 3×10^6 cells/mL.</p> <p>We recommend using a 125- or 250-mL flask containing 30 or 60 mL of medium, respectively.</p> | | |

|  Package Contents | <table border="1"> <thead> <tr> <th>Catalog Number</th> <th>Size</th> </tr> </thead> <tbody> <tr> <td>16447-100</td> <td>1.0 mL</td> </tr> <tr> <td>16447-500</td> <td>15.0 mL</td> </tr> <tr> <td>16447-750</td> <td>10 × 15.0 mL</td> </tr> </tbody> </table> | Catalog Number | Size | 16447-100 | 1.0 mL | 16447-500 | 15.0 mL | 16447-750 | 10 × 15.0 mL |
|---|---|----------------|------|-----------|--------|-----------|---------|-----------|--------------|
| Catalog Number | Size | | | | | | | | |
| 16447-100 | 1.0 mL | | | | | | | | |
| 16447-500 | 15.0 mL | | | | | | | | |
| 16447-750 | 10 × 15.0 mL | | | | | | | | |
|  Storage Conditions | <ul style="list-style-type: none"> Store at 4°C. Do not freeze. | | | | | | | | |
|  Required Materials | <ul style="list-style-type: none"> FreeStyle™ 293-F Cells, FreeStyle™ CHO-S™ Cells, or DG44 Cells FreeStyle™ 293 Expression Medium, FreeStyle™ CHO Expression Medium, or DG44 Medium Erlenmeyer flasks with vented caps Orbital shaker in temperature and CO₂ controlled incubator Plasmid DNA OptiPRO™ SFM | | | | | | | | |
|  Timing | <p>Cell Preparation: 1 day Transfection: 10–20 minutes</p> | | | | | | | | |
|  Selection Guide | <p>Protein Expression Systems Go online to view related products.</p> | | | | | | | | |
|  Product Description | <ul style="list-style-type: none"> FreeStyle™ MAX Reagent is a proprietary, animal origin-free formulation for transfecting plasmid DNA into eukaryotic cells, which can be easily scaled up to produce large amounts of recombinant proteins. This transfection reagent is formulated specifically for use with FreeStyle™ 293-F, FreeStyle™ CHO-S™, and DG44 cells. | | | | | | | | |
|  Important Guidelines | <ul style="list-style-type: none"> DNA-FreeStyle™ MAX complexes must be made in OptiPRO™ SFM and can be added directly to cells in culture medium. Cultivate FreeStyle™ 293-F and FreeStyle™ CHO-S™ Cells, or DG44 Cells, in a humidified, 37°C, 8% CO₂ environment in suspension on an orbital shaker. | | | | | | | | |
|  Online Resources | <p>Visit our product page for additional information and protocols. For support, visit www.thermofisher.com/support.</p> | | | | | | | | |

Protocol Outline

- Culture cells at least three passages after thawing.
- Prepare and add DNA-lipid complexes to cells.
- Incubate cells for 1–7 days.
- Harvest.

Transfection Protocol

- See page 2 to view a typical procedure for transfecting FreeStyle™ 293-F and FreeStyle™ CHO-S™ Cells for protein expression.
- See page 3 to view a typical procedure for transfecting DG44 cells to generate stable cell lines.

Transfection Conditions for FreeStyle™ Cells

Final transfection volume: 30 mL

Number of cells to transfect: 3×10^7

Amount of plasmid DNA: 37.5 µg








Amount of FreeStyle™ MAX Reagent: 37.5 µL

Scaling Up or Down Transfections

Limited Product Warranty and Disclaimer Details









Transfecting FreeStyle™ 293-F or FreeStyle™ CHO-S™ Cells

Use the following protocol to transfect suspension cells. All amounts are given on a per-flask basis for 30-mL cultures in 125-mL shake flasks.

| Timeline | | Steps | Procedure Details | | | | | | | | |
|-------------|---------------------------|---|--|--|--|-------------|-----------------------|-------------------------|------|---------------------------|---------|
| Day -1 | 1 |  | Expand cells | | | | | | | | |
| | 2 |  | Count cells and determine viability | | | | | | | | |
| | 3 |  | Seed cells in flask | | | | | | | | |
| Day 0 | 4 |  | Prepare DNA-lipid complexes | | | | | | | | |
| | 5 |  | Add DNA-lipid complex to cells | | | | | | | | |
| | 6 |  | Incubate | | | | | | | | |
| Days 1-7 | 7 |  | Harvest cells or media | | | | | | | | |
| | | | <p>For each 30-mL transfection, you will need 3×10^7 cells in 30 mL of FreeStyle™ 293 Expression Medium or FreeStyle™ CHO Expression Medium.</p> <p>For FreeStyle™ 293-F Cells: One day prior to transfection, passage at $6-7 \times 10^5$ cells/mL; shake at 120–135 rpm.</p> <p>For FreeStyle™ CHO-S™ Cells: One day prior to transfection, passage at $5-6 \times 10^5$ cells/mL; shake at 120–135 rpm.</p> <p>Use the trypan blue dye exclusion method to determine cell viability and clumping in a small aliquot of cells. Use an automated cell counter or a hemocytometer to determine cell counts. On the day of transfection, your cells should have a density of $1.2-1.5 \times 10^6$ cells/mL at >95% viability.</p> <p>Dilute cells to 1×10^6 cells/mL. You will need 3×10^7 cells for each 30-mL transfection.</p> <p>Use fresh, pre-warmed FreeStyle™ 293 Expression Medium or FreeStyle™ CHO Expression Medium to a total volume of 30 mL for each 30-mL transfection.</p> <p>Prepare DNA-lipid complexes as follows:</p> <ol style="list-style-type: none"> Dilute 37.5 µg of plasmid DNA in OptiPRO™ SFM reduced serum medium to a total volume of 0.6 mL. Mix gently. Dilute 37.5 µL of FreeStyle™ MAX Reagent in OptiPRO™ SFM reduced serum medium to a total volume of 0.6 mL. Mix gently and incubate for 5 minutes at room temperature. Incubation times longer than five minutes may result in decreased activity. After the 5-minute incubation, add the diluted DNA to the diluted reagent to obtain a total volume of 1.3 mL. Mix gently. Incubate for 20–30 minutes at room temperature to allow the DNA-lipid complexes to form. <p>Add 1.2 mL of complex to each cell suspension flask. Each flask should have a total volume of 30 mL, and contain approximately 1×10^6 viable cells/mL.</p> <p>To the negative control flask, add 2 mL of reduced serum medium instead of complex.</p> <table border="1"> <thead> <tr> <th>Temperature</th> <th>Humidified Atmosphere</th> <th>Orbital Shaker Platform</th> </tr> </thead> <tbody> <tr> <td>37°C</td> <td>8% CO₂ in air</td> <td>125 rpm</td> </tr> </tbody> </table> <p>Assay for recombinant protein expression. Perform this step 1–7 days post-transfection. Harvest media instead of cells if recombinant protein is secreted.</p> | | | Temperature | Humidified Atmosphere | Orbital Shaker Platform | 37°C | 8% CO ₂ in air | 125 rpm |
| Temperature | Humidified Atmosphere | Orbital Shaker Platform | | | | | | | | | |
| 37°C | 8% CO ₂ in air | 125 rpm | | | | | | | | | |

Transfecting DG44 Cells to Generate Stable Cell Lines

Use this procedure to transfect linearized DNA into DG44 cells. All amounts are on a per-flask basis for 30-mL cultures in 125-mL shake flasks.

| Timeline | | Steps | Procedure Details | | |
|----------|---|------------------------------------|---|---|---|
| Day 0 | 1  | Prepare and culture the DG44 cells | a. Passage the cells at 3×10^5 cell/mL. b. Shake at 130–135 rpm at 37°C, 8% CO ₂ . c. Culture in CD DG44 Medium (Cat. No. 12610-010) with 8 mM L-glutamine (Cat. No. 25030-081) and 18 mL/L of 10% Pluronic™ F-68 (Cat. No. 24040-032). | | |
| Day 1 | 2  | Passage the DG44 cells again | Passage cells again at 3×10^5 cell/mL. | | |
| | 3  | Prepare the cells | Count the cells. Cell viability should be >95%. In each flask, add 1.5×10^7 cells in a total volume of 30 mL CD DG44 Medium. | | |
| | 4  | Combine lipid and linearized DNA | Gently invert the tube to mix the reagent. Then, add 18 µg of linearized DNA and 15 µg of FreeStyle™ MAX Reagent into 1.2 mL of OptiPRO™ SFM (at room temperature), and gently invert to mix. | | |
| Day 2 | 5  | Incubate the DNA-lipid mixture | Incubate for 10 minutes at room temperature, but no longer than 20 minutes. | | |
| | 6  | Add DNA-lipid mixture to cells | Slowly add 1.2 mL of mixture into the 125-mL flask containing the cells while slowly swirling the flask. | | |
| | 7  | Incubate | Temperature 37°C | Humidified Atmosphere 8% CO ₂ in air | Orbital Shaker Platform 130–135 rpm |
| Day 4 | 8  | Place cells on a selective medium | Place cells on a selective medium (for example, CD OptiCHO™ Medium, Cat. No. 12681-011). | | |