INSTRUCTIONS

Binding and Elution Buffers



0411.6

Number	Description
21001	Protein A IgG Binding Buffer, 1L, pH 8.0; contains EDTA as a preservative
21007	Protein A IgG Binding Buffer, 3.75L, pH 8.0; contains EDTA as a preservative
21011	Protein G IgG Binding Buffer, 3.75L, pH 5.0; with 0.02% sodium azide
21019	Protein G IgG Binding Buffer, 1L, pH 5.0; with 0.02% sodium azide
54200	Protein A/G IgG Binding Buffer, 240mL pH 8.0; contains EDTA as a preservative
21004	IgG Elution Buffer, 1L, pH 2.8
21009	IgG Elution Buffer, 3.75L, pH 2.8
21028	IgG Elution Buffer, 1L, pH 2.0
	Note: The binding buffers are processed through a 0.2µm filter. The elution buffer contains primary amines and is purged with nitrogen to exclude oxygen.
	Storage: Upon receipt store buffers at 4°C. Buffers are shipped at ambient temperature.

Introduction

The Thermo Scientific Binding and Elution Buffers are optimized to provide the highest efficiency of IgG binding and recovery when used with various immobilized immunoglobulin-binding proteins.

Note: The Protein A IgG Binding Buffer precipitates in the presence of SDS.

Procedure for Purifying IgG

A. Additional Materials Required

• Immobilized Protein A, Protein G, Protein A/G, or Protein L packed into a gravity-flow column

Note: While the same elution buffer is effective for all of these immobilized proteins, a different binding buffer is required for optimal binding with each. For Protein L, use phosphate-buffered saline for binding (Product No. 28372).

- Neutralization Buffer: Prepare 1mL of high-ionic strength alkaline buffer such as 1M phosphate or 1M Tris, (pH 7.5-9)
- Thermo Scientific Slide-A-Lyzer Dialysis Cassette or Zeba Spin Desalting Columns (Product No. 43230) for buffer exchange

B. Procedure

- 1. Equilibrate the buffers and column to room temperature.
- 2. Dilute serum samples, ascites fluid, or tissue culture supernatant at least 1:1 with binding buffer to ensure the ionic strength and pH are maintained for optimal binding.

Note: If plasma is used, the sample may appear hazy after adding the binding buffer caused by lipoprotein precipitation. For optimal antibody recoveries, centrifuge the diluted sample at $10,000 \times g$ for 20 minutes and apply the supernatant to the equilibrated column.

- 3. Equilibrate the column with five resin-bed volumes of binding buffer.
- 4. Apply the diluted sample to the column and allow it to flow completely into the resin.
- 5. Wash the column with 5-10 resin-bed volumes of binding buffer.



- 6. Elute antibodies with five resin-bed volumes of elution buffer and collect 0.5-1mL fractions. Immediately adjust eluted fractions to physiologic pH by adding 100μL of Neutralization Buffer to 1mL of eluate. Monitor elution by measuring the absorbance at 280nm or by protein assay such as Thermo Scientific BCA Protein Assay Kit (Product No. 23225).
- 7. Pool the eluted IgG fractions that contain the highest absorbance values. Use the purified antibodies directly for SDS-PAGE or exchange the buffer to a system compatible with the specific downstream application.
- 8. Regenerate the affinity column by applying eight resin-bed volumes of Elution Buffer and allowing it to flow through.
- 9. For storage, wash column with five resin-bed volumes of water containing 0.02% sodium azide. When ~3mL of solution remains above the top disc, replace bottom cap and then the top cap on the column. Store column upright at 4°C.

Elution of Antibody/Antigen Interactions

The IgG Elution Buffer is an effective universal buffer for disrupting antigen:antibody interactions. The most effective elution conditions (pH, ionic strength, chaotrope, or denaturant) for a particular interaction, however, depends on the specific composition of ionic, hydrophobic and hydrogen bonds involved. All elution buffers cause various levels of antibody or antigen functionality loss, limiting the number of times an affinity support can be reused. For more information, see Tech Tip #27: Optimize elution conditions for immunoaffinity purification, which is available from the Pierce website.

Note: This example procedure assumes that an antigen has been immobilized to a beaded resin and is being used in a column format to purify antibody from serum or culture supernatant. Use the same procedure for the converse purification scheme (i.e., immobilized antibody column to purify an antigen) by simply switching all references to antigen and antibody.

- 1. Equilibrate buffers and column of immobilized antigen to the same temperature (room temperature or 4°C).
- 2. Prepare antibody sample for binding. Dilute concentrated samples such as serum and ascites fluid with an equal volume of a physiological buffer (i.e., PBS or TBS). Adjust cell culture supernatant and other dilute samples to pH 7.2-7.4 with TBS or PBS. Centrifuge cloudy samples and use only the clear supernatant.
- 3. Wash and prepare the antigen column by adding five resin-bed volumes of buffer and allowing it to flow through. Discard the flow-through storage buffer.
- 4. Add the prepared antibody sample to the antigen column and allow it to flow through. If desired, stop the column occasionally or control the flow rate to ensure adequate binding time. Save the flow-through non-bound sample.
- 5. Wash column by adding 5-10 resin-bed volumes of buffer and allowing it to flow through. If desired, analyze flow-through fractions to determine if washing is complete.
- 6. Elute the purified antibody by adding 5-10 resin-bed volumes of Elution Buffer and collecting the flow-through in several small fractions. Immediately adjust eluted fractions to physiologic pH by adding 100µL of Neutralization Buffer (i.e., 1M phosphate or 1M Tris; pH 7.5-9) per 1mL of eluate. Fractions confirmed to contain antibody may be pooled, and the other fractions may be discarded.
- 7. Dialyze or desalt to exchange the purified antibody fractions into a buffer suitable for storage and analysis.

Related Thermo Scientific Products

28372	BupH TM Phosphate Buffered Saline Packs, 40 packs
69576	Slide-A-Lyzer™ MINI Dialysis Device Kit, 10K MWCO, 0.1mL
66382	Slide-A-Lyzer Dialysis Cassette Kit, 10K MWCO, 3mL
89889	Zeba™ Spin Desalting Columns, 7K MWCO, 2mL, 5/pkg
89891	Zeba Spin Desalting Columns, 7K MWCO, 5mL, 5/pkg
44894	AminoLink [™] Plus Immobilization Kit
44999	SulfoLink TM Immobilization Kit for Peptides
44899	CarboxyLink TM Immobilization Kit

General References

Bigbee, W.L., *et al.* (1983). Monoclonal antibodies specific for the M- and N-forms of human glycophorin A. *Mol Immunol* **20**:1353-62. Ey, P.L., *et al.* (1978). Isolation of pure IgG₁, IgG_{2a}, and IgG_{2b} immunoglobulins from mouse serum using protein A-Sepharose. *Immunochem* **15**:429-36.



Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at <u>www.thermoscientific.com/pierce</u>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.