

## **Random Primers DNA Labeling System**

Cat. No.:	18187-013	Size:	30 Assays		
Lot No.:		Store at non-fros	-20°C in a st-free freezer	°C in a ree freezer	
COMPONEN	<u>rs</u> :		<u>PART NO.</u>	<u>AMOUNT</u>	
Random Prime 33 mM 2- oligodeox	ers Buffer Mixture: [0.67 M HEPES, 0.17 M Tris-HCl, 1 mercaptoethanol, 1.33 mg/ml BSA, 18 OD <sub>260</sub> units/ml yribonucleotide primers (hexamers), pH 6.8]	7 mM MgC	Cl <sub>2</sub> , Y01101	500 µl	
dATP Solution	e: 0.5 mM dATP in 1 mM Tris-HCl (pH 7.5)		Y01102	100 µl	
dCTP Solution	: 0.5 mM dCTP in 1 mM Tris-HCl (pH 7.5)		Y01103	100 µl	
dGTP Solution	: 0.5 mM dGTP in 1 mM Tris-HCl (pH 7.5)		Y01104	100 µl	
dTTP Solution	: 0.5 mM dTTP in 1 mM Tris-HCl (pH 7.5)		Y01105	100 µl	
<u>Klenow Fragm</u> Klenow fr 100 mM F	tent (Large Fragment of DNA Polymerase I): 3 U/μl agment in 50 mM potassium phosphate buffer (pH 7.0), KCl, 1 mM DTT, 50% (v/v) glycerol		Y01106	40 µl	
Stop Buffer: 0	0.5 M EDTA, pH 8.0		50690	500 µl	
<u>Control DNA</u> : 10 mM Tr	5 ng/μl pBR322 DNA/ <i>Rsa</i> I fragments in is-HCl (pH 7.4), 5 mM NaCl, 0.1 mM EDTA		50981	30 µl	
Distilled Wate	<u>r</u> :		50837	1.25 ml	

## QUALITY CONTROL:

When 25 ng of the control DNA provided in this system is labeled using the standard Random Primers labeling conditions, a specific activity  $\ge 1 \times 10^9$  cpm per µg of template DNA was obtained within 1 h at 25°C. Typical results are shown at right.



This product is distibuted for laboratory research use only. CAUTION: Not for diagnostic use. The safety and efficacy of this product in diagnostic or other clinical uses has not been established.

For technical questions about this product, call the Life Technologies TECH-LINE<sup>SM</sup> [U.S.A. (800) 828-6686].

## STANDARD LABELING PROTOCOL:

The following instructions are given for a standard assay using 25 ng of DNA.

- 1. Denature 25 ng of DNA dissolved in 5-20 µl of distilled water in a microcentrifuge tube by heating for 5 min in a boiling water bath, then immediately cool on ice.
- 2. Perform the following additions on ice:
  - $2\,\mu l\,\,\,dATP$  solution
  - 2 µl dGTP solution
  - 2 µl dTTP solution
  - 15 µl Random Primers Buffer Mixture
  - 5  $\mu$ l (approximately 50  $\mu$ Ci) [ $\alpha$ -<sup>32</sup>P]dCTP, 3000 Ci/mmol, 10  $\mu$ Ci/ $\mu$ l
  - Distilled water to a total volume of 49 µl

Mix briefly.

- 3. Add 1 µl Klenow Fragment, mix gently but thoroughly, and centrifuge briefly.
- 4. Incubate at 25°C for 1 h (see note 4).
- 5. Add 5 µl Stop Buffer.
- 6. Dilute a 2 μl-aliquot of the mixture with 498 μl of distilled water (provided by the user's laboratory) or with 498 μl of TE buffer [10 mM Tris-HCl, (pH 8), 1 mM Na<sub>2</sub> EDTA].
- Spot a 5-µl aliquot of this dilution on a glass fiber filter disk (Whatman<sup>®</sup> GF/C or equivalent). Wash the filter disk three times with 50 ml of ice cold 10% (w/v) TCA containing 1% (w/v) of sodium pyrophosphate, and once with 50 ml of 95% ethanol at room temperature.
- 8. Dry the filter under a heat lamp, and determine the precipitable radioactivity by liquid scintillation counting. A result of X cpm in this test corresponds to a total of 2750 × cpm incorporated in the entire incubation mixture.

For confirmation of <u>total radioactivity</u> in the incubation mixture, a second 5- $\mu$ l aliquot of the above dilution (step 6) can be spotted on another glass fiber filter disk, dried under the heat lamp (<u>without</u> intervening TCA wash), and counted in a liquid scintillation counter. A result of X cpm in this test corresponds to 2750 × cpm of <u>total radioactivity</u> in the entire reaction mixture.

NOTES:

- 1. The Random Primers DNA Labeling System contains the deoxynucleoside triphosphates in four separate solutions. This allows the system to be used with any of the four  $[\alpha^{-32}P]$ -labeled nucleoside triphosphates by suitably modifying the above protocol, which is given for the case of  $[\alpha^{-32}P]$ -labeled dCTP. The user may also elect to mix three of the nucleoside triphosphates in a cocktail in a 1:1:1 ratio, if a particular protocol is to be consistently used.
- 2. The Random Primers DNA Labeling System contains a qualified control DNA which may be used to monitor the performance of the system. The incorporation curve shown on the other side was obtained with this control DNA.
- 3. DNA probes labeled by this procedure can be used in blot hybridizations without removing the unincorporated nucleoside triphosphates (2). If desired, the labeled probe can be separated from unincorporated nucleotides by repeated ethanol precipitation (add 1/2 volume of 7.5 M ammonium acetate and 2 volumes of ethanol, repeat once.) or using the CONCERT<sup>TM</sup> PCR Purification System (Cat. Series 11458).
- 4. Incubation times longer than 1 h may give higher specific activity (see figure on other side).
- 5. For optimal results, supercoiled DNA should be linearized or alkali denatured prior to random promer labeling:

-Effect of Substrate DNA Configuration/Denaturation on Labeling:

Substrate	Sp.Act.
Linear pUC 19 (denatured by boiling)	$2.06 \times 10^9$ cpm/µg
Supercoiled pUC 19 (denatured by boiling)	$0.12 \times 10^9$ cpm/µg
Supercoiled pUC 19 (alkali dematured))	$1.30 \times 10^9$ cpm/µg

## **REFERENCES:**

- 1. Feinberg, A. P. and Vogelstein, B. (1983) Anal Biochem. 132, 6.
- 2. Feinberg, A. P. and Vogelstein, B. (1984) Anal Biochem. 137, 266.
- 3. Matathias, A.; and Komro, C. (1990) *Focus*<sup>®</sup> 12:2, 53.

Whatman<sup>®</sup> is a registered trademark of Whatman International, Ltd.