---

# CERTEST FOB+Transferrin + Calprotectin+ Lactoferrin

---

ONE STEP FOB, Transferrin,  
Calprotectin and Lactoferrin  
COMBO CARD TEST

---

CERTEST BIOTEC S.L.

---

For information purposes only

For information purposes only

# CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette

## ENGLISH

### INTENDED USE

CerTest FOB+Transferrin+Calprotectin+Lactoferrin one step combo card test is a coloured chromatographic immunoassay for the simultaneous semi-qualitative detection of human haemoglobin (hHb), human transferrin (hTf), human calprotectin (hCp) and human lactoferrin (hLf) in stool samples.

CerTest FOB+Transferrin+Calprotectin+Lactoferrin combo card test offers a simple, highly sensitive and non-invasive screening assay to make a presumptive diagnosis of gastrointestinal bleeding and inflammatory disorder caused by colorectal cancer.

### INTRODUCTION

Colorectal cancer is the second leading cause of illness and death in the Western world. The screening using fecal occult blood tests is based on the bleeding detection, as important symptom of colonic neoplasm, either at an early-stage cancer or large adenomatous polyps. Colorectal cancer is also associated with local acute inflammatory reaction being visualized, in some cases, by white cell neutrophil scanning.

Haemoglobin is the iron-containing oxygen-transport protein in the red blood cells of all vertebrates that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases.

Transferrin is a blood-derived component that may be leaked into gastrointestinal tract and then discharged with the faeces in gastrointestinal bleeding diseases. Transferrin is stable in faeces and a good marker to detect loss of blood from the upper and lower intestine (gastrointestinal bleeding).

When gastrointestinal blood is lost, the stool will contain a combination of intact or nearly intact haemoglobin, intact heme, and heme-derived porphyrins in amounts that depend on the site, the amount of bleeding and the transit time through the gut. Immunochemical tests detect intact or nearly intact human haemoglobin, being a very specific technique to detect the loss of blood from the lower intestine, because blood from lower sites is less degraded during transit. The haemoglobin is unstable in faeces, leading to false negative results. The detection of fecal transferrin, which is more stable than haemoglobin, provides an alternative way of diagnosing the disease in the upper digestive tract.

Calprotectin is a neutrophil cytosolic protein with antimicrobial properties, which is present at increased concentration in stool during bowel inflammation. The stability of the protein to degradation keeps it stable in faeces for up to 7 days at room temperature, making it an ideal analyte. Calprotectin is released by activation of leukocytes, giving increased levels in plasma, cerebral spinal fluid, synovial fluid, urine or stools as a consequence of disease in the relevant organ(s). Calprotectin inhibits zinc-dependent enzyme systems, killing microbes and inducing apoptosis in normal and cancer cells. In the presence of calcium, calprotectin is a remarkably resistant to proteolytic degradation, and being stable in stools kept at room temperature for 7 days.

Lactoferrin is a glycoprotein component of neutrophil secondary granules, a primary component of the acute inflammatory response released from fecal leukocytes. This protein is resistant to proteolysis in the faeces and may serve as a marker of inflammation in the intestine. The major cause of fecal neutrophils in patients with chronic diarrhea is chronic inflammatory bowel disease of the colon (i.e., Crohn's disease and ulcerative colitis). Lactoferrin has also been studied as a predictor of infection with invasive enteropathogens in children with diarrhea. Bacterial inflammatory diarrhea may be caused by *Shigella*, *Salmonella*, *Escherichia coli enteroinvasive (EIEC)* and *enteroaggregative (EAEC)*, *Campylobacter* and *Clostridium difficile*.

The hCp and hLf markers are useful for the differentiation between functional (for example: irritable bowel syndrome) and organic disease (for example inflammatory bowel disease).

Gastrointestinal inflammatory disorders also include Crohn's disease (CD) and ulcerative colitis (UC). Suspicion is raised in patients with persistent ( $\geq 4$  weeks) or recurrent ( $\geq 2$  episodes in six months) abdominal pain and diarrhea. Additionally, rectal bleeding, weight loss or anemia increase the probability of the condition. Endoscopic evaluation with histopathological sampling is generally considered indispensable to diagnosis patients with gastrointestinal inflammation disorder suspicion.

### TEST PRINCIPLE

CerTest FOB+Transferrin+Calprotectin+Lactoferrin is based on the principle of a semi-quantitative immunochromatographic assay for the determination of human haemoglobin, human transferrin, human calprotectin and human lactoferrin in stool samples.

**Strip A** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human haemoglobin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with label solution (mouse monoclonal antibodies anti-human haemoglobin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

**Strip B** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human transferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with label solution (mouse monoclonal antibodies anti-human transferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

## CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette

**Strip C** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human calprotectin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with label solution (mouse monoclonal antibodies anti-human calprotectin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

**Strip D** consists of a nitrocellulose membrane pre-coated with mouse monoclonal antibodies on the test line (T), in the results window, against human lactoferrin and with rabbit polyclonal antibodies, on the control line (C), against a specific protein. The label/sample absorbent pad is sprayed with label solution (mouse monoclonal antibodies anti-human lactoferrin) conjugated to red polystyrene latex and control label solution (specific binding protein) conjugated to green polystyrene latex, forming coloured conjugate complexes.

If the sample is human haemoglobin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human haemoglobin monoclonal antibodies-red polystyrene microspheres) in the strip A, if the sample is human transferrin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human transferrin monoclonal antibodies-red polystyrene microspheres) in strip B, if the sample is human calprotectin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human calprotectin monoclonal antibodies-red polystyrene microspheres) in strip C and if the sample is human lactoferrin positive, the antigens of the diluted sample react with the red-coloured conjugate complex (anti-human lactoferrin monoclonal antibodies-red polystyrene microspheres) in strip D, which were previously pre-dried on the absorbent pad. The mixture then moves upward on the membrane by capillary action. As the sample flows through the test membrane, the binding conjugate complexes migrate. The anti-human haemoglobin antibodies present on the membrane of strip A (test line), the anti-human transferrin antibodies present on the membrane of strip B (test line), the anti-human calprotectin antibodies present on the membrane of strip C (test line) and the anti-human lactoferrin antibodies present on the membrane of strip D (test line) capture the coloured conjugate and the red line will be visible in the strips. These bands are used to interpret the result.

If the sample is negative, there is no human haemoglobin, human transferrin, human calprotectin and human lactoferrin presence and yet, the antigens may be present in a concentration lower than the cut-off value, for which the reaction will not take place with any red-coloured conjugate complex. The anti-human haemoglobin, anti-human transferrin, anti-human calprotectin and anti-human lactoferrin antibodies present on the membranes (test lines) will not capture the antigen-red-coloured conjugate complex (not formed), for which the red lines will not appear.

Whether the sample is positive or not, in the four strips, the mixture continues to move across the membranes to the immobilized specific antibodies placed in the control lines. The anti-specific protein antibodies present on the four membranes will capture control green-conjugate complex and the four control lines will always appear. The presence of these green lines serve as: 1) verification that sufficient volume is added, 2) that proper flow is obtained and 3) an internal control for the reagents.

### STORAGE AND STABILITY

Store as packaged in the sealed pouch at 2-30°C. The test is stable until the expiration date marked on its sealed pouch. The test must remain in the sealed pouch until use. Do not freeze.

### PRECAUTIONS

- For professional *in vitro* diagnostic use.
- Do not use after expiration date.
- All the specimens should be considered potentially hazardous and handled in the same manner as an infectious agent. A new test must be used for each sample to avoid contamination errors.
- The tests should be discarded in a proper biohazard container after testing.
- Reagents contain preservatives. Avoid any contact with skin or mucous membrane. Consult safety data sheet, available on request.
- Components provide in the kit are approved for use in the CerTest FOB+Transferrin+Calprotectin+Lactoferrin combo card test. Do not use any other commercial kit component.
- The presence of yellow lines in the results window (control line zone and test line zone), before using the test, is completely normal and does not imply failure of the test functionality.
- Follow Good Laboratory Practices, wear protective clothing, use disposal gloves, goggles and mask. Do not eat, drink or smoke in the working area.


### SPECIMEN COLLECTION AND PREPARATION

Stool samples should be collected in clean containers. The samples can be stored in the refrigerator (2-8°C) for 7 days prior to testing. For longer storage, maximum 6 months, the specimen must be kept frozen at -20°C. In this case, the sample will be totally thawed and brought to room temperature before testing. Homogenise stool sample as thoroughly as possible prior to preparation.

## Specimen preparation (see illustration):

1. Take out the cap 1 of the stool collection tube and use the stick to pick up sufficient sample quantity. Introduce the stick once into 4 different parts of the stool sample (1), making sure that, in each insertion, the stick's screw is covered with the sample to take the appropriate amount of stool. Then, introduce the sample in the stool collection tube (2). For liquid samples, add approx. 20µL in the stool collection tube using a micropipette.
2. Close the tube with the diluent and stool sample (2). Shake the tube in order to assure good sample dispersion (3). The stool collection tube with dilute sample can be stored until 2 days at room temperature or until 5 days in the refrigerator (2-8°C) prior to testing.

**Precautions:** Patients should not collect samples during their menstrual period if they have bleeding haemorrhoids, blood in urine or if they have strained during bowel movement.

## MATERIALS

### MATERIALS PROVIDED

- CerTest FOB+Transferrin+Calprotectin+Lactoferrin combo card tests
- Instructions for use
- Stool collection tubes with diluent

### MATERIALS REQUIRED BUT NO PROVIDED

- Specimen collection container
- Disposable gloves
- Timer

## TEST PROCEDURE



Allow tests, stool samples and controls to reach room temperature (15-30°C) prior to testing. Do not open pouches until the performance of the assay.

1. Proceed to shake the stool collection tube in order to assure good sample dispersion.
2. Remove the CerTest FOB+Transferrin+Calprotectin+Lactoferrin combo card test from its sealed bag just before using it.
3. Take the stool collection tube, open the cap 2 (4) and, using each time the same tube, dispense 3 drops in the circular window marked with the letter A (5), 3 drops in the circular window marked with the letter B (6), 3 drops in the circular window marked with the letter C (7) and 3 drops in the circular window marked with the letter D (8). Avoid adding solid particles with the liquid.
4. **Read the results at 10 minutes.** Do not read the test result later than 10 minutes.

If the test does not run due to solid particles, stir the sample added in the sample window (S) with the stick. If it doesn't work, dispense a drop of extraction diluent until seeing the liquid running through the reaction zone.


# CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
 Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette


## FOB Strip A-procedure

Add 3 drops in A-circular window


## Transferrin Strip B-procedure

Add 3 drops in B-circular window


## Calprotectin Strip C-procedure

Add 3 drops in C-circular window


## Lactoferrin Strip D-procedure


Add 3 drops in D-circular window


## INTERPRETATION OF RESULTS (please refer to the illustration below)

A: Green→Negative hHb  
 B: Green→Negative hTf  
 C: Green→Negative hCp  
 D: Green→Negative hLf


A: Green/Red→Positive hHb  
 B: Green/Red→Positive hTf  
 C: Green/Red→Positive hCp  
 D: Green/Red→Positive hLf


A: Green/Red→Positive hHb  
 B: Green/Red→Positive hTf  
 C: Green/Red→Positive hCp  
 D: Green→Negative hLf

# CERTEST FOB+Transferrin+Calprotectin+Lactoferrin


One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette


A: Green/Red→Positive hHb  
B: Green/Red→Positive hTf  
C: Green→Negative hCp  
D: Green/Red→Positive hLf


A: Green/Red→Positive hHb  
B: Green→Negative hTf  
C: Green/Red→Positive hCp  
D: Green/Red→Positive hLf


A: Green→Negative hHb  
B: Green/Red→Positive hTf  
C: Green/Red→Positive hCp  
D: Green/Red→Positive hLf


A: Green/Red→Positive hHb  
B: Green/Red→Positive hTf  
C: Green→Negative hCp  
D: Green→Negative hLf


A: Green→Negative hHb  
B: Green→Negative hTf  
C: Green/Red→Positive hCp  
D: Green/Red→Positive hLf


A: Green→Negative hHb  
B: Green/Red→Positive hTf  
C: Green→Negative hCp  
D: Green/Red→Positive hLf


A: Green→Negative hHb  
B: Green/Red→Positive hTf  
C: Green/Red→Positive hCp  
D: Green→Negative hLf


A: Green/Red Positive hHb  
B: Green→Negative hTf  
C: Green/Red→Positive hCp  
D: Green→Negative hLf


A: Green/Red→Positive hHb  
B: Green→Negative hTf  
C: Green→Negative hCp  
D: Green/Red→Positive hLf


A: Green/Red→Positive hHb  
B: Green→Negative hTf  
C: Green→Negative hCp  
D: Green→Negative hLf

A: Green→Negative hHb  
B: Green/Red→Positive hTf  
C: Green→Negative hCp  
D: Green→Negative hLf

A: Green→Negative hHb  
B: Green→Negative hTf  
C: Green/Red→Positive hCp  
D: Green→Negative hLf

A: Green→Negative hHb  
B: Green→Negative hTf  
C: Green→Negative hCp  
D: Green/Red→Positive hLf

	A (hHb)	B (hTf)	C (hCp)	D (hLf)	Interpretation of the results
1.	- GREEN	- GREEN	- GREEN	- GREEN	No occult blood presence (no haemoglobin, transferrin presence→no gastrointestinal bleeding). No active gastrointestinal inflammation disorder (no calprotectin, lactoferrin presence).
2.	+ GREEN-RED	+ GREEN-RED	+ GREEN-RED	+ GREEN-RED	Patient stool sample contains <b>haemoglobin, transferrin, calprotectin</b> and <b>lactoferrin</b> markers, which might mean gastrointestinal bleeding and inflammation (organic disease) caused by Ulcerative Colitis, Crohn's disease, some types of carcinomas (colorectal cancer), ulcers or bacterial infections ( <i>Shigella</i> , <i>Salmonella</i> , <i>Campylobacter</i> , <i>E. coli</i> (EIEC or EAEC) or <i>C. difficile</i> ).

## CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

*One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette*

3.	+	+	+	-	Patient stool sample contains <b>haemoglobin, transferrin</b> and <b>calprotectin</b> markers, which might mean gastrointestinal bleeding and inflammation (organic disease) caused by CD, UC or some carcinoma (e.g. colorectal cancer).
	GREEN-RED	GREEN-RED	GREEN-RED	GREEN	
4.	+	+	-	+	Patient sample contains <b>haemoglobin, transferrin</b> and <b>lactoferrin</b> markers, which might mean gastrointestinal bleeding and inflammation (organic disease) caused by CD or UC, some type of carcinoma (e.g. colorectal cancer) and with persistent diarrhea: bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC) or C. difficile</i> ).
	GREEN-RED	GREEN-RED	GREEN	GREEN-RED	
5.	+	-	+	+	Patient stool sample contains <b>haemoglobin, calprotectin</b> and <b>lactoferrin</b> markers, which might mean lower gastrointestinal bleeding and inflammation (organic disease) caused by CD, UC, some carcinoma (e.g. colorectal cancer) and with persistent diarrhoea: bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC) or C. difficile</i> ).
	GREEN-RED	GREEN	GREEN-RED	GREEN-RED	
6.	-	+	+	+	Patient stool sample contains <b>transferrin, calprotectin</b> and <b>lactoferrin</b> markers, which might mean upper gastrointestinal bleeding and inflammation (organic disease) caused by CD (no UC), some carcinoma (e.g. gastric cancer) and with persistent diarrhoea: bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC) or C. difficile</i> ).
	GREEN	GREEN-RED	GREEN-RED	GREEN-RED	
7.	+	+	-	-	Patient stool sample contains <b>haemoglobin</b> and <b>transferrin</b> markers, which might mean gastrointestinal bleeding problems caused by some bacterial infections or some carcinomas without inflammation or some traumas.
	GREEN-RED	GREEN-RED	GREEN	GREEN	
8.	-	-	+	+	Patient stool sample contains <b>calprotectin</b> and <b>lactoferrin</b> markers, which might mean organic gastrointestinal inflammation caused by ulcerative colitis, Crohn's disease, some types of carcinomas or bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC), or C. difficile</i> ).
	GREEN	GREEN	GREEN-RED	GREEN-RED	
9.	-	+	-	+	Patient stool sample contains <b>transferrin</b> and <b>lactoferrin</b> markers, which might mean upper gastrointestinal bleeding and inflammation (organic disease) caused by CD (no UC), some carcinoma (e.g. gastric cancer) and with persistent diarrhea bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC) or C. difficile</i> ).
	GREEN	GREEN-RED	GREEN	GREEN-RED	
10.	-	+	+	-	Patient stool sample contains <b>transferrin</b> and <b>calprotectin</b> markers, which might mean upper gastrointestinal bleeding and inflammation (organic disease) caused by CD or some carcinoma (e.g. gastric cancer).
	GREEN	GREEN-RED	GREEN-RED	GREEN	
11.	+	-	+	-	Patient stool sample contains <b>haemoglobin</b> and <b>calprotectin</b> markers, which might mean lower gastrointestinal bleeding and inflammation (organic disease) caused by CD, UC or some carcinoma (e.g. colorectal cancer).
	GREEN-RED	GREEN	GREEN-RED	GREEN	
12.	+	-	-	+	Patient stool sample contains <b>haemoglobin</b> and <b>lactoferrin</b> markers, which might mean lower gastrointestinal bleeding and inflammation (organic disease) caused by CD, UC, some carcinoma (e.g. colorectal cancer) and with persistent diarrhea: bacterial enteropathies ( <i>Shigella, Salmonella, Campylobacter, E. coli (EIEC or EAEC) or C. difficile</i> ).
	GREEN-RED	GREEN	GREEN	GREEN-RED	
13.	+	-	-	-	Patient stool sample contains <b>haemoglobin</b> marker, which might mean lower gastrointestinal bleeding disease (e.g. early stage colorectal cancer). The concentration of human transferrin in blood is 100 times less than human haemoglobin, not too much blood present in faeces.
	GREEN-RED	GREEN	GREEN	GREEN	
14.	-	+	-	-	Patient stool sample contains <b>transferrin</b> marker, which might mean upper gastrointestinal bleeding (e.g. caused by gastric cancer). Human haemoglobin probably degraded in the gastrointestinal tract.
	GREEN	GREEN-RED	GREEN	GREEN	
15.	-	-	+	-	Patient stool sample contains <b>calprotectin</b> marker, which might mean some gastrointestinal inflammation (e.g. Crohn's disease, ulcerative colitis, etc.), without gastrointestinal bleeding and colorectal cancer presence.
	GREEN	GREEN	GREEN-RED	GREEN	
16.	-	-	-	+	Patient stool sample contains <b>lactoferrin</b> marker, which might mean some gastrointestinal inflammation (bacterial enteropathy), without gastrointestinal bleeding presence.
	GREEN	GREEN	GREEN	GREEN-RED	
17.	Any other result				Invalid result either A, B, C or D, we recommend repeating the assay using the same sample with another test.

## CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette

Whatever result detailed in the above-mentioned table (2 to 16) should be followed up with additional confirmatory diagnostic procedures to determine the exact cause and source of the blood or inflammation markers in the patient sample.

**INVALID:** Total absence of any control coloured lines (GREEN) regardless the appearance or not of the test lines (RED). Insufficient specimen volume, incorrect procedural techniques or deterioration of the reagents are mostly the main reasons for control lines failure. Review the procedure and repeat the assay with a new test. If the symptoms or situation still persist, discontinue using the test kit and contact your local distributor.

### NOTES ON THE INTERPRETATION OF RESULTS

The intensity of the red coloured bands in the test lines (T) in the results windows will vary depending on the concentration of human haemoglobin, human transferrin, human calprotectin and human lactoferrin present in the specimen.

### QUALITY CONTROL

Internal procedural controls are included in the test. The green lines appearing in the control lines (C) in the results windows are internal controls, which confirm sufficient specimen volume and correct procedural technique.

### LIMITATIONS

1. The test must be carried out within 2 hours after opening the sealed bag.
2. An excess of stool sample could cause wrong results (brown bands appear). Dilute the sample with the diluent and repeat the test.
3. The intensity of test line may vary from very strong at high human haemoglobin/transferrin/calprotectin/lactoferrin concentration to faint when the human haemoglobin/transferrin/calprotectin/lactoferrin concentration is close to the cut-off value of the test.
4. CerTest FOB+Transferrin+Calprotectin+Lactoferrin should be used only with samples from human faeces. The use of other samples has not been established. The quality of the test depends on the quality of the sample; proper faecal specimens must be obtained.
5. Positive results, in haemoglobin and/or transferrin, determine the presence of human haemoglobin and/or human transferrin in faecal samples; nevertheless, it can be due to several causes, besides colorectal bleeding, such as hemorrhoids, blood in urine or stomach irritations. A positive result should be followed up with additional diagnostic procedures to determine the exact cause and source of the blood in the stool.
6. Positive results, in calprotectin and/or lactoferrin, determine the presence of human calprotectin or/and human lactoferrin in faecal samples; nevertheless, it can be due to several causes (Inflammatory bowel disease, colorectal cancer and some enteropathies). A positive result should be followed up with additional diagnostic invasive procedures, a colonoscopy and a biopsy in order to confirm the diagnosis and to establish the inflammation extent.
7. A negative result, in haemoglobin and/or transferrin, is not meaningful because of it is possible the haemoglobin or transferrin concentration in the stool sample is lower than cut-off values. Negative results do not exclude bleeding, as some polyps and colorectal cancers may bleed intermittently or not during certain stages of the disease. Moreover, blood may not uniformly distributed in stool samples.
8. A negative result, in calprotectin and/or lactoferrin, is not meaningful because of it is possible the human calprotectin or human lactoferrin concentration in the stool sample is lower than cut-off values. If symptoms or situation still persist, calprotectin and lactoferrin determination should be carried out invasive techniques. Negative results do not exclude inflammation, some diseases such as celiac sprue and microscopic colitis polyps that involve mainly mononuclear inflammation.
9. Stool from patients with non-steroidal anti-inflammatory drug (NSAID) treatment could show positive calprotectin result.
10. Lactoferrin is a component of breast milk; the test will be positive in breast fed children and should not be used to evaluate neonates receiving breast milk.
11. Neonatal faecal calprotectin levels have been reported higher than normal in children with a median of 167µg/g (range 22-860µg/g).
12. Stool from patients with active inflammatory bowel diseases that usually involve significant neutrophilic inflammation of the intestine, such as Crohn's disease and ulcerative colitis, would be positive for faecal calprotectin (the level of human calprotectin in Crohn's disease is higher than in ulcerative colitis) and would be positive for faecal lactoferrin. CerTest FOB+Transferrin+Calprotectin+Lactoferrin could be sensitive for this diagnosis in patients with chronic diarrhea.
13. Patients should not collect samples during their menstrual period, if they have bleeding hemorrhoids, blood in urine or if they have strained during bowel movement.
14. Bloody stool samples can contain components that may cause non-specific reactions in the test. Every bloody stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.
15. Mucous stool samples may cause non-specific reactions in the test. Every mucous stool sample whose result is positive should be followed up with other techniques of diagnosis to confirm the result.

## EXPECTED VALUES

Colorectal cancer affects both men and women of all racial and ethnic groups, being most often found in people aged 50 years or older. Regarding men, colorectal cancer is the third most common cancer after prostate and lung cancers. Regarding women, colorectal cancer is the third most common cancer after breast and lung cancers.

Some studies established equal or higher 5.1µg hHb/g faeces for haemoglobin and 0.4µg hTf/g faeces for transferrin as cut-off values to allow detect adult patients with gastrointestinal bleeding that will require to diagnosis additional diagnostic invasive procedures.

Some studies established equal or higher 50µg hCp/g faeces and 10µg hLf/g faeces as cut-off values to allow detect adult patients with gastrointestinal inflammatory pathology that will require to diagnosis additional diagnostic invasive procedures.

## PERFORMANCE CHARACTERISTICS

### Analytical sensitivity

Cut-off values of CerTest FOB+Transferrin+Calprotectin+Lactoferrin are 50ng/mL (5.1µg hHb/g faeces) for human haemoglobin, 4ng/mL (0.4µg hTf/g faeces) for human transferrin, 500ng/mL (50µg hCp/g faeces) for human calprotectin and 100ng/mL (10µg hLf/g faeces) for human lactoferrin.

### Clinical sensitivity and specificity (Strip A and Strip B)

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest FOB+Transferrin+Calprotectin+Lactoferrin, CerTest) and commercial ELISA tests (Human Hemoglobin ELISA kit, ICL, and Human Transferrin ELISA kit, ICL). The results were as follows:

#### Data for cut-off value 5 µg hHb/g

		ELISA test: Human Hemoglobin ELISA kit		
		+	-	Total
CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip A: FOB)	+	38	11	49
	-	0	49	49
	Total	38	60	98

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip A: FOB) vs Human Hemoglobin ELISA kit		
	Mean Value	95% confidence interval
<b>Sensitivity</b>	100.0%	90.7 – 100.0%
<b>Specificity</b>	81.7%	69.6 – 90.5%
<b>PPV</b>	77.6%	63.4 – 88.2%
<b>NPV</b>	100.0%	92.7 – 100.0%

Specificity for a cut-off value of 1 µg hHb/g: 93.5% (mean value); 82.1-98.6% (95% confidence interval).

#### Data for cut-off value 0.4 µg hTf/g

		ELISA test: Human Transferrin ELISA kit		
		+	-	Total
CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip B: Transferrin)	+	38	3	41
	-	12	47	59
	Total	50	50	100

## CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip B: Transferrin) vs Human Transferrin ELISA kit		
	Mean Value	95% CI (confidence interval)
<b>Sensitivity</b>	76.0%	61.8 – 86.9%
<b>Specificity</b>	94.0%	83.5 – 98.7%
<b>PPV</b>	92.7%	80.1 – 98.5%
<b>NPV</b>	79.7%	67.2 – 98.0%

Sensitivity for a cut-off value of 2 µg hTf/g: 81.1% (mean value); 64.8-92.0% (95% confidence interval).

Specificity for a cut-off value of 0.2 µg hTf/g: 96.3% (mean value); 81.0-99.9% (95% confidence interval).

### Clinical sensitivity and specificity (Strip C, Strip D)

An evaluation, with faecal samples, was performed comparing the results obtained by an immunochromatographic test (CerTest FOB+Transferrin+Calprotectin+Lactoferrin, CerTest), and other commercial immunoassays (Calprest®, Eurospital) and other commercial ELISA test (Human Lactoferrin ELISA kit, ICL). The results were as follows:

#### Data for cut-off value 50 µg hCp/g

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip C: Calprotectin)	ELISA test: Calprest®			Total
	+	-		
+	31	7		38
-	1	41		42
Total	32	48		80

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip C: Calprotectin) vs Calprest®		
	Mean Value	95% CI (confidence interval)
<b>Sensitivity</b>	96.9%	83.8 – 99.9%
<b>Specificity</b>	85.4%	72.2 – 93.9%
<b>PPV</b>	81.6%	65.7 - 92.3%
<b>NPV</b>	97.6%	87.4 – 99.9%

Sensitivity for a cut-off value of 100 µg hCp/g: 100.0% (mean value); 85.2-100.0% (95% confidence interval).

Specificity for a cut-off value of 25 µg hCp/g: 94.9% (mean value); 82.7-99.4% (95% confidence interval).

#### Data for cut-off value 10 µg hLf/g

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip D: Lactoferrin)	ELISA test: Human Lactoferrin ELISA kit			Total
	+	-		
+	27	13		40
-	3	47		50
Total	30	60		90

## CERTEST FOB+Transferrin+Calprotectin+Lactoferrin

One Step test to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin combo card format/  
Prueba combo de un solo paso para detección de hemoglobina, transferrina, calprotectina y lactoferrina humanas en formato cassette

CerTest FOB+Transferrin+Calprotectin+Lactoferrin (Strip D: Lactoferrin) vs Human Lactoferrin ELISA kit		
	Mean Value	95% CI (confidence interval)
<b>Sensitivity</b>	90.0%	73.5 – 97.9%
<b>Specificity</b>	78.3%	65.8 – 87.9%
<b>PPV</b>	67.5%	50.9 – 81.4%
<b>NPV</b>	94.0%	83.5 – 98.7%

Specificity for a cut-off value of 2.5 µg hLf/g: 100.0% (mean value); 83.9-100.0% (95% confidence interval).

The results showed a high sensitivity and specificity to detect human haemoglobin, human transferrin, human calprotectin and human lactoferrin using CerTest FOB+Transferrin+Calprotectin+Lactoferrin.

### Cross reactivity

An evaluation was performed to determine the cross reactivity of CerTest FOB+Transferrin+Calprotectin+Lactoferrin; no cross reactivity against other faecal markers occasionally present in faeces:

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## Clinical sensitivity and specificity (Strip C, Strip D)

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1. Towler BP, Irwig L, Glasziou P, Weller D, Kewenter J. Screening for colorectal cancer using the faecal occult blood test, Hemoccult. Cochrane Database Syst Rev. 2000;(2): CD001216.
2. Ransohoff DF and Lang CA. Screening for colorectal cancer with the Fecal Occult Blood Test: a background paper. Ann Intern Med. 1997; 126: 811-822.
3. Ransohoff DF and Lang CA. Suggested technique for Faecal Occult Blood testing and interpretation in colorectal cancer screening. Ann Intern Med. 1997; 126: 808-810.
4. Chien-Hua Chiang et al., "A Comparative Study of Three Fecal Occult Blood Test in Upper Gastrointestinal Bleeding". Kaohsiung J. Med. Sci. May 2006, 22; 223-228.
5. Virtanen et al., "Higher concentrations of serum transferrin receptor in children than in adults", Am. J. Clin. Nutr. 1999; 69: 256-60.
6. Angriman I. et al. Enzymes in feces: Useful markers of chronic inflammatory bowel disease. Clinica Chimica Acta 381 Feb 2007, p. 63-68.
7. Quail, M.A. et al. Fecal Calprotectin Complements Routine Laboratory Investigations in Diagnosing Childhood Inflammatory Bowel Disease. Inflamm Bowel Dis, Vol 15 No 5; May 2009, p. 756-759.
8. Gaya D.R., et al. Faecal calprotectin in the assessment of Crohn's disease activity. Q J Med 2005, Vol 98, May 2005, p. 435-441.
9. Langhorst, M.D. et al. Noninvasive Markers in the Assessment of Intestinal Inflammation in Inflammatory Bowel Diseases: Performance of Fecal Lactoferrin, Calprotectin and PMN-Elastase, CRP, and Clinical Indices. Am. J. Gastroenterol. 2008; Vol 103, p. 162-169.

## SYMBOLS FOR IVD COMPONENTS AND REAGENTS/SÍMBOLOS PARA REACTIVOS Y PRODUCTOS PARA DIAGNÓSTICO IN VITRO

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