

# DyLight® Amine-Reactive Dyes

2032.11

Number	Description
46427	DyLight 350 NHS Ester, $5 \times 65 \mu g$
46401	DyLight 405 NHS Ester, $5 \times 50 \mu g$
46403	DyLight 488 NHS Ester, $5 \times 50 \mu g$
62263	DyLight 550 NHS Ester, $5 \times 50 \mu g$
46413	DyLight 594 NHS Ester, $5 \times 65 \mu g$
46417	DyLight 633 NHS Ester, $5 \times 50 \mu g$
62266	DyLight 650 NHS Ester, $5 \times 50 \mu g$
46419	DyLight 680 NHS Ester, $5 \times 50 \mu g$
62279	DyLight 755 NHS Ester, $5 \times 50 \mu g$
46422	DyLight 800 NHS Ester, $5 \times 50 \mu g$

**Note:** Each vial is sufficient to label 0.25-1mg of protein (> 45kDa).

**Storage:** Upon receipt store DyLight NHS Esters at -20°C. Products are shipped with an ice pack. Store all dyes in the foil pouch with desiccant to protect from light and moisture.

# Introduction

The Thermo Scientific DyLight Amine-Reactive Dyes have absorption spectra ranging from 350nm to 770nm (Table 1) and are packaged in a convenient amount suitable for one labeling reaction. These reagents fluoresce over a broad pH range, are more intense than Alexa Fluor® or Cy® Dyes in many applications and match the output wavelengths of common fluorescence instrumentation. Additionally, the single-use packaging and water solubility of the DyLight Reagents allows protein samples to be added directly to the reagent vial for high dye-to-protein ratio conjugations without precipitation.

The amine-reactive dyes contain *N*-hydroxysuccinimide (NHS) esters, the most commonly used reactive group for labeling proteins. NHS esters react with primary amines, forming a stable, covalent amide bond and releasing the NHS group.

Table 1. Properties of the DyLight NHS-Ester Dyes.

DyLight Dye	Ex/Em*	ε†	MW (g/mol)	Spectrally Similar Dyes
350	353 / 432	15,000	874	Alexa Fluor 350, AMCA
405	400 / 420	30,000	793	Alexa Fluor 405
488	493 / 518	70,000	1011	Alexa Fluor 488, Cy2
550	562 / 576	150,000	1040	Alexa Fluor 555, Cy3
594	593 / 618	80,000	1078	Alexa Fluor 594, Texas Red
633	638 / 658	170,000	1066	Alexa Fluor 633
650	652 / 672	250,000	1066	Alexa Fluor 647, Cy5
680	682 / 715	140,000	950	Alexa Fluor 680, Cy5.5
755	755 / 776	220,000	1092	Alexa Fluor 750
800	770 / 794	270,000	1050	IRDye 800

<sup>\*</sup> Excitation and emission maxima in nanometers †Molar extinction coefficient (M<sup>-1</sup> cm<sup>-1</sup>)



# Important Product Information

- NHS ester-activated fluorophores are moisture-sensitive. Store product in the original pouch at -20°C. Avoid moisture condensation onto the product by equilibrating the vial to room temperature before opening. Prepare these labeling reagents immediately before use. Do not store NHS-ester reagents prepared in aqueous solutions.
- Use the following fluorescent imagers:
  - 350 dye: UV argon-ion laser at 351-363nm
  - 405 dye: Spectral line of the blue diode laser
  - 488 dye: Green (526) laser
  - 550 and 594 dyes: Green (532) laser
  - 633 and 650 dyes: Red (633) laser
  - 680, 755 and 800 dyes: laser- and filter-based instruments that emit in the 700nm and 800nm region of the spectrum, respectively; these dyes are well-suited for the 700 and 800 channels of the LI-COR Odyssey<sup>®</sup> and the LI-COR Aerius<sup>TM</sup> Infrared Imaging Systems.
- Low concentrations of sodium azide (≤ 3mM or 0.02%) or thimerosal (≤ 0.02mM or 0.01%) will not significantly interfere with protein labeling; however, 20-50% glycerol will reduce labeling efficiency.
- To remove excess non-reacted DyLight Dye, use a dialysis membrane with a molecular-weight cutoff ≥ 10K.

# **Procedure for Protein Labeling**

The following is an example application for the DyLight Amine-Reactive Dyes. Specific applications will require optimization.

#### A. Protein Preparation

The optimal labeling buffer is 0.05M sodium borate buffer at pH 8.5 (Product No. 28384). When labeling with DyLight 594 NHS-Ester, prepare the protein in phosphate-buffered saline to avoid precipitation. Buffers that contain primary amines (e.g., Tris or glycine) will interfere because they react with the NHS-ester moiety. Dissolve protein directly in the labeling buffer. For each labeling reaction, use  $100\text{-}500\mu\text{L}$  of purified protein sample at 1-2.5mg/mL. If the protein is already in a buffer, perform a buffer exchange into the labeling buffer by dialysis or gel filtration.

**Note:** The following buffers may be substituted for borate buffer: 0.1M sodium phosphate, 0.15M NaCl at pH 7.2-7.5 (e.g., Thermo Scientific BupH Phosphate Buffered Saline Packs, Product No. 28372) or 0.1M sodium carbonate at pH 8.3-9.0.

#### **B.** Labeling Reaction

**Note:** The DyLight NHS-Ester reagents are moisture-sensitive. Store the reagent in the original container at -20°C with desiccant.

- 1. Transfer the protein solution to the vial containing the dye. Mix well by vortexing up-and-down several times and incubate at room temperature for 1 hour.
- 2. Remove excess dye reagent from the sample using the Thermo Scientific Dye Removal Columns (Product No. 22858) or a dialysis membrane with a molecular-weight cutoff ≥ 10K.

**Note:** The non-reacted dye must be completely removed for optimal results and accurate determination of the dye-to-protein ratio. For best results, use the Dye Removal Columns or dialyze for ~4 hours using three dialysis buffer changes. Gel filtration, such as with desalting columns, is not typically as effective as dialysis.

3. Store labeled protein protected from light at 4°C for up to one month. Alternatively, store labeled protein in single-use aliquots at -20°C.

#### C. Calculate the Degree of Labeling

1. Dilute a small amount of labeled, purified protein in PBS. Using a 1cm path length cuvette, measure the absorbance at 280nm and the  $A_{max}$  of the specific dye (Table 2).

**Table 2.** Properties of the DyLight Dyes.

DyLight Dye	A <sub>max</sub> *	ε†	CF‡
350	353	15,000	0.144
405	405	30,000	0.564
488	493	70,000	0.147
550	557	150,000	0.081
594	595	80,000	0.585
633	627	170,000	0.110
650	655	250,000	0.037
680	684	140,000	0.128
755	755	220,000	0.030
800	777	270,000	0.045

<sup>\*</sup> Excitation wavelength in nanometers – note that upon protein conjugation the absorption maximum shifts to the right of the spectra

2. Calculate protein concentration as follows:

Protein concentration (M) = 
$$\frac{[A_{280} - (A_{max} \times CF)]}{\epsilon_{protein}} \times dilution factor$$

- $\epsilon_{protein}$  = protein molar extinction coefficient (e.g., the molar extinction coefficient of IgG is ~210,000 M<sup>-1</sup> cm<sup>-1</sup>)
- CF = Correction factor =  $\frac{A_{280} \text{ of the fluor}}{A_{\text{max}} \text{ of the fluor}}$  (see Table 2)
- 3. Calculate the degree of labeling:

Moles dye per mole protein = 
$$\frac{A_{max} \text{ of the labeled protein} \times \text{dilution factor}}{\epsilon_{fluor} \times \text{protein concentration (M)}}$$

•  $\varepsilon_{\text{dve}} = \text{See Table 2}$ 

# Example calculations for DyLight 550 NHS Ester-conjugated antibody:

- Dilution factor = 40
- $A_{280} = 0.177$
- $A_{max}$  at 557nm = 0.411

Protein concentration (M) = 
$$\frac{[0.177 - (0.411 \times 0.081)]}{210,000} \times 40 = 0.00002737 \text{ M}$$

Moles dye per mole protein = 
$$\frac{0.177 \times 40}{150,000 \times 0.00002737} = 2.9$$

# **Troubleshooting**

Problem	Cause	Solution
The application in which the dye-labeled	The protein was not labeled	Before troubleshooting, determine if the protein is labeled by calculating the $A_{max}$ : $A_{280}$ ratio; determine
protein was used was		this ratio after thorough desalting or dialysis
unsuccessful		<b>Note:</b> For dye-labeled antibodies the $A_{max}$ : $A_{280}$ ratio
		should be > 1.
The protein was not	Conjugation Buffer contained	Use a conjugation buffer free of primary amines such
labeled	primary amines (e.g., Tris or glycine)	as borate, carbonate or PBS
	that interfered with the reaction	
	The NHS ester has hydrolyzed and is	Prepare labeling reagent immediately before use – do
	non-reactive	not store NHS-ester reagents in aqueous solutions

<sup>†</sup>Molar extinction coefficient (M<sup>-1</sup> cm<sup>-1</sup>) at A<sub>max</sub> ‡Correction factor (A<sub>280</sub>/A<sub>max</sub>)



# **Additional Information**

Visit our website for additional information including the following items:

- Tech Tip #43: Protein stability and storage
- Tech Tip #3: Determine reactivity of NHS-ester biotinylation and crosslinking reagents
- Tech Tip #30: Modify and label oligonucleotide 5' phosphate groups

# **Related Thermo Scientific Products**

46426	DyLight 350 NHS Ester, 1mg
46400	DyLight 405 NHS Ester, 1mg
46402	DyLight 488 NHS Ester, 1mg
62262	DyLight 550 NHS Ester, 1mg
46412	DyLight 594 NHS Ester, 1mg
46414	DyLight 633 NHS Ester, 1mg
62265	DyLight 650 NHS Ester, 1mg
46418	DyLight 680 NHS Ester, 1mg
62278	DyLight 755 NHS Ester, 1mg
46421	DyLight 800 NHS Ester, 1mg
69576	Slide-A-Lyzer <sup>®</sup> MINI Dialysis Unit Kit, for 10-100μL samples, 10 units plus float
66382	Slide-A-Lyzer Dialysis Cassette Kit, 10K MWCO, for 0.5-3mL samples, 10 units, buoys and syringes
66807	Slide-A-Lyzer Dialysis Cassette Kit, 10K MWCO, for 3-12mL samples, 10 units, buoys and syringes

Cy is a trademark of GE Healthcare Limited UK.

Alexa Fluor is a trademark of Life Technologies, Inc.

Odyssey and Aerius are trademarks of LI-COR Biosciences.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at <a href="https://www.thermoscientific.com/pierce">www.thermoscientific.com/pierce</a>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.