INSTRUCTIONS

HiPPRTM Detergent Removal Resin



88306 88307 88305 2300.3 Description Number 88305 HiPPR Detergent Removal Spin Column Kit, 5mL of resin supplied as 50% slurry (e.g., 5mL of settled resin is equivalent to 10mL of 50% slurry) in 0.15M sodium chloride and 0.05% sodium azide, with two Spin Column Accessory Packs, 54 columns 88306 HiPPR Detergent Removal Spin Columns, 24 columns, each column contains ~100µL of resin supplied as 25% slurry (e.g., 100µL of settled resin is equivalent to 400µL of 25% slurry) in 0.15M sodium chloride and 0.05% sodium azide HiPPR Detergent Removal 96-Well Filter Plate Kit, each well contains ~100µL of resin supplied as 88307 50% slurry (e.g., 100µL of settled resin is equivalent to 400µL of 25% slurry) in 0.15M sodium chloride and 0.05% sodium azide

Storage: Upon receipt store at 4°C. Products shipped at ambient temperature.

Table of Contents


Introduction	1
Procedure Summary	2
Important Product Information	
Additional Materials Required	2
Detergent Removal Procedure using Spin Columns	
Detergent Removal Procedure using 96-well Filter Plate	
Troubleshooting	4
Additional Information	
Related Thermo Scientific Products	
Cited References	

Introduction

The Thermo Scientific HiPPR (<u>High Protein and Peptide Recovery</u>) Detergent Removal Resin removes commonly used ionic, nonionic and zwitterionic detergents from samples with low protein/peptide concentrations (1 to 100μ g/mL). Detergents and surfactants are used for solubilizing, stabilizing and disaggregating proteins and protein complexes; however, detergents interfere with downstream analyses such as ELISA, isoelectric focusing and mass spectrometry (MS). The HiPPR Detergent Removal Resin removes detergents with > 95% efficiency at concentrations of 0.5-1%, enabling recovery and characterization of low-abundant proteins. The columns and 96-well filter plates are pre-loaded with detergent-removal resin for multiple-sample processing of 25 to 100μ L.



Procedure Summary



Important Product Information

- The sample must be added directly onto the resin to remove detergents with efficient protein/peptide recovery.
- To remove the column's bottom plug, pull out the plug in one continuous motion without twisting.
- Use one part sample to one part settled resin (e.g., 50μL of sample to 50μL of settled resin) for removing SDS, sodium deoxycholate, TritonTM X-100, Triton X-114 and NP-40. Please refer to Table 1 and 2 for the maximum removable detergent concentrations in the spin-column and plate formats, respectively.
- Use one part sample to two parts settled resin (e.g., 50µL of sample to 100µL of settled resin) for removing CHAPS, octyl glucoside and octyl thioglucoside. Please refer to Table 1 and 2 for the maximum removable detergent concentrations in the spin-column and plate formats, respectively.

Detergent	Maximum Removable Concentration (%)
SDS	1
Sodium deoxycholate	1
Triton X-100	1
Triton X-114	0.5
NP-40	0.75
CHAPS*	1
Octyl glucoside*	1
Octyl thioglucoside*	1
*Use one part comple to two	parts sattlad rasin

Table 1. Maximum removable detergent concentrationsusing the Thermo Scientific HiPPR Detergent RemovalSpin Column.

*Use one part sample to two parts settled resin

Table 2. Maximum removable detergent concentrationsusing the Thermo Scientific HiPPR Detergent Removal96-well Filter Plate.

Detergent	Maximum Removable Concentration (%)	
SDS	0.4	
Sodium deoxycholate	0.5	
Triton X-100	0.25	
NP-40	0.125	
CHAPS*	1	
Octyl glucoside*	1	
Octyl thioglucoside*	1	
*Maximum processable san	nle volume is 50µI	

⁴Maximum processable sample volume is 50µL



Detergent Removal Procedure using Spin Columns

Materials Required

- Variable-speed bench-top microcentrifuge with a minimum speed of $1500 \times g$
- 2mL microcentrifuge collection tubes, multi-channel pipettor and tips
- Equilibration buffer: use a buffer such as AMBIC, PBS, MES, MOPS, Tris or carbonate-bicarbonate at pH 5 to 10. Do not use buffers containing organic solvents.

A. Method

- 1. Place an empty spin column into a 2mL collection tube. If using a pre-dispensed HiPPR Detergent Removal Column (Product No. 88306), remove the bottom plug, place the spin column into a 2mL collection tube and proceed to step 5.
- 2. Gently swirl the bottle of HiPPR Detergent Removal Resin to obtain an even suspension.
- 3. Using a cut pipette tip, add an appropriate amount of resin slurry into the spin column (See Important Product Information).
- 4. Centrifuge at $1500 \times g$ for 1 minute to remove storage buffer. Discard the flow-through.
- 5. Add the appropriate volume of equilibration buffer (one part buffer to one part settled resin), centrifuge at $1500 \times g$ for 1 minute and discard the flow-through. Repeat this step two additional times.
- 6. Insert bottom plug and place column into a new collection tube.
- 7. Slowly apply sample directly onto the top of the compacted resin bed. Cap the column, vortex gently to mix the sample with the resin and incubate at room temperature for 10 minutes.
- 8. Slightly loosen the cap, remove the bottom plug and place the column in a collection tube. Centrifuge at $1500 \times g$ for 2 minutes to collect the detergent-free sample. Discard the used column.

Detergent Removal Procedure using a 96-well Filter Plate

Materials Required

- Centrifuge with a minimum speed of $1000 \times g$ and rotor and carriers capable of handling staked plates with a height of 4.4cm.
- Equilibration buffer: use a buffer such as AMBIC, PBS, MES, MOPS, Tris or carbonate-bicarbonate at pH 5 to 10. Do not use buffers containing organic solvents.

A. Method

- 1. Equilibrate the HiPPR Detergent Removal 96-well Filter Plate to room temperature.
- 2. Carefully remove the sealing material from the bottom of the plate. Place the detergent removal plate on top of a wash plate.
- 3. Remove the sealing material from the top of the plate.
- 4. Place the assembly into a centrifuge with a 96-well plate-carrier rotor and centrifuge at $1000 \times g$ for 2 minutes to remove the storage buffer. Discard the flow-through and replace the detergent-removal plate on top of the wash plate.
- 5. Add 200μ L of equilibration buffer to each well and centrifuge at $1000 \times g$ for 2 minutes. Discard the flow-through. Repeat this step two additional times.
- 6. Stack the detergent-removal plate on top of a sample-collection plate (blue), aligning the alphanumeric indices on the plate.
- 7. Slowly apply sample $(25-100\mu L)$ to the center of the resin bed. To expel the entire sample from the pipette tip, carefully touch the top of the resin. Cover the plate with a new plate sealer and incubate for 10 minutes at room temperature.

Note: To process > 96 samples, evenly divide the samples between two plates. To process \leq 96 samples, balance the centrifuge by using an unprocessed plate with the top and bottom seal in place.

8. Centrifuge the plate assembly at $1000 \times g$ for 2 minutes to collect the detergent-free sample. Discard the detergent removal plate or reserve it for future balancing purposes.



Troubleshooting

Problem	Possible Cause	Solution	
Sample or buffer does not flow through	Centrifugation problem	Ensure the centrifuge is working properly	
the resin or plate	Bottom closure was not removed	Ensure bottom closure is removed	
Detergent is not completely removed	Improper sample loading	Load sample directly onto the top of the resin bed	
	Insufficient amount of resin used	d Use the appropriate resin amount (See Important Product Information) or slightly increase the resin amount (e.g., increase from 1 part sample:2 parts settled resin to 1 part sample:1.5 parts settled resin)	
	Exceeded the maximum removable detergent concentration	Ensure the detergent concentration is appropriate or dilute sample before processing (See Table 1 or 2)	
	Buffer contained organic solvents	Remove organic solvents by dialysis before detergent removal	

Additional Information

A. Typical detergent removal efficiency and protein recovery using Thermo Scientific HiPPR Detergent Removal Resin.[†]

	Sample	Detergent	Protein
Detergent	Volume (µL)	Removal (%)	Recovery (%)
	25	> 99	98
SDS	50	> 99	97
(1%)	100	> 99	100
	200	> 99	100
	25	98	82
Triton X-100 (1%)	50	95	86
	100	98	86
	200	98	93
	25	95	90
NP-40	50	96	94
(0.75%)	100	97	91
	200	97	97
CHAPS (1%)	25	95	64
	50	97	70
	100	98	78
	200	98	75

[†]A volume of 0.1mL of BSA (1.5µg) with detergent in 0.15M sodium chloride and 0.05% sodium azide was mixed with 0.1mL of detergent removal resin (0.2mL for CHAPS removal). Residual SDS was measured using Stains-All dye.¹ Triton X-100 and NP-40 were measured by absorbance at 275nm (protein absorbance was subtracted); CHAPS, sodium deoxycholate, octyl glucoside and octyl thioglucoside were measured by colorimetric method using concentrated sulfuric acid and phenol.² The Thermo Scientific Micro BCA Protein Assay was used for protein determination.

(815) 968-0747






Figure 1. LC-MS/MS analysis of BSA tryptic peptides before and after detergent removal. Base peak LC chromatograms of BSA tryptic peptides. BSA (100µg/mL) in 50mM ammonium bicarbonate buffer, pH 8.0 was reduced, alkylated and digested with trypsin for 3 hours at 37°C (enzyme-to-protein ratio, 1:50) in the presence of 0.5% detergents. To prepare the sample containing 0.5% SDS, the detergent was added to the sample following tryptic digestion. The 0.1mL digested sample containing the detergent was processed with 0.1mL of Thermo Scientific HiPPR Detergent Removal Resin as shown in the protocol. Control samples (shown above as Unprocessed) were not processed with the detergent removal resin. Samples (3pmol) were loaded directly onto a C18 column and subjected to LC-MS/MS analysis using a Thermo Scientific LTQ Mass Spectrometer.

Related Thermo Scientific Products

- 87776 Pierce[™] Detergent Removal Spin Columns, 125μL, 25 columns
- 87777 Pierce Detergent Removal Spin Columns, 0.5mL, 25 columns
- 87778 Pierce Detergent Removal Spin Columns, 2mL, 5 columns
- 87779 Pierce Detergent Removal Spin Columns, 4mL, 5 columns
- 87780 Pierce Detergent Removal Resin, 10mL
- **87782 Pierce C18 Tips, 10µL bed,** 96/pkg
- **87784 Pierce C18 Tips, 100µL bed,** 96/pkg
- 89870 Pierce C18 Spin Columns, 25/pkg
- 88302 Pierce Graphite Spin Columns, 30/pkg
- 23235 Micro BCA Protein Assay Kit



Cited References

- 1. Rusconi, E., *et al.* (2001). Quantitation of sodium dodecyl sulfate in microliter-volume biochemical samples by visible light spectroscopy. *Anal Biochem* **295**:31-7.
- 2. Urbani, A. and Warne, T. (2005). A colorimetric determination for flycosidic and bile-salt based detergents: applications in membrane protein research. *Anal Biochem* **336**:117-24.

Triton is a trademark of The Dow Chemical Company.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at <u>www.thermoscientific.com/pierce</u>. For a faxed copy, call 800-874-3723 or contact your local distributor. © 2012 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.