

NAb™ Spin Kits, 1mL for Antibody Purification

89978 89979 89980 89981

1941.4

Number	Description
89978	NAb Protein A Plus Spin Kit, 1mL Binding Capacity: ≥ 34mg human IgG per column
89979	NAb Protein G Spin Kit, 1mL Binding Capacity 11-15mg human IgG per column
89980	NAb Protein A/G Spin Kit, 1mL Binding Capacity: ≥ 7mg human IgG per column
89981	NAb Protein L Spin Kit, 1mL Binding Capacity: 5-10mg human IgG per column

Kit Contents:

NAb Spin Columns, 2 each, columns are supplied with top and bottom caps. Each column contains a 1mL resin bed of crosslinked 6% beaded agarose in 0.02% sodium azide

Binding Buffer, 1 pack (100mM phosphate, 150mM sodium chloride; pH 7.2 when dissolved in 500mL of ultrapure water)

IgG Elution Buffer, 240mL, pH 2.8

Neutralization Buffer, 12mL, 1M Tris•HCl, pH 8.5

Storage: Upon receipt store kit at 4°C. These products are shipped at ambient temperature.

Introduction

The Thermo Scientific™ NAb™ Spin Kits are convenient for rapid, small-scale affinity purification of antibodies from a variety of sample types. Each pre-filled microcentrifuge spin column of the immobilized protein resin enables quick purification of 1-13mg of IgG from 0.5-2mL of serum or other sample. The actual amount of IgG purified varies depending upon the sample type and the specific spin column used. Also included in each kit are buffers and a streamlined protocol for purifying at least two antibody samples.

Proteins A, G and L are different bacterial proteins that bind with high specificity to mammalian immunoglobulins. Immobilized forms of these proteins have been widely used for affinity purification of antibodies from serum, ascites fluid and hybridoma culture supernatant samples. The particular species and class of antibody to be purified determines which of these immobilized protein resins is most appropriate. The following paragraphs provide general guidelines for making this choice; consult our catalog or website for a more detailed description and table of antibody-binding characteristics for Proteins A, G, A/G and L.

Proteins A and G bind to many of the same species and subclasses of IgG, although they have particular differences in affinity and binding capacity. Protein A is generally preferred for affinity purification of rabbit, pig, dog and cat IgG. Protein G has better binding capacity for a broader range of mouse and human IgG subclasses (IgG₁, IgG₂, etc.). Protein A/G is a recombinant fusion protein that includes the IgG-binding domains of both Protein A and Protein G. Therefore, Protein A/G is ideal for binding the broadest range of IgG subclasses from rabbit, mouse, human and other mammalian samples.

Protein L binds to certain immunoglobulin kappa light chains. Because kappa light chains occur in members of all classes of immunoglobulin (i.e., IgG, IgM, IgA, IgE and IgD), Protein L can purify these different classes of antibody. However, only those antibodies within each class that possess the appropriate kappa light chains will bind. Generally, empirical testing is required to determine if Protein L is effective for purifying a particular antibody.

Additional Materials Required

- 15mL collection tubes
- Table-top centrifuge set to $1000 \times g$ for all centrifugation steps
- Storage solution: 0.02% sodium azide in phosphate-buffered saline (PBS)

Spin Purification Protocol

Note: Typically, the immobilized protein column may be used up to 10 times without significant loss in binding capacity.

1. Equilibrate column and buffers to room temperature. Set centrifuge to $1000 \times g$.
2. Prepare sample for purification by diluting to 2mL in Binding Buffer (maximum volume).
3. Loosen top cap on spin column and snap off bottom closure. Place column in a 15mL collection tube and centrifuge for 1 minute to remove storage solution. Discard the flow-through.
4. Equilibrate column by adding 2mL of Binding Buffer. Centrifuge for 1 minute and discard the flow-through. Repeat.
5. Cap bottom of column with the included rubber cap. Apply sample to column and tightly cap the top. Incubate at room temperature with end-over-end mixing for 10 minutes.
6. Loosen top cap and remove bottom cap. Place column in a new 15mL collection tube and centrifuge for 1 minute. Save the flow-through.

Note: This first collection tube contains the nonbound sample components and can be analyzed to assess binding efficiency and capacity.

7. Transfer column to a new 15mL collection tube. Wash column by adding 2mL of Binding Buffer and centrifuging for 1 minute. Repeat for a total of three washes.
8. Add 100 μ L of Neutralization Buffer to three 15mL collection tubes and place the spin column into one of the tubes.
9. Add 1 ml Elution Buffer to the column and centrifuge for 1 minute. Transfer the spin column to another collection tube that contains Neutralization Buffer, saving the collected solution as the first elution fraction. Repeat this step two times to obtain three fractions.
10. Determine which fraction(s) contain the purified antibody by measuring the relative absorbance of each fraction at 280nm. If required for downstream applications, exchange the buffer using Thermo Scientific™ Zeba™ Spin Desalting Columns or Slide-A-Lyzer™ Dialysis Cassettes (see Related Thermo Scientific Products Section).
11. To regenerate the column for storage or re-use, add 3 ml of Elution Buffer and centrifuge for 1 minute. Repeat. Wash column with 3 ml of PBS to remove elution buffer. Add 3 ml of storage solution and store column at 4°C. Typically, the immobilized protein column may be used up to 10 times without significant loss in binding capacity, although the actual number of effective usages may vary.

Gravity-flow Purification Protocol

1. Equilibrate column and buffers to room temperature.
2. Dilute sample 1:1 with Binding Buffer.
3. Gently tap column on bench top to dislodge any resin that may be in the top cap. Gently snap off bottom closure and remove the top cap. Place column in a 15mL collection tube and allow storage solution to drain.
4. Equilibrate column by adding 5mL of Binding Buffer and allow the solution to drain.
5. Apply the diluted sample to the column. For best results, add a sample volume that is less than 80% of the column's antibody-binding capacity. Collect the flow-through.

Note: If the sample contains more IgG than can bind to the column, the flow-through will contain the excess antibody. By saving the flow-through, nonbound antibody can be recovered and analyzed.

6. Wash column with 15mL of Binding Buffer.

Note: If desired, verify that all nonbound proteins are removed from the column by collecting separate 2mL fractions as the column drains. Measure the absorbance of each fraction at 280nm. The last fractions should have an absorbance similar to the Binding Buffer.

7. Add 100 μ L Neutralization Buffer to five collection tubes. Elute antibodies with 5mL of Elution Buffer, collecting 1mL fractions in each of the buffer-containing tubes.
8. Measure the relative absorbance of each fraction at 280nm and pool fractions as desired. If required for downstream applications, exchange the buffer using Zeba Spin Desalting Columns or Slide-A-Lyzer Dialysis Cassettes (see Related Thermo Scientific Products Section).
9. Regenerate column by adding 8mL of Elution Buffer and allow solution to flow through the column.
10. Store column by adding 5 ml of storage solution. When approximately 3mL remain in the column, cap bottom and secure top cap. Store column at 4°C.

Troubleshooting

Problem	Possible Cause	Solution
No protein detected in any elution fractions by absorbance at 280nm or general protein staining of electrophoresed sample	Sample devoid of any antibody species or isotype that binds to the immobilized protein used (e.g., no antibodies in sample contain kappa light chains when using Immobilized Protein L)	Ensure by other means, such as an ELISA or isotyping kit, that the sample contains IgG-type antibody (see Related Thermo Scientific Products)
Considerable antibody purified, but no specific antibody of interest detected	Antibody of interest is at low concentration or has low binding affinity for the immobilized protein relative to other immunoglobulins in the sample	Use serum-free medium for cell supernatant samples
		Affinity-purify the antibody using the specific antigen coupled to a support such as Thermo Scientific AminoLink Plus Kit (see Related Products)
Antibody of interest purified, but it is denatured (as determined by lack of function in downstream assay)	Antibody is sensitive to low-pH Elution Buffer	Try Gentle Ag/Ab Elution Buffer (see Related Thermo Scientific Products)
	Downstream application is sensitive to neutralized Elution Buffer	Desalt or dialyze eluted sample into an application-compatible buffer

Additional Information Available on Our Website

- Tech Tip #34: Binding characteristics of antibody-binding proteins
- Tech Tip #43: Protein stability and storage

Related Thermo Scientific Products

89896	Pierce Centrifuge Columns, 2mL, 25 units
89897	Pierce Centrifuge Columns, 5mL, 25 units
21001	Protein A IgG Binding Buffer, 1L
21019	Protein G IgG Binding Buffer, 1L
54200	Protein A/G IgG Binding Buffer, 240mL
21027	Gentle Ag/Ab Elution Buffer, 500mL
37501	Monoclonal Antibody Isotyping Kit I (HRP/ABTS)
44894	AminoLink™ Plus Immobilization Kit
66385	Slide-A-Lyzer Dialysis Cassette Kit
66528	Slide-A-Lyzer Concentrating Solution, 200mL
89893	Zeba Spin Desalting Columns, 10mL, 5 columns

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2013 Thermo Fisher Scientific Inc. All rights reserved. All (other) trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.