

Zip Alexa Fluor™ Rapid Antibody Labeling Kits

Pub. No. MAN0017653 Rev. A.0

For detailed protocol, visit the product page at www.thermofisher.com.

Product	Cat. No.	Amount
Zip Alexa Fluor™ 488 Rapid Antibody Labeling Kit	Z11233	1 kit
Zip Alexa Fluor™ 555 Rapid Antibody Labeling Kit	Z11234	1 kit
Zip Alexa Fluor™ 647 Rapid Antibody Labeling Kit	Z11235	1 kit
Store desiccated and protected from light in the freezer (-5°C to -30°C). When stored as directed, the product is stable for at least 6 months.		

1. Prepare a 1 M sodium bicarbonate solution by re-suspending Component A in 1 mL of deionized water.
2. Adjust the antibody concentration to 1.0 mg/mL. To 100 μ L antibody, add 10 μ L of the 1 M sodium bicarbonate solution (from step 1).
3. Add the antibody solution (from step 2) to Component B. Invert or pipet up and down to mix, then incubate for 15 minutes at room temperature.
4. Your conjugate is now ready to use.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT.

Important Licensing Information: These products may be covered by one or more Limited Use Label Licenses. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

 **Manufacturer:** Life Technologies Corporation | 29851 Willow Creek Road | Eugene, OR 97402

© 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/support
thermofisher.com

ThermoFisher
SCIENTIFIC