

# Pierce Direct IP Kit

26148

2121.8

Number	Description
26148	<p><b>Pierce Direct IP Kit</b>, contains sufficient reagents to perform 50 reactions using 10<math>\mu</math>L of immobilized antibody support</p> <p><b>Kit Contents:</b></p> <p><b>AminoLink Plus Coupling Resin</b>, 2mL of settled resin supplied as a 50% slurry (e.g., 100<math>\mu</math>L of 50% slurry is equivalent to 50<math>\mu</math>L of settled resin)</p> <p><b>20X Coupling Buffer</b>, 25mL, when diluted results in 0.01M sodium phosphate, 0.15M sodium chloride; pH 7.2</p> <p><b>Sodium Cyanoborohydride Solution (5M)</b>, 0.5mL</p> <p><b>AminoLink Quenching Buffer</b>, 50mL, 1M Tris•HCl</p> <p><b>Wash Solution</b>, 60mL, 1M NaCl</p> <p><b>IP Lysis/Wash Buffer</b>, 2 <math>\times</math> 50mL, 0.025M Tris, 0.15M NaCl, 0.001M EDTA, 1% NP-40, 5% glycerol; pH 7.4</p> <p><b>100X Conditioning Buffer</b>, 5mL, neutral pH buffer</p> <p><b>20X Tris-Buffered Saline</b>, 25mL, when diluted results in 0.025M Tris, 0.15M NaCl; pH 7.2</p> <p><b>Elution Buffer</b>, 50mL, pH 2.8, contains primary amine</p> <p><b>Non-Reducing Lane Marker Sample Buffer (5X)</b>, 5mL, 0.3M Tris•HCl, 5% SDS, 50% glycerol, lane marker tracking dye; pH 6.8</p> <p><b>Pierce Spin Columns – Screw Cap</b>, 100 columns, includes accessories</p> <p><b>Microcentrifuge Collection Tubes</b>, 2mL, 100 each</p> <p><b>Microcentrifuge Sample Tubes</b>, 1.5mL, 50 each</p> <p><b>Pierce Control Agarose Resin (crosslinked 4% beaded agarose)</b>, 2mL of settled resin supplied as a 50% slurry (e.g., 100<math>\mu</math>L of 50% slurry is equivalent to 50<math>\mu</math>L of settled resin)</p> <p><b>Storage:</b> Upon receipt store at 4°C. Product shipped at ambient temperature.</p>

## Introduction

The Thermo Scientific™ Pierce™ Direct IP Kit enables highly effective and efficient antigen immunoprecipitations by directly immobilizing purified antibodies onto an agarose support. Immobilizing the antibody provides faster and easier immunoprecipitations, enables reuse of the immobilized antibody, and results in purified antigen free from antibody contamination. Immunoprecipitation is achieved using less than 10 $\mu$ g of antibody and a short coupling protocol. After the antibody is coupled to the Thermo Scientific™ AminoLink™ Plus Resin, the antigen sample is incubated with the immobilized antibody to form the immune complex. The complex is washed to remove non-bound material, and a low pH elution buffer is used to dissociate the bound antigen from the antibody.

In contrast to traditional methods, the direct IP kit uses an amine-reactive support that does not contain Protein A or Protein G and does not need a crosslinker for covalent immobilization. Furthermore, this method couples any primary amine-containing molecule, unlike methods requiring an antibody species and subclass that binds strongly to Protein A or Protein G. The kit includes optimized buffers for high antigen yield and efficient spin columns and collection tubes, which shortens the protocol by minimizing handling and mixing.

---

## Important Product Information

- Amines (e.g., Tris, glycine) in the antibody solution will compete for coupling sites on the resin. Remove amines before coupling using Zeba™ Spin Desalting Columns or Slide-A-Lyzer™ Dialysis Cassettes.
- Gelatin or carrier proteins in the antibody solution will compete for coupling sites on the resin. Remove gelatin and carrier proteins using the Pierce Antibody Clean-up Kit (Product No. 44600) or by performing Protein A/G purification (Product No. 20423) and dialysis.
- Perform antibody coupling at room temperature. Perform cell lysis and antigen IP at 4°C.
- Perform all resin centrifugation steps for 30-60 seconds at low speed (i.e., 1000 × g). Centrifuging at speeds greater than 5000 × g may cause the resin to clump and make resuspending difficult.
- When centrifuging spin columns, the flow-through volume should not exceed 600µL when using a 2mL collection tube and 300µL when using a 1.5mL collection tube. Exceeding these volumes may result in back pressure in the column and incomplete washing or elution.
- IP Lysis/Wash Buffer has been tested on representative cell types including but not limited to the following cell lines: HeLa, Jurkat, A431, A549, MOPC, NIH 3T3 and U2OS. Typically, 10<sup>6</sup> HeLa cells yields ~10mg of cell pellet and ~3µg/µL (or 300µg) when lysed with 100µL of IP Lysis/Wash Buffer.
- For best results, add Halt™ Protease (Product No. 78429) and Phosphatase (Product No. 78420) Inhibitor Cocktails to minimize degradation and dephosphorylation of cell lysate proteins.
- The IP Lysis/Wash buffer is compatible with the Pierce BCA Protein Assay (Product No. 23225).
- Proper controls are vital for identifying relevant interactions. The supplied Pierce Control Agarose Resin is composed of the same support material as the AminoLink Plus Coupling Resin but is not amine-reactive. This resin provides an excellent negative control.
- The Pierce Spin Columns package includes spin columns, screw caps, column plugs, Luer-Lok™ Adapter Caps, large frits and a large frit tool. The large frit is not needed for the standard IP protocol. When scaling-up the IP reaction (i.e., > 200µL of resin), the large frit can be inserted into the column to facilitate washing. The Luer-Lok Caps have a flip top that may be used during wash steps. Use the screw caps for sealing the spin columns during storage (see the Additional Information Section).

## Additional Material Required

- Microcentrifuge collection tubes, 2mL

## Procedure for the Pierce Direct IP Kit

### A. Coupling of Antibody to AminoLink Plus Coupling Resin

**Note:** The following protocol is for coupling 2-10µg of affinity-purified antibody in a solution that is free of amines (e.g., Tris, glycine) and carrier proteins (see the Important Product Information Section). This protocol can be scaled up as needed; see the Additional Information Section at the end of the protocol for suggested antibody and resin volumes.

1. Equilibrate the AminoLink Plus Coupling Resin and reagents to room temperature.
2. Prepare 2mL of 1X Coupling Buffer for each IP reaction by diluting the 20X Coupling Buffer with ultrapure water.
3. Gently swirl the bottle of AminoLink Plus Coupling Resin to obtain an even suspension. Using a wide-bore or cut pipette tip, add 20µL of the resin slurry into a Pierce Spin Column. Place column into a microcentrifuge tube and centrifuge at 1000 × g for 1 minute. Discard the flow-through.
4. Wash resin twice by adding 200µL of 1X Coupling Buffer, centrifuge and discard the flow-through.
5. Gently tap the bottom of the spin column on a paper towel to remove any excess liquid and insert the bottom plug.
6. Prepare 2-10µg of affinity-purified antibody for coupling by adjusting the volume to 200µL, using sufficient ultrapure water and 20X Coupling Buffer to produce 1X Coupling Buffer. For example, for 20µL of antibody add 10µL of 20X Coupling Buffer and 170µL of water. Add the ultrapure water, 20X Coupling Buffer and affinity-purified antibody directly to the resin in the spin column.

7. In a fume hood, add 3µL of the Sodium Cyanoborohydride Solution (5M) for every 200µL of reaction volume.  
**Note:** Sodium cyanoborohydride is highly toxic. Wear gloves and use caution when handling.
8. Attach the screw cap to the column and incubate on a rotator or mixer at room temperature for 90-120 minutes, ensuring that the slurry remains suspended during incubation.
9. Remove and retain the bottom plug and loosen the screw cap. Place the spin column into a collection tube and centrifuge. Save the flow-through to verify antibody coupling.
10. Remove the screw cap, add 200µL of 1X Coupling Buffer, centrifuge and discard the flow-through. Repeat this step once.
11. Add 200µL of AminoLink Quenching Buffer to the column, centrifuge and discard the flow-through.
12. Tap the bottom of the column on a paper towel to remove excess liquid and insert the bottom plug. Add 200µL of AminoLink Quenching Buffer to the resin.
13. In a fume hood, add 3µL of Sodium Cyanoborohydride Solution and attach the screw cap. Incubate for 15 minutes with gentle shaking or end-over-end mixing.
14. Remove plug and loosen the screw cap. Place spin column in a collection tube, centrifuge and discard the flow-through.
15. Remove screw cap, wash the resin twice with 200µL of 1X Coupling Buffer, centrifuging after each wash.
16. Wash the resin six times with 150µL of Wash Solution, centrifuging after each wash.
17. Either proceed to immunoprecipitation or, if storing the resin, proceed to the next step.
18. Wash the resin twice with 200µL of 1X Coupling Buffer, centrifuging after each wash.
19. Tap the column on a paper towel to remove excess liquid and place plug in bottom of spin column. Add 200µL of 1X Coupling Buffer, attach the screw cap and store column at 4°C. For long-term storage, add sodium azide to a final concentration of 0.02%.

## B. Mammalian Cell Lysis

### Protocol I: Lysis of Cell Monolayer (Adherent) Cultures

1. Carefully remove (decant) culture medium from cells.
2. Wash the cells once with 1X Coupling Buffer.
3. Add ice-cold IP Lysis/Wash Buffer (Table 1) to the cells. Incubate on ice for 5 minutes with periodic mixing.

**Table 1. Suggested volume of IP Lysis/Wash Buffer to use for different standard culture plates.**

<u>Plate Size/Surface Area</u>	<u>Volume of IP Lysis/Wash Buffer</u>
100 × 100mm	500µL-1mL
100 × 60mm	250-500µL
6-well plate	200-400µL per well
24-well plate	100-200µL per well

4. Transfer the lysate to a microcentrifuge tube and centrifuge at ~ 13,000 × g for 10 minutes to pellet the cell debris.
5. Transfer supernatant to a new tube for protein concentration determination and further analysis.

### Protocol II: Lysis of Cell Suspension Cultures

1. Centrifuge the cell suspension at 1000 × g for 5 minutes to pellet the cells. Discard the supernatant.
2. Wash cells once by suspending the cell pellet in 1X Coupling Buffer. Centrifuge at 1000 × g for 5 minutes to pellet cells.
3. Add ice-cold IP Lysis/Wash Buffer to the cell pellet. Use 500µL of IP Lysis/Wash Buffer per 50mg of wet cell pellet (i.e., 10:1 v/w). If using a large amount of cells, first add 10% the final volume of IP Lysis/Wash Buffer to the cell pellet and pipette the mixture up and down to mix. Add the remaining volume of buffer to the cell suspension.
4. Incubate lysate on ice for 5 minutes with periodic mixing. Remove cell debris by centrifugation at ~ 13,000 × g for 10 minutes.
5. Transfer supernatant to a new tube to determine the concentration and further analysis.

### C. Pre-clear lysate using the Control Agarose Resin

1. For 1mg of lysate, add 80 $\mu$ L of the Pierce Control Agarose Resin slurry (40 $\mu$ L of settled resin) into a spin column.
2. Centrifuge column to remove storage buffer.
3. Add 100 $\mu$ L of 1X Coupling Buffer to the column, centrifuge and discard the flow-through.
4. Add 1mg of lysate to the column containing the resin and incubate at 4°C for 30 minutes to 1 hour with gentle end-over-end mixing.
5. Centrifuge column at 1000  $\times$  g for 1 minute. Discard the column containing the resin and save the flow-through, which will be added to the immobilized antibody.

### D. Antigen Immunoprecipitation General Protocol

**Note:** The amount of sample needed and the incubation time are dependent upon each specific antibody-antigen system and may require optimization for maximum yield.

1. Remove the bottom plug and loosen the screw cap of the spin column containing the antibody-coupled resin. Place the spin column in a collection tube and centrifuge to remove storage buffer. Discard the flow-through.
2. Remove screw cap and place column into a collection tube. Wash resin twice with 200 $\mu$ L of ice-cold IP Lysis/Wash Buffer. Discard the flow-through after each wash.
3. Tap bottom of the spin column on a paper towel to remove excess liquid. Replace bottom plug.
4. Dilute the cell extract in IP Lysis/Wash Buffer. The recommended sample volume in the spin column is 300-600 $\mu$ L. The suggested amount of total protein per IP reaction is 500-1000 $\mu$ g, as determined by the Pierce BCA Protein Assay.
5. Add the sample to the antibody-coupled resin in the spin column. Attach the screw cap and incubate column with gentle end-over-end mixing or shaking for 1 hour to overnight at 4°C.
6. Remove bottom plug, loosen the screw cap and place the column in a collection tube. Centrifuge column and save the flow-through. Do not discard flow-through until confirming that the IP was successful.
7. Remove the screw cap, place the column into a new collection tube, add 200 $\mu$ L of IP Lysis/Wash Buffer and centrifuge.  
**Note:** An alternative wash buffer, 20X Tris-Buffered Saline (20X TBS), is supplied if a detergent-free wash is required. Dilute buffer to 1X before use.
8. Wash the sample three times with 200 $\mu$ L IP Lysis/Wash Buffer and centrifuge after each wash.
9. Wash the sample once with 100 $\mu$ L of 1X Conditioning Buffer (diluted from 100X using ultrapure water).

### E. Antigen Elution

**Note:** To neutralize the low pH of the Elution Buffer (e.g., for downstream enzymatic or functional assays), add 8 $\mu$ L of 1M Tris, pH 9.5 to the collection tube, which will neutralize the pH upon centrifugation (Step E3). Alternatively, use a neutral pH elution buffer (i.e., Gentle Elution Buffer, Product No. 21027).

1. Place the spin column into a new collection tube, add 25 $\mu$ L of Elution Buffer and centrifuge.
2. Keep the column in the tube and add 75 $\mu$ L of Elution Buffer. Incubate for 10 minutes at room temperature. The column does not need to be closed or mixed.  
**Note:** For a more concentrated eluate, less Elution Buffer may be used; however, overall yield might be reduced.
3. Centrifuge the tube and collect the flow-through. Analyze the eluate for presence of antigen. Perform additional elutions (i.e., Steps E1-E3) as needed. Analyze each eluate separately to ensure that the antigen has completely eluted.
4. To preserve activity of the antibody-coupled resin, immediately proceed to Section F, Resin Regeneration and Storage.

### F. Resin Regeneration and Storage

1. Add 100 $\mu$ L of 1X Coupling Buffer to the spin column, centrifuge and discard the flow-through. Repeat this step once.

2. Replace the bottom plug on the spin column. Add 200µL of 1X Coupling Buffer to spin column. Replace screw cap. Wrap the bottom of the tube with laboratory film to prevent resin from drying. For long-term storage (i.e., > 2 weeks) add sodium azide at a final concentration of 0.02%.

### G. Sample Preparation for SDS-PAGE Analysis

1. Equilibrate the Non-Reducing Lane Marker Sample Buffer (5X) to room temperature. Gently mix the sample buffer by inverting 5-10 times. For a reducing gel, add 1M DTT to a final concentration of 100mM in the 5X Sample Buffer.
2. Add Non-Reducing Lane Marker Sample Buffer (5X) to sample to make a 1X final solution (i.e., add 5µL of Non-Reducing Lane Marker Sample Buffer to 20µL of sample).
3. Heat the sample at 95-100°C for ~5 minutes. Allow the sample to cool to room temperature before applying to the gel.

### Troubleshooting

Problem	Possible Cause	Solution
Antibody detected with the eluted antigen	Non-coupled antibody was not removed sufficiently with the Wash Solution during the coupling procedure	Wash the antibody-coupled resin with Elution Buffer until no additional antibody elutes from the resin, as determined by protein assay or measuring the absorbance at 280nm
	The antibody-coupled resin was treated with a reducing agent (i.e., DTT or β-mercaptoethanol) during the IP or elution steps, which reduced the antibody and eluted antibody fragments or subunits that were not covalently linked to the resin	Use buffers that do not contain reducing agents
Antigen did not immunoprecipitate	Sample does not contain enough antigen to detect	Verify protein expression and/or lysis efficiency of the lysate by SDS-PAGE or Western blot
	Antibody did not couple to the resin	Make sure the antibody solution does not contain amines or carrier proteins; verify the antibody coupling by monitoring the flow-through and wash fractions (i.e., measure the absorbance at 280 nm or analyze by SDS-PAGE)
	Component in the IP Lysis/Wash Buffer interfered with antibody-antigen binding	Perform the IP and washes using 1X Tris-Buffered Saline
Antigen did not elute	The antibody-antigen interaction was not disrupted by the Elution Buffer	Optimize the elution conditions (see Tech Tip #27 on our website)
		Elute the antigen as follows: add 100µL of 1X non-reducing SDS sample buffer to the column and incubate at 100°C for 5-10 minutes – keep the spin column in the tube while heating and do not plug or cap the column*
Antigen is non-functional in the downstream application	Antigen is sensitive to low pH and has become inactive during the elution process	Repeat the IP and use a high-salt, neutral pH elution buffer, such as the Gentle Elution Buffer (Product No. 21027)

\*After heating the antibody-coupled resin with SDS sample buffer, the resin cannot be reused and must be discarded.

## Additional Information

### A. Large-scale Antibody Coupling

The Pierce Direct IP Protocol can be scaled up as needed. Guidelines for scaling up the amounts of antibody and resin are listed in Table 2. Adjust volumes for coupling reagents and washes accordingly.

**Table 2. Amount of Thermo Scientific AminoLink Plus Coupling Resin and antibody to use for larger scale IPs.**

<u>Antibody Amount (µg)</u>	<u>Resin Slurry Amount (µL)</u>	<u>Reaction Volume (µL)</u>
25-100	100	200
50-200	200	300
100-400	400	400


### B. Visit the website for additional information including the following:

- Tech Tip #27: Optimize elution conditions for immunoaffinity purification
- Tech Tip #43: Protein stability and storage
- Tech Tip #40: Convert between times gravity ( $\times g$ ) and centrifuge rotor speed (RPM)


### C. Pierce Spin Columns

Pierce Spin Columns can hold up to 900µL. Columns can be placed in 1.5mL or 2mL microcentrifuge tubes or used with a Luer-lok Adapter (see Figure 1) for processing samples with a syringe. When using a syringe, sample size and wash volumes are only limited by the volume capacity of the syringe. For small volumes of resin, use columns with only the small, pre-inserted frit (Figure 2A). For applications requiring more than 100µL of resin, the large frit may be used at either the top or bottom (Figure 2B and 2C). Resins may be used repeatedly when the resin is in-between the small and large frit (Figure 2B).

- To remove a frit from a column, use an unfolded paper clip and insert the wire through the column tip and push the frit.
- To insert a frit, place the frit inside the column and use the frit tool to push the frit into position.
- To remove the top frit from an already-packed column with a top and bottom frit, use an unfolded paper clip to tip the top frit up. The top frit can be then removed with tweezers.



**Figure 1.** Schematic of the Pierce Spin Column contents.



**Figure 2.** Three formats for frit placement.

## Related Thermo Scientific Products

78428	<b>Halt Phosphatase Inhibitor Single-Use Cocktail (100X)</b> , 100µL $\times$ 24 microtubes
78430	<b>Halt Protease Inhibitor Single-Use Cocktail (100X)</b> , 24 $\times$ 100µL
28348	<b>20X Phosphate Buffered Saline (Coupling Buffer)</b> , 500mL
69705	<b>Pierce Spin Columns – Screw Cap</b> , 25/pkg
69720	<b>Pierce Microcentrifuge Tubes, 2mL</b> , 72/pkg
69715	<b>Pierce Microcentrifuge Tubes, 1.5mL</b> , 72/pkg
21027	<b>Gentle Ag/Ab Elution Buffer, pH 6.6</b> , 500mL
21004	<b>IgG Elution Buffer</b> , low-pH elution buffer for general protein affinity purifications, 1L

---

Products are warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation"). No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than Buyer. Any model or sample furnished to Buyer is merely illustrative of the general type and quality of goods and does not represent that any Product will conform to such model or sample.

NO OTHER WARRANTIES, EXPRESS OR IMPLIED, ARE GRANTED, INCLUDING WITHOUT LIMITATION, IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR ANY PARTICULAR PURPOSE, OR NON INFRINGEMENT. BUYER'S EXCLUSIVE REMEDY FOR NON-CONFORMING PRODUCTS DURING THE WARRANTY PERIOD IS LIMITED TO REPAIR, REPLACEMENT OF OR REFUND FOR THE NON-CONFORMING PRODUCT(S) AT SELLER'S SOLE OPTION. THERE IS NO OBLIGATION TO REPAIR, REPLACE OR REFUND FOR PRODUCTS AS THE RESULT OF (I) ACCIDENT, DISASTER OR EVENT OF FORCE MAJEURE, (II) MISUSE, FAULT OR NEGLIGENCE OF OR BY BUYER, (III) USE OF THE PRODUCTS IN A MANNER FOR WHICH THEY WERE NOT DESIGNED, OR (IV) IMPROPER STORAGE AND HANDLING OF THE PRODUCTS.

Unless otherwise expressly stated on the Product or in the documentation accompanying the Product, the Product is intended for research only and is not to be used for any other purpose, including without limitation, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses, or any type of consumption by or application to humans or animals.

Current product instructions are available at [www.thermoscientific.com/pierce](http://www.thermoscientific.com/pierce). For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2014 Thermo Fisher Scientific Inc. All rights reserved. Luer-Lok is a trademark of Becton Dickinson and Company. Unless otherwise indicated, all other trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.