invitrogen

QUICK REFERENCE

Qubit[™] 4 Fluorometer System Verification Assay Kit

Catalog No. Q33237

Pub. No. MAN0018153

Rev. A.0

For detailed protocol, visit the product page at www.thermofisher.com.

Material	Amount	Storage
Blank Reagent (Component A)	10 mL	Room Temperature
Green Fluorescence Reagent (Component B)	10 mL	Room Temperature Protect from light
Far Red Fluorescence Reagent (Component C)	10 mL	

When stored as directed, kit components are stable for at least 12 months.

S

- 1. Set up three Qubit™ assay tubes (Cat. No. Q32856).
- 2. Label the tube lids: 1-3.
- 3. Add 200 µL of Blank Reagent to Tube 1.
- 4. Add 200 µL of Green Fluorescence Reagent to Tube 2.
- 5. Add 200 µL of Far Red Fluorescence Reagent to Tube 3.
- On the Qubit™ 4 Fluorometer, select Settings ➤ System Verification Test.
- Following the on-screen instructions, read the tubes in the following order: 1–3.

The information in this guide is subject to change without notice.

DISCLAIMER: TO THE EXTENT ALLOWED BY LAW, LIFE TECHNOLOGIES AND/OR ITS AFFILIATE(S) WILL NOT BE LIABLE FOR SPECIAL, INCIDENTAL, INDIRECT, PUNITIVE, MULTIPLE, OR CONSEQUENTIAL DAMAGES IN CONNECTION WITH OR ARISING FROM THIS DOCUMENT, INCLUDING YOUR USE OF IT. Important Licensing Information: These products may be covered by one or more Limited Use Label

Licenses. By use of these products you accept the terms and conditions of all applicable Limited Use Label Licenses.

- Manufacturer: Life Technologies Corporation | 29851 Willow Creek Road | Eugene, OR 97402
- © 2018 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified.

For support visit thermofisher.com/support thermofisher.com

