

Immobilized Reductant Columns

77701

0275.5

Number	Description
77701	<p>Immobilized Reductant Column, pre-packed with 2mL Immobilized Reductant Resin, plus column accessory pack containing two white tips</p> <p>Support: 6% crosslinked beaded agarose</p> <p>Capacity: 20-40μmol -SH groups per column per use</p> <p>Storage: Upon receipt store at 4°C. Product is shipped at ambient temperature.</p>

Introduction

Thermo Scientific Immobilized Reductant Columns allow quick and convenient reduction of disulfide bonds in peptide and protein samples. Peptides are reduced efficiently in the few minutes required for the sample to pass through the gravity-flow column. This reduction is convenient because tedious and time-consuming gel filtration or HPLC steps are not required for subsequent separation of reducing agent from the sample. Instead, the proprietary reducing agent remains attached to the beaded agarose support and does not contaminate the recovered sample. This feature is especially important when treating peptides, whose small size prevents them from being desalted effectively by standard gel filtration.

Sulfhydryl (-SH) groups that occur in the side chain of cysteine residues in peptides and proteins are excellent targets for conjugation reactions. Applications for such conjugation reactions include crosslinking peptides to carrier proteins for use as immunogens and immobilization to solid supports for affinity purification of antibodies or binding partners. However, cysteine-containing peptides tend to dimerize in solution through formation of disulfide bonds, making the sulfur atoms unavailable for conjugation; disulfide bonds must be reduced to make sulfhydryls available for conjugation.

Although reduction of disulfide bonds in peptides and proteins can be performed with either dithiothreitol (DTT) or 2-mercaptoethanol (2-ME), these reducing agents must be removed by gel filtration before sulfhydryl-directed conjugation reactions can be performed. This extra step to remove the reducing agent can be time-consuming (allowing disulfide bonds to form again) and may result in loss or dilution of the sample peptide or protein.

Thermo Scientific Immobilized Reductant Columns circumvent these processing difficulties to enable rapid and efficient preparation of reduced samples for applications that cannot also contain reducing agent. Each column may be regenerated four times, for a total of five uses.

Additional Materials Required

- Equilibration Buffer #1: The optimal buffer for the reducing reaction is 0.1M sodium phosphate, 1mM EDTA, pH 8.0. Other buffers, such as phosphate-buffered saline (PBS), Tris-buffered saline (TBS) or HEPES may also be used.
- Equilibration Buffer #2 (optional): 0.1M sodium phosphate, 1mM EDTA, 6M guanidine•HCl, pH 8.0. This or similar buffer containing denaturant will help to expose protein disulfide bonds to the reducing agent and so make complete reduction of proteins more efficient.
- Dithiothreitol (DTT) Activation Solution: See related Pierce Products. For each use of a Reductant Column, a 10mM DTT solution in Equilibration Buffer #1 will be needed (15mg DTT dissolved in 10mL Equilibration Buffer #1). This solution must be prepared immediately before use (Step 4 of procedure).
- Protein or peptide sample: Dissolve 5-10mg peptide or 5-6mg protein in 1mL of Equilibration Buffer #1 or #2.

Procedure for Peptide or Protein Disulfide Reduction

1. Equilibrate Immobilized Reductant Column and prepared Equilibration Buffers (hereafter “Buffers”) to room temperature before use. Perform all steps of this procedure at room temperature.
2. Remove top cap from the Immobilized Reductant Column and pour off storage solution, which contains 0.02% sodium azide.
Note: When uncapping a column, always remove the top cap first to avoid drawing air into the resin bed. If the top porous disc in the column is dislodged, push it down to the resin bed surface (but avoid compressing the resin bed itself).
3. Twist off bottom column tab and stand column upright.
4. Equilibrate/wash the Immobilized Reductant by adding 5mL of Buffer #1 to the column and allowing it to drain through.
5. Prepare DTT Activation Solution: Dissolve 15mg of DTT in 10mL Buffer #1 (results in 10mM DTT).
6. Activate the column by applying 10mL of DTT Activation Solution to the column and allowing it to drain through.
7. To remove the activating DTT, wash the column with 10mL of Buffer #1. (If Buffer #2 was chosen and used for the sample, wash column with only 5mL of Buffer #1, followed by 5mL of Buffer #2.)
8. Apply 1mL of peptide or protein solution to the column. Allow the sample to completely enter the resin bed. (Flow will stop automatically when the liquid drains down to the top porous disc.)
9. Cap top and bottom of column (white bottom tips supplied) and incubate for 60 minutes for protein samples; no incubation is needed for peptide samples.
10. Recover reduced peptide or protein sample from the column by applying 9mL of the appropriate Buffer and collecting separate 1mL fractions as they emerge from the column tip. To ensure that the volume of each fraction is exactly 1mL, apply only 1mL of buffer at a time and collect the entire volume that emerges until column flow stops; then change collection tubes and apply the next 1mL of buffer. This method utilizes the stop-flow feature of the top porous disc.
11. Identify the fractions that contain the peptide or protein (now reduced) by determining which ones have peak absorbance at 280nm (usually the first two fractions). Be aware that some peptides do not absorb significantly at 280nm and cannot be assayed by this method. Alternatively, identify fractions that contain the reduced peptide or protein by specifically measuring for sulfhydryl groups using the Ellman’s Reagent (Product No. 22582, see Related Thermo Scientific Products).

Column Regeneration and Storage

Reuse the Immobilized Reductant Column by simply repeating the procedure from the beginning; the DTT activation step will regenerate the column. Four cycles of regeneration (5 total uses) are possible with each Immobilized Reductant Column. For prolonged storage of used columns, first wash them with phosphate buffer or water containing 0.02% sodium azide, then replace top and bottom caps and store the columns upright at 4°C.

Related Thermo Scientific Products

20291	Dithiothreitol (DTT), No-Weigh™ Format, 7.7mg × 48 microtubes
20290	Dithiothreitol (DTT), 5g
22582	Ellman’s Reagent, 5g
44899	CarboxyLink™ Immobilization Kit, 2mL
53155	UltraLink® Iodoacetyl Resin, 10mL
21901	EZ-Link® Maleimide-PEG ₂ -Biotin, 50mg
77712	Immobilized TCEP Disulfide Reducing Gel, 5mL

Product References

- Stimmel, J.B., *et al.* (2000). Site-specific conjugation on serine-cysteine variant monoclonal antibodies. *J Biol Chem* **275(39)**:30445-50.
- Corse, E. and Machamer, C.E. (2000). Infectious bronchitis virus E protein is targeted to the golgi complex and directs release of virus-like particles. *J Virol* **74(9)**: 4319-26.
- Segura-Totten, M, *et al.* (2002). Barrier-to-autointegration factor: major roles in chromatin decondensation and nuclear assembly. *J Cell Biol* **158**: 475-85.
- Al-Hallaq, R.A., *et al.* (2002). Association of NR3A with the N-methyl-D-aspartate receptor NR1 and NR2 subunits. *Mol Pharmacol* **62(5)**: 1119-27.

This product ("Product") is warranted to operate or perform substantially in conformance with published Product specifications in effect at the time of sale, as set forth in the Product documentation, specifications and/or accompanying package inserts ("Documentation") and to be free from defects in material and workmanship. Unless otherwise expressly authorized in writing, Products are supplied for research use only. No claim of suitability for use in applications regulated by FDA is made. The warranty provided herein is valid only when used by properly trained individuals. Unless otherwise stated in the Documentation, this warranty is limited to one year from date of shipment when the Product is subjected to normal, proper and intended usage. This warranty does not extend to anyone other than the original purchaser of the Product ("Buyer").

No other warranties, express or implied, are granted, including without limitation, implied warranties of merchantability, fitness for any particular purpose, or non infringement. Buyer's exclusive remedy for non-conforming Products during the warranty period is limited to replacement of or refund for the non-conforming Product(s).

There is no obligation to replace Products as the result of (i) accident, disaster or event of force majeure, (ii) misuse, fault or negligence of or by Buyer, (iii) use of the Products in a manner for which they were not designed, or (iv) improper storage and handling of the Products.

Current product instructions are available at www.thermoscientific.com/pierce. For a faxed copy, call 800-874-3723 or contact your local distributor.

© 2011 Thermo Fisher Scientific Inc. All rights reserved. Unless otherwise indicated, all trademarks are property of Thermo Fisher Scientific Inc. and its subsidiaries. Printed in the USA.